

Recommended Handling & Plate Erasing Tips

Phosphor Plate – Handling Tips

- Avoid touching the phosphor side of the plate with long fingernails.
- When positioning a plate you may bend and flex accordingly. However, **DO NOT FOLD OR CREASE.**
- When using Snap-O-Ray positioning devices, be sure to file sharp teeth down to prevent damaging of the plate. Also, you can fill the teeth on the holder with dental composite to prevent damage.
- Never drag the plate across any surface with the phosphor side facing down, i.e. do not use a view box to erase plates.
- Clean plates with 100% cotton gauze and anhydrous isopropyl alcohol.

Phosphor Plate – Erasing Tips

- Prevent or minimize any contact with the phosphor side of the plate to reduce surface scratches.
- Only use [Air Techniques Phosphor Plate Eraser](#), to obtain an evenly distributed light source for proper erasure and no scratching.

