

# ASSISTINA 301 plus

maintenance system



## CONSERVES

lowest operating cost  
of any other handpiece  
maintenance device

## CLEANS

flushes internal air/  
water coolant channels  
with cleaning liquid

## ROTATES

helps remove debris  
by rotating gears  
during purge cycle

# An Easy, Consistent Way To Extend The Life Of Your Handpieces

## What It Offers You:

- **value**—at less than 2 cents per cycle, pays for itself after about 8,500 cycles
- **time savings**—automatically timed operation takes only 35 seconds a cycle
- **convenience**—simple to install and easy to use
- **compatibility**—includes universal adapter that accepts many handpieces\*
- **less cleanup**—filters exhaust air and keeps aerosols contained

## Why It Works:

**The Assistina provides consistent maintenance—a key to extending the life of your handpieces, motors, and attachments.**

- uses a special cleaning liquid to flush internal air/water coolant lines, removing debris and preventing scale buildup
- provides consistent lubrication with automatic oil dispensing system
- spins gears and turbines during purge cycle to help expel debris and excess oil

## How It Works:

**The Assistina cleans and lubricates internal components by using compressed air, maintenance solutions, and a few simple steps:**

- connect the instrument and close the cover
- push the start button for 2 seconds to begin automated delivery of service oil and cleaning liquid
- in 35 seconds, your instrument is properly maintained and prepared for sterilization



\* may require additional adapters

For more information about A-dec|W&H products, please contact your authorized A-dec dealer, or call us today at **1.800.547.1883**.

## Specifications:

<b>cycle time:</b>	approx. 35 seconds	<b>rqrd. air pressure:</b>	65 psi/min; 140 psi/max	<b>height/width/depth:</b>	8.8" x 7.5" x 16.3" (223 x 190 x 415 mm)
<b>capacity:</b>	250 ml	<b>connection:</b>	compressed air connection with quick-swivel coupling for air hose	<b>weight:</b>	5.95 lb. (2.70 kg)
service oil	(enough for about 3,500 cycles)				
cleaning liquid	(enough for about 2,500 cycles)				

A-dec Inc.  
2601 Crestview Drive  
Newberg, OR 97132 USA  
Phone: 1.800.547.1883, 1.503.538.7478  
Fax: 1.503.538.0276  
a-dec.com

The products shown in this literature are based upon the latest product information at the time of printing. A-dec reserves the right to make changes at any time without notice. Please contact your authorized A-dec dealer for the most current product information.

© 2010 A-dec® Inc. All rights reserved. 85.6058.00 Rev B/ PT / 5M / 05-10