

## CUSTOMER INFORMATION

Company Name

Doctor Name

Phone

Email

Address

City

State / Province

Zip Code / Postal Code

## DEALER INFORMATION

Dealer Name

Rep Name

Phone

Email

## IT INFORMATION

Name

Phone

Email

## FORM SUBMISSION CONTACT

Name

Company Name

Date \_\_\_\_\_

### ADDITIONAL COMMENTS

**Questions? Call Technical Support at:**  
**1-800-769-2909**

Email completed form to:  
**techsupport@gendex.com**

**SUBMIT**

REQUIREMENTS		RESULTS	COMMENTS
<b>Fixing Hardware</b> Required fixing hardware depends on the wall and floor material and it is not delivered with the unit. The fixing hardware, as well as wall and floor materials must endure a pull-out force of 1,124 pounds per force. For further details please refer to installation manual.		Meets Specifications Doesn't Meet Specifications	
<b>Weight Limit</b> Make sure that the floor where the unit is to be installed can support this weight. The floor should be able to carry minimum 2200 lbs/11 ft².		Meets Specifications Doesn't Meet Specifications	
<b>Installation Space Requirements</b> See page 3 for Installation Space Requirements and Diagrams.		Meets Specifications Doesn't Meet Specifications	
<b>Room Preparation</b> A radiation shielding plan is required prior to the install of any CBCT unit. Local regulatory requirements must be met. Minimum doorway widths of 28" on all doorways between the building's entrance and the installation site are recommended in order for delivery of the unit. Contact Prophysics at 800-835-3615.		Meets Specifications Doesn't Meet Specifications	
<b>Electrical Specifications</b> 100-120, 220-240VAC, 50/60Hz +/- 10% (10A@230VAC, 15A@110VAC) dedicated power supply (isolated circuit should be powering the unit and acquisition computer). Power conditioner recommended on unit.		Meets Specifications Doesn't Meet Specifications	
<b>Ethernet Connections</b> Direct connection required between unit and acquisition PC with a Cat5e or Cat6 Ethernet cable.		Meets Specifications Doesn't Meet Specifications	
<b>2D Image Acquisition and Viewing – Minimum PC Requirements</b> Windows 7 (32 or 64-bit) 2.5 GHz dual core 3 GB RAM 500 GB hard drive 2 Gb Ethernet 1000Base-T Network cards		Meets Specifications Doesn't Meet Specifications	
<b>3D Image Viewing Only on Additional Workstations – Minimum PC Requirements</b> Windows XP Professional (32-bit) 2.0 GHz dual core 4 GB RAM ATI Radeon HD 4650 or NVIDIA GeForce 260 (256 MB memory) 3 GB free disk space Gb Ethernet network interface		Meets Specifications Doesn't Meet Specifications	
<b>High Speed Internet Connection</b> Broadband internet connection to the workstation is mandatory.		Meets Specifications Doesn't Meet Specifications	
<b>Imaging Software</b>	Software installed:	Version:	Compatible Incompatible
<b>Practice Management Software</b>	Software installed:	Version:	Compatible Incompatible
<b>Other Gendex Products</b>	List all Gendex products in use:		

## INSTALLATION SPACE REQUIREMENTS

### PAN/3D ONLY:

Depth: 66.9" minimum  
 Width: 59" minimum  
 Height: 94.9" minimum / 84.25" to 96" adjustable maximum

### WITH CEPH ARM:

Depth: 66.9" minimum  
 Width: 98.4" minimum  
 Height: 94.9" minimum / 84.25" to 96" adjustable maximum

