

Dealer Installation Checklist for ScanX

Pre-Installation

| Done? | Computer Provider / Technician | |
|--------------------------|--------------------------------|---|
| <input type="checkbox"/> | Dealer / Technician | Ensure stable area and proper location for CPU and ScanX are available. |
| <input type="checkbox"/> | Technician | Ensure that CPU is equipped with a USB 1.1 port or higher. Computers with ECP/EPP will NOT work. |
| <input type="checkbox"/> | Technician | Verify appropriate and approved operating systems for ScanX and Apporved Software integration |
| <input type="checkbox"/> | Technician | Verify that any required peripherals or security devices are installed and operating without conflict. |
| <input type="checkbox"/> | Technician | Verify that computer is assembled, EPP port configuration completed, CPU fully functioning, awaiting software installation. |
| Done? | Software Provider / Technician | |
| <input type="checkbox"/> | Technician | Ensures current and approved integration module purchased. |
| <input type="checkbox"/> | Technician | Installs software and necessary drivers. If ImageFX, drivers for ScanX on CD in Accessory Pack. |
| <input type="checkbox"/> | Technician | Installs software while logged in as the Administrator. Other users with Administrative Privileges only, may not allow for complete driver install. |

Installation Day

| Done? | Dealer | Notes |
|--------------------------|------------------|---|
| <input type="checkbox"/> | Dealer | Dealer arranges for Hardware / Software, Equipment / High Tech Rep, and PSA Rep involved to be on site during final installation. |
| <input type="checkbox"/> | Dealer | Install Software from Approved Software Vendor List. |
| <input type="checkbox"/> | Dealer | Install ScanX Drivers & Utilities Disc supplied with the ScanX located in the accessory pack. |
| <input type="checkbox"/> | Dealer | Power brick connected to ScanX and power source. |
| <input type="checkbox"/> | Dealer | Parallel cable connected to ScanX and corresponding computer. |
| <input type="checkbox"/> | Dealer | Power to ScanX verified. |
| <input type="checkbox"/> | Technician | Computer functions and activates corresponding, approved software title. |
| <input type="checkbox"/> | Technician | Known software values (brightness/contrast values, histogram, filters, etc.) are applied to corresponding software. |
| <input type="checkbox"/> | Technician | Software initializes and locates ScanX. |
| <input type="checkbox"/> | Technician | X-Ray equipment is inspected. Exposure adjustments made, intensifying screens removed. |
| <input type="checkbox"/> | PSA | Image acquisition procedure verified. |
| <input type="checkbox"/> | PSA | ScanX accepts, processes, and produces images. |
| <input type="checkbox"/> | PSA / Technician | Fine tuning and adjustments to images are made (if necessary) and saved for future images. |