

EzRay Control Panel Exposure Switch

EzRay Control Panel Exposure Switch



1. To enter Admin mode, press and hold the settings and save buttons shown in red circles in the picture above until menu is displaying.

EzRay Control Panel Exposure Switch



2. Use the down arrow button to go to #6 option (exposure set), press the settings button shown in top red circle to enter Panel Expo S/W. Toggle to OFF setting.
3. Once OFF displays, press and hold Save button shown in bottom red circle for 5 seconds or until an audible sound confirmation.

EzRay Control Panel Exposure Switch



❖ For States that require Simultaneous Exposure Button Press:

Follow instruction to step 2. From the sub menu, use the down arrow key to get to option #4(Press Type). Toggle to BOTH setting then press save button for 5 seconds or wait for an audible confirmation.

