

Monthly Internal Distilled Water Tank Cleaning - Bravo 21V

1. If you fill your Bravo with distilled water from an external bottle, turn the power off on the Bravo before you begin.
2. Arrange an empty container on the floor near the sterilizer and insert the free end of the drain tube into the container.
3. Insert the other end of tube into the quick-connect port marked “Service” positioned on the front left of the chamber as shown in the figure below.
4. Allow the tank to empty completely, then disconnect the tube.
5. Prepare a 1 Gallon mixture of 90% distilled water and 10% isopropyl alcohol (70%-100% isopropyl alcohol), then pour it into the distilled water tank from the top of the machine. If you can’t access the top distilled water fill port on the machine, use the manual fill bottle with quick connect tube attached (picture below). You will need to use 2 full bottles of the 90/10 distilled water alcohol mixture to completely fill the chamber.
6. Let the solution sit for 30 minutes. If possible, agitate Bravo to move solution around inside tank.
7. Repeat step 2 and 3 to completely empty the internal tank again.
8. Fill the tank with distilled water. If you are using an external bottle to autofill, turn the power back on and the Bravo will fill with distilled water from the external bottle. If you get an error code E040, turn machine off then back on.
9. Run one empty (no load) cycle.

