CLEANSTREAM SYSTEM

Makes it easy to clean the vacuum system efficiently and reliably without polluting our environment.



1295 Walt Whitman Road, Melville, NY 11747 / www.airtechniques.com Toll free in the U.S. and Canada: 800-Air-Tech / Fax 516-433-7684 Western Facility: 800-822-2899 / Fax 909-898-7646

CleanStream™

Evacuation System Cleaner



The solution to a clean vacuum system and a healthy environment.



CleanStream offers you a comprehensive solution to maintain your vacuum system for prolonged life and uninterrupted work flow without harming our environment.



CleanStream™

Organic and inorganic debris (blood, saliva, amalgam) in the tubing of evacuation systems promote the growth of disease causing bacteria, offensive odors, and can even result in system failures.

Every dental vacuum system should be regularly cleaned to maintain peak performance. CleanStream offers you a reliable and effective system to keep your vacuum system running without system failures.

CleanStream Evacuation System Cleaner

- CleanStream, a highly effective evacuation system cleaner reliably removes and prevents debris in the evacuation lines.
- Special non-foaming agents protect and extend the life of your vacuum pump.
- Use CleanStream for all wet and dry evacuation systems as well as systems with an Amalgam Separator.
 Unlike other cleaners CleanStream does not cause leaching of mercury from amalgam.



CleanStream Cleaner PN 57850

- CleanStream Evacuation System Cleaner is environmentally friendly and completely biodegradable.
- The concentrate is diluted with water. For a 1 liter solution of CleanStream, fill the cap of the 2.5 I CleanStream bottle to the line (represents 20ml) with CleanStream concentrate and pour into the Dispenser. For a 2 liter solution of CleanStream repeat this step.

CleanStream Dispenser System

 Easy to use: Add concentrate, fill with water and attach hoses. The CleanStream Dispenser System will automatically aspirate the solution through suction hoses.

 CleanStream Evacuation System Cleaner is aspirated into a mist of very fine particles.
 No stress is added to your vacuum pump.

 The CleanStream Dispenser System holds 2-Liters of solution, sufficient for the cleaning of two units. Connection adapters are available for all suction hoses.

 Durably constructed, the Dispenser
 System is unaffected by chemical and mechanical influences.

CleanStream Dispenser System PN 57665 Includes Dispenser System with fittings and instructions.

CleanStream Starter Kit PN57660 (not shown)
Includes a 120 ml bottle of CleanStream concentrate
(sufficient for 6 applications) and a CleanStream Dispenser System with



Easy-to-use System





The CleanStream System is very easy to use. Just fill the CleanStream Dispenser with water, add one or 2 caps of CleanStream and connect hoses.

