

Forest Dental Catalog 2016

Affordable Excellence + Designer Friendly®



3900 Patient Chair Powder Coat Paint Color: Bronze

9078 LED Light Powder Coat Paint Color: Bronze

Standard IC Control Powder Coat Paint Color: Bronze | Plastic Finish: Laser Leather

Assistant's Instrumentation Powder Coat Paint Color: Bronze | Plastic Finish: Laser Leather

Plush Duo Color Upholstery: Ultraleather® Doe & Lagoon



Let's Dance.



3900 Patient Chair Powder Coat Paint Color: Charcoal
Plush Narrow Back Duo Color Upholstery: Ultraleather® Raven Wing & Sugar
9080 Halogen Light Powder Coat Paint Color: Charcoal
Euro Style Control Powder Coat Paint Color: Charcoal | Plastic Finish: Greywood





UPDATE

YOUR WAY

Updating your practice over a period of time, as you do with your home, following a plan that works with your budget, is one of your most important long-term financial considerations. You'll minimize the financial impact while improving curb appeal and the enjoyment of patients and staff.

THE FOREST WAY

Forest customers are quick to point out that their Forest Treatment Rooms are a key internal marketing component. They create a fun atmosphere that speaks volumes about our customers' aesthetics and attention to detail. Insist on being shown the beauty of Forest.



Hank Barton, President of Forest Dental

My career selling dental equipment started more than 40 years ago. For the first 25, I repped for the largest dental manufacturer in North America. At trade shows, I'd often ask dentists what kind of chair they had. After all, dental chairs are the cornerstone of their treatment room. Many times they did not recall the brand. "It's blue and shaped like this or that." Then 13 years ago, I bought Forest Dental and decided it was time for a new approach.

From the outset, our goal at Forest has been to build great looking core equipment that a dentist could be proud of, unlike anything I'd sold before. Paint colors other than the standard white soon became a major factor in our becoming an industry growth leader. Then we added three carbon fiber look trims and a chrome option, offering a complete aesthetic package. We believe our approach helps customers reflect their personality, style, and branding. We call it Designer Friendly®.

Our other big growth factor is our commitment to Affordable Excellence. In order to give our customers the best possible value for their investment, we mirror many of the expensive features of big brand names with even more precision of fit and finish, while selling for less.

Insist on being shown the beauty of Forest. I know you will be glad you did.

Thank you!

Hank Barton



Control Heads

Materials & Benefits

Standard IC Control Head

Material: aluminum hinged lid / 1/8" ABS plastic chassis with super tough DR acrylic protective coating

Benefit: hinge allows quick access to controls & ancillary products / lid & chassis are resistant to UV, impact & compatible surface cleaners / aluminum lid is strong & lightweight

Material: electroless nickel plated brass air & water T's, routing valves, pneumatic control & distribution blocks

Benefit: lifetime warranty on all parts "under the hood" (on all forest control heads) / trouble-free handpiece operation

Material: aluminum handpiece holder bar with glass filled nylon injection molded handpiece holders

Benefit: allows for up to four handpiece positions, a syringe & ancillary touchpad at your desired angle & spacing

Material: hinged aluminum control head lid

Benefit: available in nine two-coat powder coat paint colors

Material: plastic chassis

Benefit: available in three carbon fiber look & two solid color finishes

Euro Control Head

Material: whip arms rotate on nylon wheels with aluminum & brass components

Benefit: feather-touch whip movement / ultra light handpiece feel / long, trouble free life

Material: control head swivels at its center point

Benefit: 240-degrees of control head rotation allows for beside- or over-the-patient placement

Material: forward mounted aluminum whip arm system

Benefit: allows for greater reach of handpiece & syringe expanding control head placement options

Material: plastic control head lid

Benefit: available in three carbon fiber looks & two solid color finishes

Feature: autoclavable silicone handpiece pad available in black & white

Benefit: blends with Designer Friendly® choices

Value Control Head

Material: PVC pinch valve control blocks

Benefit: ultra-reliable handpiece functions (drive air & air/water coolant) with a long trouble-free life

Material: ABS plastic control cover w/ tough DR acrylic protective coating

Benefit: resistant to UV, impact & compatible surface cleaners / easy removal for quick access to controls & ancillary products

Material: aluminum chassis base plate

Benefit: strong support for handpiece holder bar, positioning handles & control lid

Material: plastic control head lid

Benefit: available in three carbon fiber looks & two solid color finishes

Material: chrome plated or powder coated handpiece holder bar

Benefit: rich look enhancing a value based control head


Standard IC

Please use graph below to determine number of ancillaries and options possible for the Standard IC.

**Standard Features**

Automatic control for three handpieces

asepsis 3-way air/water syringe with tubing, 5' handpiece tubing, 4-hole Midwest, asepsis style, 5' pneumatic master control block with diaphragm master on/off individual handpiece pressure adjustment individual water coolant adjustment non-retraction water coolant air/water coolant adjustment handpiece flush handpiece exhaust filter/oil collector handle brake release function

% Standard IC Control Head Space Absorbed per Option

Please note:

If your ancillary and option choices exceed 100%, please call customer service for additional information.

* Uses existing handpiece position

★ Standard on SE and SI deliveries

• Standard on 1500PI delivery

+ 3900 Chair touchpad cable standard on PI and SI deliveries

◆ Two or more require purchase of 0010-310 ancillary shelf (see page 8)

Ordering Ancillaries & Options

Delivery Units

Euro

Please use graph below to determine number of ancillaries and options possible for the Euro.




Standard Features

Automatic control for four handpieces

asepsis 3-way air/water syringe with tubing, 3' handpiece tubing, 4-hole Midwest, asepsis style, 3' pneumatic master control block with diaphragm master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air/water coolant adjustment
handpiece flush
handpiece exhaust filter/oil collector
dual handle brake release function
swivel instrument tray

% Euro Control Head Space Absorbed per Option



Please note:

If your ancillary and option choices exceed 100%, please call customer service for additional information.


† Uses existing handpiece position

★ Standard on SE and SI deliveries

+ 3900 Chair touchpad cable standard on PI and SI deliveries

Value


Please use graph below to determine number of ancillaries and options possible for the Value.

**Standard Features**

Automatic control for three handpieces

asepsis 3-way air/water syringe with tubing, 5' handpiece tubing, 4-hole Midwest, asepsis style, 5' individual pinch valve control blocks master on/off individual handpiece pressure adjustment individual water coolant adjustment non-retraction water coolant air coolant adjustment handpiece flush handpiece exhaust filter/oil collector handle brake release function

% Value Control Head Space Absorbed per Option



Please note:

If your ancillary and option choices exceed 100%, please call customer service for additional information.

* Uses existing handpiece position

★ Standard on SE and SI deliveries

+ 3900 Chair touchpad cable standard on PI and SI deliveries

Control Head Accessories

Ancillary Products

Bien-Air Optima Int Electric HP Control & HP Control with Endo



Factory Installed	BIEN-AIR-MCX	\$ 2295
Factory Installed	BIEN-AIR-MX2	4595

Please call 800-433-2436 for warranty and technical support of the Handpiece
Bien-Air Optima Int Electric Handpiece Control (BIEN-AIR-MCX)
Uses existing handpiece position
Attachments purchased separately from Bien-air
Highest torque output on the market
Swivel tubing reduces arm fatigue
Quick-disconnect, autoclavable motor with type E connection
(compliant with ISO 3964)
High and low speed applications using one attachment
MX motor speed from 100 to 40,000 rpm
24 volt power supply
Bien-Air warranty 3 years from date of purchase

Bien-Air Optima Int Electric Handpiece Control with endo (BIEN-AIR-MX2)
Uses existing handpiece position
Same as above with additional features
Auto-reverse and torque control for endo procedures
40 memory position with pre-programmed procedures

Satelec SP



Factory Installed	SATELEC-S	\$ 1845
Factory Installed	SATELEC-S-LED	2320

Please call 800-289-6367 for warranty and technical support of the Scaler
Uses existing handpiece position
Satelec SP Newtron LED lit Piezo Ultrasonic Scaler with power supply
Color coding system ensures proper tip and power setting
Three tips: #1, #10Z, H3 (Intro Scaling Kit)
Tips available for periodontics, endodontics, scaling, restorative, and cavity preparation
Cruise Control® automatically adjusts frequency amplitude
On/off controlled by rheostat
First ever autoclavable LED lit handpiece – 100,000 LUX
Newtron technology means power, control and reliability
Satelec warranty 1 year / tips 30 days from date of purchase

Scaler



Factory Installed	CAVITRON	\$ 3250
-------------------	----------	---------

Please call 800-800-7202 for warranty and technical support of the Scaler
Uses existing handpiece positionDentsply Cavitron Magnetostrictive Ultrasonic scaler
Blue Zone® for patient comfort
SPS® (Sustained Performance System) constant balance between
scaling efficiency and patient comfort
Handpiece with 330 degree swivel and lavage control (reduces cord drag)
Autoclavable Steri-mate® (infection control)
PowerLINE 1000 30K insert
Frequency 30 kHz
Cavitron warranty 2 years from date of purchase

HP Illumination System, 5 Hole Adapter with Krypton Lamp

Factory Installed 1700-021-CS \$ 695

HP Illumination System, 6 Pin

Factory Installed 1700-019-CS 695



Uses existing handpiece position
 Dual channel outputs individually adjustable from 2VDC to 4VDC
 Accommodates one additional handpiece tubing
 Illumination delay (adjustable between 0-25 seconds) can be set individually for 2 HP
 Automatic activation of HP illumination with the dental unit foot control
 Simple installation into any new or existing dental unit
 Indicator LEDs confirms correct operation of the control module
 Powers all brands of optic handpieces and air-motors
 ISO standard handpiece terminal with stainless steel nut
 Asepsis style handpiece tubing, (4' to 5' installed)
 Power transformer, output 12VAC, 1.5 Amps
 Installation hardware and instructions

HP Illumination Tubing, 5 Hole Adapter with Krypton Lamp

Factory Installed 0017-615-CS \$ 265



Used to replace power optic handpiece tubing
 ISO standard power optic handpiece terminal with stainless steel nut
 5 hole adapter with krypton lamp
 Asepsis style handpiece tubing, (4' to 5' factory installed)

HP Illumination Tubing, 6 Pin

Factory Installed 0017-603-CS \$ 190



Used to replace power optic handpiece tubing
 ISO standard power optic handpiece terminal with stainless steel adapter
 and nut (5-hole Midwest)
 Asepsis style handpiece tubing, (4' to 5' factory installed)

Control Head Accessories

Ancillary Products

Chair Control, Touchpad, Unit-Mounted

Additional Chair Control, Touchpad, Unit-Mounted

Factory Installed

3914-998-CS

N/C

Factory Installed

3914-999-CS

\$ 500

4 programmable preset positions



Assistant's Hygiene Instrumentation

Factory Installed

1408-CS

\$ 320

Canister with ½" drain

Mounting bracket

Drain tubing

Autoclavable standard HVE with tubing, 5'

Autoclavable saliva ejector with tubing, 5'

Standard vacuum holder

Universal holder

Installed on Standard IC and Value control heads (not available on Euro)



Ancillary Shelf

Factory Installed

0010-310

\$ 95

Accommodates additional ancillary products

Designed for easy access to installed devices

Powder coated aluminum

Installed on Standard IC control head



Seven Packages Configured by Forest

To create a Forest package, first start with the delivery system mounting method & control head of your choice. The delivery system is what differentiates one package from another. Purchase the exact package shown or add/delete items to make it your own.

Four Pivot Chair Mount Packages
pivot mounted deliveries with three control head choices.

Three Fixed Chair Mount with Sidebox Packages
sidebox mounted deliveries with three control head choices.

Chair Mount Packages

[Create a Package](#)





WORD OF MOUTH

We feel good in a warm and friendly environment. We tell our friends about a restaurant with great food at a reasonable price in the perfect setting. The same is true with dentistry. We've seen that within a year of upgrading, practices experience up to a 40% increase in new patients.

3900 Patient Chair Powder Coat Paint Color: Bronze

9078 LED Light Powder Coat Paint Color: Bronze


Standard IC Control Powder Coat Paint Color: Bronze | Plastic Finish: Laser Leather

Assistant's Instrumentation Powder Coat Paint Color: Bronze | Plastic Finish: Laser Leather

Plush Duo Color Upholstery: Ultraleather® Doe and Lagoon


Chair Mounting Options

Pivot Chair Mount



Pivot Chair Mounts easily move left to right on large roller bearings to accommodate left and right-handed practitioners. Unlike a Fixed Chair Mount delivery system, the mounting arm comes from underneath the patient rather than over the top, making it less obtrusive and more patient-friendly. Pivot Assistant's Instrumentation, Operatory Lights and Monitor Mounts are popular options to support your style of practice.

Fixed Chair Mount With Sidebox



Fixed Chair Mounts With Sidebox allow for a Cuspidor, Assistant's Instrumentation a Monitor and Light to mount on a sidebox, which is large enough to accommodate a host of ancillary support products, finish the aluminum side covers in one of our nine powder coat paint colors and choose a complementary carbon fiber trim to finish the top cap. If you use a cuspidor, select our beautiful high-gloss ceramic black for a customized and beautiful patient-friendly delivery system.

Pivot Chair Mount Packages

Package P1 | Standard IC Control Head

STANDARD FEATURES:

DELIVERY SYSTEM 4485PI

Automatic control for three handpieces

asepsis 3-way air/water syringe with tubing, 5' handpiece tubing, 4-hole Midwest, asepsis style, 5' pneumatic master control block with diaphragm master on/off

individual handpiece pressure adjustment

individual water coolant adjustment

non-retraction water coolant

air/water coolant adjustment

handpiece flush

handpiece exhaust filter/oil collector

dual handle brake release function

self-contained clean water system with quick disconnect, single bottle, 1½ liter

wet/dry foot control

Flexible control arm with pneumatic brake

vertical travel 28" / horizontal travel 62"

Cable for unit-mounted touchpad, 3900 chair

Internal umbilical, 5' with 3' external from chair to utility center

Contoured utility center

air pressure regulator with filter

master air shut-off valve

Chair adapter, pivot

ASSISTANT'S INSTRUMENTATION,

REAR PIVOT 1095 (option not shown)

Central vacuum accessory package

central vacuum canister with disposable solids collector

asepsis 3-way air/water syringe with tubing, 5'

autoclavable standard HVE lever valve with tubing, 5'

autoclavable saliva ejector lever valve with tubing, 5'

Internal umbilical, 5' with 3' external

Rear chair pivot mount

CUSPIDOR, REAR PIVOT 1580 (option not shown)

One-piece ceramic bowl with touch lever water activator

Automatic cup fill

Timed bowl rinse

Swing away or over patient's lap

Water pressure regulator with filter

Master water shut-off valve

Internal umbilical, 5' with 3' external

Rear chair pivot mount

CUSPIDOR, REAR PIVOT 1590 (option not shown)

Same as 1580 above, with additional features:

Swivel assistant's quad holder bar

Assistant's instrumentation

central vacuum canister with disposable solids collector

asepsis 3-way air/water syringe with tubing, 5'

autoclavable standard HVE lever valve with tubing, 5'

autoclavable saliva ejector lever valve with tubing, 5'

LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching temperature setting 4200K to 5500K

Shadow free illumination

CHAIR, HYDRAULIC 3900

Ultraleather® upholstery with wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad with 4 programmable preset positions

Seat height: minimum 17" maximum 29"



Popular Equipment Choices



ASSISTANT'S
INSTRUMENTATION 1095



CUSPIDOR 1580



CUSPIDOR 1590



LED LIGHT
9078



MONITOR MOUNT
1300RM



OFFICE COMFORT
STOOL 6198



DOCTOR'S COMFORT
STOOL 6176



ASSISTANT'S COMFORT
STOOL 6187

		FACTORY OPTIONS	PART NO	RETAIL	PRICE
Base PACKAGE	\$16,440	DELIVERY SYSTEM WITH STANDARD IC 4485PI* \$ 6160			
4485PI Delivery System		Handpiece illumination system			
3900 Chair		One handpiece (5-hole), accommodates 1 add'l tubing..... 1700-021-CS 695			
9080 Light		Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS 265			
1112-644 Light Post		One handpiece (6-pin), accommodates 1 add'l tubing..... 1700-019-CS 695			
		Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS 190			
As shown PACKAGE	\$16,775	Ancillaries (use existing handpiece position) (see page 6)			
Additionally equipped with:		Write in selection here.....			
1110-468 Swivel Tray, Arm-Mounted		Write in selection here.....			
3914-998-CS Chair Control, Touchpad, Unit-Mounted		Additional handpiece position, wet 1100-189-CS 405			
		Additional handpiece position, air drive only 1100-190-CS 195			
		Assistant's hygiene instrumentation (control head mount) 1408-CS 400			
		Quick disconnect with flow control 1/4" (for water) 0016-454-CS 100			
		Quick disconnect with flow control 3/8" (for air) 0016-456-CS 100			
		Swivel tray, arm-mounted (over control head) 1110-468* 335			
		Chrome accents 9000-998-CS 135			
		Additional umbilical (per foot) 1115-999-CS 30			
		Routing valve kit with water shutoff & regulator/factory installation 1115-315-CS 235			
		Routing valve kit for use with cuspidor/factory installation 1115-316-CS 60			
		ASSISTANT'S INSTRUMENTATION, REAR PIVOT ○ 1095 1495			
		or CUSPIDOR, REAR PIVOT ○ 1580 2950			
		or CUSPIDOR WITH ASSISTANT'S INSTRUMENTATION,			
		REAR PIVOT ○ 1590 3350			
		Quick disconnect with flow control 1/4" (specify air or water) (1095) 0016-455-CS 100			
		Additional autoclavable standard HVE			
		with tubing, 5' & autoclavable oral dry cup 1400-113-CS 85			
		Self-contained clean water system with quick disconnect,			
		single bottle, 1 1/2 liter (1095) 1700-319-CS 265			
		Touchpad, unit-mounted (1095) for 3900 chair 1700-088-CS 500			
		Tray assembly (1095) (see page 6) 1110-680 560			
		Dual articulating quad holder assembly, standard HVE (1095) 0024-227-CS 235			
		Deluxe cuspidor unit-mounted activation (1580/1590) 1100-480-CS 200			
		LIGHT, HALOGEN (mounting options, see pages 30-33) ○ 9080* 2090			
		or LED ○ 9078 3490			
		Light post, pivot (does not include chair adaptor) 1112-644* 500			
		• Monitor post (for top post mount) (does not include chair adaptor) 1331-644 500			
		Monitor mount assembly (monitor screen not included) (see page 34)			
		Ring mount (used with light post) 1300RM 990			
		Top post mount (used with monitor post) 1300TPM 990			
		CHAIR, HYDRAULIC ○ 3900* 7690			
		One chair control included (specify type below) 3914-998-CS* N/C			
		○ Footswitch control			
		○ Touchpad, chair mounted			
		○ Touchpad, unit-mounted			
		Additional chair control (limit 2 on chair, limit 2 on unit) 3914-999-CS 500			
		○ Footswitch control			
		○ Touchpad, chair mounted			
		○ Touchpad, unit-mounted			
		Additional chair options (See page 25 for pricing, write selection here)			
		STOOL, DR'S, ECONOMY/DELUXE/COMFORT 6156C/6159C/6176 485/825/825			
		STOOL, ASST'S, ECONOMY/DELUXE 6157C/6160C 625/925			
		STOOL, ASST'S WITH BACK, DELUXE/COMFORT 6161C/6187 1100/1100			
		STOOL, OFFICE, COMFORT 6198 975			
		TOTAL			

* ITEM included in AS SHOWN PACKAGE price

• Choose only one post when sold with unit

■ Often overlooked

Equipment dimensions, see page 66

Shipping weight estimate:

Delivery unit & chair 450 lb / 204 kg

Base package 580 lb / 264 kg

As shown package 630 lb / 286 kg

Pivot Chair Mount Packages

Package P2 | Euro Control Head

STANDARD FEATURES:

DELIVERY SYSTEM 4987PI

Automatic control for four handpieces

- asepsis 3-way air/water syringe with tubing, 3'
- handpiece tubing, 4-hole Midwest, asepsis style, 3'
- pneumatic master control block with diaphragm master on/off
- individual handpiece pressure adjustment
- individual water coolant adjustment
- non-retraction water coolant
- air/water coolant adjustment
- handpiece flush
- handpiece exhaust filter/oil collector
- dual handle brake release function
- swivel instrument tray
- self-contained clean water system with quick disconnect, single bottle, 1½ liter
- wet/dry foot control

Flexible control arm with pneumatic brake
vertical travel 28" / horizontal travel 62"

Cable for unit-mounted touchpad, 3900 chair

Internal umbilical, 5' with 3' external from chair to utility center
Contoured utility center

- air pressure regulator with filter
- master air shut-off valve

Chair adapter, pivot

CUSPIDOR, REAR PIVOT 1590

Same as 1580 below, with additional features:

Swivel assistant's quad holder bar

Assistant's instrumentation

- central vacuum canister with disposable solids collector
- asepsis 3-way air/water syringe with tubing, 5'
- autoclavable standard HVE lever valve with tubing, 5'
- autoclavable saliva ejector lever valve with tubing, 5'

CUSPIDOR, REAR PIVOT 1580 (option not shown)

One-piece ceramic bowl with touch lever water activator

Automatic cup fill

Timed bowl rinse

Swing away or over patient's lap

Water pressure regulator with filter

Master water shut-off valve

Internal umbilical, 5' with 3' external

Rear chair pivot mount

ASSISTANT'S INSTRUMENTATION, REAR PIVOT 1095 (option not shown)

Central vacuum accessory package

- central vacuum canister with disposable solids collector
- asepsis 3-way air/water syringe with tubing, 5'
- autoclavable standard HVE lever valve with tubing, 5'
- autoclavable saliva ejector lever valve with tubing, 5'

Internal umbilical, 5' with 3' external

Rear chair pivot mount

LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching temperature setting 4200K to 5500K

Shadow free illumination

CHAIR, HYDRAULIC 3900

Ultraleather® upholstery with wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad with 4 programmable preset positions

Seat height: minimum 17" maximum 29"



Popular Equipment Choices



ASSISTANT'S
INSTRUMENTATION 1095



CUSPIDOR 1580



ADD'L HVE POSITION
WITH ORAL DRY CUP
1400-113-CS



LED LIGHT
9078



DOCTOR'S COMFORT
STOOL 6176



ASSISTANT'S COMFORT
STOOL 6187



DOCTOR'S DELUXE
STOOL 6159C



ASSISTANT'S DELUXE
STOOL 6161C

		FACTORY OPTIONS	PART NO	RETAIL	PRICE
Base PACKAGE	\$17,440	DELIVERY SYSTEM WITH EURO Handpiece illumination system One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS Add'l handpiece illumination tubing (5-hole) 0017-615-CS One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS Add'l handpiece illumination tubing (6-pin) 0017-603-CS	4987PI*	\$ 7160	
4987PI Delivery System 3900 Chair 9080 Light 1112-644 Light Post					
As shown PACKAGE	\$21,780	Ancillaries use existing handpiece position (see page 6) Write in selection here..... Write in selection here..... Additional handpiece position, wet 1100-189-CS Additional handpiece position, air drive only 1100-190-CS Quick disconnect with flow control 1/4" (for water) 0016-454-CS Quick disconnect with flow control 3/8" (for air) 0016-456-CS Chrome accents 9000-998-CS Additional umbilical (per foot) 1115-999-CS Routing valve kit with water shutoff & regulator/factory installation .. 1115-315-CS Routing valve kit for use with cuspidor/factory installation 1115-316-CS ASSISTANT'S INSTRUMENTATION, REAR PIVOT ○ 1095 or CUSPIDOR, REAR PIVOT ○ 1580 or CUSPIDOR WITH ASSISTANT'S INSTRUMENTATION, REAR PIVOT . ○ 1590* Quick disconnect with flow control 1/4" (specify air or water) (1095). 0016-455-CS Additional autoclavable standard HVE with tubing, 5' & autoclavable oral dry cup 1400-113-CS Self-contained clean water system with quick disconnect, single bottle, 1 1/2 liter (1095) 1700-319-CS Touchpad, unit-mounted (1095) for 3900 chair 1700-088-CS Tray assembly (1095) (see page 6!) 1110-680 Dual articulating quad holder assembly, standard HVE (1095) 0024-227-CS Deluxe cuspidor unit-mounted activation (1580/1590) 1100-480-CS LIGHT, HALOGEN (mounting options, see pages 30-33) ○ 9080* or LED ○ 9078 Light post, pivot (does not include chair adaptor) 1112-644* • Monitor post (for top post mount) (does not include chair adaptor) 1331-644 Monitor mount assembly (monitor screen not included) (see page 34) Ring mount (used with light post) 1300RM* Top post mount (used with monitor post) 1300TPM CHAIR, HYDRAULIC ○ 3900* ● One chair control included (specify type below) 3914-998-CS* ○ Footswitch control ○ Touchpad, chair mounted ○ Touchpad, unit-mounted ● Additional chair control (limit 2 on chair, limit 2 on unit) 3914-999-CS ○ Footswitch control ○ Touchpad, chair mounted ○ Touchpad, unit-mounted Additional chair options (See page 25 for pricing, write selection here) STOOL, DR'S, ECONOMY/DELUXE/COMFORT 6156C/6159C/6176 STOOL, ASST'S, ECONOMY/DELUXE 6157C/6160C STOOL, ASST'S WITH BACK, DELUXE/COMFORT 6161C/6187 STOOL, OFFICE, COMFORT 6198 TOTAL	485/825/825		
Additionally equipped with: 1590 Deluxe Cuspidor 1300RM Monitor Assembly 3914-998-CS Chair Control, Chair Mounted					
DESIGNER FRIENDLY® SURFACE TREATMENTS					
STEP 1: Select unit plastic finishes:					
<input type="radio"/> Greywood <input type="radio"/> Laser Leather <input type="radio"/> Lemans <input type="radio"/> Forest Grey <input type="radio"/> White					
STEP 2: Select powder coat paint for unit & light:					
<input type="radio"/> Arctic Silver <input type="radio"/> Bronze <input type="radio"/> Charcoal <input type="radio"/> Copper <input type="radio"/> Forest Grey <input type="radio"/> Midnight <input type="radio"/> Pewter <input type="radio"/> Shadow <input type="radio"/> White					
STEP 3: Select standard powder coat paint for chair:					
<input type="radio"/> Forest Grey <input type="radio"/> Shadow					
OR upgrade chair base/back powder coat paint (add \$2,250):					
<input type="radio"/> Arctic Silver <input type="radio"/> Bronze <input type="radio"/> Charcoal <input type="radio"/> Copper <input type="radio"/> Midnight <input type="radio"/> Pewter <input type="radio"/> White					
STEP 4: Select Ultraleather® upholstery color choices for chair & deluxe /economy stools:					
<input type="radio"/> Apricot <input type="radio"/> Baltic <input type="radio"/> Cave <input type="radio"/> Charcoal <input type="radio"/> Diplomat Blue <input type="radio"/> Doe <input type="radio"/> Dove Grey <input type="radio"/> Eucalyptus <input type="radio"/> Lemon <input type="radio"/> Papyrus <input type="radio"/> Parrot <input type="radio"/> Raven Wing <input type="radio"/> Red <input type="radio"/> Riviera <input type="radio"/> Stone Pearlized colors: <input type="radio"/> Abalone <input type="radio"/> Midnight <input type="radio"/> Pewter					
OR select Special Color Ultraleather®					
<input type="radio"/> Chair 3907-999-SP add \$280 <input type="radio"/> Deluxe/Economy Stool 6007-999-SP add \$100					
STEP 5: Upholstery Selections:					
Plush Selections:					
<input type="radio"/> Standard Back, Wide/Thin 3907-PL-WB N/C <input type="radio"/> Narrow Back 3907-PL-NB N/C <input type="radio"/> Narrow Back with Slings 3907-PL-NS add \$400					
Plush Additional Options:					
<input type="radio"/> Contrast/Deviated Stitching 3907-999-DS add \$280 <input type="radio"/> Dual Color 3907-999-DC add \$280 <input type="radio"/> Dual Color with Contrast/Deviated Stitching Call Customer Service for information (1.800.423.3555) <input type="radio"/> Massage Option 3900-100 add \$540					
Asepsis Selections:					
<input type="radio"/> Standard Back, Wide/Thin 3907-AS-WB N/C <input type="radio"/> Narrow Back 3907-AS-NB N/C <input type="radio"/> Narrow Back with Slings 3907-AS-NS add \$400					
Asepsis Additional Options:					
<input type="radio"/> Massage Option 3900-100 add \$540					
Memory Comfort Selections:					
<input type="radio"/> Wide/Thin 3907-MC-WB \$800 <input type="radio"/> Dual Color 3907-MC-DC \$1080					
Equipment dimensions, see page 66					
Shipping weight estimate:					
Delivery unit & chair 450 lb / 204 kg					
Base package 580 lb / 264 kg					
As shown package 630 lb / 286 kg					

* ITEM included in AS SHOWN PACKAGE price

• Choose only one post when sold with unit

● Often overlooked

Pivot Chair Mount Packages

Package P3 | Value Control Head

STANDARD FEATURES:

DELIVERY SYSTEM 4195PI

Automatic control for three handpieces

asepsis 3-way air/water syringe with tubing, 5'

handpiece tubing, 4-hole Midwest, asepsis style, 5'

individual pinch valve control blocks

master on/off

individual handpiece pressure adjustment

individual water coolant adjustment

non-retraction water coolant

air coolant adjustment

handpiece flush

handpiece exhaust filter/oil collector

dual handle brake release function

self-contained clean water system with quick disconnect,
single bottle, 1½ liter

wet/dry foot control

Flexible control arm with pneumatic brake

vertical travel 28" / horizontal travel 62"

Cable for unit-mounted touchpad, 3900 chair

Internal umbilical, 5' with 3' external from chair to utility center

Contoured utility center

air pressure regulator with filter

master air shut-off valve

Chair adapter, pivot

ASSISTANT'S INSTRUMENTATION, REAR PIVOT 1095

Central vacuum accessory package

central vacuum canister with disposable solids collector

asepsis 3-way air/water syringe with tubing, 5'

autoclavable standard HVE lever valve with tubing, 5'

autoclavable saliva ejector lever valve with tubing, 5'

Internal umbilical, 5' with 3' external

Rear chair pivot mount

CUSPIDOR, REAR PIVOT 1580 (option not shown)

One-piece ceramic bowl with touch lever water activator

Automatic cup fill

Timed bowl rinse

Swing away or over patient's lap

Water pressure regulator with filter

Master water shut-off valve

Internal umbilical, 5' with 3' external

Rear chair pivot mount

CUSPIDOR, REAR PIVOT 1590 (option not shown)

Same as 1580 above, with additional features:

Swivel assistant's quad holder bar

Assistant's instrumentation

central vacuum canister with disposable solids collector

asepsis 3-way air/water syringe with tubing, 5'

autoclavable standard HVE lever valve with tubing, 5'

autoclavable saliva ejector lever valve with tubing, 5'

LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching temperature setting 4200K to 5500K

Shadow free illumination

CHAIR, HYDRAULIC 3900

Ultraleather® upholstery with wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad with 4 programmable preset positions

Seat height: minimum 17" maximum 29"



Popular Equipment Choices



CUSPIDOR 1580



CUSPIDOR 1590



ADD'L HVE POSITION
WITH ORAL DRY CUP
1400-113-CS



LED LIGHT
9078



MONITOR MOUNT
1300RM



DOCTOR'S ECONOMY
STOOL 6156C



ASST'S ECONOMY
STOOL 6157C



DOCTOR'S COMFORT
STOOL 6176

		FACTORY OPTIONS	PART NO	RETAIL	PRICE
Base PACKAGE	\$15,240	DELIVERY SYSTEM WITH VALUE Handpiece illumination system One handpiece (5-hole), accommodates 1 add'l tubing..... Add'l handpiece illumination tubing (5-hole)..... One handpiece (6-pin), accommodates 1 add'l tubing..... Add'l handpiece illumination tubing (6-pin).....	4195PI*	\$ 4960	
4195PI Delivery System		1700-021-CS	695		
3900 Chair		0017-615-CS	265		
9080 Light		1700-019-CS	695		
1112-644 Light Post		0017-603-CS	190		
As shown PACKAGE	\$17,070	• Ancillaries use existing handpiece position (see page 6) Write in selection here..... Write in selection here..... Additional handpiece position, wet	1100-189-CS	405	
Additionally equipped with:		Additional handpiece position, air drive only	1100-190-CS	195	
1110-468 Swivel Tray, Arm-Mounted		Assistant's hygiene instrumentation (control head mount)	1408-CS	400	
1095 Assistant's Instrumentation		Quick disconnect with flow control 1/4" (for water)	0016-454-CS	100	
3914-998-CS Chair Control, Chair Mounted		Quick disconnect with flow control 3/8" (for air)	0016-456-CS	100	
DESIGNER FRIENDLY® SURFACE TREATMENTS					
STEP 1: Select unit plastic finishes:		Swivel tray, arm-mounted (over control head)	1110-468*	335	
<input type="radio"/> Greywood <input type="radio"/> Laser Leather <input type="radio"/> Lemans		Chrome accents	9000-998-CS	135	
<input type="radio"/> Forest Grey <input type="radio"/> White		Additional umbilical (per foot).....	1115-999-CS	30	
STEP 2: Select powder coat paint for unit & light:		Routing valve kit with water shutoff & regulator/factory installation	1115-315-CS	235	
<input type="radio"/> Arctic Silver <input type="radio"/> Bronze <input type="radio"/> Charcoal		Routing valve kit for use with cuspidor/factory installation.....	1115-316-CS	60	
<input type="radio"/> Copper <input type="radio"/> Forest Grey <input type="radio"/> Midnight		ASSISTANT'S INSTRUMENTATION, REAR PIVOT	1095*	1495	
<input type="radio"/> Pewter <input type="radio"/> Shadow <input type="radio"/> White		OR CUSPIDOR, REAR PIVOT	1580	2950	
STEP 3: Select standard powder coat paint for chair:		OR CUSPIDOR WITH ASSISTANT'S INSTRUMENTATION, REAR PIVOT	1590	3350	
<input type="radio"/> Forest Grey <input type="radio"/> Shadow		Quick disconnect with flow control 1/4" (specify air or water) (1095)	0016-455-CS	100	
OR upgrade chair base/back powder coat paint (add \$2,250):		Additional autoclavable standard HVE with tubing, 5' & autoclavable oral dry cup.....	1400-113-CS	85	
<input type="radio"/> Arctic Silver <input type="radio"/> Bronze <input type="radio"/> Charcoal		Self-contained clean water system with quick disconnect, single bottle, 1 1/2 liter (1095).....	1700-319-CS	265	
<input type="radio"/> Copper <input type="radio"/> Midnight <input type="radio"/> Pewter <input type="radio"/> White		Touchpad, unit-mounted (1095) for 3900 chair.....	1700-088-CS	500	
STEP 4: Select Ultraleather® upholstery color choices for chair & deluxe /economy stools:		Tray assembly (1095) (see page 6).....	1110-680	560	
<input type="radio"/> Apricot <input type="radio"/> Baltic <input type="radio"/> Cave <input type="radio"/> Charcoal		Dual articulating quad holder assembly, standard HVE (1095).....	0024-227-CS	235	
<input type="radio"/> Diplomat Blue <input type="radio"/> Doe <input type="radio"/> Dove Grey		Deluxe cuspidor unit-mounted activation (1580/1590).....	1100-480-CS	200	
<input type="radio"/> Eucalyptus <input type="radio"/> Lemon <input type="radio"/> Papyrus <input type="radio"/> Parrot		LIGHT, HALOGEN (mounting options, see pages 30-33).....	9080*	2090	
<input type="radio"/> Raven Wing <input type="radio"/> Red <input type="radio"/> Riviera <input type="radio"/> Stone		OR LED	9078	3490	
Pearlized colors: <input type="radio"/> Abalone <input type="radio"/> Midnight <input type="radio"/> Pewter		• Light post, pivot (does not include chair adaptor)	1112-644*	500	
OR select Special Color Ultraleather®		• Monitor post (for top post mount) (does not include chair adaptor).....	1331-644	500	
<input type="radio"/> Chair <input type="radio"/> 3907-999-SP add \$280		Monitor mount assembly (monitor screen not included) (see page 34).....			
<input type="radio"/> Deluxe/Economy Stool <input type="radio"/> 6007-999-SP add \$100		Ring mount (used with light post)	1300RM	990	
STEP 5: Upholstery Selections:		Top post mount (used with monitor post)	1300TPM	990	
Plush Selections:		CHAIR, HYDRAULIC	3900*	7690	
<input type="radio"/> Standard Back, Wide/Thin	3907-PL-WB N/C	• One chair control included (specify type below)	3914-998-CS*	N/C	
<input type="radio"/> Narrow Back	3907-PL-NB N/C	<input type="radio"/> Footswitch control			
<input type="radio"/> Narrow Back with Slings	3907-PL-NS add \$400	<input type="radio"/> Touchpad, chair mounted			
Plush Additional Options:		<input type="radio"/> Touchpad, unit-mounted			
<input type="radio"/> Contrast/Deviated Stitching	3907-999-DS add \$280	• Additional chair control (limit 2 on chair, limit 2 on unit)	3914-999-CS	500	
<input type="radio"/> Dual Color	3907-999-DC add \$280	<input type="radio"/> Footswitch control			
<input type="radio"/> Dual Color with		<input type="radio"/> Touchpad, chair mounted			
Contrast/Deviated Stitching	Call Customer Service for information (1.800.423.3555)	<input type="radio"/> Touchpad, unit-mounted			
<input type="radio"/> Massage Option	3900-100 add \$540	Additional chair options (See page 25 for pricing, write selection here)			
Asepsis Selections:		STOOL, DR'S, ECONOMY/DELUXE/COMFORT	6156C/6159C/6176	485/825/825	
<input type="radio"/> Standard Back, Wide/Thin	3907-AS-WB N/C	STOOL, ASST'S, ECONOMY/DELUXE	6157C/6160C	625/925	
<input type="radio"/> Narrow Back	3907-AS-NB N/C	STOOL, ASST'S WITH BACK, DELUXE/COMFORT	6161C/6187	1100/1100	
<input type="radio"/> Narrow Back with Slings	3907-AS-NS add \$400	STOOL, OFFICE, COMFORT	6198	975	
Asepsis Additional Options:		* ITEM included in AS SHOWN PACKAGE price			
<input type="radio"/> Massage Option	3900-100 add \$540	• Choose only one post when sold with unit			
Memory Comfort Selections:		• Often overlooked			
<input type="radio"/> Wide/Thin	3907-MC-WB \$800	TOTAL			
<input type="radio"/> Dual Color	3907-MC-DC \$1080				

Equipment dimensions, see page 66

Shipping weight estimate:

Delivery unit & chair	445 lb / 203 kg
Base package	556 lb / 257 kg
As shown package	585 lb / 266 kg

Pivot Chair Mount Packages

Package P4 | Rear Standard IC Control Head

STANDARD FEATURES:

DELIVERY SYSTEM 1500PI

Automatic control for three handpieces

central vacuum canister with disposable solids collector

asepsis 3-way air/water syringe with tubing, 5'

autoclavable standard HVE lever valve with tubing, 5'

autoclavable saliva ejector lever valve with tubing, 5'

handpiece tubing, 4-hole Midwest, asepsis style, 5'

pneumatic master control block with diaphragm

master on/off

individual handpiece pressure adjustment

individual water coolant adjustment

non-retraction water coolant

air coolant adjustment

handpiece flush

handpiece exhaust filter/oil collector

self-contained clean water system with quick disconnect,

single bottle, 1½ liter

wet/dry foot control

Cable for unit-mounted touchpad, 3900 chair

Internal umbilical, 5' with 3' external from chair to utility center

Contoured utility center

air pressure regulator with filter

master air shut-off valve

Rear chair pivot mount

LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching temperature setting 4200K to 5500K

Shadow free illumination

CHAIR, HYDRAULIC 3900

Ultraleather® upholstery with wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad with 4 programmable preset positions

Seat height: minimum 17" maximum 29"



Popular Equipment Choices



SWIVEL TRAY,
ARM MOUNTED
1110-731



MONITOR MOUNT
1300RM



ADD'L HVE POSITION
WITH ORAL DRY CUP
1400-113-CS



LED LIGHT
9078



OFFICE COMFORT
STOOL 6198



DOCTOR'S COMFORT
STOOL 6176



DOCTOR'S DELUXE
STOOL 6159C



DOCTOR'S ECONOMY
STOOL 6156C

	FACTORY OPTIONS	PART NO	RETAIL	PRICE
Base PACKAGE	\$15,940	DELIVERY SYSTEM REAR PIVOT WITH STANDARD IC 1500PI*	\$ 5095	
1500PI Delivery System		Handpiece illumination system		
3900 Chair		One handpiece (5-hole), accommodates 1 add'l tubing..... 1700-021-CS	695	
9080 Light		Add'l handpiece illumination tubing (5-hole) 0017-615-CS	265	
1112-644-3 Light Post		One handpiece (6-pin), accommodates 1 add'l tubing..... 1700-019-CS	695	
		Add'l handpiece illumination tubing (6-pin) 0017-603-CS	190	
As shown PACKAGE	\$16,090	Ancillaries use existing handpiece position (see page 6)		
0024-054 Assistant's Quad Holder Bar		Write in selection here.....		
3914-998-CS Chair Control, Chair Mounted		Write in selection here.....		
DESIGNER FRIENDLY® SURFACE TREATMENTS				
STEP 1: Select unit plastic finishes:				
<input type="radio"/> Greywood <input type="radio"/> Laser Leather <input type="radio"/> Lemans				
<input type="radio"/> Forest Grey <input type="radio"/> White				
STEP 2: Select powder coat paint for unit & light:				
<input type="radio"/> Arctic Silver <input type="radio"/> Bronze <input type="radio"/> Charcoal				
<input type="radio"/> Copper <input type="radio"/> Forest Grey <input type="radio"/> Midnight				
<input type="radio"/> Pewter <input type="radio"/> Shadow <input type="radio"/> White				
STEP 3: Select standard powder coat paint for chair:				
<input type="radio"/> Forest Grey <input type="radio"/> Shadow				
OR upgrade chair base/back powder coat paint (add \$2,250):				
<input type="radio"/> Arctic Silver <input type="radio"/> Bronze <input type="radio"/> Charcoal				
<input type="radio"/> Copper <input type="radio"/> Midnight <input type="radio"/> Pewter <input type="radio"/> White				
STEP 4: Select Ultraleather® upholstery color choices for chair & deluxe /economy stools:				
<input type="radio"/> Apricot <input type="radio"/> Baltic <input type="radio"/> Cave <input type="radio"/> Charcoal				
<input type="radio"/> Diplomat Blue <input type="radio"/> Doe <input type="radio"/> Dove Grey				
<input type="radio"/> Eucalyptus <input type="radio"/> Lemon <input type="radio"/> Papyrus <input type="radio"/> Parrot				
<input type="radio"/> Raven Wing <input type="radio"/> Red <input type="radio"/> Riviera <input type="radio"/> Stone				
Pearlized colors: <input type="radio"/> Abalone <input type="radio"/> Midnight <input type="radio"/> Pewter				
OR select Special Color Ultraleather®				
<input type="radio"/> Chair 3907-999-SP add \$280				
<input type="radio"/> Deluxe/Economy Stool 6007-999-SP add \$100				
STEP 5: Upholstery Selections:				
Plush Selections:				
<input type="radio"/> Standard Back, Wide/Thin 3907-PL-WB N/C				
<input type="radio"/> Narrow Back 3907-PL-NB N/C				
<input type="radio"/> Narrow Back with Slings 3907-PL-NS add \$400				
Plush Additional Options:				
<input type="radio"/> Contrast/Deviated Stitching 3907-999-DS add \$280				
<input type="radio"/> Dual Color 3907-999-DC add \$280				
<input type="radio"/> Dual Color with Contrast/Deviated Stitching Call Customer Service for information (1.800.423.3555)				
<input type="radio"/> Massage Option 3900-100 add \$540				
Asepsis Selections:				
<input type="radio"/> Standard Back, Wide/Thin 3907-AS-WB N/C				
<input type="radio"/> Narrow Back 3907-AS-NB N/C				
<input type="radio"/> Narrow Back with Slings 3907-AS-NS add \$400				
Asepsis Additional Options:				
<input type="radio"/> Massage Option 3900-100 add \$540				
Memory Comfort Selections:				
<input type="radio"/> Wide/Thin 3907-MC-WB \$800				
<input type="radio"/> Dual Color 3907-MC-DC \$1080				
Equipment dimensions, see page 66				
Shipping weight estimate:				
Delivery unit & chair 445 lb / 203 kg				
Base package 533 lb / 242 kg				
As shown package 538 lb / 245 kg				

* ITEM included in AS SHOWN PACKAGE price

Often overlooked

TOTAL

Fixed Chair Mount With Sidebox Packages

Package S1 | Standard IC Control Head

STANDARD FEATURES:

DELIVERY SYSTEM 4485SI

Automatic control for three handpieces

asepsis 3-way air/water syringe with tubing, 5' handpiece tubing, 4-hole Midwest, asepsis style, 5' pneumatic master control block with diaphragm master on/off

individual handpiece pressure adjustment
individual water coolant adjustment

non-retraction water coolant

air coolant adjustment

handpiece flush

handpiece exhaust filter/oil collector
single handle brake release function
wet/dry foot control

Flexible control arm with pneumatic brake
vertical travel 28" / horizontal travel 62"

Sidebox with self-contained clean water system
with quick disconnect, single bottle, 1½ liter

Water quick disconnect ¼" with flow control

Air quick disconnect ⅜" with flow control

Toggle valve for city/bottle water routing

Telescoping arm with swivel quad holder bar

Central vacuum accessory package

central vacuum canister with disposable solids collector

asepsis 3-way air/water syringe with tubing, 5'
autoclavable standard HVE lever valve with tubing, 5'

autoclavable saliva ejector lever valve with tubing, 5'

Cable for unit-mounted touchpad, 3900 chair

Internal umbilical, 5' with 3' external from chair to utility center

Contoured utility center

air/water pressure regulators with filter

master air/water shut-off valves

Chair adapter, fixed

CUSPIDOR, DELUXE 1577 (option not shown)

One-piece ceramic bowl with touch lever water activator

Automatic cup fill

Timed bowl rinse

CUSPIDOR, BASIC 1679 (option not shown)

Removable ceramic bowl

Manual cup fill

Timed bowl rinse

LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching temperature setting 4200K to 5500K

Shadow free illumination

CHAIR, HYDRAULIC 3900

Ultraleather® upholstery with wide back for elbow support
Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad with 4 programmable preset positions

Seat height: minimum 17" maximum 29"



Popular Equipment Choices



CUSPIDOR 1577
(Sidebox comes with delivery)



CUSPIDOR 1679
(Sidebox comes with delivery)



MONITOR MOUNT
1300RM



LED LIGHT
9078



DOCTOR'S COMFORT
STOOL 6176



ASSISTANT'S COMFORT
STOOL 6187



DOCTOR'S DELUXE
STOOL 6159C



ASSISTANT'S DELUXE
STOOL 6161C

	FACTORY OPTIONS	PART NO	RETAIL	PRICE
Base PACKAGE	\$17,420	DELIVERY SYSTEM WITH STANDARD IC	4485SI*	\$ 7440
4485SI Delivery System	Handpiece illumination system			
3900 Chair	One handpiece (5-hole), accommodates 1 add'l tubing	1700-021-CS	695	
9080 Light	Add'l handpiece illumination tubing (5-hole)	0017-615-CS	265	
1110-900 Light Post	One handpiece (6-pin), accommodates 1 add'l tubing	1700-019-CS	695	
	Add'l handpiece illumination tubing (6-pin)	0017-603-CS	190	
As shown PACKAGE	\$17,755	Ancillaries use existing handpiece position (see page 6)		
Additionally equipped with:	Write in selection here.....			
1110-468 Swivel Tray, Arm-Mounted	Write in selection here.....			
3914-998-CS Chair Control, Touchpad, Unit-Mounted	Additional handpiece position, wet	1100-189-CS	405	
	Additional handpiece position	1100-190-CS	195	
	Swivel tray, arm-mounted option (over control head)	1110-468*	335	
	Less sidebox telescoping arm and assistant's instrumentation	1400-000-CS	-350	
	Chrome accents	9000-998-CS	135	
	Additional autoclavable standard HVE			
	with tubing, 5' & autoclavable oral dry cup	1400-113-CS	85	
	Additional umbilical (per foot)	1115-999-CS	30	
	Self-contained dual bottle clean water system			
	with quick disconnect in lieu of single	1700-365-CS	150	
	CUSPIDOR (sold only with delivery)			
	DELUXE	○ 1577	1450	
	or BASIC	○ 1679	995	
	Deluxe cuspidor unit-mounted activation (1577)	1100-480-CS	200	
	LIGHT, HALOGEN (mounting options, see pages 30-33)	○ 9080*	2090	
	or LED	○ 9078	3490	
	• Light post, 44"	1110-900*	200	
	• Monitor post (for top post mount)	1331-998	200	
	Monitor mount assembly (monitor screen not included) (see page 34)			
	Ring mount (used with light post)	1300RM	990	
	Top post mount (used with monitor post)	1300TPM	990	
	CHAIR, HYDRAULIC	○ 3900*	7690	
	• One chair control included (specify type below)	3914-998-CS*	N/C	
	○ Footswitch control			
	○ Touchpad, chair mounted			
	○ Touchpad, unit-mounted			
	• Additional chair control (limit 2 on chair, limit 2 on unit)	3914-999-CS	500	
	○ Footswitch control			
	○ Touchpad, chair mounted			
	○ Touchpad, unit-mounted			
	Additional chair options (See page 25 for pricing, write selection here)			
	STOOL, DR'S, ECONOMY/DELUXE/COMFORT	6156C/6159C/6176	485/825/825	
	STOOL, ASST'S, ECONOMY/DELUXE	6157C/6160C	625/925	
	STOOL, ASST'S WITH BACK, DELUXE/COMFORT	6161C/6187	1100/1100	
	STOOL, OFFICE, COMFORT	6198	975	
	TOTAL			
* ITEM included in AS SHOWN PACKAGE price				
• Choose only one post when sold with unit				
• Often overlooked				
Equipment dimensions, see page 67				
Shipping weight estimate:				
Delivery unit & chair	440 lb / 202 kg			
Base package	556 lb / 253 kg			
As shown package	561 lb / 255 kg			

Fixed Chair Mount With Sidebox Packages

Package S2 | Euro Control Head

STANDARD FEATURES:

DELIVERY SYSTEM 4987SI

Automatic control for four handpieces

asepsis 3-way air/water syringe with tubing, 3'
handpiece tubing, 4-hole Midwest, asepsis style, 3'
pneumatic master control block with diaphragm
master on/off

individual handpiece pressure adjustment

individual water coolant adjustment

non-retraction water coolant

air coolant adjustment

handpiece flush

handpiece exhaust filter/oil collector

dual handle brake release function

swivel instrument tray

wet/dry foot control

Flexible control arm with pneumatic brake

vertical travel 28" / horizontal travel 62"

Sidebox with self-contained clean water system

with quick disconnect, single bottle, 1½ liter

Water quick disconnect ¼" with flow control

Air quick disconnect ⅜" with flow control

Toggle valve for city/bottle water routing

Telescoping arm with swivel quad holder bar

Central vacuum accessory package

central vacuum canister with disposable solids collector

asepsis 3-way air/water syringe with tubing, 5'

autoclavable standard HVE lever valve with tubing, 5'

autoclavable saliva ejector lever valve with tubing, 5'

Cable for unit-mounted touchpad, 3900 chair

Internal umbilical, 5' with 3' external from chair to utility center

Contoured utility center

air/water pressure regulators with filter

master air/water shut-off valves

Chair adapter, fixed

CUSPIDOR, DELUXE 1577 (option not shown)

One-piece ceramic bowl with touch lever water activator

Automatic cup fill

Timed bowl rinse

CUSPIDOR, BASIC 1679 (option not shown)

Removable ceramic bowl

Manual cup fill

Timed bowl rinse

LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching

temperature setting 4200K to 5500K

Shadow free illumination

CHAIR, HYDRAULIC 3900

Ultraleather® upholstery with wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad with 4 programmable preset positions

Seat height: minimum 17" maximum 29"



Popular Equipment Choices



CUSPIDOR 1577
(Sidebox comes with delivery)



CUSPIDOR 1679
(Sidebox comes with delivery)



MONITOR MOUNT
1300RM



LED LIGHT
9078



DOCTOR'S COMFORT
STOOL 6176



ASSISTANT'S COMFORT
STOOL 6187



DOCTOR'S DELUXE
STOOL 6159C



ASSISTANT'S DELUXE
STOOL 6161C

		FACTORY OPTIONS	PART NO	RETAIL	PRICE
Base PACKAGE	\$18,420	DELIVERY SYSTEM WITH EURO 4987SI* \$ 8440			
4987SI Delivery System		Handpiece illumination system			
3900 Chair		One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS	695		
9080 Light		Add'l handpiece illumination tubing (5-hole) 0017-615-CS	265		
1110-900 Light Post		One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS	695		
		Add'l handpiece illumination tubing (6-pin) 0017-603-CS	190		
As shown PACKAGE	\$18,420	■ Ancillaries use existing handpiece position (see page 6)			
3914-998-CS Chair Control, Chair Mounted		Write in selection here.....			
		Write in selection here.....			
		Additional handpiece position, wet 1100-189-CS	405		
		Additional handpiece position, air drive only 1100-190-CS	195		
		Less sidebox telescoping arm and assistant's instrumentation 1400-000-CS	-350		
		Chrome accents 9000-998-CS	135		
		Additional autoclavable standard HVE			
		with tubing, 5' & autoclavable oral dry cup 1400-113-CS	85		
		Additional umbilical (per foot) 1115-999-CS	30		
		Self-contained dual bottle clean water system			
		with quick disconnect in lieu of single 1700-365-CS	150		
		CUSPIDOR (sold only with delivery)			
		DELUXE ○ 1577	1450		
		or BASIC ○ 1679	995		
		Deluxe cuspidor unit-mounted activation (1577) 1100-480-CS	200		
		LIGHT, HALOGEN (mounting options, see pages 30-33) ○ 9080*	2090		
		or LED ○ 9078	3490		
		• Light post, 44" 1110-900*	200		
		• Monitor post (for top post mount) 1331-998	200		
		Monitor mount assembly (monitor screen not included) (see page 34)			
		Ring mount (used with light post) 1300RM	990		
		Top post mount (used with monitor post) 1300TPM	990		
		CHAIR, HYDRAULIC 3900*	7690		
		● One chair control included (specify type below) 3914-998-CS*	N/C		
		○ Footswitch control			
		○ Touchpad, chair mounted			
		○ Touchpad, unit-mounted			
		● Additional chair control (limit 2 on chair, limit 2 on unit) 3914-999-CS	500		
		○ Footswitch control			
		○ Touchpad, chair mounted			
		○ Touchpad, unit-mounted			
		Additional chair options (See page 25 for pricing, write selection here)			
		STOOL, DR'S, ECONOMY/DELUXE/COMFORT 6156C/6159C/6176	485/825/825		
		STOOL, ASST'S, ECONOMY/DELUXE 6157C/6160C	625/925		
		STOOL, ASST'S WITH BACK, DELUXE/COMFORT 6161C/6187	1100/1100		
		STOOL, OFFICE, COMFORT 6198	975		
		* ITEM included in AS SHOWN PACKAGE price			
		• Choose only one post when sold with unit			
		● Often overlooked			
		TOTAL			
Equipment dimensions, see page 67					
Shipping weight estimate:					
Delivery unit & chair	445 lb / 203 kg				
Base package	551 lb / 251 kg				
As shown package	551 lb / 251 kg				

Fixed Chair Mount With Sidebox Packages

Package S3 | Value Control Head

STANDARD FEATURES:

DELIVERY SYSTEM 4195SI

Automatic control for three handpieces

asepsis 3-way air/water syringe with tubing, 5'

handpiece tubing, 4-hole Midwest, asepsis style, 5'

individual pinch valve control blocks

master on/off

individual handpiece pressure adjustment

individual water coolant adjustment

non-retraction water coolant adjustment

air coolant adjustment

handpiece flush

handpiece exhaust filter/oil collector

single handle brake release function

wet/dry foot control

Flexible control arm with pneumatic brake

vertical travel 28" / horizontal travel 62"

Sidebox with self-contained clean water system

with quick disconnect, single bottle, 1½ liter

Water quick disconnect ¼" with flow control

Air quick disconnect ⅜" with flow control

Toggle valve for city/bottle water routing

Telescoping arm with swivel quad holder bar

Central vacuum accessory package

central vacuum canister with disposable solids collector

asepsis 3-way air/water syringe with tubing, 5'

autoclavable standard HVE lever valve with tubing, 5'

autoclavable saliva ejector lever valve with tubing, 5'

Cable for unit-mounted touchpad, 3900 chair

Internal umbilical, 5' with 3' external from chair to utility center

Contoured utility center

air pressure regulator with filter

master air shut-off valve

Chair adapter, fixed

CUSPIDOR, BASIC 1679 (option not shown)

Removable ceramic bowl

Manual cup fill

Timed bowl rinse

CUSPIDOR, DELUXE 1577 (option not shown)

One-piece ceramic bowl with touch lever water activator

Automatic cup fill

Timed bowl rinse

LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching temperature setting 4200K to 5500K

Shadow free illumination

CHAIR, HYDRAULIC 3900

Ultraleather® upholstery with wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad with 4 programmable preset positions

Seat height: minimum 17" maximum 29"



Popular Equipment Choices



CUSPIDOR 1577
(Sidebox comes with delivery)



CUSPIDOR 1679
(Sidebox comes with delivery)



ADD'L HVE POSITION
WITH ORAL DRY CUP
1400-113-CS



LED LIGHT
9078



MONITOR MOUNT
1300RM



DOCTOR'S ECONOMY
STOOL 6156C



ASST'S ECONOMY
STOOL 6157C



DOCTOR'S COMFORT
STOOL 6176

	FACTORY OPTIONS	PART NO	RETAIL	PRICE
Base PACKAGE	\$16,220	DELIVERY SYSTEM WITH VALUE	4195SI*	\$ 6240
4195SI Delivery System		Handpiece illumination system		
3900 Chair		One handpiece (5-hole), accommodates 1 add'l tubing.....	1700-021-CS	695
9080 Light		Add'l handpiece illumination tubing (5-hole)	0017-615-CS	265
1110-900 Light Post		One handpiece (6-pin), accommodates 1 add'l tubing.....	1700-019-CS	695
		Add'l handpiece illumination tubing (6-pin)	0017-603-CS	190
As shown PACKAGE	\$16,220	► Ancillaries use existing handpiece position (see page 6)		
Additionally equipped with:		Write in selection here.....		
3914-998-CS Chair Control, Touchpad, Unit-Mounted		Write in selection here.....		
DESIGNER FRIENDLY® SURFACE TREATMENTS				
STEP 1: Select unit plastic finishes:				
<input type="radio"/> Greywood <input type="radio"/> Laser Leather <input type="radio"/> Lemans				
<input type="radio"/> Forest Grey <input type="radio"/> White				
STEP 2: Select powder coat paint for unit & light:				
<input type="radio"/> Arctic Silver <input type="radio"/> Bronze <input type="radio"/> Charcoal				
<input type="radio"/> Copper <input type="radio"/> Forest Grey <input type="radio"/> Midnight				
<input type="radio"/> Pewter <input type="radio"/> Shadow <input type="radio"/> White				
STEP 3: Select standard powder coat paint for chair:				
<input type="radio"/> Forest Grey <input type="radio"/> Shadow				
OR upgrade chair base/back powder coat paint (add \$2,250):				
<input type="radio"/> Arctic Silver <input type="radio"/> Bronze <input type="radio"/> Charcoal				
<input type="radio"/> Copper <input type="radio"/> Midnight <input type="radio"/> Pewter <input type="radio"/> White				
STEP 4: Select Ultraleather® upholstery color choices for chair & deluxe /economy stools:				
<input type="radio"/> Apricot <input type="radio"/> Baltic <input type="radio"/> Cave <input type="radio"/> Charcoal				
<input type="radio"/> Diplomat Blue <input type="radio"/> Doe <input type="radio"/> Dove Grey				
<input type="radio"/> Eucalyptus <input type="radio"/> Lemon <input type="radio"/> Papyrus <input type="radio"/> Parrot				
<input type="radio"/> Raven Wing <input type="radio"/> Red <input type="radio"/> Riviera <input type="radio"/> Stone				
Pearlized colors: <input type="radio"/> Abalone <input type="radio"/> Midnight <input type="radio"/> Pewter				
OR select Special Color Ultraleather®				
<input type="radio"/> Chair 3907-999-SP add \$280				
<input type="radio"/> Deluxe/Economy Stool 6007-999-SP add \$100				
STEP 5: Upholstery Selections:				
Plush Selections:				
<input type="radio"/> Standard Back, Wide/Thin 3907-PL-WB N/C				
<input type="radio"/> Narrow Back 3907-PL-NB N/C				
<input type="radio"/> Narrow Back with Slings 3907-PL-NS add \$400				
Plush Additional Options:				
<input type="radio"/> Contrast/Deviated Stitching 3907-999-DS add \$280				
<input type="radio"/> Dual Color 3907-999-DC add \$280				
<input type="radio"/> Dual Color with Contrast/Deviated Stitching Call Customer Service for information (1.800.423.3555)				
<input type="radio"/> Massage Option 3900-100 add \$540				
Asepsis Selections:				
<input type="radio"/> Standard Back, Wide/Thin 3907-AS-WB N/C				
<input type="radio"/> Narrow Back 3907-AS-NB N/C				
<input type="radio"/> Narrow Back with Slings 3907-AS-NS add \$400				
Asepsis Additional Options:				
<input type="radio"/> Massage Option 3900-100 add \$540				
Memory Comfort Selections:				
<input type="radio"/> Wide/Thin 3907-MC-WB \$800				
<input type="radio"/> Dual Color 3907-MC-DC \$1080				
Equipment dimensions, see page 67				
Shipping weight estimate:				
Delivery unit & chair	430 lb / 200 kg			
Base package	547 lb / 249 kg			
As shown package	547 lb / 249 kg			

- * ITEM included in AS SHOWN PACKAGE price
- Choose only one post when sold with unit


► Often overlooked

TOTAL

Chairs & Stools

Materials & Benefits





New doesn't have to mean plastic. While our competitors continue to move toward more and more plastic parts, Forest uses metal whenever possible so you get real value from your investment.

3900 Patient Chair

Material: cast aluminum chair back, chair base & utility center base

Benefit: strong, impervious to fading or cracking

Material: pressure-formed plastic chair & utility center covers

Benefit: warp, impact & stain resistant

Material: aluminum hydraulic manifold & speed adjustment screws

Benefit: independent speed adjustment of base & back both up & down

Material: optional soft-touch back massage with integrated compressor

Benefit: saves central compressor capacity for delivery systems

Material: ultraleather® upholstery standard

Benefit: withstands heavy use & cleaning with glove-soft leather feel

Material: forest grey or shadow powder coat standard on 3900 chair / seven optional Designer Friendly® colors

Benefit: complements treatment room interior design color palette

6176/6187/6198 Comfort Stool

Material: tension woven material with rocking back

Benefit: eliminates pressure points, massages & supports the seat & back

Material: all-surface 22mm polyurethane casters

Benefit: smooth rolling eliminates squeaks / long-life

3900 Patient Chair Powder Coat Paint Color: Charcoal
Plush Upholstery: Ultraleather® Raven Wing w/ Red Contrast Stitching

3900 Chair Accessories

Knee Pillow

3907-317 \$ 380

Knee support ideal for patients with joint discomfort and/or long procedures

Zero gravity posture, eliminates pressure on knees

Specify upholstery color

Dual color upholstery option (upcharge)

3907-317-DC 75

Contrast/Deviated stitching option (upcharge)

3907-317-DS 75

Chair Control, Footswitch**Additional Chair Control, Footswitch**

Factory Installed

3914-998-CS N/C

Factory Installed 3914-999-CS \$ 500

After Market 3914-010-SH 500

4 programmable preset positions

3900 Remote Mount Touchpad with 25' Cable

4 programmable preset positions

Fixed Horizontal Mount

1700-079 \$ 625

White 1700-079-G 625

Forest Grey 1700-079-C 625

Designer Friendly® (specify color)



Swing Horizontal Mount

1700-089-HM \$ 625

White 1700-089-HM-G 625

Forest Grey 1700-089-HM-C 625

Designer Friendly® (specify color)



Swing Vertical Mount

1700-089-VM \$ 625

White 1700-089-VM-G 625

Forest Grey 1700-089-VM-C 625

Designer Friendly® (specify color)

Contoured Utility Center Housing & Cover

Included with PI and SI chair mounted packages

Measurements: 8 x 17 3/16 x 17 3/16"

3900-101 \$ 420

Factory option: 3900-101-G 420

Forest Grey: 3900-101-SH 420

Shadow 3900-101-C 520

Designer Friendly® (specify color)

Chairs + Stools

3900 Hydraulic Patient Chair




Plush



Asepsis

MODEL 3900 \$7690



STANDARD FEATURES

HYDRAULIC CHAIR 3900

Ultraleather® upholstery with wide/thin back for elbow support
Plastic toe guard, 25"
Double articulating pinch-style headrest
Thin back with full lumbar support
Right/left patient access
Swing-down locking armrests
Hydraulic cantilever mechanism
Emergency shut-off
30 degree base rotation with lock
Multifunction touchpad with 4 programmable preset positions, accommodates 3 additional controls
Seat height: minimum 17" maximum 29"

FACTORY OPTIONS

PART NO	RETAIL	PRICE
CHAIR (Forest Grey or Shadow).....	3900	\$ 7690
CHAIR (Designer Friendly®).....	3900C	9940
One chair control included (specify type below).....	3914-998-CS*	N/C
<input type="radio"/> Footswitch control		
<input type="radio"/> Touchpad, chair mounted		
Additional chair control (limit 2 on chair, limit 2 on unit)	3914-999-CS	500
<input type="radio"/> Footswitch control		
<input type="radio"/> Touchpad, chair mounted		
Chair adapter, fixed (for external umbilical)	1115-195	430
Chair adapter, fixed (for internal umbilical)	1115-530	565
Chair adapter, pivot (for internal umbilical)	1115-528	565
Contoured utility center housing & cover	3900-101	420
Contoured utility center housing & cover (Designer Friendly®)	3900-101-C	520
UPHOLSTERY SELECTIONS:		
Write in selection here (see page 27).....		
Write in selection here (see page 27).....		
Knee pillow (see page 25).....	3907-317	380
Soft-touch back massage option (not available with Memory Comfort)	3900-100	540
Plastic toe guard 33" in lieu of 25".....	3907-028-CS	20
LIGHT POSTS (for unit mount lights/post mount monitor assemblies see page 30) Fixed (chair adaptor sold separately)		
44" light post without slot, chair mount	1110-794	200
44" light post with slot for monitor cord.....	1110-050	200
Pivot (includes chair adaptor)		
Pivot light post, 3900 (light not included)	1112-644-3	1065
Pivot monitor post, 3900 (monitor mount assembly not included)	1331-644-3	1065
Monitor mount assembly, top post mount (see page 34)	1300TPM	990
Monitor mount assembly, ring mount (see page 34)	1300RM	990

TOTAL

Often overlooked

DESIGNER FRIENDLY® SURFACE TREATMENTS

Powder coat paint for chair:

Standard colors: Forest Grey Shadow

OR upgrade chair base/back powder coat paint (add \$2,250):

<input type="radio"/> Arctic Silver	<input type="radio"/> Bronze	<input type="radio"/> Charcoal
<input type="radio"/> Copper	<input type="radio"/> Midnight	<input type="radio"/> Pewter
		<input type="radio"/> White

ULTRALEATHER® UPHOLSTERY

Ultraleather® upholstery color choices for chair:

<input type="radio"/> Apricot	<input type="radio"/> Baltic	<input type="radio"/> Cave	<input type="radio"/> Charcoal
<input type="radio"/> Diplomat Blue	<input type="radio"/> Doe	<input type="radio"/> Dove Grey	<input type="radio"/> Eucalyptus
<input type="radio"/> Lemon	<input type="radio"/> Papyrus	<input type="radio"/> Parrot	
<input type="radio"/> Raven Wing	<input type="radio"/> Red	<input type="radio"/> Riviera	<input type="radio"/> Stone
Pearlized colors:	<input type="radio"/> Abalone	<input type="radio"/> Midnight	<input type="radio"/> Pewter
<input type="radio"/> Special color (specify _____)	3907-999-SP	\$280	

Upholstery selections

(see page 27 for pricing)

Shipping weight estimate:

335 lb / 152 kg



3900 Chair Upholstery Options



Memory Comfort



Narrow Back or Narrow Back with Slings

STANDARD FEATURES	FACTORY OPTIONS	PART NO	RETAIL	PRICE
MEMORY COMFORT UPHOLSTERY Designed for all day/ultra-long procedures or greater general patient comfort 3 layers of foam topped with 1.5 inches of memory foam for extra comfort Side-bolsters on seat cushion support legs Memory foam headrest	UPHOLSTERY SELECTIONS: Plush selections: Standard back, wide/thin 3907-PL-WB N/C Narrow back 3907-PL-NB N/C Narrow back with slings 3907-PL-NS \$400			
NARROW BACK Perfect for Pediatric or Orthodontic dentistry Unrestricted patient access	Plush additional options: Contrast/deviated stitching 3907-999-DS \$280 Dual color 3907-999-DC \$280 Dual color with contrast/deviated stitching Call customer service for information (1.800.423.3555) Massage option 3900-100 \$540			
NARROW BACK WITH SLINGS <i>Same as above, with additional features:</i> Slings provide alternative for additional patient elbow support Slings create a barrier between patient and practitioner	Asepsis selections: Standard back, wide/thin 3907-AS-WB N/C Narrow Back 3907-AS-NB N/C Narrow Back with Slings 3907-AS-NS \$400			
	Asepsis additional options: Massage option 3900-100 \$540			
	Memory comfort selections: Wide/thin 3907-MC-WB \$800 Dual color 3907-MC-DC \$1080			
			TOTAL	

Shipping weight estimate:
335 lb / 152 kg

Chairs + Stools

Economy & Deluxe Stools



Model 6156C \$485



Model 6157C \$625



Model 6159C \$825



Model 6160C \$925



Model 6161C \$1100

STANDARD FEATURES

Ultraleather® upholstery
Injection molded foam
5-leg fully polished base
Twin wheel hard floor casters
Accommodates maximum of 300 lbs.

DOCTOR'S ECONOMY 6156C

Backrest, non-adjusting
Round seat with single lever height adjustment
Pneumatic cylinder height adjustment, 18" - 23"

ASSISTANT'S ECONOMY 6157C

Round seat with single lever height adjustment
Torso support, non-ratcheting, with height adjustment
Pneumatic cylinder height adjustment, 21½" - 28"
Foot ring affixed to base

DOCTOR'S DELUXE 6159C

Backrest with height adjustment
Contoured seat with triple lever adjustment
Pneumatic cylinder height adjustment, 18½" - 23"

ASSISTANT'S DELUXE 6160C

Contoured seat with single lever adjustment
Torso support with ratcheting and height adjustment
Pneumatic cylinder height adjustment, 21 ½" - 28"
Foot ring with adjustment

ASSISTANT'S DELUXE WITH BACK 6161C

Backrest with height adjustment
Contoured seat with triple lever adjustment
Torso support with ratcheting and height adjustment
Pneumatic cylinder height adjustment, 21 ½" - 28"
Foot ring with adjustment

FACTORY OPTIONS

STOOL	MODEL CHOICE	PRICE:
Vinyl in lieu of Ultraleather®	6007-999-CS*	\$ - 100
Special color upholstery (upcharge)	6007-999-SP	100
CYLINDER CHOICES:		
Doctor's short cylinder, 16" - 21" (6156C/6159C)	6000-628-CS	N/C
Doctor's tall cylinder, 21½" - 28" (6156C/6159C) foot ring with support highly recommended	6000-604-CS	N/C
Assistant's short cylinder, 18" - 23" (6157C/6160C/6161C)	6000-605-CS	N/C
Foot ring with support (6156C/6159C)	6000-602	115
TOTAL		

Often overlooked

Shipping weight estimate:
6156C/6159C 35 lb / 15 kg
6161C 45 lb / 21 kg
6157C/6160C 40 lb / 18 kg

* Specify vinyl manufacturer and color (not all colors may be available. Call Forest for ordering details 1-800-423-3555.)

ULTRALEATHER® UPHOLSTERY

- | | | | | |
|-------------------------------------|-------------------------------|----------------------------------|----------------------------------|--------------------------------|
| Ultraleather®: | <input type="radio"/> Apricot | <input type="radio"/> Baltic | <input type="radio"/> Cave | <input type="radio"/> Charcoal |
| <input type="radio"/> Diplomat Blue | <input type="radio"/> Doe | <input type="radio"/> Dove Grey | <input type="radio"/> Eucalyptus | <input type="radio"/> Lemon |
| <input type="radio"/> Papyrus | <input type="radio"/> Parrot | <input type="radio"/> Raven Wing | <input type="radio"/> Red | <input type="radio"/> Riviera |
| <input type="radio"/> Stone | Pearlized colors: | <input type="radio"/> Abalone | <input type="radio"/> Midnight | <input type="radio"/> Pewter |

OR select Special Color Ultraleather®

Deluxe/Economy Stool

6007-999-SP add \$100



Model 6176 \$825



Model 6187 \$1100



Model 6198 \$975

STANDARD FEATURES

Air circulating weave construction
Active tilt® seat (actively moves with user)
Tilt seat lever function
Contoured seat
5-leg base
Twin wheel hard floor casters
Accommodates maximum of 300 lbs.

FACTORY OPTIONS

STOOL	MODEL CHOICE	PRICE :
► CYLINDER CHOICES:		
Doctor's/Office short cylinder, 16" – 21" (6176/6198)	6000-628-CS	N/C
Doctor's/Office tall cylinder, 21½" – 28" (6176/6198) foot ring with support highly recommended.....	6000-604-CS	N/C
Assistant's short cylinder, 18" – 23" (6187).....	6000-605-CS	N/C
Foot ring with support (6176/6198)	6000-602	115

TOTAL _____

DOCTOR'S AND OFFICE COMFORT 6176 / 6198

Triple lever adjustment - independent back and seat
Backrest with angle & height adjustment
Fully adjustable office arms (6198)
Pneumatic cylinder height adjustment, 19 ¼" – 24 ½"

► Often overlooked

ASSISTANT'S COMFORT 6187

Triple lever adjustment - independent back and seat
Backrest with angle & height adjustment
Fully adjustable ratcheting assistant arm
Foot ring with height adjustment
Pneumatic cylinder height adjustment, 23" – 30 ½"

Shipping weight estimate:

6176 48 lb / 23 kg
6187 60 lb / 28 kg
6198 54 lb / 24 kg

Lights & Monitor Mounts


Materials & Benefits



Powder Coat Paint Colors




Forest Grey



Shadow




White



Midnight




Arctic Silver




Pewter



Copper



Charcoal



Bronze

9080/9090 Halogen

Material: reflector material embedded into glass / lexan light shield

Benefit: reflector will not peel or crack / shield is scratch resistant

Material: oil impregnated bronze bushing / 3-axis light head joint

Benefit: no light drift / 3-axis light keeps light pattern parallel with the long axis of mouth while freeing up space above patient's head

9078/9079 LED

Material: ten 1-watt LED lamps

Benefit: creates a light pattern with feathered edges reducing eye fatigue / consumes approximately 75 percent less power than halogen lights / no heat / 22 year projected LED life span

Material: cast aluminum light head (9080/9090/9078/9079)

Benefit: allows for nine Designer Friendly® powder coat paint colors

9078 LED Light Powder Coat Paint Color: Charcoal



MODEL 9080 \$2090



MODEL 9078 \$3490



Optional: LED light core white cosmetic whitening accelerator ancillary 9072-061



Optional: LED light anti-polymerization filter and mirror 9072-060

STANDARD FEATURES**HALOGEN LIGHT 9080**

Automatic on/off activation with override
3-Axis light head movement
Additional bulb
Color temperature, 4,200 to 5,500 Kelvin
Designer Friendly® powder coat paint light head

LED LIGHT 9078

Automatic on/off activation with override
3-Axis light head movement
Easy release autoclavable handles
Light source 10 LED lamps
Color temperature, 4,200 to 6,000 Kelvin
Designer Friendly® powder coat paint light head

FACTORY OPTIONS

LIGHT	MODEL CHOICE	PRICE:	_____
Replacement bulb (9080)	0014-300	\$ 35	_____
Anti-polymerization filter and mirror (9078)	9072-060	310	_____
Core white cosmetic whitening accelerator ancillary (9078)	9072-061	1160	_____
Short rigid arm 20" in lieu of standard 32"	1110-227-CS	N/C	_____
LIGHT POSTS:			
44" with slot for sidebox access	1110-900	200	_____
44" without slot, chair mount	1110-794	200	_____
35" with post clamp	1110-820	260	_____
Pivot (mounts on Forest 3900 chair only)	1112-644	500	_____
Pivot (includes 3900 cast chair adapter)	1112-644-3	1065	_____
LIGHT POSTS, WITH MONITOR CORD SLOT:			
44" with dual slot for sidebox access	1110-170	200	_____
44" with single slot, chair mount	1110-050	200	_____
35" with single slot & post clamp	1110-045	260	_____

TOTAL




Shipping weight estimate:
35 lb / 16 kg

DESIGNER FRIENDLY® SURFACE TREATMENTS

- Powder coat paint finishes: Arctic Silver Bronze Charcoal Copper
 Forest Grey Midnight Pewter Shadow White

Lights + Monitor Mounts

Track Mount Lights



STANDARD FEATURES

LIGHT + TRACK ASSEMBLY

Pre-wired through flex arm
Track/Trolley assembly
Vertical light tube with 6 inch horizontal side offset for clearance of optional monitor
Flex arm & light head assembly
Trolley mounting plate and shroud
Drop post (specify mounting height)

MONITOR MOUNT ASSEMBLY

Pre-wired from track through flex arm
Coax cable
Low voltage wiring (2 conductor 18 awg wire)
HDMI shielded cable
Lower positioning handle

Monitor Requirements:
12 volt with power supply
weight: up to 18 lbs max, width of 19"

Shipping weight estimate:

Light/track 115 lb / 53 kg

1305T 140 lb / 64 kg (includes mount/track light weight)



FACTORY OPTIONS

LIGHT / TRACK.....	MODEL CHOICE_____	PRICE :	_____
Monitor mount assembly	1305T	\$ 2320	_____
Replacement bulb (9080)	0014-300	35	_____
Anti-polymerization filter and mirror (9078)	9072-060	310	_____
Core white cosmetic cosmetic whitening accelerator ancillary (9078) ...	9072-061	1160	_____
S Video, cable pre-installed.....	0014-282	40	_____
RCA cables (3 conductor), pre-installed	0014-283	40	_____

PART NO RETAIL

PRICE

TOTAL

Often overlooked

CHOOSE YOUR STANDARD MOUNTING OPTION

Mounting height:	Light post length:	Monitor post length:
<input type="radio"/> 8'	8.25"	20.25"
<input type="radio"/> 8' 6"	14.25"	26.25"
<input type="radio"/> 9'	20.25"	32.25"
<input type="radio"/> 9' 6"	26.25"	38.25"
<input type="radio"/> 10'	32.25"	44.25"
<input type="radio"/> 10' 6"	38.25"	50.25"
<input type="radio"/> 11'	44.25"	56.25"
<input type="radio"/> _____ Custom mounting height	<input type="radio"/> _____ Custom post length	<input type="radio"/> _____ Custom post length

DESIGNER FRIENDLY® SURFACE TREATMENTS

Powder coat paint finishes: Arctic Silver Bronze Charcoal Copper
 Forest Grey Midnight Pewter Shadow White



MODEL 9080-C \$2690
HALOGEN LIGHT



MODEL 9078-C \$4090
LED LIGHT

STANDARD FEATURES

LIGHTS

Ceiling mounting plate and shroud with drop post
Specify mounting height

HALOGEN LIGHT 9080-C

Automatic on/off activation with override
3-Axis light head movement
Additional bulb
Color temperature, 4,200 to 5,500 Kelvin
Designer Friendly® powder coat paint light head

LED LIGHT 9078-C

Automatic on/off activation with override
3-Axis light head movement
Easy release autoclavable handles
Light source 10 LED lamps
Color temperature, 4,200 to 6,000 Kelvin
Designer Friendly® powder coat paint light head

FACTORY OPTIONS

LIGHT	MODEL CHOICE	PRICE :	
Replacement bulb (9080)	0014-300	\$ 35	
Anti-polymerization filter and mirror (9078)	9072-060	310	
Core white cosmetic whitening accelerator ancillary (9078)	9072-061	1160	
			TOTAL

Often overlooked

CHOOSE YOUR STANDARD MOUNTING OPTION

Mounting height:

- 8' 12"
- 8' 6" 18"
- 9' 24"
- 9' 6" 30"
- 10' 36"
- 10' 6" 42"
- 11' 48"
- _____ Custom mounting height

Post length:

- _____ Custom post length

DESIGNER FRIENDLY® SURFACE TREATMENTS

Shipping weight estimate:
50 lb / 23 kg



- | | | | | |
|-----------------------------|-------------------------------------|--------------------------------|--------------------------------|------------------------------|
| Powder coat paint finishes: | <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze | <input type="radio"/> Charcoal | <input type="radio"/> Copper |
| | <input type="radio"/> Forest Grey | <input type="radio"/> Midnight | <input type="radio"/> Pewter | <input type="radio"/> Shadow |
| | | | | <input type="radio"/> White |

Lights + Monitor Mounts

Cabinet + Wall Mount Lights



**MODEL 9079C \$4590
(STANDARD) 9079C-T (TALL)
LED CABINET**

**MODEL 9090C \$3190
(STANDARD) 9090C-T (TALL)
HALOGEN CABINET**

**MODEL 9090W \$3190
(STANDARD) 9090W-T (TALL)
HALOGEN WALL**

**MODEL 9079W \$4590
(STANDARD) 9079W-T (TALL)
LED WALL**

STANDARD FEATURES

LIGHTS

Cabinet mount, internal vertical mount with post
Wall mount, vertical mounting plate
Standard rigid arm
Tall rigid arm (additional 1' drop)

HALOGEN LIGHT 9090

Automatic on/off activation with override
3-Axis light head movement
Additional bulb
Color temperature, 4,200 to 5,500 Kelvin
Designer Friendly® powder coat paint light head

LED LIGHT 9079

Automatic on/off activation with override
3-Axis light head movement
Easy release autoclavable handles
Light source 10 LED lamps
Color temperature, 4,200 to 6,000 Kelvin
Designer Friendly® powder coat paint light head

FACTORY OPTIONS

LIGHT	MODEL CHOICE	PRICE :
Replacement bulb (9090)	0014-300	\$ 35
Anti-polymerization filter and mirror (9079)	9072-060	310
Core white cosmetic cosmetic whitening accelerator ancillary (9079) ...	9072-061	1160
Cabinet mounting options: (for dimensions see page 70)		
Horizontal top cabinet mount in lieu of internal mount.....	1110-881-CS	N/C
Vertical exterior cabinet mount in lieu of internal mount	1110-886-CS	N/C

TOTAL

► ** SPECIFY CABINET/WALL MOUNTING HEIGHT **

► Often overlooked


Shipping weight estimate:

Cabinet 60 lb / 28 kg
Wall 65 lb / 30 kg



DESIGNER FRIENDLY® SURFACE TREATMENTS

Powder coat paint finishes: Arctic Silver Bronze Charcoal Copper
 Forest Grey Midnight Pewter Shadow White

**MODEL 1300RM RING POST MOUNT \$990***Monitor screens and chair-mounted posts not included***MODEL 1300TPM TOP POST MOUNT \$990****MODEL 1300-025-CS \$410****ORGANIZER BIN MOUNT***Used with rear delivery systems and carts with storage organizer bin***STANDARD FEATURES****FACTORY OPTIONS****PART NO****RETAIL****PRICE****MONITOR MOUNT ASSEMBLY,
RING POST MOUNT 1300RM**

ICW Quicklink monitor mount
Rex positioning handles
22" ring mount rigid arm

**MONITOR MOUNT ASSEMBLY,
TOP POST MOUNT 1300TPM**

ICW Quicklink monitor mount
Rex positioning handles
22" top post mount rigid arm

ORGANIZER BIN MOUNT 1300-025-CS

Forest's unique locking tilt mechanism with knob
 $\frac{1}{2}$ " pin

MONITOR MOUNT ASSEMBLY MODEL CHOICE**PRICE:****► LIGHT POST (1300RM):**

44" with dual slot, fixed sidebox mount, with Forest light bushing 1110-170
35" with slot & post clamp, fixed mount, with Forest light bushing 1110-045
Pivot light posts & requirements (call to order)

► MONITOR POST (1300TPM):

30" with dual slot, fixed sidebox mount, with oilite bushing 1331-998
30" with slot & post clamp, fixed mount, with oilite bushing 1331-997
Pivot monitor posts & requirements (call to order)

PRICE:

TOTAL

► Often overlooked**Monitor Requirements:**

12 volt with power supply
weight: 8 $\frac{1}{4}$ to 13 lbs max, width of 19"

Shipping weight estimate:

1300RM, 1300TPM 15 lb / 7 kg

1300-025-CS 5 lb / 2 kg

**DESIGNER FRIENDLY® SURFACE TREATMENTS**

Powder coat paint finishes:

 Arctic Silver Bronze Charcoal Copper Forest Grey Midnight Pewter Shadow White



Rear Mount Side Delivery

Materials & Benefits



FUN FACTOR

More and more people are realizing FUN is the most important WHY in life. Forest embraces fun. Whether it's red stitching on black upholstery or choosing the colors of your favorite sports team, all you've got to do is set your imagination free.

7000 Rear Mount Delivery Series

Material: north america's only all metal control head chassis

Benefit: resistant to UV & cleaners

Material: hinged solid surface lid

Benefit: long lasting surface gives easy access to ancillary products

Material: optional exclusive brass handpiece lockout valves

Benefit: allows assistant to remove a handpiece from its holder to change a burr while dentist continues working

Material: four steel mounting configurations

Benefit: wall, under cabinet, beneath counter & cart

Material: secondary aluminum work surface

Benefit: dedicated work surface for burr block, electric handpiece attachment, mirror, explorer

Material: steel height adjust mechanism


Benefit: adjustable work surface & handpiece delivery arm height

Material: nine designer friendly® powder coat paint colors including white & forest grey

Benefit: dental equipment with a WOW factor


7010 Beneath Counter Control Powder Coat Paint Color: Pewter | Solid Surface: Fossil
3900 Patient Chair Powder Coat Paint Color: Shadow
Plush Duo Color Upholstery: Ultraleather® Waxed Canvas and Tusk
9078 LED Light Powder Coat Paint Color: Pewter
Ring Mount Monitor Arm Powder Coat Paint Color: Pewter

M&CC PFR8144



Made for each other

From the ground up, Forest delivery systems and M&CC cabinets are designed for easy integration.



Shown with optional storage organizer tray and assistants holder bar.

see pages 36-37 for ordering information

Good thinking


Low profile $\frac{1}{2}$ " hardened steel mount leaves room for your stool casters.

Start with a strong foundation


- Your choice of $\frac{3}{4}$ " hardened steel $11'' \times 14\frac{1}{2}''$ (pre-installed on the PFR8144) or optional $10'' \times 24''$ mounting plate that bolts to the bottom of your counter top.
- The $6'' \times 6''$ leveling plate is next, then a set of needle roller bearings with friction adjustment.
- A super strong $2\frac{1}{2}$ " steel tube scissor arm assembly.
- Needle roller bearings are used where the two arm sections pivot.

Rear Mount Side Delivery Systems

7000 Series



*Shown with optional storage organizer and swivel tray
1110-179-CS and 1110-909-CS*



*Shown with optional
assistants holder 0024-156*

*Shown with optional
handpiece lockouts 0010-748*

STANDARD FEATURES

DUO SWING DELIVERY SYSTEM 7000 SERIES

Duo control

- automatic control for three handpieces
- asepsis 3-way air water syringe with tubing, 4'
- handpiece tubing, 4-hole Midwest, asepsis style, 4'
- pneumatic master control block with diaphragm
- master on/off
- individual handpiece pressure adjustment
- individual water coolant adjustment
- non-retraction water coolant
- air coolant adjustment
- handpiece flush
- handpiece exhaust filter/oil collector
- self-contained clean water system with quick disconnect,
single bottle, 1½ liter
- wet/dry foot control

Height adjustable work surface, 5"

"Burr Buddy" tray holder and covered Zirc tray

Central vacuum accessory package

- central vacuum canister with disposable solids collector
- asepsis 3-way air/water syringe with tubing, 5'
- autoclavable standard HVE lever valve with tubing, 5'
- autoclavable saliva ejector lever valve with tubing, 5'

Solid work surface, 19" x 24"

Universal heavy duty mounting plate

- 11" x 14 ½", beneath counter (7010BC)

Laminate wall mount panel 18" x 24" (7010WM)

Standard utility center housing and cover (7010MC)

Air pressure regulator with filter


Master air shut-off valve

Umbilical, 7'

Mounting plate dimensions, see pages 71-72

Shipping weight estimate:

7010BC, UC, WM 130 lb / 59 kg



MODEL 7010BC \$6190

BENEATH COUNTER MOUNT

*Universal heavy-duty mounting plate 11" x 14 1/2"
included (mounts to under counter)*

Specify cabinet manufacturer

MODEL 7010UC \$6190

UNDER CABINET FLOOR MOUNT

MODEL 7010WM \$6190

WALL MOUNT

► DESIGNER FRIENDLY® SURFACE TREATMENTS

- Powder coat paint finishes: Arctic Silver Bronze
 Charcoal Copper Forest Grey
 Midnight Pewter Shadow
 White

Solid work surface colors:

Corian®: Antarctica

Gibraltar®: Arctic Melange Baja Melange

Meganite®: Cinder Granite Goldrush Gemstone
 Moonstone Platinum Granite
 Rocky Road Granite Sanibel Granite
 Silver Mist Starry Starry Night

Laminate work surface colors in lieu of solid surface 0023-776-LAM
deduct —\$325

Formica®: Folkstone Oyster Gray Surf
 White Blanc

Laminate wall mount panel (Model WM)

Formica®: Folkstone Oyster Gray Surf
 White Blanc

FACTORY OPTIONS

DELIVERY SYSTEM	MODEL CHOICE	PRICE :
Handpiece Illumination System		
One handpiece (5-hole), accommodates 1 add'l tubing.....	1700-021-CS	\$ 695
Add'l handpiece illumination tubing (5-hole)	0017-615-CS	265
One handpiece (6-pin), accommodates 1 add'l tubing.....	1700-019-CS	695
Add'l handpiece illumination tubing (6-pin)	0017-603-CS	190
► Ancillaries use existing handpiece position (see page 6)		
Write in selection here		
Write in selection here		
Additional handpiece position, wet	1100-189-CS	405
Additional handpiece position, air drive only.....	1100-190-CS	195
Lockout feature for 3 handpieces	0010-748	150
Additional standard autoclavable HVE		
with tubing, 5' & autoclavable oral dry cup	1400-113-CS	85
Additional assistant's quad holder bar	0024-156	150
Quick disconnect with flow control 1/4" (for water)	0016-454-CS	100
Quick disconnect with flow control 3/8" (for air)	0016-456-CS	100
Storage organizer with bins & frame.....	1110-179-CS	420
► Swivel tray	1110-909-CS	335
► Touchpad for 3900 chair (25 ft. Cable).....	1700-081-1-CS	500
Chrome accents	9000-998-CS	135
Monitor mount assembly (monitor screen not included).....	1300-025-CS	395
Standard utility center housing & cover (BC/UC/WM).....	1115-020	130
In floor recessed utility center housing & cover.....	1115-305	270
Additional umbilical (per foot).....	1115-999-CS	30
Routing valve kit with water shutoff & regulator/factory installation..	1115-315-CS	235
Universal mounting plate, 10" x 24" in lieu of 11" x 14 1/2"	0010-084-CS	125
Standard color laminate work surface in lieu of solid surface	0023-776-LAM	-325
Special color solid or laminate work surface.....	0023-776-SSV / 0023-776-X	(call for price)

TOTAL

► Often overlooked

Cabinet mounting configuration see page 35

Rear Mount Side Delivery Systems

7000 Series Mobile Cart



Shown with optional storage organizer and swivel tray
1110-179-CS and 1110-909-CS



MODEL 7010MC \$5590

Shown with optional handpiece lockouts 0010-748
and assistants holder 0024-156

STANDARD FEATURES

DUO SWING CART

Automatic control for three handpieces

asepsis 3-way air water syringe with tubing, 4' handpiece tubing, 4-hole Midwest, asepsis style, 4' pneumatic master control block with diaphragm master on/off

individual handpiece pressure adjustment

individual water coolant adjustment

non-retraction water coolant

air coolant adjustment

handpiece flush

handpiece exhaust filter/oil collector

self-contained clean water system with quick disconnect, single bottle, 1½ liter

wet/dry foot control

"Burr Buddy" sliding tray holder and covered Zirc tray

Pneumatic work surface easy lift/close mechanism

Central vacuum accessory package

canister with disposable solids collector

asepsis 3-way air/water syringe with tubing, 5'

autoclavable standard HVE lever valve with tubing, 5'

autoclavable saliva ejector lever valve with tubing, 5'

Solid work surface, 19" x 24"

Adjustable height frame, 24" x 46"

Twin wheel hard floor casters

External umbilical, 6'

Standard utility center housing and cover

Air pressure regulator with filter

Master air shut-off valve

FACTORY OPTIONS

CART..... 7000MC \$ 5950

Handpiece Illumination System

One handpiece (5-hole), accommodates 1 add'l tubing..... 1700-021-CS

695

Add'l handpiece illumination tubing (5-hole) 0017-615-CS

265

One handpiece (6-pin), accommodates 1 add'l tubing..... 1700-019-CS

695

Add'l handpiece illumination tubing (6-pin) 0017-603-CS

190

► Ancillaries use existing handpiece position (see page 6)

Write in selection here.....

Write in selection here.....

Additional handpiece position, wet 1100-189-CS

405

Additional handpiece position, air drive only 1100-190-CS

195

Lockout feature for 3 handpieces 0010-748

150

Additional standard autoclavable HVE

with tubing, 5' & autoclavable oral dry cup..... 1400-113-CS

85

Additional assistant's quad holder bar 0024-156

150

Quick disconnect with flow control ¼" (for water) 0016-454-CS

100

Quick disconnect with flow control ¾" (for air) 0016-456-CS

100

► Storage organizer with bins & frame 1110-179-CS

420

► Swivel tray 1110-909-CS

335

Touchpad for 3900 chair (25 ft. Cable) 1700-081-1-CS

500

Chrome accents 9000-998-CS

135

Less standard utility center housing & cover 1115-020-CS

-50

Additional umbilical (per foot) 1115-999-CS

30

Monitor assembly (monitor screen not included) 1300-025-CS

395

Routing valve kit with water shutoff & regulator/factory installation... 1115-315-CS

235

Standard color laminate work surface in lieu of solid surface 0023-776-LAM

-325

Special color solid or laminate work surface 0023-776-SSV / 0023-776-X

(call for price)

TOTAL

► Often overlooked

Solid work surface colors:

Corian®: Antarctica

Baja Melange

Gibraltar®: Arctic Melange

Goldrush Gemstone

Meganite®: Cinder Granite

Moonstone

Platinum Granite

Rocky Road Granite

Sanibel Granite

Silver Mist

Starry Starry Night

Laminate work surface colors in lieu of solid surface 0023-776-LAM

deduct —\$325

Formica®: Folkstone

Oyster Gray

Surf

White Blanc

Shipping weight estimate: 119 lb / 54 kg

DESIGNER FRIENDLY® SURFACE TREATMENTS

Powder coat paint finishes:

Arctic Silver

Bronze

Charcoal

Copper

Shadow

White



Forest Simplified

Ordering Made Easy



SIMPLIFIED

FOREST

Too many options can sometimes be overwhelming. That's why this product offering has just one color choice. You can make upgrading a treatment room or matching existing equipment easy.



Forest Grey

Delivery System, Chair Mounted page 40

Fixed Mount with Sidebox, External Umbilical
4485SE Standard IC / 4987SE Euro / 4195SE Value

Fixed Mount without Sidebox, External Umbilical
4485FE Standard IC / 4987FE Euro / 4195FE Value

Cuspidors / Assistant's Instrumentation page 44

Post Mount, External Umbilical
1636 / 1637 Basic Cuspidor
1075 Vacuum Telescoping Arm

Rear Mount Side Delivery Systems page 45

Cabinetry Mounting Configurations
M&CC PFR8144 / M&CC DR Cabinet

Work Surface Rear Delivery
Work Surface Only (3386), Assistant's Instrumentation Only (3486)
Beneath Counter (BC), Under Cabinet (UC), Wall Mount (WM)

Work Surface Cart
Work Surface Only (3386MC), Assistant's Instrumentation Only (5386), Duo (5886)

Standard IC Cart
Doctor's (5350), Duo (5355)

Standard IC Side Delivery
Cabinet (3235), Wall Mount (3230), Beneath Counter (2400)

Rear Beneath Counter Mount
Duo with Vertical Mount Panel (2184), Assistant's (1284)

Panel Mount
3 handpiece (2500), 2 handpiece (2520)

MT Compact Control
3 handpiece (MT600), 2 handpiece (MT610)


Factory Options + Parts page 56

General Information page 66

Dimensions
Ordering Guide
Terms + Warranty
Professional Sales Associates Contact
Index

Delivery Systems, Chair Mounted

Fixed Mount with Sidebox, External Umbilical



MODEL
PRICE

4485SE STANDARD IC
\$7440

4987SE EURO
\$8440

4195SE VALUE
\$6240

PACKAGE OFFERINGS SEE PAGES 18-23

STANDARD FEATURES


PNEUMATIC
3 HP
SYRINGE

PNEUMATIC
4 HP
SYRINGE

PINCH
3 HP
SYRINGE

delivery systems can be mounted on any 2" fixed post mount brand of chair

CHOOSE YOUR
CONTROL HEAD



OPTIONAL EQUIPMENT: CHOOSE ONE:

- 1577 Add \$1450
Cuspidor, Deluxe
(Sidebox comes with delivery)



- 1679 Add \$995
Cuspidor, Basic
(Sidebox comes with delivery)



- 9080 Add \$2090
Light, Halogen
1110-900 Add \$ 200
Light post



- 9078 Add \$3490
Light, LED
1110-900 Add \$ 200
Light post



STANDARD FEATURES	FACTORY OPTIONS	PART NO	RETAIL	PRICE
DELIVERY SYSTEM	DELIVERY SYSTEM	MODEL CHOICE	PRICE:	
Choice of Control Head	Handpiece illumination system			
Automatic control for handpieces	One handpiece (5-hole), accommodates 1 add'l tubing.....	1700-021-CS	\$ 695	
asepsis 3-way air/water syringe with tubing, 5' (Euro tubing 3')	Add'l handpiece illumination tubing (5-hole)	0017-615-CS	265	
handpiece tubing, 4-hole Midwest, asepsis style, 5' (Euro tubing 3')	One handpiece (6-pin), accommodates 1 add'l tubing	1700-019-CS	695	
master control block (with pneumatic or pinch valves)	Add'l handpiece illumination tubing (6-pin)	0017-603-CS	190	
master on/off	● Ancillaries use existing handpiece position (see pages 4-5)			
individual handpiece pressure adjustment	Write in selection here.....			
individual water coolant adjustment	Write in selection here.....			
non-retraction water coolant	Additional handpiece position, wet (4195/4485)	1100-189-CS	405	
air coolant adjustment	Additional handpiece position, air drive only (4195/4485).....	1100-190-CS	195	
handpiece flush	Swivel tray, arm-mounted (euro standard)	1110-468	335	
handpiece exhaust filter/oil collector	Less telescoping arm and assistant's instrumentation.....	1400-000-CS	-350	
single handle brake release function	Chrome accents	9000-998-CS	135	
dual handle brake release function (Euro)	Additional autoclavable standard HVE			
swivel instrument tray (Euro)	with tubing, 5' & autoclavable oral dry cup.....	1400-113-CS	85	
wet/dry foot control	Additional umbilical (per foot).....	1115-999-CS	30	
Flexible control arm with pneumatic brake	Self-contained dual bottle clean water system			
vertical travel 28" / horizontal travel 62"	with quick disconnect in lieu of single	1700-365-CS	150	
Sidebox with self-contained clean water system with quick disconnect, single bottle, 1½ liter	Without chair adapter	1115-099-CS	-150	
electrical outlet/UL	CUSPIDOR, (sold only with delivery)			
Water quick disconnect 1/4" with flow control	● DELUXE ○ 1577		1450	
Air quick disconnect 3/8" with flow control	or BASIC ○ 1679		995	
Toggle valve for city/bottle water routing	Deluxe cuspidor unit mount activation (1577)	1100-480-CS	200	
Telescoping arm with swivel quad holder bar	LIGHT, HALOGEN (mounting options, see pages 30-33)	○ 9080	2090	
Central vacuum accessory package	or LED ○ 9078		3490	
central vacuum canister with disposable solids collector	● LIGHT POST (specify light make/model)			
asepsis 3-way air/water syringe with tubing, 5'	Light post, 44" (Forest).....	1110-900	200	
autoclavable standard HVE lever valve with tubing, 5'	Light post, 44" (Pelton & Crane).....	1110-890	200	
autoclavable saliva ejector lever valve with tubing, 5'	● MONITOR MOUNT ASSEMBLY (monitor screen not included) (see page 34)			
External umbilical, 5'	Ring mount (used with light post).....	1300RM	990	
Deluxe utility center	Top post mount (used with monitor post)	1300TPM	990	
air pressure regulator with filter	● Monitor post (for top post mount).....	1331-998	200	
master air shut-off valve				
● Chair adapter, fixed (specify chair make/model)				
DELUXE CUSPIDOR 1577		TOTAL		
One-piece ceramic bowl with touch lever water activator				
Automatic cup fill				
Timed bowl rinse				
BASIC CUSPIDOR 1679				
Removable ceramic bowl				
Manual cup fill				
Timed bowl rinse				

Shipping weight estimate:

4485SE Delivery Unit 105 lb / 48 kg

4987SE Delivery Unit 110 lb / 50 kg

4195SE Delivery Unit 100 lb / 46 kg

DESIGNER FRIENDLY® SURFACE TREATMENTS

Unit plastic finishes: Forest Grey White

Powder coat paint finishes: Forest Grey White

Delivery Systems, Chair Mounted

Fixed Mount without Sidebox, External Umbilical



MODEL	4485FE STANDARD IC	4987FE EURO	4195FE VALUE
PRICE	\$5090	\$6170	\$3890
PACKAGE PRICE WITH 3900 CHAIR, 9080 LIGHT AND 1110-820 LIGHT POST			
	\$15,130	\$16,210	\$13,930
STANDARD FEATURES	PNEUMATIC 3 HP SYRINGE	PNEUMATIC 4 HP SYRINGE	PINCH 3 HP SYRINGE

delivery systems can be mounted on any 2" fixed post mount brand of chair



3900 Add \$7690
Chair, Hydraulic



OPTIONAL EQUIPMENT:

- 1075FE Add \$880
Assistant's telescoping arm with instrumentation



- 1636FE Add \$1250
Cuspidor, Basic



- 1637FE Add \$1720
Cuspidor, Basic with assistant's instrumentation



- 9080 Add \$2090
Light, Halogen
1110-820 Add \$ 260
Light post



- 9078 Add \$3490
Light, LED
1110-820 Add \$260
Light post



STANDARD FEATURES	FACTORY OPTIONS	PART NO	RETAIL	PRICE
DELIVERY SYSTEM	DELIVERY SYSTEM MODEL CHOICE * Handpiece illumination system			
Choice of Control Head	One handpiece (5-hole), accommodates 1 add'l tubing..... 1700-021-CS	\$ 695		
Automatic control for handpieces	Add'l handpiece illumination tubing (5-hole) 0017-615-CS	265		
asepsis 3-way air/water syringe with tubing, 5' (Euro tubing 3')	One handpiece (6-pin), accommodates 1 add'l tubing..... 1700-019-CS	695		
handpiece tubing, 4-hole Midwest, asepsis style, 5' (Euro tubing 3')	Add'l handpiece illumination tubing (6-pin) 0017-603-CS	190		
master control block (with pneumatic or pinch valves)	■ Ancillaries use existing handpiece position (see pages 4-5)			
master on/off	Write in selection here.....			
individual handpiece pressure adjustment	Write in selection here.....			
individual water coolant adjustment	Additional handpiece position, wet (4195/4485) 1100-189-CS	405		
non-retraction water coolant	Additional handpiece position, air drive only (4195/4485) 1100-190-CS	195		
air coolant adjustment	Assistant's Hygiene Instrumentation			
handpiece flush	(control head mount) (4195/4485) 1408-CS	320		
handpiece exhaust filter/oil collector	Swivel tray, arm-mounted (standard for Euro) 1110-468	335		
single handle brake release function	Chrome accents 9000-998-CS	135		
self-contained clean water system with quick disconnect, single bottle, 1½ liter	Quick disconnect with flow control, ¼" (specify air or water) 0016-455-CS	100		
wet/dry foot control	Additional umbilical (per foot) 1115-999-CS	30		
Flexible control arm with pneumatic brake	Routing valve kit with water shutoff & regulator/factory installation 1115-315-CS	235		
vertical travel 28" / horizontal travel 62"	Routing valve kit for use with cuspidor/factory installation 1115-316-CS	60		
External umbilical, 5'	Without chair adapter 1115-099-CS	-150		
*Deluxe utility center	ASSISTANT'S INSTRUMENTATION 1075FE	880		
air pressure regulator with filter	Additional autoclavable standard HVE			
master air shut-off valve	with tubing, 5' & autoclavable oral dry cup 1400-113-CS	85		
Chair adapter, fixed (specify chair make/model)	CUSPIDOR, BASIC ○ 1636FE	1250		
	with assistant's instrumentation ○ 1637FE	1720		
	Additional autoclavable standard HVE			
	with tubing, 5' & autoclavable oral dry cup (1637FE) 1400-113-CS	85		
	LIGHT, HALOGEN (mounting options, see pages 30-33) ○ 9080*	2090		
	or LED ○ 9078	3490		
	■ LIGHT POST with 2" post clamp (specify light make/model)			
	Light post, 35" (Forest) 1110-820*	260		
	Light post, 35" (Pelton & Crane) 1110-822	260		
	MONITOR MOUNT ASSEMBLY (monitor screen not included) (see page 34)			
	Ring mount (used with light post) 1300RM	990		
	Top post mount (used with monitor post) 1300TPM	990		
	• Monitor post (for top post mount) 1331-997	260		
	CHAIR, HYDRAULIC ○ 3900*	7690		
	■ One chair control included (specify type below) 3914-998-CS*	N/C		
	○ Footswitch control			
	○ Touchpad, chair mounted			
	○ Touchpad, unit-mounted			
	■ Additional chair control (limit 2 on chair, limit 2 on unit) 3914-999-CS	500		
	○ Footswitch control			
	○ Touchpad, chair mounted			
	○ Touchpad, unit-mounted			
	Additional chair options (See page 25 for pricing, write selection here) _____			
	STOOL, DR'S, ECONOMY/DELUXE/COMFORT 6156C/6159C/6176	485/825/825		
	STOOL, ASST'S, ECONOMY/DELUXE 6157C/6160C	625/925		
	STOOL, ASST'S WITH BACK, DELUXE/COMFORT 6161C/6187	1100/1100		
	STOOL, OFFICE, COMFORT 6198	975		
	TOTAL _____			
	• Choose only one post when sold with unit			
	* Item included in package price			
	■ Often overlooked			
	DESIGNER FRIENDLY® SURFACE TREATMENTS			
Shipping weight estimate:	Unit plastic finishes: <input type="radio"/> Forest Grey <input type="radio"/> White			
4485FE Delivery Unit	95 lb / 44 kg			
4987FE Delivery Unit	100 lb / 46 kg			
4195FE Delivery Unit	90 lb / 42 kg			
Powder coat paint finishes:	<input type="radio"/> Forest Grey <input type="radio"/> White			
Package plastic and powder coat paint finishes:	<input type="radio"/> Forest Grey			

* When optional cuspidor added, utility center housing and cover change from deluxe to standard for two umbilical accommodations.

Shipping weight estimate:

4485FE Delivery Unit 95 lb / 44 kg
 4987FE Delivery Unit 100 lb / 46 kg
 4195FE Delivery Unit 90 lb / 42 kg

Cuspidors / Assistant's Instrumentation

Fixed Chair Post Mount, External Umbilical



MODEL 1636 \$1250



MODEL 1637 \$1720



MODEL 1075 \$880

STANDARD FEATURES

CUSPIDOR, BASIC 1636

Gravity drain cuspidor with swing away arm
timed bowl rinse

manual cup fill

removable ceramic bowl

Post mount adapter, 2"

External umbilical, 5'

Does not include post or chair adapter

CUSPIDOR, BASIC 1637

Same as 1636 above, with additional features:

Individual assistant's instrumentation holders

Central vacuum accessory package

central vacuum canister with disposable solids collector
asepsis 3-way air/water syringe with tubing, 5'
autoclavable standard HVE lever valve with tubing, 5'
autoclavable saliva ejector lever valve with tubing, 5'

Does not include post or chair adapter

POST MOUNT VACUUM TELESCOPING ARM 1075

ARM 1075

Post mount adapter, 2"

Central vacuum accessory package

post mounted central vacuum canister
with disposable solids collector

asepsis 3-way air/water syringe with tubing, 5'
autoclavable standard HVE lever valve with tubing, 5'
autoclavable saliva ejector lever valve with tubing, 5'

External umbilical, 8'

Does not include post or chair adapter

FACTORY OPTIONS

CUSPIDOR, BASIC MODEL CHOICE _____

Additional autoclavable standard HVE

with tubing, 5' & autoclavable oral dry cup (1637) 1400-113-CS 85

Quick disconnect with flow control, 1/4" (specify air or water) 0016-455-CS 100

Standard utility center housing & cover 1115-020 130

Standard utility center * 1115-001 475

Additional umbilical (per foot) 1115-999-CS 30

MOUNTING POSTS 2" diameter:

Control post with cap, 18" 1110-487 110

LIGHT POSTS: 44" 2" diameter

Forest model 1110-794 200

Pelton & Crane 1110-785 200

ASSISTANT'S INSTRUMENTATION 1075 880

Additional standard autoclavable HVE

with tubing, 5' & autoclavable oral dry cup 1400-113-CS 85

Assistant's air shut-off & regulator utility center w/o housing & cover 1115-008 225

Standard utility center housing & cover 1115-020 130

Additional umbilical (per foot) 1115-998-CS 20

PART NO

RETAIL

PRICE

PRICE:

Often overlooked

Shipping weight estimate: * All optional utility centers
1636 17 lb / 8 kg include: housing and cover, air &
1637 20 lb / 9 kg water pressure regulator with filter
1075 15 lb / 7 kg master air and water shut-off valve,
unless otherwise stated.

DESIGNER FRIENDLY® SURFACE TREATMENTS

Unit plastic finishes: Forest Grey White

Powder coat paint finishes: Forest Grey White


Cabinetry Mounting Configurations

Start with a strong foundation


- Your choice of $\frac{3}{4}$ " hardened steel 11"x 14 $\frac{1}{2}$ " (pre-installed on the PFR8144) or optional 10"x 24" mounting plate that bolts to the bottom of your counter top.
- The 6"x 6" leveling plate is next, then a set of needle roller bearings with friction adjustment.
- A super strong 2 $\frac{1}{2}$ " steel tube scissor arm assembly.
- Needle roller bearings are used where the two arm sections pivot.

see page 46-47 for ordering information

3386UC
3486UC




3386BC
3486BC




shown with optional storage organizer and tray

M&CC PFR8144




M&CC DR CABINET



see page 50 for ordering information

3235



OPTIONAL

STANDARD


2400

see page 51 for ordering information



MT600BS/MT610BS


see page 54 for ordering information



1284

1284 shown with optional clean water system

see page 52 for ordering information



2184



Rear Mount Side Delivery Systems

3000 Series



Shown with optional storage organizer and swivel tray
1110-179-CS and 1110-909-CS



MODEL 3386
Photo not shown

MODEL 3486

STANDARD FEATURES

ASSISTANT'S SWING DELIVERY SYSTEM 3486 SERIES

Central vacuum accessory package
arm-mounted central vacuum canister
w/disposable solids collector
asepsis 3-way air/water syringe w/tubing, 5'
autoclavable standard HVE lever valve w/tubing, 5'
autoclavable saliva ejector lever valve w/tubing, 5'
self-contained clean water system w/quick disconnect,
single bottle, 1½ liter
Assistant's quad holder assembly
Laminate work surface, 19" x 24"
Beneath counter rigid arm assembly heavy duty mounting
plate, 11" x 14½" (3386/3486BC)
Under cabinet floor mount delivery post (3386/3486UC)
Laminate wall mount panel, 18" x 24" (3386/3486WM)
Air pressure regulator with filter
Master shut-off valve
Umbilical, 7'

WORK SURFACE SYSTEM, NO CONTROLS 3386 SERIES

Laminate work surface, 19" x 24"
Beneath counter rigid arm assembly heavy duty mounting
plate, 11" x 14½" (3386/3486BC)
Under cabinet floor mount delivery post (3386/3486UC)
Laminate wall mount panel, 18" x 24" (3386/3486WM)

Mounting plate dimensions, see pages 73-74

Shipping weight estimate:
3386 100 lb / 46 kg
3486 115 lb / 53 kg



MODEL 3386BC \$1950

MODEL 3486BC \$2750

BENEATH COUNTER MOUNT

Specify cabinet manufacturer

MODEL 3386UC \$2390

MODEL 3486UC \$3190

UNDER CABINET FLOOR MOUNT

MODEL 3386WM \$1950

MODEL 3486WM \$2750

WALL MOUNT

► DESIGNER FRIENDLY® SURFACE TREATMENTS

Unit plastic finishes: Forest Grey

Powder coat paint finishes: Forest Grey

Laminate work surface colors

Formica®: Folkstone Oyster Gray Surf
 White Blanc

Laminate wall mount panel (WM)

Formica®: Folkstone Oyster Gray Surf
 White Blanc

Solid work surface upgrade colors 0023-003-SS \$325

Corian®: Antarctica

Gibraltar®: Arctic Melange Baja Melange

Meganite®: Cinder Granite Goldrush Gemstone
 Moonstone Platinum Granite
 Rocky Road Granite Sanibel Granite
 Silver Mist Starry Starry Night

FACTORY OPTIONS

DELIVERY SYSTEM	MODEL CHOICE	PRICE :
Additional standard autoclavable HVE with tubing, 5' & autoclavable oral dry cup (3486)	1400-113-CS	85
Additional assistant's quad holder bar	0024-149	150
Quick disconnect with flow control 1/4" (for water) (3486)	0016-454-CS	100
Quick disconnect with flow control 3/8" (for air) (3486)	0016-456-CS	100
► Storage organizer with bins & frame	1110-179-CS	420
► Swivel tray	1110-909-CS	335
Touchpad for 3900 chair (25 ft. Cable)	1700-081-1-CS	500
Standard utility center housing & cover (3486)	1115-020	130
Additional umbilical (per foot) (3486)	1115-999-CS	30
Routing valve kit with water shutoff & regulator/factory installation (3486)	1115-315-CS	235
Rectangular work surface, 17" x 19" in lieu of semi-round	0023-262-CS	N/C
Standard color solid work surface upgrade	0023-003-SS	325
Special color laminate or solid work surface upgrade...0023-003-X/0023-003-SSV	(call for price)	

TOTAL

► Often overlooked


Cabinet mounting configuration see page 45

Rear Mount Side Delivery Systems

Mobile Cart



Shown with optional storage organizer and swivel tray
1110-179-CS and 1110-909-CS



MODEL 3386MC \$2390
photo not shown

MODEL 5386 \$3190

MODEL 5886 \$5190

STANDARD FEATURES

ASSISTANT'S CART W/INSTRUMENTATION 5386

Central vacuum accessory package, swing mount arm-mounted central vacuum canister w/disposable solids collector asepsis 3-way air/water syringe w/tubing, 5' autoclavable standard HVE lever valve w/tubing, 5' autoclavable saliva ejector lever valve w/tubing, 5' self-contained clean water system w/quick disconnect, single bottle 1½ liter Assistant's holder quad assembly Lamine work surface, 19" x 24" Adjustable height frame, 30" to 42" Twin wheel hard floor casters Electrical outlet/UL External umbilical, 7' Assistant's utility center with housing and cover air pressure regulator with filter master shut-off valve

DUO SWING CART 5886

Same as above, with additional features:

Standard IC Control Head automatic control for three handpieces asepsis 3-way air/water syringe w/tubing, 4' handpiece tubing, 4-hole Midwest, asepsis style, 4' pneumatic master control block w/diaphragm master on/off individual handpiece pressure adjustment individual water coolant adjustment non-retraction water coolant air coolant adjustment handpiece flush handpiece exhaust filter/oil collector wet/dry foot control Standard utility center housing and cover Air pressure regulator with filter Master air shut-off valve

CART, NO CONTROLS 3386MC

Lamine work surface, 19" x 24" Adjustable height frame, 32" to 42"

Shipping weight:

3386 100 lb / 46 kg

5386 115 lb / 53 kg

5886 120 lb / 55 kg

FACTORY OPTIONS

CART

Handpiece Illumination System (5886)

MODEL CHOICE	PRICE :
One handpiece (5-hole), accommodates 1 add'l tubing.....	1700-021-CS \$ 695
Add'l handpiece illumination tubing (5-hole)	0017-615-CS 265
One handpiece (6-pin), accommodates 1 add'l tubing	1700-019-CS 695
Add'l handpiece illumination tubing (6-pin)	0017-603-CS 190

Additional standard autoclavable HVE

with tubing, 5' & autoclavable oral dry cup (5386/5886).....	1400-113-CS 85
Additional assistant's quad holder bar	0024-149 150
Quick disconnect with flow control ¼" (for water) (5386/5886).....	0016-454-CS 100
Quick disconnect with flow control ¾" (for air) (5386/5886)	0016-456-CS 100

Storage organizer with bins & frame	1110-179-CS 420
Swivel tray	1110-909-CS 335
Additional umbilical (per foot) (5386/5886).....	1115-999-CS 30
Routing valve kit with water shutoff & regulator/	

factory installation (5386/5886).....	1115-315-CS 235
Standard color solid work surface upgrade.....	0023-003-SS 325
Rectangular work surface, 17" x 19" in lieu of semi-round	0023-262-CS N/C
Special color laminate or solid work surface upgrade..	0023-003-X / 0023-003-SSV (call for price)

TOTAL

Often overlooked

DESIGNER FRIENDLY® SURFACE TREATMENTS

Unit plastic finishes: Forest Grey

Powder coat paint finishes: Forest Grey

Laminate work surface Formica®: Surf Folkstone Oyster Gray White Blanc

Solid work surface upgrade colors 0023-003-SS \$325


Corian®: Antarctica

Gibraltar®: Arctic Melange Baja Melange

Meganite®: Cinder Granite Goldrush Gemstone Moonstone Platinum Granite

Rocky Road Granite Sanibel Granite Silver Mist Starry Starry Night


Mobile Cart



MODEL 5350 \$3590



*Shown with optional swivel tray,
arm mounted 1110-731*



MODEL 5355 \$3910

STANDARD FEATURES

DOCTOR STANDARD IC CART 5350

Standard IC Control Head

- automatic control for three handpieces
- asepsis 3-way air/water syringe with tubing, 4' handpiece tubing, 4-hole Midwest, asepsis style, 4' pneumatic master control block with diaphragm master on/off
- individual handpiece pressure adjustment
- individual water coolant adjustment
- non-retraction water coolant
- air coolant adjustment
- handpiece flush
- handpiece exhaust filter/oil collector
- self-contained clean water system with quick disconnect, single bottle 1½ liter
- wet/dry foot control

Adjustable height frame, 26" to 35"

Twin wheel hard floor casters

External 7' umbilical

Standard utility center housing and cover

Air pressure regulator with filter

Master air shut-off valve

DUO STANDARD IC CART 5355

Same as above, with additional features:

Central vacuum canister with disposable solids collector
autoclavable standard HVE lever valve with tubing, 5'
autoclavable saliva ejector lever valve with tubing, 5'

FACTORY OPTIONS

CART	MODEL CHOICE	PRICE :
Handpiece Illumination System		
One handpiece (5-hole), accommodates 1 add'l tubing.....	1700-021-CS	\$ 695
Add'l handpiece illumination tubing (5-hole).....	0017-615-CS	265
One handpiece (6-pin), accommodates 1 add'l tubing.....	1700-019-CS	695
Add'l handpiece illumination tubing (6-pin)	0017-603-CS	190
Additional standard autoclavable HVE		
with tubing, 5' & autoclavable oral dry cup (5355).....	1400-113-CS	85
Quick disconnect with flow control, ¼" (specify air or water).....	0016-455-CS	100
Swivel tray, arm-mounted.....	1110-731	335
Foot control support shelf, clamps to 2" base	0023-665	180
Additional umbilical (per foot).....	1115-999-CS	30
Routing valve kit with water shutoff & regulator/factory installation ...	1115-315-CS	235
		TOTAL

Often overlooked

Shipping weight estimate:

5350 50 lb / 23 kg

5355 55 lb / 25 kg



DESIGNER FRIENDLY® SURFACE TREATMENTS

Unit plastic finishes: Forest Grey

Powder coat paint finishes: Forest Grey

Rear Mount Side Delivery Systems

Side Cabinet & Wall Mount



MODEL 3235 \$4230
12" rigid arm cabinet mount



MODEL 3230 \$4360
22" rigid arm wall mount

STANDARD FEATURES

DOCTOR STANDARD IC DELIVERY SYSTEM 3235/3230

Standard IC control head
automatic control for three handpieces
asepsis 3-way air/water syringe with tubing, 5' handpiece tubing, 4-hole Midwest, asepsis style, 5' pneumatic master control block with diaphragm master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air coolant adjustment
handpiece flush
handpiece exhaust filter/oil collector
single handle brake release function
self-contained clean water system with quick disconnect, single bottle, 1½ liter
wet/dry foot control
Flexible control arm with pneumatic brake (3235)
cabinet vertical mount 12" rigid arm & 18" flex arm
Flexible control arm with pneumatic brake (3230)
wall vertical mount 22" rigid arm & 18" flex arm
External umbilical, 5'
Standard utility center housing and cover (3230)
Air pressure regulator with filter
Master air shut-off valve

FACTORY OPTIONS

DELIVERY SYSTEM	MODEL CHOICE	PRICE :
Handpiece Illumination System		
One handpiece (5-hole), accommodates 1 add'l tubing	1700-021-CS	\$ 695
Add'l handpiece illumination tubing (5-hole)	0017-615-CS	265
One handpiece (6-pin), accommodates 1 add'l tubing	1700-019-CS	695
Add'l handpiece illumination tubing (6-pin)	0017-603-CS	190
Assistant's hygiene instrumentation (control head mount)	1408-CS	320
Quick disconnect with flow control, ¼" (specify air or water)	0016-455-CS	100
Swivel tray, arm-mounted (over control head)	1110-468	335
Standard utility center housing & cover (3235)	1115-020	130
Additional umbilical (per foot)	1115-999-CS	30
Routing valve kit with water shutoff & regulator/factory installation ..	1115-315-CS	235
Laminate wall mount panel, 18" x 24"	1110-498	285
Ring mount rigid arm in lieu of top post mount		
22" without vertical mount	1110-536-RM	N/C
26" without vertical mount	1110-537-RM	N/C
26" rigid with 24" flex without vertical mount	1110-771-RM	N/C

TOTAL

► Often overlooked

For cabinet mounting configuration see page 45
Mounting plate dimensions, see page 75

Shipping weight estimate:
3235, 3230 50 lb / 23 kg



DESIGNER FRIENDLY® SURFACE TREATMENTS

Unit plastic finishes: Forest Grey

Powder coat paint finishes: Forest Grey

Rear Mount Side Delivery Systems

Rear or Side Beneath Counter Mount



MODEL 2400 \$2900

Beneath counter mount with swivel included
Countertop/work surface not included



Optional beneath counter slide 1700-140



Optional folding arm assembly 1110-115-CS

STANDARD FEATURES

DOCTOR STANDARD IC DELIVERY SYSTEM 2400
Beneath counter swivel mount (work surface not included)
Standard IC control head
automatic control for three handpieces
asepsis 3-way air/water syringe with tubing, 5' handpiece tubing, 4-hole Midwest, asepsis style, 5' pneumatic master control block with diaphragm master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air coolant adjustment
handpiece flush
handpiece exhaust filter/oil collector
self-contained clean water system with quick disconnect, single bottle, 1½ liter
wet/dry foot control
External braided umbilical, 7'
Air pressure regulator with filter
Master air shut-off valve

FACTORY OPTIONS

DELIVERY SYSTEM	MODEL CHOICE	2400	PRICE:
Handpiece Illumination System			
One handpiece (5-hole), accommodates 1 add'l tubing	1700-021-CS	\$ 695	
Add'l handpiece illumination tubing (5-hole)	0017-615-CS	265	
One handpiece (6-pin), accommodates 1 add'l tubing	1700-019-CS	695	
Add'l handpiece illumination tubing (6-pin)	0017-603-CS	190	
Assistant's hygiene instrumentation (control head mount)	1408-CS	320	
Quick disconnect with flow control, ¼" (specify air or water)	0016-455-CS	100	
Slide mount with swivel, 7½" extension	1700-140	280	
Folding arm assembly, 8" x 12" with wall mount	1110-115-CS	335	
Standard utility center housing & cover	1115-020	130	
Routing valve kit with water shutoff & regulator/factory installation ..	1115-315-CS	235	

TOTAL

Often overlooked

Cabinet mounting configuration, see page 45
Mounting plate dimensions, see page 75

Shipping weight estimate:
2400 28 lb / 13kg



DESIGNER FRIENDLY® SURFACE TREATMENTS

Unit plastic finishes: Forest Grey

Powder coat paint finishes: Forest Grey

Rear Mount Side Delivery Systems

Rear Beneath Counter Mount



MODEL 1284 \$1190

Shown with optional clean water system
1700-312 not included in price



MODEL 2184 \$3710

Mounts under counter or floating work surface

STANDARD FEATURES

HORIZONTAL ASSISTANT'S DELIVERY SYSTEM 1284

Central vacuum accessory package, swing mount
arm-mounted central vacuum canister
with disposable solids collector
asepsis 3-way air/water syringe with tubing, 5'
autoclavable universal HVE lever valve with tubing, 5'
autoclavable saliva ejector lever valve with tubing, 5'
Horizontal surface mount
External umbilical, 7'

HORIZONTAL DUO DELIVERY SYSTEM 2184

Same as above, with additional features:

Automatic panel control for three handpieces
asepsis 3-way air/water syringe with tubing, 5'
handpiece tubing, 4-hole Midwest, asepsis style, 5'
pneumatic master control block with diaphragm
master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air coolant adjustment
handpiece flush
wet/dry foot control
Articulating left/right under counter handpiece positioning Often overlooked
Self-contained clean water system with quick disconnect,
single bottle, 1½ liter
External umbilical, 7'

FACTORY OPTIONS

DELIVERY SYSTEM	MODEL CHOICE _____	PRICE: _____
Handpiece Illumination System (2184)		
One handpiece (5-hole), accommodates 1 add'l tubing	1700-021-CS	\$ 695
Add'l handpiece illumination tubing (5-hole)	0017-615-CS	265
One handpiece (6-pin), accommodates 1 add'l tubing	1700-019-CS	695
Add'l handpiece illumination tubing (6-pin)	0017-603-CS	190
Additional universal autoclavable HVE position with tubing, 5' & autoclavable oral dry cup	1400-113-CS	85
Slide mount, 7½" extension	1700-142	280
Touchpad for 3900 chair (25 ft. Cable) (2184)	1700-081-I-CS	500
Standard utility center housing & cover	1115-020	130
Self-contained clean water system with quick disconnect, single bottle (1284)	1700-312	265
Assistant's air shut-off & regulator utility center without housing & cover (1284)	1115-008	225
Quick disconnect for air or water	dealer installed (see page 63)	
Routing valve kit for city water	dealer installed (see page 63)	

TOTAL _____

For work surface and mount order 3386 from pages 46-48

Mounting plate dimensions, see page 76

Panel cut-out, see page 76

Shipping weight estimate:

1284 12 lb / 5 kg

2184 30 lb / 14 kg



DESIGNER FRIENDLY® SURFACE TREATMENTS

Unit plastic finishes: Forest Grey

Powder coat paint finishes: Forest Grey



MODEL 2500 \$2395



Optional assistant's hygiene instrumentation 1408-CS



MODEL 2520 \$1990

STANDARD FEATURES

DOCTOR CONTROL 2500

Automatic control for three handpieces
asepsis 3-way air/water syringe with tubing, 7'
handpiece tubing, 4-hole Midwest, asepsis style, 7'
pneumatic master control block with diaphragm
master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air coolant adjustment
self-contained clean water system with quick disconnect,
single bottle, 1½ liter
wet/dry foot control
Remote mount holder assembly with control tubing, 6'
External umbilical, 6'
Air pressure regulator with filter
Master air shut-off valve

DOCTOR CONTROL 2520

Automatic control for two handpieces
asepsis 3-way air/water syringe with tubing, 7'
handpiece tubing, 4-hole Midwest, asepsis style, 7'
pneumatic master control block with diaphragm
master on/off
individual handpiece pressure adjustment
non-retraction water coolant
air coolant adjustment
self-contained clean water system with quick
disconnect, single bottle, 1½ liter
wet/dry foot control
Remote mount holder assembly with control tubing, 6'
External umbilical, 6'
Air pressure regulator with filter
Master air shut-off valve

FACTORY OPTIONS

DELIVERY SYSTEM MODEL CHOICE _____

Handpiece Illumination System

One handpiece (5-hole), accommodates 1 add'l tubing.....	1700-021-CS	\$ 695
Add'l handpiece illumination tubing (5-hole)	0017-615-CS	265
One handpiece (6-pin), accommodates 1 add'l tubing.....	1700-019-CS	695
Add'l handpiece illumination tubing (6-pin)	0017-603-CS	190
Standard utility center housing & cover.....	1115-020	130
Assistant's hygiene instrumentation with extended 10 ½" holder bar (Drain Tubing 12') standard autoclavable HVE, saliva ejector with 2 holders & canister... 1408-CS	320	
Quick disconnect for air or water	dealer installed (see page 63)	
Routing valve kit for city water	dealer installed (see page 63)	

TOTAL

Often overlooked

Wall/cabinet cut-out dimensions, see page 76



DESIGNER FRIENDLY® SURFACE TREATMENTS

Unit plastic finishes:

Forest Grey

Powder coat paint finishes:

Forest Grey


Shipping weight estimate:
2500 18 lb / 8 kg
2520 15 lb / 7 kg

Rear Mount Side Delivery Systems

Rear or Side Mount Compact Controls



MODEL MT600 \$2095
Work surface not included



MODEL MT610 \$1990
Work surface not included

STANDARD FEATURES

DOCTOR CONTROL WITH 3 HP

Top Shelf Mount MT600TS
Beneath Shelf Mount MT600BS
with Swivel

DOCTOR CONTROL WITH 2 HP

Top Shelf Mount MT610TS
Beneath Shelf Mount MT610BS
with Swivel

All units include the following features:

Automatic control for three handpieces, or
Automatic control for two handpieces
asepsis 3-way air/water syringe with tubing, 4'
handpiece tubing, 4-hole Midwest, asepsis style, 4'
pneumatic master control block with diaphragm
master on/off
individual handpiece pressure adjustment
non-retraction water coolant
air coolant adjustment
self-contained clean water system with quick
disconnect, single bottle, 1½ liter
wet/dry type foot control
External umbilical, 12'
Air pressure regulator with filter
Master air shut-off valve

FACTORY OPTIONS

DELIVERY SYSTEM	MODEL CHOICE	PRICE :
Handpiece Illumination System		
One handpiece (5-hole), accommodates 1 add'l tubing.....	1700-021-CS	\$ 695
Add'l handpiece illumination tubing (5-hole)	0017-615-CS	265
One handpiece (6-pin), accommodates 1 add'l tubing.....	1700-019-CS	695
Add'l handpiece illumination tubing (6-pin)	0017-603-CS	190
Slide mount with swivel, 7 ½" extension (order BS).....	1700-140	280
Folding arm assembly, 8" x 12" with wall mount (order TS).....	1110-115-CS	335
Standard utility center housing & cover.....	1115-020	130
Assistant's hygiene instrumentation with extended 10 ½" holder bar (Drain Tubing 12')		
standard autoclavable HVE, saliva ejector with 2 holders & canister... 1408-CS		320
► Quick disconnect for air or water	dealer installed (see page 63)	
► Routing valve kit for city water	dealer installed (see page 63)	

TOTAL

► Often overlooked

Mounting plate dimensions, see page 77

Shipping weight estimate:
MT600, MT610 20 lb / 10 kg



DESIGNER FRIENDLY® SURFACE TREATMENTS

Unit plastic finishes: Forest Grey

Powder coat paint finishes: Forest Grey

Factory Options + Parts

Chair Adapters, Monitor Posts + Light Posts

PART NO	RETAIL
---------	--------

STANDARD CHAIR ADAPTERS - for 2" post

3900 fixed accommodates external umbilical.....	1115-195	\$ 430
3900 fixed sidebox accommodates internal umbilical.....	1115-530	565
3900 pivot accommodates internal umbilical.....	1115-528	565

MONITOR POSTS - 2" diameter with oilite bushing (used with 1300TPM)

30" with dual slot, fixed sidebox mount	1331-998	\$ 200
30" with single slot, fixed (includes post clamp)	1331-997	260

LIGHT POSTS - 2" with light bushing

24" fixed (includes post clamp)

Forest	1110-518	\$ 260
Pelton & Crane	1110-517	260

35" fixed (includes post clamp)

Forest	1110-820	260
Pelton & Crane	1110-822	260

35" fixed with slot (for monitor cord) (includes post clamp)

Forest	1110-045	260
--------------	----------	-----

44" fixed

Forest	1110-794	200
Pelton & Crane	1110-785	200

44" fixed with slot (for monitor cord)

Forest	1110-050	200
--------------	----------	-----

44" fixed with slot (sidebox option)

Forest	1110-900	200
Pelton & Crane	1110-890	200

44" fixed with dual slot (for monitor cord & sidebox access)

Forest	1110-170	200
--------------	----------	-----


PART NO	RETAIL
---------	--------



250 WATT (UNIVERSAL POWER SUPPLY)

- Factory option 0014-386 \$ 435
- Required when using two or more ancillary devices, or when using a Forest Operatory light and any Forest ancillary product (excluding the 3900 soft-touch massager power source)
 - Provides power for MULTIPLE electrical accessories such as: fiber-optics, heated syringe tubing, operatory lights, scalers, curing lights, and electric micro motors
 - Size: 7.5" - 5.5" - 2.5". Easily fits in all Forest floor boxes
 - Quick disconnect connectors for all of Forest ancillary products
 - Available in both 230 volt and 110 volt
 - Auto-resettable Klixon switch: inbuilt overcurrent and over temperature protection, eliminating the need for manually changing fuses or circuit breakers

Please note: built in air-electric switch, power source turns on/off with pilot air only



HANDPIECE ILLUMINATION SYSTEM, 5 HOLE ADAPTER WITH KRYPTON LAMP

Factory option	1700-021-CS	\$ 695
Sterling tubing	1700-052-STG	695
Shadow tubing	1700-052-SHW	695
Black tubing	1700-052-BLK	695
HANDPIECE ILLUMINATION SYSTEM, 6 PIN		
Factory option	1700-019-CS	695
Sterling tubing	1700-051-STG	695
Shadow tubing	1700-051-SHW	695
Black tubing	1700-051-BLK	695

- Dual channel outputs individually adjustable from 2VDC to 4VDC
- Accommodates one additional handpiece tubing
- Illumination delay (adjustable between 0-25 seconds) can be set individually for two handpieces
- Automatic activation of handpiece illumination with the dental unit foot control
- Simple installation into any new or existing dental unit
- Indicator LEDs confirms correct operation of the control module
- Powers all brands of optic handpieces and air-motors
- ISO standard handpiece terminal with stainless steel nut
- Asepsis style handpiece tubing, 7' long (4' to 5' installed)
- Power transformer, output 12VAC, 1.5 Amps
- Installation hardware and instructions



HANDPIECE ILLUMINATION TUBING, 5 HOLE ADAPTER WITH KRYPTON LAMP

Factory option	0017-615-CS	\$ 265
Sterling	0017-615-STG	265
Shadow	0017-615-SHW	265
Black	0017-615-BLK	265

- Used to replace power optic handpiece tubing
- ISO standard power optic handpiece terminal with stainless steel nut
 - 5 hole adapter with krypton lamp
 - Asepsis style handpiece tubing, 7' (4' to 5' factory installed)



HANDPIECE ILLUMINATION TUBING, 6 PIN

Factory option	0017-603-CS	\$ 190
Sterling	0017-603-STG	190
Shadow	0017-603-SHW	190
Black	0017-603-BLK	190

- Used to replace power optic handpiece tubing
- ISO standard power optic handpiece terminal with stainless steel adapter and nut (5-hole Midwest)
 - Asepsis style handpiece tubing, 7' (4' to 5' factory installed)

Factory Options + Parts

PART NO RETAIL




KRYPTON LAMP	0014-067	\$ 80
For use in handpiece illumination system, 3.5VAC, 720mA		
Used with 0017-692, 0017-592 and 0017-675 tubing		

ADAPTER WITH KRYPTON LAMP	0016-472	80
For use in handpiece illumination system, 3.5VAC, 720mA		
Used with 1700-019 or 1700-052 system and 0017-615 tubing		



HANDPIECE TUBING, ASEPSIS STYLE (VINYL) 4-Hole, straight midwest terminal		
Sterling	0017-810	\$ 35
Shadow	0017-510	35
Black	0017-398	35




STANDARD DELUXE AIR/WATER SYRINGE WITH QUICK-RELEASE TIP	1109-010	\$ 85
with sheathed tubing assembly, 5'		
Sterling	1107-755	135
Shadow	1107-755-SH	135
Black	1107-755-BLK	135
STRAIGHT DELUXE AIR/WATER SYRINGE WITH QUICK-RELEASE TIP	1109-025	105
with sheathed tubing assembly, 4'		
Sterling	110-585-STG	155
Shadow	1109-585-SHW	155
Black	1109-585-BLK	155
STANDARD DELUXE ONE BUTTON SYRINGE WITH QUICK-RELEASE TIP	1109-035	100
with sheathed tubing assembly, 5'		
Sterling	1109-485-STG	175
Shadow	1109-485-SHW	175
Black	1109-485-BLK	175



SYRINGE TUBING ASSEMBLY (POLYURETHANE), SHEATHED		
Sterling	0017-865	\$ 30
Shadow	0017-565	30
Black	0017-368	30



HEATED SYRINGE TUBING ASSEMBLY, ASEPSIS STYLE		
Sterling	0017-122	\$ 130
Shadow	0017-522	130
Black	0017-113	130

	PART NO	RETAIL
	DELUXE STANDARD HVE VACUUM VALVE Autoclavable with metal lever	1100-342 \$ 45
	DELUXE UNIVERSAL VACUUM VALVE Autoclavable with metal lever	1100-343 \$ 45
	DELUXE SALIVA EJECTOR VALVE Autoclavable with metal lever with grey tip..... with black tip.....	1100-341 \$ 40 1100-341-BLK 40
	ADDITIONAL HVE with tubing 5' and autoclavable oral dry cup Standard vacuum valve Factory option..... Sterling..... Shadow..... Black..... Universal vacuum valve Sterling..... Shadow..... Black..... ORAL DRY CUP..... Autoclavable polysulfone (beige). Fits all vacuum valves with 7/16" dia. tips.	1400-113-CS \$ 85 1400-113-G 85 1400-113-SH 85 1400-113-BLK 85 1400-119-G 85 1400-119-SH 85 1400-119-BLK 85 1400-010 15
	SMOOTH VACUUM TUBING, 1/2" ID Factory option..... Sterling..... Shadow Black	0017-933-CS \$ 4.25 / foot 0017-933 4.25 / foot .0017-533 4.25 / foot .0017-333 4.25 / foot
	SMOOTH VACUUM TUBING, 1/2" HEX ID Sterling..... Shadow Black	0017-932 3.25 / foot 0017-532 3.25 / foot 0017-332 3.25 / foot
	SALIVA EJECTOR TUBING (VINYL) Sterling..... Shadow Black	0017-930 \$.80 / foot 0017-530 .80 / foot 0017-330 .80 / foot

Factory Options + Parts

PART NO RETAIL




QUAD HOLDER ASSEMBLY, UNIVERSAL

White	0024-114	\$ 75
Black	0024-114-BLK	75



QUAD HOLDER ASSEMBLY, DOUBLE HVE

White	0024-127	75
Black	0024-127-BLK	75



DUAL ARTICULATING QUAD HOLDER ASSEMBLY

Factory option	0024-227-CS	\$ 235
4 standard, 4 universal, white	0024-227	310
4 standard, 4 universal, forest grey	0024-227-G	310
4 standard, 4 universal, Designer Friendly® (specify color)	0024-227-C	310
Factory option	0024-224-CS	235
8 universal, white	0024-224	310
8 universal, forest grey	0024-224-G	310
8 universal, Designer Friendly® (specify color)	0024-224-C	310



WET/DRY FOOT CONTROL

Forest Grey	1105-174-G	\$ 170
Designer Friendly® (specify color)	1105-174-C	170



CONTOURED UTILITY CENTER HOUSING AND COVER (8 x 17 3/16 x 17 3/16")

Factory option	3900-101	\$ 420
Forest Grey	3900-101-G	420
Shadow	3900-101-SH	420
Designer Friendly® (specify color)	3900-101-C	520

COMPLETE CONTOURED UTILITY CENTER

• Air and water pressure regulator with filter		
• Master air and water shut-off valves		
Factory option	3900-102	765
Forest Grey	3900-102-G	765
Shadow	3900-102-SH	765
Designer Friendly® (specify color)	3900-102-C	865



DELUXE UTILITY CENTER HOUSING AND COVER (7 x 13 5/16 x 12 11/16")

Rear umbilical access		
Forest Grey	1115-007-G	225
Shadow	1115-007-SH	225
Black	1115-007-BLK	225
White	1115-007-WH	225
Right umbilical access		
Forest Grey	1115-009-G	225
Shadow	1115-009-SH	225
Black	1115-009-BLK	225
White	1115-009-WH	225

COMPLETE DELUXE UTILITY CENTER

Right umbilical access		
Air and water pressure regulator with filter		
Master air and water shut-off valves		
Forest Grey	1115-005-G	640
Shadow	1115-005-SH	640
Black	1115-005-BLK	640
White	1115-005-WH	640

PART NO	RETAIL
---------	--------



STANDARD UTILITY CENTER HOUSING AND COVER (4 1/2 x 9 x 14")

Factory option	1115-020	\$ 130
Forest Grey	1115-020-G	130
Designer Friendly® (specify color).....	1115-020-C	130

COMPLETE STANDARD UTILITY CENTER

Air and water pressure regulator with filter		
Master air and water shut-off valves		
Factory option	1115-001	475
Forest Grey	1115-001-G	475
Designer Friendly® (specify color).....	1115-001-C	475

RECESSED UTILITY CENTER HOUSING AND COVER.....	1115-305	\$ 270
--	----------	--------



TRAY HOLDER ASSEMBLY, use with 1095 rear pivot assistant's instrumentation


Factory option	1110-680	\$ 560
Forest Grey	1110-680-G	560
Designer Friendly® (specify color).....	1110-680-C	560

Includes:

- Pivoting powder coated tray holder, 9 3/4" x 13 1/2"
- Clamp assembly adapts to 1 1/2" horizontal or vertical post

SLIDE WITH SWIVEL 7 5/8" EXTENSION, use with 1284BS, 2184VM, 2400MT's

Factory option (2400/MT's).....	1700-140	\$ 280
Forest Grey (2400/MT's).....	1700-140-G	280
Factory option (1284BS, 2184VM).....	1700-142	280



FOLDING ARM ASSEMBLY, use with 2400, MT's

Factory option	1110-115-CS	\$ 335
Forest Grey	1110-115-G	335

Includes:

- 8" folding arm with 1/2" pin and umbilical retainer
- 12" folding arm with 1/2" pin and umbilical retainer
- 4" x 4" wall mounting plate which accepts 1/2" pin
- Control mount with 1/2" pin



Factory Options + Parts

PART NO RETAIL



UNIVERSAL HEAVY DUTY MOUNTING PLATE

10" x 24" beneath counter.....	0010-084-BLK	\$ 320
11" x 14 1/2" beneath counter	0010-086-BLK	195
10" x 24" in lieu of 11" x 14 1/2".....	0010-084-CS	125



FLOATING WORK SURFACE ARM ASSEMBLY, 20"

Forest Grey	0023-815-G	730
Designer Friendly® (specify color).....	0023-815-C	730



SUPPORT SHELF

Factory option	0023-665	\$ 180
Forest Grey	0023-665-G	180

Provides adequate support for use of foot control while allowing cart mobility.



26" RIGID WITH 24" FLEXIBLE ARM



In lieu of top post mount factory option	1110-771-RM	N/C
22" rigid with 18" flexible arm in lieu of top post mount, factory option	1110-536-RM	N/C
26" rigid with 18" flexible arm in lieu of top post mount, factory option	1110-537-RM	N/C

Designed for use with 2" diameter post.



AIR AND WATER ACCESSORY PANEL Forest Grey	1300-030-G	\$ 165
Designer Friendly® (specify color).....	1300-030-C	165

- 1/4" female air quick disconnect
- Toggled shut off for air quick disconnect
- 1/4" female water quick disconnect
- Flow adjustment for water quick disconnect
- 4" x 3" flush mount panel

	PART NO	RETAIL
	QUICK DISCONNECTS	
flow control	1/4" female and male QD with flow control, factory installed	0016-454-CS \$ 100
	1/4" female QD with flow control, black / grey / shadow / white	0016-454-BLK / -G / -SH / -WH 100
	1/4" female and male QD with flow control, arm and control head mount, factory installed.....	0016-455-CS 100
	1/4" female QD with flow control, arm and control head mount, black / grey	0016-455-BLK / -G 100
	3/8" female and male QD with flow control, factory installed.....	0016-456-CS 100
	3/8" female QD with flow control, black / grey / shadow / white	0016-456-BLK / -G / -SH / -WH 100
	ROUTING VALVE KITS — Check with local authorities about installation requirements	
water shut-off and regulator	Factory installed with water shut-off and regulator.....	1115-315-CS \$ 235
	Factory installed without water shut-off and regulator.....	1115-316-CS 60
	Dealer installed includes; toggle valve, bracket, tubing, water shut-off and regulator black / grey / shadow / white	1115-310-BLK / 1115-310-G / 1115-310-SH / 1115-310-WH 225
	Dealer installed includes; toggle valve, bracket, and tubing black / grey / shadow / white	1115-312-BLK / 1115-312-G / 1115-312-SH / 1115-312-WH 60
	HIGH DENSITY BOTTLE	
0010-464/364	1 1/2 liter	0010-464 \$ 30
1100-431	1 liter	0010-364 20
1100-423	1 1/2 liter with quick disconnect adapter	1100-431 65
	1 1/2 liter with quick disconnect upgrade kit	1100-423 100
	For use with clean water system	
	• Opaque bottle	
	• High Density Polyethylene (HDPE)	
	AUTOMATIC HANDPIECE CONTROL BLOCK DIAPHRAGM	
	Tri-block	0022-020 \$ 3.25
	Quad-block	0022-024 7.35
	N/C lifetime warranty. Serial number and purchase order required. Freight and labor paid by dealer.	
	FOREST ALLEN WRENCH KIT.....	0008-019 \$ 18
	Includes one each of the following sizes: .050, 1/16, 5/64, 3/32, 7/64, 1/8, 9/64, 5/32, 3/16, 7/32, 1/4, 5/16	

Factory Options + Parts



	PART NO	RETAIL
	HOLDER for Automatic handpiece, 1/2" holder bar Forest Grey/Sterling 0010-836-G Shadow 0010-836-SH Black 00010-836-BLK White 0010-836-WH	\$70 70 70 70
	HOLDER, NOTCHED for automatic handpiece, fixed holder bar Forest Grey/Sterling 0010-835-G Shadow 0010-835-SH Black 0010-835-BLK White 0010-835-WH	\$70 70 70 70
	HOLDER, ELECTRIC automatic handpiece/universal vacuum/saliva ejector, normally open Forest Grey/Sterling 0010-856-G Shadow 0010-856-SH Black 0010-856-BLK White 0010-856-WH	\$75 75 75 75
	HOLDER, ELECTRIC automatic handpiece/universal vacuum/saliva ejector, normally closed Forest Grey/Sterling 0010-854-G Shadow 0010-854-SH Black 0010-854-BLK White 0010-854-WH	75 75 75 75
	HOLDER, NOTCHED ELECTRIC automatic handpiece/universal vacuum/saliva ejector, normally open, fixed holder bar Forest Grey/Sterling 0010-857-G Shadow 0010-857-SH Black 0010-857-BLK White 0010-857-WH	\$75 75 75 75
	HOLDER, NOTCHED ELECTRIC automatic handpiece/universal vacuum/saliva ejector, normally closed, fixed holder bar Forest Grey/Sterling 0010-855-G Shadow 0010-855-SH Black 0010-855-BLK White 0010-855-WH	75 75 75 75
	HOLDER for universal vacuum/syringe/saliva ejector Forest Grey/Sterling 0010-886-G Shadow 0010-886-SH Black 0010-886-BLK White 0010-886-WH	\$35 35 35 35
	HOLDER, NOTCHED for universal vacuum/syringe/saliva ejector, fixed holder bar Forest Grey/Sterling 0010-885-G Shadow 0010-885-SH Black 0010-885-BLK White 0010-885-WH	\$35 35 35 35
	HOLDER for standard vacuum Forest Grey/Sterling 0010-887-G Shadow 0010-887-SH Black 0010-887-BLK White 0010-887-WH	\$35 35 35 35

	PART NO	RETAIL
	HOLDER, NOTCHED for standard vacuum, fixed holder bar Forest Grey/Sterling 0010-888-G Shadow 0010-888-SH Black 0010-888-BLK White 0010-888-WH	\$ 35
	HOLDER for automatic standard vacuum Forest Grey/Sterling 0010-837-G Shadow 0010-837-SH Black 0010-837-BLK White 0010-837-WH	\$ 70
	HOLDER, NOTCHED for automatic standard vacuum, for fixed holder bar Forest Grey/Sterling 0010-838-G Shadow 0010-838-SH Black 0010-838-BLK White 0010-838-WH	\$ 70
	HOLDER for automatic accessory/universal vacuum Forest Grey/Sterling 0010-846-G Shadow 0010-846-SH Black 0010-846-BLK White 0010-846-WH	\$ 70
	HOLDER, NOTCHED for automatic accessory/universal vacuum, for fixed holder bar Forest Grey/Sterling 0010-847-G Shadow 0010-847-SH Black 0010-847-BLK White 0010-847-WH	\$ 70
	HOLDER, ELECTRIC for automatic standard vacuum, normally open Forest Grey/Sterling 0010-858-G Shadow 0010-858-SH Black 0010-858-BLK White 0010-858-WH	\$ 75
	HOLDER, ELECTRIC for automatic standard vacuum, normally closed Forest Grey/Sterling 0010-868-G Shadow 0010-868-SH Black 0010-868-BLK White 0010-868-WH	\$ 75
	HOLDER, NOTCHED ELECTRIC for automatic standard vacuum, normally open for fixed holder bar Forest Grey/Sterling 0010-859-G Shadow 0010-859-SH Black 0010-859-BLK White 0010-859-WH	\$ 75
	HOLDER, NOTCHED ELECTRIC for automatic standard vacuum, normally closed for fixed holder bar Forest Grey/Sterling 0010-869-G Shadow 0010-869-SH Black 0010-869-BLK White 0010-869-WH	\$ 75



General Information

Dimensions

PIVOT DELIVERY SYSTEMS




From tip of cuspidor assistant's holder to end of pivot light post is 120"














From tip of assistant's holder to end of pivot light post is 111"

1500PI



* shown with optional tray and asst's holder bar




FIXED WITH SIDEBOX DELIVERY SYSTEMS**FIXED WITHOUT SIDEBOX DELIVERY SYSTEMS****1636/1637 BASIC**

General Information


Dimensions

3900 CHAIR




UTILITY CENTERS


Contoured




Deluxe




Standard




Recessed




6176/6198 DOCTOR'S/OFFICE COMFORT STOOL




6187 ASSISTANT'S COMFORT STOOL




6159C DOCTOR'S DELUXE STOOL




6160C/6161C ASSISTANT'S DELUXE STOOL








6156C DOCTOR'S ECONOMY STOOL



6157C/6158C ASSISTANT'S ECONOMY STOOL




CHAIR-MOUNTED LIGHT**1110-227-CS
SHORT RIGID ARM****CEILING-MOUNTED LIGHT****TRACK MOUNT LIGHT/MONITOR**


General Information

Dimensions


STANDARD WALL & CABINET MOUNTED LIGHT





TALL WALL & CABINET MOUNTED LIGHT



CABINET MOUNT




OPTIONAL CABINET MOUNT





Vertical Exterior


Horizontal Top

1300RM











1300TPM



7010BC BENEATH COUNTER


- * work surface $\frac{9}{16}$ " thick (semi-round and rectangular)
- * shown with optional storage organizer and tray
- * shown with additional optional asst's holder bar
- * height adjustable work surface, 5", 7010

**7010WM MOUNTING PLATE****7010WM WALL MOUNT****7010BC MOUNTING PLATES****0010-086****0010-084-CS OPTION****7010WM LEVELING/
MOUNTING BOARD****1300-025-CS MONITOR
ASSEMBLY OPTION**

General Information


Dimensions

- * work surface $\frac{9}{16}$ " thick (semi-round and rectangular)
- * shown with optional storage organizer and tray
- * shown with additional optional asst's holder bar
- * height adjustable work surface, 5", 7010




semi-round work
surface standard
for 7010


1300-025-CS MONITOR ASSEMBLY OPTION




7010UC UNDER CABINET FLOOR MOUNT




7010UC UNDER CABINET FLOOR MOUNT LEVELING/MOUNTING PLATE







7010MC DUO CART



3386/3486BC & WM REAR BENEATH COUNTER MOUNT


3386 does not include assistants instrumentation

3386/3486BC HEAVY DUTY MOUNTING PLATE**3386/3486WM WALL MOUNT HUB****3386/3486WM LEVELING/MOUNTING BOARD****3386/3486BC & WM MOUNTING ARM**


* work surface $9/16"$ thick (semi-round and rectangular)

* shown with optional storage organizer and tray

* shown with additional optional asst's holder bar




semi-round work surface standard for 3386/3486

0023-262-CS RECTANGULAR WORK SURFACE OPTION FOR 3386/3486

General Information


Dimensions

3386/3486UC UNDER CABINET FLOOR MOUNT ASSISTANT'S




3386 does not include assistants instrumentation

- * work surface $\frac{9}{16}$ " thick (semi-round and rectangular)
- * shown with optional storage organizer and tray
- * shown with additional optional asst's holder bar




semi-round work surface standard for 3386/3486


3386/3486UC UNDER CABINET FLOOR MOUNT LEVELING/MOUNTING PLATE



0023-262-CS RECTANGULAR WORK SURFACE OPTION FOR 3386/3486








5386 & 3386MC ASSISTANT'S CART







3386 does not include assistants instrumentation

5886 DUO CART



5350/5355 STANDARD IC CART**5350/5355 WORK SURFACE****3235 CABINET SIDE DELIVERY****3230 WALL SIDE DELIVERY**


Minimum cabinet opening 18" deep x 30" wide

**2400 REAR MOUNT
BENEATH COUNTER****2400 LEVELING/
MOUNTING PLATE****1110-115-CS FOLDING ARM
SYSTEM OPTION****1700-140 SLIDE MOUNT
OPTION**


General Information

Dimensions


1284 REAR, BENEATH WORK SURFACE




1284/2184 LEVELING/MOUNTING PLATE




2184 REAR, BENEATH WORK SURFACE DUO UNIT




2184 WALL CUT-OUT




1700-142 SLIDE MOUNT OPTION, 1284/2184




2500 PANEL MOUNT




2500 WALL CUT-OUT










2520 PANEL MOUNT



2520 WALL CUT-OUT



MT600/MT610 COMPACT CONTROLS**MT600/610BS LEVELING/MOUNTING PLATE****1700-140 SLIDE MOUNT WITH SWIVEL OPTION, MT'S****1110-115-CS FOLDING ARM SYSTEM OPTION****1110-115-CS LEVELING/MOUNTING PLATE****0023-815 FLOATING WORK SURFACE ARM ASSEMBLY****CABINET MOUNTING PLATE**

General Information

Ordering Guide

ORDERING: When ordering Forest products, the following information should be included:

QUANTITY	ACCESSORIES OR OPTIONS
MODEL/PART #	DATE REQUIRED
DESCRIPTION	PURCHASE ORDER NUMBER
COLOR	METHOD OF SHIPPING

Deviation from our standard colors and Designer Friendly® options are subject to an upcharge.

STANDARD LAMINATE WORK SURFACE COLORS

(special laminates: call customer service for pricing and quantity discounts.)

FORMICA®:	Folkstone 927-58	Oyster Gray 929-58	Surf 923-58	White Blanc 949-58
-----------	------------------	--------------------	-------------	--------------------

STANDARD SOLID WORK SURFACE COLORS

(Special non-stock solid work surface colors: call customer service for pricing and quantity discounts.)

CORIAN®:	Antarctica
----------	------------

GIBRALTAR®:	Arctic Melange	Baja Melange
-------------	----------------	--------------

MEGANITE®:	Cinder Granite	Goldrush Gemstone	Moonstone	Platinum Granite
	Rocky Road Granite	Sanibel Granite	Silver Mist	Starry Starry Night

STANDARD UPHOLSTERY COLORS (Pearlized colors tend to wear similarly to soft leather)

*Ultraleather®:	Apricot	Baltic	Cave	Charcoal
	Diplomat Blue	Doe	Dove Grey	Eucalyptus
	Lemon	Papyrus	Parrot	Raven Wing
	Red	Riviera	Stone	
Pearlized Colors:	Abalone	Midnight	Pewter	

DESIGNER FRIENDLY® PAINT COLORS

	Arctic Silver	Bronze	Charcoal	Copper
	Forest Grey	Midnight	Pewter	Shadow
	White			

DESIGNER FRIENDLY® PLASTIC COLORS

Solid Colors:	Forest Grey	White	
Carbon Fiber Patterns:	Greywood	Lemans	Laser Leather

TUBING COLORS (no upcharge - colors configured by Forest)

Black	Shadow	Sterling
-------	--------	----------

STOOL PLASTIC COLOR (colors configured by Forest)

Black	Beige
-------	-------

CHROME ACCENTS (\$135 upcharge)

CUSTOM COLORS DEVIATING FROM STANDARD SUBJECT TO A 10%-15% UPCHARGE AND EXTENDED DELIVERY TIME

*ULTRALEATHER® IS A TRADEMARK OF ULTRAFABRICS® LLC.

TERMS

All prices unless otherwise listed are F.O.B. Hillsboro, Oregon U.S.A. Payments are net 30 days of invoice date. A service charge of 1-1/2% per month (annual rate 18%) will be charged on past due accounts, \$1.00 minimum. All prices are subject to change without notice.

RETURN MERCHANDISE

Authorized Forest dealers wishing to return merchandise for credit consideration must:

Obtain prior factory authorization.

Provide basic reason for return.

Provide serial number, purchase order or invoice number of original order.

Provide a purchase order for warranty items when necessary. Credit will be issued by Forest upon determination of warranty coverage in accordance with the stated Warranty Policy.

Return product within a period of 90 days—if product is not returned within time period, authorization becomes invalid.

Forest will approve return and credit amount after evaluation of product. If product is used or modified, the reason for return will be used to determine if replacement product or credit will be allowed.

Handling Charges:

25% or \$15.00 minimum restock fee will be charged on returned equipment approved for credit.

10% or \$15.00 minimum restock fee will be charged on new and unopened parts/components.

Returned items must be shipped freight prepaid.

WARRANTY

FOREST warrants its products to be free from defects in material and workmanship only. No other warranties are expressed or implied. This warranty shall extend for five years unless otherwise stated. Written notice of breach must be given to FOREST within this period. Buyer's remedy for breach of this warranty is limited to repair parts or replacement of items by FOREST. Factory repairs are covered for a period of 90 days. No claim for labor or consequential damages or freight will be allowed. FOREST is not responsible for ensuring equipment or installation conformity to local code or regulatory requirements. Check with local code or regulatory requirements before purchase and installation of all FOREST product.

Warranty is void if:

Items are carelessly used or improperly maintained or installed.

Damage results from using cleaning, disinfecting or sterilization chemicals and processes.

Damage is the result of freight damage and as such claims must be filed with the carrier by dealer.

Warranty term exceptions cover items that are normally expected service items and items not covered by Forest's warranty:

Handpiece illumination tubings, heated syringe tubings & assemblies and light bulbs are warranted for six months.

Power supplies, syringe buttons, HVE & SE valves are warranted for one year.

Water bottles, gas springs and handpiece tubings are warranted for two years.

Air & Water tubing connectors (T's), Distribution Block, Master Control Block, Flow Valves, 7000 Chassis are warranted for lifetime of dental unit (maximum of 20 years).

Bien-Air handpiece motors, Satelec scaler and Cavitron are covered under manufacturer's warranty. Two year upholstery warranty, excluding wear and tear. Installation of ancillary devices are warranted for two years.

CANCELLATION POLICY

Standard orders:

- Fax cancellation request required 10 working days prior to promised ship date, otherwise a 25% restock fee will be charged.
- Expedited ship date requests are non-cancellable.

CUSTOM ORDERS:

- Custom configurations are those that deviate from Forest's published product offering:

Orders are non-cancellable and non-returnable.

May be subject to an additional "Deviation" charge and/or a non-refundable engineering charge.

- Custom powder coat, custom plastic surface treatments and/or custom treatments such as multiple color combinations or custom colors may result in an upcharge of 10-15% per model and are non-cancellable and non-returnable.

General Information

Professional Sales Associates (PSA)

MILT BAISLEY

Cell: 203-645-7431
Fax: 203-248-0427
Email: mbaisley@profsales.com
Upstate NY-Albany, Syracuse, Newburgh

MICHEL CABANIE

Cell: 514-793-3368
Fax: 514-221-3368
Email: mcabanie@profsales.com
Quebec, Nova Scotia, New Brunswick,
Ottawa, Canada

STEVE CARLIN

Cell: 952-491-1223
Email: scarlin@profsales.com
MN, SD, ND, IA, MB Canada

CHARLIE CLARKE

859-268-8856
Cell: 859-509-6298
Fax: 859-263-0302
Email: cclarke@profsales.com
KY, Southern OH, IN, Western WV

JIM DOMITRZ

Cell: 617-921-4580
EFax: 603-218-7400
Email: jdomitrz@profsales.com
CT, MA, ME, NH, RI, VT

SKIP ELLENDER

Cell: 617-894-1224
Email: sellender@profsales.com
AZ, NM, Southern WY, CO, El Paso, TX

RICH ELLIS

Cell: 301-502-7424
EFax: 888-814-8630
Email: rellis@profsales.com
MD, VA, Southern WV, DC

BILL HIGHLAND

Cell: 810-623-7297
Fax: 517-540-9690
Email: bhighland@profsales.com
MI, OH, IN

MATT HIGHLAND

Cell: 901-606-3633
Fax: 901-756-0564
Email: mhightland@profsales.com
New Orleans, MS, Memphis, Nashville,
Little Rock

MATTHEW LANGDON

Cell: 516-316-5913
Fax: 516-706-2844
Email: mlangdon@profsales.com
Greater NYC, Westchester & Long Island,
Northern & Central NJ

MILTON LORRING

Cell: 216-396-2341
Fax: 216-803-9910
Email: mlorryng@profsales.com
Central and Northern OH, Western PA; Buffalo
& Rochester NY, Northern WV

JOHN MILLER

Cell: 503-758-9895
Email: jmiller@profsales.com
OR, WA, Northern ID, MT, Northern WY, AK

CHARLIE MARTIN

Cell: 330-421-4959
Fax: 330-769-2853
Email: cmartin@profsales.com
Central Regional Manager

PAT MOYLAN

Cell: 816-729-2900
Fax: 208-730-2759
Email: pmoylan@profsales.com
KS, NE, MO, Southern IL

PAUL MYERS

Cell: 214-704-5100
Fax: 253-369-1308
Email: pmyers@profsales.com
Northern TX, OK, Fayetteville AR,
Shreveport, LA

KEITH O'BERRY

Cell: 336-210-1429
Fax: 336-855-0184
Email: koberry@profsales.com
NC, SC, Savannah, GA

TOM O'BERRY

Cell: 404-307-4931
Email: toberry@profsales.com
GA (except Savannah & Valdosta), AL, East TN,
Pensacola, FL

RICK ORME

Cell: 916-837-3839
Fax: 916-984-8407
Email: rorme@profsales.com
Northern CA, NV (except L.V.) HI

SAM PAOLA

Cell: 310-569-1178
Fax: 206-770-7477
Email: spaola@profsales.com
Southern CA, Southern ID, UT, Las Vegas

MICHAEL PARK

281-251-9155
Cell: 281-702-3992
Fax: 281-251-9125
Email: mpark@profsales.com
Southern TX

BILL REINHARDT

847-382-2364
Cell: 847-840-3330
Email: breinhardt@profsales.com
Northern & Central IL, WI

AARON RENNER

Cell: 519-319-7457
Fax: 519-637-2466
Email: arenner@profsales.com
Ontario, Canada (except Ottawa)

ROBERT RODRIGUEZ

Cell: 305-458-6538
Fax: 305-723-5909
Email: rrodriguez@profsales.com
FL, Valdosta GA, Puerto Rico

DON ROGERS, III

Cell: 717-201-7322
Fax: 866-764-5980
Email: drogers@profsales.com
Eastern PA, Southern NJ, DE

SCOTT RUMBLE

Cell: 403-801-5390
Email: srumble@profsales.com
BC, AB, SK Canada

JIM SCHLOSSER

Cell: 909-215-3277
Email: jschlosser@profsales.com
Southern CA, Southern ID, UT, Las Vegas

VICTOR YOON

Cell: 778-688-8429
Email: vyoon@profsales.com
BC, AB, SK Canada

MODEL NUMBER + PAGE NUMBER INDEX

1075.....	44	3230	50, 75	5386	56, 74
1075FE.....	42	3235	50, 75	5886.....	56, 74
1095.....	10, 12, 14, 66	3386	45, 46, 47, 48, 73, 74	6156C	28, 68
1284.....	52, 76	3486	45, 46, 47, 48, 73, 74	6157C.....	28, 68
1300.....	12, 34, 70	3900	26, 27, 68	6159C	28, 68
1305T.....	31	4195FE.....	42, 67	6160C	28, 68
1500PI.....	16, 66	4195PI.....	14, 66	6161C.....	28, 68
1577.....	18, 20, 22, 40, 67	4195SE.....	40, 67	6176.....	29, 68
1580.....	10, 12, 14, 66	4195SI.....	22, 67	6187.....	29, 68
1590.....	10, 12, 14, 66	4485FE	42, 67	6198	29, 68
1636.....	44, 67	4485PI	10, 66	7010.....	35, 36, 37, 38, 71, 72
1636FE.....	42, 67	4485SE	40, 67	9078.....	30, 31, 32, 69
1637.....	44, 67	4485SI.....	18, 67	9079.....	33, 70
1637FE.....	42, 67	4987FE	42, 67	9080.....	30, 31, 32, 69
1679.....	18, 20, 22, 40, 67	4987PI	12, 66	9090.....	33, 70
2184.....	52, 76	4987SE	40, 67	MT600.....	54, 77
2400.....	51, 75	4987SI	20, 67	MT610.....	54, 77
2500.....	53, 76	5350	49, 75		
2520.....	53, 76	5355	49, 75		

FACTORY OPTION + PART NUMBER INDEX

0010-084	62	1110-731	16, 49	1700-089	25
0010-086	62	1110-771-RM	62	1700-140	51, 61, 75, 77
0010-310.....	8	1110-909-CS	36, 38, 46, 48	1700-142	61, 76
0014-386	57	1115-001.....	61	3900-101.....	25, 60
0016-454	63	1115-007	60	3900-102	60
0016-455	63	1115-020	61	3907-AS-NB.....	27
0016-456	63	1115-305.....	61	3907-AS-NS.....	26, 27
0017-603	7, 57	1115-315	63	3907-AS-WB	26
0017-615	7, 57	1115-316	63	3907-MC-WB.....	27
0023-665	62	1300-025	34	3907-PL-WB.....	26
0023-815	62, 77	1300-030	62	3907-317	25
0024-054	16	1400-113	59	3914-010-SH.....	25
0024-056	36	1408-CS	8, 53	3914-998-CS.....	8, 25
0024-224	60	1700-019-CS.....	7	3914-999-CS.....	8, 25
0024-227	60	1700-021-CS.....	7	9072-060	30
1110-115	51, 61, 75, 77	1700-051	57	9072-061	30
1110-179-CS	36, 38, 46, 48	1700-052	57	BIEN-AIR.....	6
1110-468	10, 14, 18	1700-140	51, 61, 75, 77	CAVITRON.....	6
1110-536-RM.....	62	1700-142	61, 76	SATELEC.....	6
1110-537-RM.....	62	1700-312	52		
1110-680	61	1700-079	25		

Surface Treatments

Forest is the only dental equipment with a WOW factor. It really is a kick watching a dentist react to Forest equipment when they first see it. A common reaction is “wow, that is beautiful.”

Carbon Fiber

Control Heads

- Standard IC - chassis
- Euro - lid
- Value - lid

Sidebox - cover

Basic Cuspidor - cover

Assistant's Instrumentation - vac canister cover

Powder Coat Paint

Dentists and interior designers love being able to put something besides white into their beautiful patient friendly decors.

Control Heads

- Standard IC - lid, brake handles & holder bar
- Euro - chassis, whip arms, brake handles & swivel tray
- Value - chassis, brake handles & holder bar

Sidebox - housing

Basic Cuspidor - housing & holder bar

Assistant's Instrumentation - arm

Ultraleather® (standard)

What if you could only order a white car with your sole color being your upholstery? You could live with it, but if given the choice, why would you? Once you see a Forest chair, delivery system and light package with color coordinated paint and upholstery in a well-designed treatment room, you wonder why it's taken so long for a manufacturer to offer dental equipment with a full range of surface treatments.

3900 Chair


- plush upholstery
- asepsis upholstery
- knee pillow
- deluxe and economy stools

Solid Surface

- 7000 Series rear mount, side delivery & mobile cart (standard)
- 3386/3486 rear mount side delivery (option)
- 5386/5886 mobile cart (option)

Laminate

- 7000 Series rear mount, side delivery & mobile cart (option)
- 3386/3486 rear mount side delivery (standard)
- 5386/5886 mobile cart (standard)



Raven Wing



Cave



Stone



Doe



Papyrus

Ultraleather®

Charcoal



Diplomat Blue



Dove Grey



Eucalyptus



Riviera



Pearl Midnight



Pearl Pewter



Pearl Abalone

Customize Your Treatment Room ~ Call for Details

Forest will:

- Powder coat paint your drawer and door pulls, trash receptacles, or other metal surfaces in your treatment room to match your delivery system paint color
- Source custom powder coat paint colors to match your provided sample
- Special order any Ultraleather® upholstery color shown on the Ultrafabrics® website (*from original or pearlized palette*)
- Silkscreen your practice logo/name on your control head
- Embroider your logo/name on your upholstery (*under chair toeboard cover*)



Baltic



Parrot



Lemon




Red



Apricot



Laser Leather



Lemans




Greywood

Carbon Fiber



Forest Grey




Shadow




White



Arctic Silver




Pewter




Midnight



Copper




Charcoal



Bronze

Powder Coat


Formica®



Surf




Folkstone



Rocky Road Granite




Silver Mist



Platinum Granite




Oyster Grey




White Blanc



Moonstone



Sanibel Granite



Goldrush Gemstone

Gibraltar®




Arctic Melange




Baja Melange


Corian®



Antarctica



Cinder Granite



Starry Starry Night

3	Ordering Ancillaries and Options 6 Control Head Accessories
10	Pivot Chair Mount Packages
18	Fixed Chair Mount w/ Sidebox Packages
25	Chair + Stools 25 Chair Accessories 26 3900 Hydraulic Patient Chair 28 Stools
30	Lights + Monitor Mounts 30 Operatory Lights 34 Monitor Mounts
35	Rear Mount Side Delivery Systems 35 Cabinetry Mounting Configurations 36 7000 Series
40	Delivery System, Chair Mounted 40 Fixed Mount with Sidebox, External Umbilical 42 Fixed Mount w/o Sidebox, External Umbilical 44 Post Mount Cuspidors & Assistant's Instrumentation
45	Rear Mount Side Delivery Systems 45 Cabinetry Mounting Configurations 46 Work Surface with & without Assistant's 48 Work Surface Cart with & without Delivery 49 Standard IC Carts 50 Standard IC Side Delivery 51 Standard IC Rear Beneath Counter Delivery 52 Beneath Counter Mount Duo & Assistant's 53 Panel Mount 54 MT Compact Control
56	Factory Options + Parts
66	General Information 66 Dimensions 78 Ordering Guide 79 Terms & Warranty 80 Professional Sales Associates Contact Information 81 Index

Forest Dental

Phone

- 800 . 423 . 3555
- 503 . 640 . 3012

Fax

- 503 . 693 . 9715

Digital

- sales@forestdental.com
- forestdental.com

Address

- 6200 ne cherry drive hillsboro oregon 97124 ISO 13485



* product videos available online at forestdental.com

