



# 2014 KaVo Collections

instruments & small equipment catalog



KaVo. Dental Excellence.

# KaVo Innovation

Over a Century of Dental Excellence

Innovation and quality are traditions that have kept KaVo at the forefront of the dental industry for more than 100 years. The company was founded in 1909 by Alois Kaltenbach in Berlin-Steglitz. It became known as Kaltenbach & Voigt (KaVo) when Richard Voigt joined the company in 1919 and is now a trusted partner for dentists and dental technicians worldwide.

KaVo's continued commitment to advancing dental technology for the betterment of your practice drives us every day. It is innovation and our passion for quality that has fueled our steady growth and world-wide recognition. You can see the outcome in our revolutionary new products, such as the COMFORTdrive, DIAGNOdent, MASTERtorque, and innovative programs such as TryKaVo.com. We believe these products, and our entire portfolio, are the best in the industry and we encourage you to experience KaVo for yourself.






In 1909 KaVo is founded by Alois Kaltenbach and in 1919 Richard Voigt joins the company.



In 1928, the first sterilizable handpiece is created by KaVo.



In 1958 KaVo launches its first turbine called Borden-Air 60.




In 1977 the KaVo MULTIflex is the first quick-release coupling.



In 1982 KaVo integrates the first light conductors into turbines.



In 2000 KaVo introduces the revolutionary DIAGNOdent, a caries detection aid.



In 2005, KaVo introduces GENTLEsilence LUX 6500 B – one of the quietest turbines in the world.



In 2008 KaVo introduces the GENTLEmini LUX 4500 B, a mini head, mighty power handpiece. Also launched is the GENTLEray 980 Premium Soft Tissue Diode Laser, the first diode laser with micro-pulsing power.



In 2009 KaVo celebrates 100 years of serving the dental industry, joining Pelton & Crane with 110 years of providing dental excellence. KaVo introduces the COMFORTdrive – KaVo's most innovative high-speed.



In 2010 KaVo introduces TRYKAVO.COM, the only risk free trial in the industry. KaVo also launches MULTI LED Illumination, an LED bulb that replaces your existing halogen lamp.



In 2012 KaVo introduces INTRAmatic LUX KL703 Short Motor, one of the shortest and lightest motors.



In 2014 KaVo introduces the MASTERtorque M8900 L, one of world's most powerful and advanced high-speed turbines.





## Electric Collections

Electric Collections Overview . . . . .	6-9
Master Series - COMFORTdrive . . . . .	10
Master Series - ELECTROtorque TLC Systems . . . . .	11
Master Series - GENTLEpower Attachments . . . . .	12
Master Series - INTRA LUX Heads . . . . .	13-14
Expert Series - ELECTROtorque Plus Systems . . . . .	18-19
Expert Series - EXPERTmatic Attachments . . . . .	19-20
Smart Series - INTRAmatic ES Systems . . . . .	21

## Air-Driven Collections

Air Driven Collections Overview . . . . .	22-25
Master Series - MASTERtorque & MASTERtorque Mini . . . . .	26
Expert Series - EXPERTtorque . . . . .	27
Smart Series - SMARTtorque . . . . .	28-29
INTRAflex Air Motor System . . . . .	30
MULTIflex Couplers . . . . .	31
INTRAmatic E Air Motor System . . . . .	32
Low-Speed Air Motors . . . . .	33

## KaVo Parts

Turbine Cartridges . . . . .	37
Turbine Back Cups And Tools . . . . .	38
Accessories for High-Speed Handpieces . . . . .	39

## Bundle Program

KaVo Bundle Programs . . . . .	40-43
--------------------------------	-------

## Maintenance Systems

QUATTROcare Plus . . . . .	44-48
Spray / Accessories . . . . .	49

## Hygiene

Hygiene 360° Overview . . . . .	50-51
KaVo SONICflex System . . . . .	52
SONICflex Scalers . . . . .	53
SONICflex - Hygiene & Periodontology . . . . .	54-55
SONICflex - Endodontics . . . . .	56
SONICflex - Minimally invasive caries therapy . . . . .	57-58
SONICflex - Inlay/Onlay preparation . . . . .	59
SONICflex - Surgery . . . . .	60
Accessories for SONICflex . . . . .	61
PROPHYflex System and Accessories . . . . .	62
Hygiene Polishing Instruments . . . . .	63
DIAGNODent Caries Detection . . . . .	64-65
RONDOflex plus . . . . .	66
Accessories for RONDOflex plus . . . . .	67
CORONAflex 2005 . . . . .	68

## Surgical Systems

Surgical Systems Overview . . . . .	70
INTRAsurg Systems . . . . .	71
Surgical Handpieces and Accessories . . . . .	72-74

## Soft Tissue Laser Accessories

GENTLEray 980 Premium Accessories . . . . .	75
---	----

## In-Office Laboratory

In-Office Laboratory . . . . .	76
--------------------------------	----

## Illuminating Systems

Illuminating Systems . . . . .	77
--------------------------------	----

## Water Treatment

Water Treatment System . . . . .	78
----------------------------------	----

## Product Index

Product Index . . . . .	79
-------------------------	----

## Maintenance Information

Maintenance Information . . . . .	80
-----------------------------------	----

## General Information

Terms of Sale / Registered Trademarks . . . . .	81
Warranty . . . . .	82
KaVo's Industry Award Recognition . . . . .	83


## Legend



NOTE: The prices listed in the catalog are retail prices effective January 1st, 2014. Prices are subject to change without notice. Possession of list price does not constitute an offer to sell.




# Electric Collections



## Master Series

The Pinnacle of KaVo Innovation.

### INTRA LUX Heads



Color Key  
● Blue = 1:1   ● Green = Reduction   ● Yellow = Specialty

## Expert Series

Superior KaVo Performance.

### EXPERTmatic Attachments

#### Optic



#### Non-Optic



## GENTLEpower Attachments



## ELECTROtorque Systems



ELECTROtorque TLC Torque Limit Control LUX with KL703 Brushless Electric Motor Kit  
1003.5840



COMFORTdrive Drive Full Operatory Package  
1007.4682

## Smart Series

Reliable and Affordable KaVo Technology.

### ELECTROtorque Systems



ELECTROtorque plus LUX with  
KL703 Brushless Electric Motor Kit  
1002.5015

To be used with  
air or electric motor  
with "E"  
type connection



INTRAmatic 23ES  
1004.5380




INTRAmatic 80ES  
1004.4375



INTRAmatic 19ES  
1004.4503

To be used with  
any air or electric  
motor with "E"  
type connection  
(181L, 181H,  
ELECTROtorque plus,  
ELECTROtorque TLC)






High-Speed

## Master Series

The Pinnacle of  
KaVo Innovation.

Featuring the COMFORTdrive  
200 XDR, KaVo's most  
innovative high-speed ever.  
Also featuring the award  
winning GENTLEpower 25LPR  
high-speed, 20LP attachment  
and INTRA LUX heads for a  
wide range of applications.



High-Speed

## Expert Series

Superior KaVo Performance.

Featuring the new  
EXPERTmatic E25L high-  
speed. Fast and precise  
working thanks to an even  
smaller head, snap-on  
facility and ceramic bearings.





### High-Speed

#### Smart Series

Reliable and Affordable  
KaVo Technology.

The KaVo Smart Series high-speed instruments are redesigned with a new surface finish for a perfect grip and no rough edges. Seamless stainless steel profile means there's nowhere for debris to accumulate and the push-button chuck is standard.



### Low-Speed

#### Smart Series

Reliable and Affordable  
KaVo Technology.

The KaVo Smart Series low-speed instruments can be thermally disinfected and sterilized in autoclaves at temperatures of up to 275°F and can also be washed in the thermodisinfector. Durable hygienic coating with textured surface for easy grip.



# Master Series

## The Pinnacle of KaVo Innovation



townie  
choice  
awards.2013

### COMFORTdrive 200 XDR Handpiece

High-speed contra-angle handpiece features precise performance, superior comfort, and allows for seamless integration into dental units. Speed range 30,000 to 200,000 rpm. Integrated micro-motor with 30W power, Cellular optic, Ceramic ball bearings, Exchangeable micro-filter. Only connects to the KaVo COMFORTbase coupling.

Part No. 1007.3570

\$1,699.00



### COMFORTdrive System



The COMFORTdrive system consists of the COMFORTdrive 200 XDR handpiece, the COMFORTtronic control box and the COMFORTbase coupling. The 200 XDR high-speed features a 2½ year warranty, to guarantee worry-free performance when maintained with KaVo QUATTROcare. Standard warranty is 2 years. The COMFORTbase and the COMFORTtronic carry a 3 year warranty.

Part No. 1007.3656

\$3,349.00

### COMFORTdrive Full Operatory Package



The Full Operatory Package includes three COMFORTdrive 200 XDR handpiece, one COMFORTtronic control box and one COMFORTbase coupling. The 200 XDR high-speed features a 2½ year warranty to guarantee worry-free performance when maintained with KaVo QUATTROcare. Standard warranty is 2 years. The COMFORTbase and the COMFORTtronic carry a 3 year warranty.

Part No. 1007.4682

\$5,549.00

### Parts and Accessories for COMFORTdrive

Part No.	Price	Description
1005.1707	\$92.00	QUATTROcare COMFORTdrive Adapter - Maintain your COMFORTdrive handpiece on your QUATTROcare handpiece maintenance system
1005.3154	\$20.00	COMFORTdrive Nozzle For maintenance of COMFORTdrive 200 XDR handpiece using the KaVo spray
1005.1773	\$15.00	Lamp / Bulb change tool
1007.5372	\$121.00	KaVo MULTI LED
1002.2928	\$98.00	High - pressure bulb (pack of 2)
1005.0327	\$3.25	COMFORTdrive O-Rings (set of 3)
1002.0271	\$20.00	Water filter
1002.0321	\$2.00	Filter tool
1007.9742	\$86.00	External Water Filter
1007.9736	\$11.00	Water Filter Element
1005.0169	\$2,436.00	COMFORTtronic Unit
1005.0120	\$384.00	COMFORTtronic Transformer
1005.0008	\$399.00	COMFORTtronic Tubing
1006.5370	\$11.00	COMFORTbase Cleaning Cap



WWW.KAVO.COM

## ELECTROtorque TLC Torque Limit Control LUX Brushless Motor Kit

Air/Electric conversion kit for any delivery unit. To be used with KaVo GENTLEpower LUX LP or other "E" type attachments.

townie  
choice  
awards 2011

### Includes:

- ELECTROtorque TLC module with 6.5' cable (1003.3490)
- ELECTROtorque TLC transformer with connector (1003.8542)
- Electric tubing (1003.5876)
- KL703 brushless electric motor (sterilizable)
- Mounting bracket with hardware (1002.0622)

### System Advantages:

- Brushless electric motor for ALL your high-speed, low-speed and torque control endodontic needs.
- No more searching and set up of a separate endo motor and no more wasted counter space – endo operation is available with the push of a button.
- Performs exact rotary endo procedures safely with torque control, auto-reverse and cellular optic advantages.
- Remote mountable faceplate for convenient mounting options and easy to access controls.
- KL703 is one of the world's shortest and lightest brushless electric motors (no brushes to exchange means less maintenance and downtime)
- Motor speed range 100-40,000 rpm (with various attachments 2-200,000 rpm).
- Endodontic torque 0.4 - 8.0 Ncm. Motor torque 3 Ncm.
- 3 year warranty (includes motor, tubing and control box).

### Also Available:

- Introduction manual (1003.5925) • Service manual (1003.8441)
- Quick reference guide (1003.8440)

### Technical:

- W: 5-1/8", D: 8-1/2", H: 3-1/2"
- Motor weight: 0.30 lb. or 0.14 kg
- Total weight: 2.30 lb. of 1.04 kg
- Electric motor for surgical procedures refer to Surgical Systems
- "KaVo INTRAsurg series" on page 70

Part No. 1003.5840      \$3,349.00

### Parts and Accessories for ELECTROtorque TLC

Part No.	Price	Description
1003.5710	\$509.00	ELECTROtorque TLC remote mounting face plate connection kit
1003.8537	\$56.00	3' cable for above kit
1008.2671	\$422.00	Tubing
1003.8542	\$347.00	Transformer 4881
1001.7115	\$5.50	Fuse
0200.6120	\$16.00	Motor O-ring
1003.5822	\$2.00	Large O-ring (base of motor)
1003.5656	\$5.50	Large O-ring (under bulb sheath)
1007.0150	\$2,393.00	INTRAmatic LUX KL703
1007.0146	\$23.00	KL703 Motor Bulb Ring
1002.2928	\$98.00	Light Bulb (pack of 2)
1007.9742	\$86.00	External Water Filter
1007.9736	\$11.00	Water Filter Element
1007.8474	\$102.00	KL703 motor light bulb

Note: KL703 will not work with ET TLC Serial # SN13651 and earlier.



WWW.TRY  
KAVO.COM



## GENTLEpower LUX 25LPR High-Speed Attachment

- Increasing 1:5
- Ceramic bearings
- Speed range up to 200,000 rpm
- Cellular optic
- Internal air/water coolant
- Friction grip push button auto chuck
- 3-port spray
- Works on electric or air motors that have a universal "E" type connection
- Plasmatec hygienic coating
- Improved angulation for easy access

Part No. 1007.3580

\$1,699.00

25LPR

Part No. 1002.0271

\$20.00

Water filter for 25LPA/25LPR

Part No. 1002.0321

\$2.00

Tool for filter



WWW.TRY  
KAVO.COM



## GENTLEpower LUX 20LP Contra Angle

- Transmission 1:1
- Speed range up to 40,000 rpm
- Cellular optic
- Internal air/water coolant

- Works on electric or air motors that have a universal "E" type connection
- Plasmatec hygienic coating
- Improved angulation

Part No. 1001.7453

\$920.00



## GENTLEpower LUX 7LP Contra Angle

- Reduction 2.7:1
- Speed range up to 15,000 rpm
- Cellular optic
- Internal air/water coolant

- Works on electric or air motors that have a universal "E" type connection
- Plasmatec hygienic coating
- Improved angulation

Part No. 1001.7452

\$1,120.00



## GENTLEpower LUX 29LP Contra Angle

- Reduction 7.4:1
- Speed range up to 5,500 rpm
- Cellular optic
- Internal air/water coolant

- Works on electric or air motors that have a universal "E" type connection
- Plasmatec hygienic coating
- Improved angulation

Part No. 1001.7454

\$1,250.00



## GENTLEpower LUX 10LP Straight Handpiece

- Transmission 1:1
- Speed range up to 40,000 rpm
- Works on electric or air motors that have a universal "E" type connection
- Plasmatec hygienic coating

- Cellular optic
- Internal air/water coolant
- For straight handpiece burs 2.334 to 2.35mm and contra angle burs 2.35mm (burstop must be used, and is included with purchase)

Part No. 1001.5660

\$1,325.00

Part No. 0524.0892

\$18.00 Burstop



\*\*AXIS Dental Toll Free 800.355.5063



### L68 B Latch Head – Operative

- Transmission 1:1
- For latch burs
- Push button auto chuck
- Blue dot
- Replaces 68LU

Part No. 1008.1834 \$565.00

### L52 Mini Latch Head

- Transmission 1:1
- Ideal for Mini Latch Burs (16mm) available from Axis Dental\*\*
- Push button Auto chuck
- Ideal for treatment in narrow areas or when treating children

Part No. 1008.1825 \$565.00

### L80 Friction Grip Head – Operative

- Transmission 1:1
- For friction grip burs
- Push button auto chuck
- Blue dot
- Replaces 80LD

Part No. 1004.4456 \$565.00

### L67 Reduction Latch Head – Operative

- Reduction 2:1
- Push button auto chuck
- For latch burs
- Green dot
- Replaces 67LH

Part No. 1008.1832 \$585.00

### L66 B Reduction Latch Head – Endodontics

- Reduction 3:1
- Push button auto chuck
- For latch type endo files or latch burs
- 360° rotation
- For use with KaVo 20LP and ELECTROtorque TLC Torque Limit Control
- Green dot
- replaces 66LU

Part No. 1008.1831 \$585.00

### L62 Reduction Latch Head – Endodontics

- Reduction 10:1
- Push button auto chuck
- For latch type endo files or latch burs
- 360° rotation
- Green dot
- Replaces 62LDN

Part No. 1008.1830 \$585.00

### L53 Oscillating Head – Endodontics

- The drive speed is converted into an oscillating movement of about 90 degrees
- Used with latch type endo files
- Max Speed 3,000 rpm, use with KaVo 29LP shank
- Now also for polishing attachments (2.35mm diameter)

Part No. 1008.1826 \$550.00



### L3 Reduction Head – Endodontics

- Speed max 3,000 rpm
- The drive speed is converted into an oscillating movement of about 90°
- To be used with hand NiTi files
- Push-button chuck
- Use with KaVo 29LP shank
- Replaces 3LD

Part No. 1008.1822

\$550.00



### L3 Y Reduction Head – Endodontics

- Speed max. 3,000 rpm
- The drive speed is converted into a stroke movement of about 0.4 mm
- To be used with hand NiTi files
- Push-button chuck
- Use with KaVo 29LP shank
- Replaces 3LDY

Part No. 1008.1823

\$550.00



### L61 G Orthodontics / Periodontics / Endodontics Head

- To be used with Intensiv files (Proxoshape) available through Axis Dental\*\*
- Free rotation
- For interproximal cleaning, polishing and preparation
- Buffered stroke (protects soft tissue)
- Elevation 0.8 mm
- Includes file ejector tool (0410.0634)
- Yellow dot
- Replaces 61LG

Part No. 1008.1828

\$550.00

Part No. 0410.0634

\$16.00 Ejector Tool for L61 G and L61 R

\*\*AXIS Dental Toll Free 800.355.5063



### L61 R Orthodontics / Periodontics / Endodontics Head

- To be used with Intensiv files available through Axis Dental\*\*
- 36 individual locking positions
- For interproximal polishing, stripping and preparation
- 10,000 rpm optimal drive speed
- 20,000 rpm maximum drive speed
- Buffered stroke (protects soft tissue)
- Elevation 0.8mm
- Includes file ejector tool (0410.0634)
- Recommended use with 7LP (on electric motor); 20LP/CH or INTRAFlex LUX 2307LN
- Yellow dot
- Replaces 61LRG

Part No. 1008.1829

\$550.00

Part No. 0410.0634

\$16.00 Ejector Tool for L61 G and L61 R

\*\*AXIS Dental Toll Free 800.355.5063



### L31 Reduction Prophy Head

- Reduction 2:1
- For snap on or screw in cups or brushes
- Green dot
- Replaces 31LR

Part No. 1008.1824

\$305.00

Part No. 0549.0602

\$20.00 Screw for L31



### L22 Micro Head – Surgical

- Transmission ratio 1:1
- With round drill size 2, also fits KaVo round
- Drill size Ø or U-inverted cone burs sizes 2 and 0
- Microhead with extremely small dimensions for microsurgical procedures such as apicectomies (through bone window)
- Maximum operating speed 20,000 rpm

Part No. 1008.1835

\$550.00

Part No. 1003.1523

\$144.00 Special Screwdrive



# TryKaVo.com

## Try the Best Handpieces, Risk Free



Step 1.  
**CHOOSE**



Choose a handpiece from  
our collections listed

Step 2.  
**TRY**



Try a handpiece risk free  
for 5 business days

Step 3.  
**RETURN**



After 5 day trial, return the  
handpiece. Like it, buy it.  
Contact 1-877-746-3995



Go to **TRYKAVO.COM** to start  
your free and easy trial today!





## KaVo Chuck Strength

KaVo handpieces feature a patented chuck system with exceptional chuck strength. The chuck strength will deliver bur concentricity for precise cutting and safety.

# EXPERT SERIES

Superior KaVo Performance



## ELECTROtorque plus LUX Motor Kit

Air/Electric conversion kit for any delivery unit. To be used with KaVo GENTLEpower LUX LP or other "E" type attachments.

### System Advantages:

Superior power & torque. Unique brushless motor (no brushes to exchange - no downtime). Quieter than brush motors. 3-year warranty (includes motor, tubing and control box). Two programmable speed settings. Water spray can be adjusted at motor. Motor speed range 100–40,000 rpm (with various attachments 2–200,000 rpm). Motor torque 3 Ncm. Integrated light and coolant. Forward/reverse speed adjustment. Small, easy to install control box. Unit can be calibrated to existing air pressure for optimized performance. Easy handpiece maintenance. Motor is autoclavable. Includes KL703 – one of the shortest and lightest electric motors with exclusive SAFEdrive and MULTI LED Light. Optional handpiece stand (0411.9501).

### Technical:

- Width: 5"
- Depth: 6"
- Height: 2"
- Motor weight: 99 grams, 2.97 oz
- Total weight: 29oz or 880g

- Electric motor for surgical procedures refer to Surgical Systems

Part No. 1002.5015      \$2,549.00

### Includes:

- ELECTROtorque plus module with cable
- ELECTROtorque plus transformer with connector (1003.8542)
- Electric tubing (1008.2671)
- Mounting bracket with hardware (1001.5113)
- KL703 Short Motor (1007.0150)
- SAFEdrive

Note: KL703 will not work with ET Plus Serial # SN18796 and earlier.

### Parts and Accessories for ELECTROtorque plus

Part No.	Price	Description
1007.0150	\$2,393.00	KL703 Short Motor
1001.9951	\$192.00	Flow meter
1008.2671	\$422.00	Tubing for ELECTROtorque plus box
1001.5113	\$39.00	Mounting bracket
1001.8748	\$2,675.00	Control unit only
0692.6891	\$49.00	Transformer power cord
0200.6120	\$16.00	O-ring
0200.6323	\$2.00	Large O-ring (base of motor)
0200.6308	\$5.50	Large O-ring (under bulb sheath)
1003.8542	\$347.00	Transformer 4881
1007.0146	\$23.00	KL703 Motor Ring
1007.8474	\$102.00	LED Bulb for KL703 Motor Only
1007.9742	\$86.00	External Water Filter
1007.9736	\$11.00	Water Filter Element
1001.9081	\$503.00	Tubing for KL 702 Motor

### Parts and Accessories for ELECTROtorque

Part No.	Price	Description
0905.8095	\$199.00	Remote panel kit. Includes 6' remote cable (0906.0753) & remote panel (0906.0754)
0201.8135	\$4.00	Water valve spring
0200.6120	\$16.00	Replacement small O-ring. 4 required for motors predating manufacturing date of 2/2008. After that date only 3 required.
0200.6042	\$4.00	Water valve small O-ring
0200.6038	\$4.00	Water valve large O-ring
0676.0901	\$16.00	Water valve ball bearing
0535.1661	\$49.00	Bulb holder
1002.2928	\$98.00	Light Bulbs (pack of 2)

## Parts and Accessories for K200 / 196 Motors

Part No.	Price	Description
1000.6717	\$2,300.00	K200 INTRA LUX Motor
0516.5850	\$749.00	Electric tubing for 196/200 motor. 6.8' length. Use with ELECTROtorque
0200.6120	\$16.00	Replacement small O-ring. 4 required for motors predating manufacturing date of 2/2008. After that date only 3 required
0676.1852	\$75.00	Sleeve for 196 Motor
0676.1922	\$55.00	196 motor water control ring
0905.3035	\$95.00	Carbon brush kit for 196 Motor. Includes 2 brushes (0676.0881) & 1 tool (0411.0041)
0676.0881	\$55.00	Carbon brushes for 196 Motor. 2 in box.
0411.0041	\$25.00	Tool for 196 Motor to change carbon brushes
1002.2928	\$98.00	Light bulb for 196, K200, KL700, KL701, K702 and 181L Motors. Pack of 2 bulbs.

## Expert Series - EXPERTmatic Attachments

**NEW**townie  
choice  
awards 2013**EXPERTmatic E25L**

- Speed Increasing 1:5
- Ceramic bearings
- Speed range up to 200,000 rpm
- Cellular optic
- Internal air/water coolant, Friction grip push button auto-chuck
- 3-port spray, EXPERTgrip, Universal "E" type connection

Part No. 1007.5550 \$1,369.00

Part No. 1002.0271 \$20.00 Water filter

Part No. 1002.0321 \$2.00 Tool for filter

**NEW**townie  
choice  
awards 2013**EXPERTmatic E25C**

- Speed Increasing 1:5
- Ceramic bearings
- Speed range up to 200,000 rpm
- Internal air/water coolant, Friction grip push button auto-chuck
- 3-port spray, EXPERTgrip, Universal "E" type connection

Part No. 1007.5551 \$1,169.00

Part No. 1002.0271 \$20.00 Water filter

Part No. 1002.0321 \$2.00 Tool for filter

**NEW**townie  
choice  
awards 2013**EXPERTmatic E20L**

- Transmission 1:1
- Cellular optic
- For latch type burs
- Speed range up to 40,000 rpm (electric motor)
- Speed range up to 20,000 rpm (air motor)
- Internal air/water coolant
- Push button auto chuck
- 1-port spray
- EXPERTgrip
- Universal "E" type connection

Part No. 1007.5540 \$1,299.00

**NEW**townie  
choice  
awards 2013**EXPERTmatic E20C**

- Transmission 1:1
- For latch type burs
- Speed range up to 40,000 rpm (electric motor)
- Speed range up to 20,000 rpm (air motor)
- Internal air/water coolant
- Push button auto chuck
- 1-port spray
- EXPERTgrip
- Universal "E" type connection

Part No. 1007.5541 \$919.00





townie  
choice  
awards 2013



### EXPERTmatic E15L

- Reduction 5.4:1
- Cellular optic
- For latch type burs
- Speed range up to 7,400 rpm (electric motor)
- Speed range up to 3,700 rpm (air motor)
- Internal air/water coolant
- Push button auto chuck
- 1-port spray
- EXPERTgrip
- Universal "E" type connection

Part No. 1007.5530      \$1,135.00



townie  
choice  
awards 2013



### EXPERTmatic E15C

- Reduction 5.4:1
- For latch type burs
- Speed range up to 7,400 rpm (electric motor)
- Speed range up to 3,700 rpm (air motor)
- Internal air/water coolant
- Push button auto chuck
- 1-port spray
- EXPERTgrip
- Universal "E" type connection

Part No. 1007.5531      \$919.00



townie  
choice  
awards 2013



### EXPERTmatic E61C

- EVA
- 10,000 rpm optimal drive speed (20,000 rpm max. drive speed)
- To be used with Intensiv files (ORTHOstrips or Proxoshapes)
- 36 individual locking positions
- For interproximal polishing
- Stripping and preparation
- Buffered stroke (protects soft tissue)
- Elevation 0.8 mm
- Internal air/water coolant
- 1-port spray
- EXPERTgrip
- Universal "E" type connection

Part No. 1008.4999      \$999.00  
Part No. 0410.0634      \$16.00 Ejector Tool for L61 G and L61 R



townie  
choice  
awards 2013



### EXPERTmatic E31C

- Reduction 5.4:1
- Ideal speed less than 5000 rpm
- Speed range up to 7,400 rpm (electric motor)
- Speed range up to 3,700 rpm (air motor)
- EXPERTgrip
- Universal "E" type connection

Part No. 1008.4998      \$1,075.00



townie  
choice  
awards 2013



### EXPERTmatic E10C

- Transmission 1:1
- For straight handpiece burs
- Speed range up to 40,000 rpm (electric motor)
- Speed range up to 20,000 rpm (air motor)
- Internal air/water coolant
- 1-port spray
- EXPERTgrip
- Universal "E" type connection

Part No. 1007.5560      \$789.00

# SMART SERIES

## Reliable and Affordable KaVo Technology



### 23ES Contra Angle Attachment

- Transmission 1:3
- Speed range up to 132,000 rpm
- 1-port spray
- Push button
- 100° head angle
- For friction grip burs
- For use with any "E" style connection
- Air or Electric motor
- Durable hygienic coating with textured surface for easy grip

Part No. 1004.5380 \$963.00

Part No. 0065.5188 \$16.00 External waterline



### 80ES Contra Angle Attachment

- Transmission 1:1
- Speed range up to 40,000 rpm
- 1-port spray
- Push button
- 100° head angle
- For friction grip burs
- For use with any "E" style connection
- Air or Electric motor
- Durable hygienic coating with textured surface for easy grip

Part No. 1004.4375 \$600.00

Part No. 0065.5188 \$16.00 External waterline




### 19ES Hygiene Handpiece

- Reduction 7.4:1
- Speed range up to 5,500 rpm
- For disposable prophy heads only
- For use with any "E" style connection Air or Electric Motor
- Durable hygienic coating with textured surface for easy grip

Part No. 1004.4503 \$499.00



# Air Collections



## Master Series

The Pinnacle of KaVo Innovation.



## Expert Series

Superior KaVo Performance.



## Smart Series

Reliable and Affordable KaVo Technology.



\*Non-optic fixed back handpiece does not use MULTIflex coupler

NEW

NEW



1390 MULTIflex LUX Coupler  
0553.1390



460LED MULTIflex LED Coupler  
1007.3201



460LE MULTIflex LUX Coupler  
1001.7599



465LED MULTIflex LED Coupler  
1007.3197

465LRN MULTIflex LUX Coupler  
0553.1550




457 MULTIflex Coupler  
0553.1210



Upgrade 460LE and 465LRN  
KaVo MULTI LED Bulb  
1007.5372






High-Speed

## Master Series

The Pinnacle of  
KaVo Innovation.

Featuring the new  
MASTERtorque handpiece  
– the best just got better  
with unique Direct Stop  
Technology for improved  
bur control and reduced  
drawback of contamination.  
Now with even higher power  
and same low noise level for  
the highest treatment quality  
and patient comfort.



High-Speed

## Expert Series

Superior KaVo Performance.

Featuring the new  
EXPERTtorque E679 L standard  
head and E675 L mini head  
handpieces with new innovative  
design and EXPERTgrip.





### High-Speed

#### Smart Series

Reliable and Affordable  
KaVo Technology.

Featuring the new  
SMARTtorque S619 L  
standard head and S615 L  
mini head handpiece with  
new ergonomic design and  
legendary KaVo durability.



### Low-Speed

#### Smart Series

Reliable and Affordable  
KaVo Technology.

The KaVo Smart Series  
INTRAmatic E Air Motor  
low-speed handpieces  
have a number of  
straightforward options  
designed around  
preferences and intended  
clinical application.  
Simply choose the right  
handpiece for the job, and  
the rest is easy.



# MASTER SERIES

The Pinnacle of KaVo Innovation

## MASTERtorque Handpieces



### MASTERtorque M8900 L

- Replaces 6500BR
- New: Most powerful (23 watts) KaVo high-speed with ceramic bearing
- New: Direct Stop Technology
- DST stops the bur in one second, increasing efficiency and improving safety during treatment
- DST reduces drawback of particles into the ball bearings for improved durability and hygiene
- Quietest KaVo air high-speed at 57 decibels
- Field replaceable water filter
- Improved 4-port spray for better visibility
- Plasmatec coating
- New angles for improved access
- Push button auto chuck
- 3rd generation cellular optics with 5-year warranty
- Automatic drive air adjustment - Use at minimum drive air pressure 38 psi (2.6 bar)
- MULTIflex LUX connection
- 370,000 rpm



### MASTERtorque Mini M4500 L

- Replaces 4500BR
- Mini Head (17 watts) KaVo high-speed with ceramic bearing
- Quietest KaVo mini head air high-speed at 59 decibels
- Field replaceable water filter
- 3-port spray
- Plasmatec coating
- New angles for improved access
- Push button auto chuck
- 3rd generation cellular optics with 5-year warranty
- Automatic drive air adjustment - Use at minimum drive air pressure 38 psi (2.6 bar)
- MULTIflex LUX connection
- 440,000 rpm

Part No. 1008.7800

\$1,550.00

Part No. 1002.0271

\$20.00 Water filter for M8900

Part No. 1002.0321

\$2.00 Tool for filter



Part No. 1007.1700

\$1,550.00

Part No. 1000.4823

\$33.00 Water filter for M4500L



\*2.5 years when maintained with QUATTROcare Plus.



# EXPERT SERIES

Superior KaVo Performance.

## EXPERTtorque Handpieces



townie  
choice  
awards 2013



townie  
choice  
awards 2013



### EXPERTtorque LUX E679 L

- Ceramic bearing
- Standard head (18 watt)
- MULTIflex LUX connection
- Push button auto w/ chuck
- 4-port spray
- Replaces 659BR
- 3rd generation cellular optics with 5-year warranty
- 36 psi (2.5 bar)
- 380,000 rpm

Part No. 1006.9200 \$1,120.00

Part No. 1002.0271 \$20.00 Water filter

Part No. 1002.0321 \$2.00 Tool for filter

### EXPERTtorque Mini LUX E675 L

- Ceramic bearing
- Miniature head (15 watts)
- MULTIflex LUX connection
- Push button auto-chuck
- 3-port spray
- Replaces 635BR
- 3rd generation cellular optics with 5-year warranty
- 33 psi (2.3 bar)
- 420,000 rpm

Part No. 1006.6400 \$1,120.00

Part No. 1006.0800 \$1,150.00 EXPERTtorque E675 LMW Midwest Stylus Coupler

Part No. 1002.0271 \$20.00 Water filter

Part No. 1002.0321 \$2.00 Tool for filter



KaVo. Dental Excellence.

# SMART SERIES

## Reliable and Affordable KaVo Technology

### SMARTtorque Handpieces



townie  
choice  
awards 2013



#### SMARTtorque LUX S619 L

- Fiber optic standard head (16 watts)
- Optimum balance between speed and torque
- MULTIflex LUX connection
- Push-button auto chuck
- 4-port spray
- 40 psi (2.8 bar)
- 400,000 rpm
- 1-year warranty

Part No. 1008.1641      \$715.00



townie  
choice  
awards 2013



#### SMARTtorque mini LUX S615 L

- Fiber optic miniature head (13 watts)
- Optimum balance between speed and torque
- MULTIflex LUX connection
- Push-button auto chuck
- 1-port spray
- 450,000 rpm
- 40 psi (2.8 bar)
- 1-year warranty

Part No. 1008.1643      \$715.00



townie  
choice  
awards 2013



#### SMARTtorque S619C

- Non-optic standard head (16 watts)
- Optimum balance between speed and torque
- MULTIflex connection
- Push-button auto chuck
- 4-port spray
- 40 psi (2.8 bar)
- 400,000 rpm
- 1-year warranty

Part No. 1008.1642      \$635.00

## SMARTtorque Handpieces



townie  
choice  
awards 2013



townie  
choice  
awards 2013

### SMARTtorque S609C

- Non-optic standard head (16 watts)
- Optimum balance between speed and torque
- Standard 4-hole connection
- Push-button auto chuck
- 4-port spray
- 40 psi (2.8 bar)
- 400,000 rpm
- 1-year warranty

Part No. 1008.1644      \$529.00

### SMARTtorque mini S605C

- Non-optic miniature head (13 watts)
- Optimum balance between speed and torque
- Standard 4-hole connection
- Push-button auto chuck
- 1-port spray
- 40 psi (2.8 bar)
- 400,000 rpm
- 1-year warranty

Part No. 1008.1645      \$529.00



# INTRAflex Air Motor System

townie  
choice  
awards 2009

135°C

## Cellular Optic INTRAflex LUX 2 Air Motor System 2325LN with Integrated Motor and Fixed Head

- Increase 1:4
- Speed range up to 140,000 rpm
- Cellular optic
- Internal air/water coolant
- MULTIflex LUX connection
- Friction grip push button auto chuck
- 3-port spray
- Use at minimum drive air pressure 43 psi (2.9 bar)
- 1-year warranty

Part No. 0535.6210

\$1,999.00

townie  
choice  
awards 2009

135°C

## Cellular Optic INTRAflex LUX 2 Air Motor System 2320LN Contra Angle with Integrated Motor

- Transmission 1:1
- Speed range up to 35,000 rpm
- Cellular optic
- Internal air/water coolant
- MULTIflex LUX connection
- Use at minimum drive air pressure 43 psi (2.9 bar)
- Fits INTRA LUX Heads
- 1-year warranty

Part No. 0535.6180

\$1,499.00

townie  
choice  
awards 2009

135°C

## Cellular Optic INTRAflex LUX 2 Air Motor System 2307LN Contra Angle with Integrated Motor

- Reduction 2.7:1
- Speed range up to 13,000 rpm
- Cellular optic
- Internal air/water coolant
- MULTIflex LUX connection
- Use at minimum drive air pressure 43 psi (2.9 bar)
- Fits INTRA LUX Heads
- 1-year warranty

Part No. 0535.6190

\$1,599.00

townie  
choice  
awards 2009

135°C

## Non-Optic INTRAflex 2 Air Motor System 2313LN Straight Handpiece w/ Integrated Motor

- Reduction 2.7:1
- Speed range up to 13,000 rpm
- Non optic
- Internal air/water coolant
- MULTIflex connection
- For handpiece and contra angle burs
- Comes with bur stop
- Maintain with chuck in closed position 43 PSI
- 1-year warranty

Part No. 0535.6220

\$1,299.00



NEW

**460LED MULTIflex LUX Coupler**

6-pin fiber optic coupler with MULTI LED bulb in coupler. Anti-retraction valve. (Same as 465LRN but without water spray regulator.) Use with LCM Illuminator and tubing. Not sterilizable.

Part No. 1007.3201 \$365.00



NEW

**465LED MULTIflex LUX Coupler**

6-pin fiber optic coupler with MULTI LED bulb in coupler. Spray regulator. Anti-retraction valve. Use with LCM Illuminator and tubing. Not sterilizable.

Part No. 1007.3197 \$465.00



NEW

**KaVo MULTI LED Bulb**

KaVo MULTI LED for all your 6-pin MULTIflex couplers. KaVo's innovative solution is an LED that can easily replace your existing halogen lamp. MULTI LED offers efficient and superior light, with color temperature similar to daylight (5,500 K and 25,000 LUX) at the bur tip. Operating voltage range 3.0 - 3.6 V DC. Nominal voltage 3.4 V DC. Current < 120 mA.

Part No. 1007.5372 \$121.00

WWW.TRY  
KAVO.COM**460LE MULTIflex LUX Coupler**

- 6-pin fiber optic coupler with light bulb in coupler
- Anti-retraction valve (Same as 465LRN but without water spray regulator.)
- Use with LCM Illuminator and tubing
- Not sterilizable

Part No. 1001.7599 \$315.00

**465LRN MULTIflex LUX Coupler**

- 6-pin connection
- Spray regulator
- Anti-retraction valve
- Use with LCM Illuminator and tubing
- Light bulb in coupler
- Not sterilizable

Part No. 0553.1550 \$405.00

WWW.TRY  
KAVO.COM

Fits product 465LRN / LED and 460LED

**1390 MULTIflex LUX Coupler**

- 5-hole connection
- Use with KaVo standard 5-hole fiber optic tubing or any standard 5-hole tubing
- Anti-retraction valve
- Use with LS2 Illuminator and tubing
- Not sterilizable

Part No. 0553.1390 \$354.00

**457 MULTIflex Coupler**

- 4-hole connection
- Non-optic
- Not sterilizable

Part No. 0553.1210 \$260.00

For Coupler parts see page 39. 2 pack halogen light bulb part #1002.2928. Retail \$98.00.



# INTRAmatic E Air Motor System



townie  
choice  
awards 2013



## 20E Low-Speed Handpiece

- Transmission 1:1
- Speed range up to 20,000 rpm
- Push button auto chuck for latch burs
- 100° head for superior access
- Durable hygienic coating with textured surface for easy grip
- Use with 181K Air Motor

Part No. 1003.1107 \$599.00



townie  
choice  
awards 2013



## 19E Hygiene Handpiece

- Reduction 7.4:1
- Speed range up to 2,700 rpm
- For disposable prophy heads only
- Durable hygienic coating with textured surface for easy grip
- Use with 181K Air Motor

Part No. 1004.6403 \$499.00



townie  
choice  
awards 2013



## 10E Low-Speed Handpiece

- Transmission 1:1
- Speed range up to 20,000 rpm
- For use with handpiece burs
- Durable hygienic coating with textured surface for easy grip
- Use with 181K Air Motor

Part No. 1003.6977 \$529.00



townie  
choice  
awards 2013



## 80E Low-Speed Handpiece

- Transmission 1:1
- Speed range up to 20,000 rpm
- Push button auto chuck for friction grip burs 1.6mm
- 100° head for superior access
- Durable hygienic coating with textured surface for easy grip
- Use with 181K Air Motor

Part No. 1004.6401 \$599.00



townie  
choice  
awards 2013



## 14E INTRAmatic Contra Angle

- Reduction 4:1 and for latch burs
- Max. 5,000 rpm
- Push button chuck
- 100° head angle for superior access
- Use with 181K Air motor

Part No. 1003.1108 \$545.00



townie  
choice  
awards 2013



## 31E Low-Speed Handpiece

- Reduction 4:1
- Speed range up to 5,000 rpm
- For use with snap-on / screw-in cups and brushes
- Extremely small head with 100° for superior access
- Durable hygienic coating with textured surface for easy grip
- Use with 181K Air Motor
- Comes with Prophy cup holder

Part No. 1003.6976 \$599.00

Part No. 1004.0596 \$3.50 Holder Tool

Master Series

**181L INTRA LUX Air Motor**dental  
townie  
choice  
awards 2007

Works with E attachments

- Optic
- Speed range up to 20,000 rpm
- Forward/reverse speed adjustment
- Internal air/water coolant
- 6-pin connection

- Fits all attachments optic or non-optic with Universal "E" type connection
- Anti-retraction system
- Use with KaVo LCM
- 31-38 PSI
- 1-year warranty

Part No. 0535.5210

\$1,899.00

Expert Series

**181H INTRA Air Motor**dental  
townie  
choice  
awards 2007

Works with E attachments

- Non-optic
- Speed range up to 20,000 rpm
- Forward/reverse speed adjustment
- Internal air/water coolant
- 4-hole connection

- Fits all non optic attachments with Universal "E" type connection
- Anti-retraction system
- Sterilizable
- 31-38 PSI
- 1-year warranty

Part No. 0535.5200

\$1,649.00

Smart Series

**181K Air Motor**REALITY  
Five Star Award  
REALITY'S  
CHOICES

Works with INTRAmatic E attachments

- Speed range up to 20,000 rpm
- 4-hole connection
- Non optic
- Forward/reverse speed adjustment
- Sterilizable

- 36-43 PSI
- No water
- No chip air
- 1-year warranty

Part No. 0535.5260

\$755.00

**INTRAmatic Air Motor Accessories**

Part No.	Price	Description
0200.6120	\$16.00	O-ring for 181L, 181H, 196, K200, 700KL, 701KL, 702KL Motors Requires 4 O-rings
0553.1872	\$11.00	181H, 181K gasket 4-hole for 181H Motor & non optic MULTIflex Coupler 465RN, 457, 459 & non optic high-speeds
1002.2928	\$98.00	Light bulbs for 181L, 196, K 200, 701KL, 702KL (pack of 2) 3.0-3.3 Volts      No warranty on bulbs
0246.5019	\$9.25	Set screw for 181H
0535.4242	\$55.00	Sheath for 181L
0535.1661	\$49.00	Bulb holder for 181L



## Genuine KaVo Turbines

Genuine KaVo turbines are individually balanced and electronically tested to meet the highest concentricity requirements. Properly balanced turbines leads to lower noise during use.



Kavo Dental Excellence

# KaVo Repair Parts

Quality in the Finest Details

When repairing equipment and handpieces, choose KaVo Original Parts. Original parts are manufactured to the same standards as your original purchase, which means better performance, durability, and safety. Alternate choices may not have all the innovation that KaVo puts into its products, like computer balanced turbines and two-piece carbide chuck sleeves. Restore your KaVo equipment to the standards by which it was originally manufactured.



## Turbine Cartridges



	<b>Product</b>	<b>Part No.</b>	<b>Price</b>	<b>Description</b>
	6500BR	1003.6883	\$400.00	1 auto chuck turbine 1 back cap tool
	4500BR	1006.5535	\$400.00	1 auto chuck turbine 1 back cap tool
	E679L	1007.9314	\$320.00	1 auto chuck turbine 1 back cap tool
	E675L	1007.9458	\$320.00	1 auto chuck turbine 1 back cap tool
	S619 / 609	1008.4165	\$215.00	1 auto chuck turbine 1 back cap tool
	S615 / 605	1008.4166	\$215.00	1 auto chuck turbine 1 back cap tool
	6500B	1003.9663	\$415.00	1 auto chuck turbine 1 back cap tool
Single Packs	4500B	1006.5711	\$415.00	1 auto chuck turbine 1 back cap tool
	6000B	1001.3889	\$445.00	1 auto chuck turbine 1 back cap tool
	659B	1004.5633	\$434.00	1 auto chuck turbine 1 back cap tool
	659BR	1000.2080	\$382.00	1 auto chuck turbine 1 back cap tool
	635B	0553.8591	\$414.00	1 auto chuck turbine 1 back cap tool
	635BR	0553.8671	\$414.00	1 auto chuck turbine 1 back cap tool
	647B, 649B	0553.8551	\$443.00	1 auto chuck turbine 1 back cap tool
	636CP, 646B	0553.8931	\$223.00	1 auto chuck turbine 1 back cap tool
	607C	0980.6649	\$254.00	1 auto chuck turbine 1 back cap tool
	640B, 630B	0553.6681	\$462.00	1 auto chuck turbine 1 back cap tool
Triple Packs	642B, 632, 643	0553.6901	\$462.00	1 auto chuck turbine 1 back cap tool
	630	0553.6671	\$417.10	1 friction grip turbine 1 back cap tool
	633, 632	0553.5551	\$474.00	1 Jacobs turbine 1 back cap tool
	625, 630	0553.5561	\$474.00	1 Jacobs turbine 1 back cap tool
	647B, 649B	0906.2179	\$1,258.95	3 turbines 1 back cap tool
	640B, 630B	0906.0300	\$1,258.95	3 turbines 1 back cap tool
	642B, 632B	0906.0310	\$1,258.95	3 turbines 1 back cap tool



## Accessories for High-Speed Handpieces

### Back Caps



Product	Part No.	Price	Description
6500B, 6500BR	1003.5525	\$179.00	Push button auto chuck
4500B, 4500BR	1006.1079	\$93.00	Push button auto chuck
E679 L	1007.8930	\$93.00	Push button auto chuck
E675 L, E675 LMW	1007.7051	\$93.00	Push button auto chuck
S619 L, S615C, S609C	1007.0854	\$63.00	Push button auto chuck
S615 L, S605C	1007.7277	\$63.00	Push button auto chuck
6000B	1000.4770	\$179.00	Push button auto chuck
659B	1000.2076	\$179.00	Push button auto chuck
635B	0553.8491	\$179.00	Push button auto chuck
647B, 649B, 646B, 636CP	0553.8061	\$179.00	Push button auto chuck
632, 642	0553.5081	\$179.00	Push button auto chuck
630B	0553.4041	\$179.00	Push button auto chuck
643B	0553.7631	\$179.00	Push button auto chuck
640B	0553.6511	\$179.00	Push button auto chuck
640FG	0553.8742	\$141.00	Friction grip
630	0553.1782	\$141.00	Friction grip

### Back Cap Removal Tools



Handpiece	Part No.	Price	Type
635 B/635 BR	0411.3063	\$30.00	Push button auto chuck
647B, 649B, 646B, 6000B, 6500B/BR, 636CP, E679L	0411.3053	\$18.00	Push button auto chuck
640, 642, 643, 630, 632	0411.1861	\$49.00	Push button auto chuck
4500B/BR, 5000B, E675L	1008.6133	\$25.00	Push button auto chuck
625C, 630	0411.0831	\$56.00	Friction grip
634 F4, 632, 642	0411.1701	\$56.00	Friction grip
S619 / 609	0411.2361	\$162.00	Push button auto chuck
S615 / 605	1006.3384	\$125.00	Push button auto chuck

**Coupler Bulbs**

Sterilizable. 175 hour of use.

Handpiece(s)	Type
465LRN	2 - MULTIflex
459AN	LUX bulbs
181L, 460LE	Voltage 3.0 - 3.3DC. 2.5 watts

Part No. 1002.2928	\$98.00	No warranty on bulbs
Part No. 1007.5372	\$121.00	No warranty on bulbs

**Coupler Tools**

Handpieces	Part No.	Price	Includes
all models	0411.1563	\$11.00	Coupler tightening tool
all models	0905.8088	\$416.00	Coupler straightening tool

**Accessories**

Couplers	Part No.	Price	Includes
465LRN / LED, 459AN, 460LE / LED	0553.9482	\$28.00	Screw-on bulb sheath Must order O-ring 0200.6033
all models	0200.6033	\$5.50	O-ring for bulb sheath
all models	1005.1033	\$22.00	Set of 5 O-rings: 3 large, 2 small
457, 459, 465RN	0553.1872	\$11.00	4-hole coupler gasket non optic high-speed handpieces, motors
1390	0553.4762	\$20.00	5-hole coupler gasket
465LRN / LED, 459AN, 460LE / LED	0553.5262	\$20.00	6-pin coupler gasket

**Bur Tools**

Handpiece(s)	Part No.	Price	Type
632, 633, 642	0411.3871	\$55.00	Friction grip (black)
630, 640, 625C	0411.3461	\$55.00	Friction grip (blue)
80G, 23D	0411.2501	\$49.00	Friction Grip (tan)
8A	0411.3521	\$196.00	For 8A Surgical Attachment

**Spray Repair Kits**

Handpiece(s)	Part No.	Price
640, 630 (A, B or C)	0410.0600	\$218.00
647 (A, B or C)	0410.0650	\$218.00
649	0410.0660	\$218.00
Ratchet	0411.2321	\$42.00
Water Filter, 6000B, 6500B, 6500BR, 4500BR, 4500B	1000.4823	\$33.00


**Accessories**

Item(s)	Part No.	Price
O-ring for High-Speed Handpiece	0200.6172	\$3.75
Cleaning needle	0410.0921	\$7.50
Cleaning brush	0229.3005	\$14.00



# KaVo Bundle Programs

Purchasing large equipment can be daunting, that's why KaVo provides handpiece bundles – a simple, easy choice to outfit operatories with new KaVo handpieces and hygiene equipment. Bundles are the best deal on handpieces when you purchase any dental delivery unit.



Package Part #	Part #	Retail Price	Quantity	Hygiene Non-Optic Package	Total Savings
1008.3301	1007.1606	\$3,200.00	1	SONICflex 2008	\$500
	1006.9926		1	PROPHYflex 3 Violet	
	0553.1210		1	457 Coupler	

Package Part #	Part #	Retail Price	Quantity	Hygiene Optic Package	Total Savings
1008.3302	1007.1604	\$3,300.00	1	SONICflex 2008 L	\$650
	1006.9926		1	PROPHYflex 3 Violet	
	1001.7599		1	460LE	

Package Part #	Part #	Retail Price	Quantity	Hygiene 360 Package	Total Savings
1008.6390	1003.2278	\$7,500.00	1	PROPHYwiz	\$700
	1006.9926		1	PROPHYflex 3 Violet	
	1007.1604		1	SONICflex quick 2008L set	
	1004.3400		1	DIAGNOdent pen	
	1004.8400		1	DIAGNOdent display	
	1001.7599		1	MULTIflex coupler 460 LE	

Part #	Retail Price	Quantity	Hygiene Package Buy	Get	Savings up to
1003.2278	\$819.00	1	PROPHYwiz	ONE Box of Kerr Perfect Pearl prophy angles FREE and ONE Box of Kerr mint prophy paste cups FREE All boxes contain 200 ct	\$700
1003.9882	\$1,715.00	1	PROPHYwiz 3-pack	THREE Boxes of Kerr Perfect Pearl prophy angles FREE and THREE Boxes of Kerr mint prophy paste cups FREE All boxes contain 200 ct	

Part #	Pricing	Quantity	Hygiene Package	Savings up to
1007.8889	\$1,600.00	1	181K Motor	\$650
		3	19E Hygiene Handpiece	

Purchase of dental unit is not required to take advantage of hygiene bundle pricing.




## Handpiece Bundle Program\*

	Part #	Pricing	Quantity	COMFORTdrive / High Speed	Total Savings**		
C1	1007.8108	\$3,855.00	3	COMFORTdrive 200 XDR	\$1,200	- Electric	
E1	1005.7646	\$3,855.00	3	GENTLEpower LUX 25LPR	\$1,200	- Hygiene	
E2	1010.7113	\$2,900.00	3	EXPERTmatic E25L	\$1,200	- Maintenance	
						- Air	
	Part #	Pricing	Quantity	Endo Premium / Endo Basic	Quantity	Low Speed / Endo Package	Total Savings**
E3	1007.1232	\$6,800.00	3	GENTLEpower LUX 25LPR	3	GENTLEpower LUX 20LP	\$2,750
					2	L68 B Latch Head - Operative	
					1	L66 B Reduction Latch Head - Endodontics	
E4	1005.7647	\$4,700.00	3	GENTLEpower LUX 25LPR	1	GENTLEpower LUX 20LP	\$1,900
					1	L66 B Reduction Latch Head - Endodontics	
	Part #	Pricing	Quantity	High Speed Package	Quantity	Low Speed Package	Total Savings**
E5	1007.1231	\$5,750.00	3	GENTLEpower LUX 25LPR	2	GENTLEpower LUX 20LP	\$2,300
					2	L68 B Latch Head - Operative	
E6	1007.8887	\$2,000.00			2	GENTLEpower LUX 20LP	\$950
					2	L68 B Latch Head - Operative	
	Part #	Pricing	Quantity	Expert Electric Premium	Quantity	Low Speed Package	Total Savings**
E7	1010.7114	\$5,525.00	3	EXPERTmatic E25L	2	EXPERTmatic E20L	\$2,300
					1	EXPERTmatic E15L	
	Part #	Pricing	Quantity	Expert Electric Basic	Quantity	Low Speed Package	Total Savings**
E8	1010.7115	\$4,725.00	3	EXPERTmatic E25L	2	EXPERTmatic E20L	\$1,950
	Part #	Pricing	Quantity	Expert Electric Low Speeds	Quantity	Low Speed Package	Total Savings**
E9	1010.7116	\$2,630.00	2	EXPERTmatic E20L	2	EXPERTmatic E15L	\$1,100
	Part #	Pricing		Add-ons			
D1 <sup>†</sup>	1004.3400	\$2,689.00		Diagnostic Package - DIAGNOdent pen Caries Detection Aid			
D2 <sup>†</sup>	0574.0500	\$2,689.00		Diagnostic Package - DIAGNOdent Classic Caries Detection Aid			
M1 <sup>†</sup>	1008.3805	\$3,198.00		Maintenance Package - QUATTROcare Plus			

\* Prices valid only with purchase of any manufacturing dental unit.

\*\* Total savings over individual instruments full retail price (rounded off).

† Purchase of dental unit is not required to take advantage of this bundle pricing.



6-Pin	Part #	Pricing	Quantity	High Speed 6-Pin Package	Quantity	Low Speed Package	Total Savings**
A1	1007.1233	\$5,100.00	1	MASTERtorque M8900 L	2	20E Low-Speed Handpiece	\$2,350
			2	MASTERtorque Mini M4500 L	1	10E Low-Speed Handpiece	
			1	MULTIflex LUX 460LED Coupler	1	181K Motor	
A2	1005.7648	\$3,600.00	3	MASTERtorque M8900 L			\$1,400
			1	MULTIflex LUX 460LED Coupler			
A3	1007.1226	\$3,600.00	3	MASTERtorque Mini M4500 L			\$1,400
			1	MULTIflex LUX 460LED Coupler			
A4	1010.7117	\$2,600.00	3	EXPERTtorque E679L			\$1,100
			1	MULTIflex LUX 460LED Coupler			
A5	1010.7118	\$2,600.00	3	EXPERTtorque E675L			\$1,100
			1	MULTIflex LUX 460LED Coupler			
5-Hole	Part #	Pricing	Quantity	High Speed 5-Hole Package	Quantity	Low Speed Package	Total Savings**
A6	1007.1234	\$5,100.00	1	MASTERtorque M8900 L	2	20E Low-Speed Handpiece	\$2,350
			2	MASTERtorque Mini M4500 L	1	10E Low-Speed Handpiece	
			1	MULTIflex LUX 1390 Coupler	1	181K Motor	
A7	1006.8293	\$3,600.00	3	MASTERtorque M8900 L			\$1,400
			1	MULTIflex LUX 1390 Coupler			
A8	1007.1227	\$3,600.00	3	MASTERtorque Mini M4500 L			\$1,400
			1	MULTIflex LUX 1390 Coupler			
A9	1010.7119	\$2,600.00	3	EXPERTtorque E679L			\$1,100
			1	MULTIflex LUX 1390 Coupler			
A10	1010.7120	\$2,600.00	3	EXPERTtorque E675L			\$1,100
			1	MULTIflex LUX 1390 Coupler			
4-Hole	Part #	Pricing	Quantity	Low-Speed 4-Hole Package	Quantity	Low Speed Package	Total Savings**
A11	1007.8888	\$1,700.00			2	20E Low-Speed Handpiece	\$750
					1	10E Low-Speed Handpiece	
					1	181K Motor	

\* Prices valid only with purchase of any manufacturing dental unit.

\*\* Total savings over individual instruments full retail price (rounded off).



## QUATTROcare Plus

QUATTROcare Plus Maintenance System is designed to clean, lubricate and purge your instruments with the simple push of a button. Chuck maintenance is now an easy step included with the QUATTROcare Plus.



# KaVo QUATTR0care Plus

Protect Your Investments

Proper cleaning of rotating instruments is essential before sterilization. If debris remains inside of the instrument, it can cause handpiece failure. Proper care means longer handpiece life, reduced repair costs, less frustration and elimination of costly downtime.



## The QUATTROCare Plus Advantage

### Easy to use

The user-friendly QUATTROCare Plus system is designed to clean, lubricate and purge your instruments with the simple push of a button. Loading and operation is easy as well as the automatic care and maintenance procedure.

### Safe controlled cleaning

Reduce errors with visual and audio signals during the completely controlled cleaning and purging procedure. With a front-loading enclosed maintenance chamber, the entire process is hygienic and thus eliminates mishandling or insufficient care and maintenance.

### Economical investment

Using the QUATTROCare Plus system means less wear, less repairs, less staff down time, smoother operation, and increased longevity of your handpieces – what an excellent return on investment! The QUATTROCare Plus can be used for virtually all manufacturers' high-speeds, low-speed contra-angles, straight handpieces, air driven scalers, and air motors.

### Fast cleaning cycle

Up to four (4) handpieces or air motors are cleaned, lubricated and purged in the same cycle. The care and maintenance of one complete set of handpieces in the same cycle is possible which means there's no need to change adapters between high-speed and low-speed attachments. Flexible programming allows you to choose and set total cycle time from 20-45 seconds per handpiece. Handpieces can be sterilized immediately after maintenance cycle; there is no need for separate purging of the handpiece.

### Chuck care maintenance

New chuck care system provides easy maintenance of handpiece chuck.

### Easy As...



1

Position the instruments...



2

... close the door ...




3

... press Start and everything is DONE.



4

Chuck care



# QUATTROcare Plus

## QUATTROcare Plus

Automatic handpiece maintenance and purging system. Correct maintenance of dental handpieces is critical to ensure top performance and longevity. KaVo QUATTROcare Plus is equipped with two MULTIflex adapters, one INTRAmatic (E-Type) adapter and one COMFORTdrive. Unit also includes one can of QUATTROcare Plus spray.

### System advantages:

- Up to four (4) handpieces or air motors are cleaned, lubricated and purged in the same cycle
- Flexible programming allows setting of total cycle time from 20-45 seconds per handpiece
- Can be used for virtually all manufacturers' handpieces, air driven scalers and air motors
- Unique spray formula expands inside each handpiece, creating a highly effective cleaning foam
- QUATTROcare Plus reduces mishandling or insufficient care and maintenance
- Completely controlled cleaning and purging procedure
- Chuck care system

### Technical:

- 120V, 50/60Hz, 20 Amps electrical connection

- Simply connects to 1/4" airline (58 to 87 psi air pressure ideal (max 145 psi) built in automatic pressure regulator controls flow)

### Dimensions:

- Height: 15", Width: 13-9/16", Depth: 11-1/2"
- Weight: 22 lbs

- Width with door ajar (right side - to exchange spray can): 16-15/16". Depth with flap open (front): 18-7/8"

Part No. 1008.3805

\$3,198.00



NEW



## Accessories for QUATTROcare

Part No.	Price	Description
0411.7991	\$96.00	MULTIflex adapter for KaVo Cellular Optic High-Speeds, SONICflex Air Scaler, INTRAflex Low-Speeds
0411.7981	\$96.00	E-Type adapter for INTRAmatic Low-Speed and surgical attachments
1005.1707	\$92.00	Adapter for COMFORTdrive
0411.7941	\$135.00	Intra Head adapter
0411.5532	\$71.00	Bumper sheath for INTRA head adapter
1000.8459	\$154.00	Adapter for Midwest "Shorty" attachment
1001.7860	\$154.00	Adapter for Bien Air High-Speeds
1001.7888	\$154.00	Adapter for W&H High-Speeds
1000.8786	\$154.00	Adapter for NSK High-Speeds
1000.1003	\$154.00	Adapter for SIRONA Electric T1 Instruments
0411.5412	\$14.00	U-clip for door, 2 needed
0200.6098	\$3.75	O-ring, Green, Small for front, 2 each
0200.6177	\$3.75	O-ring, Green, Large for front, 2 each
0200.6089	\$3.75	Low-speed O-ring
0324.4301	\$9.25	Blue hose for unit air hook up , 1/4"
0692.6891	\$49.00	Power Cord
0220.0977	\$43.00	Male Connector
0223.2786	\$2.00	Fuse for QUATTROcare
0411.4472	\$7.50	Union Nut
0411.7603	\$27.00	Adapter for chuck maintenance

**KaVo QUATTROcare plus Spray\* (6 cans)**


- 6 cans
- For cleaning and care of instruments / air motors in the KaVo QUATTROcare
- Designed specifically to withstand sterilization
- Relubrication after sterilization is unnecessary, eliminating the risk of recontaminating sterilized handpieces
- 500 ml

Part No. 1005.4524 \$355.00

**KaVo Spray\* (6 cans)**


- 6 cans
- Designed specifically to withstand sterilization
- Relubrication after sterilization is unnecessary, eliminating the risk of recontaminating sterilized handpieces
- 500 ml
- Nozzles sold separately

Part No. 0411.9660 \$369.00

**MULTIflex Nozzle**

- For maintenance of all KaVo MULTIflex Handpieces, INTRAFlex LUX Motors, SONICflex Scalers and the 181 M Air Motor

Part No. 0411.9921 \$20.00

**Quick Nozzle**

- For all INTRA and INTRAFlex LUX Heads
- 4-hole fixed back handpieces
- All KaVo air motors (except 181M), all handpieces with Universal "E" type connection and surgical handpieces

Part No. 0411.9931 \$20.00

Part No. 1005.3154 \$20.00 COMFORTdrive Nozzle

**KaVo CleanOil 2025 and Cleanpac**

- Keeps hands clean during cleaning and lubricating handpieces
  - Cleanpac traps excess lubricant while allowing view for flushing
  - Clean Pac (pack of 10)
- Part No. 0411.9691 \$50.00 (pack of 10)  
 Part No. 0411.7570 \$334.00 CleanOil 2025 2.5 L (84.5 fl. oz)

**Accessories**

Part No.	Price	Description
0411.9501	\$181.00	2151 Handpiece stand sterilizable stand for handpieces, attachments, motors and heads. Includes 3 inserts for MULTIflex instruments.
0411.9902	\$16.00	Insert for MULTIflex instruments
0411.9862	\$65.00	Drip Pads. 100 pads for handpiece stand Absorbs excess lubricant and/or condensation.

\*Hazmat Fees - refer to page 81 for Terms of Sale




# HYGIENE 360°

Everything Easy. Everything KaVo.



KaVo Hygiene 360° is a complete solution which helps you provide your patients with great dental care. Start the 360° cycle with the SONICflex Scaler for scaling, then PROPHYflex for tooth cleaning, PROPHYwiz for tooth polishing, and DIAGNOdent for caries and sub-gingival calculus detection. Complete your hygiene workflow on the Pelton & Crane Hygienist Stool.





### Step 1 Scaling

**Difficult stains?**

**The problem is solved so easily.**

Both thorough and gentle at the same time.

Your hygiene treatment center takes difficult cases in its stride. Tenacious calculus requires thorough, yet gentle removal. One of the most powerful scalers in its class - The KaVo SONICflex Quick 2008L easily conquers tough challenges.



KaVo SONICflex quick  
2008L



### Step 2 Tooth Cleaning

**Just conjure up a smile for the face.**

Your patients see success immediately.

The next step is about preventing caries, periodontitis and gingivitis. The main causes of these oral indicators, such as biofilm, deposits from smoking and discolorations from coffee and tea can be easily removed with the KaVo 2018 PROPHYflex 3 air polishing device.



KaVo 2018 PROPHYflex 3



### Step 3 Tooth Polishing

**You've got the skills and technology to maintain healthy dental enamel.**

Uncompromising outcome, means everything works great.

For a great finish to every professional tooth cleaning, the final polish and seal of the enamel is essential. KaVo offers the right tool for all indications, in each price segment.



KaVo PROPHYwiz



### Step 4 Caries Detection

**Your aid in detection of caries.**

The DIAGNOdent pen will show you where the decay lies.

It works quickly and reliably. The laser fluorescence detector within the DIAGNOdent pen, is a precise method for identifying fissure caries, proximal caries and periodontitis (with the Perio probe). A virtually unbeatable method with a 90% success rate. The wireless display ensures additional comfort, with a large status and measurement display.



KaVo DIAGNOdent pen +  
Display



### Step 5 Workflow Integration

**Your comfort is of the utmost importance.**

The Pelton & Crane hygienist stool.

With adjustable levers at your finger tips proper positioning is easy. Stay seated and on task during patient appointments.



Hygienist Stool



# KaVo SONICflex System

Variety Developed to a Fine Point

The right fit for every use. The wide range of options offered by the SONICflex will soon make it one of your most-used instruments. Choose from the right tip for your needs from the wide range of types. Over 50 tips for every application.





## SONICflex 2008/2008 L LUX Air Scaler



The new SONICflex 2008 has single turn tip installation and removal. Cellular optic glass rod (SONICflex 2008L) and Non-Fiber Optic (SONICflex 2008). Air driven. For removal of calculus and a wide range of other applications.

- Air/water coolant exits at tip
- Hygienic grip
- Three different performance levels
- Standard with three tips (Universal, Sickle, Perio) and tip removal tool
- MULTIflex connection - optic (2008 L), non-optic (2008)
- Automatic drive air regulator
- Operates at 36psi
- (2.5 bar)
- Set water at 15 psi
- Only compatible with KaVo SONICflex "A" tips

Part No. 1007.1604	\$1,999.00	SONICflex quick 2008 L Set
Part No. 1005.9310	\$1,645.00	SONICflex quick 2008 L (w/o tips)
Part No. 1007.1606	\$1,830.00	SONICflex quick 2008 Set
Part No. 1005.9311	\$1,545.00	SONICflex quick 2008 (w/o tips)



## SONICflex 2003/2003 L Air Scaler



Cellular optic glass rod (SONICflex 2003L) and Non-Fiber optic (SONICflex 2003). Air driven. For removal of calculus and a wide range of other applications.

- Hygienic grip
- Three different performance levels
- 50% more powerful than previous 2000L/N
- Standard with three tips (Universal, Sickle, Perio) and tip removal tool
- MULTIflex connection - optic (2003 L), non-optic (2003)
- Automatic drive air regulator
- Operates at 36psi
- (2.5 bar)
- Set water at 15 psi

Part No. 1000.8333	\$1,999.00	SONICflex LUX 2003L Set
Part No. 1000.4243	\$1,645.00	SONICflex 2003L (w/o tips)
Part No. 1000.8334	\$1,830.00	SONICflex 2003 Set
Part No. 1000.4246	\$1,545.00	SONICflex 2003 (w/o tips)



## Tip Removal Tool and Key for Adapter

For changing tips. Sterilizable. Torque feature prevents over tightening of tips. Can be used as protective cover when SONICflex is attached to delivery system.

Part No. 1000.4887	\$68.00	Tip Removal Tool
Part No. 0411.0892	\$18.00	Key
Part No. 0410.0911	\$15.00	Cleaning Tip



## Hygiene and periodontology

## SONICflex Scaler Tips\*



A variety of tips for tenacious calculus removal. Clean shallow pockets tips for effective and gentle calculus removal, hygiene.

Tip Kit (No. 5, 6, 7)		Part No.	Price	Description	Tip No.
		1000.4983	\$480.00	Three tips for supra- and subgingival calculus removal	5, 6, 7
		1005.9328	\$480.00	SONICflex scaler Tip Set A	5A, 6A, 7A
SONICflex Scaler Tips Universal		0571.5171	\$171.00	For calculus removal on all tooth surfaces	5
		1005.8949	\$171.00		5A
		0571.5181	\$171.00	For pedo use and areas of difficult access	6
		1005.8950	\$171.00		6A
		0571.5191	\$171.00	For subgingival scaling	7
		1005.8951	\$171.00		7A
		0571.5371	\$171.00	For subgingival Calculus	8
		1006.1953	\$171.00		8A

## SONICflex Scaler Tips\*



Wide spectrum vibrating brush. For gentle cleaning and polishing.

SONICflex Clean Sets and Tips		Part No.	Price	Description	Tip No.
		0571.0401	\$137.00	Clean Insert tip (w/o brush)	48
		1006.1982	\$137.00	Clean Insert tip (w/o brush)	48A
		0571.0001	\$273.00	Clean Set	(incl. Tip 48, plus 5 units of 6 brushes)
		1006.2026	\$273.00	Clean Set	(incl. Tip 48A, plus 5 units of 6 brushes)
SONICflex Clean Brushes		0571.0402	\$16.00	SONICflex Clean Brush	1
Clean Brushes fit into SONICflex tip 48		0571.0412	\$16.00	SONICflex Clean Brush	2
		0571.0422	\$16.00	SONICflex Clean Brush	3
		0571.0432	\$16.00	SONICflex Clean Brush	4
		0571.0472	\$16.00	SONICflex Clean Brush	5
		0571.0482	\$16.00	SONICflex Clean Brush	6

\*All tips also work with SONICflex 2000L and 2000N.

Note: Tips ending in suffix 'A' are compatible with SONICflex 2008 / SONICflex 2008L only.

## Hygiene and periodontology (cont.)

## SONICflex Perio



Different tips for minimally invasive, gentle and efficient cleaning of periodontal pockets up to 9 mm deep.

## SONICflex Perio tips



	Part No.	Price	Description	Tip No.
60	0571.7402	\$172.00	Straight - For anterior region	60
61	1006.1934	\$172.00		60A
62	0571.7412	\$172.00	Left - for furcations, molar regions	61
	1006.1935	\$172.00		61A
	0571.7422	\$172.00	Right- for furcations, molar regions	62
	1006.1936	\$172.00		62A

## SONICflex Perio accessories



Coupling for the supply of therapeutic fluids such as sterile saline solution  
(Part No. : 1000.6296) Retail Price \$61.00

## SONICflex Implant



Effective removal of sub-gingival plaque and calculus. No abrasion, no roughening. Disposable cleaning attachment to prevent the transmission of diseases; a very favorable price/performance ratio.

## SONICflex Implant tips



	Part No.	Price	Description	Tip No.
48	1003.8167	\$221.00	Incl. Tip 48, plus 10 cleaning pins	48
	1006.2027	\$221.00	Incl. Tip 48, plus 10 cleaning pins	48A
	0571.0401	\$137.00	Clean tip	48
	1006.1982	\$137.00	Clean tip	48A
	1003.8168	\$192.00	SONICflex Implant Refill Add-on kit with 30 pins	

## SONICflex Root planer




For easily processing of difficult-to-access or narrow furcations

## SONICflex Root planer insert – kit

	Part No.	Price	Description	Tip No.
	0571.0301	\$708.00	Root planer Tips	24, 25, 26, 27
	1006.2022	\$708.00	Root planer Tip Set A	24A, 25A, 26A, 27A

## SONICflex Root planer tips



	Part No.	Price	Description	Tip No.
24	0571.5621	\$201.00	Small Bud, Left - hard to reach furcations	24
25	1006.1957	\$201.00	Furcations	24A
26	0571.5631	\$201.00	Small Bud, Right - hard to reach furcations	25
27	1006.1959	\$201.00	Furcations	25A
	0571.5641	\$201.00	Small Bud, Universal	26
	1006.1961	\$201.00		26A
	0571.5651	\$201.00	Large Bud, Perio	27
	1006.1963	\$201.00		27A

\*All tips also work with SONICflex 2000L and 2000N.

Note: Tips ending in suffix 'A' are compatible with SONICflex 2008 / SONICflex 2008L only.

## Endodontics

## SONICflex Endo



SONICflex endo – a special line of tips for probing for canals.  
Removing soft materials during revision. Fashioning canal inlets.  
Shaping the pulp chamber for coronal root preparation

SONICflex Endo Tip Set A		Part No.	Price	Description	Tip No.
		1000.5828	\$600.00	SONICflex endo Tip Set	66, 67, 68, 69, 70
		1006.2031	\$600.00	SONICflex endo Tip Set A	66A, 67A, 68A, 69A, 70A
SONICflex Endo tips		1000.5825	\$192.00	Large bud, eval of access cavities	66
	66	1006.1992	\$192.00		66A
	67	1000.5822	\$192.00	Tip angle 125° expansion of larger channels	67
	68	1006.1994	\$192.00		67A
	69	1000.5823	\$192.00	Tip angle 117° Expansion of larger channels	68
	70	1006.1998	\$192.00		68A
	66	1000.5827	\$192.00	Removal of residual root canal fillings	69
	67	1006.1998	\$192.00		69A
	68	1000.5821	\$192.00	Tip angle 117° – for large channels	70
	69	1006.2000	\$192.00		70A

SONICflex Retro		Part No.	Price	Description	Tip No.
		1006.2034	\$924.00	Retro Tip Set	56, 57, 16, 17, 20, 21
		0571.0381	\$924.00	Retro Tip Set A	56A, 57A, 16A, 17A, 20A, 21A
SONICflex Retro tips		0571.5541	\$204.00	Cylindrical, small, left, Prep of channel	16
	16	1006.2036	\$204.00		16A
	17	0571.5581	\$204.00	Cylindrical, small, right, Prep of channel	17
	20	1006.2035	\$204.00		17A
	21	0571.5521	\$204.00	T-shape, small, left, Prep of the undercut	20
	55	1006.2037	\$204.00		20A
	56	0571.5561	\$205.00	T-shape, small, right, Prep of the undercut	21
	57	1006.2038	\$204.00		21A
	55	0571.7342	\$204.00	Anterior tooth	55
	56	1006.2039	\$204.00		55A
	56	0571.7322	\$204.00	Isthmus, left	56
	57	1006.2040	\$204.00		56A
	57	0571.7332	\$204.00	Isthmus, right	57
		1006.2041	\$204.00		57A

\*All tips also work with SONICflex 2000L and 2000N.

Note: Tips ending in suffix 'A' are compatible with SONICflex 2008 / SONICflex 2008L only.

## Minimally invasive caries therapy

## SONICflex Microinvasive



Special tip system for small, difficult-to-access lesions in the proximal region. Tips diamond-coated on one side allow direct access to initial lesions from an occlusal or buccal direction. Maintenance of healthy tooth substance.

	Part No.	Price	Description	Tip No.
SONICsys Micro kit	0571.0421	\$959.00	SONICsys Micro Tip Set	28, 29, 30, 31, 32, 33
SONICsys Micro Tip Set A	1006.2023	\$959.00	SONICsys Micro Tip Set	28A, 2A9, 30A, 31A, 32A, 33A
SONICflex SONICsys Microinvasive tip	0571.6741	\$164.00	Mesial Tip	28
	1006.1965	\$164.00		28A
	0571.6731	\$164.00	Distal Tip	29
	1006.1967	\$164.00		29A
	0571.6811	\$164.00	Small hemisphere, mesial	30
	1006.1969	\$164.00		30A
	0571.6801	\$164.00	Small hemisphere, distal	31
	1006.1971	\$164.00		31A
	0571.6831	\$164.00	Large hemisphere, mesial	32
	1006.1973	\$164.00		32A
	0571.6821	\$164.00	Large hemisphere, distal	33
	1006.1975	\$164.00		33A

## Accessories for SONICsys Micro Tip



1000.4887 \$68.00 For changing tips. Sterilizable. Torque feature prevents over tightening of tips. Can be used as protective cover when SONICflex is attached to delivery system.

0411.9101 \$95.00 Sterilization cassette

## SONICflex Bevel



Tips diamond-coated on one side for the precise bevelling of prepared small lesions. Optimization of adhesion. Minimization of the marginal gap.

	Part No.	Price	Description	Tip No.
SONICflex Bevel Tip Set A	1006.2030	\$249.00	Bevel Tip Set	58A, 59A
SONICflex Bevel Tips	0571.7392	\$105.00	Mesial Tip	58
	1006.1988	\$105.00		58A
	0571.7382	\$105.00	Distal Tip	59
	1006.1990	\$105.00		59A

\*All tips also work with SONICflex 2000L and 2000N.

Note: Tips ending in suffix 'A' are compatible with SONICflex 2008 / SONICflex 2008L only.



## Minimally invasive caries therapy (cont.)

## SONICflex Seal



Preparation of fissures. The perfect complement to the KaVo DIAGNOdent.

## SONICflex Seal Tip



45

	Part No.	Price	Description	Tip No.
	1000.8323	\$95.00	Three tips for supra- and subgingival calculus removal.	45
	1007.1506	\$95.00		45A

## SONICflex Cariex D + TC



Minimally invasive removal of carious dentine. Cariex TC: Removal of large carious lesions. Cariex D: for small carious lesions with the ability to shape the cavity


## SONICflex Cari-ex D Tip Set A

	Part No.	Price	Description	Tip No.
	1006.2025	\$235.00	SONICflex cari-ex D Tip Set A	42A, 43A

## SONICflex Cari-ex TC Tip Set A

	1006.2032	\$235.00	SONICflex cari-ex TC Tip Set A	71A, 72A
--	-----------	----------	--------------------------------	----------

## SONICflex Cari-ex D + TC tips




	1000.7105	\$154.00	Ball shaped for excavation 0.8mm	42
	1006.1980	\$154.00		42A
	1000.7167	\$154.00	Ball shaped for excavation 1.2mm	43
	1006.1981	\$154.00		43A
	1000.7362	\$154.00	Ball shaped for excavation 1.0mm	71
	1006.2002	\$154.00		71A
	1000.7363	\$154.00	Ball shaped for excavation 1.4mm	72
	1006.2004	\$154.00		72A

\*All tips also work with SONICflex 2000L and 2000N.

Note: Tips ending in suffix 'A' are compatible with SONICflex 2008 / SONICflex 2008L only.

## Inlay / onlay preparation

## SONICflex prep gold



Preparation tips diamond-coated on one side for finishing the edges of cavities, specially developed for gold inlays and partial crowns. The attachment geometry for gold inlays and partial.

## SONICflex prep gold Tip Set A




Part No.	Price	Description	Tip No.
1006.2028	\$380.00	SONICflex Prep Gold Tip Set A	49A, 50A

## SONICflex prep gold tips



0571.7212	\$193.00	Mesial Tip	49
1006.1983	\$193.00		49A
0571.7222	\$193.00	Distal Tip	50
1006.1984	\$193.00		50A

## SONICflex prep ceram



For ceramic inlays/onlays. For finishing the margins of cavities; directly transfers the geometry of the tips to the tooth substance. No more undercuts with laboratory-fabricated inlays and onlays.

## SONICflex prep ceram Tip Set A




Part No.	Price	Description	Tip No.
1006.2029	\$380.00	SONICflex Prep Ceram Tip Set A	51A 52A

## SONICflex prep ceram tips



0571.7262	\$193.00	Mesial Tip	51
1006.1985	\$193.00	Mesial Tip	51A
0571.7272	\$193.00	Distal Tip	52
1006.1986	\$193.00		52A

## SONICflex prep CAD/CAM



Ideal shaping of the cavity margins for CAD/CAM systems such as KaVo EVEREST and CEREC.

## SONICflex prep CAD-CAM insert-kit

Part No.	Price	Description	Tip No.
1002.1988	\$390.00	SONICflex CAD/CAM Tip Set	34, 35

## SONICflex prep CAD-CAM Tip Set A

1006.2024	\$390.00	SONICflex CAD/CAM Tip Set A (incl. tips 34A, 35A)	34A, 35A
-----------	----------	---	----------

## SONICflex retro tips



1002.1984	\$209.00	Mesial Tip	34
1006.1977	\$209.00		34A
1002.1986	\$209.00	Distal Tip	35
1006.1979	\$209.00		35A

## SONICflex cem



Gentle and precise insertion of inlays. Reduction of the marginal gap.

## SONICflex cem tips



Part No.	Price	Description	Tip No.
0571.5431	\$133.00	Cem tip with white cap	12
1006.2021	\$133.00		12A
0571.7142	\$15.00	Cem Inlay top	

\*All tips also work with SONICflex 2000L and 2000N.

Note: Tips ending in suffix 'A' are compatible with SONICflex 2008 / SONICflex 2008L only.



## Surgery

## SONICflex Bone



Can be used in many areas of bone treatment (osteotomies, augmentation, extractions and sinus lift).

		Part No.	Price	Description	Tip No.
SONICflex Bone Set		1004.3874	\$1,094.00	SONICflex Bone Set	80, 81, 83, 84, 85, 86, 87
SONICflex Bone Set A		1006.2033	\$1,094.00	SONICflex Bone Set A	80A, 81A, 83A, 84A, 85A, 86A, 87A
SONICflex Bone tips		1004.3875	\$187.00	Square cut for bone-splitting, distraction	80
	80	1006.2006	\$170.00		80A
	81	1004.3876	\$170.00	Large ball head for sinus lift	81
	81	1006.2008	\$170.00		81A
	82	1004.3877	\$187.00	Sagittal Saw	82
	82	1006.2010	\$170.00		82A
	83	1004.3878	\$187.00	Axial Saw	83
	83	1006.2012	\$170.00		83A
	84	1004.3879	\$170.00	Plate-shaped, for releasing Schneider's membrane in a sinus lift	84
	84	1006.2014	\$170.00		84A
	85	1006.0645	\$170.00	Plate-shaped, for releasing Schneider's membrane in a sinus lift	85
	85	1007.1624	\$170.00		85A
	86	1006.0741	\$170.00	Tip for obtaining greater amount of bone sinus lift	86
	86	1007.1625	\$170.00		86A
	87	1006.1275	\$170.00	SONICflex bone tip 87 plus 5 saw blades	87
	87	1007.1626	\$170.00	SONICflex bone tip 87A plus 5 saw blades	87A

\*All tips also work with SONICflex 2000L and 2000N.

Note: Tips ending in suffix 'A' are compatible with SONICflex 2008 / SONICflex 2008L only.

## SONICflex Accessories

## SONICflex Angle



Angled tips with diamond coated working surfaces (the ends are uncoated) for micro preparations (eg. tunnel preps) in difficult to access areas. Marginal ridge can be left intact.

## SONICflex Angle Kit



Part No.	Price	Description	Tip No.
0571.0341	\$240.00	mesial, distal	53, 54

## Accessories for SONICflex 2000 only

Part No.	Price	Tip Number	Description
1000.4887	\$68.00	Tip tool	For removing tips
0571.5641	\$201.00	No. 26	Root planer
0571.7382	\$105.00	No. 59	Distal tip
0571.7392	\$105.00	No. 58	Mesial tip



# PROPHYflex System and Accessories

## PROPHYflex 3 Black



WWW.TRY  
KAVO.COM

## PROPHYflex 3 Violet



NEW

## PROPHYpearls



NEW

## PROPHYflex Cleaning Powder



## Tips for PROPHYflex 3



Focused spray



Wide spray



## Silicone Reservoir Cover



1000.2678 \$16.00 1 cover, for storing powder canisters

## Powder Screw



0573.0412 \$36.00 For 100% powder flow Ceramic, Blue or Silver

0573.0002 \$35.00 For 70% powder flow Plastic, White or Bronze

## Cleaning Tools



0573.0321 \$25.00 Long Cleaning Brush

0573.6052 \$6.50 Nozzle Needle (small) for cleaning tip

## Parts / Accessories



0573.6152 \$36.00 Plastic powder reservoir, sterilizable

0573.6072 \$9.25 Black seal for reservoir

0200.6019 \$5.50 O-ring for tip (2 required)

Air cleaning system for removal of stains and bacterial plaque, cleaning of enamel surfaces prior to fissure sealing. The PROPHYflex 3 features a shorter, more ergonomic handle; a lighter plastic reservoir and a new improved design. MULTIflex connection, sterilizable up to 275°F (135°C), 44 psi. Water 7.2 - 10 psi.

### PROPHYflex 3 box includes:

- PROPHYflex 3
- 2 Plastic reservoirs
- 5 Packets of powder / PROPHYpearls
- 1 Large cleaning drill
- 2 Small nozzle needles
- 1 Restrictor screw
- 2 Tips (one small, one large)

Part No.	Price	Description
1000.4672	\$1,648.00	PROPHYflex 3 Black
1006.9926	\$1,648.00	PROPHYflex 3 Violet

Calcium carbonate powder eliminates salty taste. Smooth, soft & gentle to the tooth. Allows for efficient and faster cleaning. No clogging in the unit. Can be used with any manufacturers' air polishing system.

Part No.	Price	Description
1010.1798	\$115.00	PROPHYpearls Neutral (4 bottles x 8.8 oz, 250 g)
1010.1826	\$148.00	PROPHYpearls Neutral (80 sticks x 0.5 oz, 15 g)
1010.1828	\$160.00	PROPHYpearls mint (80 sticks x 0.5 oz, 15 g)
1010.1829	\$160.00	PROPHYpearls peach (80 sticks x 0.5 oz, 15 g)
1010.1830	\$160.00	PROPHYpearls orange (80 sticks x 0.5 oz, 15 g)
1010.1831	\$160.00	PROPHYpearls berry (80 sticks x 0.5 oz, 15 g)

Part No.	Price	Description
1007.0014	\$185.00	Prophy orange powder. 1 box of 80 packs
1007.0015	\$185.00	Prophy berry powder. 1 box of 80 packs
1007.0016	\$185.00	Prophy cherry powder. 1 box of 80 packs
1007.0017	\$185.00	Prophy mint powder. 1 box of 80 packs
1009.3732	\$98.00	PROPHYflex Perio Powder (4 bottles x 3.5 oz, 10g)

Part No.	Price	Description
0573.0151	\$198.00	Focused spray
0573.0181	\$198.00	Wide spray
1000.2678	\$16.00	1 cover, for storing powder canisters
0573.0412	\$36.00	For 100% powder flow Ceramic, Blue or Silver
0573.0002	\$35.00	For 70% powder flow Plastic, White or Bronze
0573.0321	\$25.00	Long Cleaning Brush
0573.6052	\$6.50	Nozzle Needle (small) for cleaning tip
0573.6152	\$36.00	Plastic powder reservoir, sterilizable
0573.6072	\$9.25	Black seal for reservoir
0200.6019	\$5.50	O-ring for tip (2 required)



### PROPHYwiz 181P Hygiene Handpiece

- Ultra light handpiece, minimizes hand fatigue
- Ergonomic and well-balanced
- Whisper quiet, smoothest running
- Autoclavable to 275°F (135°C)
- 3,100 rpm
- For metal Doriot attachments or disposable prophy angles
- 360° rotation at center
- 36-43 PSI
- 1-year warranty
- 2.7 oz

Part No. 1003.2278

\$819.00

Part No. 1003.9882

\$1,715.00 PROPHYwiz 3 pack

### 19E / 19ES Hygiene Handpiece

- Reduction 7.4:1
- Speed range up to 5,000 rpm
- For disposable prophy heads only
- Use with 181K Air Motor (19E)
- Use with universal E-connection (19ES)

Part No. 1004.6403

\$499.00 19E

Part No. 1004.4503

\$499.00 19ES

### 181K Air Motor

- Speed range up to 20,000 rpm
- 4-hole connection
- Non-optic
- Forward/reverse speed adjustment
- Sterilizable
- 36-43 PSI
- 1-year warranty
- 4.8 oz

Part No. 0535.5260

\$755.00

### 31E Low-Speed Handpiece

- Reduction 4:1
- Speed range up to 5,000 rpm
- For use with snap-on / screw-in cups and brushes
- Extremely small head with 100° for superior access
- Durable hygienic coating with textured surface for easy grip
- Use with 181K Air Motor
- Comes with Prophy cup holder

Part No. 1003.6976

\$599.00



## Caries Detection | DIAGNOdent pen and Display



Shown: DIAGNOdent pen and Display Unit 2191.  
Units sold separately.



### DIAGNOdent pen

High Tech Laser Fluorescence Caries Detection Aid which measures changes in tooth structure caused by caries. It is especially useful for pit and fissure areas - even when the outer tooth surface seems to be intact. 90% accurate in detecting incipient caries. Only caries detection device backed by numerous clinical studies. Allows monitoring over time. Easy to use. Portable, handheld, pen-like device. Moves easily from room to room. Ideal for offices that already use the DIAGNOdent Classic and would like to have the added mobility of the hand-held DIAGNOdent pen. Weight 5 oz Length: 8-3/4". Display sold separately.

#### Includes:

- DIAGNOdent pen unit
- 2 Sterilizable tips for occlusal surfaces
- Sterilization cassette
- Sterilizable calibration standard
- 1 standard AA battery
- Clinical guidelines and operation manual
- Laminated quick reference guide

Part No. 1004.3400 \$2,689.00



### Display Unit 2191

With the new DIAGNOdent pen Display, we have enhanced the patient experience in caries diagnosis. The remote display for DIAGNOdent pen uses infrared signals, and the values are shown simultaneously on the DIAGNOdent pen and on the remote display. The display also serves as a holder for the DIAGNOdent pen when not in use. Pen sold separately.

Part No. 1004.8400 \$750.00



### DIAGNOdent pen Perio Tip

DIAGNOdent now goes beyond caries detection to detect subgingival calculus for the treatment of periodontal disease. Sterilizable.

Perio Tip for detection of subgingival calculus in periodontal pockets.

Part No. 1004.0370 \$392.00

### DIAGNOdent pen Accessories



Sterilization and storage cassette



Diagnodent pen Barriers (250 Barriers)

Part No.	Price	Description
1002.6967	\$156.00	Tip for occlusal surfaces.
1002.7003	\$77.00	Sterilizable grip sleeve
1002.7020	\$110.00	Removable "C" standard calibration disk
1002.7011	\$125.00	Sterilization storage cassette without calibration standard
1007.1480	\$30.00	DIAGNOdent pen Barriers (250 pack)
0906.8088	\$22.00	Patient Education Flyers (100 pack)

**DIAGNOdent Classic**

High Tech Laser Fluorescence Caries Detection Aid which measures changes in tooth structure caused by caries. It is especially useful for pit and fissure area – even when the outer tooth surface seems to be intact. 90% accurate in detecting incipient caries. Only caries detection device backed by numerous clinical studies. Allows monitoring over time. Easy to use.

## Includes:

- DIAGNOdent Classic unit
- Sterilization Cassette
- Handpiece with Tubing
- 2 sterilizable tips - type A for occlusal surfaces
- Rechargeable Battery
- Battery Charger
- Sterilizable calibration standard
- Clinical guidelines and operation manual
- Laminated quick reference guide

Part No. 0574.0500

\$2,689.00

**DIAGNOdent Classic Perio Tip**

DIAGNOdent now goes beyond caries detection to detect subgingival calculus for the treatment of periodontal disease. Sterilizable.

Perio tip for detection of subgingival calculus in periodontal pockets.

Part No. 1004.1640

\$392.00

**DIAGNOdent Classic Accessories**

A Tip



B Tip



Snap in Ceramic standard



Magnetized hand-piece holder



DIAGNOdent Classic Barriers (250 Barriers)



A Tip kit



Sterilization cassette without standard

Part No.	Price	Description
0574.1311	\$161.00	Type A Tip for occlusal surfaces
0574.1321	\$161.00	Type B Tip for flat surfaces
1002.8569	\$410.00	Type A Tip Kit Includes: 4 Type A Tips, Sterilization cassette without standard
0574.1271	\$52.00	Magnetized handpiece holder
1003.5528	\$65.00	Snap in Ceramic Standard "C" (for DIAGNOdents with serial numbers higher than 1500 000)
0574.1692	\$36.00	Snap in Ceramic Standard "B" (for DIAGNOdents with serial numbers smaller than 1500 000)
1000.1533	\$129.00	Sterilization cassette w/o standard
0574.1341	\$65.00	Mounting bracket (to mount DIAGNOdent to delivery system)
1007.1479	\$30.00	DIAGNOdent Classic Barriers (250 pack)
1007.5158	\$85.00	Rechargeable Battery
1007.5160	\$85.00	Battery Charger
1007.7209	\$125.00	Rechargeable battery plus charger (Replaces part # 0574.1261)
0574.1211	\$93.00	DIAGNOdent Classic Sheath
0906.8088	\$22.00	Patient Education Flyers (100 pack)



# RONDOflex plus

## RONDOflex plus 360 with water spray



**REALITY**  
Four Star Award  
★★★



### Tips for RONDOflex plus 360



### Tips for RONDOflex 2013\*



The RONDOflex plus is a portable, easy-to-use abrasion handpiece that adds water spray to air abrasion technology. The RONDOflex plus is an efficient, "snap-on" alternative to your existing air abrasion unit or small investment to get you started. Ideal adjunct to the DIAGNOdent Caries Detection Aid. MULTIflex connection, sterilizable up to 275°F (135°C), 46-55 psi.

#### Indications

- Preparation of small carious lesions
- Opening and extension of fissures

#### Surface roughening of preparations and restorations

- Removes residual cement from crowns and bridges

#### Includes:

- RONDOflex plus
- 2 Tips (Both 110° tips)
- 1 Reservoir 27 micron
- 1 Reservoir 50 micron
- 1 - 27 micron powder (75g)
- 1 - 50 micron powder (75g)
- Small nozzle needle
- Large cleaning drill
- 2 Silicone reservoir covers

Part No. 1002.2179      \$2,850.00

Part No.	Price	Tip Angulation	Size
1002.9176	\$417.00	110°	0.018" small
1002.6251	\$417.00	110°	0.025" wide
1002.9182	\$417.00	90°	0.018" small
1002.9179	\$417.00	90°	0.025" wide
1001.3254	\$245.00	120°	0.018" small
1000.5698	\$249.00	120°	0.025" wide
1001.3253	\$245.00	90°	0.018" small
1000.5697	\$249.00	90°	0.025" wide

\* Discontinued model without water spray.

	Part No.	Price	Description
Reservoir for RONDOflex	1003.1235	\$33.00	27 micron reservoir
		\$46.00	50 micron reservoir
Reservoir Seal	0573.6072	\$9.25	Black seal for reservoir
			
27 Micron aluminum oxide powder 1000g (2.2 lbs.). Recommended for intraoral use.	1000.5957	\$326.00	6 pack
			
50 Micron aluminum oxide powder 1000g (2.2 lbs.). Recommended for extraoral use.	1000.5956	\$280.00	6 pack
			
Silicone reservoir cover	1000.2678	\$16.00	1 cover. For storing powder canisters.
			
Powder screw	0573.0412	\$36.00	For normal flow powder (100%) Ceramic
			
Cleaning Tools	0573.0321	\$25.00	Long cleaning brush, (like for the prophyflex)
			
	0573.6052	\$6.50	Nozzle needle (small) for cleaning tip
			
Tip Tools	1001.1206	\$19.00	For changing tips on ROND
			



## CORONAflex 2005 Crown and Bridge Remover



- Pneumatically operated instrument to remove crowns and bridges
- In most cases the restoration can be removed undamaged
- Excellent tool if patient needs recementation of a loose bridge anchor, endodontic treatment, etc
- Comes with detailed instructions; easy to use
- MULTIflex connection
- Operates between 35-60 psi
- 1 bar = 14.50 psi

## Includes:

- CORONAflex 2005
- Crown forceps
- 2 small clamps
- 2 large clamps
- Loop holder
- 2 Loops

Part No. 0579.1000 \$2,050.00

## CORONAflex Tools




Part No.	Price	Description
0579.0111	\$942.00	Crown forceps

## Clamps



0579.0392	\$29.00	Clamp Size: Small
0579.0402	\$29.00	Large

## Loops holder and Loops



0579.0412	\$140.00	Loop holder
0579.0422	\$33.00	Loop





# TryKaVo.com

## Try the Best Handpieces, Risk Free



Step 1.

### CHOOSE



Choose a handpiece from  
our collections listed

Step 2.

### TRY



Try a handpiece risk free  
for 5 business days

Step 3.

### RETURN



After 5 day trial, return the  
handpiece. Like it, buy it.  
Contact 1-877-746-3995



Go to **TRYKAVO.COM** to start  
your free and easy trial today!



# KaVo Surgical Systems

Driven to Precision

The KaVo INTRAsurg 300 and INTRAsurg 300 plus, are precise implant and oral surgery systems. They are simple, intuitive and easy to use. The units can be used for all surgical applications and offer intelligent power and precise torque. The KaVo Innovation story continues with INTRAsurg 300 plus.



**INTRAsurg 300 plus**

Handpiece not included with INTRAsurg 300/300 plus



Unique surgical unit with light. Compact stand-alone unit for implantology, oral and maxillofacial and periodontal surgery. Ideal for beginners in implantology – works with any implant system! Easy to customize to doctor's needs.

**System Advantages:**

- Displays actual speed and torque.
- 25,000 LUX light output
- Automatically recognizes optic surgical attachment and gear ratio
- Memory chip to remember max. torque used during each procedure step
- 1 program with max. of 6 steps – in each step user can program desired parameters: gear ratio of attachment, max. speed, max. torque, amount of coolant flow, motor fwd/reverse

**Technical:**

- Voltage: 100V to 240V
- Brushless Optic Motor: INTRA-SL 550
- Motor Speed Range: 300 – 40,000 rpm
- Max Torque: 300 to 12,000 rpm – 5.5 Ncm; 12,000 – 40,000 rpm – 5.5 – 2.0 Ncm
- Pump Delivery: 20 – 110 ml/min.
- Height: 16-1/16" with water bottle
- Height: 5" without water bottle
- Width: 15-1/2" with water bottle
- Depth: 13" (allows plug / cables)
- Weight: 15.55 lbs.

Part No. 1001.3439      \$8,999.00

**INTRAsurg 300**

Handpiece not included with INTRAsurg 300/300 plus

Compact stand-alone unit for implantology, oral and maxillofacial and periodontal surgery. Ideal for beginners in implantology – works with any implant system! Easy to customize to doctor's needs.

**System Advantages:**

- 1 program with max. of 6 steps – in each step user can program desired parameters: gear ratio of attachment, max. speed, max. torque, amount of coolant flow, motor fwd/reverse
- Displays actual speed and torque
- Autoclavable brushless motor INTRA-S 550 with sterilizable tubing
- Easy to operate – only five buttons
- Multifunctional foot control can be used to control unit hands-free
- Integrated irrigant pump with automatic rinse cycle for saline removal

**Technical:**

- Voltage: 100V to 240V
- Brushless Motor: INTRA-S 550
- Motor Speed Range: 300 – 40,000 rpm
- Max Torque: 300 to 12,000 rpm – 5.5 Ncm; 12,000 – 40,000 rpm – 5.5 – 2.0 Ncm
- Pump Delivery: 20 – 110 ml/min.
- Height: 16-1/16" with water bottle
- Height: 5" without water bottle
- Width: 15-1/2" with water bottle
- Depth: 13" (allows plug / cables)
- Weight: 15.55 lbs

Part No. 1001.3440      \$6,499.00





## CL 09 Cellular Optic Contra Angle Attachment

- Transmission 9:1
- Speed range up to 4,444 rpm
- Ideal to use with INTRA CL3 head (3:1) (speed range 11 to 1,500 rpm)
- Extremely strong reduction gears
- For use with INTRAsurg 300 plus only
- 1 green ring

Part No. 1001.9258 \$1,442.00

## C 09 Contra Angle Attachment

- Transmission 9:1
- Speed range up to 4,444 rpm
- Ideal to use with INTRA C3 head (3:1) (speed range 11 to 1,500 rpm)
- Easily dismantled for cleaning
- Extremely strong reduction gears
- 1 green ring

Part No. 1001.2278 \$888.00

## Cellular Optic Attachments for Oral Surgery and Implantology

### CL 10 Straight Attachment

- 25,000 LUX glass rod
- For handpiece burs with standard length of 44.5mm
- With external spray attachment
- Strengthened drill chuck
- Speed range up to 40,000 rpm
- For use with INTRAsurg 300 plus only
- 1 green ring

Part No. 1002.0911 \$2,199.00

## Cellular Optic Attachments for Oral Surgery and Implantology

### CL 3-09 Contra Angle Attachment

- 25,000 LUX glass rod
- Head and shank can be dismantled
- Reduction 27:1
- For latch burs
- Speed range 11 to 1500 rpm
- Small head with internal cooling system
- External cooling medium connection attached to shank
- For use with INTRAsurg 300 plus only

Part No. 1002.3693 \$2,472.00

## CL 3 Reduction Latch Head

- For use with CL09 cellular optic Contra Angle
- For internally cooled burs
- Reduction 3:1
- Internal/external coolant
- Small head diameter
- Max torque 43 Ncm
- For use with INTRAsurg 300 plus only

Part No. 1001.8993 \$970.00  
Part No. 0540.3711 \$56.00 Clip for Water filter

## INTRAsurg 300 / 300 plus Accessories and Replacement Parts

Part No.	Price	Description
1001.4968	\$5,125.00	Steri-set fiber optic surgery motor INTRA-SL 550 Includes: INTRA-SL 550 fiber optic motor and internal irrigant tubing, motor support
1001.3421	\$3,986.00	INTRA-SL 550 fiber optic surgery motor only
1002.2928	\$98.00	INTRA-SL 550 bulb. Autoclavable bulb for INTRA-SL 550 fiber optic surgical motor
1001.4953	\$3,374.00	Steri-set non-optic surgery motor INTRA-S 550. Includes: INTRA-S 550 non-optic motor and internal irrigant tubing, motor support
1000.8072	\$2,008.00	INTRA-S 550 non-optic surgery motor nly
1001.3347	\$123.00	Sterilizable pump tubing (Bag to handpiece tubing)
1001.9902	\$203.00	Sterilizable tube fittings. 10 each Single use
1001.2651	\$1,523.00	Autoclavable tubing (Complete set)



## INTRAsapt 905 Accessories

Part No.	Price	Description
0641.7301	\$298.00	Water flush assembly
0516.7870	\$1,287.00	Micromotor tubing. Light beige.
0516.3011	\$111.00	Connected pump tubing. Grey.
0065.5389	\$16.00	Inner tube - 2x1. Sold by meter.
1001.3436	\$2.00	Handpiece tubing. Sterilizable. Clear.
0516.5512	\$19.00	Intermediate piece
1002.3569	\$151.00	INTRAsurg calibration



## 3620N Contra Angle Attachment

- Transmission 1:1
- Speed range up to 40,000 rpm
- Easily dismantled for cleaning
- 1 blue ring

Part No. 0534.5650 \$718.00

## 3624N Contra Angle Attachment

- Reduction 4:1
- Speed range up to 10,000 rpm
- Easily dismantled for cleaning
- 1 green ring

Part No. 0534.5660 \$925.00

## 3630N Contra Angle Attachment

- Reduction 30:1
- Speed range up to 1,333 rpm
- Easily dismantled for cleaning
- 2 green rings

Part No. 0534.5670 \$1,129.00

## C3 Reduction Latch Head for Internally Cooled Burs

- Reduction 3:1
- Internal/external coolant
- Small head diameter
- Ideal to use with C09 contra-angle attachment (pg. 72)
- Max torque 43 Ncm
- Green dot

Part No. 1001.4503 \$723.00

# Surgical Handpieces and Accessories



## 67RIC Reduction Latch Head for Internally Cooled Burs

- Reduction 4:1
- Internal/external coolant
- Green dot
- Reinforced gearing
- Max torque 55 Ncm

Part No. 0540.8340 \$670.00



## 3610N1 Straight Attachment

- Transmission 1:1
- Speed range up to 40,000 rpm
- Easily dismantled for cleaning
- Reinforced chuck
- Maintain with closed chuck without bur
- 1 blue ring

Part No. 0524.5600 \$1,424.00



## 3610N2 Surgical Straight Attachment

- Handpiece 1:1; For handpiece burs with special length 62.5 mm, with spray attachment and with locking device for alternative use of standard handpiece burs (44.5 mm)
- Can be easily dismantled
- Speed (max): 40,000 rpm Torque (max): 10 Ncm

Part No. 0524.5610 \$972.00



## 3610N3 Surgical Straight Attachment

- Handpiece 1:1
- For handpiece burs with special length 70.0 mm with spray attachment
- Can be easily dismantled and features strengthened bur chuck
- Speed (max): 40,000 rpm Torque (max): 10 Ncm

Part No. 0524.5620 \$1,648.00

## INTRA Surgical Head Accessories

Part No.	Price	Description
0540.1971	\$57.00	Internal spray line for 67IC
0593.0361	\$66.00	Y-Distributor for 67IC/67RIC/C3/CLO3
1001.3436	\$2.00	External spray tubing for heads with external coolant. Clear sterilizable tubing.
0540.3711	\$56.00	Spray tube for 67RIC
0065.5279	\$35.00	Irrigation tubing to attach to the exterior of the electric motor tubing (INTRAsurg 500). Can be cut to desired length. Sold by meter.

## Attachment Accessories for ENT Surgery

Part No.	Price	Description
0549.0020	\$302.00	Surgical attachment for use with 8A contra angle attachment. Accepts regular length handpiece burs. Bur tool sold separately.
0411.3521	\$196.00	Bur tool for surgical attachment
0524.1531	\$55.00	Spray collar (3554, 3555, 3556)
0524.3062	\$29.00	Bur stop (3555, 3556)
0410.1963	\$23.00	J Hook for bur stop
0410.0941	\$16.00	Nozzle wire

**GENTLEray 982 Laser Whitening Handpiece Accessories**

Sterilizable handpiece expands chairside treatment options; creates a uniform parallel laser beam profile (spot diameter of approx. 5 mm) that accelerates the whitening gel without causing unnecessary heating. Available as option for both GENTLEray 980, Classic and Premium lasers.

Part No. 1005.2895      \$1,531.00

**GENTLEray 980 Accessories**

Part No.	Price	Description
1005.2894	\$317.00	Laser Safety goggles Value Line
1005.2896	\$517.00	Whitening Fiber
1006.0031	\$208.00	Straight Cannula Premium 3/pk
1006.0907	\$173.00	30° Cannula Premium 3/pk
1006.0910	\$173.00	70° Cannula Premium 3/pk
0065.5279	\$35.00	Handpiece tubing per meter
1006.0500	\$6.50	Instrument Tubing
1004.3926	\$75.00	Extra Keys
0236.2394	\$109.00	Interlock
1006.8257	\$610.00	981 Premium water cooling handpiece set
1006.3976	\$111.00	Locking mechanism for foot pedal

**GENTLEray 980 Accessories (cont.)**

Part No.	Price	Description
1006.4428	\$458.00	300 micron laser fiber, 3M (red)
1006.5431	\$458.00	200 micron laser fiber, 3M (blue)
1005.2887	\$49.00	Straight Cannula Classic 3/pk
1005.2888	\$49.00	30° Bended Cannula Classic 3/pk
1005.2889	\$49.00	70° Bended Cannula Classic 3/pk
1004.9800	\$508.00	981 Classic Laser Handpiece Set
1005.3955	\$379.00	Laser safety glasses XC with insert for glasses
1005.3609	\$75.00	Ceramic scissors
0236.2439	\$6.50	Laser warning sign
1005.2885	\$1,701.00	Foot switch



## In-Office Laboratory

### K-ERGOgrip Installation



Part No.	Price	Description
1005.4207	\$2,300.00	K-ERGOgrip/K-Control TLC-Install. EWL 4956/4944 foot control
1005.4206	\$2,300.00	K-ERGOgrip/K-Control TLC-Install Knee EWL 4955/4944
1005.4204	\$2,300.00	K-ERGOgrip/K-Control TLC-Install. EWL 4957/4944 Bench

### K-POWERgrip Installation



1007.2082	\$2,300.00	K-POWERgrip/K-Control TLC Knee Control Installation 4956/4941
1007.2084	\$2,300.00	K-POWERgrip/K-Control TLC Foot Control Installation 4955/4941
1007.2083	\$2,300.00	K-POWERgrip/K-Control TLC On-bench Control Installation 4957/4941

### K5 Plus Installation



1009.4244	\$1,535.00	K5 Plus/K-Control TLC Knee Control Installation 4956/4911
1009.4245	\$1,535.00	K5 Plus/K-Control TLC Foot Control Installation 4955/4911
1009.4243	\$1,535.00	K5 Plus/K-Control TLC On-bench Control Installation 4957/4911

### SMARTair Mobile

Portable unit with 2 filter bags.



Part No.	Price	Description
0658.2160	\$155.00	Pack of 5 filter bags
1001.8196	\$1,997.00	

Please note:

This is only a limited selection of our wide range of laboratory products.

For complete in-office laboratory set-up solutions, call KaVo at 888.ASK.KAVO for pricing.

## LCM Illuminator



- Compact bulb-in-coupler system
- For foot pedal or hanger installation
- Operates up to four handpieces
- Includes module, transformer, remote intensity control, and installation kit
- LCM tubing and MULTIflex LUX (465LRN / 465LED or 459AN, 460LE / 460LED) couplers ordered separately

Part No. 0905.8400 \$626.00

## LS2 Illuminator



- Compact remote fiber optic light source for hanger or foot pedal installation
- Operates up to three handpieces
- High/low intensity switch
- Installation kit included
- LS2 tubing and MULTIflex LUX (1390) coupler ordered separately

Part No. 0905.8500 \$920.00

## LCM Tubing—6-Pin Type Tubing

Part No.	Price	Color	Description
0905.4113	\$398.00	Gray	Straight 11' LCM tubing

## LS2 Tubing—Standard 5-Hole Type Tubing

Part No.	Price	Color	Description
0905.5793	\$485.00	Gray	Straight 7' x 9'
0905.5243	\$485.00	Gray	Straight 12' x 14'

## Parts and Accessories for LCM

Part No.	Price	Description
1006.6225	\$148.00	Transformer (120 Volt AC)
0905.8405	\$47.00	Low voltage cable (15 Volt AC)
0905.8406	\$96.00	LCM remote intensity control
0905.8440	\$43.00	LCM installation kit
1002.2928	\$98.00	LCM bulb (2 pack). (3 to 3.3 Volt DC)
1007.5372	\$121.00	MULTI LED bulb (3.0 – 3.6 V DC. Current < 120 mA)
0905.8411	\$29.00	Voltage clip
0905.8403	\$499.00	LCM module only (3.6 Volt DC)

## Parts and Accessories for LS2

Part No.	Price	Description
0905.8510	\$74.00	Halogen 12 Volts / 100 Watt bulb
0905.8540	\$67.00	Installation kit
0905.8501	\$565.00	LS2 module only
0905.8520	\$55.00	Connector cable assembly
0513.0812	\$35.00	Coupler nut (protects crimping)
0905.8550	\$59.00	Time delay switch kit





VistaClear™ 1000 Treatment System  
(1 Operatory)



VistaClear™ 1000-47-C Treatment System  
(7 Operatories)

## Pelton & Crane VistaClear™ System

VistaClear provides improved water quality and reduces bacterial contamination in dental unit waterlines used for irrigation, cooling, lubrication and scaling procedures. These patented, award-winning systems also remove disinfection residuals from municipal water supplies while reducing hard water scaling potential to help prevent plugging of dental unit ports and attached appliances.

VistaClear Systems save time and money since they completely eliminate the need for dental bottles. And, no dangerous chemicals, ozone, iodine or other potentially allergic or corrosive agents are used with patients. No maintenance or attention is required daily or between patients like other products. All that's necessary is a simple quarterly maintenance procedure and changing the filter elements about once per year.

Multiple operatory systems suitable for treating up to 7 operatories per zone are available. All VistaClear Systems are EPA listed as pest control devices.

Part No.	Description
2961	VistaClear 1000 Central Treatment System (1 Operatory)
2963	VistaClear 1000-23-C Central Treatment System (2-3 Operatories)
2964	VistaClear 1000-34-C Central Treatment System (4 Operatories)
2966	VistaClear 1000-46-C Central Treatment System (5-6 Operatories)
2967	VistaClear 1000-47-C Central Treatment System (7 Operatories)



VistaPure™

## Pelton & Crane VistaPure™

VistaPure provides a convenient and inexpensive supply of high purity, distilled-quality water for labs, dental and medical offices. The system produces two grades of water — distilled-quality water for autoclaves and ultra filtered non-corrosive water for bottles used with dental delivery systems. For complete selection of Pelton & Crane products visit [www.pelton.net](http://www.pelton.net) or call 800.659.6560.

Type	Part No.	Page No.	Retail
200XDR	1007.3570	10	\$1,699.00
25LPR	1007.3580	12	\$1,699.00
7LP	1001.7452	12	\$1,120.00
10LP	1001.5660	12	\$1,325.00
20LP	1001.7453	12	\$920.00
29LP	1001.7454	12	\$1,250.00
L68B	1008.1834	13	\$565.00
L80	1004.4456	13	\$565.00
L52	1008.1825	13	\$565.00
L62	1008.1830	13	\$585.00
L66B	1008.1831	13	\$585.00
L67	1008.1832	13	\$585.00
L53	1008.1826	13	\$550.00
L3	1008.1822	14	\$550.00
L3Y	1008.1823	14	\$550.00
L31	1008.1824	14	\$305.00
L61G	1008.1828	14	\$550.00
L61R	1008.1829	14	\$550.00
L22	1008.1835	14	\$550.00
E25L	1007.5550	19	\$1,369.00
E25C	1007.5551	19	\$1,169.00
E20L	1007.5540	19	\$1,299.00
E20C	1007.5541	19	\$919.00
E10C	1007.5560	20	\$789.00
E15L	1007.5530	20	\$1,135.00
E15C	1007.5531	20	\$919.00
E61C	1008.4999	20	\$999.00.00
E31C	1008.4998	20	\$1,075.00
19ES	1004.4503	21	\$499.00
23ES	1004.5380	21	\$963.00
80ES	1004.4375	21	\$600.00
M4500L	1007.1700	26	\$1,550.00
M8900L	1008.7800	26	\$1,550.00
E675L	1006.6400	27	\$1,120.00
E679L	1006.9200	27	\$1,120.00
E675LMW	1006.0800	27	\$1,150.00
S619L	1008.1641	28	\$715.00
S619C	1008.1642	28	\$635.00
S615L	1008.1643	28	\$715.00
S609C	1008.1644	28	\$529.00
S605C	1008.1645	28	\$529.00

Type	Part No.	Page No.	Retail
2307LN	0535.6190	30	\$1,599.00
2313LN	0535.6220	30	\$1,299.00
2320LN	0535.6180	30	\$1,499.00
2325LN	0535.6210	30	\$1,999.00
457	0553.1210	31	\$260.00
460LE	1001.7599	31	\$315.00
460LED	1007.3201	31	\$365.00
465LED	1007.3197	31	\$465.00
465LRN	0553.1550	31	\$405.00
1390	0553.1390	31	\$354.00
10E	1003.6977	32	\$529.00
19E	1004.6403	32	\$499.00
20E	1003.1107	32	\$599.00
31E	1003.6976	32	\$599.00
14E	1003.1108	32	\$545.00
80E	1004.6401	32	\$599.00
181L	0535.5210	33	\$1,899.00
181H	0535.5200	33	\$1,649.00
181K	0535.5260	33	\$755.00
QUATTROcare Plus	1008.3805	48	\$3,198.00
2003L	1000.8333	53	\$1,999.00
2003	1000.8334	53	\$1,830.00
2008	1007.1606	53	\$1,830.00
2008L	1007.1604	53	\$1,999.00
181P	1003.2278	63	\$819.00
CORONAFlex	0579.1000	68	\$2,050.00
CL3-09	1002.3693	72	\$2,472.00
CL3	1001.8993	72	\$970.00
C09	1001.2278	72	\$888.00
CL09	1001.9258	72	\$1,442.00
CL10	1002.0911	72	\$2,199.00
3620N	0534.5650	73	\$718.00
3624N	0534.5660	73	\$925.00
3630N	0534.5670	73	\$1,129.00
C3	1001.4503	73	\$723.00
67RIC	0540.8340	74	\$670.00
3610N1	0524.5600	74	\$1,424.00
3610N2	0524.5610	74	\$972.00
3610N3	0524.5620	74	\$1,648



## Maintenance Information



Sample maintenance page

For questions contact customer service at  
1-888-ASK KAVO or visit  
[www.kavousa.com](http://www.kavousa.com)

Maintenance Brochure Chart	Brochure No.	Description
	906.8045	KaVo E-type low-speed attachments with removable heads
	906.8049	Surgical handpieces
	906.8047	Intra air motors 181L, 181H, 181K, PROPHYwiz 181P
	906.8045	E-type low-speed attachments
	906.8048	SONICflex 2008, 2000, 2000L/N and 2003L/N scaler
	906.8043	INTRAfex LUX 2 handpieces
	906.8046	INTRA heads
	906.8044	Speed increasing attachment INTRAmatic LUX 25LCA / 25LCR / LUX 25LHA / GENTLEpower LUX 25LPA / 25LPR and e-type low-speed attachments with fixed heads

KaVo Repair Fax: Do you have a filled out repair form complete with a purchase order number from your dealer.  
Fax the completed form to **704-927-0818**.

## Terms of Sale

### Pricing:

Prices are subject to change without notice. Possession of list price does not constitute an offer to sell.

### Terms:

Payments are due 30 days from date of invoice. A 1% per month service charge (12% annually) will be applied to accounts past due. KaVo reserves the right to change this percentage based on inflation rates.

### Minimum Order:

Orders totaling less than \$50.00 wholesale cannot be economically processed and will be subject to a \$7.00 service charge. Repairs are exempt from the minimum order requirement.

### Shipment:

All shipments are FOB Factory unless otherwise stated and will be sent the most economical way unless otherwise specified. Shortage claims must be made within 10 days from receipt of goods. Damaged goods must be reported to the carrier within 24 hours and damage claims must be filed by the consignee.

### Drop Shipment Fees:

All shipments of product from our warehouse that bypass the normal distribution channel will be subject to a drop shipment fee and freight charges.

### Returns:

Prior approval must be obtained before merchandise is returned for credit or exchange. Merchandise that is less than nine (9) months from date of sale and still manufactured will be accepted for credit. Shipment of returned merchandise is customer's responsibility. Merchandise out of warranty may be returned for repair, not for credit. KaVo will not accept returns on products with retail price less than \$25.00.

Merchandise will be credited at the lower of market or purchase price. Merchandise returned in the original, unopened container and accompanied by the invoice number and date of purchase will be subject to a minimum 15% (\$5.00 minimum) restocking charge in addition to any charges for reconditioning or repackaging. Used products are not acceptable for return.

### Product Changes:

KaVo Dental Corporation reserves the right to make changes in product design without the obligation to accept previously manufactured products for alteration, exchange or credit.

### Warranty:

KaVo Dental Corporation warrants its products to the original purchaser against defective workmanship and materials under normal procedures of installation, use and service within the dental profession for a period of 90 days from the date of sale. Selected KaVo Dental Corporation products contain further written warranty terms. Warranty period begins no later than dealer invoice date.

KaVo dealers will sell products within a geographic location that the dealer can adequately service by making calls on current and potential customers.

### Training Policy:

Dealer shall employ and maintain an adequate, knowledgeable, full time sales staff who regularly contact users and potential users within the territory. Dealers will utilize service and installation staff who are technically trained in the installation of products within the territory.

Dealer's staff shall possess and maintain adequate knowledge concerning the products and their uses, features and benefits. Dealer will cooperate with KaVo in conducting meetings for training of dealer personnel for the above mentioned reason.

### Hazmat Fees:

Hazardous transport fees by carriers will be passed along to customer. Contact Customer Service for economic shipping solutions on all hazardous goods.

## Registered Trademarks

The following product names are registered trademarks used by KaVo group and are therefore protected by the trademark license:

BELLAtorque	GENTLEforce	INTRAsurg	PROPHYflex
BELLAtorque LUX	GENTLEpower	KAVO Care	PROPHYpearls
COMFORTbase	GENTLEpulse	KAVO spray	PROPHYwiz
COMFORTdrive	GENTLEray	KaVo SPRAYrotor	QUATTROcare
COMFORTtronic	GENTLEmini	LUX 2	RONDOflex
Companion	GENTLESilence	LUX 3	ROTAspray
CORONAFlex	INTRA	MASTERtorque	SAFEdrive
DiaLUX	INTRAccompact	MINItech	SMARTair
DIAGNOdent	INTRAflex	Mira	SMARTtorque
ELECTROtorque	INTRAflex LUX	MULTIflex	SONICflex
EXPERTgrip	INTRALUX	MULTIflex LUX	SONICflex LUX
EXPERTmatic	INTRAmatic	MULTI LED	SONICSYS
EXPERTtorque	INTRAmatic LUX	Opti-Torque	
EVEREST	INTRAsapt	PHYSIO	



# Warranty

## Warranty

For the periods of time specified below, KaVo guarantees that the following products, when furnished by KaVo in new and unused condition, will not fail during normal service due to defects in material and workmanship that are not due to apparent abuse, misuse or accident:

Cellular Optic Glass Rod in all High-Speed, Low-Speed & SONICflex	5 years	MULTIflex Couplers	1 year
Smart Series (646B, S619L, S619C, S615L, S609C, S605C)	1 year	Expert Series (E679L, E675LMW, E675L)	2.5 years*
COMFORTdrive 200 XDR	2.5 years*	Master Series (M8900 L and M4500 L)	2.5 years*
COMFORTtronic and COMFORTbase	3 years	PROPHYflex 3	1 year
DIAGNOdent / DIAGNOdent pen (except tips, tubing and display)	1 year	PROPHYwiz	1 year
ELECTROtorque TLC / ELECTROtorque plus (control box, tubing, motor)	3 year	QUATTROcare / QUATTROcare Plus	1 year
GENTLEray 980 (Classic or Premium)	2 years	Repairs on Low-Speed Heads, Attachments and Speed Increasing Attachments (if repaired at an authorized KaVo repair center)	1 year
INTRAfex LUX	1 year	RONDOflex plus	1 year
INTRAsurg and all surgical handpieces	1 year	SONICflex / SONICflex LUX	1 year
KaVo LCM, KaVo/Vicon Light Sources	1 year	Tubings (all types)	1 year
Low-Speed Heads and Attachments, Speed Increasing Attachments	2.5 years*	Turbines (if handpiece repaired at an authorized KaVo repair center or repaired in-office)	1 year
Low-Speed Motors	1 year		

Additional Warranty terms apply. See product warranty card.

\* 2.5 year warranty if maintained with KaVo QUATTROcare Plus, if not maintained with QUATTROcare Plus the standard warranty of 2 years applies.

## Dental Advisor Awards



COMFORTdrive  
QUATTROcare Plus



635B Handpiece  
CORONAFlex  
ELECTROtorque plus  
PROPHYwiz  
QUATTROcare  
GENTLEsilence 6500BR



## Townie Awards



Air Driven High-Speed Handpieces - winner 2013, 2012, 2011, 2010, 2009, 2008, 2007 and 2006

ELECTROtorque TLC - winner 2011, 2010, 2009, 2008, 2007 and 2006

ELECTROtorque plus - winner 2011, 2010, 2009, 2008, 2007, 2006, 2005, 2004 and 2003

DIAGNOdent pen - winner 2013, 2012, 2011, 2010, 2009, 2008, 2007 and 2006

DIAGNOdent - winner 2013, 2012, 2011, 2010, 2009, 2008, 2007, 2006, 2005, 2004 and 2003

Cellular Optic High-Speed Handpieces - winner 2010, 2005, 2004 and 2003

Air Driven Low-Speed Systems - winner 2005 and 2004

Electric Driven Low-Speed Handpieces - winner 2013, 2012 and 2010

QUATTROcare - winner 2005, 2004 and 2003

PROPHYwiz - winner 2004

SONICflex LUX - winner 2004

GENTLEpower 25 LPA / 25 LPR - winner 2013, 2012, 2011, 2010, 2008 and 2009

## Reality Awards



GENTLEsilence 6500B  
PROPHYflex 3  
DIAGNOdent  
K1 System  
SONICflex LUX 2003L  
CORONAFlex  
ELECTROtorque plus  
GENTLEpower 25LPA  
INTRAFlex LUX  
635B Handpiece  
QUATTROcare  
COMFORTdrive



PROPHYwiz 181P  
DIAGNOdent pen  
RONDOflex plus



DIAGNOdent pen



ELECTROtorque plus  
GENTLEpower 25LPA

## Dentistry Today Award



DIAGNOdent pen

## Dental Products Report



COMFORTdrive





# TryKaVo.com

## Try the Best Handpieces, Risk Free



Step 1.

### CHOOSE



Choose a handpiece from  
our collections listed

Step 2.

### TRY



Try a handpiece risk free  
for 5 business days

Step 3.

### RETURN



After 5 day trial, return the  
handpiece. Like it, buy it.  
Contact 1-877-746-3995



Go to **TRYKAVO.COM** to start  
your free and easy trial today!

KaVo Dental Corporation  
11729 Fruehauf Drive  
Charlotte, NC 28273  
888.ASK.KAVO (888.275.5286)  
Fax 704.927.0816  
[www.kavousa.com](http://www.kavousa.com)



KaVo. Dental Excellence.