


# Parts Book




# Table of Contents

Table of Contents	2
Table of Contents	3
Statim 2000 Parts Front View	4
Statim 2000 Parts Rear View	5
Statim 2000 Parts Bottom View	6
Statim 2000 Parts Top View No Cover	7
Statim 2000 Revision 5XX & Higher Software Alex Parts	8
Statim 2000 Individual Parts Pictures	9-16
Statim 5000 Parts Front View	17
Statim 5000 Parts Rear View	18
Statim 5000 Parts Top View No Cover	19
Statim 5000 Revision 5XX & Higher Software Alex Parts	20
Statim 5000 Individual Parts Pictures	21-27
Statim 7000 Front View	28
Statim 7000 Front Accessories & Rear View	29
Statim 7000 Front & Side View No Cover	30
Statim 7000 Rear View No Cover	31
Statim 7000 Cassette	32
Statim 7000 Condenser, Waste & Overflow Bottles	33
Statim 900 Parts Front View	34
Statim 900 Parts Rear View	35
Statim 900 Parts Left Side View	36
Statim 900 Parts Right Side View	37
Statim 900 Parts Top View	38
Statim 900 Individual Parts Pictures	39-42
Bravo Accessories	43
Bravo Parts Table 1	44
Bravo Parts List for Table 1	45
Bravo Parts Table 2	46
Bravo Parts List for Table 2	47
Bravo Parts Table 3	48
Bravo Parts List for Table 3 & Bravo Software	49
Bravo Parts Table 4	50
Bravo Parts List for Table 4	51
Bravo Parts Table 4A & 4B	52
Bravo Parts List for Tables 4A & 4B	53
Bravo Parts Table 5	54
Bravo Parts List for Table 5	55
Bravo Parts Table 6	56
Bravo Parts List for Table 6	57
Bravo Parts Table 7	58
Bravo Parts List for Table 7	59
Bravo Parts PCB Fuses	60

# Table of Contents


Quantim Parts Front View Door Closed	61
Quantim Parts Front View Door Open	62
Quantim Parts Rear View with Cover	63
Quantim Parts Right Side No Cover	64
Quantim Parts Right Side Front View	65
Quantim Parts Rear View No Cover	66
Quantim Parts Left Side Front & Back No Cover	67
Quantim Parts Left Side Top Front No Cover	68
Quantim Parts Bottom View	69
Hydrim C51W Parts Front View Door Closed	70
Hydrim C51W Parts Front View Door Open	71
Hydrim C51W Parts Right Side	72
Hydrim C51W Right Side View & Inside Chamber with Dryer	73
Hydrim C51W Parts Bottom	74
Hydrim C51W Spare Parts Update 06-TSB-230	75
Hydrim C51W w/Touchscreen	76-78
Hydrim L110W Parts Front View Door Closed	79
Hydrim L110W Parts Front View Door Open	80
Hydrim L110W Parts Left Side No Cover	81
Hydrim L110W Parts Right Side No Cover	82
Hydrim L110W Parts Inlet Solenoids & Controllers	83
Hydrim L110W Preventative Maintenance Kit	84
Hydrim L110W w/Glass Door Front View	85
Hydrim L110W w/Glass Door Chamber w/Trolleys	86
Hydrim L110W w/Glass Door Chamber	87
Hydrim L110W w/Glass Door Behind Front Kickplate	88
Hydrim L110W w/Glass Door Control Boards	89
Hydrim L110W w/Glass Door Right Side Top View	90
Hydrim L110W w/Glass Door Right Side Bottom View	91
Hydrim L110W w/Glass Door Rear View	92
Hydrim L110W w/Glass Door Top View	93
Statmatic Parts	94-98
Statmatic Couplers	99
Aquastat Boiler	100
Aquastat Chassis	101
Aquastat Reservoir	102
Aquastat Individual Parts Pictures	103-106
VistaBrite	107-109

## Statim 2000 Front View




Item #	Part #	Description
1	01-101611S	Cover without Facia
2	01-100564S	Cover Complete for Model 101102
2	01-106121S	Cover Complete for Model 01-121101
3	01-100563S	Facia Complete for Model 101102
3	01-106122S	Facia Complete for Model 01-121101
3	01-103518S	Facia without LCD for Model 101102
3	01-106124S	Facia without LCD for Model 01-121101
4	01-100008A	Cassette Complete
5	01-104381S	LCD includes Decal (located inside Facia)
6	01-106950S	LCD Decal Only
7	01-101783S	Reservoir Cap & Filter
7	01-109300S	Filter Only, Reservoir (not shown under Reservoir Cap)
8	01-101707S	Keypad (old style not shown)
8	01-109680S	Keypad (pictured above)
9	01-100842S	Facia Gasket (not shown located inside Facia)
10	01-106675S	Label Kit Qty 10
11	01-103856S	Isoplate (not shown, not field replaceable requires alignment jig)

## Statim 2000 Rear View




Item #	Part #	Description
1	01-100564S	Cover Complete (Model #101102)
1	01-106121S	Cover Complete (Model #01-121101)
2	01-102119S	Biological Filter (not on all units)
3	01-104284S	Biological Filter Bracket (behind Biological Filter)
4	01-100207S	Air Filter (not shown behind cover on Air Compressor)
5	01-107817S	Fitting 3/8" for Biological Filter Tubes (Pkg. of 2)
6	01-100782S	Push In Fitting
7	01-106087S	Line Filter
8	01-100573S	Power Switch
9	01-101783S	Reservoir Cap & Filter
9	01-109300S	Filter Only, Reservoir (not shown under Reservoir Cap)

## Statim 2000 Bottom View



Item #	Part #	Description
1	01-100780S	Bumper
2	01-101553S	Leveler
2	01-104180S	Leveler Repair Kit (includes Leveler and Insert for Leveler)
3	01-108152S	Drain Tubing Kit
4	01-104343S	Plug Drain Tubing
5	74-100335A	Drain Tube 1"


## Statim 2000 Revision 4XX & Lower Software Top View No Cover




Item #	Part #	Description
1	01-100684S	Bubble Level
2	01-100561S	Armature (not field replaceable requires alignment jig)
3	01-100574S	Check Valve
4	01-100560S	Compressor
5	01-106787S	Pressure Relief Valve
6	01-104461S	Thermal Fuse (not shown, under Boiler)
7	01-108979S	Boiler (software must be 5XX or higher to replace in field)
8	01-100571S	Microswitch (not shown, behind Boiler)
9	01-100570S	Probe Bracket (see note)
10	01-100782S	Push In Fitting
11	01-100557S	Solenoid
12	01-103571S	Water Quality Sensor
13	01-100568S	Pump
14	01-106087S	Line Filter
15	01-100573S	Power Switch
16	01-108541S	Power Cable Filter/Switch
17	01-106349S	PCB Power Cable
18	01-110868S	Reservoir with Water Quality Sensor
18	01-100565S	Reservoir with Float Switch (not shown)
18	01-106459S	Reservoir with Float Switch & Sensor (not shown)
19	01-109365S	PCB Rev. 7 (see note)

**Note: This part is not field replaceable requires complete calibration**

## Statim 2000 Revision 5XX & Higher Software Alex Parts



Item #	Part #	Description
1	01-108979S	Boiler (software must be 5XX or higher to replace in field)
2	01-108983S	Validation Thermocouple 2000



Item #	Part #	Description
1	01-109365S	PCB Rev. 7 (not shown) (see note)
2	01-108981S	Temperature Adapter Board for Alex
3	01-108985S	PCB Adapter Board Revision 5XX Software Alex 2000
4	01-108976S	Microprocessor Kit Revision 5XX
4	01-109385S	Microprocessor Kit Revision 6XX for Rev. 7 PCB

**Note: This part is not field replaceable requires complete calibration**

## Statim 2000 Parts



Air Filter  
#01-100207S



Alex Kit 2XX Software  
#01-108996S



Alex Kit 4XX Software  
#01-108993S



Alex Kit PCB Rev 3.0 &  
4.2 #01-109716S



Allen Key 9/64  
#01-103469S



Armature  
#01-100561S



Biological Filter  
#01-102119S



Biological Filter Bracket  
#01-104284S



Boiler Alex  
#01-108979S



Bubble Level  
#01-100684S



Bumper  
#01-100780S



Cassette Complete  
#01-100008A

## Statim 2000 Parts



Cassette Lid  
#01-100834S



Cassette Seal & Lube Kit  
#01-100028S



Cassette Tray Only  
#01-103475S



Cassette Tray with Rack  
#01-100271A



Chassis Kit (Frame)  
#01-108388S



Check Valve  
#01-100574S



Compressor  
#01-100560S



Conductivity Meter  
#01-103139S



Control Box  
#103141S



Copper Tubing  
Validation Thermocouple  
#01-109128S



Cover Complete  
Model 01-121101  
#01-106121S



Cover Complete  
Model 101102  
#01-100564S

## Statim 2000 Parts



Cover without Facia  
#01-101611S



Drain Tube 1"  
#74-100335A



Drain Tubing Kit  
#01-108152S




Exhaust Tubing  
#01-100204S



Facia Complete  
Model 01-121101  
#01-106122S



Facia Complete  
Model 101102  
#01-100563S



Facia Gasket  
#01-100842S



Facia without LCD  
Model 01-121101  
#01-106124S



Facia without LCD  
Model 101102  
#01-103518S



Filter Reservoir  
#01-109300S



Fitting 3/8 for Biological  
Filter Tubes  
#01-107817S



Fuse 2A  
#01-103498S


## Statim 2000 Parts




**Fuse 15A**  
**#01-103472S**




**Fuse .25A**  
**#01-103497S**



**Isoplate**  
**#01-103856S**




**Keypad Model 01-121101**  
**#01-109680S**



**Keypad Model 101102**  
**#01-101707S**



**Label Kit**  
**#01-106675S**



**LCD**  
**#01-104381S**



**LCD Decal**  
**#01-106950S**



**Leveler**  
**#01-101553S**



**Leveler Repair Kit**  
**#01-104180S**



**Line Filter**  
**#01-106087S**




**Microprocessor 4XX**  
**#01-108742S**

## Statim 2000 Parts



**Microprocessor 5XX**  
**#01-108976S**



**Microswitch**  
**#01-100571S**



**Operator's Manual**  
**#01-106938S**



**Packaging with Inserts**  
**#01-108348S**




**PCB 4XX**  
**#01-108741S**



**PCB 5XX**  
**#01-109043S**



**PCB Adapter Board 4XX**  
**#01-108775S**



**PCB Adapter Board 5XX**  
**#01-108985S**



**PCB Power Cable**  
**#01-106349S**



**Plug Drain Tubing**  
**#01-104343S**



**Power Cable Filter/Switch**  
**#01-108541S**



**Power Cord**  
**#01-100572S**

## Statim 2000 Parts




**Power Cord 2000**  
**#01-101647S**



**Power Switch**  
**#01-100573S**



**Pressure Relief Valve (PRV)**  
**#01-106787S**



**Probe Bracket**  
**#01-100570S**



**Probe Bracket Gasket**  
**#01-100865S**



**Pump**  
**#01-100568S**




**Pump Bracket**  
**#01-108363S**




**Pump Tester Bottle**  
**#01-100713S**



**Push In Fitting**  
**#01-100782S**



**Rack Unwrapped**  
**#01-103945S**



**Rack Wrapped (Mesh)**  
**#01-106653**




**Reservoir Cap with Filter**  
**#01-101783S**

## Statim 2000 Parts



**Reservoir with Float**  
**#01-100565S**



**Reservoir with  
Water Quality Sensor**  
**#01-110868S**



**Solenoid Coil**  
**#01-104303S**



**Solenoid Plunger Wrench**  
**#01-103471S**




**Solenoid Valve Complete**  
**#01-100557S**



**Solenoid Valve Repair Kit**  
**#01-100998S**



**Strain Relief**  
**#53-100219A**



**Temperature Adapter  
Board**  
**#01-108981S**



**Test Strip Indicators**  
**#SCI134**



**Thermal Fuse**  
**#01-104461S**



**Validation Thermocouple**  
**#01-108983S**



**Waste Bottle Complete**  
**#01-100812S**

## Statim 2000 Parts



**Waste Bottle Fitting**  
**#01-100735S**




**Waste Bottle Only**  
**#01-100724S**




**Water Quality Sensor**  
**#01-103571S**

## Statim 5000 Front View




<b>Item #</b>	<b>Part #</b>	<b>Description</b>
1	01-101783S	Reservoir Cap & Filter
1	01-109300S	Filter Only, Reservoir (not shown under Reservoir Cap)
2	01-104300S	Cover Complete
3	01-104299S	Facia Complete
4	01-101650S	Facia Gasket (not shown located inside Facia)
5	01-101613S	Cassette Complete
6	01-103857S	Isoplate (not shown, not field replaceable requires alignment jig)
7	01-104382S	LCD (located inside Facia)
8	01-109671S	Keypad
9	01-101625S	Cover without Facia
10	01-101658S	Cassette Handles (Lid & Tray)

## Statim 5000 Rear View



Item #	Part #	Description
1	01-101755S	Push In Fitting
2	01-101652S	Air Compressor Filter (not on all units)
3	01-104284S	Biological Filter Bracket
4	01-102119S	Biological Filter
5	01-103993S	Power Socket (older units)
5	01-106087S	Line Filter (newer units)
6	01-100573S	Power Switch


## Statim 5000 Parts Top View No Cover




Item #	Part #	Description
1	01-100684S	Bubble Level
2	01-101623S	Armature (not field replaceable requires alignment jig)
3	01-106787S	Pressure Relief Valve
4	01-108980S	Boiler (software must be 5XX or higher to replace in field)
5	01-104462S	Thermal Fuse (not shown, under boiler)
6	01-101628S	Solenoid
7	01-101627S	Check Valve
8	01-101621S	Microswitch
9	01-101618S	Probe Bracket (see note)
10	01-101652S	Air Compressor Filter (not on all units)
11	01-102119S	Biological Filter
12	01-101761S	Pump
13	01-101619S	Compressor
14	01-103571S	Water Quality Sensor
15	01-110869S	Reservoir with Water Quality Sensor
15	01-101622S	Reservoir with Float Switch & Sensor (not shown)
16	01-109395S	PCB Rev. 7 for 5000 no printer USA (see note)
16	01-109823S	PCB Rev. 7 for 5000 with printer USA (see note)

**Note: This part is not field replaceable requires complete calibration**

## Statim 5000 Revision 5XX & Higher Software Alex Parts



Item #	Part #	Description
1	01-108980S	Boiler (software must be 5XX or higher to replace in field)
2	01-108984S	Validation Thermocouple 5000



Item #	Part #	Description
1	01-109395S	PCB Rev. 7 for 5000 no printer USA (see note)
1	01-109823S	PCB Rev. 7 for 5000 with printer USA (see note)
2	01-108981S	Temperature Adapter Board for Alex
3	01-108987S	PCB Adapter Board Alex 5000 5XX & Higher Software
4	01-108978S	Microprocessor Kit Alex 5000 Revision 5XX & Higher Software

**Note: This part is not field replaceable requires complete calibration**

## Statim 5000 Parts



Air Compressor Filter  
#01-101652S



Air Compressor  
Filter Holder  
#01-104329S



Alex Kit 2XX  
#01-108998S



Alex Kit 4XX  
#01-108995S



Armature  
#01-101623S



Boiler Alex  
#01-108980S



Bubble Level  
#01-100684S




Bumper  
#01-100780S



Cassette Complete  
#01-101613S



Cassette Extended  
#01-104104



Cassette Handles  
#01-101658S



Cassette Handles Extended  
#01-104472S

## Statim 5000 Parts



Cassette Lid  
#01-101757S



Cassette Seal & Lube Kit  
#01-101649S



Cassette Tray with Rack  
#01-101614S



Chassis Kit (Frame)  
#01-108389S



Check Valve  
#01-101627S



Compressor  
#01-101619S



Conductivity Meter  
#01-103139S



Control Box  
#01-103141S



Copper Tubing  
Validation Thermocouple  
#01-109129S



Cover Complete  
#01-104300S




Cover without Facia  
#01-101625S



Drain Tube 1'  
74-100335A

## Statim 5000 Parts



Drying Plates  
#01-103935



Exhaust Tubing  
#01-100204S



Facia Complete  
#01-104299S



Facia Gasket  
#01-101650S



Filter, Reservoir  
#01-109300S



Fuse 2A  
#01-103498S



Fuse 15A  
#01-103472S




Fuse .25A  
#01-103497S



Isoplate  
#01-103857S



Keypad  
#01-109671S



LCD  
#01-104382S



Leveler  
#01-101553S

## Statim 5000 Parts



**Leveler Repair Kit  
#01-104180S**




**Line Filter  
#01-106087S**



**Microprocessor 4XX  
#01-108746S**



**Microprocessor 5XX  
#01-108978S**



**Microswitch  
#01-101621S**



**Operator's Manual  
#01-106938S**



**Packaging with Inserts  
#01-108327S**



**PCB 4XX  
#01-108744S**



**PCB 5XX  
#01-109045S**



**PCB Adapter Board 4XX  
#01-108777S**



**PCB Adapter Board 5XX  
#01-108987S**



**PCB Power Cable  
#01-106350S**

## Statim 5000 Parts



**Plug Drain Tubing**  
**#01-104343S**



**Power Cable Filter/Switch**  
**#01-108541S**



**Power Cord**  
**#01-101647S**



**Power Socket**  
**#01-103993S**



**Power Switch**  
**#01-100573S**




**Pressure Relief Valve**  
**#01-106787S**



**Printer 5XX> Software**  
**#01-210000**




**Printer Case with Damper**  
**#01-103470S**



**Printer Door**  
**#01-102159S**



**Printer Paper (Pkg. 10)**  
**#01-101657S**



**Printer Paper Holder**  
**#01-101973S**



**Printer Spring/Latch**  
**#01-101971S**

## Statim 5000 Parts




**Probe Bracket**  
**#01-101618S**



**Probe Bracket Gasket**  
**#01-100865S**



**Pump**  
**#01-101761S**



**Pump Recovery Kit**  
**#01-104313S**




**Pump Tester Bottle**  
**#01-100713S**



**Push In Fitting**  
**#01-101755S**



**Rack Unwrapped Mesh**  
**#01-101709S**



**Rack Wrapped**  
**#01-102054S**



**Reservoir w/Sensor**  
**#01-110869S**



**Reservoir Cap with Filter**  
**#01-101783S**



**Solenoid Coil**  
**#01-104303S**



**Solenoid Plunger Wrench**  
**01-103471S**

## Statim 5000 Parts



**Solenoid Valve Complete**  
**01-101628S**



**Solenoid Valve Repair Kit**  
**01-100998S**




**Temperature Adapter  
Board**  
**01-108981S**



**Test Strip Indicators**  
**SCI134**



**Thermal Fuse**  
**01-104462S**



**Validation Thermocouple**  
**#01-108984S**



**Waste Bottle Complete**  
**01-100812S**



**Waste Bottle Fitting**  
**01-100735S**



**Waste Bottle Only**  
**01-100724S**







**Water Quality Sensor**  
**01-103571S**

## Statim 7000 Front View







Item #	Part #	Description
1	01-110298S	Cover Complete
2	01-110523S	Cover Facia Front Cassette
3	01-110520S	Cover Reservoir Compartment
4	01-110522S	Cover Facia Front Left
5	01-110521S	Cover Logo SciCan
6	01-110526S	Cover Keypad Complete
7	01-110525S	Cover LCD Only
8	01-110524S	Cover Facia Gasket (not shown located inside cover)
9	01-110293S	Reservoir Complete
10	01-110328S	Reservoir Lid
11	01-110512S	Leveling Foot
12	01-110288S	Cassette Complete

## Statim 7000 Front Accessories & Rear View





Item #	Part #	Description
1	SCWF1	Water Filter Cartridge
	SCWF6	Water Filter Cartridge (Package of 6)
	SCWFBP	Water Bypass Cartridge (not shown)
2	01-110293S	Reservoir Complete
3	01-102119S	Biological Filter
4	01-110503S	Push In Fitting
5	01-101652S	Air Compressor Filter
6	01-104329S	Air Compressor Filter Holder
7	01-104343S	Plug Drain Tubing
8	01-110505S	Line Filter (Plug Receptacle)
9	01-110504S	Power Switch
10	01-110327S	Reservoir Plunger Kit w/Large Outer O-ring 7000

## Statim 7000 Parts









Item #	Part #	Description
1	01-102119S	Biological Filter
2	01-110472S	Bio Filter Holder
3	01-110496S	Armature
4	01-110498S	Armature Isoplate
5	01-110493S	Reservoir Inlet
6	01-110488S	PCB RFID Antenna Assembly
7	01-110502S	Armature Microswitch
8	01-110515S	Thermocouple Validation
9	01-110471S	Pump, Water
10	01-110497S	Armature Guide Rails (Pkg. Of 2)
11	01-110490S	Fan
12	01-110482S	PCB Power Board Auxiliary Heater
13	01-110716S	Fuse 10A 250V (Pkg. Of 2)
14	01-103498S	Fuse 2A
15	01-110477S	Auxiliary Heater Triac (Solid State Relay)
16	01-109816S	Fuse 0.8A
17	01-110483S	PCB RFID Board
18	01-110480S	PCB 7000 ( <b>not field replaceable</b> )

## Statim 7000 Rear View No Cover






Item #	Part #	Description
1	01-110474S	Auxiliary Heater
2	01-110527S	Solenoid Chamber (Valve 2)
3	01-110513S	Thermocouple Inlet/Outlet
4	01-110528S	Solenoid Validation (Valve 1)
5	01-110514S	Thermocouple Outlet ( <b>not field replaceable</b> )
6	01-110472S	Boiler 7000
7	01-110510S	Check Valve
8	01-110471S	Pump, Water
9	01-110478S	Compressor

## Statim 7000 Cassette




Item #	Part #	Description
1	01-110288S	Cassette Complete
2	01-110290S	Cassette Lid
3	01-110289S	Cassette Tray
4	01-110291S	Cassette Inlet Coupling
5	01-110292S	Cassette Outlet Coupling
6	01-110296S	Cassette Coupling Seal Kit (Seals for both Couplings included)
7	01-110294S	Cassette Perforated Rack
8	01-110330S	Cassette Tray Handle
9	01-110297S	Cassette Tray Exhaust Duct
10	01-110825S	Cassette Steam Deflector
11	01-106848S	Filter Exhaust Duct Insert
12	01-110824S	Cassette Steam Deflector Clip
13	01-110295S	Cassette Seal Kit
14	01-110329S	Cassette Lid Handle
Not Shown	01-111590S	Cassette O-ring Lubricator Kit

## Statim 7000 Condenser, Waste & Overflow Bottles



Item #	Part #	Description
1	01-110506S	Bottle Waste Complete
2	01-100735S	Waste Bottle Fitting
3	01-111497S	Bottle Tubing Condenser to Waste (Includes Quick Disconnect)
4	01-111660S	Bottle Condenser Fitting 3/8
5	01-111395S	Bottle Condenser
6	01-110508S	Bottle Waste Condenser Fitting Quick Disconnect
7	01-110300S	Bottle Waste Overflow Kit

## Statim 900 Parts Front View




Item #	Part #	Description
1	01-101554S	Bubble Level
2	01-103527S	Keypad
3	01-103508S	Facia
4	01-103636S	Wire Basket W/Insert
5	01-103608S	Sterilization Chamber
6	01-101560S	Drawer Front
7	01-101562S	Drawer
8	01-103559S	Reservoir Cap
9	01-101558S	Cover

## Statim 900 Parts Rear View




Item #	Part #	Description
1	01-104343S	Plug Drain Tubing
2	01-100573S	Power Switch
3	01-106087S	Line Filter
4	01-100782S	Push In Fitting

## Statim 900 Parts Left Side View




Item #	Part #	Description
1	01-103508S	Facia
2	01-101542S	Reservoir
3	01-103608S	Sterilization Chamber
4	01-101562S	Drawer
5	01-110242S	Boiler (software must be 5XX or higher to replace in field)
6	01-102038S	Pump
7	01-101550S	Solenoid
8	01-100782S	Push In Fitting

## Statim 900 Parts Right Side View



Item #	Part #	Description
1	01-103508S	Facia
2	01-101542S	Reservoir
3	01-103608S	Sterilization Chamber
4	01-101562S	Drawer
5	01-101552S	Microswitch
6	01-101552S	Microswitch

## Statim 900 Parts Top View



Item #	Part #	Description
1	01-103508S	Facia
2	01-101542S	Reservoir
3	01-110242S	Boiler (software must be 5XX or higher to replace in field)
4	01-103085S	Pressure Relief Valve
5	01-101550S	Solenoid
6	01-100782S	Push In Fitting
7	01-110462S	PCB Rev. 7 (see note)
8	01-102038S	Pump
9	01-102061S	Motor Control PCB
10	01-104461S	Thermal Fuse (not shown, under boiler)
11	01-101545S	Motor

**Note: This part is not field replaceable requires complete calibration**

## Statim 900 Parts



**Basket**  
**#01-103636S**



**Boiler Alex**  
**#01-110242S**



**Bubble Level**  
**#01-101554S**



**Cap Replacement**  
**#01-101582S**




**Chamber Sterilization**  
**#01-103608S**



**Cover**  
**#01-101558S**



**Drawer**  
**#01-101562S**



**Drawer Front**  
**#01-101560S**



**Facia**  
**#01-103508S**



**Filter Insert Exhaust**  
**#01-106848S**



**Fuse .25A**  
**#01-103497S**



**Fuses 15A**  
**#01-103472S**

## Statim 900 Parts



Fuse 2A  
#01-103498S



Keypad  
#103527S




LCD  
#01-103138S



Leveler  
#01-101553S




Leveler Repair Kit  
#01-104180S



Manual Op  
#01-106847S



Microswitch  
#01-101552S



Motor  
#01-101545S



Motor Control PCB  
#01-102061S



Plug Drain Tubing  
#01-104343S



Power Cord  
#01-101647S



Power Switch  
#01-100573S

## Statim 900 Parts



**Pressure Relief Valve**  
**#01-103085S**



**Pump**  
**#01-102038S**




**Push-in Fitting**  
**#01-100782S**



**Reservoir**  
**#01-101542S**



**Reservoir Cap**  
**#01-103559S**



**Seal (Exhaust Only)**  
**#01-103144S**



**Seal (Pressure Only)**  
**#01-103143S**



**Seal Assembly**  
**(Complete w/Plate)**  
**#01-106898S**



**Seal Frame**  
**#01-102064S**



**Seal Kit**  
**(Exhaust & Pressure)**  
**#01-101543S**



**Seal Plate Flat**  
**#01-102065S**



**Seal Screw**  
**#01-102063S**

## Statim 900 Parts



**Solenoid Coil**  
**#01-104303S**



**Solenoid Valve Complete**  
**#01-101550S**



**Solenoid Valve Repair Kit**  
**#01-100998S**



**Thermal Fuse**  
**#01-104461S**




**Waste Bottle Complete**  
**#01-100812S**




**Water Filter In-Line**  
**#01-106637S**

# Bravo Accessories



Item #	Part #	Description
1	C1XP076000Y	Bravo Tray for 17 & 17V
1	C1BG345000Y	Bravo Tray for 21V
2	C1BP583000Y	Bravo Rack for 17 & 17V
2	C1BG534000Y	Bravo Rack for 21V
3	STXX0080000	Bravo Tray Handle
4	A0XP0010000	Bravo Fill Bottle 2500 ml
5	110000003W0	Bravo Drain Tubing w/Fitting
6	47200010000	Bravo Biological Filter
7	47200030000	Bravo Chamber Filter Screen Metal
8	STXX0250000	Bravo Printer Paper Thermal Roll
9	261000009MN	Bravo Elbow Barbed Fitting 1/8"
	SD-306	Bravo Operator's Manual (not shown)


# Bravo Parts Table 1



# Parts List for Table 1

<b>Item #</b>	<b>Table #</b>	<b>Part #</b>	<b>Description</b>
1	1	C1XP076000Y	Tray for Bravo 17 & 17V
1	1	C1BG345000Y	Tray for Bravo 21 & 21V
2	1	C1BP583000Y	Rack for Chamber Bravo 17 & 17V
2	1	C1BG534000Y	Rack for Chamber Bravo 21 & 21V
3	1	STXX0080000	Tray Extractor Handle Bravo
4	1	A0XP0010000	Reservoir Fill Bottle 2500 ml Bravo
5	1	110000003W0	Drain Tubing w/Fitting Bravo
6	1	47200010000	Biological Filter Bravo
7	1	47200030000	Water Filter Screen Bravo
9	1	C1BH0060001	Cover for Bravo 17 & 17V
9	1	C1BG4680001	Cover for Bravo 21 & 21V
10	1	C3BP1850000	Door Plastic Cover Bravo
11	1	C5BP5500000	Printer Thermal Bravo
11a	1	STXX0080000	Printer Paper Thermal Roll Bravo
12	1	C3BP1840000	Printer Door Bravo
13	1	230B04L20K0	Printer Door Pin 4 x 20 Bravo
14	1	42000070000	Power Switch Bravo
15	1	C3BP1830000	Facia Front Plastic Bravo
16	1	C5BP1230000	LCD Bravo
17	1	C6JP0010000	Keypad Bravo
18	1	C6JM0010000	Label Bravo 17
18	1	C6JP0020000	Label Bravo 17V
18	1	C6JG0010000	Label Bravo 21V
19	1	47200010000	Biological Filter Bravo
20	1	C1BH005000M	Support Bar Left Bravo 17 & 17V
20	1	C1BG317000M	Support Bar Left Bravo 21 & 21V
21	1	C1BP152000M	Support Bar Right Bravo 17 & 17V
21	1	C1BG318000M	Support Bar Right Bravo 21 & 21V


## Bravo Parts Table 2



## Parts List for Table 2

<b>Item #</b>	<b>Table #</b>	<b>Part #</b>	<b>Description</b>
22	2	C4BP2740000	Door Thermo insulating panel Bravo
23	2	C0BP260000P	Door Front Metal Bravo
24	2	C7BP2640000	Door Spring Bravo
25	2	C0BP060000Z	Door Fork Bravo
26	2	49000050000	Door Bushing Bravo
27	2	C0BP0500000	Door Fork Pin Bravo
28	2	203C06L30K0	Door Center Screw Bravo
29	2	C0BP079000N	Door Plate Adjustment Screw Bravo
30	2	C4BP5690000	Door Plate Thermo Insulating Bravo
31	2	A1BP027000Y	Door Seal Plate Bravo
32	2	48000050000	Door Gasket Bravo
33	2	C0BP2720000	Door Hinge Pin Bravo
34	2	A1BP5520000	Door Hinge
35	2	C2BP0010000	Chamber Bravo 17 & 17V
35	2	C2BG5300000	Chamber Bravo 21 & 21V
36	2	C0BP0430000	Door Locking System Pin Bravo
37	2	A1BP2710000	Door Locking System Complete Bravo
38	2	40300010000	Door Locking System Motor
39	2	232A03L25K0	Door Locking System Pin 3 x 25 Bravo
40	2	C0BP0450000	Door Locking System Motor Coupling Bravo
41	2	C0BP0460000	Door Locking System Intermediate Coupling Bravo
42	2	C0BP0590000	Door Locking System Closing Hook Bravo
43	2	43300010000	Door Locking System Hook Microswitch (open) Bravo
44	2	43300010000	Door Locking System Door Microswitch Bravo
45	2	43300040000	Door Locking System Hook Microswitch (close) Bravo
46	2	28600160000	Temperature Sensor Ferrule Bravo
47	2	2860014000N	Temperature Sensor Nut
48	2	43000190000	Temperature Sensor PT1000 – PT1 Bravo
49	2	41000140000	Heating Element Bravo
50	2	43000180000	Temperature Sensor PT1000 – PT3 Bravo
51	2	43000140000	Safety Thermostat TG400 Bravo
52	2	C4BP2140100	Thermo Insulating Cover for Chamber Bravo 17 & 17V
52	2	C4BG3280000	Thermo Insulating Cover for Chamber Bravo 21 & 21V
53	2	C4BP5570000	Thermo Insulating Ring for Chamber Bravo


# Bravo Parts Table 3



# Parts List for Table 3

<b>Item #</b>	<b>Table #</b>	<b>Part #</b>	<b>Description</b>
<b>54</b>	3	A0BP2900000	Reservoir Assembly Complete Bravo
<b>54a</b>	3	C3BP1700000	Reservoir and Waste Water Tank Bravo
<b>54b</b>	3	43200060000	Float Switch Waste Water Tank Bravo
<b>54c</b>	3	43200050000	Float Switch Water Reservoir Bravo
<b>55</b>	3	C1BP175000M	Support Bar Reservoir Bravo 17 & 17V
<b>55</b>	3	C1BG320000M	Support Bar Reservoir Bravo 21 & 21V
<b>56</b>	3	A9BP5860000	Transformer Complete Bravo 120V
<b>56</b>	3	A9BP5970000	Transformer Complete Bravo 220V
<b>56a</b>	3	41200140000	Transformer Only Bravo 120V
<b>56a</b>	3	41200080000	Transformer Only Bravo 220V
<b>56b</b>	3	C1BP161000M	Transformer Support Plate Bravo
<b>57</b>	3	43100060000	Pressure Switch Bravo
<b>58</b>	3	C1BP538000M	Pressure Switch Support Bravo
<b>59</b>	3	260000009MN	Pressure Switch/Transducer Fitting Bravo
<b>60</b>	3	43100090000	Pressure Transducer Bravo
<b>61</b>	3	C5BP5630000	Basic Board (GAM) Bravo 120V/220V
<b>61c</b>	3	C5BP4350000	PCB Filter Board 120V Bravo
<b>61c</b>	3	C5BP1420000	PCB Filter Board 220V Bravo
<b>61d</b>	3	C5BP1430000	PCB Printer Board Bravo
Not	Shown	SB17V120V	Bravo Software – 120V, 17V
Not	Shown	SB17V220V	Bravo Software – 220V, 17V
Not	Shown	SB17120V	Bravo Software – 120V, 17
Not	Shown	SB17220V	Bravo Software – 220V, 17
Not	Shown	SB21V120V	Bravo Software – 120V, 21V
Not	Shown	SB21V220V	Bravo Software – 220V, 21V


# Bravo Parts Table 4




# Parts List for Table 4

<b>Item #</b>	<b>Table #</b>	<b>Part #</b>	<b>Description</b>
63	2	28600290000	Filter Fitting Chamber Front Bravo
65	2	47200030000	Filter Screen Chamber Front LP1 Bravo
67	4	48300010000	Gasket Rubber Bravo
68	4	260000006V0	Fitting CPC Female Water Tap Chamber Front Bravo
69	4	A2BP2380000	Printer Plug RS232 9 pin w/Cable Bravo
70	4	C3BP5610000	Plate ABS Bravo
71	4	25600000400	Foot Front PVC Bravo
72	4	25600001000	Foot Rear Plastic Bravo
73	4	55V702016000	Plate Ground Bravo
74	4	C1BP147000M	Support Electronic Angle Bravo
75	4	C1BP269000M	Support Front Top Angle Bar Bravo
76	4	C1BP270000M	Support Front Side Angle Bar Bravo
77	4	C1BP150000M	Chassis Bravo 17 & 17V
77	4	A1BG350000M	Chassis Bravo 21 & 21V
78	4	A0BP2860000	Water Pump Complete w/Solenoid Bravo
78a	4	40000220000	Water Pump Only Bravo
78b	4	40100370000	Valve EV6 for Water Pump Bravo
79	4	261000009MN	Fitting Elbow for Water Pump Solenoid Bravo
80	4	47200020000	Water Filter Plastic for Water Pump Bravo
81	4	111000006W0	Water Pump Adhesive Mousse Bravo
82	4	25600000800	Water Pump Mounting Bracket Foot Bravo
83	4	40000050000	Water Pump Mounting Bracket Bravo
84	4	A9BM5510000	Steam Generator Complete Bravo 120V
84	4	A9BP5250000	Steam Generator Complete Bravo 220V
84a	4	C1BP136000M	Steam Generator Support Bravo
84b	4	C0BP1310000	Steam Generator Support Spacer Insulator Bravo
84c	4	C4BP2360000	Steam Generator Thermo Insulating Panel End Bravo
84d	4	C4BP5700000	Steam Generator Thermo Insulating Cover Bravo
84e	4	C0BP548000P	Steam Generator Lower Part Bravo
84f	4	48100080000	Steam Generator O-ring Bravo
84g	4	C0BP547000P	Steam Generator Upper Part Bravo
85	4	41000170000	Steam Generator Heating Element Bravo 115V 850W
85	4	41000130000	Steam Generator Heating Element Bravo 220V 1000W
86	4	43000180000	Temperature Sensor PT1000 – PT2 Bravo
87	4	2860014000N	Steam Generator Temperature Sensor Nut Bravo
88	4	28600160000	Steam Generator Temperature Sensor Ferrule Bravo
89	4	43000150000	Steam Generator Safety Thermostat TG400 Bravo
90	4A	A9BM5520000	Vacuum Pump Small Complete Bravo 120V
90	4A	A9BM5270000	Vacuum Pump Small Complete Bravo 220V
90	4B	A9BP5870000	Vacuum Pump Large Complete Bravo 120V
90	4B	A9BP5270000	Vacuum Pump Large Complete Bravo 220V

## Bravo Parts Table 4A



## Bravo Parts Table 4B




## Parts List for Table 4A

<b>Item #</b>	<b>Table #</b>	<b>Part #</b>	<b>Description</b>
<b>90</b>	4A	A9BM5520000	Vacuum Pump Small Complete Bravo 120V
<b>90</b>	4A	A9BM5270000	Vacuum Pump Small Complete Bravo 220V
<b>90b</b>	4A	40000150000	Vacuum Pump Small Central Gasket Bravo
<b>90c</b>	4A	40000250000	Vacuum Pump Small, Small Membrane Bravo
<b>90d</b>	4A	48100150000	Vacuum Pump Small O-ring Bravo
<b>90e</b>	4A	40000180000	Vacuum Pump Small Shutter Bravo
<b>90f</b>	4A	40000120000	Vacuum Pump Small Head Complete Bravo

## Parts List for Table 4B

<b>Item #</b>	<b>Table #</b>	<b>Part #</b>	<b>Description</b>
<b>90</b>	4B	A9BP5870000	Vacuum Pump Large Complete Bravo 120V
<b>90</b>	4B	A9BP5270000	Vacuum Pump Large Complete Bravo 220V
<b>90a</b>	4B	2860034000N	Vacuum Pump Large Elbow Fitting Bravo
<b>90b</b>	4B	40000240000	Vacuum Pump Large Big Membrane with Screw Bravo
<b>90c</b>	4B	40000250000	Vacuum Pump Small Membrane Bravo
<b>90d</b>	4B	48100150000	Vacuum Pump Large Head O-ring Bravo









# Bravo Parts Table 5



# Parts List for Table 5

<b>Item #</b>	<b>Table #</b>	<b>Part #</b>	<b>Description</b>
<b>91</b>	5	C1BP151000M	Chassis Rear Bravo
<b>92</b>	5	A0BS0520000	Pressure Relief Valve Complete Bravo
<b>92a</b>	5	C1BS020000P	Pressure Relief Valve Support Bravo
<b>92b</b>	5	47000020000	Pressure Relief Valve Only Bravo
<b>93</b>	5	C1BP6220000	Heat Exchanger Complete Bravo
<b>93a</b>	5	40400040000	Heat Exchanger Protection Grid Bravo
<b>93b</b>	5	C1BP6220000	Heat Exchanger Only Bravo
<b>93c</b>	5	C1BP159000M	Heat Exchanger Support Bravo
<b>93d</b>	5	40400030000	Heat Exchanger Fan
<b>94</b>	5	A0BP2850000	Valve EV5 Complete Bravo
<b>95</b>	5	55V432011000	Valve EV5 Support Bravo
<b>96</b>	5	A0BP2820000	Valve EV1 Complete Bravo
<b>97</b>	5	C1BP241000M	Valve EV1 Support Bravo
<b>98</b>	5	A0BP2830000	Valve EV2/EV3 Complete Bravo
<b>99</b>	5	A0BP2840000	Valve EV4 Complete Bravo
<b>100</b>	5	C1BP173000M	Valve EV2/EV3/EV4 Support Bravo
<b>101</b>	5	A0BP2870000	Water Pump Automatic Filling Complete Bravo
<b>102</b>	5	47200020000	Water Filter Plastic for Water Pump Bravo
<b>103</b>	5	41700260000	Fuse Holder Bravo
<b>104</b>	5	41700330000	Fuse 15A 250V
<b>105</b>	5	261000014V0	Tee Fitting Bravo
<b>106</b>	5	A2BP2630000	Level Sensor Cable External Bravo
<b>107</b>	5	27700003500	Strain Relief Bravo
<b>108</b>	5	A2BP5890000	Power Cord Bravo 115V
<b>108</b>	5	A2BP6140000	Power Cord Bravo 220V
<b>109</b>	6	261000009MN	Valve EV1Elbow Fitting Bravo


# Bravo Parts Table 6



# Parts List for Table 6

<b>Item #</b>	<b>Table #</b>	<b>Part #</b>	<b>Description</b>
<b>94</b>	6	A0BP2850000	Valve EV5 Complete Bravo
<b>94a</b>	6	40100240000	Valve EV5 3 ways Fitting Bravo
<b>95</b>	6	55V432011000	Valve EV5 Support Bravo
<b>96</b>	6	A0BP2820000	Valve EV1 Complete Bravo
<b>96a</b>	6	40100240000	Valve EV1 3 ways Fitting Bravo
<b>97</b>	6	C1BP241000M	Valve EV1 Support Bravo
<b>98</b>	6	A0BP2830000	Valve EV2/EV3 Complete Bravo
<b>98a</b>	6	40100370000	Valve EV2 2 ways Fitting Bravo
<b>98b</b>	6	40100230000	Valve EV3 2 ways Fitting Bravo
<b>99</b>	6	A0BP2840000	Valve EV4 Complete Bravo
<b>99a</b>	6	40100260000	Valve EV4 3 ways Fitting Bravo
<b>100</b>	6	C1BP173000M	Valve EV2/EV3/EV4 Support Bravo
<b>109</b>	6	261000009MN	Valve EV1 Elbow Fitting Bravo
<b>110</b>	6	263000002MN	Valve EV1 Tee Fitting Female/Female/Female Bravo
<b>111</b>	6	2860034000N	Valve EV1/EV3 Pipe Fitting Bravo
<b>113</b>	6	263000003MN	Valve EV3 Tee Fitting Female/Male/Female Bravo
<b>114</b>	6	265000009MN	Valve EV3 Coupling Bravo
<b>115</b>	6	262000003MN	Valve EV2/EV4/EV5 Elbow Fitting Female/Male Bravo
<b>116</b>	6	2860034000N	Valve EV4/EV5 Elbow Fitting Bravo
<b>117</b>	6	261000004MN	Valve EV4/EV5 Pipe Fitting Bravo
<b>118</b>	6	47000010000	Valve EV5 One Way Valve Connector Bravo


# Bravo Parts Table 7




# Parts List for Table 7

<b>Item #</b>	<b>Table #</b>	<b>Part #</b>	<b>Description</b>
A1	7	A2BP2220000	Cable Microswitch Bravo
A2	7	A2BP4010000	Cable Safety Thermostat Bravo
A3	7	A2BG3940000	Cable Main Switch Bravo
A4	7	A2BM3390000	Cable Vacuum Pump (small) Bravo
A4	7	A2BP4030000	Cable Vacuum Pump (large) Bravo
A5	7	A2BG3520000	Cable Motor Reducer Bravo
A6	7	A2BP2270000	Cable Water Pump Bravo
A7	7	A2BP2280000	Cable Electro Valve Bravo
A8	7	A2BP4070000	Cable Filter Bravo
A9	7	A2BP4820000	Cable Ground Chamber Bravo
A10	7	A2BP4830000	Cable Ground Steam Generator Bravo
A11	7	A2BP4810000	Cable Ground Rear Frame Bravo
A12	7	A2BP4700000	Printer Cable Power Bravo
A13	7	A2BP4690000	Printer Cable Signal Bravo
P1	7	110000005W0	Silicon Tube 6 x 12 Black Bravo
P2	7	110000011W0	Silicon Tube 4 x 7 Transparent Bravo
P3	7	110000003W0	Silicon Tube 6 x 10 Transparent Bravo
P4	7	110000002V0	Teflon Tube 4 X 2.5 Bravo

## Bravo Parts PCB Fuses




## Quantim Front View Door Closed



Item #	Part #	Description
1	289124	Door Molding Kit
2	279159	Door Microswitch 2 Wire (Not Shown Under Door Molding)
2	279494	Door Microswitch 3 Wire (Not Shown Under Door Molding)
3	289162	Decal Quantim 16 (New Style no Liquids Cycle Not Pictured)
4	289125	Printer Door with Hinges


**Note:** Door Microswitch 2 Wire, wires are connected to switch terminals.  
Door Microswitch 3 Wire, switch terminals are connected to door PCB.

## Quantim Front View Door Open




Item #	Part #	Description
1	279011	Door Gasket (Includes Door Plate Seals)
2	279490	Acorn Nut (Pkg 2) (Pictured Above)
2	279491	Dome Nut (Pkg 2) (Not Shown in Picture)
3	279107	O-ring Kit Entry Port White (Pkg 2)
3	279467	Bolt & O-ring Entry Port
4	279527	Tubing & Fitting Kit
5	279102	Drain Tap
6	289132	Rack for Chamber 3 Tray (Wire)
7	309077	Air Filter
8	289006	Temperature Probe Chamber

## Quantum Rear View with Cover




<b>Item #</b>	<b>Part #</b>	<b>Description</b>
1	319003	Cover w/Feet
2	319001	Power Cord
3	319014	Fan Assembly Rear
4	279469	Cover Circular
5	279503	Waste Bottle Complete Quantim 16 (Bottle Not Shown)

## Quantum Right Side No Cover




Item #	Part #	Description
1	289125	Printer Door with Hinges
2	279537	Front Bezel
3	279170	Door Opening Solenoid
4	309039	Reservoir Single Use (Quantum with Waste Bottle)
4	279117	Reservoir Multi Use (Quantum no Waste Bottle)
5	309077	Air Filter & Clip
6	299051	Air Filter Nozzle
7	289134	Silencer Vacuum Pump
8	289118	Vacuum Pump

## Quantum Right Side Front View




Item #	Part #	Description
1	279225	Display Module
2	279171	Door Solenoid Spring Kit (2 Springs)
3	279170	Door Opening Solenoid

## Quantum Rear View No Cover




Item #	Part #	Description
1	309039	Reservoir Single Use (Quantum with Waste Bottle)
1	279117	Reservoir Multi Use (Quantum no Waste Bottle)
2	309077	Air Filter & Clip
3	319015	Fuse 15A
4	319001	Power Cord
5	319006	Control Module

## Quantum Left Side Lower Front No Cover




Item #	Part #	Description
1	289011	Temperature Probe Boiler
2	279155	Insulator Heating Element
3	289136	Transformer
4	279544	Boiler Thermistor

## Quantum Left Side Lower Back No Cover




Item #	Part #	Description
1	279157	ASME Safety Valve***
2	279155	Insulator Heating Element
<b>***Note: Boiler must be removed to replace ASME valve. A new boiler to vessel seal (#319022) will need to be installed when reinstalling boiler.</b>		

## Quantim Left Side Top Front No Cover




Item #	Part #	Description
1	289006	Temperature Probe Chamber
2	289015	Door Opening Spring
3	289028	Door Spring Pillars (Qty 2)

## Quantum Bottom View



<b>Item #</b>	<b>Part #</b>	<b>Description</b>
1	319019	Valve Block Kit
2	279176	Water Filter
3	319022	Seal Boiler to Vessel (Not Shown in Picture)
4	319024	Boiler w/PRV
5	279266	Thermal Cutout
6	279203	Thermal Cutout Mounting Kit
7	279155	Insulator Heating Element
8	289136	Transformer
9	289117	Heater Pad (Not Shown in Picture)

## Hydrim C51W Parts Front View Door Closed




Item #	Part #	Description
1	01-107974S	Door Facia
2	01-108790S	Door Actuator Kit
3	01-108917S	Facia Control Panel
3	01-109147S	Facia Control Panel w/Dryer (HD01FK0001 & newer)
4	01-107935S	Chemical Door w/Hinge
Not Pictured	01-108924S	Screw Kit (Cover Screws)

## Hydrim C51W Parts Front View Door Open




Item #	Part #	Description
1	01-107805S	Wash Arm Upper
2	01-108154S	Plug Upper Arm
3	01-107786S	Seal Door
4	01-107804S	Wash Arm Lower
5	01-107806S	Drain Screen
6	01-107812S	Door Latch
7	01-107935S	Chemical Door New Style
Not Pictured	01-108090S	Seal Kit Door Bottom
Not Pictured	01-107807S	Drain Filter


## Hydrim C51W Parts Right Side




**HD01DF9999 & Older**



**HD01DG0001 to HD01FJ9999**





**Old Style Dosing Reservoir**



**View behind Right Side Panel**


<b>Item #</b>	<b>Part #</b>	<b>Description</b>
1	01-108030S	Quick Connect Female
2	01-108031S	Tube Chemical (HD01DF9999 & older)
3	01-108122S	Drip Tray
4	01-109058S	Dosing Reservoir Upgrade Kit (Old Style Reservoir not available must upgrade)
5	01-107790S	Dosing Pump
6	01-108700S	Fuse Holders (Pkg. 3)
7	01-108699S	Tube Chemical (HD01DG0001 & newer)
8	01-108253S	Dosing Reservoir New Style
9	01-109909S	Tube Dosing Pump Internal
10	01-108351S	Float Dosing Reservoir
11	01-107938S	Switch Pressure Assembly
12	01-107803S	Switch Overflow Chamber
13	01-107937S	Switch Overflow Bottom
Not Pictured	01-107800S	Fuse 10 Amp (Pkg. 10)
Not Pictured	01-107799S	Fuse 2 Amp (Pkg. 5)

**Hydrim C51W Right Side View & Inside Chamber with Dryer  
HD01FK0001 & Newer**



Item #	Part #	Description
1	01-109143S	Dryer Tubing
2	01-109142S	Dryer Assembly
3	01-109145S	Vent Assembly Dryer
4	01-109144S	Dryer Fitting

## Hydrim C51W Parts Bottom



Item #	Part #	Description
1	01-107794S	Pump Recirculation Lower Arm
2	01-108701S	Tubing (blue) Chamber
3	01-107797S	Pump Drain
4	01-107808S	Water Heater
5	01-107815S	Valve Hot Water
6	01-107801S	Valve Cold Water
7	01-107809S	Thermostat Heater NC Sensor
8	01-107795S	Pump Recirculation Upper Arm

**Note:** See next page for upper and lower pump impellers and upper pump Y fittings with stainless steel barbs.

## Hydrim C51w Spare Parts Update

06-TSB-230

Please be advised of the new Hydrim C51w spare parts.

The Impellers (items # 4&5) have been set-up separately due to cost concerns of changing circulation pumps. The Y-fittings have been upgraded and are now made of stainless steel. The part numbers for the items have remained the same.



1. Stainless Steel Fitting For Y's  
p/n 01-109590S (Package of 2)

2. Y-Fitting, Upper Arm Pump  
p/n 01-109303S

3. Y- Fitting, Upper Arm  
p/n 01-108702S



4. Impeller - Upper  
p/n 01-109447S

5. Impeller - lower  
p/n 01-109448S

## Hydrim C51W w/Touchscreen Parts



01-107786S  
Seal Door, J



01-107788S  
Inlet Hose N.A., J



01-107789S  
Drain Hose, J



01-107790S  
Dosing Pump, J



01-107791S  
Valve Salt Regeneration, J



01-107794S  
Pump Recirculation Lower Arm, J



01-110259S  
Pump Upper Arm w/o  
Fitting, J



01-107797S  
Drain Pump N.A. 230V, 60Hz, J



01-107799S  
2A Fuse, J



01-107800S  
10A Fuse, J



01-107801S  
Inlet Valve Cold Water, J



01-107803S  
Switch Overflow Chamber, J



01-107804S  
Lower Wash Arm, J



01-107806S  
Screen Drain, J



01-107807S  
Filter Drain, J



01-107808S  
Water Heater, J



01-107810S  
Switch Pressure Heater, J



01-107812S  
Clip Door, J



01-107815S  
Inlet Valve Hot Water, J



01-107938S  
Switch, Pressure Assembly  
(Chamber Full Switch Assembly)

## Hydrim C51W w/Touchscreen Parts



01-107802S  
Switch Full Chamber, J



01-107975S  
Adjustable Feet, J



01-108030S  
Quick Connect Female, J



01-108699S  
Tube, chemical



01-108121S  
Cap Quick Disconnect, J



01-108122S  
Drip Tray, J



01-108253S  
Reservoir



01-108309S  
Kickplate, J



01-108310S  
Mesh Drain, J



01-108351S  
Float Dosing Reservoir, J



01-108700S  
Fuse Holder, J



01-108795S  
Plastic Trim Edge, J



01-108797S  
Support Bracket 1-st Pump, J



01-109143S  
Dryer Tubing, J



01-109144S  
Dryer Fitting, J



01-109145S  
Dryer Vent Assy, J



01-109776S  
Top Arm Fitting Hydrim, S



01-110204S  
Tubing, Upper Arm Fitting, C51W, J

## Hydrim C51W w/Touchscreen Parts



01-109780S  
Door Solenoid C51WD, S




01-110202S  
Fascia Controller w/  
Label, C51W, J



01-109782S  
Power Supply Board, S



01-110206S  
Logic Board Hydrim, C51W, J



01-109783S  
I/O Board Hydrim, S



01-109785S  
Touch Display Hydrim, S



01-110203S  
Single Temperature  
Sensor C51W, J



01-109789S  
Rope, Door Latch, C51WD, S



01-109790S  
Upper Spray Arm  
Hydrim, S

SciCan

01-110201S  
Fascia Door w. Label, C51W, J



01-110200S  
Detergent Door w. Label  
C51W, J



01-107809S  
NC Sensor



01-110205S  
Cover Top, C51W, J



01-107973S  
Packaging Hydrim  
C51W, S




01-109834S  
Seals, Door Bottom Hydrim, S




01-109835S  
Dryer Assy Hydrim, S

## Hydrim L110W Solid Door Parts Front View Door Closed




Item #	Part#	Description
1	50750.10	Cover Top
2	13790.06	Decal Front (Set of 2)
3	13459.00	Power Switch
3	13471.01	Button ON/OFF Silver
4	16945.00	Handle Door Lock
5	16979.01	Cover Front Door Complete
6	16648.01	Kick Plate Front
7	16695.00	Cover Right Side
8	16662.06	Door Cleaning Solution
9	16663.06	Cover Slave Controller (No Decal)
10	13471.02	Button Program Silver
11	13476.25	Front Facia Complete with Decal
Not Pictured	16696.00	Cover Left Side
Not Pictured	16656.00	Cover Back
Not Pictured	13633.00	Screw Kit (Cover Screws)

## Hydrim L110W Solid Door Parts Front View Door Open




Item #	Part#	Description
1	13492.00	Wash Arm Middle
2	13456.00	Seal Peripheral Door Seal
3	13638.00	Lid for Salt Container
4	13481.00	Wash Arm Bottom
5	13611.06	Door Inner Stainless Steel
6	13466.00	Filter Bottom Flat
7	13465.00	Filter Bottom Round
8	13629.00	Pipe Water Feed to Wash Arm

## Hydrim L110W Solid Door Parts Left Side No Cover





Item #	Part#	Description
1	13462.00	Side Bag
2	13455.00	Drain Tubing
3	13614.00	Solenoid Valve Drain
4	13615.00	Solenoid Valve Regeneration
5	13447.01	Emitter Housing Red
Not Pictured	13447.00	Emitter Housing Clear

## Hydrim L110W Solid Door Parts Right Side No Cover




Item #	Part #	Description
1	16590.01	Dryer Unit Complete
2	13319.10	Dryer Air Filter
3	18466.00	Dryer Air Turbine
4	13416.01	Check Valve Dryer 1"
5	01-111568S	Flow Switch w/Bracket
6	13599.00	Foot Front
7	21624.02	Dosing Pump Clamp External
8	13634.00	Dosing Pump
9	21970.00	Dosing Pump Hose
10	21027.00	Tubing Exhaust Air


## Hydrim L110W Solid Door Parts Inlet Solenoids & Controllers



Item #	Part#	Description
1	N/A	Valve Solenoid Inlet***
<b>***Note: If Solenoid Inlet Valve is defective install Preventative Maintenance Kit #01-109707S see next page.</b>		




Item #	Part#	Description
1	16945.00	Handle Door Lock
2	13475.03	Controller Master



Item #	Part#	Description
1	13475.04	Controller Slave

# Hydrim L110W Solid Door Preventative Maintenance Kit

Part #01-109707S



Item #	Description
1	Inlet tubing's for hot and cold water
2	Tubing assembly with sensor from inlet valve block to sidebag
3	Water inlet valve block
4	Vent tubing for top of sidebag
5	Tubing's from flow switch to dosing pump and dosing pump to chamber
6	Flow switch
7	Flow switch mounting bracket
8	Wiring harness for new valve block
9	Plate for water inlet valve block
10	Screws to mount valve block plate
11	Hardware to mount flow switch
12	Tie wraps for wiring harness

## Hydrim L110W w/Glass Door Front View




Item #	Part #	Description
1	01-111469S	Cover Top
2	01-111484S	Screw 3mm Hex Head M5x16 Back Cover/Controller Door
3	01-110219S	Decal Kit
4	01-112024S	Rocker (Power) Switch
5	01-111470S	Door Cleaning Solution
6	01-111468S	Cover Right Side
7	01-111485S	Screw Kickplate M5x12 Slot Round Head
8	01-111476S	Kickplate Front
9	01-111466S	Cover Left Side
10	01-111483S	Screw 2.5mm Hex Head M4x8 Top Cover
Not Shown	01-111467S	Cover Rear
Not Shown	01-111477S	Kickplate Rear

## Hydrim L110W w/Glass Door Chamber with Trolleys




Item #	Part #	Description
1	01-111491S	Trolley Upper
2	01-109963S	5 Cassette Rack Full Size
3	01-111490S	Trolley Lower

## Hydrim L110W w/Glass Door Chamber




Item #	Part #	Description
1	01-111559S	Door Seal Kit
2	01-109790S	Wash Arm Upper & Lower Metal
3	01-111495S	Wash Arm Middle Complete
4	01-109144S	Dryer Fitting
5	13448.22	Salt Cap
6	01-111589S	Mesh Filter Wash Chamber

## Hydrim L110W w/ Glass Door Behind Front Kickplate Left Side




Item #	Part #	Description
1	22815.03	Switch Safety (behind plate)
2	01-112080S	Door Latch Cord Assembly
3	01-111702S	Sensor, Chamber Temperature
4	01-111489S	Temperature Sensor PT-100

## Behind Front Kickplate Right Side




Item #	Part #	Description
1	01-111471S	Dosing Pump
Not Shown	01-109909S	Tube Dosing Pump Internal

## Hydrim L110W w/Glass Door Control Boards




Item #	Part #	Description
1	01-109783S	I/O Board
2	01-110217S	Logic Board
3	01-109785S	Touch Display

## Hydrim L110W Right Side Top View




Item #	Part #	Description
1	01-109782S	Power Supply Board
2	01-111408S	Switch Full (Level Control)
3	01-109142S	Dryer Assembly
4	01-111474S	Flow Switch

## Hydrim L110W w/Glass Door Right Side Bottom View




Item #	Part #	Description
1	01-111474S	Flow Switch
2	01-111580S	Tension Door Spring (1/2 diameter)
2	18156.00	Door Tension Spring (large diameter)
3	01-111412S	Drain Pump
4	01-111494S	Valve Single Hot Water
5	01-111493S	Valve Double Cold Water
6	01-111473S	EMI Filter Kit
Not Shown	01-112026S	Check Valve Drain

## Hydrim L110W w/Glass Door Rear View





Item #	Part #	Description
1	01-111409S	Switch Pressure
2	01-111585S	Breather
3	01-109142S	Dryer Assembly
4	01-111492S	Drain Tubing
5	01-107788S	Water Inlet Hose
6	01-111480S	Power Cord
7	01-111465S	Circulation Pump

## Hydrim L110W w/Glass Door Top View




Item #	Part #	Description
1	01-111478S	Nut Top Cover for 2.5mm Screw
2	01-111783S	Door Latch Assembly
3	01-109782S	Power Supply Board
4	01-111408S	Switch Full (Level Control)
5	01-109142S	Dryer Assembly
6	01-111585S	Breather
7	01-111409S	Switch Pressure
8	01-112080S	Door Latch Cord Assembly

## Statmatic Parts





<b>Item #</b>	<b>Part #</b>	<b>Description</b>
1	S30153	Cover, Top with keypad
2	S30123	Keypad
3	S30141	Spray Can Black Plastic Adapter (no brass sliding lock ring)
4	S30143	Door Collet





<b>Item #</b>	<b>Part #</b>	<b>Description</b>
5	S30113	Power Inlet Assembly
6	S30112	Regulator, Air
7	S30102	Absorption pad (6 large pads & 2 bottom pads under large pad)
8	S30120	Seal Kit Front Door

## Statmatic Parts





<b>Item #</b>	<b>Part #</b>	<b>Description</b>
9	S30103	Hinge (for front door, no nuts or screws)
10	S30108	Rubber Foot





<b>Item #</b>	<b>Part #</b>	<b>Description</b>
11	S30101-CX	Coupling Multiflex/Connex (for SciCan & Kavo high speeds)
12	S30101-E	Coupling E-Type (for SciCan & Kavo low speed contra angles)
13	S30109	Chuck Adapter Complete
14	S30118	O-ring Intramatic (o-ring set for E-type coupler)
15	S30144	Clamping Nut
16	S30150	Recognition Flap
17	S30116	O-ring Set Multiflex/Connex Coupling

## Statmatic Parts




Item #	Part #	Description
18	S30115	Fuse Drawer
19	S30117	Fuse 1A





Item #	Part #	Description
20	S30107	Comp. Air Cylinder (air solenoid for solution bottle)
21	S30105	Solenoid Valve Block
22	S30130	Power Cable (from power inlet to power supply board)
23	S30129	Cable PE (ground wire for power inlet assembly)
24	S30132	Tubing Blue (2" piece of blue air line tubing)
25	S30121	Control Board
26	S30119	Power Supply Board
27	S30151	Coupler Receptacle with 5 Hole Configuration
27	S30140	Connection Bush (Coupler receptacle with 4 hole configuration)
28	S30106	Microswitch Assembly (3 Switches with Wiring)
29	S30125	Tubing (clear 10" tubing)
30	S30110	Hose Set (complete tubing assembly including "Y" fittings)
31	S30139	Hex Nut 12 X 1 (hex-nut which holds the coupler receptacles)

## Statmatic Parts



<b>Item #</b>	<b>Part #</b>	<b>Description</b>
32	S30107	Air Solenoid for Solution Bottle
33	S30127	Tubing Labeled (4" yellow tubing labeled "oil")
34	S30145	Sensor Collet Complete (with oil tubes, optical sensor & wiring)
34	S30148	Sensor Collet (with oil tubes & brass elbow fitting)
35	S30104	Fluid Sensor (optical fluid sensor)
36	S30137	Clamp Bugel (wire fastener the holds cable to fluid sensor)
37	S30138	Cable Photo Sensor (cable from fluid sensor to main PCB)
38	S30126	Tubing Labeled (10" yellow tubing labeled "oil")
39	S30133	Tubing Switch (door magnet sensor and cable)
40	S30131	Cable SNT/Elect. (cable from power supply board to main PCB)

## Statmatic Parts



Item #	Part #	Description
41	S30124	Spring (for microswitch activating)
42	S30142	Guide Set (plunger assembly for solution bottle)

## Statmatic Couplers



**S30101-BA**  
BienAir



**S30101-CX**  
Multiflex/Connex



**S30101-E**  
E-type



**S30101-KVH**  
Kavo Intra head



**S30101-MW**  
Midwest Slow speed



**S30101-NSK**  
Coupling NSK



**S30101-ST1I**  
Sirona T1 Slow speed  
Attachments




**S30101-ST1T**  
Sirona T1 High speed




**S30101-W&H**  
Coupling W&H

## Aquastat Boiler




Item #	Part #	Description
1	P492-2A	Boiling Chamber Cap (No Gasket)
2	P167-72	Boiling Chamber Cap Gasket
3	KP05EC001	Cap for Baffle
4	KP03WS010	Baffle without O-ring
Not Shown	KP32GA162	Gasket for Baffle
Not Shown	KP32GA163	Baffle O-ring
5	P186-1X	Boiler Vessel with Heater Only
6	KP28TS091	Thermostat
6	KP28TS098	Thermostat Round (not shown)
7	KP64CO001U	Control Board
8	P338-26	Strix Connector Boiler
9	KP44WP005	Boiler Base Bottom with Grey Feet
10	P128-134	Spring
11	P81-174	Button
12	P118-96	Retainer
13	P19-283	Rubber Feet Boiler Bottom Base

## Aquastat Chassis



Item #	Part #	Description
17	P8-34	Vent Piece for Top Cover
18	P72-8	Top Piece
19	KP76PL001	Fan Blade Clear
20	KP63MT019	Fan Motor 110V
21	P450-42	Clip Thermostat
22	P185-97	Thermostat Condensing Coil
23	P3-192	Elbow
24	P26-93	Hose without Vent Hole
25	P84-126	Condensing Coil
26	P26-93D	Hose with Vent Hole
28	KP03WS013	Spout
30	P338-26	Strix Connector Boiler
30	P314-26	Otter Connector Boiler
31	P16-78	Strix Connector Chassis
31	P13-78	Otter Connector Chassis
33	P118-14	Power Socket Main Frame

## Aquastat Reservoir



Item #	Part #	Description
35	P160-89M	Reservoir Cap
36	W50211	Water Jug (clear w/spigot)
37	P96-104	Faucet for Reservoir
38	W8010	Carbon Filters (6 Per Pack)(see note)

**Note: Carbon filters should be used for making drinking water only.**

## Aquastat Parts



**Baffle Gasket**  
**#KP32GA162**



**Baffle Nut**  
**#KP05EC001**



**Baffle O-ring**  
**#KP32GA163**



**Baffle without O-ring**  
**#KP03WS010**



**Boiler Base Bottom**  
**# KP44WP005**



**Boiling Chamber Lid**  
**(no gasket)**  
**#P492-2A**



**Boiling Chamber New  
Style  
Oval Connector**  
**#P206-1H**



**Boiling Vessel  
with Heater Only**  
**#P186-1X**



**Button**  
**#P81-174**



**Carbon Filters 6 per Pack**  
**#W8010**



**Chassis Kit**  
**#P35-197**




**Clip Thermostat**  
**#P450-42**

## Aquastat Parts



**Condensing Coil**  
**#P84-126**



**Control Board**  
**#KP64CO001U**



**Cover Top without  
Vent Piece**  
**#P72-8**



**Distiller Cleaner &  
Descaler**  
**#D1-500200S**



**Distiller Wiring Kit**  
**#P376-90K**



**Elbow**  
**#P3-192**



**Fan Blade Clear**  
**#KP76PL001**




**Faucet for Reservoir**  
**#P96-104**



**Feet Gray**  
**#P19-283**



**Gasket for Boiling  
Chamber Lid**  
**#P167-72**



**Hose with Vent Hole**  
**#P26-93D**



**Hose without Vent Hole**  
**#P26-93**

## Aquastat Parts



**Motor, 110V**  
**#KP63MT019**



**Nut**  
**#P340-81**



**Otter Connector Boiler**  
**#P314-26**



**Otter Connector Chassis**  
**#P13-78**



**Owners Manual**  
**#L5446**



**Packaging**  
**#J3366K**



**Power Cord**  
**#P313-74**



**Power Socket Main  
Frame**  
**#P118-14**



**Replacement Reservoir**  
**#W50211**



**Reservoir Cap**  
**#P160-89M**



**Retainer**  
**#P118-96**



**Screw**  
**#KP13MS570**

## Aquastat Parts



**Spacer**  
**#KP12SN254**



**Spout**  
**#KP03WS013**




**Spring**  
**#P128-134**



**Strix Connector Boiler**  
**#P338-26**



**Strix Connector Chassis**  
**#P16-78**



**Thermostat Otter**  
**#KP28TS091**



**Thermostat Strix**  
**#P197-97**




**Thermostat Condensing  
Coil**  
**#P185-97**



**Vent Piece for Top Cover**  
**#P8-34**

## VistaBrite Part Identification




Name	Part #	Description
Supply Hose Connectors	S4231	Garden hose adaptor – Female hose x ¼" tube
<b>VistaCheck™ Backflow Preventer</b>	VC250	Backflow check valve - In-line ¼" FNPT
System Check Valve	S6265P	In-line check valve – ¼" tube
System Regulator	S9130	Automatic shut-off valve - Hydraulic
Temperature Gauge	S6435	Temp gauge – ¼" MNPT back mount (0-250 F/C)
Filter Sump Brackets	S9170	Filter mounting bracket
Filter Sump Heads	S9171	Filter head
Inlet Pressure Gauge	S6431	Pressure gauge – Back mount ¼" MNPT (0-100 PSI)
Storage Tank Pressure Gauge	S6431	Pressure gauge – Back mount ¼" MNPT (0-100 PSI)
System Blending Manifold	S4310	System manifold ONLY (does not include fittings)
<b>Carbon Prefilter</b>	<b>R4151</b>	<b>Replacement carbon prefilter (every 6 months)</b>
<b>Hyper Prefilter</b>	<b>R4155</b>	<b>Replacement hyperfiltration filter (approx 2 yrs)*</b>
Drain Line Flow Restrictor	S9165	Waste water flow restrictor – ¼" tube (1500 ml/min)
Sump Drain Elbow	S6121P	Male adaptor elbow – 1/8" NPT x ¼" tube
Drain Saddle	S9175	Drain saddle – ¼" tube
Valved Hose Wye	S9180	Garden hose valved wye
Male-Male Hose Adaptor	S9181	Garden hose male to male adaptor
Female Hose Adaptor	S4230	Garden hose adaptor – Female hose x 3/8" tube
Storage Tank (Standard 7 GAL)	S9111	Pressure storage tank w/ stand – 7 gallon nominal
Storage Tank (Optional 14 GAL)	S9113	Pressure storage tank w/ stand – 14 gallon nominal
Tank Valve Assembly	*****	(see top-down view on following page)

\* Hyperfiltration filter should be changed every 2 years, or when TDS reduction is less than 70%

## VistaBrite Part Identification




(Storage Tank Top-Down View)



Name	Part #	Description
Storage Tank (Standard 7 GAL)	S9111	Pressure storage tank w/ stand – 7 gallon nominal
Storage Tank (Optional 14 GAL)	S9113	Pressure storage tank w/ stand – 14 gallon nominal
Tank Valve	S6160P	Elbow shut-off valve – 1/4" tube x 1/4" FNPT
1/4" x 3/8" Tube Union	S6310P	Union – 1/4" tube x 3/8" tube
1/4" Cross	S6175P	Cross – 1/1" tube
Test Cock Valve	S6155P	Valve – In-line – 1/4" tube

## VistaBrite Part Identification

(continued)



Item #	Part Number	Part Description
1	S6112P	Swivel Plug-In Elbow – ¼" tube x ¼" stem
2	S6110P	Elbow Union – ¼" tube x ¼" tube
3	S9130	Automatic Shut-Off Valve – Hydraulic / (Regulator)
4	S6145P	Valve – Straight Male – ¼" MNPT x ¼" tube
5	S6570	Wall Mount "C-Clip" – 2"
6	S6122	Male Adapter Straight – ¼" NPT x ¼" tube
7	S6120P	Male Adapter Elbow – ¼" NPT x ¼" tube
8	S9172	Filter Head Self-Tapping Screws
9	S6171P	Male Adapter Elbow – 3/8" NPT x ¼" tube
10	S6123P	Male Adapter Straight – 3/8" NPT x ¼" tube
11	S6170P	Male Adapter Straight – 3/8" NPT x 3/8" tube