

2016

KaVo COLLECTIONS INSTRUMENTS & SMALL EQUIPMENT CATALOG



KaVo. Dental Excellence.

KaVo INNOVATION

Over a Century of Dental Excellence

Innovation and quality are traditions that have kept KaVo at the forefront of the dental industry for more than 100 years. The company was founded in 1909 by Alois Kaltenbach in Berlin-Steglitz. It became known as Kaltenbach & Voigt (KaVo) when Richard Voigt joined the company in 1919 and is now a trusted partner for dentists and dental technicians worldwide.

KaVo's continued commitment to advancing dental technology for the betterment of your practice drives us every day. It is innovation and our passion for quality that has fueled our steady growth and world-wide recognition. You can see the outcome in our revolutionary new products, such as the COMFORTdrive™, MASTERtorque™, MASTERmatic™, MASTERsurg™, and innovative programs such as TryKaVo.com. We believe these products, and our entire portfolio, are the best in the industry and we encourage you to experience KaVo for yourself.






In 1909 KaVo is founded by Alois Kaltenbach and in 1919 Richard Voigt joins the company.



In 1928, the first sterilizable handpiece is created by KaVo.



In 1958 KaVo launches its first turbine called Borden-Air 60.



In 1977 the KaVo MULTIflex™ is the first quick-release coupling.



In 1982 KaVo integrates the first light conductors into turbines.



In 2000 KaVo introduces the revolutionary DIAGNOdent™, a caries detection aid.



In 2005, KaVo introduces GENTLEsilence™ LUX 6500 B – one of the quietest turbines in the world.



In 2008 KaVo introduces the GENTLEmini™ LUX 4500 B, a mini head, mighty power handpiece. Also launched is the GENTLEray™ 980 Premium Soft Tissue Diode Laser, the first diode laser with micro-pulsing power.



In 2009 KaVo celebrates 100 years of serving the dental industry, joining Pelton & Crane with 110 years of providing dental excellence. KaVo introduces the COMFORTdrive - KaVo's most innovative High-speed.



In 2012 KaVo introduces INTRAmatic LUX™ KL703 Short Motor, one of the shortest and lightest motors.



In 2014 KaVo introduces the MASTERtorque™ M8900 L, one of world's most powerful and advanced High-speed turbines. KaVo also introduces the new surgical portfolio featuring the MASTERsurg and EXPERTsurg – simple and intuitive surgical motors, redefining the best.



In 2015, KaVo introduces the MASTERtorque M8700 L with the advanced technology of the M8900 L in a mini head. Also launched is the MASTERmatic™ Series, featuring significant improvements in head size without compromising power and performance.





Electric Collections

Electric Collections Overview	6-9
Master Series - COMFORTdrive	10
Master Series - ELECTROtorque™ TLC System.	11
Master Series - MASTERmatic™ Attachments	12
Master Series - INTRA LUX™ Heads	14-15
Expert Series - ELECTROtorque plus System	18-19
Expert Series - EXPERTmatic™ Attachments	19-20
Smart Series - INTRAmatic™ ES Systems	21

Air-Driven Collections

Air Driven Collections Overview	22-25
Master Series - MASTERtorque	26
Expert Series - EXPERTtorque™	27
Smart Series - SMARTtorque™	28-29
INTRAFlex™ Air Motor System	30
MULTIflex Couplers	31
INTRAmatic E Air Systems	32
Low-speed Air Motors	33

KaVo Parts

Turbine Cartridges	37
Turbince Back Cups And Tools.	38
Accessories for High-speed Handpieces.	39

Bundle Program

KaVo Bundle Programs.	40-43
-------------------------------	-------

Maintenance Systems

QUATTROcare™ Plus	44-48
Spray / Accessories.	49

Hygiene

Hygiene 360° Overview	50-51
KaVo SONICflex™ System	52
SONICflex Scalers	53
SONICflex - Hygiene & Periodontology	54-55
SONICflex - Endodontics	56
SONICflex - Minimally invasive caries therapy	57-58
SONICflex - Inlay/Onlay preparation	59
PROPHYflex System and Accessories	62
Hygiene Polishing Instruments	63
DIAGNODent Caries Detection.	64-65
RONDOflex plus	66
Accessories for RONDOflex plus.	67
CORONAFlex 2005	68

Surgical Systems

Surgical Systems Overview	70
MASTERsurg™ LUX Wireless.	71
EXPERTsurg™ LUX	71
Surgical Handpieces and Accessories	72-74

Soft Tissue Laser Accessories

GENTLEray 980 Premium Accessories	75
---	----

In-Office Laboratory

In-Office Laboratory	76
--------------------------------	----

Illuminating Systems

Illuminating Systems.	77
-------------------------------	----

Succession Chart

KaVo Instruments and Their Successors.	78
--	----

Product Index

Product Index	79
-------------------------	----

Maintenance Information

Maintenance Information	80
-----------------------------------	----

General Information

Terms of Sale / Registered Trademarks	81
Warranty	82
KaVo's Industry Award Recognition.	83

Legend



NOTE: The prices listed in the catalog are retail prices effective January 1st, 2016. Prices are subject to change without notice. Possession of list price does not constitute an offer to sell.




ELECTRIC COLLECTIONS



Master Series

The Pinnacle of KaVo Innovation.

INTRA LUX Heads



Color Key
● Blue = 1:1 ● Green = Reduction ● Yellow = Specialty

MASTERmatic Attachments



ELECTROtorque Systems



ELECTROtorque TLC Torque Limit Control LUX with KL703 Brushless Electric Motor Kit
1003.5840



COMFORTdrive Full Operatory Package
1007.4682

Expert Series

Superior KaVo Performance.

EXPERTmatic Attachments



Non-Optic



ELECTROtorque Systems




ELECTROtorque plus LUX with
KL703 Brushless Electric Motor Kit
1002.5015

To be used with
air or electric motor
with "E"
type connection

Smart Series


Reliable and Affordable KaVo Technology.

INTRAmatic™ Attachments



To be used with
any air or electric
motor with "E"
type connection
(181L, 181H,
ELECTROtorque plus,
ELECTROtorque TLC)






High-speed

Master Series

The Pinnacle of
KaVo Innovation.

Featuring the new
MASTERmatic LUX M25 L
and MASTERmatic LUX
M05 L Mini High-speeds
with KaVo's smallest head
size yet on an electric
attachment. Also featuring
the COMFORTdrive 200 XDR
– KaVo's most innovative
High-speed ever.



High-speed

Expert Series

Superior KaVo Performance.

Featuring the new
EXPERTmatic E25L High-
speed. Fast and precise
working thanks to an
even smaller head, snap-
on facility and ceramic
bearings.





High-speed

Smart Series

Reliable and Affordable
KaVo Technology.

The KaVo Smart Series High-speed instruments are redesigned with a new surface finish for a perfect grip and no rough edges. Seamless stainless steel profile means there's nowhere for debris to accumulate and the push-button chuck is standard.



Low-speed

Smart Series

Reliable and Affordable
KaVo Technology.

The KaVo Smart Series Low-speed instruments can be thermally disinfected and sterilized in autoclaves at temperatures of up to 275°F and can also be washed in the thermodisinfector. Durable hygienic coating with textured surface for easy grip.



MASTER SERIES

The Pinnacle of KaVo Innovation



COMFORTdrive 200 XDR Handpiece

KaVo's most innovative High-speed handpiece features precise performance, superior comfort, and allows for seamless integration into dental units. Speed range 30,000 to 200,000 rpm. Integrated micro-motor with 30W power, cellular optic, ceramic ball bearings, and exchangeable micro-filter. Only connects to the KaVo COMFORTbase coupling.

Part No. 1007.3570

\$1,699.00



COMFORTdrive System



The COMFORTdrive system consists of the COMFORTdrive 200 XDR handpiece, the COMFORTtronic control box and the COMFORTbase coupling. The 200 XDR High-speed features a 2½ year warranty, to guarantee worry-free performance when maintained with KaVo QUATTROcare. Standard warranty is 2 years. The COMFORTbase and the COMFORTtronic carry a 3 year warranty.

Part No. 1007.3656

\$3,499.00

COMFORTdrive Full Operatory Package



The Full Operatory Package includes three COMFORTdrive 200 XDR handpieces, one COMFORTtronic control box and one COMFORTbase coupling. The 200 XDR High-speed features a 2½ year warranty to guarantee worry-free performance when maintained with KaVo QUATTROcare. Standard warranty is 2 years. The COMFORTbase and the COMFORTtronic carry a 3 year warranty.

Part No. 1007.4682

\$5,699.00

Parts and Accessories for COMFORTdrive

Part No.	Price	Description
1005.1707	\$99.00	COMFORTdrive Adapter for maintenance of your 200 XDR handpiece using KaVo's QUATTROcare Plus handpiece maintenance system
1005.3154	\$22.00	COMFORTdrive Nozzle for maintenance of your 200 XDR handpiece using KaVo spray
1005.1773	\$16.00	Lamp / Bulb change tool
1007.5372	\$125.00	KaVo MULTI LED Bulb
1002.2928	\$105.00	Halogen Light Bulb (pack of 2)
1005.0327	\$4.00	COMFORTdrive O-Rings (set of 3)
1002.0271	\$23.00	Water filter
1002.0321	\$2.75	Filter tool
1007.9742	\$91.00	External Water Filter
1007.9736	\$14.00	Water Filter Element
1005.0169	\$2,525.00	COMFORTtronic Unit
1005.0120	\$399.00	COMFORTtronic Transformer
1005.0008	\$429.00	COMFORTtronic Tubing
1006.5370	\$14.00	COMFORTtronic Rinsing Cap
3001.9302	\$50.00	COMFORTbase Cleaning Cap



ELECTROtorque TLC Torque Limit Control LUX Brushless Motor Kit

Air/Electric conversion kit for any delivery unit. To be used with KaVo MASTERmatic LUX or other "E" type attachments.



townie
choice
awards 2011



System Advantages:

- Brushless electric motor for ALL your High-speed, Low-speed and torque control endodontic needs.
- No more searching and set up of a separate endo motor and no more wasted counter space – endo operation is available with the push of a button.
- Performs exact rotary endo procedures safely with torque control, auto-reverse and cellular optic advantages.
- Remote mountable faceplate for convenient mounting options and easy to access controls.
- KL703 is one of the world's shortest and lightest brushless electric motors (no brushes to exchange means less maintenance and downtime).
- Motor speed range 100-40,000 rpm (with various attachments 2-200,000 rpm).
- Endodontic torque 0.4 - 8.0 Ncm. Motor torque 3 Ncm.
- 3-year warranty (includes motor, tubing and control box).

Includes:

- ELECTROtorque TLC module with 6.5' cable (1003.3490)
- ELECTROtorque TLC transformer with connector (1003.8542)
- Electric tubing (1008.2671)
- KL703 brushless electric motor (sterilizable)
- Mounting bracket with hardware (1002.0622)

Technical:

- W: 5-1/8", D: 8-1/2", H: 3-1/2"
- Motor weight: 0.30 lb. or 0.14 kg
- Total weight: 2.30 lb. of 1.04 kg

- Electric motor for surgical procedures refer to Surgical Systems on page 70

Also Available:

- Introduction manual (1003.5925)
- Quick reference guide (1003.8440)
- Service manual (1003.8441)

Part No. 1003.5840 \$3,489.00

Parts and Accessories for ELECTROtorque TLC

Part No.	Price	Description
1003.5710	\$529.00	ELECTROtorque TLC remote mounting face plate connection kit
1003.5707	\$53.00	Remote Cable for Face Plate
1003.8537	\$63.00	3' Cable for ELECTROtorque TLC LUX Brushless Motor Kit
1008.2671	\$449.00	Tubing
1003.8542	\$365.00	Transformer 4881
1001.7115	\$7.25	Fuse
0200.6120	\$17.25	Motor O-ring
1003.5822	\$2.50	Large O-ring (base of motor)
1003.5656	\$6.75	Large O-ring (under bulb sheath)
1007.0150	\$2,479.00	INTRAmatic LUX KL703
1007.0146	\$26.00	KL703 Motor Bulb Ring
1007.9742	\$91.00	External Water Filter
1007.9736	\$14.00	Water Filter Element
1007.8474	\$116.00	KL703 motor light bulb

Note: KL703 will not work with ET TLC Serial # SN13651 and earlier.

Master Series | MASTERmatic Attachments



MASTERmatic LUX M25 L High-speed Attachment MASTERmatic LUX M05 L Mini High-speed Attachment

- New smaller head size
- Transmission 1:5
- Triple gear transmission for reduced vibration
- Ceramic bearings
- Speed range up to 200,000 rpm
- Cellular optic
- Internal air/water coolant
- Friction grip push button auto chuck
- 3-port spray M25 L (1-port spray M05 L)
- Now with field exchangeable spray filter for consistent spray performance
- Works on electric or air motors that have a universal "E" type connection
- Plasmatec hygienic coating
- Improved angulation for easy access
- 100° head angle for superior access

Part No. 1009.3630	\$1,669.00	M25 L
Part No. 1009.3640	\$1,669.00	M05 L Mini
Part No. 1002.0271	\$23.00	Water filter
Part No. 1002.0321	\$2.75	Filter tool



MASTERmatic LUX M20 L Contra Angle

- Transmission 1:1
- Speed range up to 40,000 rpm
- Cellular optic
- Internal air/water coolant
- Now with field exchangeable spray filter for consistent spray performance
- Works on electric or air motors that have a universal "E" type connection
- Plasmatec hygienic coating
- Improved angulation

Part No. 1009.3620	\$960.00	M20 L
Part No. 1002.0271	\$23.00	Water filter
Part No. 1002.0321	\$2.75	Filter tool



MASTERmatic LUX M07 L Contra Angle

- Reduction 2.7:1
- Speed range up to 15,000 rpm
- Cellular optic
- Internal air/water coolant
- Now with field exchangeable spray filter for consistent spray performance
- Works on electric or air motors that have a universal "E" type connection
- Plasmatec hygienic coating
- Improved angulation

Part No. 1009.3610	\$1,175.00	M07 L
Part No. 1002.0271	\$23.00	Water filter
Part No. 1002.0321	\$2.75	Filter tool



MASTERmatic LUX M29 L Contra Angle

- Reduction 7.4:1
- Speed range up to 5,500 rpm
- Cellular optic
- Internal air/water coolant
- Now with field exchangeable spray filter for consistent spray performance
- Works on electric or air motors that have a universal "E" type connection
- Plasmatec hygienic coating
- Improved angulation

Part No. 1009.3580	\$1,319.00	M29 L
Part No. 1002.0271	\$23.00	Water filter
Part No. 1002.0321	\$2.75	Filter tool



MASTERmatic LUX M10 L Straight Handpiece

- Transmission 1:1
- Speed range up to 40,000 rpm
- Works on electric or air motors that have a universal "E" type connection
- Plasmatec hygienic coating
- Cellular optic
- Internal air/water coolant
- For straight handpiece burs 2.334 to 2.35mm and contra angle burs 2.35mm (burstop must be used, and is included with purchase)

Part No. 1009.3570	\$1,385.00	M10 L
Part No. 0524.0892	\$21.00	Burstop

*2.5 years when maintained with QUATTROcare Plus






NEW



2.5
YEAR
WARRANTY



L68 B Latch Head – Operative

- Transmission 1:1
- For latch burs
- Push button auto chuck
- Replaces 68LU

Part No. 1008.1834

\$595.00



2.5
YEAR
WARRANTY



L52 Mini Latch Head

- Transmission 1:1
- Ideal for Mini Latch Burs (16mm) available from Kerr Dental**
- Push button Auto chuck
- Ideal for treatment in narrow areas or when treating children

Part No. 1008.1825

\$595.00



**Kerr Dental 800-537-7233



2.5
YEAR
WARRANTY



L80 Friction Grip Head – Operative

- Transmission 1:1
- For friction grip burs
- Push button auto chuck
- Replaces 80LD

Part No. 1004.4456

\$595.00



2.5
YEAR
WARRANTY



L67 Reduction Latch Head – Operative

- Reduction 2:1
- Push button auto chuck
- For latch burs
- Replaces 67LH

Part No. 1008.1832


\$615.00



NEW



2.5
YEAR
WARRANTY



L66 B Reduction Latch Head – Endodontics

- Reduction 3:1
- Push button auto chuck
- For latch type endo files or latch burs
- 360° rotation
- Recommended use with M20 L and ELECTROtorque TLC Torque Limit Control
- Replaces 66LU

Part No. 1008.1831

\$615.00



2.5
YEAR
WARRANTY



L62 Reduction Latch Head – Endodontics

- Reduction 10:1
- Push button auto chuck
- For latch type endo files or latch burs
- 360° rotation
- Replaces 62LDN

Part No. 1008.1830

\$615.00



2.5
YEAR
WARRANTY



L53 Oscillating Head – Endodontics

- The drive speed is converted into an oscillating movement about 90°
- Used with latch type endo files
- Max Speed 3,000 rpm, recommended use with M29 L shank
- Now also for polishing attachments (2.35mm diameter)

Part No. 1008.1826

\$575.00

NEW**L3 Reduction Head - Endodontics**

- Speed max 3,000 rpm
- The drive speed is converted into an oscillating movement about 90°
- To be used with hand NiTi files
- Push-button chuck
- Use with M29 L shank
- Replaces 3LD

Part No. 1008.1822

\$575.00

NEW**L3 Y Reduction Head - Endodontics**

- Speed max. 3,000 rpm
- The drive speed is converted into a stroke movement of about 0.4 mm
- To be used with hand NiTi files
- Push-button chuck
- Use with M29 L shank
- Replaces 3LDY

Part No. 1008.1823

\$575.00

NEW**L61 G Orthodontics / Periodontics / Endodontics Head**

- To be used with Intensiv files (Proxoshape) available through Kerr Dental**
- Free rotation
- For interproximal cleaning, polishing and preparation
- Buffered stroke (protects soft tissue)
- Elevation 0.8 mm
- Includes file ejector tool (0410.0634)
- Recommended use with M07 L (on electric motor); 20LP/CH, M20 L or INTRAFlex LUX 2307LN
- Replaces 61LG

Part No. 1008.1828

\$575.00

Part No. 0410.0634

\$19.00 Ejector Tool for L61 G and L61 R

NEW**L61 R Orthodontics / Periodontics / Endodontics Head**

- To be used with Intensiv files available through Kerr Dental**
- 36 individual locking positions
- For interproximal polishing, stripping and preparation
- 10,000 rpm optimal drive speed
- 20,000 rpm maximum drive speed
- Buffered stroke (protects soft tissue)
- Elevation 0.8mm
- Includes file ejector tool (0410.0634)
- Recommended use with M07 L (on electric motor); 20LP/CH, M20 L or INTRAFlex LUX 2307LN
- Replaces 61LR

Part No. 1008.1829

\$575.00

Part No. 0410.0634

\$19.00 Ejector Tool for L61 G and L61 R

NEW**L31 Reduction Prophy Head**

- Reduction 2:1
- For snap on or screw in cups or brushes
- Replaces 31LR

Part No. 1008.1824

\$320.00

Part No. 0549.0602

\$23.00 Screw for L31

NEW**L22 Micro Head - Surgical**

- Transmission ratio 1:1
- With round drill size 2, also fits KaVo round
- Drill size Ø or U-inverted cone burs sizes 2 and 0
- Microhead with extremely small dimensions for microsurgical procedures such as apicectomies (through bone window)
- Maximum operating speed 20,000 rpm

Part No. 1008.1835

\$575.00

Part No. 1003.1523

\$162.00 Special Screwdriver



KaVo Chuck Strength

KaVo handpieces feature a patented chuck system with exceptional chuck strength. The chuck strength will deliver bur concentricity for precise cutting and safety.

EXPERT SERIES

Superior KaVo Performance



ELECTROtorque plus LUX Motor Kit

Air/Electric conversion kit for any delivery unit. To be used with KaVo MASTERmatic LUX, EXPERTmatic or other "E" type attachments.

System Advantages:

Superior power & torque. Unique brushless motor (no brushes to exchange - no downtime). Quieter than brush motors. 3-year warranty includes motor, tubing and control box. Two programmable speed settings. Water spray can be adjusted at motor. Motor speed range 100–40,000 rpm (with various attachments 2–200,000 rpm). Motor torque 3 Ncm. Integrated light and coolant. Forward/reverse speed adjustment. Small, easy to install control box. Unit can be calibrated to existing air pressure for optimized performance. Motor is autoclavable. Includes KL703 – one of the shortest and lightest electric motors with exclusive SAFEdrive and MULTI LED Light. Optional handpiece stand (0411.9501).

Technical:

- Width: 5" Depth: 6" Height: 2"
- Motor weight: 99 grams, 2.97 oz
- Total weight: 29oz or 880g

- Electric motor for surgical procedures refer to Surgical Systems on page 70

Part No. 1002.5015 \$2,629.00

Includes:

- ELECTROtorque plus module with cable
- ELECTROtorque plus transformer with connector (1003.8542)
- Electric tubing (1008.2671)
- Mounting bracket with hardware (1001.5113)
- KL703 Short Motor (1007.0150)

Note: KL703 will not work with ET Plus Serial # SN18796 and earlier.

Parts and Accessories for ELECTROtorque plus

Part No.	Price	Description
1007.0150	\$2,479.00	INTRAmatic LUX KL703
1001.9951	\$216.00	Flow meter
1008.2671	\$449.00	Tubing for ELECTROtorque plus box
1001.5113	\$45.00	Mounting bracket with hardware
1001.8748	\$2,795.00	ELECTROtorque plus control unit
0692.6891	\$56.00	Transformer power cord
0200.6120	\$17.25	Motor O-ring
0200.6323	\$2.75	Large O-ring (base of motor)
0200.6308	\$7.25	Large O-ring (under bulb sheath)
1003.8542	\$365.00	Transformer 4881
1007.0146	\$26.00	KL703 Motor Bulb Ring
1007.8474	\$116.00	KL703 Motor Light Bulb
1007.9742	\$91.00	External Water Filter
1007.9736	\$14.00	Water Filter Element
1001.9081	\$535.00	Tubing for KL 700 and 701 motor series

Parts and Accessories for ELECTROtorque

Part No.	Price	Description
0905.8095	\$215.00	Remote panel kit. Includes 6' remote cable (0906.0753) & remote panel (0906.0754)
0201.8135	\$5.25	Water valve spring
0200.6120	\$17.25	Motor O-ring, 4 required for motors predating manufacturing date of 2/2008. After that date only 3 required.
0200.6042	\$5.25	Water valve small O-ring
0200.6038	\$5.25	Water valve large O-ring
0676.0901	\$19.00	Water valve ball bearing
0535.1661	\$53.00	Bulb holder
1002.2928	\$105.00	Halogen Light Bulb (pack of 2)

Parts and Accessories for ELECTROtorque with K200 / 196 Motors

Part No.	Price	Description
1000.6717	\$2,447.00	K200 INTRA LUX Motor
0516.5850	\$840.00	Electric tubing for 196/200 motor. 6.8' length. Use with ELECTROtorque
0200.6120	\$17.25	Motor O-ring. 4 required for motors predating manufacturing date of 2/2008. After that date only 3 required
0676.1852	\$85.00	Sleeve for 196 Motor
0676.1922	\$62.00	196 motor water control ring
0905.3035	\$108.00	Carbon brush kit for 196 Motor. Includes 2 brushes (0676.0881) & 1 tool (0411.0041)
1002.2928	\$105.00	Halogen Light Bulb (pack of 2) for 196, K200, KL700, KL701, K702 Motors

Expert Series - EXPERTmatic Attachments



EXPERTmatic E25L

- Speed Increasing 1:5
- Ceramic bearings
- Speed range up to 200,000 rpm
- Cellular optic
- Internal air/water coolant, Friction grip push button auto-chuck
- 3-port spray, EXPERTgrip, Universal "E" type connection

Part No. 1007.5550 \$1,429.00
 Part No. 1002.0271 \$23.00 Water filter
 Part No. 1002.0321 \$2.75 Filter tool



EXPERTmatic E25C

- Speed Increasing 1:5
- Ceramic bearings
- Speed range up to 200,000 rpm
- Internal air/water coolant, Friction grip push button auto-chuck
- 3-port spray, EXPERTgrip, Universal "E" type connection

Part No. 1007.5551 \$1,225.00
 Part No. 1002.0271 \$23.00 Water filter
 Part No. 1002.0321 \$2.75 Filter tool



EXPERTmatic E20L

- Transmission 1:1
- Cellular optic
- For latch type burs
- Speed range up to 40,000 rpm (electric motor)
- Speed range up to 20,000 rpm (air motor)
- Internal air/water coolant
- Push button auto chuck
- 1-port spray
- EXPERTgrip
- Universal "E" type connection

Part No. 1007.5540 \$1,359.00



EXPERTmatic E20C

- Transmission 1:1
- For latch type burs
- Speed range up to 40,000 rpm (electric motor)
- Speed range up to 20,000 rpm (air motor)
- Internal air/water coolant
- Push button auto chuck
- 1-port spray
- EXPERTgrip
- Universal "E" type connection

Part No. 1007.5541 \$969.00



townie
choice
awards 2015

EXPERTmatic E15L

- Reduction 5.4:1
- Cellular optic
- For latch type burs
- Speed range up to 7,400 rpm (electric motor)
- Speed range up to 3,700 rpm (air motor)
- Internal air/water coolant
- Push button auto chuck
- 1-port spray
- EXPERTgrip
- Universal "E" type connection

Part No. 1007.5530 \$1,185.00

townie
choice
awards 2015

EXPERTmatic E15C

- Reduction 5.4:1
- For latch type burs
- Speed range up to 7,400 rpm (electric motor)
- Speed range up to 3,700 rpm (air motor)
- Internal air/water coolant
- Push button auto chuck
- 1-port spray
- EXPERTgrip
- Universal "E" type connection

Part No. 1007.5531 \$969.00

townie
choice
awards 2015

EXPERTmatic E61C

- EVA
- 10,000 rpm optimal drive speed (20,000 rpm max. drive speed)
- To be used with Intensiv files (ORTHOstrips or Proxoshapes)
- 36 individual locking positions
- For interproximal polishing
- Stripping and preparation
- Buffered stroke (protects soft tissue)
- Elevation 0.8 mm
- Internal air/water coolant
- 1-port spray
- EXPERTgrip
- Universal "E" type connection

Part No. 1008.4999 \$1,035.00
Part No. 0410.0634 \$19.00 Ejector Tool for E61Ctownie
choice
awards 2015

EXPERTmatic E31C

- Reduction 5.4:1
- Ideal speed less than 5000 rpm
- Speed range up to 7,400 rpm (electric motor)
- Speed range up to 3,700 rpm (air motor)
- EXPERTgrip
- Universal "E" type connection

Part No. 1008.4998 \$1,125.00

townie
choice
awards 2015

EXPERTmatic E10C

- Transmission 1:1
- For straight handpiece burs
- Speed range up to 40,000 rpm (electric motor)
- Speed range up to 20,000 rpm (air motor)
- Internal air/water coolant
- 1-port spray
- EXPERTgrip
- Universal "E" type connection

Part No. 1007.5560 \$819.00

SMART SERIES

Reliable and Affordable KaVo Technology

Smart Series - INTRAmatic Attachments



23ES Contra Angle Attachment

- Transmission 1:3
- Speed range up to 132,000 rpm
- 1-port spray
- Push button
- 100° head angle
- For friction grip burs
- For use with any "E" style connection
- Air or Electric motor
- Durable hygienic coating with textured surface for easy grip

Part No. 1004.5380 \$999.00

Part No. 0065.5188 \$19.00 External waterline

80ES Contra Angle Attachment

- Transmission 1:1
- Speed range up to 40,000 rpm
- 1-port spray
- Push button
- 100° head angle
- For friction grip burs
- For use with any "E" style connection
- Air or Electric motor
- Durable hygienic coating with textured surface for easy grip

Part No. 1004.4375 \$619.00

Part No. 0065.5188 \$19.00 External waterline

19ES Hygiene Handpiece

- Reduction 7.4:1
- Speed range up to 5,500 rpm
- For disposable prophy heads only
- For use with any "E" style connection Air or Electric Motor
- Durable hygienic coating with textured surface for easy grip

Part No. 1004.4503 \$519.00

Part No. 0201.0116 \$21.00 Screw for Prophy Angle



AIR COLLECTIONS



Master Series

The Pinnacle of KaVo Innovation.



MASTERtorque M8900 L
1008.7800

NEW



MASTERtorque Mini M8700 L
3001.0000

NEW

Expert Series

Superior KaVo Performance.



EXPERTtorque LUX E679 L
1006.9200



EXPERTtorque Mini LUX E675 L
1006.6400

Smart Series

Reliable and Affordable KaVo Technology.



SMARTtorque LUX S619 L
1008.1641




SMARTtorque Mini LUX S615 L
1008.1643




SMARTtorque S609 C
1008.1644




SMARTtorque Mini S605 C
1008.1645



SURGtorque LUX S459 L
3000.4808



SURGtorque S459 C
3000.5063



1390 MULTIflex LUX Coupler
0553.1390



460LED MULTIflex LUX Coupler
1007.3201



460LE MULTIflex LUX Coupler
1001.7599



465LED MULTIflex LUX Coupler
1007.3197

465LRN MULTIflex LUX Coupler
0553.1550




457 MULTIflex Coupler
0553.1210



Upgrade 460LE and 465LRN
KaVo MULTI LED Bulb
1007.5372

*Non-optic fixed back handpiece do not use MULTIflex coupler






High-speed

Master Series

The Pinnacle of
KaVo Innovation.

Featuring the new
MASTERtorque series
handpieces with unique
Direct Stop Technology™
for improved bur control
and reduced drawback of
contamination. Now with
even higher power and
same low noise level for the
highest treatment quality
and patient comfort.



High-speed

Expert Series

Superior KaVo Performance.

Featuring the EXPERTtorque
E679 L standard head and
E675 L mini head handpieces
with new innovative design and
EXPERTgrip.





High-speed

Smart Series

Reliable and Affordable
KaVo Technology.

Featuring the SMARTtorque
S619 L standard head and
S615 L mini head handpiece
with new ergonomic
design and legendary KaVo
durability.



Low-speed

Smart Series

Reliable and Affordable
KaVo Technology.

The KaVo Smart Series
INTRAmatic E Air Motor
Low-speed handpieces
have a number of
straightforward options
designed around
preferences and intended
clinical application.
Simply choose the right
handpiece for the job,
and the rest is easy.



MASTER SERIES

The Pinnacle of KaVo Innovation

MASTERtorque Handpieces



townie
choice
awards 2015



TRY KAVO NEW

MASTERtorque M8900 L

- New: Most powerful (23 watts) KaVo High-speed with ceramic bearing
- New: Direct Stop Technology
- DST stops the bur in one second, increasing efficiency and improving safety during treatment
- DST reduces drawback of particles into the ball bearings for improved durability and hygiene
- Quietest KaVo air High-speed at 57 decibels
- Field replaceable water filter
- Improved 4-port spray for better visibility
- Plasmatec coating
- New angles for improved access
- Push button auto chuck
- 3rd generation cellular optics with 5-year warranty
- Automatic drive air adjustment - Use at minimum drive air pressure 38 psi (2.6 bar)
- MULTIflex LUX connection
- 370,000 rpm
- Replaces 6500BR

Part No. 1008.7800 \$1,599.00

Part No. 1002.0271 \$23.00 Water filter

Part No. 1002.0321 \$2.75 Filter tool

MASTERtorque Mini M8700 L

- Mini Head (19 watts) KaVo High-speed with ceramic bearing
- Quietest KaVo mini head air High-speed at 59 decibels
- New: Direct Stop Technology
- DST stops the bur in one second, increasing efficiency and improving safety during treatment
- DST reduces drawback of particles into the ball bearings for improved durability and hygiene
- Field replaceable water filter
- 3-port spray
- Plasmatec coating
- New angles for improved access
- Push button auto chuck
- 3rd generation cellular optics with 5-year warranty
- Automatic drive air adjustment- Use at minimum drive air pressure 38 psi (2.6 bar)
- MULTIflex LUX connection
- 440,000 rpm

Part No. 3001.0000 \$1,599.00

Part No. 1000.4823 \$37.00 Water filter



townie
choice
awards 2015



TRY KAVO NEW

*2.5 years when maintained with QUATTROcare Plus.



EXPERT SERIES

Superior KaVo Performance.

EXPERTtorque Handpieces



townie
choice
awards 2015



EXPERTtorque LUX E679 L

- Ceramic bearing
- Standard head (18 watt)
- MULTIflex LUX connection
- Push button auto-chuck
- 4-port spray
- Replaces 659BR
- 3rd generation cellular optics with 5-year warranty
- 36 psi (2.5 bar)
- 380,000 rpm

Part No. 1006.9200 \$1,195.00

Part No. 1002.0271 \$23.00 Water filter

Part No. 1002.0321 \$2.75 Filter tool



townie
choice
awards 2015



EXPERTtorque Mini LUX E675 L

- Ceramic bearing
- Miniature head (15 watts)
- MULTIflex LUX connection
- Push button auto-chuck
- 3-port spray
- Replaces 635BR
- 3rd generation cellular optics with 5-year warranty
- 33 psi (2.3 bar)
- 420,000 rpm

Part No. 1006.6400 \$1,195.00

Part No. 1006.0800 \$1,195.00 EXPERTtorque LUX E675 LMW Midwest Stylus Coupler

Part No. 1002.0271 \$23.00 Water filter

Part No. 1002.0321 \$2.75 Filter tool



SMART SERIES

Reliable and Affordable KaVo Technology

SMARTtorque Handpieces



townie
choice
awards 2015



SMARTtorque LUX S619 L

- Fiber optic standard head (16 watts)
- Optimum balance between speed and torque
- MULTIflex connection
- Push-button auto chuck
- 4-port spray
- 40 psi (2.8 bar)
- 400,000 rpm
- 1-year warranty

Part No. 1008.1641

\$749.00



townie
choice
awards 2015



SMARTtorque Mini LUX S615 L

- Fiber optic miniature head (13 watts)
- Optimum balance between speed and torque
- MULTIflex connection
- Push-button auto chuck
- 1-port spray
- 450,000 rpm
- 40 psi (2.8 bar)
- 1-year warranty

Part No. 1008.1643

\$749.00



townie
choice
awards 2015



SMARTtorque S619 C

- Non-optic standard head (16 watts)
- Optimum balance between speed and torque
- MULTIflex connection
- Push-button auto chuck
- 4-port spray
- 40 psi (2.8 bar)
- 400,000 rpm
- 1-year warranty

Part No. 1008.1642

\$669.00



SMARTtorque Handpieces



townie
choice
awards 2015



townie
choice
awards 2015

SMARTtorque S609 C

- Non-optic standard head (16 watts)
- Optimum balance between speed and torque
- Standard 4-hole connection
- Push-button auto chuck
- 4-port spray
- 40 psi (2.8 bar)
- 400,000 rpm
- 1-year warranty

Part No. 1008.1644 \$569.00

SMARTtorque Mini S605 C

- Non-optic miniature head (13 watts)
- Optimum balance between speed and torque
- Standard 4-hole connection
- Push-button auto chuck
- 1-port spray
- 40 psi (2.8 bar)
- 400,000 rpm
- 1-year warranty

Part No. 1008.1645 \$569.00

SURGtorque Handpieces



TRY
KAVO
NEW

SURGtorque LUX S459 L

- Unique head angle for improved access to molar region (45°)
- Fiber optic, Powerful 19W turbine for increased efficiency
- MULTIflex connection
- Ceramic Bearings
- Exhaust ports direct excess exhaust away from the cutting area
- Single port water with no spray
- Push button auto chuck
- Accepts standard and surgical length burs
- 1 year warranty
- 380,000 rpm

Part No. 3000.4808 \$995.00



TRY
KAVO
NEW

SURGtorque S459 C

- Unique head angle for improved access to molar region (45°)
- Non-Fiber optic, Powerful 19W turbine for increased efficiency
- Standard 4-hole connection
- Ceramic Bearings
- Exhaust ports direct air away from the cutting area
- Single port water with no spray
- Push button auto chuck
- Accepts standard and surgical length burs
- 1 year warranty
- 420,000 rpm

Part No. 3000.5063 \$845.00

INTRAflex Air Motor System

This gear-driven system is ergonomic and compact with high torque and it fits on KaVo's unique MULTIflex system. The slow speed air handpiece's lighting system uses fused glass rods versus fiber bundle technology, resulting in no degradation of light quality after repeated sterilization.



townie
choice
awards 2009



Cellular Optic INTRAflex LUX 2 Air Motor System 2325LN with Integrated Motor and Fixed Head

- Transmission 1:4 speed-increasing
- Speed range up to 140,000 rpm
- Cellular optic
- Internal air/water coolant
- MULTIflex LUX connection
- Friction grip push button auto chuck
- 3-port spray
- Use at minimum drive air pressure 43 psi (2.9 bar)
- 1-year warranty

Part No. 0535.6210

\$1,999.00



townie
choice
awards 2009



Cellular Optic INTRAflex LUX 2 Air Motor System 2320LN Contra Angle with Integrated Motor

- Transmission 1:1
- Speed range up to 35,000 rpm
- Cellular optic
- Internal air/water coolant
- MULTIflex LUX connection
- Use at minimum drive air pressure 43 psi (2.9 bar)
- Fits INTRA LUX Heads
- 1-year warranty

Part No. 0535.6180

\$1,599.00



townie
choice
awards 2009



Cellular Optic INTRAflex LUX 2 Air Motor System 2307LN Contra Angle with Integrated Motor

- Reduction 2.7:1
- Speed range up to 13,000 rpm
- Cellular optic
- Internal air/water coolant
- MULTIflex LUX connection
- Use at minimum drive air pressure 43 psi (2.9 bar)
- Fits INTRA LUX Heads
- 1-year warranty

Part No. 0535.6190

\$1,699.00



townie
choice
awards 2009



Non-Optic INTRAflex 2 Air Motor System 2313LN Straight Handpiece w/ Integrated Motor

- Reduction 2.7:1
- Speed range up to 13,000 rpm
- Non optic
- Internal air/water coolant
- MULTIflex connection
- For handpiece and contra angle burs
- Comes with bur stop
- Maintain with chuck in closed position 43 PSI
- 1-year warranty

Part No. 0535.6220

\$1,359.00

* See page 14 for INTRA LUX heads



460LED MULTIflex LUX Coupler

- 6-pin fiber optic coupler with MULTI LED bulb in coupler
- Anti-retraction valve (Same as 465LRN but without water spray regulator)
- Use with LCM Illuminator and tubing
- Not sterilizable

Part No. 1007.3201 \$372.00



465LED MULTIflex LUX Coupler

- 6-pin fiber optic coupler with MULTI LED bulb in coupler
- Spray regulator
- Anti-retraction valve (Same as 465LRN but without water spray regulator)
- Use with LCM Illuminator and tubing
- Not sterilizable

Part No. 1007.3197 \$472.00



KaVo MULTI LED Bulb

KaVo MULTI LED for all your 6-pin MULTIflex couplers. KaVo's innovative solution is an LED that can easily replace your existing halogen lamp. MULTI LED offers efficient and superior light, with color temperature similar to daylight (5,500 K and 25,000 LUX) at the bur tip. Operating voltage range 3.0 – 3.6 V DC. Nominal voltage 3.4 V DC. Current < 120 mA.

Part No. 1007.5372 \$125.00



460LE MULTIflex LUX Coupler

- 6-pin fiber optic coupler with light bulb in coupler
- Anti-retraction valve (Same as 465LRN but without water spray regulator)
- Use with LCM Illuminator and tubing
- Not sterilizable

Part No. 1001.7599 \$323.00



465LRN MULTIflex LUX Coupler

- 6-pin connection
- Spray regulator
- Anti-retraction valve
- Use with LCM Illuminator and tubing
- Light bulb in coupler
- Not sterilizable

Part No. 0553.1550 \$415.00



Fits product 465LRN / LED and 460LED

1390 MULTIflex LUX Coupler

- Anti-retraction valve
- Use with LS2 Illuminator and tubing
- Not sterilizable

Part No. 0553.1390 \$359.00



457 MULTIflex Coupler


- 4-hole connection
- Non-optic
- Not sterilizable

Part No. 0553.1210 \$265.00

For Coupler parts, see page 39.
For light source parts, see page 77.



INTRAmatic E Air System



townie
choice
awards 2005

2.5
YEAR
WARRANTY

135°C

TRY
KAVO

20E Low-speed Handpiece

- Transmission 1:1
- Speed range up to 20,000 rpm
- Push button auto chuck for latch burs
- 100° head for superior access
- Durable hygienic coating with textured surface for easy grip
- Use with 181K Air Motor

Part No. 1003.1107

\$619.00



townie
choice
awards 2005

2.5
YEAR
WARRANTY

135°C

TRY
KAVO

19E Hygiene Handpiece

- Reduction 7.4:1
- Speed range up to 2,700 rpm
- For disposable prophy heads only
- Durable hygienic coating with textured surface for easy grip
- Use with 181K Air Motor

Part No. 1004.6403

\$519.00

Part No. 0201.0116

\$21.00 Screw for Prophy Angle



townie
choice
awards 2005

2.5
YEAR
WARRANTY

135°C

TRY
KAVO

10E Low-speed Handpiece

- Transmission 1:1
- Speed range up to 20,000 rpm
- For use with handpiece burs
- Durable hygienic coating with textured surface for easy grip
- Use with 181K Air Motor

Part No. 1003.6977

\$550.00



townie
choice
awards 2005

2.5
YEAR
WARRANTY

135°C

80E Low-speed Handpiece

- Transmission 1:1
- Speed range up to 20,000 rpm
- Push button auto chuck for friction grip burs 1.6mm
- 100° head for superior access
- Durable hygienic coating with textured surface for easy grip
- Use with 181K Air Motor

Part No. 1004.6401

\$619.00



townie
choice
awards 2005

2.5
YEAR
WARRANTY

135°C

14E Low-speed Handpiece

- Reduction 4:1 and for latch burs
- Max. 5,000 rpm
- Push button chuck
- 100° head angle for superior access
- Use with 181K Air motor

Part No. 1003.1108

\$569.00



townie
choice
awards 2005

2.5
YEAR
WARRANTY

135°C

31E Hygiene Handpiece

- Reduction 4:1
- Speed range up to 5,000 rpm
- For use with snap-on / screw-in cups and brushes
- Extremely small head with 100° head angle for superior access
- Durable hygienic coating with textured surface for easy grip
- Use with 181K Air Motor
- Comes with Prophy cup holder

Part No. 1003.6976

\$615.00

Part No. 1004.0596

\$4.25 Holder Tool

Master Series

181L INTRA LUX Air Motor

dental
townie
choice
awards2007

Works with E attachments

- Optic
- Speed range up to 20,000 rpm
- Forward/reverse speed adjustment
- Internal air/water coolant
- 6-pin connection

- Fits all attachments optic or non-optic with Universal "E" type connection
- Anti-retraction system
- Use with KaVo LCM
- 31-38 PSI
- 1-year warranty

Part No. 0535.5210

\$1,959.00

Expert Series

181H INTRA Air Motor

dental
townie
choice
awards2007

Works with E attachments

- Non-optic
- Speed range up to 20,000 rpm
- Forward/reverse speed adjustment
- Internal air/water coolant
- 4-hole connection

- Fits all non optic attachments with Universal "E" type connection
- Anti-retraction system
- Sterilizable
- 31-38 PSI
- 1-year warranty

Part No. 0535.5200

\$1,720.00

Smart Series

181K Air Motor

REALITY
Five Star Award
REALITY'S
CHOICE

Works with INTRAmatic E attachments

- Speed range up to 20,000 rpm
- 4-hole connection
- Non optic
- Forward/reverse speed adjustment
- Sterilizable

- 36-43 PSI
- No water
- No chip air
- 1-year warranty

Part No. 0535.5260

\$789.00

INTRAmatic Air Motor Accessories



Part No.	Price	Description
0200.6120	\$17.25	Motor O-ring for 181L, 181H, 196, K200, 700KL, 701KL, 702KL, KL703 Motors Requires four O-rings
0553.1872	\$14.00	4-hole gasket for 181H and 181K Motors & non-optic MULTIflex Coupler 465RN, 457, 459 & non-optic High-speeds
1002.2928	\$105.00	Halogen Light Bulb (pack of 2) for 181L, 196, K 200, 701KL, 702KL 3.0-3.3 Volts No warranty on bulbs
0246.5019	\$11.00	Set screw for 181H
0535.4242	\$62.00	Sheath for 181L
0535.1661	\$53.00	Bulb holder for 181L



Genuine KaVo Turbines

Genuine KaVo turbines are individually balanced and electronically tested to meet the highest concentricity requirements. Properly balanced turbines leads to lower noise during use.



KaVo REPAIR PARTS

Quality in the Finest Details

When repairing equipment and handpieces, choose KaVo Original Parts. Original parts are manufactured to the same standards as your original purchase, which means better performance, durability, and safety. Alternate choices may not have all the innovation that KaVo puts into its products, like computer balanced turbines and two-piece carbide chuck sleeves. Restore your KaVo equipment to the standards by which it was originally manufactured.



Turbine Cartridges



	Product	Part No.	Price	Description
	6500BR	1003.6884	\$395.00	1 auto chuck turbine
	4500BR, M4500 L	1006.4564	\$395.00	1 auto chuck turbine
	E679L	1007.9313	\$315.00	1 auto chuck turbine
	E675L	1007.9457	\$315.00	1 auto chuck turbine
	S619/609	1008.4163	\$205.00	1 auto chuck turbine
	S615/605	1008.4164	\$205.00	1 auto chuck turbine
	S459 L/C	1011.2518	\$315.00	1 auto chuck turbine
	6500B	1003.9664	\$410.00	1 auto chuck turbine
	6000B	1001.3888	\$415.00	1 auto chuck turbine
Single Packs	659B	1004.5632	\$350.00	1 auto chuck turbine
	659BR	1000.2078	\$375.00	1 auto chuck turbine
	635B	0553.8581	\$415.00	1 auto chuck turbine
	635BR	0553.8661	\$415.00	1 auto chuck turbine
	647B, 649B	0553.8541	\$465.00	1 auto chuck turbine
	636CP, 646B	0553.8911	\$225.00	1 auto chuck turbine
	607C	0980.6649	\$255.00	1 auto chuck turbine
Triple Packs	640B, 630B, 630A	0553.5881	\$445.00	1 auto chuck turbine
	642B, 632, 643	0553.6911	\$445.00	1 auto chuck turbine
	630	0553.6671	\$490.00	1 friction grip turbine 1 back cap tool
	647B, 649B	0906.2179	\$1,360.00	3 turbines 1 back cap tool
	640B, 630B	0906.0300	\$1,360.00	3 turbines 1 back cap tool
	642B, 632B	0906.0310	\$1,360.00	3 turbines 1 back cap tool



Accessories for High-speed Handpieces

Back Caps



Product	Part No.	Price	Description
6500B, 6500BR	1003.5525	\$199.00	Push button auto chuck
4500B, 4500BR	1006.1079	\$109.00	Push button auto chuck
E679 L	1007.8930	\$62.00	Push button auto chuck
E675 L, E675 LMW	1007.7051	\$62.00	Push button auto chuck
S619 L, S615C, S609C	1007.0854	\$69.00	Push button auto chuck
S615 L, S605C	1007.7277	\$69.00	Push button auto chuck
S459 L/C	1007.0854	\$69.00	Push button auto chuck
6000B	1000.4770	\$201.00	Push button auto chuck
659B	1000.2076	\$201.00	Push button auto chuck
635B	0553.8491	\$201.00	Push button auto chuck
647B, 649B, 646B, 636CP	0553.8061	\$201.00	Push button auto chuck
632, 642	0553.5081	\$201.00	Push button auto chuck
630B	0553.4041	\$201.00	Push button auto chuck
643B	0553.7631	\$201.00	Push button auto chuck
640B	0553.6511	\$201.00	Push button auto chuck
640FG	0553.8742	\$159.00	Friction grip
630	0553.1782	\$159.00	Friction grip

Back Cap Removal Tools



Handpiece	Part No.	Price	Type
635 B/635 BR	0411.3063	\$32.00	Push button auto chuck
647B, 649B, 646B, 6000B, 6500B/BR, 636CP, E679L	0411.3053	\$21.00	Push button auto chuck
640, 642, 643, 630, 632	0411.1861	\$52.00	Push button auto chuck
4500B/BR, 5000B, E675L	1008.6133	\$29.00	Push button auto chuck
625C, 630	0411.0831	\$62.00	Friction grip
634 F4, 632, 642	0411.1701	\$62.00	Friction grip
S619 / 609	0411.2361	\$183.00	Push button auto chuck
S615 / 605	1006.3384	\$129.00	Push button auto chuck

Coupler Bulbs

Couplers	Part No.	Price	Includes
465LRN, 459AN, 181L, 460LE	1002.2928	\$105.00	Halogen bulb (2 pack) Voltage 3.0 - 3.3V. 2.5 watts
460LED, 465LED	1007.5372	\$125.00	LED bulb

Sterilizable. 175 hour of use.

Coupler Tools

Handpieces	Part No.	Price	Includes
all models	0411.1563	\$14.00	Coupler tightening tool
all models	0905.8088	\$449.00	Coupler straightening tool

Accessories

Couplers	Part No.	Price	Includes
465LRN / LED, 459AN, 460LE / LED	0553.9482	\$32.00	Screw-on bulb sheath Must order O-ring 0200.6033
all models	0200.6033	\$7.25	O-ring for bulb sheath
all models	1005.1033	\$25.00	Set of 5 O-rings: 3 large, 2 small
457, 459, 465RN	0553.1872	\$14.00	4-hole coupler gasket non optic High-speed handpieces, motors
1390	0553.4762	\$23.00	5-hole coupler gasket
465LRN / LED, 459AN, 460LE / LED	0553.5262	\$23.00	6-pin coupler gasket

Bur Tools

Handpiece(s)	Part No.	Price	Type
632, 633, 642	0411.3871	\$59.00	Friction grip (black)
630, 640, 625C	0411.3461	\$59.00	Friction grip (blue)
80G, 23D	0411.2501	\$53.00	Friction Grip (tan)
8A	0411.3521	\$215.00	For 8A Surgical Attachment

Spray Repair Kits

Handpiece(s)	Part No.	Price
640, 630 (A, B or C)	0410.0600	\$239.00
647 (A, B or C)	0410.0650	\$239.00
649	0410.0660	\$239.00
Water Filter, 6000B, 6500B, 6500BR, 4500BR, 4500B	1000.4823	\$37.00


Accessories

Item(s)	Part No.	Price
O-ring for High-speed Handpiece	0200.6172	\$4.50
Cleaning needle	0410.0921	\$9.50
Cleaning brush	0229.3005	\$17.00




KaVo BUNDLE PROGRAMS

Purchasing large equipment can be daunting, that's why KaVo provides handpiece bundles – a simple, easy choice to outfit operatories with new KaVo handpieces and hygiene equipment. Bundles are the best deal on handpieces when you purchase any dental delivery unit.





Package Part #	Part #	Retail Price	Quantity	Hygiene Non-Optic Package	Total Savings
1008.3301	1007.1606	\$3,295	1	SONICflex quick 2008 Set	\$565
	1006.9926		1	PROPHYflex 3 Violet	
	0553.1210		1	457 MULTIflex Coupler	



Package Part #	Part #	Retail Price	Quantity	Hygiene Optic Package	Total Savings
1008.3302	1007.1604	\$3,395	1	SONICflex quick 2008 L Set	\$525
	1006.9926		1	PROPHYflex 3 Violet	
	1001.7599		1	460LE MULTIflex LUX Coupler	



Package Part #	Part #	Retail Price	Quantity	Hygiene 360 Package	Total Savings
1008.6390	1003.2278	\$7,725	1	PROPHYwiz 181P	\$815
	1006.9926		1	PROPHYflex 3 Violet	
	1007.1604		1	SONICflex quick 2008 L Set	
	1004.3400		1	DIAGNOdent pen	
	1004.8400		1	DIAGNOdent display	
	1001.7599		1	460LE MULTIflex LUX Coupler	



Part #	Retail Price	Quantity	Hygiene Package Buy	Get	Savings up to
1003.2278	\$865	1	PROPHYwiz 181P	TWO Boxes of Kerr Perfect Pearl prophy angles FREE All boxes contain 200 ct	\$700
1003.9882	\$1,815	1	PROPHYwiz 3-pack	SIX Boxes of Kerr Perfect Pearl prophy angles FREE All boxes contain 200 ct	



Part #	Pricing	Quantity	Hygiene Package	Savings up to
1007.8889	\$1,665	1	181K Motor	\$680
		3	19E Hygiene Handpiece	

Purchase of dental unit is not required to take advantage of hygiene bundle pricing.

Handpiece Bundle Program*

- - Electric
- - Hygiene
- - Maintenance
- - Air

	Part #	Pricing	Quantity	COMFORTdrive / High-speed	Total Savings**
C1	1007.8108	\$3,855	3	COMFORTdrive 200 XDR	\$1,225
E1	1005.7646	\$3,935	3	MASTERmatic LUX M25 L	\$1,070
E2	1011.9866	\$3,935	3	MASTERmatic LUX M05 L Mini	\$1,070
E3	1010.7113	\$2,965	3	EXPERTmatic E25L	\$1,320

	Part #	Pricing	Quantity	Endo Premium / Endo Basic	Quantity	Low-speed / Endo Package	Total Savings**
E4	1007.1232	\$7,025	3	MASTERmatic LUX M25 L	3	MASTERmatic LUX M20 L	\$2,650
					2	L68 B Latch Head – Operative	
					1	L66 B Reduction Latch Head – Endodontics	
E5	1005.7647	\$4,855	3	MASTERmatic LUX M25 L	1	MASTERmatic LUX M20 L	\$1,725
					1	L66 B Reduction Latch Head – Endodontics	

	Part #	Pricing	Quantity	High-speed Package	Quantity	Low-speed Package	Total Savings**
E6	1007.1231	\$5,925	3	MASTERmatic LUX M25 L	2	MASTERmatic LUX M20 L	\$2,280
					2	L68 B Latch Head – Operative	
E7	1011.9867	\$5,925	2	MASTERmatic LUX M25 L	2	MASTERmatic LUX M20 L	\$2,280
					2	L68 B Latch Head – Operative	
E8	1007.8887	\$2,095			2	MASTERmatic LUX M20 L	\$1,015
					2	L68 B Latch Head – Operative	

	Part #	Pricing	Quantity	Expert Electric Premium	Quantity	Low-speed Package	Total Savings**
E9	1010.7114	\$5,745	3	EXPERTmatic E25L	2	EXPERTmatic E20L	\$2,445
					1	EXPERTmatic E15L	

	Part #	Pricing	Quantity	Expert Electric Basic	Quantity	Low-speed Package	Total Savings**
E10	1010.7115	\$4,915	3	EXPERTmatic E25L	2	EXPERTmatic E20L	\$2,090

	Part #	Pricing	Add-ons		
D1 ⁺	1004.3400	\$2,849	Diagnostic Package – DIAGNODent pen Caries Detection Aid		
D2 ⁺	0574.0500	\$2,848	Diagnostic Package – DIAGNODent Classic Caries Detection Aid		
M1 ⁺	1008.3805	\$3,259	Maintenance Package – QUATTROcare Plus		

* Prices valid only with purchase of any manufacturing dental unit.

** Total savings over individual instruments full retail price (rounded off).

† Purchase of dental unit is not required to take advantage of this bundle pricing.

Handpiece Bundle Program*

6-Pin	Part #	Pricing	Quantity	High-speed 6-Pin Package	Quantity	Low-speed Package	Total Savings**
A1	1007.1233	\$5,265	1	MASTERtorque M8900 L	2	20E Low-speed Handpiece	\$2,480
			2	MASTERtorque Mini M8700 L	1	10E Low-speed Handpiece	
			1	MULTIflex LUX 460LED Coupler	1	181K Motor	
A2	1005.7648	\$3,685	3	MASTERtorque M8900 L			\$1,480
			1	MULTIflex LUX 460LED Coupler			
A3	1007.1226	\$3,685	3	MASTERtorque Mini M8700 L			\$1,480
			1	MULTIflex LUX 460LED Coupler			
A4	1010.7117	\$2,675	3	EXPERTtorque E679L			\$1,280
			1	MULTIflex LUX 460LED Coupler			
A5	1010.7118	\$2,675	3	EXPERTtorque E675L			\$1,280
			1	MULTIflex LUX 460LED Coupler			
A6	1008.6971	\$1,825	3	SMARTtorque LUX S619L			\$790
			1	MULTIflex LUX 460LED Coupler			
A7	1008.7000	\$1,825	3	SMARTtorque mini LUX S615L			\$790
			1	MULTIflex LUX 460LED Coupler			
5-Hole	Part #	Pricing	Quantity	High-speed 5-Hole Package	Quantity	Low-speed Package	Total Savings**
A8	1007.1234	\$5,265	1	MASTERtorque M8900 L	2	20E Low-speed Handpiece	\$2,470
			2	MASTERtorque Mini M8700 L	1	10E Low-speed Handpiece	
			1	MULTIflex LUX 1390 Coupler	1	181K Motor	
A9	1006.8293	\$3,685	3	MASTERtorque M8900 L			\$1,470
			1	MULTIflex LUX 1390 Coupler			
A10	1007.1227	\$3,685	3	MASTERtorque Mini M8700 L			\$1,470
			1	MULTIflex LUX 1390 Coupler			
A11	1010.7119	\$2,675	3	EXPERTtorque E679L			\$1,270
			1	MULTIflex LUX 1390 Coupler			
A12	1010.7120	\$2,675	3	EXPERTtorque E675L			\$1,270
			1	MULTIflex LUX 1390 Coupler			
A13	1008.6999	\$1,825	3	SMARTtorque LUX S619L			\$780
			1	MULTIflex LUX 1390 Coupler			
A14	1008.7001	\$1,825	3	SMARTtorque mini LUX S615L			\$780
			1	MULTIflex LUX 1390 Coupler			
4-Hole	Part #	Pricing	Quantity	Low-speed 4-Hole Package	Quantity	Low-speed Package	Total Savings**
A15	1007.8888	\$1,725			2	20E Low-speed Handpiece	\$850
					1	10E Low-speed Handpiece	
					1	181K Motor	

* Prices valid only with purchase of any manufacturing dental unit.

** Total savings over individual instruments full retail price (rounded off).



QUATTROcare Plus

QUATTROcare Plus Maintenance System is designed to clean, lubricate and purge your instruments with the simple push of a button. Chuck maintenance is now an easy step included with the QUATTROcare Plus.



KaVo QUATTROcare PLUS

Protect Your Investments

Proper cleaning of rotating instruments is essential before sterilization. If debris remains inside of the instrument, it can cause handpiece failure. Proper care means longer handpiece life, reduced repair costs, less frustration and elimination of costly downtime.



The QUATTROcare Plus Advantage

Easy to use

The user-friendly QUATTROcare Plus system is designed to clean, lubricate and purge your instruments with the simple push of a button. Loading and operation is easy as well as the automatic care and maintenance procedure.

Safe controlled cleaning

Reduce errors with visual and audio signals during the completely controlled cleaning and purging procedure. With a front-loading, enclosed maintenance chamber, the entire process is hygienic and thus eliminates insufficient care and maintenance.

Economical investment

Using the QUATTROcare Plus system means less wear, less repairs, less staff down time, smoother operation, and increased longevity of your handpieces – what an excellent return on investment! The QUATTROcare Plus can be used for virtually all manufacturers' High-speeds, Low-speed contra-angles, straight handpieces, air driven scalers, and air motors.

Fast cleaning cycle

Up to four (4) handpieces or air motors are cleaned, lubricated and purged in the same cycle. The care and maintenance of one complete set of handpieces in the same cycle is possible which means there's no need to change adapters between High-speed and Low-speed attachments. Flexible programming allows you to choose and set total cycle time from 20–45 seconds per handpiece. Handpieces can be sterilized immediately after maintenance cycle; there is no need for separate purging of the handpiece.

Chuck care maintenance

New chuck care system provides easy maintenance of handpiece chuck. A clean chuck helps improve the bur retention strength and ensures optimal longevity of the turbine components.

Easy As...



1

Position the instruments...



2

... close the door ...




3

... press Start and everything is DONE.



NEW! Chuck care



QUATTROcare Plus

QUATTROcare Plus

Automatic handpiece maintenance and purging system. Correct maintenance of dental handpieces is critical to ensure top performance and longevity. KaVo QUATTROcare Plus is equipped with two MULTIflex adapter and two INTRAmatic (E-Type) adapters. Unit also includes one can of QUATTROcare Plus spray.

System advantages:

- Up to four (4) handpieces or air motors are cleaned, lubricated and purged in the same cycle
- Flexible programming allows setting of total cycle time from 20-45 seconds per handpiece
- Can be used for virtually all manufacturers' handpieces, air driven scalers and air motors
- Unique spray formula expands inside each handpiece, creating a highly effective cleaning foam
- Completely controlled cleaning and purging procedure
- Chuck care system

Technical:

- 120V, 50/60Hz, 20 Amps electrical connection
- Air consumption: 50 mL/min or 1.765 cfm
- Simply connects to 1/4" airline (58 to 87 psi air pressure ideal (max 145 psi)
- Built in automatic pressure regulator controls flow

Dimensions:

- Height: 15", Width: 13-9/16", Depth: 11-1/2"
- Weight: 22 lbs
- Width with door ajar (right side - to exchange spray can): 16-15/16". Depth with flap open (front): 18-7/8"

Part No. 1008.3805 \$3,259.00

Accessories for QUATTROcare

Part No.	Price	Description
1009.6142	\$75.00	Adapter for MULTIflex
1009.6143	\$75.00	Adapter for E-Style
1005.1707	\$99.00	Adapter for COMFORTdrive
0200.6084	\$6.50	O-rings for COMFORTdrive adapter
0411.7941	\$139.00	Intra Head adapter
0411.5532	\$81.00	Bumper sheath for INTRA head adapter
1000.8459	\$169.00	Adapter for Midwest "Shorty" attachment
1001.7860	\$169.00	Adapter for Bien Air High-speeds
1001.7888	\$169.00	Adapter for W&H High-speeds
1000.8786	\$169.00	Adapter for NSK High-speeds
1000.1003	\$169.00	Adapter for SIRONA Electric T1 Instruments
0411.5412	\$17.00	U-clip for door, 2 needed
0200.6098	\$4.50	O-ring, Green, Small for front, 2 each
0200.6177	\$4.50	O-ring, Green, Large for front, 2 each
0200.6089	\$4.50	Low-speed O-ring
0324.4301	\$11.00	Blue hose for unit air hook up , 1/4"
0692.6891	\$56.00	Transformer power cord
0220.0977	\$49.00	Male Connector
0223.2786	\$2.75	Fuse for QUATTROcare
0411.4472	\$9.50	Union Nut
0411.7603	\$29.00	Adapter for chuck maintenance
1007.4633	\$9.25	Internal Air Filter
0411.4572	\$50.00	Catch Pan



**KaVo QUATTROcare plus Spray* (6 cans)**

- 6 cans
- For cleaning and care of instruments / air motors in the KaVo QUATTROcare
- Designed specifically to withstand sterilization
- Relubrication after sterilization is unnecessary, eliminating the risk of recontaminating sterilized handpieces
- 500 ml

Part No. 1005.4524 \$369.00

**KaVo Spray* (6 cans)**


- 6 cans
- Designed specifically to withstand sterilization
- Relubrication after sterilization is unnecessary, eliminating the risk of recontaminating sterilized handpieces
- 500 ml
- Nozzles sold separately

Part No. 0411.9660 \$389.00

**MULTIflex Nozzle**

- For maintenance of all KaVo MULTIflex Handpieces, INTRAflex LUX Motors, SONICflex Scalers and the 181 M Air Motor

Part No. 0411.9921 \$23.00

**Quick Nozzle**

- For all INTRA and INTRAflex LUX Heads
- 4-hole fixed back handpieces
- All KaVo air motors (except 181M), all handpieces with Universal "E" type connection and surgical handpieces

Part No. 0411.9931 \$23.00

Part No. 1005.3154 \$22.00 COMFORTdrive Nozzle

**KaVo Cleanpac**

- Keeps hands clean during cleaning and lubricating handpieces
- Cleanpac traps excess lubricant while allowing view for flushing
- Clean Pac (pack of 10)

Part No. 0411.9691 \$57.00 (pack of 10)



Part No.	Price	Description
0411.9501	\$189.00	2151 Handpiece stand sterilizable stand for handpieces, attachments, motors and heads. Includes 3 inserts for MULTIflex instruments.
0411.9902	\$19.00	Insert for MULTIflex instruments
0411.9862	\$72.00	Drip Pads. 100 pads for handpiece stand Absorbs excess lubricant and/or condensation.

*Hazmat Fees - refer to page 81 for Terms of Sale



HYGIENE 360°

Everything Easy. Everything KaVo.



KaVo Hygiene 360° is a complete solution which helps you provide your patients with great dental care. Start the 360° cycle with the SONICflex Scaler for scaling, then PROPHYflex for tooth cleaning, PROPHYwiz for tooth polishing, and DIAGNOdent for caries and sub-gingival calculus detection. Complete your hygiene workflow by using a stool that properly supports you ergonomically.





Step 1 Scaling

Difficult cases?

The problem is solved so easily.

Both thorough and gentle at the same time.

Your hygiene treatment center takes difficult cases in its stride. Tenacious calculus requires thorough, yet gentle removal. One of the most powerful scalers in its class - the KaVo SONICflex Quick 2008L easily conquers tough challenges.



Step 2 Tooth Cleaning

Just conjure up a smile for the face.

Your patients see results immediately.

The next step is about preventing caries, periodontitis and gingivitis. The main causes of these oral indicators, such as biofilm, deposits from smoking and discolorations from coffee and tea can be easily removed with the KaVo 2018 PROPHYflex 3 air polishing device.

KaVo 2018 PROPHYflex 3



Step 3 Tooth Polishing

You've got the skills and technology to maintain healthy dental enamel.

Uncompromising outcome, means everything works great.

For a great finish to every professional tooth cleaning, the final polish and seal of the enamel is essential. KaVo offers the right tool for all indications, in each price segment.

KaVo PROPHYwiz



Step 4 Caries Detection

Your aid in detection of caries.

The DIAGNOdent pen will show you where the decay lies.

It works quickly and reliably. The laser fluorescence detector within the DIAGNOdent pen, is a precise method for identifying fissure caries, proximal caries and periodontitis (with the Perio probe). A virtually unbeatable method with a 90% success rate. The wireless display ensures additional comfort, with a large status and measurement display.

KaVo DIAGNOdent pen + Display



Step 5 Workflow Integration

Your comfort is of the utmost importance.

With adjustable levers at your finger tips proper positioning is easy. Stay seated and on task during patient appointments.



KaVo SONICflex SYSTEM

Variety Developed to a Fine Point




The right fit for every use. The wide range of options offered by the SONICflex will soon make it one of your most-used instruments. Choose from the right tip for your needs from the wide range of types. Over 50 tips for every application.





SONICflex 2008/2008 L LUX Air Scaler



The SONICflex 2008 has single turn tip installation and removal. Air driven. For removal of calculus and a wide range of other applications.

- Air/water coolant exits at tip
- Hygienic grip
- Three different performance levels
- Standard with three tips (Universal, Sickle, Perio) and tip removal tool
- MULTIflex connection - optic (2008 L), non-optic (2008)
- Automatic drive air regulator
- (2.5 bar)
- Operates at 36psi
- Set water at 15 psi
- Only compatible with KaVo SONICflex "A" tips

Part No. 1007.1604

\$2,029.00 SONICflex quick 2008 L Set

Part No. 1005.9310

\$1,689.00 SONICflex quick 2008 L (w/o tips)

Part No. 1007.1606

\$1,899.00 SONICflex quick 2008 Set

Part No. 1005.9311

\$1,589.00 SONICflex quick 2008 (w/o tips)

SONICflex 2003/2003 L Air Scaler



Cellular optic glass rod (SONICflex 2003 L) and Non-Fiber optic (SONICflex 2003). Air driven. For removal of calculus and a wide range of other applications.

- Hygienic grip
- Three different performance levels
- 50% more powerful than previous 2000L/N
- Standard with three tips (Universal, Sickle, Perio) and tip removal tool
- MULTIflex connection - optic (2003 L), non-optic (2003)
- Automatic drive air regulator
- (2.5 bar)
- Operates at 36psi
- Set water at 15 psi

Part No. 1000.8333

\$2,029.00 SONICflex LUX 2003L Set

Part No. 1000.4243

\$1,689.00 SONICflex 2003L (w/o tips)

Part No. 1000.8334

\$1,899.00 SONICflex 2003 Set

Part No. 1000.4246

\$1,589.00 SONICflex 2003 (w/o tips)



Tip Removal Tool and Key for Adapter



For changing tips. Sterilizable. Torque feature prevents over tightening of tips. Can be used as protective cover when SONICflex is attached to delivery system.

Part No. 1000.4887

\$75.00 Tip Removal Tool

Part No. 0411.0892

\$21.00 Key

Part No. 0410.0911

\$18.00 Cleaning Tip



Hygiene and periodontology

SONICflex Scaler Tips*



A variety of tips for tenacious calculus removal. Clean shallow pockets tips for effective and gentle calculus removal, hygiene.

Tip Kit (No. 5, 6, 7)		Part No.	Price	Description	Tip No.
		1000.4983	\$495.00	Three tips for supra- and subgingival calculus removal	5, 6, 7
		1005.9328	\$495.00	SONICflex scaler Tip Set A	5A, 6A, 7A
SONICflex Scaler Tips Universal		0571.5171	\$179.00	For calculus removal on all tooth surfaces	5
		1005.8949	\$179.00		5A
		0571.5181	\$179.00	For pedo use and areas of difficult access	6
		1005.8950	\$179.00		6A
		0571.5191	\$179.00	For subgingival scaling	7
		1005.8951	\$179.00		7A
		0571.5371	\$179.00	For subgingival Calculus	8
		1006.1953	\$179.00		8A

SONICflex Scaler Tips*



Wide spectrum vibrating brush. For gentle cleaning and polishing.


SONICflex Clean Sets and Tips		Part No.	Price	Description	Tip No.
		0571.0401	\$145.00	Clean Insert tip (w/o brush)	48
		1006.1982	\$145.00	Clean Insert tip (w/o brush)	48A
		0571.0001	\$299.00	Clean Set	(incl. Tip 48, plus 5 units of 6 brushes)
		1006.2026	\$299.00	Clean Set	(incl. Tip 48A, plus 5 units of 6 brushes)
SONICflex Clean Brushes		0571.0402	\$19.00	SONICflex Clean Brush	1
Clean Brushes fit into SONICflex tip 48		0571.0412	\$19.00	SONICflex Clean Brush	2
		0571.0422	\$19.00	SONICflex Clean Brush	3
		0571.0432	\$19.00	SONICflex Clean Brush	4
		0571.0472	\$19.00	SONICflex Clean Brush	5
		0571.0482	\$19.00	SONICflex Clean Brush	6

*All tips also work with SONICflex 2000L and 2000N.

Note: Tips ending in suffix 'A' are compatible with SONICflex 2008 / SONICflex 2008L only.

Hygiene and periodontology (cont.)

SONICflex Perio



Different tips for minimally invasive, gentle and efficient cleaning of periodontal pockets up to 9 mm deep.

SONICflex Perio tips



	Part No.	Price	Description	Tip No.
60	0571.7402	\$179.00	Straight - For anterior region	60
61	1006.1934	\$179.00		60A
61	0571.7412	\$179.00	Left - for furcations, molar regions	61
61A	1006.1935	\$179.00		61A
62	0571.7422	\$179.00	Right - for furcations, molar regions	62
62A	1006.1936	\$179.00		62A

SONICflex Perio accessories



Coupling for the supply of therapeutic fluids such as sterile saline solution
(Part No. : 1000.6296) Retail Price \$68.00

SONICflex Implant



Effective removal of sub-gingival plaque and calculus. No abrasion, no roughening. Disposable cleaning attachment to prevent the transmission of diseases; a very favorable price/performance ratio.

SONICflex Implant tips



	Part No.	Price	Description	Tip No.
48	1003.8167	\$235.00	Incl. Tip 48, plus 10 cleaning pins	48
48	1006.2027	\$235.00	Incl. Tip 48, plus 10 cleaning pins	48A
48	0571.0401	\$145.00	Clean tip	48
48	1006.1982	\$145.00	Clean tip	48A
	1003.8168	\$201.00	SONICflex Implant Refill Add-on kit with 30 pins	

SONICflex Root planer




For easily processing of difficult-to-access or narrow furcations

SONICflex Root planer insert – kit

	Part No.	Price	Description	Tip No.
	0571.0301	\$739.00	Root planer Tips	24, 25, 26, 27
	1006.2022	\$739.00	Root planer Tip Set A	24A, 25A, 26A, 27A

SONICflex Root planer tips



	Part No.	Price	Description	Tip No.
24	0571.5621	\$215.00	Small Bud, Left - hard to reach furcations	24
25	1006.1957	\$215.00	Furcations	24A
26	0571.5631	\$215.00	Small Bud, Right - hard to reach furcations	25
26	1006.1959	\$215.00	Furcations	25A
26	0571.5641	\$215.00	Small Bud, Universal	26
26	1006.1961	\$215.00		26A
27	0571.5651	\$215.00	Large Bud, Perio	27
27	1006.1963	\$215.00		27A

*All tips also work with SONICflex 2000L and 2000N.

Note: Tips ending in suffix 'A' are compatible with SONICflex 2008 / SONICflex 2008L only.

Endodontics

SONICflex Endo



SONICflex endo – a special line of tips for probing for canals.
Removing soft materials during revision. Fashioning canal inlets.
Shaping the pulp chamber for coronal root preparation

SONICflex Endo Tip Set A		Part No.	Price	Description	Tip No.
		1000.5828	\$699.00	SONICflex endo Tip Set	66, 67, 68, 69, 70
		1006.2031	\$699.00	SONICflex endo Tip Set A	66A, 67A, 68A, 69A, 70A
SONICflex Endo tips		1000.5825	\$205.00	Large bud, eval of access cavities	66
	66	1006.1992	\$205.00		66A
	67	1000.5822	\$205.00	Tip angle 125° expansion of larger channels	67
	68	1006.1994	\$205.00		67A
	69	1000.5823	\$205.00	Tip angle 117° Expansion of larger channels	68
	70	1006.1996	\$205.00		68A
	66	1000.5827	\$205.00	Removal of residual root canal fillings	69
	67	1006.1998	\$205.00		69A
	68	1000.5821	\$205.00	Tip angle 117° - for large channels	70
	69	1006.2000	\$205.00		70A
SONICflex Retro		Retrograde root preparation with contra-angle, diamond-coated tips. Prevention of root damage from hairline cracks, support of tooth stability, protection of the periapical bone. Preparation along a precise axis.			
SONICflex Retro Tip Set A		Part No.	Price	Description	Tip No.
		1006.2034	\$969.00	Retro Tip Set	56, 57, 16, 17, 20, 21
		0571.0381	\$969.00	Retro Tip Set A	56A, 57A, 16A, 17A, 20A, 21A
SONICflex Retro tips		0571.5541	\$215.00	Cylindrical, small, left, Prep of channel	16
	16	1006.2036	\$215.00		16A
	17	0571.5581	\$215.00	Cylindrical, small, right, Prep of channel	17
	20	1006.2035	\$215.00		17A
	21	0571.5521	\$215.00	T-shape, small, left, Prep of the undercut	20
	55	1006.2037	\$215.00		20A
	56	0571.5561	\$215.00	T-shape, small, right, Prep of the undercut	21
	57	1006.2038	\$215.00		21A
	55	0571.7342	\$215.00	Anterior tooth	55
	56	1006.2039	\$215.00		55A
	56	0571.7322	\$215.00	Isthmus, left	56
	57	1006.2040	\$215.00		56A
	57	0571.7332	\$215.00	Isthmus, right	57
		1006.2041	\$215.00		57A

*All tips also work with SONICflex 2000L and 2000N.

Note: Tips ending in suffix 'A' are compatible with SONICflex 2008 / SONICflex 2008L only.

Minimally invasive caries therapy

SONICflex Microinvasive



Special tip system for small, difficult-to-access lesions in the proximal region. Tips diamond-coated on one side allow direct access to initial lesions from an occlusal or buccal direction. Maintenance of healthy tooth substance.

SONICflex SONICsys Microinvasive tip



Part No.	Price	Description	Tip No.
0571.6811	\$171.00	Small hemisphere, mesial	30
1006.1969	\$171.00		30A
0571.6801	\$171.00	Small hemisphere, distal	31
1006.1971	\$171.00		31A
0571.6831	\$171.00	Large hemisphere, mesial	32
1006.1973	\$171.00		32A
0571.6821	\$171.00	Large hemisphere, distal	33
1006.1975	\$171.00		33A

Accessories for SONICsys Micro Tip



1000.4887	\$75.00	For changing tips. Sterilizable. Torque feature prevents over tightening of tips. Can be used as protective cover when SONICflex is attached to delivery system.
0411.9101	\$98.00	Sterilization cassette

SONICflex Bevel



Tips diamond-coated on one side for the precise bevelling of prepared small lesions. Optimization of adhesion. Minimization of the marginal gap.

SONICflex Bevel Tip Set A

Part No.	Price	Description	Tip No.
1006.2030	\$215.00	Bevel Tip Set	58A, 59A

SONICflex Bevel Tips



0571.7392	\$109.00	Mesial Tip	58
1006.1988	\$109.00		58A
0571.7382	\$109.00	Distal Tip	59
1006.1990	\$109.00		59A

SONICflex Seal



Preparation of fissures. The perfect complement to the KaVo DIAGNOdent.

SONICflex Seal Tip



Part No.	Price	Description	Tip No.
1000.8323	\$105.00	Three tips for supra- and subgingival calculus removal.	45
1007.1506	\$105.00		45A

*All tips also work with SONICflex 2000L and 2000N.

Note: Tips ending in suffix 'A' are compatible with SONICflex 2008 / SONICflex 2008L only.




Minimally invasive caries therapy (cont.)

SONICflex Cariex D + TC



Minimally invasive removal of carious dentine. Cariex TC: Removal of large carious lesions. Cariex D: for small carious lesions with the ability to shape the cavity.

		Part No.	Price	Description	Tip No.
SONICflex Cari-ex D Tip Set		1000.7378	\$295.00	SONICflex cari-ex D Tip Set	42, 43
		1006.2025	\$295.00	SONICflex cari-ex D Tip Set A	42A, 43A
SONICflex Cari-ex TC Tip Set		1000.7379	\$295.00	SONICflex cari-ex TC Tip Set	71, 72
		1006.2025	\$295.00	SONICflex cari-ex TC Tip Set A	71A, 72A
SONICflex Cari-ex D + TC tips		1000.7105	\$165.00	Ball shaped for excavation 0.8mm	42
		1006.1980	\$165.00		42A
		1000.7167	\$165.00	Ball shaped for excavation 1.2mm	43
		1006.1981	\$165.00		43A
		1000.7362	\$165.00	Ball shaped for excavation 1.0mm	71
		1006.2002	\$165.00		71A
		1000.7363	\$165.00	Ball shaped for excavation 1.4mm	72
		1006.2004	\$165.00		72A

SONICflex Angle



Angled tips with diamond coated working surfaces (the ends are uncoated) for micro preparations (eg. tunnel preps) in difficult to access areas. Marginal ridge can be left intact.

		Part No.	Price	Description	Tip No.
SONICflex Angle Kit		0571.0341	\$245.00	Mesial, distal	53, 54

SONICflex prep gold



Preparation tips diamond-coated on one side for finishing the edges of cavities, specially developed for gold inlays and partial crowns. The attachment geometry for gold inlays and partial.

		Part No.	Price	Description	Tip No.
SONICflex prep gold Tip Set A		1006.2028	\$399.00	SONICflex Prep Gold Tip Set A	49A, 50A
SONICflex prep gold tips		0571.7212	\$202.00	Mesial Tip	49
		1006.1983	\$202.00		49A
		0571.7222	\$202.00	Distal Tip	50
		1006.1984	\$202.00		50A

*All tips also work with SONICflex 2000L and 2000N.

Note: Tips ending in suffix 'A' are compatible with SONICflex 2008 / SONICflex 2008L only.

Inlay / onlay preparation

SONICflex prep ceram



For ceramic inlays/onlays. For finishing the margins of cavities; directly transfers the geometry of the tips to the tooth substance. No more undercuts with laboratory-fabricated inlays and onlays.

SONICflex prep ceram Tip Set A



Part No.	Price	Description	Tip No.
1006.2029	\$399.00	SONICflex Prep Ceram Tip Set A	51A 52A

SONICflex prep ceram tips

51
52

0571.7262	\$202.00	Mesial Tip	51
1006.1985	\$202.00		51A
0571.7272	\$202.00	Distal Tip	52
1006.1986	\$202.00		52A

SONICflex prep CAD/CAM



Ideal shaping of the cavity margins for CAD/CAM systems such as KaVo EVEREST and CEREC.

SONICflex prep CAD-CAM insert-kit

Part No.	Price	Description	Tip No.
1002.1988	\$405.00	SONICflex CAD/CAM Insert Kit	34, 35

SONICflex prep CAD-CAM Tip Set A


1006.2024	\$405.00	SONICflex CAD/CAM Tip Set A	34A, 35A
-----------	----------	-----------------------------	----------

SONICflex retro tips

34
35

1002.1984	\$215.00	Mesial Tip	34
1006.1977	\$215.00		34A
1002.1986	\$215.00	Distal Tip	35
1006.1979	\$215.00		35A

SONICflex crown prep tips



For crown preparation. Used for finishing work on margins in crown preparations.

97
98

Part No.	Price	Description	Tip No.
1008.6383	\$215.00	Crown prep with round tip	97
1008.6384	\$215.00		97A
1008.6385	\$215.00	Crown prep with flat tip	98
1008.6386	\$215.00		98A

SONICflex cem



Gentle and precise insertion of inlays. Reduction of the marginal gap.

SONICflex cem tips



Part No.	Price	Description	Tip No.
0571.5431	\$140.00	Cem tip with white cap	12
1006.2021	\$140.00		12A
0571.7142	\$18.00	Cem Inlay top	

*All tips also work with SONICflex 2000L and 2000N.

Note: Tips ending in suffix 'A' are compatible with SONICflex 2008 / SONICflex 2008L only.



PROPHYflex 3 Black



Air cleaning system for removal of stains and bacterial plaque, cleaning of enamel surfaces prior to fissure sealing.

The PROPHYflex 3 features a shorter, more ergonomic handle; a lighter plastic reservoir and a new improved design. MULTIflex connection, sterilizable up to 275°F (135°C), 44 psi. Water 7.2 – 10 psi.

PROPHYflex 3 box includes:

- PROPHYflex 3
- 2 Plastic reservoirs
- 5 Packets of powder PROPHYpearls
- 1 Large cleaning drill
- 2 Small nozzle needles
- 1 Restrictor screw
- 2 Tips (one small, one large)

PROPHYflex 3 Violet



Part No.	Price	Description
1000.4672	\$1,699.00	PROPHYflex 3 Black
1006.9926	\$1,699.00	PROPHYflex 3 Violet



PROPHYpearls



Calcium carbonate powder eliminates salty taste. Smooth, soft & gentle to the tooth. Allows for efficient and faster cleaning. No clogging in the unit. Can be used with any manufacturers' air polishing system.

Part No.	Price	Description
1010.1798	\$119.00	PROPHYpearls Neutral (4 bottles)
1010.1826	\$150.00	PROPHYpearls Neutral (80 sticks)
1010.1828	\$150.00	PROPHYpearls Mint (80 sticks)
1010.1829	\$150.00	PROPHYpearls Peach (80 sticks)
1010.1830	\$150.00	PROPHYpearls Orange (80 sticks)
1010.1831	\$150.00	PROPHYpearls Berry (80 sticks)

PROPHYflex Clearing Powder



Part No.	Price	Description
1007.0014	\$195.00	PROPHYflex Powder Orange
1007.0015	\$195.00	PROPHYflex Powder Berry
1007.0016	\$195.00	PROPHYflex Powder Cherry
1007.0017	\$195.00	PROPHYflex Powder Mint
1009.3732	\$105.00	PROPHYflex Perio Powder



PROPHYflex System and Accessories

Tips for PROPHYflex 3



Focused spray

Wide spray

Part No.	Price	Description
0573.0151	\$209.00	Focused spray
0573.0181	\$209.00	Wide spray

Silicone Reservoir Cover



1000.2678 \$19.00 1 cover, for storing powder canisters

Powder Screw



0573.0412	\$40.00	For 100% powder flow Ceramic, Blue or Silver
0573.0002	\$40.00	For 70% powder flow Plastic, White or Bronze

Cleaning Tools



0573.0321	\$26.00	Long Cleaning Brush
0573.6052	\$8.25	Nozzle Needle (small) for cleaning tip

Parts / Accessories



0573.6152	\$38.00	Plastic powder reservoir, sterilizable
0573.6072	\$11.00	Black seal for reservoir
0200.6019	\$7.25	O-ring for tip (2 required)



**townie
choice**
awards 2005

**THE DENTAL
ADVISOR**
2005 Preferred Product



PROPHYwiz 181P Hygiene Handpiece

- Ultra light handpiece, minimizes hand fatigue
- Ergonomic and well-balanced
- Whisper quiet, smoothest running
- Autoclavable to 275°F (135°C)
- 3,100 rpm
- For metal Doriot attachments or disposable prophy angles
- 360° rotation at center
- 36-43 PSI
- 1-year warranty
- 2.7 oz

Part No. 1003.2278 \$865.00

Part No. 1003.9882 \$1,815.00 PROPHYwiz 3-pack

Part No. 0201.0116 \$21.00 Screw for Prophy Angle



**townie
choice**
awards 2005



19E / 19ES Hygiene Handpiece

- Reduction 7.4:1
- Speed range up to 5,000 rpm
- For disposable prophy heads only
- Use with 181K Air Motor (19E)
- Use with universal E-connection (19ES)

Part No. 1004.6403 \$519.00 19E

Part No. 1004.4503 \$519.00 19ES

Part No. 0201.0116 \$21.00 Screw for Prophy Angle

181K Air Motor

- Speed range up to 20,000 rpm
- 4-hole connection
- Non-optic
- Forward/reverse speed adjustment
- Sterilizable
- 36-43 PSI
- 1-year warranty
- 4.8 oz

Part No. 0535.5260 \$789.00

31E Low-speed Handpiece

- Reduction 4:1
- Speed range up to 5,000 rpm
- For use with snap-on / screw-in cups and brushes
- Extremely small head with 100° for superior access
- Durable hygienic coating with textured surface for easy grip
- Use with 181K Air Motor
- Comes with Prophy cup holder

Part No. 1003.6976 \$615.00

**townie
choice**
awards 2005



**townie
choice**
awards 2005

**THE DENTAL
ADVISOR**
2005 Preferred Product



**townie
choice**
awards 2005



Caries Detection | DIAGNOdent pen and Display



Shown: DIAGNOdent pen and Display Unit 2191.
Units sold separately.



DIAGNOdent pen

High Tech Laser Fluorescence Caries Detection Aid which measures changes in tooth structure caused by caries. It is especially useful for pit and fissure areas - even when the outer tooth surface seems to be intact. 90% accurate in detecting incipient caries. Only caries detection device backed by numerous clinical studies. Allows monitoring over time. Easy to use. Portable, handheld, pen-like device. Moves easily from room to room. Ideal for offices that already use the DIAGNOdent Classic and would like to have the added mobility of the hand-held DIAGNOdent pen. Weight 5 oz Length: 8-3/4". Display sold separately.

Includes:

- DIAGNOdent pen unit
- 2 Sterilizable tips for occlusal surfaces
- Sterilization cassette
- Sterilizable calibration standard
- 1 standard AA battery
- Clinical guidelines and operation manual
- Laminated quick reference guide

Part No. 1004.3400 \$2,849.00



DIAGNOdent display

With the new DIAGNOdent pen Display, we have enhanced the patient experience in caries diagnosis. The remote display for DIAGNOdent pen uses infrared signals, and the values are shown simultaneously on the DIAGNOdent pen and on the remote display. The display also serves as a holder for the DIAGNOdent pen when not in use. Pen sold separately.

Part No. 1004.8400 \$779.00



DIAGNOdent pen Perio Tip

DIAGNOdent now goes beyond caries detection to detect subgingival calculus for the treatment of periodontal disease. Sterilizable.

Perio Tip for detection of subgingival calculus in periodontal pockets.

Part No. 1004.0370 \$415.00

DIAGNOdent pen Accessories



Sterilization and storage cassette



Diagnodent pen Barriers (250 Barriers)

Part No.	Price	Description
1002.6967	\$165.00	Tip for occlusal surfaces.
1002.7003	\$83.00	Sterilizable grip sleeve
1002.7020	\$124.00	Removable "C" standard calibration disk
1002.7011	\$142.00	Sterilization storage cassette without calibration standard
1007.1480	\$34.00	DIAGNOdent pen Barriers (250 pack)
0906.7279	\$23.00	Patient Education Flyers (50 pack)

DIAGNOdent Classic

High Tech Laser Fluorescence Caries Detection Aid which measures changes in tooth structure caused by caries. It is especially useful for pit and fissure area – even when the outer tooth surface seems to be intact. 90% accurate in detecting incipient caries. Only caries detection device backed by numerous clinical studies. Allows monitoring over time. Easy to use.

Includes:

- DIAGNOdent Classic unit
- Sterilization Cassette
- Handpiece with Tubing
- 2 sterilizable tips - type A for occlusal surfaces
- Rechargeable Battery
- Battery Charger
- Sterilizable calibration standard
- Clinical guidelines and operation manual
- Laminated quick reference guide

Part No. 0574.0500

\$2,849.00

**DIAGNOdent Classic Perio Tip**

DIAGNOdent now goes beyond caries detection to detect subgingival calculus for the treatment of periodontal disease. Sterilizable.

Perio tip for detection of subgingival calculus in periodontal pockets.

Part No. 1004.1640

\$415.00

DIAGNOdent Classic Accessories

A Tip



B Tip



Snap in Ceramic standard



Magnetized hand-piece holder

DIAGNOdent Classic Barriers
(250 Barriers)

A Tip kit



Sterilization cassette without standard

Part No.	Price	Description
0574.1311	\$169.00	Type A Tip for occlusal surfaces
0574.1321	\$169.00	Type B Tip for flat surfaces
1002.8569	\$465.00	Type A Tip Kit Includes: 4 Type A Tips, Sterilization cassette without standard
0574.1271	\$56.00	Magnetized handpiece holder
1003.5528	\$70.00	Snap in Ceramic Standard "C" (for DIAGNOdents with serial numbers higher than 1500 000)
0574.1692	\$42.00	Snap in Ceramic Standard "B" (for DIAGNOdents with serial numbers smaller than 1500 000)
1000.1533	\$139.00	Sterilization cassette w/o standard
0574.1341	\$65.00	Mounting bracket (to mount DIAGNOdent to delivery system)
1007.1479	\$34.00	DIAGNOdent Classic Barriers (250 pack)
1007.5158	\$96.00	Rechargeable Battery
1007.5160	\$96.00	Battery Charger
1007.7209	\$149.00	Rechargeable battery plus charger (Replaces part # 0574.1261)
0574.1211	\$105.00	DIAGNOdent Classic Sheath
0906.7279	\$23.00	Patient Education Flyers (50 pack)



RONDOflex plus

RONDOflex plus 360 with water spray



REALITY
Four Star Award



Tips for RONDOflex plus 360



The RONDOflex plus is a portable, easy-to-use abrasion handpiece that adds water spray to air abrasion technology. The RONDOflex plus is an efficient, "snap-on" alternative to your existing air abrasion unit or small investment to get you started. Ideal adjunct to the DIAGNOdent Caries Detection Aid. MULTIflex connection, sterilizable up to 275°F (135°C), 46-55 psi.

Indications

- Preparation of small carious lesions
- Opening and extension of fissures

- Surface roughening of preparations and restorations
- Removes residual cement from crowns and bridges

Includes:

- RONDOflex plus
- 2 Tips (Both 110° tips)
- 1 Reservoir 27 micron
- 1 Reservoir 50 micron
- 1 - 27 micron powder (75g)
- 1 - 50 micron powder (75g)
- Small nozzle needle
- Large cleaning drill
- 2 Silicone reservoir covers

Part No. 1002.2179

\$2,979.00

Part No.	Price	Tip Angulation	Size
1002.9176	\$449.00	110°	0.018" small
1002.6251	\$449.00	110°	0.025" wide
1002.9182	\$449.00	90°	0.018" small
1002.9179	\$449.00	90°	0.025" wide
0200.6019	\$7.25		O-Rings (pack of 2)
1001.3254	\$265.00	120°	0.018" small
1000.5698	\$265.00	120°	0.025" wide
1001.3253	\$265.00	90°	0.018" small
1000.5697	\$265.00	90°	0.025" wide

* Discontinued model without water spray.

Tips for RONDOflex 2013*



	Part No.	Price	Description
Reservoir for RONDOflex	1003.1235	\$48.00	27 micron reservoir
		\$48.00	50 micron reservoir
Reservoir Seal	0573.6072	\$11.00	Black seal for reservoir
			
27 Micron aluminum oxide powder 1000g (2.2 lbs.). Recommended for intraoral use.	1000.5957	\$335.00	6 pack
			
50 Micron aluminum oxide powder 1000g (2.2 lbs.). Recommended for extraoral use.	1000.5956	\$335.00	6 pack
			
Silicone reservoir cover	1000.2678	\$19.00	1 cover. For storing powder canisters.
			
Powder screw	0573.0412	\$40.00	For normal flow powder (100%) Ceramic
			
Cleaning Tools	0573.0321	\$26.00	Long cleaning brush
			
	0573.6052	\$8.25	Nozzle needle (small) for cleaning tip
			
Tip Tools	1001.1206	\$23.00	For changing tips on RONDOflex
			



CORONAflex 2005 Crown and Bridge Remover



- Pneumatically operated instrument to remove crowns and bridges
- In most cases the restoration can be removed undamaged
- Excellent tool if patient needs recementation of a loose bridge anchor, endodontic treatment, etc
- Comes with detailed instructions; easy to use
- MULTIflex connection
- Operates between 35-60 psi
- 1 bar = 14.50 psi

Includes:

- CORONAflex 2005
- Crown forceps
- 2 small clamps
- 2 large clamps
- Loop holder
- 2 Loops

Part No. 0579.1000

\$2,179.00

CORONAflex Tools



Part No.	Price	Description
0579.0111	\$1,056.00	Crown forceps

Clamps



0579.0392	\$33.00	Camp Size: Small
0579.0402	\$33.00	Camp Size: Large

Loops holder and Loops



0579.0412	\$150.00	Loop holder
0579.0422	\$38.00	Loop



TryKaVo.com

Try the Best Products, Risk Free



Step 1.

CHOOSE



Choose a handpiece from
our collections listed

Step 2.

TRY




Try a handpiece risk free
for 5 business days

Step 3.

RETURN



After 5 day trial, return the
handpiece. Like it, buy it.
Contact 1-877-746-3995



Go to **TRYKAVO.COM** to
start your free and easy trial today!



KaVo SURGICAL SYSTEMS

Driven to Precision

The New KaVo Surgical Portfolio redefines the best. They are simple, intuitive and easy to use with large color screens and pictorial guides for step-by-step procedures.

The units can be used for all surgical applications and offer intelligent power and precise torque. KaVo redefines the best with the MASTERSurg and EXPERTsurg.



MASTERsurg LUX Wireless

Part No. 1009.1200

\$5,599.00

Includes: wireless foot control, instrument holder, INTRA LUX S600 motor, sterilizable motor cable, power cord, and 5 sterile irrigation tubes. Handpiece not included with MASTERsurg.

**More comfort features:**

- Simple, intuitive programming with exclusive glass panel touchscreen
- Choose a convenient location in your workplace thanks to wireless foot control
- Support procedure due to graphical real time display of torque and digital data

More you, every step of the way:

- 10 customizable programs to meet individual requirements
- Flexibility due to multi program mode to save settings

Maximum power and precision:

- One of the world's lightest, smallest surgical motors
- 80 Ncm torque at the instrument
- SMARTdrive technology for high torque and minimum vibration at the start
- One touch auto-calibration guarantees precise speed and torque for maximum safety

Technical Specifications:

- Maximum torque: 80 Ncm at the instrument
- Speed range: 0 - 40,000 rpm
- Power supply voltage: 100 V - 240 V
- Input frequency: 50/60 Hz
- Pump performance: 0 - 110 mL/min

Customizable Parameters:

- 10 programs, each with individually programmable treatment steps (max. 10)
- Contra-angle transmission ratio
- Maximum speed
- Maximum torque
- Pump delivery (0 - 4)
- Motor clockwise/anticlockwise
- Selection of implant position

EXPERTsurg LUX

Part No. 1008.3500

\$4,099.00

Includes: foot control, instrument holder, INTRA LUX S600 motor, sterilizable motor cable, power cord, and 5 sterile irrigation tubes. Handpiece not included with EXPERTsurg.

**Guaranteed precision, maximum safety:**

- One touch auto calibration guarantees precise speed and torque for maximum safety
- One of the World's lightest and smallest surgical micro motors
- Smooth running, low vibration

High reliability, top performance:

- 80 Ncm torque at the instrument
- SMARTdrive technology for high torque and minimum vibration at the start

Easy to use functions, saves time:

- Intuitive, straightforward programming
- Symbolic references and colour display make it simple to operate
- Easy to clean

Technical specifications:

- Input voltage: 100 V - 240 V AC
- Input frequency: 50 / 60 Hz
- Pump flow rate: 0 - 110 mL / min

6 adjustable parameters:

- Number of treatment steps (max. 10)
- Contra-angle transmission ratio
- Maximum speed
- Maximum torque
- Pump power (0 - 4)
- Motor forward/reverse



Surgical Handpieces KaVo SURGmatic instruments

NEW



KaVo SURGmatic handpiece 1:1

- S11 L (with light)
S11 C (without light)
- Transmission 1:1
- EXPERTgrip for better grip and hygiene
- For use with handpiece burs of 44.5 mm standard length and can also be used with contra-angle burs using the included backstop
- With external spray

Technical specifications (in connection with EXPERTsurg or MASTERsurg):

- Max. speed: 40,000 rpm
- Max. torque: 5.5 Ncm

S11L	Part No. 1009.1010	\$1,390.00
S11C	Part No. 1009.1005	\$1,275.00



Cellular Optic Attachments for Oral Surgery and Implantology CL 10 Straight Attachment

- 25,000 LUX glass rod
- For handpiece burs with standard length of 44.5mm
- With external spray attachment

- Strengthened drill chuck
- Speed range up to 40,000 rpm

Part No. 1002.0911	\$2,340.00
--------------------	------------



3610N1 Straight Attachment

- Transmission 1:1
- Speed range up to 40,000 rpm
- Easily dismantled for cleaning

- Reinforced chuck
- Maintain with closed chuck without bur

Part No. 0524.5600	\$1,499.00
--------------------	------------



3610N2 Surgical Straight Attachment

- Transmission 1:1
- For handpiece burs with special length 62.5 mm, with spray attachment and with locking device for alternative use of standard handpiece burs (44.5 mm)

- Can be easily dismantled
- Speed (max): 40,000 rpm Torque (max): 10 Ncm

Part No. 0524.5610	\$1,199.00
--------------------	------------



3610N3 Surgical Straight Attachment

- Transmission 1:1
- For handpiece burs with special length 70.0 mm with spray attachment

- Can be easily dismantled and features strengthened bur chuck
- Speed (max): 40,000 rpm Torque (max): 10 Ncm

Part No. 0524.5620	\$1,799.00
Part No. 0524.3062	Bur stop



NEW

KaVo SURGmatic contra-angle 20:1

- S201 L / XL (with light)
S201 C / XC (without light)
- S201 XL / XC - Hexagon clamping system and optimized head gearing and ball bearing for long service life (also compatible with burs/screw-in tools without hexagon clamping system)
- Reduction ratio 20:1
- Removable head, easy to clean
- EXPERTgrip for better grip and hygiene
- Small head with internal cooling system
- External cooling with coolant connector on the head

Technical specifications (in connection with EXPERTsurg or MASTERSurg):

- Speed: 15 - 2000 rpm
- Maximum torque: 55 Ncm

S201 L Part No. 1009.0470

\$1,559.00

S201 XL Part No. 1010.2332

\$1,559.00

S201 C Part No. 1009.1100

\$1,365.00

S201 XC Part No. 1010.2333

\$1,365.00



Cellular Optic Attachments for Oral Surgery and Implantology CL 3-09 Contra Angle Attachment

- 25,000 LUX glass rod
- Head and shank can be dismantled
- Reduction 27:1
- For latch burs
- Speed range 11 to 1500 rpm
- Small head with internal cooling system
- External cooling medium connection attached to shank

Part No. 1002.3693

\$2,599.00



CL 09 Cellular Optic Contra Angle Attachment

- Transmission 9:1
- Speed range up to 4,444 rpm
- Ideal to use with INTRA CL3 head (3:1) (speed range 11 to 1,500 rpm)
- Extremely strong reduction gears

CL 09 Part No. 1001.9258

\$1,499.00

C 09 Part No. 1001.2278

\$925.00



3620N Contra Angle Attachment

- Transmission 1:1
- Speed range up to 40,000 rpm
- Easily dismantled for cleaning

Part No. 0534.5650

\$765.00



3624N Contra Angle Attachment

- Reduction 4:1
- Speed range up to 10,000 rpm
- Easily dismantled for cleaning

Part No. 0534.5660

\$985.00



3630N Contra Angle Attachment

- Reduction 30:1
- Speed range up to 1,333 rpm
- Easily dismantled for cleaning

Part No. 0534.5670

\$1,201.00



Surgical Handpieces and Accessories



67RIC Reduction Latch Head for Internally Cooled Burs

- Reduction 4:1
- Internal/external coolant
- Reinforced gearing
- Max torque 55 Ncm

Part No. 0540.8340 \$725.00



C3 Reduction Latch Head for Internally Cooled Burs

- Reduction 3:1
- Internal/external coolant
- Small head diameter
- Ideal to use with CO9 contra-angle attachent (pg. 72)
- Max torque 43 Ncm

Part No. 1001.4503 \$765.00



CL 3 Reduction Latch Head

- For use with CL09 cellular optic Contra Angle
- For internally cooled burs
- Reduction 3:1
- Internal/external coolant
- Small head diameter
- Max torque 43 Ncm

Part No. 1001.8993 \$1,039.00

Part No. 0540.3711 \$63.00 Clip for Water filter

INTRAsurg 300 / 300 plus Accessories and Replacement Parts

Part No.	Price	Description
1001.4968	\$5,299.00	Steri-set fiber optic surgery motor INTRA-SL 550 Includes: INTRA-SL 550 fiber optic motor and internal irrigant tubing, motor support
1001.3421	\$4,299.00	INTRA-SL 550 fiber optic surgery motor only
1002.2928	\$105.00	INTRA-SL 550 bulb. Autoclavable bulb for INTRA-SL 550 fiber optic surgical motor
1001.4953	\$3,599.00	Steri-set non-optic surgery motor INTRA-S 550. Includes: INTRA-S 550 non-optic motor and internal irrigant tubing, motor support
1000.8072	\$2,199.00	INTRA-S 550 non-optic surgery motor only
1001.3347	\$139.00	Sterilizable pump tubing (Bag to handpiece tubing)
1001.9902	\$225.00	Sterilizable tube fittings. 10 each Single use
1001.2651	\$1,699.00	Autoclavable tubing (Complete set)

MASTERsurg / EXPERTsurg Accessories

Part No.	Price	Description
1009.8757	\$131.00	Set of 10 sterile tubes
1009.3411	\$79.00	Instrument holder
1008.8000	\$1,649.00	INTRA LUX S600 motor
1009.1700	\$575.00	S600 Motor Cable
1002.3377	\$85.00	Spray clip and internal tube for S201 series
0593.0361	\$72.00	Y-distributor for internally cooled burs
1001.3436	\$2.75	70 mm tubing for use with Y-distributor
1010.9727	\$1,059.00	MASTERsurg travel case
1010.8009	\$1,059.00	EXPERTsurg travel case

INTRA Surgical Head Accessories

Part No.	Price	Description
0540.1971	\$65.00	Internal spray line for 67IC
0593.0361	\$72.00	Y-Distributor for 67IC/67RIC/C3/CL03
1001.3436	\$2.75	External spray tubing for heads with external coolant. Clear sterilizable tubing.
0540.3711	\$63.00	Spray tube for 67RIC
0065.5279	\$41.00	Irrigation tubing to attach to the exterior of the electric motor tubing (INTRAsurg 500). Can be cut to desired length. Sold by meter.

GENTLEray 982 Laser Whitening Handpiece

Sterilizable handpiece expands chairside treatment options; creates a uniform parallel laser beam profile (spot diameter of approx. 5 mm) that accelerates the whitening gel without causing unnecessary heating. Available as option for both GENTLEray 980, Classic and Premium lasers.

Part No. 1005.2895 \$1,649.00

**GENTLEray 980 Accessories**

Part No.	Price	Description
1005.2894	\$325.00	Laser Safety goggles Value Line
1005.2896	\$575.00	Whitening Fiber
1006.0031	\$225.00	Straight Cannula Premium 3/pk
1006.0907	\$225.00	30° Cannula Premium 3/pk
1006.0910	\$225.00	70° Cannula Premium 3/pk
0065.5279	\$41.00	Handpiece tubing per meter
1006.0500	\$8.25	Instrument Tubing
1004.3926	\$85.00	Extra Keys
0236.2394	\$123.00	Interlock
1006.8257	\$684.00	981 Premium water cooling handpiece set
1006.3976	\$125.00	Locking mechanism for foot pedal
1006.4428	\$505.00	300 micron laser fiber, 3M (red)
1006.5431	\$505.00	200 micron laser fiber, 3M (blue)
1005.2887	\$56.00	Straight Cannula Classic 3/pk
1005.2888	\$56.00	30° Bended Cannula Classic 3/pk
1005.2889	\$56.00	70° Bended Cannula Classic 3/pk
1004.9800	\$999.00	981 Classic Laser Handpiece Set
1005.3955	\$395.00	Laser safety glasses XC with insert for glasses
1005.3609	\$85.00	Ceramic scissors
0236.2439	\$8.25	Laser warning sign
1005.2885	\$1,799.00	Foot switch

INTRAsert 905 Accessories

Part No.	Price	Description
0641.7301	\$335.00	Water flush assembly
0516.3011	\$125.00	Connected pump tubing. Grey.
0065.5389	\$19.00	Inner tube - 2x1. Sold by meter.
1001.3436	\$2.75	Handpiece tubing. Sterilizable. Clear.
0516.5512	\$22.00	Intermediate piece
1002.3569	\$170.00	INTRAsurg calibration

Attachment Accessories for ENT Surgery

Part No.	Price	Description
0549.0020	\$325.00	Surgical attachment for use with 8A contra angle attachment. Accepts regular length handpiece burs. Bur tool sold separately.
0411.3521	\$215.00	Bur tool for surgical attachment
0524.1531	\$62.00	Spray collar (3554, 3555, 3556)
0524.3062	\$33.00	Bur stop (3555, 3556)
0410.1963	\$27.00	J Hook for bur stop
0410.0941	\$19.00	Nozzle wire



In-Office Laboratory

K-ERGOgrip Installation



Part No.	Price	Description
1005.4207	\$2,459.00	K-ERGOgrip/K-Control TLC-Install. EWL 4956/4944 foot control
1005.4206	\$2,459.00	K-ERGOgrip/K-Control TLC-Install Knee EWL 4955/4944
1005.4204	\$2,459.00	K-ERGOgrip/K-Control TLC-Install. EWL 4957/4944 Bench

K-POWERgrip Installation



1007.2082	\$2,498.00	K-POWERgrip/K-Control TLC Knee Control Installation 4956/4941
1007.2084	\$2,498.00	K-POWERgrip/K-Control TLC Foot Control Installation 4955/4941
1007.2083	\$2,498.00	K-POWERgrip/K-Control TLC On-bench Control Installation 4957/4941

K5 Plus Installation



1009.4244	\$1,625.00	K5 Plus/K-Control TLC Knee Control Installation 4956/4911
1009.4245	\$1,625.00	K5 Plus/K-Control TLC Foot Control Installation 4955/4911
1009.4243	\$1,625.00	K5 Plus/K-Control TLC On-bench Control Installation 4957/4911

SMARTair Mobile

Portable unit with 2 filter bags.



Part No.	Price	Description
0658.2160	\$175.00	Pack of 5 filter bags
1001.8196	\$2,199.00	SMARTair Mobile

Please note:

This is only a limited selection of our wide range of laboratory products.

For complete in-office laboratory set-up solutions, call KaVo at 888.ASK.KAVO for pricing.

LCM Illuminator



- Compact bulb-in-coupler system
- For foot pedal or hanger installation
- Operates up to four handpieces
- Includes module, transformer, remote intensity control, and installation kit
- LCM tubing and MULTIflex LUX (465LRN / 465LED or 459AN, 460LE / 460LED) couplers ordered separately

Part No. 0905.8400

\$645.00

LS2 Illuminator



- Compact remote fiber optic light source for hanger or foot pedal installation
- Operates up to three handpieces
- High/low intensity switch
- Installation kit included
- LS2 tubing and MULTIflex LUX (1390) coupler ordered separately

Part No. 0905.8500

\$945.00

LCM Tubing—6-Pin Type Tubing

Part No.	Price	Color	Description
0905.4113	\$415.00	Gray	Straight 11' LCM tubing

LS2 Tubing—Standard 5-Hole Type Tubing

Part No.	Price	Color	Description
0905.5793	\$515.00	Gray	Straight 7' x 9'
0905.5243	\$515.00	Gray	Straight 12' x 14'

Parts and Accessories for LCM

Part No.	Price	Description
1006.6225	\$155.00	Transformer (120 Volt AC)
0905.8405	\$54.00	Low voltage cable (15 Volt AC)
0905.8406	\$109.00	LCM remote intensity control
0905.8440	\$49.00	LCM installation kit
1002.2928	\$105.00	Halogen Light Bulb (pack of 2)
1007.5372	\$125.00	KaVo MULTI LED Bulb
0905.8411	\$33.00	Voltage clip
0905.8403	\$555.00	LCM module only (3.6 Volt DC)

Parts and Accessories for LS2

Part No.	Price	Description
0905.8510	\$84.00	Halogen 12 Volts / 100 Watt bulb
0905.8540	\$76.00	Installation kit
0905.8501	\$599.00	LS2 module only
0905.8520	\$62.00	Connector cable assembly
0513.0812	\$41.00	Coupler nut (protects crimping)
0905.8550	\$67.00	Time delay switch kit



KaVo Instruments and Their Successors

GENERAL TYPE / FUNCTION		HISTORY OF KaVo INSTRUMENTS AND THEIR SUCCESSORS						CURRENT SERIES
HIGH SPEED AIR DRIVEN								CURRENT
Standard Head, MULTIflex, Optic	647B	649B	6000B	6500B	6500BR	M8900 L		MASTER series
Mini Head, MULTIflex, Optic		632B	642B	4500B	4500BR	M8700 L		MASTER series
Standard Head, MULTIflex, Optic	630B	640B	650 B	655 B	659B	E679L		EXPERT series
Mini Head, MULTIflex, Optic		639 B	642 B	645B	635 B	E675L		EXPERT series
Mini Head, Midwest Stylus coupling, Optic						E675 LMW		EXPERT series
Standard Head, MULTIflex, Optic					646B	S619 L		SMART series
Standard Head, MULTIflex, Non Optic						S619 C		SMART series
Mini Head, MULTIflex, Optic						S615 L		SMART series
Standard Head, 4-hole fixed back, Non Optic		625CD	625C	607	636CP	S609 C		SMART series
Mini Head, 4-hole fixed back, Non Optic				634 P4	634 F4	S605 C		SMART series
E-STYLE HIGH AND LOW SPEED								CURRENT
High-speed Contra angle internal spray, internal Motor, with light					200XD	200XDR		MASTER series
High-speed Contra angle internal spray, with Light			25 LP	25 LPA	25LPR	M25 L		MASTER series
Shank, with light		20 L	20 LN	20 LH	20 LP	M20 L		MASTER series
Shank, with light		7 L	7 LN	7 LH	7 LP	M07 L		MASTER series
Shank, with light		29 L	29 LN	29 LH	29 LP	M29 L		MASTER series
Straight Handpiece with Light			10 LN	10 LH	10 LP	M10 L		MASTER series
High-speed Contra angle, internal spray, with Light				24 LN	25 LCA	E25 L		EXPERT series
High-speed Contra angle, internal spray, with Light			25 LN	25 LH	25LHC	E25 L		EXPERT series
High-speed Contra angle, internal spray, no light		24 B	24 C	24CN	25CHC	E25C		EXPERT series
High-speed Contra angle, internal spray, no light	25 B	25 C	25CN	25CH	25CHC	E25 C		EXPERT series
1:1 Low-speed, internal spray, with light						E20 L		EXPERT series
1:1 Low-speed, internal spray, without light	20 B	20 C	20 CN	20 CH	2068CHC	E20 C		EXPERT series
Reducing Low-speed, internal spray, with light						E15L		EXPERT series
Shank internal spray without light	7 B	7 C	7 CN	7 CH	0767CHC	E15 C		EXPERT series
Shank internal spray without light	29 B	29 C	29 CN	29 CH	2967CHC	E15 C		EXPERT series
						E61C		EXPERT series
						E31C		EXPERT series
Handpiece internal spray without light	10 B	10 C	10 CN	10 CH	10 CHC	E10 C		EXPERT series
E-STYLE HIGH AND LOW SPEED								CURRENT
Shank without spray, Latch					20 K	20E		SMART series
Shank without spray, FG						80E		SMART series
Straight Handpiece without spray					10 K	10E		SMART series
Handpiece without spray					14 K	19E		SMART series
Shank without spray					29 K	14E		SMART series
High-speed Contra angle external spray		23	23 A	23D	23ES			SMART series
Shank external spray				7	7 A	14ES		SMART series
Shank external spray				20	20 A	20ES		SMART series
Shank external spray				29	29 A	14ES		SMART series

Type	Part No.	Page No.	Retail
200XDR	1007.3570	10	\$1,699.00
M25 L	1009.3630	12	\$1,669.00
M05 L	1009.3640	12	\$1,669.00
M07 L	1009.3610	12	\$1,175.00
M10 L	1009.3570	12	\$1,385.00
M20 L	1009.3620	12	\$960.00
M29 L	1009.3580	12	\$1,319.00
L68B	1008.1834	14	\$595.00
L80	1004.4456	14	\$595.00
L52	1008.1825	14	\$595.00
L62	1008.1830	14	\$615.00
L66B	1008.1831	14	\$615.00
L67	1008.1832	14	\$615.00
L53	1008.1826	14	\$575.00
L3	1008.1822	15	\$575.00
L3Y	1008.1823	15	\$575.00
L31	1008.1824	15	\$320.00
L61G	1008.1828	15	\$575.00
L61R	1008.1829	15	\$575.00
L22	1008.1835	15	\$575.00
E25L	1007.5550	19	\$1,429.00
E25C	1007.5551	19	\$1,225.00
E20L	1007.5540	19	\$1,359.00
E20C	1007.5541	19	\$969.00
E10C	1007.5560	20	\$819.00
E15L	1007.5530	20	\$1,185.00
E15C	1007.5531	20	\$969.00
E61C	1008.4999	20	\$1,035.00
E31C	1008.4998	20	\$1,125.00
19ES	1004.4503	21	\$519.00
23ES	1004.5380	21	\$999.00
80ES	1004.4375	21	\$619.00
M8700L	3001.0000	26	\$1,599.00
M8900L	1008.7800	26	\$1,599.00
E675L	1006.6400	27	\$1,195.00
E679L	1006.9200	27	\$1,195.00
E675LMW	1006.0800	27	\$1,195.00
S619L	1008.1641	28	\$749.00
S619C	1008.1642	28	\$669.00
S615L	1008.1643	28	\$749.00
S609C	1008.1644	29	\$569.00
S605C	1008.1645	29	\$569.00
S459L	3000.4808	29	\$995.00
S459C	3000.5063	29	\$845.00

Type	Part No.	Page No.	Retail
2307LN	0535.6190	30	\$1,699.00
2313LN	0535.6220	30	\$1,359.00
2320LN	0535.6180	30	\$1,599.00
2325LN	0535.6210	30	\$1,999.00
457	0553.1210	31	\$265.00
460LE	1001.7599	31	\$323.00
460LED	1007.3201	31	\$372.00
465LED	1007.3197	31	\$472.00
465LRN	0553.1550	31	\$415.00
1390	0553.1390	31	\$359.00
10E	1003.6977	32	\$550.00
19E	1004.6403	32	\$519.00
20E	1003.1107	32	\$619.00
31E	1003.6976	32	\$615.00
14E	1003.1108	32	\$569.00
80E	1004.6401	32	\$619.00
181L	0535.5210	33	\$1,959.00
181H	0535.5200	33	\$1,720.00
181K	0535.5260	33	\$789.00
QUATTROcare Plus	1008.3805	48	\$3,259.00
2003L	1000.8333	53	\$2,029.00
2003	1000.8334	53	\$1,899.00
2008	1007.1606	53	\$1,899.00
2008L	1007.1604	53	\$2,029.00
PROPHYflex	1009.9926	60	\$1,699.00
181P	1003.2278	63	\$865.00
CORONAflex	0579.1000	68	\$2,179.00
S11L	1009.1010	72	\$1,390.00
CL10	1002.0911	72	\$2,340.00
3610N1	0524.5600	72	\$1,499.00
3610N2	0524.5610	72	\$1,199.00
3610N3	0524.5620	72	\$1,799.00
S201L	1009.0470	73	\$1,559.00
CL3-09	1002.3693	73	\$2,599.00
CL3	1001.8993	73	\$1,039.00
CL09	1001.9258	73	\$1,499.00
C09	1001.2278	73	\$925.00
3620N	0534.5650	73	\$765.00
3624N	0534.5660	73	\$985.00
3630N	0534.5670	73	\$1,201.00
C3	1001.4503	74	\$765.00
67RIC	0540.8340	74	\$725.00





KaVo Handpiece Maintenance Guide

For the reprocessing of KaVo air driven highspeeds, Electric highspeeds, lowspeed contra angles, and SONICflex™ scalers using KaVo products and methods.

Basic Preparations



To minimize the risk of infection during reprocessing, always wear protective gloves.



Always remove the bur from the chuck system first.



Remove the handpiece from either the MULTIflex™ coupler or the motor - pull straight off - do not twist or bend. Some fixed back Handpiece require screwing off from the treatment unit tubing.



Manually clean the exterior of debris using a soft bristle brush under running tap water.



Do not submerge a Handpiece in water; do not place KaVo Handpiece in disinfectant solutions; do not clean Handpiece in ultrasonic devices.



Remove the heads from the attachment (if applicable) for separate care. Surgical instruments may be partially disassembled for cleaning and care (See instructions for use).

Internal Cleaning, Lubrication and Purge



Manual Procedure using KaVo Spray

KaVo recommends servicing the handpiece as part of the reprocessing after each use (i.e. after use in operatory and before sterilization). Always remove the bur prior to reprocessing.

NORMAL CLEANING AND LUBRICATION:

1. Place the appropriate nozzle on the KaVo spray can
2. Place the handpiece on the spray nozzle
3. Wrap the head of the handpiece in a soft cloth (either cloth or paper towel) or use KaVo's Cleanpac to contain excess lubricant.
4. Press the spray nozzle for 2 seconds or until lubricant is visibly expelled from the head of the handpiece
5. Place the Handpiece on an air source and purge the Handpiece for 30 seconds.
6. Wipe any excess oil off the outside of the handpiece
7. The Handpiece is now ready for sterilization

CHUCK CARE

KaVo recommends caring for the Handpiece chuck once per week

1. Place the Quick Nozzle on the KaVo Spray can
2. Ensure the small tip/nozzle is used
3. Place the opening for the bur on the small nozzle
4. Press the spray button for one second
5. Purge the Handpiece for 10 seconds



Sterilization



KaVo product with the sterilization symbol may be sterilized in steam autoclaves per EN13060/ISO 17665-1. KaVo MULTIflex couplers are not sterilizable.

The KaVo medical device has a maximum temperature resistance up to 138° (280.4°F). Depending on the available autoclave, select a suitable procedure from the following sterilization processes:

- Autoclave with three times initial vacuum: at least 3 minutes at 134 °C -1 °C / +4 °C (273 °F -1.6 °F / +7.4°F)

- Autoclave using the gravitation method: at least 10 minutes at 134 °C -1 °C / +4 °C (273 °F -1.6 °F / +7.4°F)
- Autoclave using the gravitation method: at least 60 minutes at 121 °C -1 °C / +4 °C (250 °F -1.6 °F / +7.4°F)
- Use in accordance with the sterilizer manufacturer's Instructions for Use

Terms of Sale

Pricing:

Prices are subject to change without notice. Possession of list price does not constitute an offer to sell.

Terms:

Payments are due 30 days from date of invoice. A 1% per month service charge (12% annually) will be applied to accounts past due. KaVo reserves the right to change this percentage based on inflation rates.

Minimum Order:

Orders totaling less than \$50.00 wholesale cannot be economically processed and will be subject to a \$7.00 service charge. Repairs are exempt from the minimum order requirement.

Shipment:

All shipments are FOB Factory unless otherwise stated and will be sent the most economical way unless otherwise specified. Shortage claims must be made within 10 days from receipt of goods. Damaged goods must be reported to the carrier within 24 hours and damage claims must be filed by the consignee.

Drop Shipment Fees:

All shipments of product from our warehouse that bypass the normal distribution channel will be subject to a drop shipment fee and freight charges.

Returns:

Prior approval must be obtained before merchandise is returned for credit or exchange. Merchandise that is less than nine (9) months from date of sale and still manufactured will be accepted for credit. Shipment of returned merchandise is customer's responsibility. Merchandise out of warranty may be returned for repair, not for credit. KaVo will not accept returns on products with retail price less than \$25.00.

Merchandise will be credited at the lower of market or purchase price. Merchandise returned in the original, unopened container and accompanied by the invoice number and date of purchase will be subject to a minimum 15% (\$5.00 minimum) restocking charge in addition to any charges for reconditioning or repackaging. Used products are not acceptable for return.

Product Changes:

KaVo Dental Corporation reserves the right to make changes in product design without the obligation to accept previously manufactured products for alteration, exchange or credit.

Warranty:

KaVo Dental Corporation warrants its products to the original purchaser against defective workmanship and materials under normal procedures of installation, use and service within the dental profession for a period of 90 days from the date of sale. Selected KaVo Dental Corporation products contain further written warranty terms. Warranty period begins no later than dealer invoice date.

KaVo dealers will sell products within a geographic location that the dealer can adequately service by making calls on current and potential customers.

Training Policy:

Dealer shall employ and maintain an adequate, knowledgeable, full time sales staff who regularly contact users and potential users within the territory. Dealers will utilize service and installation staff who are technically trained in the installation of products within the territory.

Dealer's staff shall possess and maintain adequate knowledge concerning the products and their uses, features and benefits. Dealer will cooperate with KaVo in conducting meetings for training of dealer personnel for the above mentioned reason.

Hazmat Fees:

Hazardous transport fees by carriers will be passed along to customer. Contact Customer Service for economic shipping solutions on all hazardous goods.

Registered Trademarks

The following product names are registered trademarks used by KaVo group and are therefore protected by the trademark license:

BELLAtorque	GENTLEforce	KaVo Care	PHYSIO
BELLAtorque LUX	GENTLEpower	KaVo spray	PROPHYflex
COMFORTbase	GENTLEpulse	KaVo SPRAYrotor	PROPHYpearls
COMFORTdrive	GENTLEray	LUX 2	PROPHYwiz
COMFORTtronic	GENTLEmini	LUX 3	QUATTROcare
Companion	GENTLESilence	MASTERmatic	RONDOflex
CORONAFlex	INTRA	MASTERmatic LUX	ROTAspray
DiaLUX	INTRAccompact	MASTERtorque	SAFEdrive
DIAGNOdent	INTRAflex	MASTERsurg	SMARTair
ELECTROtorque	INTRAflex LUX	MINItech	SMARTtorque
EXPERTgrip	INTRALUX	Mira	SONICflex
EXPERTmatic	INTRAmatic	MULTIflex	SONICflex LUX
EXPERTsurg	INTRAmatic LUX	MULTIflex LUX	SONICSYS
EXPERTtorque	INTRAsert	MULTI LED	SURGmatic
EVEREST	INTRAsurg	Opti-Torque	SURGtorque



Warranty

Warranty

For the periods of time specified below, KaVo guarantees that the following products, when furnished by KaVo in new and unused condition, will not fail during normal service due to defects in material and workmanship that are not due to apparent abuse, misuse or accident:

Cellular Optic Glass Rod in all High-speed, Low-speed & SONICflex	5 years	MULTIflex Couplers	1 year
COMFORTdrive 200 XDR	2.5 years*	Master Series (M8900 L and M8700 L)	2.5 years*
COMFORTtronic and COMFORTbase (tubing)	3 years	Expert Series (E679L, E675LMW, E675L)	2.5 years*
DIAGNOdent / DIAGNOdent pen (except tips, tubing, and display)	1 year	Smart Series (646B, S619L, S619C, S615L, S609C, S605C, S459L, S459C)	1 year
ELECTROtorque TLC / ELECTROtorque plus (control box, tubing, motor)	3 year	PROPHYflex 3	1 year
GENTLEray 980 (Classic or Premium)	2 years	PROPHYwiz	1 year
INTRAflex LUX	1 year	QUATTROcare / QUATTROcare Plus	1 year
MASTERsurg / EXPERTsurg and surgical handpieces	1 year	Repairs on Low-speed Heads, Attachments and Speed Increasing Attachments (if repaired at an authorized KaVo repair center)	1 year
KaVo LCM, KaVo/Vicon Light Sources	1 year	RONDOflex plus	1 year
MASTERmatic, EXPERTmatic, INTRAflex heads	2.5 years*	SONICflex / SONICflex LUX	1 year
Low-speed Motors	1 year	Tubings (all other types)	1 year
		Turbines (if handpiece repaired at an authorized KaVo repair center or repaired in-office)	1 year

Additional Warranty terms apply. See product warranty card.

* 2.5 year warranty if maintained with KaVo QUATTROcare Plus, if not maintained with QUATTROcare Plus the standard warranty of 2 years applies.

Dental Advisor Awards



COMFORTdrive
QUATTROcare Plus



635B Handpiece
CORONAflex
ELECTROtorque plus
PROPHYwiz
QUATTROcare
GENTLEsilence 6500BR



COMFORTdrive

Townie Awards



Air Driven High-speed Handpieces - winner 2015, 2013, 2012, 2011, 2010, 2009, 2008, 2007 and 2006

ELECTROtorque TLC - winner 2011, 2010, 2009, 2008, 2007 and 2006

ELECTROtorque plus - winner 2011, 2010, 2009, 2008, 2007, 2006, 2005, 2004 and 2003

DIAGNOdent pen - winner 2013, 2012, 2011, 2010, 2009, 2008, 2007 and 2006

DIAGNOdent - winner 2013, 2012, 2011, 2010, 2009, 2008, 2007, 2006, 2005, 2004 and 2003

Cellular Optic High-speed Handpieces - winner 2010, 2005, 2004 and 2003

Air Driven Low-speed Systems - winner 2005 and 2004

Electric Driven Low-speed Handpieces - winner 2015, 2013, 2012 and 2010

QUATTROcare - winner 2005, 2004 and 2003

PROPHYwiz - winner 2004

SONICflex LUX - winner 2004

GENTLEpower 25 LPA / 25 LPR - winner 2013, 2012, 2011, 2010, 2008 and 2009

MASTERmatic - winner 2015

Reality Awards



GENTLEsilence 6500B
PROPHYflex 3
DIAGNOdent
K1 System
SONICflex LUX 2003L
CORONAflex
ELECTROtorque plus
GENTLEpower 25LPA
INTRAFlex LUX
635B Handpiece
QUATTROcare
COMFORTdrive



PROPHYwiz 181P
DIAGNOdent pen
RONDOflex plus



DIAGNOdent pen



ELECTROtorque plus
GENTLEpower 25LPA

Dentistry Today Award



DIAGNOdent pen

Dental Products Report



COMFORTdrive





TryKaVo.com

Try the Best Handpieces, Risk Free



Step 1.

CHOOSE



Choose a handpiece from
our collections listed

Step 2.

TRY



Try a handpiece risk free
for 5 business days

Step 3.

RETURN



After 5 day trial, return the
handpiece. Like it, buy it.
Contact 1-877-746-3995



Go to **TRYKAVO.COM** to
start your free and easy trial today!

KaVo Dental Corporation
11729 Fruehauf Drive
Charlotte, NC 28273
888.ASK.KAVO (888.275.5286)
Fax 704.927.0816
www.kavousa.com

©2016 KaVo Dental Corporation. 906.8970 Rev 07 01.16



KaVo. Dental Excellence.