



IMAGING | UTILITY ROOM | MERCHANDISE



DENTAL
CATALOG
2015

AIR
TECHNIQUES

equipped for life®

Equipped for Life® with Air Techniques

Dear Valued Customer:

At Air Techniques, we are focused on high quality products, innovative new technology and unparalleled customer satisfaction – all to support your dental practice. We are proud that in our 50+ years of service to the dental profession, our company has attracted the best in design, engineering and customer care professionals. We are especially proud of our dealer partners who choose to support Air Techniques - companies with the best in sales, installation, and service support.

2014 witnessed Air Techniques expand its digital product line with the ProVecta S-Pan, bringing a new level of image sharpness in dental radiography. Building upon our experience with diagnostic image quality, the ProVecta S-Pan delivers a full adult panoramic X-ray in 7 seconds and offers 17 image programs including 4 specifically for children all with reduced radiation exposure. Patients and practitioners both appreciate the comfort and convenience of the touch screen control, optional voice commands, accommodating height adjustment and assured positioning with 3 laser guides. With the foundation of the digital clarity that the S-Pan delivers, Air Techniques is looking to build an extended offering in dental radiography with an emphasis on enhanced image quality and workflow with other new product introductions this year.

With respect to our practice supply/utility product lines, 2015 will bring new products that provide enhanced diagnostic and control capabilities, enabling you to better manage your dental air and vacuum needs, and will help the environment with superior air, water and waste management.

Also in 2015, we will further expand our Monarch® product portfolio, a line of products that help you better manage your patients and your team. Monarch includes products that are easier to use and cost effectively clean and disinfect surfaces; that enable cleaning of instruments prior to sterilization; that remove build-up in the lines of dental units and evacuation systems to maximize the performance of your equipment; and that care for your team, leaving skin and hands smooth, silky and healthy.

Today, more dental offices choose Air Techniques air and vac systems than all other brands combined. In North America, 9 out of 10 offices who choose flexible phosphor sensors, choose ScanX (and in Europe and Asia, more offices use this technology than wired sensors – ask us why?).

Thank you for your support and trust in Air Techniques.



Christoph Roeer
CEO



CONTENTS

Imaging

Panoramic X-ray	4
Intraoral X-ray	6
Digital Radiography Systems	8
Intraoral Cameras	12
Caries Detection Aid	16
Cam-X™ Spectra & Cam-X™ Polaris Clinical Images	18
Imaging Software	19
Film Processors	20

Utility Room

Air Compressors	24
Dry Vacuums	30
Wet Vacuums	34
Amalgam Separator	36
Utility Accessories	37
Utility Room Configuration	42
Clinical Utility Packages	43

Merchandise

Imaging Accessories	46
Chemistry	50
Processor Accessories	51
Surface Disinfectant	54
Instrument Cleaner	56
Hand Sanitizer & Hand Lotion	58
Evacuation System Cleaner	60
Water Line Cleaner	61

Warranty & Return Policy

Professional Sales Associates

62

68



“

ProVecta S-Pan's ease of use, high quality images and quick image acquisition helps to reduce stress for me, my staff and my patients.”

– Melissa Rosenfeld, D.D.S.



IMAGING

Panoramic X-ray

Intraoral X-ray

Digital Radiography

Intraoral Cameras

Caries Detection Aid

Imaging Software

Film Processors


IMAGING

ProVecta® S-Pan Panoramic X-Ray


ProVecta S-Pan
PN A7350



Optional Base Plate
PN A7355



Selecting from 20,000 image fragments to create one exceptional image



ProVecta S-Pan Panoramic X-ray: Panoramic Perfection.

Image Quality:

- S-Pan technology for sharper images and to make diagnosis easier.
- Adult Panoramic images within 7 seconds.
- CsI sensor for better image quality and reduced radiation.

User-friendly Design:


- 7" touch-display for intuitive operation.
- Small footprint, only 39" x 44" x 90".
- TWAIN compliant imaging software.

Patient Comfort:

- 17 X-ray programs to make diagnosis easier, including 4 child modes.
- 3 positioning laser guides for error-free, ergonomic face-to-face positioning.
- Long vertical path accommodates tall patients and wheelchair users.




Child OPG image



Bitewing image



Sinus image




TMJ image

PART NO.	DESCRIPTION	PRICE
A7350	ProVecta S-Pan Panoramic X-ray	\$ 42,500
A7355	Base Plate	\$ 4,000
A7395	Bite Block Covers (100)	\$ 40
A7396	Hygienic Protective Silicone for Temple Support (3 pairs)	\$ 75
A7376	Bite Block (3 piece)	\$ 100
A7390	Chin Support for Edentulous Jaws	\$ 75
A7391	Chin Support for Mandibular Joint Image	\$ 75
A7392	Chin Support for Sinus Image	\$ 75
A7375	Holder for Bite Block	\$ 75

Specifications	ProVecta S-Pan
Dimensions (inches WxDxH)	39 x 48 - 51 x 90
Weight Without Base Plate (Lbs)	231
Weight With Base Plate (Lbs)	342
Main Voltage	208 – 240 VAC
Frequency	50/60 Hz

IMAGING

ProVecta® HD Intraoral X-Ray



ProVecta HD

ProVecta HD Intraoral X-ray: Perfect in Function & Design.

Advanced ergonomics for advanced comfort.


The ProVecta HD Intraoral X-Ray unit offers a modern, slim design to fit the modern operatory. Featuring a precision of a DC generator with a 0.4mm focal spot, this intraoral x-ray is a perfect companion to any ScanX PSP scanner. And its operator friendly, including a right-angle technology tube head and an intuitive icon driven control panel.

- Modern, slim design for easy, precise positioning
- Perfectly matched to PSP image plates and sensors
- Highly reliable and long service life X-ray tubes
- mA and kV levels are adjustable
- DC generator with 0.4 mm focal spot

Choice of 3 arm lengths.

- Short: 68"
- Medium: 74"
- Long: 86"

The ProVecta HD works with wired sensors, cordless flexible phosphor sensors (PSPs) and film.



Horizontal Arm Length	A	B
17.7" (450 mm) PN A6350-15	68.4" (1738 mm)	74.1" (1883 mm)
23.6" (600 mm) PN A6350-24	74.3" (1888 mm)	80.0" (2033 mm)
35.4" (900 mm) PN A6350-33	86.1" (2188 mm)	91.9" (2333 mm)

PART NO.	DESCRIPTION	PRICE
A6350-15S	ProVecta HD Intraoral X-Ray Wall Mount Unit, 15" arm Single Stud Wall Mount	\$ 5,200
A6350-24S	ProVecta HD Intraoral X-Ray Wall Mount Unit, 24" arm Single Stud Wall Mount	\$ 5,200
A6350-33S	ProVecta HD Intraoral X-Ray Wall Mount Unit, 33" arm Single Stud Wall Mount	\$ 5,200
A6350-15D	ProVecta HD Intraoral X-Ray Wall Mount Unit, 15" arm Double Stud Wall Mount	\$ 5,200
A6350-24D	ProVecta HD Intraoral X-Ray Wall Mount Unit, 24" arm Double Stud Wall Mount	\$ 5,200
A6350-33D	ProVecta HD Intraoral X-Ray Wall Mount Unit, 33" arm Double Stud Wall Mount	\$ 5,200

Specifications	ProVecta HD (short / medium / long)
Arm Length (inches)	68 / 74 / 86
Storage Height (inches)	45
Weight (Lbs)	55 / 56 / 59
Focal Spot Size (mm)	0.4 (IEC 60336)
Power Consumption (W)	500
Communication Connection	RJ45, RJ11, and 3 Wire

IMAGING

ScanX® Digital Radiography



BW



CEPH



PAN



ScanX Duo
PN D1000F



ScanX Intraoral
PN F3600



ScanX Classic
PN F3700



PART NO.	DESCRIPTION	PRICE
D1000F	ScanX Duo for Chairside Imaging (Simultaneous Scanning of 2 Phosphor Sensors (PSP), sizes 0–4) Includes: USB cable; power supply and cord; 2 each size 2 PSP guides; 1 each size 1 PSP guide; 1 each size 0 PSP guide; 8 each size 2 PSPs; 1 box barrier envelopes (qty. 300) for size 2 PSPs; 1 transfer box for transporting exposed intraoral PSPs.	\$10,495
F3600	ScanX Intraoral for Centralized Imaging (Simultaneous Scanning of 4 Phosphor Sensors (PSP), sizes 0–4) Includes: USB cable; power supply and cord; 4 each size 2 PSP guides; 1 each size 0, 1, and 3 PSP guides; 20 each size 2 PSPs; 1 box barrier envelopes (qty. 300) for size 2 PSPs; 1 transfer box for transporting exposed intraoral PSPs.	\$14,495
F3700	ScanX Classic for Intraoral & Extraoral Imaging (Scans Intraoral Phosphor Sensors (PSP), and Extraoral PAN, CEPH and TMJ Phosphor Sensors) Includes: USB cable; power supply and cord; 4 each size 2 PSP guides; 1 each size 0, 1, and 3 PSP guides; 20 each size 2 PSPs; 1 box barrier envelopes (qty. 300) for size 2 PSPs; 1 transfer box for transporting exposed intraoral PSPs.	\$22,995
F3750	ScanX Ortho Includes: USB cable; power supply and cord; 1 each size 0, 1, 2 and 3 PSP guides.	\$22,495



Superior quality with any patient: pedo to adult.



Bite wings that are correctly positioned. No cord, no gap.



Sharp contrast to the apex.



ScanX Swift
PN G8000



"Flexible" PSP's

PART NO.	DESCRIPTION	PRICE
G8000	ScanX Swift for Chairside Imaging (Single-Track Scanning of 1 Phosphor Sensor (PSP), size 0-2) Includes: Installation disc; USB cable and power supply; 4 each of size 0 and 1, and 8 of size 2; Barrier envelopes - 100 each of size 0 and 1, and 300 of size 2; Barrier film roll; 2 Sensor wipes; 1 Scanner Cleaning Sheet.	\$8,400

Specifications	ScanX Swift	ScanX Duo	ScanX Intraoral	ScanX Classic	ScanX Ortho
Dimensions (inches WxHxD)	8.9 x 9.2 x 9.6	10 x 12 x 10.75	15.5 x 16 x 15.5	15.5 x 24.5 x 15.5	15.5 x 24.5 x 15.5
Weight (Lbs)	14.5	21.5	44	51	51
Electrical	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz


IMAGING

ScanX® Accessories



PART NO.	DESCRIPTION	PRICE
73445-0	Intraoral Phosphor Sensor (PSP); Size 0 (2 pack)	\$ 82
73445-1	Intraoral Phosphor Sensor (PSP); Size 1 (2 pack)	\$ 82
73445-2	Intraoral Phosphor Sensor (PSP); Size 2 (4 pack)	\$ 162
73445-3	Intraoral Phosphor Sensor (PSP); Size 3 (2 pack)	\$ 98
73445-4	Intraoral Phosphor Sensor (PSP); Size 4 (Single Pack)	\$ 126
73578-5	Extraoral Phosphor Sensor (PSP); Panoramic Size 5" x 12"	\$ 1,134
73578-6	Extraoral Phosphor Sensor (PSP); Panoramic Size 6" x 12" (15cm x 30cm)	\$ 1,168
73578-8	Extraoral Phosphor Sensor (PSP); Cephalometric Size 8" x 10"	\$ 1,222
73578-57	Extraoral Phosphor Sensor (PSP); TMJ Size 5" x 7"	\$ 682
73248-0	Patented, Easy-Open Barrier Envelopes; Size 0 (100 pack)	\$ 22
73248-1	Patented, Easy-Open Barrier Envelopes; Size 1 (100 pack)	\$ 22
73248-2	Patented, Easy-Open Barrier Envelopes; Size 2 (300 pack)	\$ 58
73248-2K	Patented, Easy-Open Barrier Envelopes; Size 2 (1,000 pack)	\$ 166
73248-3	Patented, Easy-Open Barrier Envelopes; Size 3 (100 pack)	\$ 30
73248-4	Patented, Easy-Open Barrier Envelopes; Size 4 (50 pack)	\$ 30

NOTE: Air Techniques will accept returned product that is unopened and in resalable condition within 90 days of the invoice date to receive full credit without a stocking fee.



Cleaning Sheets for
ScanX Phosphor Sensors



Intraoral Phosphor Sensor (PSP)
Guides



Transfer Box
PN 73470

PART NO.	DESCRIPTION	PRICE
73566-0	Intraoral Phosphor Sensor (PSP) Guide; Size 0	\$ 37
73566-1	Intraoral Phosphor Sensor (PSP) Guide; Size 1	\$ 37
73566-2	Intraoral Phosphor Sensor (PSP) Guide; Size 2	\$ 37
73566-3	Intraoral Phosphor Sensor (PSP) Guide; Size 3	\$ 37
73566-4	Intraoral Phosphor Sensor (PSP) Guide; Size 4	\$ 37
73470	Intraoral Phosphor Sensor (PSP) Transfer Box	\$ 37
B2010	Cleaning Sheets for ScanX (12 pack)	\$ 159
B2020	Cleaning Sheets for ScanX (25 pack)	\$ 271
B3010	Cleaning Sheets for ScanX Intraoral (12 pack)	\$ 128
B3020	Cleaning Sheets for ScanX Intraoral (25 pack)	\$ 220
D1010	Cleaning Sheets for ScanX Duo (12 pack)	\$ 71
D1020	Cleaning Sheets for ScanX Duo (25 pack)	\$ 136
B8910	Cleaning Wipes for Phosphor Sensors (PSP) (50 pack)	\$ 71
G8195	Barrier Film Roll (900/roll)	\$ 40

NOTE: Air Techniques will accept returned product that is unopened and in resalable condition within 90 days of the invoice date to receive full credit without a stocking fee.

IMAGING

Cam-X™ Luna HD Intraoral Camera



Cam-X Luna HD
USB Handpiece
PN G1100



MACRO



SINGLE TOOTH



ARCH



FULL MOUTH

Cam-X Luna HD Intraoral Camera:

High Definition for the Highest Diagnostic Clarity.

Building on years of experience producing Intraoral cameras, Air Techniques offers the Luna, a high definition video intraoral camera that can focus continuously from tooth surface to full patient portraits. Now dentist can use one device to take full smile reference images and macro inspection of a patient's tooth or recent procedure results.

HD Quality Imaging

- Live 720p video at 30 frames per second.
- Continuous focusing from 2mm-300mm.
- Precision 7 element lens.

Ergonomic Design

- Left and right hand use with 120 degrees accessible controls.
- Precision crafted focus ring enables one device to capture all your dental imaging needs with; patient portrait, full mouth, bridge, single tooth and tooth surface.

Versatile

- USB2 connection, allows simple connection to PCs and portability between operatories.
- Twain interface allows integration with leading practice software.
- A 16:9 aspect ratio fits modern displays.



Cam-X Luna HD Intraoral Camera PN G1000
(Laptop not included)

Cam-X Luna HD USB Handpiece in Wall Mount



Cam-X Luna HD Sheaths

PART NO.	DESCRIPTION	PRICE
G1000	Cam-X Luna HD Camera System Includes Luna HD Hand piece, 50 Sheaths, 10' Umbilical Cable, 1-Wall Mount	\$ 5,995
G1100	Cam-X Luna HD Hand Piece Only	\$ 5,695
G1072	Sheaths for Cam-X Luna HD (100 Pack)	\$ 97
G1073	Sheaths for Cam-X Luna HD (500 Pack)	\$ 388
G1151	Cam-X Luna HD - USB Handpiece Cable, 10 ft.	\$ 104
G4011	Wall Mount for Cam-X Luna HD	\$ 108
D7365	Visix Pro Imaging Software for Cam-X Luna HD	\$ 795


IMAGING

Cam-X™ Polaris Intraoral Camera



Cam-X Polaris
USB Handpiece
PN G5010

PART NO.	DESCRIPTION	PRICE
G5000	Cam-X Polaris Intraoral Camera Includes: Polaris handpiece, 50 sheaths, 10' USB umbilical cable, 1-wall mount. Compatible with Windows 7 and Windows 8, 32/64 Bit.	\$ 3,825
G5010	Cam-X Polaris Handpiece Only	\$ 3,520
G4600	Cam-X Polaris/Spectra Intraoral Imaging Combo Includes: 1 Polaris and 1 Spectra. Plus all accessories	\$ 6,630



Cam-X Polaris Intraoral Camera PN G5000
(Laptop not included)

Cam-X Polaris USB
Handpiece in Wall Mount






Cam-X Polaris Sheaths

PART NO.	DESCRIPTION	PRICE
G5111	Sheaths for Cam-X Polaris (100 pack)	\$ 97
G5110	Sheaths for Cam-X Polaris (500 pack)	\$ 388
G4011	Wall Mount for Cam-X Polaris	\$ 108
G1151	Cam-X Polaris - USB 2 Handpiece Cable, 10 ft.	\$ 104
D7365	Visix Pro Imaging Software for Cam-X Polaris	\$ 812


IMAGING

Cam-X™ Spectra Caries Detection Aid



Cam-X Spectra
USB Handpiece
PN G4010



PART NO.	DESCRIPTION	PRICE
G4000	Cam-X Spectra Caries Detection Aid Includes: Spectra handpiece, 50 sheaths, 10 positioning spacers, 10' USB umbilical, and 1 wall mount handpiece holder. Compatible with Windows 7 and Windows 8, 32/64 Bit.	\$ 4,075
G4010	Cam-X Spectra Handpiece Only	\$ 3,770
G4600	Cam-X Spectra/Polaris Intraoral Imaging Combo Includes: 1 Spectra and 1 Polaris. Plus all accessories	\$ 6,630



PART NO.	DESCRIPTION	PRICE
G4640	Cam-X Spectra – 10mm Distance Spacers (25 pack) Blocks both LED and stray ambient light. Positions Cam-X Spectra at the correct height above the tooth surface.	\$ 164
G5111	Sheaths for Cam-X Spectra (100 pack)	\$ 97
G5110	Sheaths for Cam-X Spectra (500 pack)	\$ 388
G4011	Wall Mount for Cam-X Spectra	\$ 108
G4650	Cam-X Spectra - USB 2 Handpiece Cable, 10 ft.	\$ 104
D4065	Visix Pro Imaging Software for Cam-X Spectra	\$ 812

IMAGING

Cam-X™ Polaris & Cam-X™ Spectra Clinical Images



Cam-X Polaris: High Image Quality
with Low-Glare LED Light


Cam-X Spectra: Existence &
Extent of Caries Obvious



VISIX®

Visix Pro Dental Imaging Software
PN 74500 (Monitor not included)

PART NO.	DESCRIPTION	PRICE
74580	Visix Pro Imaging Complete Includes: 5 licenses, training, 1 year support and 3-hr Onsite Training. Visix Tel# 1.866.908.4749	\$ 4,490
74500	Visix Pro Dental Imaging Software Includes: 5 licenses.	\$ 2,242
74505	Visix Pro Dental Imaging Software License Pack Add-On 5 license/pack.	\$ 1,119
74509	Visix Pro License Add-On	\$ 282
74510	Visix Pro Dental Imaging In-Office Training On-site, 3 hours.	\$ 1,070
74525	Visix Pro Software Support 1 year support contract.	\$ 1,140
74540	Visix Film Scanning/Intraoral Camera Module (Standard Visix)	\$ 1,680
D4065	Visix for Cam-X™ Spectra	\$ 812
D7365	Visix for Cam-X™ Polaris	\$ 812



Spectra Interface



Polaris Interface

TWAIN communication is the “**No Pain Solution**” for software integration with Spectra & Polaris.

- TWAIN interface allows for display of live video
- Button functionality
- Capture images directly into most TWAIN 2.1 compliant dental imaging software programs
- Seamlessly integrate Spectra and Polaris - two revenue generating products - directly into your existing software



TWAiN
Linking Images With Applications

IMAGING

A/T 2000®XR Film Processors



A/T 2000®XR
PN 45005

PART NO.	DESCRIPTION	PRICE
45005	A/T 2000®XR Includes: 1 leveling base, 1 bottle spray 2000 Plus, 4 Air Techniques cleaning films, 1 instructional DVD, 1 case chemistry (# 40270), 1 bottle of starter (#40280). Electrical: 115 Volts, 60 Hz, 8 Amps Dimensions: 18" x 15" x 25" (H x W x D) Weight: 75 Lbs. Water Pressure: 80 psi max. / 30 psi min. Water Flow: 1/2 gallon per minute (while processing film)	\$11,016
45009	A/T 2000®XR with Water Recirculation System Includes: Same as above plus 2 1/2 gallon water container. Cannot be connected to fresh water supply.	\$11,980
45110	Daylight Loader for A/T 2000®XR Allows processor to be installed outside the darkroom. Fits A/T 2000®, A/T 2000®PLUS and A/T 2000®XR Dimensions: 21 1/2" x 16 1/4" x 13 3/4" (H x W x D) Color: Grey	\$ 1,686
45115	Panoramic Flat Bottom Daylight Loader for A/T 2000® Flat bottom enables you to save space. Dimensions: 10 3/4" x 16 1/4" x 13 3/4" (H x W x D)	\$ 1,708



Peri-Pro III
(Shown with Daylight Loader)
PN 94000DL



Film Duplicator
Daylight Loader
PN 92000



Daylight Loader
PN 90090

PART NO.	DESCRIPTION	PRICE
94000	Peri-Pro III For automatic chemistry temperature control. Includes: 1 case chemistry (# 90800) Film Carriers: 3 Pedo film carriers, 2 adult anterior film carriers, 2 bite-wing film carriers, 1 occlusal film carrier, 1 water bottle Electrical: 115 VAC, 5 Amps. Dimensions: 8 1/2" x 25" x 10" (H x W x D). Weight: 22 Lbs. Processing Capacity: Accepts 8 #2 films, or carriers loaded with 0, 1, 3, or 4 films Processing Time: 5 minutes	\$ 3,758
94000DL	Peri-Pro III w/ Daylight Loader For film processing inside or outside a darkroom. Dimensions: 15 3/4" x 25" x 10" (H x W x D). Weight: 23 Lbs.	\$ 4,168
94000DUP	Peri-Pro III w/ Film Duplicator and Daylight Loader For film processing and duplication inside or outside a darkroom. Electrical: 115 VAC, 5 Amps. Dimensions: 15 3/4" x 25" x 13" (H x W x D). Weight: 24 Lbs.	\$ 4,682
92000	Film Duplicator and Daylight Loader Attaches to all Peri-Pro models.	\$ 956
90090	Daylight Loader For film processing inside or outside of a darkroom. Attaches to the cover of all Peri-Pro models.	\$ 424



“

I've relied on my AirStar air compressor every day since I purchased it in 1988. In all that time, it's never let me down a single time."

– Justin Fried, D.D.S.

UTILITY ROOM



Air Compressors

Dry Vacuums

Wet Vacuums

Amalgam Separator

Utility Accessories

Clinical Utility
Packages

UTILITY ROOM

AirStar® NEO Air Compressors



AirStar 10 NEO
PN AS10 NEO



AirStar 50 NEO
PN AS50 NEO



AirStar 30 NEO
PN AS30 NEO



AirStar 70 NEO
PN AS70 NEO

PART NO.	DESCRIPTION	PRICE
AS10 NEO	AirStar 10 NEO; Max. 2 Simultaneous Users; 115 VAC	\$ 6,375
AS21 NEO	AirStar 21 NEO; Max. 3 Simultaneous Users; 115 VAC	\$ 7,095
AS22 NEO	AirStar 22 NEO; Max. 3 Simultaneous Users; 208/230 VAC	\$ 7,095
AS30 NEO	AirStar 30 NEO; Max. 4 Simultaneous Users; 208/230 VAC	\$ 9,675
AS40 NEO	AirStar 40 NEO; Max. 5 Simultaneous Users; 208/230 VAC	\$10,565
AS50 NEO	AirStar 50 NEO; Max. 7 Simultaneous Users; 208/230 VAC	\$11,450
AS70 NEO	AirStar 70 NEO; Max. 10 Simultaneous Users; 208/230 VAC	\$15,550

♦ Recommend using Preventative Maintenance Kits found on page 28.

Specifications	AS10 NEO	AS21 NEO	AS22 NEO	AS30 NEO
Dimensions (inches WxHxD)	26.5 x 29.5 x 20			
Weight (Lbs)	175	185	185	230
Compressor Head	Single .75 HP	Single 1.5 HP	Single 1.5 HP	Twin .75 HP
Voltage Rating	115	115	208/230	208/230
Voltage (Min/Max)♦	105/125	105/125	200/250	200/250
Breaker Size	20A	30A	20A	20A
Wire Size (Gauge)	12	10	12	12
CFM (Cubic Ft./Min.) @ 80psi Rated Capacity	2.5	5.0	5.0	5.0
Tank Size (Gallons)	9.6	9.6	9.6	9.6
Max. Users	2	3	3	4

Warranty: All AirStar products feature a 5 year limited extended warranty.

Specifications	AS40 NEO	AS50 NEO	AS70 NEO
Dimensions (inches WxHxD)	30.75 x 32.5 x 21	30.75 x 32.5 x 21	48 x 32.5 x 21
Weight (Lbs)	255	285	400
Compressor Head	Single 1.5 HP & Single .75 HP	Twin 1.5 HP	Triple 1.5 HP
Voltage Rating	208/230	208/230	208/230
Voltage (Min/Max)♦	200/250	200/250	200/250
Breaker Size	20A	30A	40A
Wire Size (Gauge)	12	10	8
CFM (Cubic Ft./Min.) @ 80psi Rated Capacity	7.5	10.0	15.0
Tank Size (Gallons)	16	16	27
Max. Users	5	7	10

♦ See Page 38 on Buck/Boost Transformer if voltage is near or outside of min/max ratings.

UTILITY ROOM

AirStar® Air Compressor



AirStar 40
PN AS40



AirStar 70
PN AS70



PART NO.	DESCRIPTION	PRICE
AS10	AirStar 10; Max. 2 Simultaneous Users; 115 VAC	\$ 6,086
AS10C	AirStar 10 with Factory-Installed Sound Cover; Max. 2 Simultaneous Users; 115 VAC	\$ 7,030
AS12	AirStar 12; Max. 2 Simultaneous Users; 230 VAC	\$ 6,086
AS12C	AirStar 12 with Factory-Installed Sound Cover; Max. 2 Simultaneous Users; 230 VAC	\$ 7,030
AS21	AirStar 21; Max. 3 Simultaneous Users; 115 VAC	\$ 6,758
AS21C	AirStar 21 with Factory-Installed Sound Cover; Max. 3 Simultaneous Users; 115 VAC	\$ 7,718
AS22	AirStar 22; Max. 3 Simultaneous Users; 208/230 VAC	\$ 6,774
AS22C	AirStar 22 with Factory-Installed Sound Cover; Max. 3 Simultaneous Users; 208/230 VAC	\$ 7,780
AS30	AirStar 30; Max. 4 Simultaneous Users; 208/230 VAC	\$ 9,060
AS30C	AirStar 30 with Factory-Installed Sound Cover; Max. 4 Simultaneous Users; 208/230 VAC	\$10,296
AS40	AirStar 40; Max. 5 Simultaneous Users; 208/230 VAC	\$ 9,890
AS40C	AirStar 40 with Factory-Installed Sound Cover; Max. 5 Simultaneous Users; 208/230 VAC	\$11,678
AS50	AirStar 50; Max. 7 Simultaneous Users; 208/230 VAC	\$10,738
AS50C	AirStar 50 with Factory-Installed Sound Cover; Max. 7 Simultaneous Users; 208/230 VAC	\$12,340
AS70	AirStar 70; Max. 10 Simultaneous Users; 208/230 VAC	\$14,560

♦ Recommend using Preventative Maintenance Kits found on page 28.

AirStar® Air Compressor Specifications

Specifications	AS10	AS12	AS21	AS22
Dimensions (inches WxHxD)	25 x 28.5 x 19.75	25 x 28.5 x 19.75	29 x 30.5 x 20	29 x 30.5 x 20
Dimensions (inches WxHxD; with Sound Reducing Cover)	25 x 30 x 22.5	25 x 30 x 22.5	31 x 32 x 22.25	31 x 32 x 22.25
Weight (Lbs)	170	170	200	200
Compressor Head	Single .75 HP	Single .75 HP	Single 1.5 HP	Single 1.5 HP
Voltage Rating	115	208/230	115	208/230
Voltage (Min/Max)♦	105/125	200/250	105/125	200/250
Breaker Size	20A	20A	30A	20A
Wire Size (Gauge)	12	12	10	12
CFM (Cubic Ft./Min.) @ 80psi Rated Capacity	2.5	2.5	5.0	5.0
Tank Size (Gallons)	6	6	12	12
Max. Users	2	2	3	3

Warranty: All AirStar products feature a 5 year limited extended warranty.

Specifications	AS30	AS40	AS50	AS70
Dimensions (inches WxHxD)	29 x 30.5 x 20	32.5 x 30.5 x 20	35.5 x 33.5 x 20.5	47.75 x 35 x 21.75
Dimensions (inches WxHxD; with Sound Reducing Cover)	31 x 32 x 22.25	33.25 x 32.25 x 22.5	36.5 x 33.5 x 22.75	N/A
Weight (Lbs)	240	255	290	430
Compressor Head	Twin .75 HP	Single 1.5 HP & Single .75 HP	Twin 1.5 HP	Triple 1.5 HP
Voltage Rating	208/230	208/230	208/230	208/230
Voltage (Min/Max)♦	200/250	200/250	200/250	200/250
Breaker Size	20A	20A	30A	40A
Wire Size (Gauge)	12	12	10	8
CFM (Cubic Ft./Min.) @ 80psi Rated Capacity	5.0	7.5	10.0	15.0
Tank Size (Gallons)	12	12	20	30
Max. Users	4	5	7	10

♦ See Page 38 on Buck/Boost Transformer if voltage is near or outside of min/max ratings.

UTILITY ROOM

AirStar® Accessories



Membrane Dryer



Filter w/ Sound Reducer

PART NO.	DESCRIPTION	PRICE
PREVENTATIVE MAINTENANCE KITS		
87351	Preventative Maintenance Kit for AS10 / AS10 NEO / AS12 Includes: Qty. 1 compressor head intake filter 10 microns, Qty. 1 membrane top filter 2 microns, Qty. 1 membrane bottom filter 35 microns, Qty. 1 top cover o-ring, Qty. 1 bottom filter bowl o-ring, Qty. 1 maintenance record card.	\$ 99
87352	Preventative Maintenance Kit for AS21 / AS21 NEO / AS22 / AS22 NEO / AS30 / AS30 NEO Includes: Qty. 2 compressor head intake filter 10 microns, Qty. 1 membrane top filter 2 microns, Qty. 1 membrane bottom filter 35 microns, Qty. 1 top cover o-ring, Qty. 1 bottom filter bowl o-ring, Qty. 1 maintenance record card.	\$ 160
87355	Preventative Maintenance Kit for AS40 / AS40 NEO Includes: Qty. 3 compressor head intake filter 10 microns, Qty. 1 membrane top filter 2 microns, Qty. 1 membrane bottom filter 35 microns, Qty. 1 top cover o-ring, Qty. 1 bottom filter bowl o-ring, Qty. 1 maintenance record card.	\$ 230
87353	Preventative Maintenance Kit for AS50 / AS50 NEO Includes: Qty. 4 compressor head intake filter 10 microns, Qty. 1 membrane top filter 2 microns, Qty. 1 membrane bottom filter 35 microns, Qty. 1 top cover o-ring, Qty. 1 bottom filter bowl o-ring, Qty. 1 maintenance record card.	\$ 299
87354	Preventative Maintenance Kit for AS70 / AS70 NEO Includes: Qty. 6 compressor head intake filter 10 microns, Qty. 2 membrane top filter 2 microns, Qty. 2 membrane bottom filter 35 microns, Qty. 2 top cover o-ring, Qty. 2 bottom filter bowl o-ring, Qty. 2 maintenance record card.	\$ 428
85694	AirStar Tandem Install Kit Connect 2 compressors to work as one.	\$ 352
85695	Multiple AirStar Timer Relay Installation Kit Use to connect up to 4 AirStars together to work as 1 unit.	\$ 2,180
87168	3/8" Filter, 5 Micron Replacement, Tank Outlet	\$ 240



Sound Cover

PART NO.	DESCRIPTION	PRICE
SOUND REDUCING COVERS - Reduces operating noise by 2-4 decibels.		
85961	AirStar 10	\$ 1,133
85182	AirStar 12	\$ 1,133
85962-1M	AirStar 21	\$ 1,342
85962-2M	AirStar 22	\$ 1,342
85963M	AirStar 30	\$ 1,437
87440M	AirStar 40	\$ 1,613
89523M	AirStar 50	\$ 1,788
89574M	AirStar 70 - comes with a 10ft cord. Requires: <ul style="list-style-type: none">• Dedicated 220v line in addition to the 220v line for the compressor• Receptacle• 15amp breaker	\$ 4,168

NOTE: Sound Reducing Covers will only apply to original AirStar models 10 - 70.

REMOTE AIR INTAKE KITS - Routes clean air from a remote location to the AirStar.

85491	AS10 / AS10 NEO, AS12	\$ 338
85492	AS21 / AS21 NEO, AS22 / AS22 NEO, AS30 / AS30 NEO	\$ 362
87361	AS40 / AS40 NEO	\$ 375
85493	AS50 / AS50 NEO	\$ 388
85494	AS70 / AS70 NEO	\$ 408

UTILITY ROOM

Mojave® Dry Vacuums



Mojave Dry Vac
PN V3



As awarded by
The Chicago Athenaeum:
Museum of Architecture and Design



Mojave Dry Vac
PN V5




Mojave Dry Vac
PN 2V3



PART NO.	DESCRIPTION	PRICE
V3	Mojave V3; Max. 3 Simultaneous Users; Includes 10-gal Tank, MMC	\$10,790
V5	Mojave V5; Max. 5 Simultaneous Users; Includes 10-gal Tank, MMC	\$13,120
V7	Mojave V7; Max 7 Simultaneous Users; Includes 10-gal Tank, MMC	\$14,900
2V3	Mojave 2V3; Dual System; Max. 6 Simultaneous Users; Includes 10-gal Tank, MMC	\$18,244
2V5	Mojave 2V5; Dual System; Max. 10 Simultaneous Users; Includes 10-gal Tank, MMC	\$23,364
2V7	Mojave 2V7; Dual System; Max. 14 Simultaneous Users; Includes Dual Continuum Tank, MMC	\$27,000
3V5	Mojave 3V5; Triple System; Max. 15 Simultaneous Users; Includes Dual Continuum Tank, MMC	\$33,386
4V5	Mojave 4V5; Quad System; Max. 20 Simultaneous Users; Includes Dual Continuum Tank, MMC	\$42,836
2V3CT	Mojave 2V3 for Large Volume Applications; Dual System; Max. 6 Simultaneous Users; Includes Dual Continuum Tank, MMC	\$20,524
2V5CT	Mojave 2V5 for Large Volume Applications; Dual System; Max. 10 Simultaneous Users; Includes Dual Continuum Tank, MMC	\$25,594



Monarch CleanStream helps keep your vacuum running smoothly.
For more information turn to page 60.



NEW Easy-View
Code Window

Mojave Dry Vac
PN 2V7

Mojave Dry Vac
PN 3V5

Warranty: All Mojave products feature a 5 year limited extended warranty.

Specifications	V3	V5	V7	2V3	2V5
Pump Dimensions (inches WxHxD)	25 x 17 x 21	25 x 17 x 21	25 x 17 x 21	25 x 34** x 21	25 x 34** x 21
Tank Dimensions (inches WxHxD)	29 x 33 x 22	29 x 33 x 22			
Weight (Lbs)	Pump -145 Tank - 75	Pump -145 Tank - 75	Pump -145 Tank - 75	Pump -145(x2) Tank -75	Pump -145(x2) Tank -75
Circuit Breaker	20A	20A	20A	20A (x2)	20A (x2)
Voltage Range (Min/Max) ♦	205/240	205/240	205/240	205/240	205/240
Horsepower	1.5	1.5	2	1.5 (x2)	1.5 (x2)
Max. Users	3	5	7	6	10

Note: 1 User = 1HVE or 2SEs

Specifications	2V7	3V5	4V5	MMC
Pump Dimensions (inches WxHxD)	25 x 34** x 21	25 x 51** x 21	25 x 34** x 21 (x2)	11 x 10 x 3
Tank Dimensions (inches WxHxD)	29 x 42 x 23	29 x 42 x 23	29 x 42 x 23	N/A
Weight (Lbs)	Pump -145(x2) Tank - 150	Pump -145(x3) Tanks - 150	Pump -145(x4) Tanks - 150	13
Circuit Breaker	20A (x2)	20A (x3)	20A (x4)	15A
Voltage Range (Min/Max) ♦	205/240	205/240	205/240	108/132
Horsepower	2 (x2)	1.5 (x3)	1.5 (x4)	N/A
Max. Users	14	15	20	N/A

♦ See Page 38 on Buck/Boost Transformer if voltage is near or outside of min/max ratings.

**Pumps/Tanks Stacked

UTILITY ROOM

Mojave® Accessories



Mojave 10-gal. Tank
PN MT10



Mojave Dual
Continuum Tank
PN CT20-120V

PART NO.	DESCRIPTION	PRICE
TANKS		
MT10	10-gal Separation Tank for Mojave V3, V5, V7, 2V3, and 2V5, 120VAC	\$ 2,947
H5470	10-gal Separation Tank for Mojave V3, V5, V7, 2V3 AND 2V5, 220VAC	\$ 3,590
CT20-120V	Dual Continuum Tank for Mojave 2V7, 3V5, 4V5, 2V3CT and 2V5CT	\$ 4,208
CT20-220V	Dual Continuum Tank for Mojave 2V7, 3V5, 4V5, 2V3CT and 2V5CT	\$ 4,858
PUMPS		
V3PUMP	Mojave V3 Add-On Pump Only (For office expansion - requires Installation Kit)	\$ 9,390
V5PUMP	Mojave V5 Add-On Pump Only (For office expansion - requires Installation Kit)	\$11,820
V7PUMP	Mojave V7 Add-On Pump Only (For office expansion - requires Installation Kit)	\$14,250



Monarch CleanStream helps keep your vacuum running smoothly.
For more information turn to page 60.



Mojave Accessory Kit



Pump Accessory Kit

Mojave Dual System Installation Kit
PN MIK2

PART NO.	DESCRIPTION	PRICE
KITS		
H5231	Mojave Powered Automatic Drain Pump (115 VAC only)	\$ 1,057
H5217	Replacement Bowl, Screen & Gasket for Inlet Filter (use 1 per pump for all Mojave solutions)	\$ 56
MIK2	Mojave Dual System Installation Kit (included with dual systems; sold separately for office expansion)	\$ 408
MIK4	Mojave Triple and Quad System Installation Kit (included with triple and quad systems; sold separately for office expansion)	\$ 584

UTILITY ROOM

VacStar® Wet Vacuums



VacStar 20
PN VS20



VacStar 50
PN VS50

PART NO.	DESCRIPTION	PRICE
VS20	VacStar 20; Max. 2 Simultaneous Users; 115/230 VAC	\$ 3,188
VS40	VacStar 40; Max. 3 Simultaneous Users; 230 VAC	\$ 3,838
VS50	VacStar 50; Max. 4 Simultaneous Users; 230 VAC	\$ 6,112
VS50H	VacStar 50 with Factory-Installed Hydromiser; Max. 4 Simultaneous Users; 230 VAC	\$ 7,000
VS80	VacStar 80; Max. 7 Simultaneous Users; 230 VAC	\$ 6,888
VS80H	VacStar 80 with factory-installed hydromiser; Max. 7 Simultaneous Users; 230 VAC	\$ 7,598
53240	Water Lockout Relay; Control Panel on VacStar	\$ 113

HYDROMISER

H-2	HydroMiser for VacStar 20 (recycles up to 75% of water)	\$ 1,212
H-4	HydroMiser for VacStar 40 (recycles up to 75% of water)	\$ 1,212
56041	Hydromiser for VacStar 50 (recycles up to 75% of water)	\$ 996
56042	Hydromiser for VacStar 80 (recycles up to 75% of water)	\$ 996
55087	Hydromiser Wall-Mount Kit	\$ 170
55540	Air/Water Separator	\$ 720
55880C	Solids Collector Replacement Kit for VacStar 20/40 (12 pack)	\$ 434
55094C	Solids Collector Replacement Kit for VacStar 50/80 (8 pack)	\$ 628
55007	VacStar Single to Dual Pump Conversion Install Kit (Pump sold separately)	\$ 857



Monarch CleanStream helps keep your vacuum running smoothly.
For more information turn to page 60.



VacStar 40
PN VS40



VacStar 80H
PN VS80H


Warranty: All VacStar products feature a 5 year limited extended warranty.

Specifications	VS20	VS40	VS50	VS50H	VS80	VS80H
Dimensions (inches WxHxD)	11 x 14 x 11	11 x 17 x 11	28 x 22 x 16	28 x 25 x 16	28 x 22 x 16	28 x 25 x 16
Weight (Lbs)	68	85	160	170	200	210
# of Motors	1	1	2	2	2	2
Voltage Rating	115/230	230	230	230	230	230
Voltage (Min/Max)♦	105/125 205/253	205/253	205/253	205/253	205/253	205/253
Circuit Breaker	30A (115V) 20A (230V)	20A	20A (x2) or 30 amp	20A (x2) or 30 amp	20A (x2) or 40 amp	20A (x2) or 40 amp
Water Flow Rate / pump (gpm)	.50	.75	.50 per pump	.12 per pump	.75 per pump	.18 per pump
Max. Users	2	3	4	4	7	7

♦ See Page 38 on Buck/Boost Transformer if voltage is near or outside of min/max ratings.

UTILITY ROOM

Acadia® Amalgam Separator



Acadia Alert
PN A1330



Acadia
Amalgam Separator
PN A1250

PART NO.	DESCRIPTION	PRICE
A1250	Acadia Amalgam Separator Supports up to ten operatories and can be installed with wet and dry vacuum systems. Installation kit included.	\$ 1,304
A1350	Acadia Amalgam Separator w/electronic Alert for Filter Cartridge Replacement	\$ 1,600
A1330	Electronic Alert for Filter Cartridge Replacement Automatically monitors your system and activates an audible alarm along with visual LED indicators when the filter reaches 90% capacity. Acadia Alert can be added to A1250, manufactured before January 1, 2012, and installed in any operatory to give you real-time notification.	\$ 388
A1390	Acadia Replacement Filter w/ Recycling Kit Includes: Prepaid shipping label and the necessary packaging for the filter to be sent to a certified recycling facility.	\$ 531
A1290	Acadia Replacement Filter Only (without recycling kit) Recommended to be changed annually, the four stage filter easily twists on and off making replacement simple.	\$ 440
A1370	Acadia Recycling Kit Only Includes: Prepaid shipping label and the necessary packaging for the filter to be sent to a certified recycling facility. Does not include replacement filter.	\$ 210
A5135	Acadia Power Supply Only to be used for Acadia Alert on Mojave models.	\$ 84



Check out Acadia's different configurations with Utility Stacking Rack, pgs 40-42.

Acadia Technical Specifications

Service	Up to 10 operatories
Liquid Flow Rate	1 liter per minute.
Certification	ISO 11143
Connection Type	3.81cm (1½") Hose barb adaptor with retaining clip
Dimensions (inches WxDxH)	20.6 x 14.7 x 26.9
Installation	Free standing or wall-mounted



Vacuum Equalizer
PN 56200

PART NO.	DESCRIPTION	PRICE
56200	Vacuum Equalizer Equilibrates run time among 2 or more pumps; alternates primary pump each day.	\$ 1,662
56230	Remote Switch Panel for Vacuum Equalizer	\$ 332

Vacuum Equalizer Specifications

Power Requirements	100 - 240 VAC, 1 Amp
Dimensions (inches LxWxH)	10 x 6.2 x 3.5
Capacity	Controls up to 4 pumps

UTILITY ROOM

Utility Accessories

Product Name	115 VAC Buck/Boost Transformer 2 kVA Rating (not supplied by Air Techniques)	208/230 VAC Buck/Boost Transformer 3.9 kVA Rating PN 67002	208/230 VAC Buck/Boost Transformer 7.8 kVA Rating PN 67000-1	Qty. Needed	Retail Price (each)
AirStar 10 / NEO	X			1	N/A
AirStar 12		X		1	N/A
AirStar 21 / NEO	X			1	N/A
AirStar 22 / NEO		X		1	\$ 744.00
AirStar 30 / NEO		X		1	\$ 744.00
AirStar 40 / NEO		X		1	\$ 744.00
AirStar 50 / NEO		X		1	\$ 744.00
AirStar 70 / NEO			X	1	\$ 1,620.00
VacStar 20	X			1	N/A
VacStar 20		X		1	\$ 744.00
VacStar 40		X		1	\$ 744.00
VacStar 50/50H (1 panel breaker)		X		1	\$ 744.00
VacStar50/50H (2 panel breakers)		X		2	\$ 744.00
VacStar80/80H (1 panel breaker)			X	1	\$ 1,620.00
VacStar80/80H (2 panel breakers)		X		2	\$ 744.00
Mojave V3		X		1	\$ 744.00
Mojave V5		X		1	\$ 744.00
Mojave V7		X		1	\$ 744.00
Mojave 2V3/2V3CT		X		2	\$ 744.00
Mojave 2V5/2V5CT		X		2	\$ 744.00
Mojave 2V7		X		2	\$ 744.00
Mojave 3V5		X		3	\$ 744.00
Mojave 4V5		X		4	\$ 744.00

Buck/Boost Transformers are needed if voltage is below or above the following: <i>Air Techniques' stocked Buck/Boost Transformers will decrease or increase voltage by 10%.</i>	AirStar 10/NEO and AirStar 21/NEO	105-125 VAC
	AirStar 12 / NEO, 22 - 70 / NEO	200-250 VAC
	VS20 (Wired to 115 VAC)	105-125 VAC
	VS20 (Wired to 208-230 VAC)	205-240 VAC
	VS40, VS50, VS50H, VS80, VS80H	205-240 VAC
	V3, V5, V7, 2V3/2V3CT, 2V5/2V5CT, 2V7, 3V5, 4V5	205-240 VAC



PN 53251



PN 53250



PN 53133

Remote Water
Control Valve
PN 53020T

PART NO.	DESCRIPTION	PRICE
53111	Remote Control Panel 24V 1 Switch Needs 2 gang electrical box*	\$ 240
53251	Remote Control Panel 24V 2 Switches Needs 2 gang electrical box*	\$ 298
53250	Remote Control Panel 24V 3 Switches Needs 2 gang electrical box*	\$ 360
53133	Remote Control Panel 24V 4 Switches Needs 2 gang electrical box*	\$ 394
53240	Two Pump Water Lockout Relay Kit for New style 24V Panel	\$ 113
53245	Low Voltage Transformer Assembly; Water Systems, 80VA for 53170 & 53173	\$ 475
53020	Remote Water Control Valve w/ Filter & Transformer Pipe Size: 3/4" Solenoid Voltage: 24V	\$ 1,742
53020-1	Remote Water Control Valve w/ Filter Pipe Size: 3/4" Solenoid Voltage: 115V	\$ 1,524
53021	Remote Water Control Valve w/ Transformer, w/o Filter Pipe Size: 3/4" Solenoid Voltage: 24V	\$ 1,276
53021-1	Remote Water Control Valve w/o Filter Pipe Size: 3/4" Solenoid Voltage: 115V	\$ 964
53170	Remote Water Control Valve w/ Transformer, w/ Filter Pipe Size: 1" Solenoid Voltage: 24V	\$ 1,904
53171	Remote Water Control Valve w/ Filter Pipe Size: 1" Solenoid Voltage: 115V	\$ 1,440
53173	Remote Water Control Valve w/ Transformer, w/o Filter Pipe Size: 1" Solenoid Voltage: 24V	\$ 1,576
53174	Remote Water Control Valve w/o Filter Pipe Size: 1" Solenoid Voltage: 115V	\$ 1,256
53150	Low Voltage Transformer Assembly for 53020 & 53021	\$ 424

*Air, Vac, Water designation rings will be used to customize each panel.
For multiple switch models, additional rings are included.

UTILITY ROOM

Utility Accessories



AS30 with VS50H & Acadia using Utility Stack Rack Mini



AS70 with VS80H & Acadia using Utility Stack Rack XL

Stack Rack Features:

- Saves square footage in your mechanical room.
- Freestanding stack rack.
- Three height adjustments for upper shelf.
- Compressor mounted over vacuum pump or vice versa.

PART NO.	DESCRIPTION	PRICE
UTILITY STACKING RACKS		
56600	Utility Stack Rack Mini Used to stack AirStar 30, VS50/50H, and smaller in any configuration. Stand Dimensions: 38.25 x 26.50 x 42.50 (highest option) (inches LxWxH) 38.25 x 26.50 x 37.50 (middle option) 38.25 x 26.50 x 30.50 (lowest option)	\$ 1,168
56650	Utility Stack Rack XL Used to stack AirStar 50/70, V3, V5, 2V3, 2V5, VS80/80H in any configuration. Stand Dimensions: 54.25 x 26.50 x 42.50 (highest option) (inches LxWxH) 54.25 x 26.50 x 37.50 (middle option) 54.25 x 26.50 x 30.50 (lowest option)	\$ 1,362



AS22C with VS40 with Air/Water Separator & Acadia using Utility Stack Rack Mini



Mojave 2V5 Pump with AS70 using Utility Rack XL


PART NO.	DESCRIPTION	PRICE
OPTIONAL ACCESSORIES		
56635	Drip Tray Accessory Kit for Utility Stack Rack Mini Used when mounting vacuum pump over compressor.	\$ 367
56655	Drip Tray Accessory Kit for Utility Stack Rack XL Used when mounting vacuum pump over compressor.	\$ 488
54145	Wall Mount Anchor Bracket Kit Used when building codes require stack to be secured to a wall.	\$ 139
56645	VacStar 20/40 Top Shelf Mount Kit Used when single VS20/VS40 will be mounted on the top shelf over the compressor.	\$ 68


*See pg. 42 for Utility Stack Rack Diagram

UTILITY ROOM

Utility Room Configuration

2V5 with AS70 using Utility Stack Rack XL





Clinical - Mojave Dry Vacuum
PN 2V15

PART NO.	DESCRIPTION	PRICE
V15	One V15 Dry Vacuum Pump and One CT12 Tank	\$21,606
2V15	Two V15 Dry Vacuum Pumps and One CT12 Tank	\$33,986
3V15	Three V15 Dry Vacuum Pumps, One CT12 and One CT34 Tank	\$55,760
4V15	Four V15 Dry Vacuum Pumps, One CT12 and One CT34 Tank	\$68,284

Specifications	V15	2V15	3V15	4V15
Pump Dimensions (inches WxHxD)	32.5 x 23.5 x 26	32.5 x 47** x 26	32.5 x 47** x 26 + 32.5 x 23.5 x 26	32.5 x 47** x 26 (x2)
Tank Dimensions (inches WxHxD)	25 x 42 x 23	25 x 42 x 23	50 x 42** x 23	50 x 42*** x 23
Weight (Lbs)	Pump – 250 Tank – 175	Pump – 250(x2) Tank – 175	Pump – 250(x3) Tank – 175(x2)	Pump – 250(x4) Tank – 175(x2)
Circuit Breaker	40A	40A(x2)	40A(x3)	40A(x4)
Voltage Range (Min/Max)	205/240	205/240	205/240	205/240
Horsepower	6.2	12.4	18.6	24.8
Max. Users	15	30	45	60

Note: 1 User = 1HVE or 2SEs
**Pumps Stacked
***Tanks Side-by-Side



“

*Monarch CleanStream Evacuation System
Cleaner keeps my VacStar 50H running
like a true workhorse.”*

– Robert Branca, D.D.S.



MERCHANDISE

Imaging Accessories

Chemistry

Processor Accessories

Surface Disinfectant

Instrument Cleaner

Hand Sanitizer & Hand Lotion

Evacuation System Cleaner

Water Line Cleaner

MERCHANDISE

Imaging Accessories



Base Plate



Hygienic Protective Silicone
Temple Supports



Bite Block w/ Holder



Mandibular Joint




Edentulous Jaws



Sinus Image

PART NO.	DESCRIPTION	PRICE
A7355	Base Plate	\$ 4,000
A7395	Bite Block Covers (100)	\$ 40
A7396	Hygienic Protective Silicone for Temple Support (3 pairs)	\$ 75
A7376	Bite Block (3 piece)	\$ 100
A7390	Chin Support for Edentulous Jaws	\$ 75
A7391	Chin Support for Mandibular Joint Image	\$ 75
A7392	Chin Support for Sinus Image	\$ 75
A7375	Holder for Bite Block	\$ 75



PART NO.	DESCRIPTION	PRICE
73445-0	Intraoral Phosphor Sensor (PSP); Size 0 (2 pack)	\$ 82
73445-1	Intraoral Phosphor Sensor (PSP); Size 1 (2 pack)	\$ 82
73445-2	Intraoral Phosphor Sensor (PSP); Size 2 (4 pack)	\$ 162
73445-3	Intraoral Phosphor Sensor (PSP); Size 3 (2 pack)	\$ 98
73445-4	Intraoral Phosphor Sensor (PSP); Size 4 (Single Pack)	\$ 126
73578-5	Extraoral Phosphor Sensor (PSP); Panoramic Size 5" x 12"	\$ 1,134
73578-6	Extraoral Phosphor Sensor (PSP); Panoramic Size 6" x 12" (15cm x 30cm)	\$ 1,168
73578-8	Extraoral Phosphor Sensor (PSP); Cephalometric Size 8" x 10"	\$ 1,222
73578-57	Extraoral Phosphor Sensor (PSP); TMJ Size 5" x 7"	\$ 682

NOTE: Air Techniques will accept returned product that is unopened and in resalable condition within 90 days of the invoice date to receive full credit without a stocking fee.

MERCHANDISE

Imaging Accessories



Barrier Envelopes



Phosphor Sensor
(PSP) Guides



Transfer Box
PN 73470

PART NO.	DESCRIPTION	PRICE
73470	Intraoral Phosphor Sensor (PSP) Transfer Box	\$ 37
73248-0	Patented, Easy-Open Barrier Envelopes; Size 0 (100 pack)	\$ 22
73248-1	Patented, Easy-Open Barrier Envelopes; Size 1 (100 pack)	\$ 22
73248-2	Patented, Easy-Open Barrier Envelopes; Size 2 (300 pack)	\$ 58
73248-2K	Patented, Easy-Open Barrier Envelopes; Size 2 (1,000 pack)	\$ 166
73248-3	Patented, Easy-Open Barrier Envelopes; Size 3 (100 pack)	\$ 30
73248-4	Patented, Easy-Open Barrier Envelopes; Size 4 (50 pack)	\$ 30
73566-0	Intraoral Phosphor Sensor (PSP) Guide; Size 0	\$ 37
73566-1	Intraoral Phosphor Sensor (PSP) Guide; Size 1	\$ 37
73566-2	Intraoral Phosphor Sensor (PSP) Guide; Size 2	\$ 37
73566-3	Intraoral Phosphor Sensor (PSP) Guide; Size 3	\$ 37
73566-4	Intraoral Phosphor Sensor (PSP) Guide; Size 4	\$ 37
G8195	Barrier Roll (900/roll)	\$ 40

NOTE: Air Techniques will accept returned product that is unopened and in resalable condition within 90 days of the invoice date to receive full credit without a stocking fee.

Imaging Accessories



PSP Cleaning Wipes
PN B8910



Cleaning Sheets for ScanX
Phosphor Sensors



Cam-X Spectra 10mm
Distance Spacers
PN G4640



ScanX PSP Cleaning Wipes



Sheaths for Cam-X™ Spectra/Polaris

PART NO.	DESCRIPTION	PRICE
B8910	Cleaning Wipes for Phosphor Sensors (PSP) (50 pack)	\$ 71
G4011	Wall Mount for Cam-X Spectra / Polaris	\$ 108
G1072	Sheaths for Cam-X Luna HD (100 pack)	\$ 97
G1073	Sheaths for Cam-X Luna HD (500 pack)	\$ 388
G5111	Sheaths for Cam-X Spectra / Polaris (100 pack)	\$ 97
G5110	Sheaths for Cam-X Spectra / Polaris (500 pack)	\$ 388
G1151	Cam-X Luna HD / Polaris - 10 ft. USB Cable	\$ 104
G4650	Cam-X Spectra - 10 ft. USB Cable	\$ 104
G4640	Cam-X Spectra – 10mm Distance Spacers (25 pack) Blocks both LED and stray ambient light. Positions Spectra at the correct height above the tooth surface.	\$ 164
B2010	Cleaning Sheets for ScanX Classic (12 pack)	\$ 159
B2020	Cleaning Sheets for ScanX Classic (25 pack)	\$ 271
B3010	Cleaning Sheets for ScanX Intraoral (12 pack)	\$ 128
B3020	Cleaning Sheets for ScanX Intraoral (25 pack)	\$ 220
D1010	Cleaning Sheets for ScanX Duo (12 pack)	\$ 71
D1020	Cleaning Sheets for ScanX Duo (25 pack)	\$ 136

MERCHANDISE

Chemistry

Peri-Pro Chemistry PN 90800	Air Techniques Developer & Fixer PN 40270	Film Carriers (#0 - #4) PN 90980, 90985, 90990, 90831	Air Techniques Starter PN 40280

PART NO.	DESCRIPTION	PRICE
90980	Pedo Film Carrier #0 (Box of 5)	\$ 60
90985	Adult Anterior Film Carrier #1 (Box of 5)	\$ 60
90990	Bite-Wing Film Carrier #3 (Box of 5)	\$ 60
90831	Occlusal Film Carrier #4 (Box of 3)	\$ 60
90800	Peri-Pro Chemistry Case includes: (3) one-quart bottles of Developer (3) one-quart bottles of Fixer	\$ 53
90801	Peri-Pro Developer (1 Quart)	\$ 15
90802	Peri-Pro Fixer (1 Quart)	\$ 15
40270	Air Techniques Developer & Fixer For roller processors only. Case includes 4 bottles: (2) 1-gal. Developer, (2) 1-gal. Fixer, 2-color coded sponges	\$ 98
40280	Air Techniques Chemistry Starter For use with Air Techniques Developer. Case includes: (6) 12-ounce bottles	\$ 75
44170	Chemistry Recovery Containers For capture and recovery of developer and fixer. Includes: (2) 2½-gal. containers For A/T 2000®, A/T 2000® PLUS, A/T 2000®XR	\$ 122
44160	Chemistry Recovery Containers For capture and recovery of developer and fixer. Includes: (2) 5-gal. containers For A/T 2000®, A/T 2000® PLUS, A/T 2000®XR	\$ 142
40140	Cleaning Film For automatic X-ray film processors and other roller transport processors. Package of 12 films.	\$ 72
40240	A/T 2000 Developer - Case includes: (4) 1-gal. bottles	\$ 105
40245	A/T 2000 Fixer - Case includes: (4) 1-gal. bottles	\$ 105



Peri-Pro Transport Cleaning Kit
PN 43975



Formula 2000 Plus
Cleanser Twin-Pack PN 43945



Spray 2000 Plus
PN 43965

PART NO.	DESCRIPTION	PRICE
43975	Peri-Pro Transport Cleaning Kit Includes: All items needed to chemically clean a Peri-Pro Transport along with (1) Cleaning Tub, (1) Brush, and (1) Twin-Pack of Formula 2000 Plus.	\$ 145
43945	Formula 2000 PLUS Cleanser Twin-Pack Use for roller-type processors, manual processing tanks & Peri-Pro transports. Includes: (2) Bottles of Formula 2000 Plus.	\$ 67
43965	Spray 2000 PLUS* Kit Includes: (1) liter Spray 2000 Plus, (1) Spray Attachment, (4) Cleaning Webs (non-abrasive scrubbing pads for cleaning roller transports).	\$ 53

* For weekly cleaning of A/T 2000 roller transports only; not for Peri-Pro.



Monarch®



*Monarch products are not sold by our Professional Sales Associates (PSA) team.

IMAGING | UTILITY ROOM | MERCANDISE



Reach for Monarch® by Air Techniques, a simply smart way to clean and disinfect your office, your instruments and your equipment.

Monarch's complete, color-coded line of products helps make it easy to use the right product for each job.

- **Green for Surfaces**
- **Blue for Instrument Cleaning**
- **Pink to protect Skin and Hands**
- **Yellow for your Equipment**

Packaging and accessories that make life easier. Unit dose packs to simplify filling baths. Push-button dispenser for easy refills. Safety trays that really protect against instrument sticks. Dispensers for no-fuss, no-muss use.

But mostly, the focus is on you. Effective formulas with your patients and your team in mind. Packaging that meets your environmental demands. Efficient usage that keeps costs down and increases your profitability.



MERCHANDISE

Surface Disinfectant



Monarch Surface
Disinfectant Spray
PN H6110

Monarch Surface
Disinfectant Dispenser Refill
PN H6110R-3

Monarch Surface Disinfectant provides broad spectrum kill of organisms in 1 minute. Our unique spray & wipe application means that no pre-cleaning is necessary saving time and money.

Monarch Surface Disinfectant Spray

- Broad spectrum kill of key organisms.
- Free of aldehydes and phenols.
- Compatible with hard surfaces in the dental office. Non-staining. Dries without residue.
- No strong chemical odor.

*Monarch products are not sold by our Professional Sales Associates (PSA) team.



Monarch Surface Disinfectant
Premoistened Wipes
PN H6171



Hand to Wipe Comparison

Monarch Surface Disinfectant Wipes

- Kills TB, MRSA, HIV, HBV, Polio and 20 other pathogens in just 1 minute.
- Ethanol-based and bleach-free.
- Kills odor-causing bacteria.
- Large 7" x 9" size.

PART NO.	DESCRIPTION	PRICE (Individual)
H6110	Monarch Surface Disinfectant Spray 710 mL / 24 oz. (10 bottles per shipper case)	\$ 13.50
H6110R-3	Monarch Surface Disinfectant Dispenser Refill 3 L / 0.8 gal. Dispenser. (4 dispenser refills per shipper case)	\$ 32.25
H6171	Monarch Surface Disinfectant Premoistened Wipes 1 Canister with Wipes (100 ea. 7" x 9")	\$ 12.95
H6182	Dry Wipes Intro Kit Includes 2 wipe canisters and 4 wipe refills (100 ea. 7" x 9")	\$ 29.95
H6183	Dry Wipes Refill Kit Includes 12 wipe refills (100 ea. 7" x 9")	\$ 49.95

MERCHANDISE

Instrument Cleaner



Monarch Enzymatic Cleaner
(3 qts. / Box of 20 Unit Dose Pouches)
PN H6310

Monarch Enzymatic Cleaner is an effective multi-enzymatic cleaning concentrate specially formulated for use in ultrasonic and holding baths.

- 4 synergistic enzymes for complete instrument cleaning prior to sterilization:
 - Proteases** — Help break up, dissolve and remove proteins and protein stains from instruments.
 - Amylases** — Break down and remove carbohydrates, carbon compounds and resistant starch residue.
 - Lipases** — Assist in the removal of fatty materials and oily stains.
 - Cellulases** — Break down cellulose, plant structures and fibers.
- Leaves instruments sparkly clean.
- Convenient unit dose packaging; no measuring and dissolves instantly.
- So effective, it can be used for up to two days (or until cloudy).
- Deodorizes plaster traps.

*Monarch products are not sold by our Professional Sales Associates (PSA) team.



Monarch Safety Tray
PN H6150

Monarch Safety Tray

- Facilitates easy and safe instrument processing.
- Removable perforated tray allows for easy rinsing, draining and pick-up of instruments.
- Hinged cover protects staff from bath-generated aerosols.
- Constructed of heavy gauge, high impact plastic.
- Holds 3 quarts; fill lines designate 1-quart increments.

PART NO.	DESCRIPTION	PRICE (Individual)
H6310	Monarch Enzymatic Cleaner for Standard (3 qts.) Ultrasonic Cleaners 20 x 50 mL / 1.7 oz. unit dose pouches each makes 3 qts. (10 boxes of 20 pouches each per shipper case)	\$ 32.65
H6150	Monarch Safety Tray 3 Quart Bath.	\$ 86.95

MERCHANDISE

Hand Sanitizer



Monarch Hydrating Instant
Hand Sanitizer (Large)
PN H6500



Monarch Hydrating Instant
Hand Sanitizer (Travel Size)
PN H6510



1
Rub hands palm to palm, including wrists.



2
Rub back of hands with interlaced fingers.



3
Palm to palm with fingers interlaced.



4
Interlock fingers to opposing palms.



5
Rotate left thumb clasped in right hand and reverse.



6
Rub backward and forward with clasped fingers into both palms.

Technique adapted from Berufsgenossenschaft für Gesundheitsdienst und Wohlfahrtspflege - BGW Essential skin protection; First published 02/2005; Revision 05/2010

Monarch Hydrating Instant Hand Sanitizer

- Waterless, alcohol-based gel for use in between patients or when hands are not visibly soiled.
- One pump is sufficient to cover hands.
- Smooth application and silky feel.
- Microbead polymer formula with Aloe and Vitamin B.
- Skin-enriching emollients prevent dryness during repeated use.
- Refreshing apple scent.

PART NO. DESCRIPTION

PRICE (Individual)

H6500	Monarch Hydrating Instant Hand Sanitizer (Large) 450 mL / 15.2 oz. Pump Bottle. (6 bottles per shipper case)	\$ 12.95
H6510	Monarch Hydrating Instant Hand Sanitizer (Travel Size) 60 mL / 2 oz. Pump Bottle. (10 bottles per shipper case)	\$ 3.95

*Monarch products are not sold by our Professional Sales Associates (PSA) team.



Monarch Revitalizing
Hand Lotion (Large)
PN H6550



Monarch Revitalizing
Hand Lotion (Travel Size)
PN H6560



1
Spread lotion
evenly over both
hands.



2
Massage
completely
into skin.



3
Including cuticles
and fingertips.



4
Don't forget the
wrists!

Technique adapted from Berufsgenossenschaft für Gesundheitsdienst und Wohlfahrtspflege - BGW Essential skin protection; First published 02/2005; Revision 05/2010

Monarch Revitalizing Hand Lotion

- Enriched with Aloe, Vitamin E and B complex.
- Latex compatible (petroleum free).
- Leaves hands smooth and silky; never oily.
- Special wax formulation makes donning gloves easier.
- Residual glove powder will not adhere to the skin.
- Long-lasting moisturizing effect, even after several hand washes.
- Clean, orange scent.
- Free of parabens and sulfates.

PART NO. DESCRIPTION

PRICE (Individual)

H6550	Monarch Revitalizing Hand Lotion (Large) 450 mL / 15.2 oz. Pump Bottle. (6 bottles per shipper case)	\$ 34.95
H6560	Monarch Revitalizing Hand Lotion (Travel Size) 60 mL / 2 oz. Pump Bottle. (10 bottles per shipper case)	\$ 7.95

MERCHANDISE

Evacuation System Cleaner



Monarch CleanStream
Evacuation System Cleaner
PN 57850



Monarch CleanStream Dispenser
PN 57665

Directions/Replacement
Parts Card Insert



Monarch CleanStream Starter Kit
PN 57660


Monarch CleanStream® designed for use with Air Techniques and all vacuum systems, wet or dry.

- Highly effective, non-foaming concentrated formula.
- Removes line build-up, increasing vacuum performance and suction.
- Cleans and deodorizes evacuation lines.
- Fresh scent.
- Non-corrosive and biodegradable.
- 8.5 pH makes CleanStream safe for use in all amalgam separators.

*Monarch products are not sold by our Professional Sales Associates (PSA) team.



Monarch Lines
PN H6340



Build-up
on internal
walls in just
3 weeks.

Clean,
smooth
walls. No
build-up.



Monarch Lines is the perfect solution to keep water lines free of deposits and build-ups including minerals, lime scale, hard water, salt and similar contaminants.

- For daily use in water bottle, weekly cleaning, and intermittent treatment to clear deposits from waterlines.
- Ready to use; no mixing or diluting is required.
- Ideal for use in your ultrasonic scaler, helps keep air-water spray clean.
- Decreases bad odors and tastes that occur during hygiene procedures.
- Does not affect bond strength.

PART NO.	DESCRIPTION	PRICE (Individual)
H6340	Monarch Lines 10 x 500 mL / 16.9 oz. bottles. (Case Price; not sold individually)	\$ 122.50
57850	Monarch CleanStream Evacuation System Cleaner 2.5 liter bottle (125 applications). (4 bottles per shipper case)	\$ 69.95
57665	Monarch CleanStream Dispenser System Includes: 1 Dispenser system with accessories.	\$ 67.95
57660	Monarch CleanStream Starter Kit Includes: 1 Dispenser System with all accessories and 12 applications of CleanStream Evacuation System Cleaner.	\$ 72.50

WARRANTY & RETURN POLICY

All parts and component returns and replacement of equipment under warranty require a Return Materials Authorization (RMA).

Warranty credit will not be applied to product returns that exhibit damage due to shipping, misuse, careless handling or repairs by unauthorized service personnel. Credit, or partial credit, will not be issued until product/parts have been received and evaluated. Warranty is void if product is installed or serviced by anyone other than authorized Air Techniques dealer service personnel.

Dealer labor, shipping & handling charges are not covered by this warranty.

Air Techniques will issue partial credit without a restocking fee if product is returned in a resalable condition, unopened packaging, and received within 90 days of invoice date. Product received older than 90 days will be charged a 15% restocking fee. Any product older than 12 months from invoice date and/or any products with a tamper-resistant seal (including electrical products, PSP's, barrier envelopes, sheaths, etc.) that have been opened, cannot be returned. Monarch products must be returned within 90 days in order to receive credit.

Air Techniques' equipment is warranted to be free from defects in material and workmanship from the date of installation for a period of:

Product	Full System	Specific Components
Acadia	12 months	n/a
A/T2000 XR	12 months	24 months: electronic modules
Peri-Pro	12 months	n/a
ScanX Classic	24 months	n/a
ScanX Intraoral	24 months	n/a
ScanX Duo	24 months	n/a
ScanX Swift	24 months	n/a
ProVecta S-Pan	24 months	n/a
ProVecta HD	24 months	n/a
Cam-X Spectra	36 months	n/a
Cam-X Luna HD	36 months	n/a
Cam-X Polaris	36 months	n/a
Mojave	24 months	60 months: all pumps and motors
AirStar	24 months	60 months: all compressor heads and motors
AirStar NEO	24 months	60 months: all compressor heads and motors
Vacstar	24 months	60 months: all pumps and motors
Hydromiser	24 months	n/a

Contact Information

For general and immediate questions please call:

Toll Free in the US and Canada: 1- 800-AIRTECH (247-8324)

Monday through Friday, 8:00am – 8:00pm (EST)

Disclaimer

1. Items returned without an RMA, or included with other products for which an RMA has been issued, may be returned to the customer at the discretion of Air Techniques.
2. Advanced replacement will be provided and invoiced at full price.
3. For Shipped Items - All freight costs are paid for by dealers, unless part of an approved promotion/program.
All shipping damage costs should be reviewed and addressed with carrier or customer/dealer as needed.
4. For all non-warranty orders, there is a \$25 minimum order size for US dealers and \$40 for Canadian dealers.

NOTES

NOTES

We Serve You NATIONWIDE



Victor Yoon
Cell: 778.688.8429
vyoon@profsales.com



Scott Rumble
Cell: 403.801.5390
srumble@profsales.com



Rick Orme
Cell: 916.837.3839
rorme@profsales.com



John Miller
Cell: 503.758.9895
jmiller@profsales.com



Jim Schlosser
Cell: 909.215.3277
jschlosser@profsales.com



Pat Moylan
Cell: 816.729.2900
pmoylan@profsales.com



Sam Paola
Cell: 310.569.1178
spaola@profsales.com




Paul Myers
Cell: 214.704.5100
pmyers@profsales.com



Mel Seagle
Cell: 949.678.4550
mseagle@profsales.com



Skip Ellender
Cell: 617.894.1224
sellender@profsales.com




Michael Park
Cell: 281.702.3992
mpark@profsales.com



Steve Carlin
Cell: 952.491.1223
scarlin@profsales.com

Professional Sales Associates

At Your Service



Aaron Renner
Cell: 519.319.7457
arenner@profsales.com



Michel Cabanie
Cell: 514.793.3368
mcabanie@profsales.com



Jim Domitrz
Cell: 617.921.4580
jdomitrz@profsales.com



Don Rogers, III
Cell: 717.201.7332
drogers@profsales.com



Rich Ellis
Cell: 301.502.7424
rellis@profsales.com



Matt Langdon
Cell: 516.316.5913
mlangdon@profsales.com



Milt Baisley
Cell: 203.645.7431
mbaisley@profsales.com



Tom O'Berry
Cell: 404.307.4931
toberry@profsales.com



Charlie Clarke
Cell: 859.509.6298
cclarke@profsales.com



Bill Highland
Cell: 810.623.7297
bhighland@profsales.com



Bill Reinhardt
Cell: 847.840.3330
breinhardt@profsales.com



Matt Highland
Cell: 901.606.3633
mhigland@profsales.com



Keith O'Berry
Cell: 336.210.1429
koberry@profsales.com



Robert Rodriguez
Cell: 305.458.6538
rrodriguez@profsales.com



For over 50 years, **Air Techniques** has been a leading innovator and manufacturer of dental products. Our priority is ensuring complete satisfaction by manufacturing reliable products and providing excellent customer and technical support. Whether the need is digital imaging, utility room equipment or merchandise, Air Techniques can provide the solution via our network of authorized professional dealers.

Proudly designed, tested and manufactured in the United States, our products are helping dental professionals take their practices to the next level.

Air Techniques' products are distributed through authorized dealers only.

Corporate Headquarters

1295 Walt Whitman Road | Melville, NY 11747-3062 | Phone 800.247.8324 | Fax 888.247.8481

Western Facility

291 Bonnie Lane, Suite 101 | Corona, CA 92880-2804 | Phone 800.247.8324 | Fax 951.898.7646

www.airtechniques.com



**AIR
TECHNIQUES**
equipped for life®