


CHAPMAN-HUFFMAN, INC.

FOOTSWITCH AND CABLE ASSEMBLY - 4 PIN

12-060-00



- A footswitch assembly heavy enough to stay where you put it, no chasing it around with your foot.
- Same light foot pressure to operate as original equipment switches
- Square 4 pin male plug as illustrated, with 8 1/2 foot cable
- Two position switch
- Fits a variety of DENTSPLY products
- DENTSPLY cross ref. # 60473

FOOTSWITCH AND CABLE ASSEMBLY - 3 PIN


12-070-00



- A footswitch assembly heavy enough to stay where you put it, no chasing it around with your foot.
- Same light foot pressure to operate as original equipment switches
- Rectangular 3 pin male plug as illustrated, with 8 1/2 foot cable
- Single position switch
- Fits DENTSPLY products
- DENTSPLY cross ref. # 61187

FOOTSWITCH AND CABLE ASSEMBLY - BOBCAT TYPE


12-080-00



- A footswitch assembly heavy enough to stay where you put it, no chasing it around with your foot.
- Same light foot pressure to operate as original equipment switch
- Round 3 pin male plug as illustrated, with 8 1/2 foot cable
- Single position switch
- Fits DENTSPLY BOBCAT
- DENTSPLY cross ref. # 77885

FOOTSWITCH AND CABLE ASSEMBLY

12-085-00




- A footswitch assembly heavy enough to stay where you put it, no chasing it around with your foot.
- Same light foot pressure to operate as original equipment switch
- Round 4 pin male plug as illustrated, with 8 1/2 foot cable
- Two position switch
- Fits DENTSPLY PROPHY-JET 122
- DENTSPLY cross ref. # 565138003

CHAPMAN-HUFFMAN, INC.

FOOTSWITCH AND CABLE ASSEMBLY - SPS TYPE

12-090-00




- A footswitch assembly heavy enough to stay where you put it, no chasing it around with your foot.
- Same light foot pressure to operate as original equipment switch
- Round 4 pin male plug as illustrated with 8 1/2 foot cable
- Two position switch
- Fits DENTSPLY CAVITRON SPS or CAVI-JET with SPS Technology
- DENTSPLY cross ref. # 565138002

FOOTSWITCH AND CABLE ASSEMBLY - CAVITRON SELECT STYLE

12-088-00


NEW!



- A footswitch assembly heavy enough to stay where you put it, no chasing it around with your foot.
- Same light foot pressure to operate as original equipment switch
- Fits DENTSPLY CAVITRON **SELECT**
- Single position
- Round 1 pin male plug as illustrated, with 8 foot long cable
- DENTSPLY cross ref. # 565162002

FOOT CONTROL ASSEMBLY - CAVITRON SELECT SPS TYPE

12-095-00



- A footswitch assembly heavy enough to stay where you put it, no chasing it around with your foot.
- Same light foot pressure to operate as original equipment switch
- Round pin male plug as illustrated with 8 foot cable
- Two position switch
- Fits DENTSPLY CAVITRON **Select SPS**
- DENTSPLY cross ref. # 565162005

FOOT CONTROLS

800-225-4621

CHAPMAN-HUFFMAN, INC.

ECONOMY DISC TYPE FOOT CONTROL



Gray

21-005-00

Sterling

21-405-00

- Same features as 21-010-00 standard foot control
- 6" diameter base plate
- Supplied with 7' of 2 hole asepsis tubing

STANDARD DISC TYPE FOOT CONTROL



Gray

21-010-00

Sterling

21-410-00

- Variable speed
- Disc-type foot control for use with air turbine controls
- Pressure may be applied at any point
- 6" in diameter and 1 1/2" high
- Supplied with 7' of 2 hole asepsis tubing

CHIP BLOWER FOOT CONTROL



Gray

21-015-00

Sterling

21-415-00

- Same features as 21-010-00 standard foot control
- Push button chip blower valve
- Supplied with 7' of 3 hole tubing

WET/DRY FOOT CONTROL



Gray

21-020-00

Sterling

21-420-00

- Same features as 21-010-00 standard foot control
- On-Off toggle for coolant water
- Supplied with 7' of 3 hole tubing

WET/DRY WITH CHIP BLOWER FOOT CONTROL



Gray

21-025-00

Sterling

21-425-00

- Same features as 21-010-00 standard foot control
- Push button chip blower valve
- On-Off toggle for coolant water
- Supplied with 7' of 4 hole asepsis tubing

SCALER TYPE FOOT CONTROL



Gray

21-030-00

Sterling

21-430-00

- Same features as 21-010-00 standard foot control
- Push button valve for ON-OFF of ultrasonic scaler
- Air electric switch included
- Supplied with 7' of 3 hole tubing

CHAPMAN-HUFFMAN, INC.

RELAY TYPE FOOT CONTROLS

STANDARD DISC TYPE FOOT CONTROL WITH RELAY



Gray

21-210-00

Surf 6

21-310-00

- Disc type foot control with relay for use with air turbine controls
- Pressure may be applied at any point
- 6" in diameter and 1 1/2" high
- Supplied with 7' of ADEC style asepsis 5 hole color and dash coded internal tubings
- ADEC cross ref. # 38.1700.00 Surf 6

CHIP BLOWER FOOT CONTROL WITH RELAY



Gray

21-215-00

Surf 6

21-315-00

- Same features as 21-210-00 Standard Disc Type Foot Control with relay
- Push button chip blower valve
- Supplied with 7' of ADEC style asepsis 5 hole color and dash coded internal tubings

WET/DRY FOOT CONTROL WITH RELAY



Gray

21-220-00

Surf 6

21-320-00

- Same features as 21-210-00 Standard Disc Type Foot Control with relay
- On-Off toggle for coolant water
- Supplied with 7' of ADEC style asepsis 5 hole color and dash coded internal tubings
- ADEC cross ref. # 38.1603.00

WET/DRY CHIP BLOWER FOOT CONTROL WITH RELAY



Gray

21-225-00

Surf 6

21-325-00

- Same features as 21-210-00 Standard Disc Type Foot Control with relay
- Push button chip blower valve
- On-Off toggle for coolant water
- Supplied with 7' of ADEC style asepsis 5 hole color and dash coded internal tubings
- ADEC cross ref. # 38.1601.00 Surf 6

FOOT CONTROLS

800-225-4621

CHAPMAN-HUFFMAN, INC.

SCALER TYPE FOOT CONTROL WITH RELAY



Gray
Surf 6

21-230-00
21-330-00

- Same features as 21-210-00 Standard Disc Type Foot Control with relay
- Push button valve for ON-OFF of Ultrasonic Scaler
- Air electric switch included
- Supplied with 7' of ADEC style asepsis 5 hole color and dash coded internal tubings

FOOT CONTROL RING KIT



95-180-00

- Replacement for all **Non-relay Type** Chapman-Huffman Foot Controls
- Gray plastic

RELAY TYPE FOOT CONTROL REPLACEMENT RING KIT



95-280-00

- Replacement for all **Relay Type** Chapman-Huffman Foot Controls
- Black plastic
- Includes four mounting screws
- Also fits some ADEC foot controls cross ref. # 38.0237.00

WATER TOGGLE TOWER ASSEMBLY



Gray
Sterling
Surf 6

21-160-00
21-460-00
21-360-00

- Fits all Chapman-Huffman Foot Controls with **Water Toggle Tower Assembly**
- 3-way valve

UNIVERSAL PILOT VALVE ASSEMBLY



19-035-00

- Normally closed 3-way
- Replacement pilot valve for Water Toggle Tower Assemblies 21-160-00, 21-460-00 and 21-360-00
- Replacement pilot valve for Scaler Tower Assemblies 21-267-00, 21-467-00 and 21-367-00
- (Also fits Automatic Handpiece Holders 25-010-00, 25-011-00, 25-012-00, 25-070-00 and 25-080-00)

CHAPMAN-HUFFMAN, INC.

CHIP BLOWER TOWER ASSEMBLY N/C 2-WAY VALVE



Gray	21-265-00
Sterling	21-465-00
Surf 6	21-365-00

- Fits all Chapman-Huffman Foot Controls with **Chip Blower**
- Replaces Obsolete # 21-165-00 beg. June 2005

PILOT VALVE N/C 2-WAY



19-435-00

- Fits Chapman-Huffman Foot Control **Chip Blower Tower Assemblies**
- Normally Closed 2-way valve

SCALER TOWER ASSEMBLY



Gray	21-267-00
Sterling	21-467-00
Surf 6	21-367-00

- Fits all Chapman-Huffman **Scaler** Foot Controls
- 3-Way valve

FOOT CONTROL TOGGLE WITH PIN



Gray	21-161-00
Sterling	21-461-00
Surf 6	21-361-00

- Fits all Chapman-Huffman Foot Controls with Water Toggle Tower Assembly

FOOT CONTROL POPPET



Pkg of 10	95-181-00
------------------	------------------

- Fits Chapman-Huffman Foot Controls 21-005-00, 21-010-00, 21-015-00, 21-020-00, 21-025-00, 21-030-00, 21-210-00, 21-215-00, 21-220-00, 21-225-00, 21-230-00, 21-310-00, 21-315-00, 21-320-00, 21-325-00, 21-330-00, 21-405-00, 21-410-00, 21-415-00, 21-420-00, 21-425-00, and 21-430-00
- Also fits ADEC Foot Control 38.0050.00
- ADEC cross ref. # 22.0060.00

800-225-4621

CHAPMAN-HUFFMAN, INC.

FOOT CONTROL REPAIR KIT

95-185-00



- For repair of Chapman-Huffman foot controls.
- Economy Foot Control 21-005-00 and 21-405-00

FOOT CONTROL REPAIR KIT

95-186-00



- For repair of Chapman-Huffman foot controls.
- Standard Foot Controls 21-010-00 and 21-410-00
- Chip Blower Foot Controls 21-015-00 & 21-415-00
- Wet/Dry Foot Controls 21-020-00 and 21-420-00
- Wet/Dry with Chip Blower Foot Controls 21-025-00 and 21-425-00
- Scaler Type Foot Controls 21-030-00 & 21-430-00

RELAY TYPE FOOT CONTROL VALVE ASSEMBLY

21-240-00



NEW!

- Fits all Chapman-Huffman **Relay** Type Foot Controls
- Mounting hardware and fittings included
- Replaces ADEC 38.0607.01 Foot Control II Valve and 38.0659.00 Foot Control I Valve
- **Not for ADEC 500 Foot Control**
- ADEC cross ref. # 38.0761.00

RELAY TYPE FOOT CONTROL REPAIR KIT

95-287-00



NEW!

- Fits all Chapman-Huffman **Relay** Type Foot Controls
- For use with ADEC Foot Control III Valve 38.0761.00
- Also fits ADEC 500 Foot Control
- ADEC cross ref. # 90.0593.00