NEW - DENTAL LIGHT REFLECTOR

16-800-00



- Fits PELTON & CRANE LFI and LFI+
- PELTON & CRANE cross ref. # 022438
- New manufactured reflector
- · Four year warranty on optics coating

NEW - DENTAL LIGHT REFLECTOR

16-805-00



- Fits PELTON & CRANE LFII
- PELTON & CRANE cross ref. # 3321887
- New manufactured reflector
- Four year warranty on optics coating

NEW - DENTAL LIGHT REFLECTOR

16-810-00



Retainer Ring 16-716-00 (For retaining lens splash shield)

- Fits models with handles
- RITTER J STARLITE cross ref. # 148774
- RITTER K SUPER STARLITE cross ref. # 148774
- HENRY SCHEIN cross ref. # 148774
- · New manufactured reflector
- Four year warranty on optics coating

NEW - DENTAL LIGHT REFLECTOR

16-815-00



- Fits models without handles
- KNIGHT L cross ref. # 152615
- MARUS cross ref. # 52R412
- HENRY SCHEIN cross ref. # 52R412
- · New manufactured reflector
- · Four year warranty on optics coating

CCESSORIES



- Fits PELTON & CRANE LFI & LFI+
- Cross ref. # 006777

LENS SPLASH SHIELD

16-705-00



- Fits PELTON & CRANE LFII
- Two Piece set
- Cross ref. # 3321598

LENS SPLASH SHIELD - RITTER

16-710-00



• Fits RITTER "J" STARLITE & RITTER "K" Super STARLITE

LENS SPLASH SHIELD - SCHEIN / MARUS

16-715-00



- Fits SCHEIN/MARUS Lights
- Cross ref. # 52R430

RETAINER RING

16-716-00



• Fits RITTER, SCHEIN, and MARUS older lights for retaining lens splash shield

LENS SPLASH SHIELD

16-725-00



• Fits HEALTHCO Aristocrat Light

LENS SPLASH SHIELD



• Fits BELMONT Lights 046, 047, 048 and 049

MECHANICAL SPRING UNIT MOUNT ASSEMBLY

16-660-00

16-650-00

16-750-00

16-720-00



• Fits PELTON & CRANE LFII unit mount lights and LFWII Wall Mount Light

• PELTON & CRANE cross ref. # 94 39 738

MECHANICAL SPRING CEILING MOUNT ASSEMBLY



 Fits PELTON & CRANE ceiling mount lights LFCII, LFCII-D, LFTII, LFTII-D and LFWII Wall mount lights

• PELTON & CRANE cross ref. # 15 32 907

GAS SPRING ADJUSTMENT TOOL



- Fits PELTON & CRANE LFII lights
- 5/16 hex
- Used to assist in removal of gas spring or its adjustment
- PELTON & CRANE cross ref. # 016749

SCREW KIT



95-068-00

- Fits PELTON & CRANE LFII Unit and Ceiling mount light arms
- Four 10-24 x 1/2 flat head socket cap screws
- PELTON & CRANE cross ref. # 017668

10-24 FLAT HEAD SOCKET CAP SCREW - BLACK



PELTON & CRANE Cross Ref. #S

090721 10-24 X 3/8 001-47-00 33 23 115 10-24 X 1/2 001-53-00 33 24 170 10-24 X 3/4 001-55-00

Fits PELTON & CRANE LFII & LFCII-D light arms

ARM ADAPTER KIT



95-075-00

- Fits PELTON & CRANE lights LFII, LFCII & LFCII-D, rear arm to front arm assembly
- Not painted, Aluminum construction
- Includes (2) set screws
- PELTON & CRANE cross ref. # (Old #) HB4-014530 (New #) 33 22 620

SNUBBING PIN KIT



14-450-00

• Fits PELTON & CRANE LFII, LFCII & LFCII-D Lights

SNUBBING PIN



14-452-00

- Fits PELTON & CRANE LFII, LFCII & LFCII-D Lights
- Nickel plated steel construction
- PELTON & CRANE cross ref. # 33 22 992/016705

WAVE WASHER



14-454-00

- Fits PELTON & CRANE LFII, LFCII & LFCII-D Lights
- PELTON & CRANE cross ref. # 33 23 008/017186

FLAT WASHER



14-456-00

- Fits PELTON & CRANE LFII, LFCII & LFCII-D Lights
- PELTON & CRANE cross ref. # 33 23 016/016189

14-435-00

ACCESSORIES

KNUCKLE HINGE SCREW

CHAPMAN-HUFFMAN, INC.



- Fits PELTON & CRANE LFII, LFCII & LFCII-D Lights
- PELTON & CRANE cross ref. # 33 23 073

REFLECTOR PIVOT BUSHING



Pkg of 2

16-625-00

95-077-00

12-100-00

20-150-00

- Fits PELTON & CRANE LFII light yoke and frame assembly
- Plastic construction
- PELTON & CRANE cross ref. # (Old) 014901 (Pkg 1) (New) 023490 (Pkg 2)

DRIFT SCREW KIT



- Fits PELTON & CRANE LF and LF+ series
- PELTON & CRANE cross ref. # 006804

LIGHT SOCKET



• Fits PELTON & CRANE LFI Light

• PELTON & CRANE cross ref. # 007028


LIGHT SOCKET



12-580-00

- Fits PELTON & CRANE LFII Light
- PELTON & CRANE cross ref. # 33 36 992

HANDPIECE OIL COLLECTOR AND MUFFLER ASSEMBLY



- · Filter element is a common gauze found in the dental office
- · Minimum back pressure
- Can mount through unit base plate so jar can be removed for cleaning & filter element replacement
- ADEC cross ref. # 24.0410.00

20-075-00

- KAVO Style handpiece Coupler
- 4 hole standard **Non-Fiber Optic** Swivel Coupler
- KAVO cross ref. # 457/553-1210

SWIVEL COUPLER



20-077-00

- KAVO Style Handpiece Coupler
- 5 hole Fiber Optic Standard Swivel
- Used with a fiber optic tubing having a fiber optic bundle to convey light to the Coupler & Handpiece
- KAVO cross ref. # 1390/553-1390

SWIVEL COUPLER



20-079-00

- KAVO Style Handpiece Coupler
- 6-Pin Bulb in Coupler Style Swivel
- · Bulb included with Coupler
- Fits Chapman-Huffman 6-Pin tubings and others
- KAVO cross ref. # 465LRN/553-1550

WASTE DROP 4 INCH



20-275-00

- This 4" diameter stainless steel waste drop fits through a 4 1/8" hole in the counter top. This is a nice transition piece between the counter top and waste basket below. The waste drop can be easily removed for cleaning.
- Stainless steel construction

WASTE DROP 6 INCH



20-276-00

- This 6" diameter stainless steel waste drop fits through a 6 1/8" hole in the counter top. This is a nice transition piece between the counter top and waste basket below. The waste drop can be easily removed for cleaning.
- Stainless steel construction

CUP FILLER ASSEMBLY




800-225-4621

50-075-00

- Fits 5 oz. cup
- Two inch post mounting
- Push and release water on-off filler tube
- Filler tube rests out of the way of cup rim
- 12 ft. supply tubing

CHEMICLAVE REPLACEMENT DOOR GASKET

20-250-00



- Fits HARVEY MDT 4000, 5000, "E" & Aquaclave 10
- MDT cross ref. # 250007801

20-251-00

- Fits MDT Chemiclave 5500, EM5500, 7000 and Aquaclave 20
- MDT cross ref. # 280000301

DOOR GASKET -1730

20-400-00



- Fits TUTTNAUER Autoclave model # 1730
- TUTTNAUER cross ref. # 02610020
- Manufactured by Chapman-Huffman to fit TUTTNAUER Autoclave

DOOR GASKET - 2340

20-401-00



CROSS SECTION VIEW

- Fits TUTTNAUER Autoclave model # 2340 units with serial # 8806 and higher
- TUTTNAUER cross ref. # 02610118
- Manufactured by Chapman-Huffman to fit TUTTNAUER Autoclave

DOOR GASKET - 2540



CROSS SECTION VIEW 20-402-00

- Fits TUTTNAUER Autoclave model # 2540
- TUTTNAUER cross ref. # 02610023
- Manufactured by Chapman-Huffman to fit TUTTNAUER Autoclave

AIR FILTER

20-070-00



- Fits PELTON & CRANE Hustler I & II
- Cross ref. # C42-005127

SALIVA EJECTOR VALVE CLEANING BRUSH



95-300-00

- 1/2" diameter brush
- Nylon wire twisted brush
- Package of 5

VACUUM VALVE CLEANING BRUSH



95-301-00

- 3/4" diameter brush
- · Nylon wire twisted brush
- Package of 5

SURGICAL TIP CLEANING BRUSH



95-302-00

- 3/16" diameter brush
- · Nylon wire twisted brush
- Package of 5

RIGHT ANGLE AIR CHUCK



For pressurizing valve core type filler stems

- Brass construction
- 1/4" FPT

AIR BLOWGUN

20-125-00

019-21-00



- Lever operated
- 1/4 FPT inlet
- 125 PSI working pressure
- Standard safety tip (venturi holes increase air volume & allow air to escape if nozzle is dead-ended)
- · Maximum of 30 PSI outlet pressure
- · Conforms to OSHA standards

AIR BLOWGUN WITH TUBING



Blowgun w/ Coiled 10 Ft. Tubing

 Neon Green
 20-126-00

 Gray
 20-127-00

Coiled 10 Ft. Tubing w/ Fittings

 Neon Green
 15-700-00

 Gray
 15-225-00

 Same features as Blowgun 20-125-00 plus 10 ft. coiled hose with 1/4 MPT swivel fitting both ends

800-225-4621

ACCESSORIES

AIR COMPRESSOR SUPPLY HOSE ASSEMBLY



3 Ft. Hose Assembly 50-003-00

6 Ft. Hose Assembly 50-006-00

- For new and existing compressor installation between air compressor and building air supply line
- Flexible 1/2" ID, 250 PSI blue high pressure air hose
- 1/2" MPT swivel couplings both ends

AIR COMPRESSOR TANK DRAIN



120 Volt Model 50-020-00

240 Volt Model 50-021-00

- · Drains moisture buildup in air compressor tank
- Fully automatic
- Adjustable drain open time from 1 sec. to 500 sec.
- Purges every time compressor cycles
- Unit is off when compressor is off
- · Valve will not foul with oil like other drain systems
- Electrical source from pressure switch
- Includes service valve, screen and installation kit
- · All metal construction
- Used on dental compressors w/ cycles of 80-175 PSI

STATIM 2000 TYPE COMPRESSOR

12-135-00



- 115 VAC 50/60 Hz 0.28A
- 2.84 psig
- Statim 2000 cross ref. # 01-100560S
- Also fits PELTON & CRANE Delta 8 & Delta 10 Autoclave

WATER BOTTLE PICK-UP TUBES



750 ml Bottle (Pkg. of 6) 95-090-00

Large Bottle (Pkg. of 6) 95-095-00

AUXILIARY CLEAN WATER SUPPLY SYSTEM



Assembly with 750 ml Bottle 50-100-00
Assembly with Large Capacity Bottle 50-200-00
750 ml Heavy Duty Bottle Only 50-060-00
Large Capacity Heavy Duty Bottle Only 50-260-00
Gasket 50-054-00

- · For antiseptic coolant or water
- · Easy to install with installation instructions
- L bracket mounts to any flat vertical surface. Order 50-030-00 2" post mount kit for mounting to post
- · On-off air supply toggle valve, with tubing
- Regulator, preset at factory
- · Water pressure gauge
- Heavy Duty Plastic Bottle is not autoclavable

IMPORTANT! AUXILIARY CLEAN WATER SUPPLY SYSTEMS ARE NOT RECOMMENDED FOR USE WITH WATER HEATERS

AUXILIARY CLEAN WATER SUPPLY SYSTEM



Assembly with Bottled Water/Air Purge Function Assembly with 750 ml Bottle 50-120-00 Assembly with Large Capacity Bottle 50-220-00 750 ml Heavy Duty Bottle Only 50-060-00 Large Capacity Heavy Duty Bottle Only 50-260-00 Gasket 50-054-00

Assembly with Bottled Water/City Water Function

Assembly with 750 ml Bottle 50-125-00

Assembly with Large Capacity Bottle 50-225-00

- For antiseptic coolant or water
- Easy to install with installation instructions
- L bracket mounts to any flat vertical surface. Order 50-030-00 2" post mount kit for mounting to post
- On-off air supply toggle valve, with tubing
- Regulator, preset at factory
- Two position routing valve allowing use of bottled water/air purge or bottled water/city water
- · Water pressure gauge
- Heavy Duty Plastic Bottle is not autoclavable

IMPORTANT! AUXILIARY CLEAN WATER SUPPLY SYSTEMS ARE NOT RECOMMENDED FOR USE WITH WATER HEATERS

DUAL WATER BOTTLE SUPPLY SYSTEM



Post Mount 50-135-00
Wall Mount 50-137-00
Under the Counter Mount 50-139-00

- If your doctor is asking for more water capacity, but the assistant has a hard time holding and threading a large bottle full of water onto the existing clean water system, this dual system will solve both of their needs.
- This clean water system utilizes two 750 ml heavy duty HDPE bottles that operate independently of each other with toggle control valves. One bottle may be pressurized and in service while the other bottle is removed for cleaning or refilling, with pressure gauge reading bottle in use. This system also provides plumbing to facilitate use of an air purge, with the selector toggle to switch from bottled water to air purge or vice versa.
- · Asepsis Housing for easy cleaning
- Sterling color
- Includes mounting hardware, supply tubing and installation instructions

IMPORTANT! AUXILIARY CLEAN WATER SUPPLY SYSTEMS ARE NOT RECOMMENDED FOR USE WITH WATER HEATERS

HEAVY DUTY PLASTIC BOTTLE



750 ml Heavy Duty 50-060-00 Large Capacity Heavy 50-260-00 Duty Bottle (Approx. 2-Liter)

- For use with Auxiliary Clean Water Supply System
- Do not exceed 50 PSI air pressure
- High Density Polyethylene (HDPE) construction, translucent for checking fluid level
- Bottle material contains no latex
- Compatible with coolant liquids used in dentistry
- Heavy wall construction, safely operated at normal dental water system pressure
- Large capacity bottle includes pick-up tube
- This bottle is not autoclavable

UNIVERSAL POST MOUNT ASSEMBLY - SURF



- Mounts vacuum canisters 30-012-00 or 30-013-00 & auxiliary clean water system 50-100-00, 50-120-00, 50-200-00 & 50-220-00 on a 2" dia. post (2) 6-32 threaded mounting holes
- Includes mounting screws
- Surf finish

800-225-4621

50-030-00

KAVO MULTIFLEX COUPLING TYPE O-RING KIT

95-226-00



- Fits KAVO Multiflex Coupling
- Recommended to replace all (5) O-rings on coupler
- NEW KAVO cross ref. #1005.1033
 All O-rings are black
 OLD KAVO cross ref. # 9053111
 (2) small O-rings were white (3) larger O-rings were black

O-RING INSTALLATION TOOL KIT FOR HANDPIECE COUPLERS

20-095-00



- Includes two separate O-ring tools to assist when installing O-rings on KAVO MULTIFLEX couplers
- Use with Chapman-Huffman O-ring Kit 95-226-00
- · Stainless steel construction
- Fits Chapman-Huffman couplers: 20-075-00, 20-077-00 and 20-079-00

MIDWEST XGT/STYLUS TYPE 5-HOLE COUPLER O-RING KIT


95-253-00



- Fits MIDWEST XGT/Stylus 5-hole coupler
- MIDWEST cross ref. # 790255R
- · Pkg. of 2 sets

MIDWEST XGT/STYLUS TYPE 6-HOLE COUPLER O-RING KIT


95-255-00



- Fits MIDWEST XGT/Stylus 6-hole coupler
- MIDWEST cross ref. # 790256
- · Pkg. of 2 sets

SWIVEL O-RING KIT

95-227-00



- Fits STAR 430 SWL handpiece swivel
- Fits TITAN S/SW Sonic Scaler Swivels
- Autoclavable
- · Green color
- Fits all two, three and four line swivels
- Package contains two large and two small O-rings
- STAR cross ref. old # 59884 new # 259884

95-229-00



 Fits both NSK and SCHEIN master Q.D. handpiece swivel couplers

CAVITRON TYPE INSERT O-RING KIT

95-228-00



- Fits CAVITRON Jet nozzle, Slimline, P Style, CAVITRON Med and Pel inserts
- Autoclavable
- · Green color
- DENTSPLY cross ref. # 62605
- Pkg. of 12 in resealable bag

CAVITRON TYPE INSERT O-RING KIT

95-245-00



- Fits FSI, FSI-SLI AND TFI inserts
- Autoclavable
- · Black color
- DENTSPLY cross ref. # 62351
- Pkg. of 12

NOTE: Vibration from dental handpiece turbines is absorbed by O-rings or quad rings. These O-rings and quad rings eventually breakdown, making the turbine sound like it has worn or defective bearings. By replacing these O-rings or quad rings, the turbine may sound normal again, this may also extend the life of the bearing.

KAVO STYLE AIR OR ELECTRIC MOTOR O-RING


95-237-00



Fits the following KAVO handpieces

- 181H Air Motor requires (3)
- 181L Air Motor requires (3)
- KL700 Electric Motor requires (4)
- 196 Series Electric Motor requires (4)
- KAVO cross ref. # 0200.6120
- Pkg. of 1

HANDPIECE TURBINE REPLACEMENT O-RINGS



 MIDWEST Quiet-Air
 95-231-00

 STAR 430
 95-232-00

 KAVO
 95-233-00

 ATHENA/CHAMPION
 95-234-00

- · For use in dental handpiece turbines
- Autoclavable
- Pkg. of 2

HANDPIECE TURBINE REPLACEMENT QUAD RINGS




MIDWEST TRADITION

95-235-00

- For use in dental handpiece turbines
- Autoclavable
- Pkg. of 2

KAVO QUATTROCARE TYPE REPLACEMENT O-RINGS



O-RING KIT



95-140-00

- Consists of convenient 20 compartment tool box size case, more than 200 O-rings in 20 sizes from 3/32" OD to 1" OD
- · O-rings packaged in resealable bags
- Use only Parker O-Ring Lube. Chapman-Huffman item 20-211-00 sold on pg. 44

ACCESSORIES


CHAPMAN-HUFFMAN, INC.

NITRILE / BUNA N REPLACEMENT O-RINGS PACKAGED IN RESEALABLE BAGS FOR O-RING KIT 95-140-00



ID	OD	WALL	Pkg. Q	ty
1/32	3/32	1/32	20	007-01Q-00
3/64	9/64	3/64	20	007-02Q-00
1/16	11/64	1/16	20	007-03Q-00
5/64	13/64	1/16	20	007-04Q-00
3/32	7/32	1/16	20	007-05Q-00
1/8	1/4	1/16	20	007-06Q-00
5/32	9/32	1/16	20	007-07Q-00
3/16	5/16	1/16	20	007-08Q-00
7/32	11/32	1/16	10	007-09X-00
1/4	3/8	1/16	10	007-10X-00
5/16	7/16	1/16	10	007-11X-00
3/8	1/2	1/16	10	007-12X-00
7/16	9/16	1/16	10	007-13X-00
1/2	5/8	1/16	10	007-14X-00
9/16	11/16	1/16	5	007-15V-00
5/8	3/4	1/16	5	007-16V-00
11/16	13/16	1/16	5	007-17V-00
3/4	7/8	1/16	5	007-18V-00
13/16	15/16	1/16	5	007-19V-00
7/8	1	1/16	5	007-20V-00
1-1/16	1-3/16	1/16	5	007-23V-00
1-1/8	1-1/4	1/16	5	007-24V-00
1/4	5/16	1/32	10	007-55X-00
.362	.568	.103	10	007-61X-00

VITON REPLACEMENT O-RINGS



1/32	3/32	1/32	10	007-71X-00
3/64	9/64	3/64	10	007-72X-00
1/16	11/64	1/16	10	007-73X-00
5/64	13/64	1/16	10	007-74X-00
1/8	1/4	1/16	10	007-76X-00
5/32	9/32	1/16	10	007-77X-00
3/16	5/16	1/16	10	007-78X-00
1/4	3/8	1/16	10	007-80X-00
5/16	7/16	1/16	10	007-81X-00
3/8	1/2	1/16	10	007-82X-00
7/16	9/16	1/16	10	007-83X-00
1/2	5/8	1/16	10	007-84X-00
1/4	5/16	1/32	10	007-85X-00

- 75 durometer
- Brown color
- Use only Parker Super O-Lube. Chapman-Huffman item 20-212-00 sold on pg. 44
- 10 O-rings per package in resealable bag


DCI and ADEC O-ring cross reference guide

Note! The below part numbers are non-standard O-rings. Standard O-rings in

BUNA-N, VITON and ETHYLENE (EP) can be found on pgs. 33 and 34

DCI CROSS REF.	ADEC CROSS REF.	MATERIAL	FRACTIONAL DIMENSION	ID X CROSS SECTION	PKG QTY	
2245	035.003.00	BUNA-N	1/32	.029 X .040	10	009-01X-00
2255		BUNA-N	3/64	.053 X .064	10	009-03X-00
2252		VITON	3/64	.053 X .064	10	009-04X-00
2220		BUNA-N	1/16	.057 X .040	10	009-05X-00
2217		BUNA-N	5/64	.075 X .035	10	009-06X-00
2266		BUNA-N	5/64	.081 X .103	10	009-08X-00
2267		VITON	5/64	.081 X .103	10	009-09X-00
2218		BUNA-N	3/32	.085 X .055	10	009-10X-00
2262		VITON	1/8	.129 X .050	10	009-12X-00
2282		VITON	1/8	.136 X .040	10	009-13X-00
2271		BUNA-N	9/64	.144 X .025	10	009-14X-00
2208		BUNA-N	3/16	.172 X .037	10	009-15X-00
2209		BUNA-N	3/16	.176 X .050	10	009-16X-00
2212		BUNA-N	3/16	.180 X .040	10	40-019X-00
2247		BUNA-N	1/4	.248 X .048	10	009-20X-00
2283		VITON	1/4	.250 X .035	10	009-21X-00
2260		BUNA-N	17/64	.270 X .070	10	009-22X-00
2285		VITON	19/64	.300 X .050	10	009-23X-00
2284		VITON	7/16	.438 X .110	10	009-25X-00
2206		BUNA-N	7/16	.447 X .052	10	009-26X-00
2253		VITON	7/16	.447 X .052	10	009-27X-00
2224		BUNA-N	15/16	.924 X .103	10	009-31X-00
2270		BUNA-N	15/16	.926 X .070	10	009-32X-00
2221		BUNA-N	1-55/64	.186 X .070	10	009-33X-00
2246		BUNA-N	2-7/64	2.109 X .139	10	009-34X-00

ETHYLENE PROPYLENE (EP) REPLACEMENT O-RINGS



- (EP) Ethylene Propylene, Superior O-ring for water applications
- •ADEC cross reference numbers listed below
- •Use only Parker Super O-Lube. Item 20-212-00 sold on pg. 44
- Autoclavable

I.D.	O.D.	WALL	ADEC Cross Ref. #	Pkg. Qty	
.029	109	.040	034.001.01	10	037-01X-00
.042	142	.050	034.002.01	10	037-02X-00
.042	142	.050	Lavender	10	037-02X-01
.056	176	.060	034.003.01	10	037-03X-00
.114	254	.070	034.006.01	10	037-06X-00
.239	.379	.070	N/A	10	037-10X-00
.364	.504	.070	034.012.01	10	037-12X-00
.426	.566	.070	034.013.01	10	037-13X-00
.489	.629	.070	034.014.02	10	037-14X-00
.739	.879	.070	034.018.02	10	037-18X-00
.801	.941	.070	034.019.01	10	037-19X-00
.070	150	.040	035.027.01	10	037-82X-00
.101	.241	.070	035.033.01	10	037-86X-00
.176	.256	.040	035.017.01	10	037-90X-00
.206	.412	.103	034.107.01	10	037-94X-00

WATER HEATER MANUAL ON/OFF



110 volt 12-655-00 220 volt 12-690-00

- ADEC Style
- Adjustable thermostat for controlling water temperature
- Manual On/Off, 110 volt or 220 volt
- Reservoir type heater with quick recovery
- Continuous warm water to syringe & handpiece
- Compact 2" Dia. X 4" High
- · Finished in noncorrosive nickel plating
- Hospital grade cord with right angle plug on 110 volt model only.
- 220 volt model comes with cord less right angle plug

NOTE: Not recommended for use with self-contained water system due to potential damage

WATER HEATER WITH AIR ACTUATED ON/OFF

110 volt 220 volt 12-660-00

12-695-00



- ADEC Style
- Adjustable thermostat for controlling water temperature
- Air actuated On/Off, 110 volt or 220 volt
- Reservoir type heater with quick recovery
- Continuous warm water to syringe & handpiece
- Compact 2" Dia. X 4" High
- Finished in noncorrosive nickel plating
- Hospital grade cord with right angle plug on 110 volt model only.
- 220 volt model comes with cord less right angle plug

NOTE: Not recommended for use with self-contained water system due to potential damage

ROUND GAUGE



- 0-160 PSI
- 1 5/8" diameter face
- 1/8 MPT or 10-32 port
- For air or water service

CCESSORIES

20-105-00

20-108-00

- 0-15 PSI (for low pressure use only)
- Larger gauge face graduations for éasier reading at low pressures
- 1 5/8" diameter face
- 1/8 MPT back mount
- For air or water

PANEL MOUNT GAUGE



20-100-00

- 0-100 PSI
- 1 5/8" diameter face
- 1/8 MPT or 10-32 port
- · For air or water service

GAUGE REPLACEMENT LENS



20-103-00

- Fits Panel Mount Gauge 20-100-00
- Push in type
- Clear polyurethane construction

RECTANGULAR GAUGE



20-005-00

- 0-60 PSI scale
- Kg/cm² 1-4 scale
- 1/2" high X 1 3/4" wide
- Panel mount
- 10-32 port
- For air service

AIR COMPRESSOR MAIN GAUGE



20-109-00

- Designed for reading main tank pressure on air compressors
- 0-300 PSI range
- 1/4 MPT bottom mount
- 2 1/8" diameter face

20-130-00





- Fits wet type vacuum pumps
- Fits the following: AIR TECHNIQUES vacuum pumps, VACSTAR 2, 20, 4, 40; Single Pumps VACSTAR 50, 50H, 80, 80H dual pumps
- 1/8 MPT back mount
- 1 5/8" diameter face
- 0-30 inches of mercury
- AIR TECHNIQUES cross ref. # 55616

VACUUM GAUGE



20-106-00

- Fits wet type vacuum pumps
- · Commonly used on single pumps
- 1/4 MPT bottom mount
- 2" diameter face
- · 0-30 inches of mercury

LIQUID FILLED VACUUM GAUGE





- When your customer has a rough running wet type vacuum pump & the extra vibration keeps destroying vacuum gauges, this is the gauge you need.
- · Liquid filled
- Commonly used for single wet type vacuum pumps
- 1/4 MPT bottom mount
- 2-3/4" diameter face
- 0-30 Inches of mercury

VACUUM GAUGE ASSEMBLY





- Fits DENTAL-EZ MC-201 and MC-202 Dynamic-Dual evacuation System (some modification to mounting hole needed on older models)
- Cross ref. # 6-4555-003
- Fits VACSTAR 5, 8; dual vacuum pumps
- 0-30 inches of mercury
- Panel mount
- 1/4 MPT back mount

BRAKE CABLE ASSEMBLY

20-015-00



- Fits ADEC Priority 1005 Chair
- ADEC cross ref. # 61.0025.00
- Includes cable assembly and nut

CASSETTE CARRIAGE CABLE ASSEMBLY

20-256-00



- Fits S.S. White panorex cassette carriage
- · Includes (3) cable assembly screws

CASSETTE CARRIAGE CABLE SCREW

95-003-00



- Fits S.S. White Panorex cassette carriage cable
- Pkg. of 12

CAVITRON HANDPIECE SLEEVE ASSEMBLY

25-036-00



 Using this sleeve on a CAVITRON handpiece will allow it to hang in universal holders 25-015-00 or 25-016-00 and others

ARMREST STOP SCREW

20-012-00



- · Stainless steel construction
- Fits DENTAL-EZ JS and VS chairs
- DENTAL-EZ cross ref. # 3800-112

STORAGE RACKS WITH ORGANIZERS

60-400-00



- · Four storage tubs and racks with mounting brackets
- · Designed for mounting to cart, cabinet or wall
- · Perfect for custom mounting
- White plastic tubs with hinged lids
- Each tub comes with one vertical divider and one half-divider for depth adjustment
- Two position, horizontal and 30° tilt for convenient access
- Mounting surface is 6 1/2" X 21"

STORAGE TUB WITH DIVIDERS

60-405-00



- Fits single rack (60-402-00)
- Tub measures (outside) 5 1/4"L X 4 1/4"W X 5 3/4"H inside 4 3/4" deep
- Includes one half and one full divider
- White plastic tub, hinged cover and dividers

SINGLE RACK

60-402-00



- Fits storage tub (60-405-00) and single rack mounting bracket (60-403-00)
- Nonslip feet
- Color white
- 5 1/4" W X 6 1/2" D

SINGLE RACK MOUNTING BRACKET

60-403-00



- For mounting single rack (60-402-00) to cart, cabinet or wall
- · Color white

STORAGE TUB DIVIDERS



- Full (3 13/16 x 4 3/4) Half (2 1/4 x 3 7/8)
- *60-407-00 60-408-00*
- Fits 60-405-00 storage tub
- · White plastic

STORAGE TUB COVER / TUB HINGE



 White Cover
 60-410-00

 Clear Cover
 60-409-00

 White Hinge Only
 60-406-00

- Fits 60-405-00 storage tub
- Plastic construction
- Note: Order hinges (60-406-00) separately