

# Automatic Series IV Arm Mounted Delivery System



- Smooth lines for easy cleaning
- Automatic control for 3 handpieces
- Flex arm system mounts on 2" post
- Individual drive air & water coolant controls
- Stainless steel tray (PN 8013)
   & holder (PN 4280)
- 9 3/4" x 13 1/2" instrument tray holder

- Precision Comfort<sup>™</sup> Syringe
- Wet/Dry foot control
- Optional storage area for high-tech accessories
- Available in left or right hand configuration
- 7 ft. of 5/8" Umbilical
- Specify tubing color & type

DIMENSIONS PAGE 29

GRAY	WHITE	DESCRIPTION	PRICE
4477	4478	Arm Mounted Automatic Control for 3 HP	\$1,955.00
8956	8923	Top Post Mount Adapter	\$150.00
	8035	Touch-Up Paint	\$20.00
		•	

# Manual Series IV Arm Mounted Delivery System



- Smooth lines for easy cleaning
- Manual control for 2 handpieces
- $\bullet$  Precision Comfort  $^{\text{TM}}$  Syringe
- Stainless steel tray & holder or top mount
- Disk type foot control
- Plastic tray (4273) & holder (4274)
- 7 ft. of 5/8" Umbilical
- Specify tubing color & type

PAGE 29

GRAY	WHITE	DESCRIPTION	PRICE
4526	4529	Arm Mounted Manual Control for 2 HP (1 Wet & 1 Dry)	\$1,730.00
4528	4530	Arm Mounted Manual Control for 2 HP (2 Wet)	\$1,795.00
8956	8923	Top Post Mount Adapter	\$150.00



# Series IV Cabinet or Wall Mounted Automatic Delivery System



4470



- Smooth lines for easy cleaning
- Individual drive air & water coolant controls
- Precision Comfort<sup>™</sup> Syringe
- Optional storage area for high-tech accessories
- Wet/Dry foot control
- 6 ft. of 5/8" Umbilical
- Specify tubing color & type



PN	DESCRIPTION	PRICE
4470	Automatic Control for 3 HP	\$935.00
4278	Wall Mount Adapter	\$40.00
8035	Touch-Up Paint (White)	\$20.00

# Series IV Cabinet or Wall Mounted Manual Delivery System



- Smooth lines for easy cleaning
- Manual control for 2 handpieces
- Precision Comfort<sup>™</sup> Syringe
- Disk type foot control
- Multiple mounting options (refer to pages 93-95)
- 6 ft. of 1/2" Umbilical
- Specify tubing color & type



PN	DESCRIPTION	PRICE
4520	Manual Control for 2 HP (1 Wet & 1 Dry)	\$670.00
4521	Manual Control for 2 HP (2 Wet)	\$700.00
8035	Touch-Up Paint (White)	\$20.00
	·	



# Series IV Cart Systems

- Sturdy all-metal control head w/ smooth lines for easy cleaning
- Automatic control for 3 handpieces or manual control for 2 handpieces
- Individual drive air & water coolant controls
- SS tray & holder
- Precision Comfort<sup>™</sup> Syringe
- Wet/Dry foot control (Auto units only)
- Deluxe clean water system
- 1/4" Water Q.D. w/ flow control
- 7 ft. of 1 1/4" Umbilical
- Utility center w/ white box
- Specify tubing color & type

#### Vacuum Instrumentation Options

- Precision Comfort<sup>™</sup> Premium vacuum pkg includes two 5105 Premium HVE's & one 5215 Premium SE
- Precision Comfort<sup>™</sup> vacuum pkg includes two 5115 standard HVE's & one 5225 standard SE







PN	DESCRIPTION	PRICE
4511	Automatic Controls for 3 HP	\$2,450.00
4540	Automatic Controls for 3 HP w/ Premium Vacuum	\$2,745.00
4542	Automatic Controls for 3 HP w/ Standard Vacuum	\$2,720.00
4240	Manual Control for 2 HP w/ Premium Vacuum	\$2,050.00
4242	Manual Control for 2 HP w/ Standard Vacuum	\$2,025.00
4245	Manual Control for 2 HP	\$1,615.00
8035	Touch-Up Paint (White)	\$20.00

			Д	Arm Mounted			Cabinet Mounted			
Series IV Handpiece Controls		WHITE	4478	4529	4530	4470	4520	4521		
		GRAY	4477	4526	4528					
Standard Features		PRICE	\$1,955	\$1,730	\$1,795	\$935	\$670	\$700		
Automatic Controls for 3 Handpieces (expandable to 4)			•			•				
Manual Control for 2 Handpieces (1 wet & 1 dry)				•			•			
Manual Control for 2 Handpieces (2 wet)					•			•		
Wet/Dry Foot Control			•			•				
Economy Foot Control				•	•		•	•		
Precision Comfort™ Syringe	•	•	•	•	•	•				
Individual Non Retracting Water Coolant			•	•	•	•	•	•		
Handpiece Purge Valve			•			•				
50" Flex Arm for 2" Post w/ Air Brake			•	•						
Cabinet Mount						•	•			
Stainless Steel Tray			•							
Oil Collector for Handpiece Exhaust			•							
Options	P/N	Price								
Wall Mount Bracket	8647	\$150.00								
40" Chair Mount Post w/ P&C Bushing	8444	\$160.00	•	•						
White Top Post Mount Adapter	8923	\$150.00	•	•	•					
Gray Top Post Mount Adapter	8956	\$150.00	•	•	•					
Chair Adapters *Refer to Chair & Light Tab for com	nplete l	isting*								
Belmont Bel 20	8493	\$300.00								
A-dec® Priority	8515	\$300.00	•	•	•					
Pelton & Crane® 1500/Marus Empire	8918	\$300.00	•	•	•					
Light Source Systems *Refer to Handpiece Accessori	es Tab	for complete	listing*							
ISO 5 hole light source system w/ bulb, transformer & 1 straight gray tubing	8744	\$665.00	•							
ISO-C 6 pin light source system w/ transformer & 1straight gray tubing	8749	\$600.00	•			•				
Tubing (expandable to 3) 110V **Specify HP Position**										
Control Options										
Add Handpiece Position	4467	\$110.00	•			•				
Delete Handpiece Position		-\$54.40	•	•	•	•	•			
Delete Syringe, Tubing & Holder		-\$105.80	•	•		•	•			
Circulating Syringe & Water Heater	3575	\$449.00	•	•		•				
Add Chip Air or Scaler Button on Foot Control		\$45.00	•			•				
Stainless Steel Tray	8013	\$39.90		•	•					
Piezo Scaler Installed	8053	\$1,085.00	•			•				
Sub Chassis for Electric Motor Handpiece Control — (Internal Size 6"x1 3/4" x 2 1/2")	4230	\$45.00	•			•				
Accessories *Refer to DCI Catalog for complete	regulat	or & water sy	stem opt <u>ic</u>	ns*						
Utility Center w/ White Box & Regulators	8329	\$398.00	•			•	•			
2-Liter Asepsis Deluxe Water System-Cabinet Mount	8184	\$220.00				•	•			
2-Liter Asepsis Deluxe Water System w/ 2" Post Clamp	8928	\$315.00	•	•	•					
White Telescoping Arm w/ Premium Vacuum Valves	5476	\$725.00								



Series IV Handpiece Controls				utomatic C	arts	Manual Carts			
			4511	4540	4542	4240	4242	4245	
Standard Features			\$2,450	\$2,745	\$2,720	\$2,050	\$2,025	\$1,615	
Automatic Controls for 3 Handpieces (expandable to 4)			•	•	•				
Manual Control for 2 Handpieces (1 wet & 1 dry)				•	•	•			
Wet/Dry Foot Control	•	•	•						
Economy Foot Control				•	•	•			
Precision Comfort™ Syringe				•	•	•	•	•	
Individual Non Retracting Water Coolant				•	•				
Handpiece Purge Valve	•	•	•						
U-Frame				•	•	•	•	•	
Stainless Steel Tray				•	•	•	•	•	
Oil Collector for Handpiece Exhaust				•	•				
Deluxe Clean Water System				•	•	•	•	•	
Utility Center w/ White Box & Regulators			•	•	•	•	•	•	
1/4" Water Quick Disconnect			•	•	•	•	•	•	
Precision Comfort™ Standard HVE & SE Vacuum Package					•		•		
Precision Comfort™ Premium HVE & SE Vacuum Package				•		•			
Light Source Systems *Refer to Handpiece Accessori	ies Tab	for complete	e listing*	•				,	
ISO 5 hole light source system w/ bulb, transformer & 1 straight gray tubing	8744	\$665.00	•	•	•				
ISO-C 6 pin light source system w/ transformer & 1 straight gray tubing	8749	\$600.00							
Tubing (expandable to 3) 110V **Specify HP Position**			•	•	•				
Control Options									
Piezo Scaler Installed	8053	\$1,085.00	•	•	•	•	•	•	
Add Handpiece Position	4467	\$110.00	•	•	•				
Delete Handpiece Position		-\$54.40	•	•	•	•	•	•	
Delete Syringe, Tubing & Holder		-\$105.80	•	•	•	•	•	•	
Circulating Syringe w/ Water Heater	3575	\$449.00	•	•	•	•	•	•	
Add Chip Air or Scaler Button on Foot Control		\$45.00	•	•	•				
Delete Utility Center & Box		-\$310.00	•	•	•	•	•	•	
Delete Clean Water System		-\$90.00	•	•		•		•	



## Alternative Arm-Mounted Delivery System

- Automatic control for 3 handpieces or manual control for 2 handpieces
- Flex arm system mounts on 2" post
- Individual drive air & water coolant controls
- Instrument tray holder

- Precision Comfort<sup>™</sup> Syringe
- Wet/Dry foot control (4376 only)
- Dimensions: Horizontal reach
   62 3/4" Fully closed: 30 3/4"L
   x 21 1/4" W
- 7 ft. of 5/8" Umbilical
- Specify tubing color & type

PN	DESCRIPTION	PRICE
4376	Arm-Mounted Automatic Control for 3 Handpieces	\$1,795.00
4128	Arm-Mounted Manual Control 4120 for 1 Wet & 1 Dry Handpiece	\$1,535.00
4129	Arm-Mounted Manual Control 4121 for 2 Wet Handpieces	\$1,615.00
8923	Top Post Mount Adapter (White)	\$150.00
8035	Touch-Up Paint (White)	\$20.00



# Alternative Cabinet or Wall Mounted Automatic Delivery System

- Individual drive air & water coolant controls
- Multiple mounting options (refer to pages 93-95)
- Precision Comfort<sup>™</sup> Syringe
- Wet/Dry foot control
- Dimensions: 12" x 8" x 5"
- 6 ft. of 1/2" Umbilical
- Specify tubing color & type



PN DESCRIPTION

4370 Alternative Automatic Control for 3 Handpieces

PRICE \$775.00



1120

# Alternative Cabinet or Wall Mounted Manual Delivery System

- Manual control for 2 handpieces
- Precision Comfort<sup>™</sup> Syringe
- Disk type foot control
- Multiple mounting options (refer to pages 93-95)
- Dimensions: 10" x 3" x 4"
- 6 ft. of 1/2" Umbilical
- Specify tubing color & type



PN	DESCRIPTION	PRICE
	Manual Control for 1 Wet & 1 Dry Handpiece Manual Control for 2 Wet Handpieces	\$525.00 \$570.00



# Alternative Cart Systems

- Automatic control for 3 handpieces or manual control for 2 handpieces
- Individual drive air & water coolant controls
- ullet Precision Comfort Syringe
- Wet/Dry foot control (Auto units only)
- 7 ft. 3/4" Umbilical
- Specify tubing color & type

#### Vacuum Instrumentation Options

Premium Package features Precision Comfort™ autoclavable premium vacuum instruments:

• Two 5105 Premium HVE's & one 5215 Premium SE

Standard Package features autoclavable Precision Comfort™Standard vacuum instruments:

Two 5115 standard HVE's & one 5225 standard SE



PN	DESCRIPTION	PRICE
4411	Automatic Control for 3 Handpieces	\$2,145.00
4440	Automatic Control for 3 Handpieces w/ Premium Vacuum Pkg	\$2,435.00
4442	Automatic Control for 3 Handpieces w/ Standard Vacuum Pkg	\$2,400.00
4145	Manual Control for 2 Handpieces	\$1,305.00
4140	Manual Control for 2 Handpieces w/ Premium Vacuum Pkg	\$1,735.00
4142	Manual Control for 2 Handpieces w/ Standard Vacuum Pkg	\$1,700.00



# **Operatory Support Cart**

- Solid 1" white laminate top w/ gray edge banding
- 19" x 19" work surface
- Offered in U-frame or H-frame configurations
- Room to mount assistants instruments & water systems
- Space to mount 4416 unit w/ 4445 holder bar & 5417 assistant's package

PN	DESCRIPTION	PRICE
4227	U-Frame Support Cart	\$865.00
4228	H-Frame Support Cart	\$870.00

DCI Alternative Controls			,	Arm Mount	ed	Cabinet Mounted			Specialty Units			
			4376	4128	4129	4370	4120	4121	4502	4102	4446	4447
Standard Features	Standard Features			\$1,535	\$1,615	\$775	\$525	\$570	\$510	\$440	\$2,645	\$2,435
Automatic Controls for 3 Handpieces (expandable to 4)			•			•						•
Manual Control for 2 Handpieces (1 wet & 1 dry)				•			•		•	•		
Manual Control for 2 Handpieces (2 wet)					•							
Wet/Dry Foot Control			•			•					•	•
Economy Foot Control												
Precision Comfort™ Syringe												
Individual Non Retracting Water Coolant			•	•	•	•		•				•
50" Flex Arm for 2" post w/ Air Brake			•	•	•							
Cabinet Mount						•	•	•	•	•	•	•
Options	P/N	Price										
Wall Mount Bracket	8647	\$150.00	•	•	•							
40" Chair mount post w/ P&C bushing	8444	\$160.00	•	•	•							
White Top Post Mount Adapter	8923	\$150.00			•							
Chair Adapters *Refer to Chair & Light Tab for com	plete li	sting*										
Belmont Bel 20	8493	\$300.00	•	•	•							
A-dec® Priority	8515	\$300.00	•	•	•							
Pelton & Crane® 1500/Marus Empire	8918	\$300.00	•	•	•							
Light Source Systems *Refer to Handpiece Accessorie	es Tab f	or complete	listing*									
ISO 5 hole light source system w/ bulb, transformer & 1 straight gray tubing	8744	\$665.00	•			•					•	•
ISO-C 6 pin light source system w/ transformer & 1 straight gray tubing	8749	\$600.00	•			•					•	•
Tubing (expandable to 3) 110V **Specify HP Position**												
Control Options												
Piezo Scaler Installed	8053	\$1,085.00	•									
Delete Handpiece Position		-\$54.40	•	•	•	•	•	•	•	•	•	•
Delete Syringe, Tubing & Holder		-\$105.80	•	•	•	•	•	•	•	•	•	•
Circulating Syringe w/ Water Heater	3575	\$449.00	•	•	•	•	•	•	•	•	•	•
Add Chip Air or Scaler Button on Foot Control		\$45.00	•			•					•	•
Accessories *Refer to DCI Catalog for complete re	egulato	or and water	system op	tions*								
Utility Center w/ Gray Box & regulators	8321	\$375.00	•	•	•	•	•	•	•	•		
2-Liter Standard Water System - Cabinet Mount	8144	\$155.00	•	•	•	•	•	•	•	•	•	•
2-Liter Deluxe Water System - Cabinet Mount	8143	\$190.00	•	•	•	•		•	•	•	•	•
750ml Standard Water System - Cabinet Mount	8145	\$150.00	•	•	•	•	•	•	•	•	•	•
750ml Deluxe Water System - Cabinet Mount	8139	\$185.00	•	•	•	•	•	•	•	•	•	•
Adapter for 2" Post, White Powder Coat	8162	\$60.00	•	•	•							



DCI Alternative Controls		Automatic Carts		Manual Carts				
			4411	4440	4442	4140	4142	4145
Standard Features		\$2,145	\$2,435	\$2,400	\$1,735	\$1,700	\$1,305	
Automatic Controls for 3 Handpieces (expandable to 4)			•	•	•			
Manual Control for 2 Handpieces (1 wet & 1 dry)						•	•	•
Wet/Dry Foot Control			•	•	•			
Economy Foot Control							•	
Precision Comfort™ Syringe							•	•
Individual Non Retracting Water Coolant			•	•	•		•	•
H-Frame			•	•	•		•	•
Stainless Steel Tray		•	•	•				
Precision Comfort™ Standard HVE & SE Vacuum Package				•				
Precision Comfort™ Premium HVE & SE Vacuum Package					•			
Light Source Systems *Refer to Handpiece Accessorie	es Tab f	or complete	listing*					
ISO 5 hole light source system w/ bulb, transformer & 1 straight gray tubing	8744	\$665.00	•	•	•			
ISO-C 6 pin light source system w/ transformer & 1 straight gray tubing	8749	\$600.00	•	•	•			
Tubing (expandable to 3) 110V **Specify HP Position**								
Control Options								
Piezo Scaler Installed	8053	\$1,085.00	•		•			
Add Handpiece Position	4467	\$110.00	•	•	•			
Delete Handpiece Position		-\$54.40	•	•	•		•	•
Delete Syringe, Tubing & Holder		-\$105.80	•	•	•	•	•	•
Circulating Syringe w/ Water Heater	3575	\$449.00	•	•	•	•	•	•
Add Chip Air or Scaler Button on Foot Control		\$45.00	•	•	•			
Accessories *Refer to DCI Catalog for complete r	egulato	or and water	system opt	tions*				
Utility Center w/ Gray Box & regulators	8321	\$375.00	•		•	•	•	•
Water System Integrated into Carts	4448	\$90.00			•		•	•







# Panel & Horizontal Surface Mount Automatic Unit for 3 Handpieces

- Individual drive air & water coolant controls
- Precision Comfort<sup>™</sup> Syringe
- Wet/Dry foot control
- Control dimensions: 8 1/2" x 1 3/4" x 5"
- Cut-out dimensions: 7 1/2" x 1 1/2"
- Horizontal system can be mounted to top or bottom of surface
- Specify tubing color & type



PN	DESCRIPTION	PRICE
4415	Panel-Mount Automatic Control	\$845.00
4416	Horizontal-Mount Automatic Control	\$1,010.00
4445	Optional Swing Arm Holder Bar	\$140.00







# Panel & Horizontal Surface Mount Manual Unit for 2 Handpieces

- Individual drive air & water coolant controls
- Precision Comfort<sup>™</sup> Syringe
- Wet/Dry foot control
- Control dimensions: 8 1/2" x 1 3/4" x 5"
- Cut-out dimensions: 7 1/2" x 1 1/2"
- Horizontal system can be mounted to top or bottom of surface
- Specify tubing color & type



PN	DESCRIPTION	PRICE
4400	Panel-Mount Manual Control	\$545.00
4401	Horizontal-Mount Manual Control	\$610.00
4445	Optional Swing Arm Holder Bar	\$140.00



#### DCI III Manual Unit

- Features smooth lines for easy cleaning
- Two HP manual selection
- Precision Comfort<sup>™</sup> Syringe
- Compact design mounts to any vertical surface or folding ridged arms
- Pneumatic controls eliminate most moving parts for better reliability
- 6 ft. 1/2" Umbilical
- Specify tubing color & type



DIMENSIONS PAGE 30

PN DESCRIPTION
4502 DCI III Manual Unit

\$510.00



#### DCI II Manual Unit

- Two HP manual selection
- Precision Comfort<sup>™</sup> Syringe
- Compact design mounts to any vertical surface or folding rigid arms
- Dimensions: 6" x 4 1/4" x 2 1/4"
- 6 ft. 1/2" Umbilical
- Specify tubing color & type

PN	DESCRIPTION	PRICE
4102	DCI II Manual Unit	\$440.00



# Air & Water Auxiliary QD Panel

- 1/4" Water guick connect w/ flow control
- Dimensions: 3" x 4"
- 3/8" Air quick connect w/ on/off switch
- Available in white or gray

PN	DESCRIPTION	PRICE
6559	Gray Auxiliary Panel	\$150.00
6569	White Auxiliary Panel	\$150.00



### Pelton & Crane Executive FEATURED ITEM Conversion Kit

- Precision Comfort™ Premium Autoclavable Vacuum Valves include two 5105 HVEs & one 5215 SE
- Automatic control for three handpieces
- Individual flow control for water coolant to each handpiece
- Individual handpiece pressure adjustments
- Includes 8323 utility center
- Automatic handpiece holders w/ manual lockout toggle
- Variable speed disc-type foot control w/ wet/dry toggle
- Precision Comfort<sup>™</sup> three-way syringe
- Dimensions: 12 1/2" x 18" x 4 1/2"
- 6 ft. Umbilical
- Specify tubing color & type

PN	DESCRIPTION	PRICE
4446	Executive Conversion Kit	\$2,645.00
4447	Executive Conversion Kit Less Utility Center	\$2,435.00



## Laboratory Handpiece Control Kit

- Requires air only
- No water connections
- Four-hole handpiece tubing 432T
- Handpiece holder & bar 5910 & 5929
- Economy foot control 6301
- Air regulator w/ gauge 7220 & 7260

PN	DESCRIPTION	PRICE
4198	Laboratory Handpiece Control Kit	\$168.00



## **Economy Self-Contained Water Systems**

- Deluxe models have a selector toggle allowing you to switch between city & bottled water
- 750 ml & 2 liter heavy duty bottles
- Standard 28 mm thread for water bottle
- Gauge 7260 shows bottle pressure
- Bracket mounts to any flat vertical surface

- Optional post-mount available
- On/off toggle releases bottle pressure
- Dimensions for 8144 & 8143: 3" x 5" x 16"
- Dimensions for 8145 & 8139: 3" x 5" x 12"







	<u>( 1937) </u>	
PN	DESCRIPTION	PRICE
8144	2-Liter Standard System	\$155.00
8143	2-Liter Deluxe System	\$190.00
8145	750-ml Standard System (includes extra bottle)	\$150.00
8139	750-ml Deluxe System (includes extra bottle)	\$185.00
8137	Adapter to Mount Water System to 2" Post, Anodized Aluminum	\$55.00
8162	Adapter to Mount Water System to 2" Post, White Powder Coat	\$60.00
8255	2-3/8" Post, White Powder Coat	\$57.00

# Asepsis Self-Contained Water Systems

- 8178
- Deluxe models have a selector toggle allowing you to switch between city & bottled water
- Models for 750-ml & 2-liter bottles
- Standard 28 mm thread for water bottle
- Pressure gauge 7265
- Mounts to vertical surfaces
- Optional post-mount available
- On/off toggle releases
- Air purge toggle bottle pressure

- Metal cover w/ white powder coat finish
- Dimensions for 8177 & 8178: 8 1/4" x 4 1/2" x 13"
- Dimensions for 8179 & 8180: 8 1/4" x 4 1/2" x 17"
- Dimensions for 8181 & 8182:
   5" x 4 1/2" x 13"
- Dimensions for 8183 & 8184:
   5" x 4 1/2" x 17"







PN	DESCRIPTION	PRICE
8177	Dual 750-ml Asepsis Standard System	\$245.00
8178	Dual 750-ml Asepsis Deluxe System	\$280.00
8179	Dual 2-liter Asepsis Standard System	\$250.00
8180	Dual 2-liter Asepsis Deluxe System	\$285.00
8181	Single 750-ml Asepsis Standard System (includes extra bottle)	\$180.00
8182	Single 750-ml Asepsis Deluxe System (includes extra bottle)	\$215.00
8183	Single 2-liter Asepsis Standard System	\$185.00
8184	Single 2-liter Asepsis Deluxe System	\$220.00
8137	Adapter to Mount Water System to 2" Post, Anodized Aluminum	\$55.00
8162	Adapter to Mount Water System to 2" Post, White Powder Coat	\$60.00
8255	Adapter to Mount Water System to 2 3/8" Post, White Powder Coa	t \$57.00



# DCI Quick-Switch Water Bottle System

- Easy installation requires no tools
- Saves time in changing bottles
- Self-locking feature ensures proper positioning when pressurized
- Automatic pressure release when removing bottle
- Patented design exclusive to DCI International

DESCRIPTION	PRICE
Kit (contains converter & adapter)	\$70.00
Converter, Manifold QS w/ O-Ring	\$46.00
Adapter, Bottle QS w/ Washer	\$32.00
Adapter, Bottle QS (pkg of 5)	\$125.00
O-Ring for Converter Manifold (pkg of 12)	\$2.80
	Kit (contains converter & adapter) Converter, Manifold QS w/ O-Ring Adapter, Bottle QS w/ Washer Adapter, Bottle QS (pkg of 5)



# Heavy-Duty Bottles for Dental Unit Water Systems

#### Standard Bottles

- Heavy wall construction, safely operated at normal water system pressure
- High Density Polyethylene (HDPE) is compatible w/ coolant liquids used in dentistry
- Can be used on most dental systems (w/ the Universal Manifold 8134 or 8119)
- Bottle dimensions: (750-ml, 3" W x 9 1/2" H),
   (1-liter, 3 1/2"W x 11" H), (2-liter, 4" W x 13 7/8" H)
- Standard 28 mm threads

PN	DESCRIPTION	PRICE
8128	750 ml Plastic Bottle w/ Cap & Pick-Up Tube	\$6.45
8669	1-Liter Plastic Bottle w/ Cap & Pick-Up Tube	\$6.75
8164	2-Liter Plastic Bottle w/ Cap & Pick-Up Tube	\$11.25



# BioFree™ Bottles for Dental Unit Water Systems

- BioFree bottles are antimicrobial which inhibits the growth of microorganisms
- One year service life

PN	DESCRIPTION	PRICE
8109	750 ml BioFree Bottle	\$26.00



### **Dental Waterline Control Tablets**

• Assists in controlling dental waterline contaminants

PN	DESCRIPTION	PRICE
6824	Citrisil Complete Shock Dental Waterline Cleaner Tablets (10 tablets)	\$32.00
6825	Citrisil Blue Dental Waterline Cleaner – 1 Liter (20 tablets)	\$40.00
6826	Citrisil Blue Dental Waterline Cleaner – 2 Liter (20 tablets)	\$32.00





# Water System Component Parts

PN	DESCRIPTION	PRICE)
0165	Barb, 1/16" for Water System Manifold	\$.80
0166	Barb, 1/8" for Water System Manifold	\$.80
8134	Universal Manifold w/ Barbs (fits most dental water systems)	\$30.00
8136	Manifold Washer	\$3.00
8146	Includes DCI 750-ml Bottle 8128 & Universal Manifold w/ Barbs 8119	\$35.00
8147	Includes DCI 2-Liter Bottle 8164 & Universal Manifold w/ Barbs 8119	\$36.50
8119	Universal Manifold, Long Self Centering Neck, w/ Barbs	\$25.50
8037	Pick-up Tube for 750 ml Bottle	\$2.90
8038	Pick-up Tube for 2 L Bottle	\$3.00



# **Chair Adapters**

PN	DESCRIPTION	PRICE
8522	A-dec® Cascade	\$300.00
8515	A-dec® Priority	\$300.00
8493*	Belmont Bel 20 Westar 2001, Westar 5000 LSM DC-1000, LSM C0H0	\$300.00
8561	Marus Pro (RH) 1300-1600 DCI Equipment (RH)	\$300.00
8918	Pelton & Crane® Spirit 1500 / 2005 Marus Empire	\$300.00
_	wder coated finish uplete list of chair adapters see pg 89	



# **Light Posts**

• 2" Diameter posts

METAL	WHITE	GRAY	LENGTH	PRICE
8478	8912	8939	24" Post	\$136.00
8218	8421	8940	30" Post	\$145.00
8454	8913	8945	40" Post	\$190.00
8453	8449	8949	48" Post	\$210.00
8450	8914	8952	60" Post	\$219.00
8475	8474	8408	Post Cap	\$43.00
PN	DESCRIPT	ION		PRICE
8476	Crush Bush	ning for Base	of Post	\$38.00
8470	P&C Bushing			\$40.00
8471	Ritter Bush	ing		\$22.00
8472	Power Cor	d (8')		\$16.00



#### **Built-In Piezo Scaler**

- Integrates into dental unit
- Broad range of power
- Small control unit
- Frequency 33.5 KHZ
- Dimensions: 1 3/8" x 2 1/2" x 1 1/8"
- Water pressure 15 55 PSI
- Includes 3 tips & holder, metal wrench & sterilization box


PN	DESCRIPTION	PRICE
8053	Built-In Piezo Scaler	\$1,085.00
8017	Universal #1 Tip	\$140.00
8018	Universal #2 Slim Tip	\$140.00
8019	#3 Paddle Tip	\$140.00
	•	





# Cabinet Mounted Assistant's Instrumentation


- Mount in vertical or horizontal position
- Solids collector 5860 uses removable screen
- HVE & SE include 6' of tubing
- Syringe includes 7' of tubing
- Rotational 3 or 4 position holder
- Syringe flow control block 3021
- Premium pkg includes Precision Comfort<sup>™</sup> Syringe 3404, 5105 Premium HVE, 5215 Premium SE
- Standard package includes 3404 Precision Comfort™ Syringe, 5115 Standard HVE, 5225 Standard SE


PN	DESCRIPTION	PRICE
5416	Premium 3 Position (Syr, HVE, SE), 11 3/4" x 4" x 5 1/2"	\$608.00
5417	Premium 4 Position (Syr, HVE, SE, open position), $13 3/4$ " x 4" x 5 $1/2$ "	\$628.00
5418	Standard 3 Position (Syr, HVE, SE), 11 3/4" x 4" x 5"	\$575.00
5419	Standard 4 Position (Syr, HVE, SE, open position), 13 $3/4$ " x $4$ " x $5 1/2$ "	\$595.00














































Panel & Horizontal Surface Mount Control 4415 & 4416



Panel & Horizontal Surface Mount Control 4400 & 4401







# Asepsis Self-Contained Water System 8178

