

HANDPIECE CONTROLS

PM 500 DOCTOR'S HANDPIECE CONTROL SYSTEM

SYMBOL	PM 500	IN	CM
Α	•	5–5/8	14
В	•	9–5/8	24
С	•	2-5/8	7
D	•	7–3/4"	20
Е	•	2–1/2	6



STANDARD FEATURES INCLUDE

- 3-way air/water syringe
- · Automatic control for three handpieces
- · Master on/off with indicator
- · Individual handpiece pressure adjustment
- Individual water coolant adjustment
- · Non-retraction water coolant on/off control
- · Air Coolant Adjustment
- · Panel mount holder bar assembly
- Handpiece tubing
 4-hole Midwest, asepsis style, straight sterling
- · Disc type foot control
- Umbilical, 12'
- Remote mount holder assembly with control tubing, 2'



PRICE 1,295.00

В

DESCRIPTION

HANDPIECE CONTROL SYSTEM

PART NO. 2001-011

PRICE 360.00 65.00

DESCRIPTION(FACTORY OPTIONS)
STANDARD UTILITY CENTER
WET/DRY TYPE FOOT CONTROL

PART NO. 1115-001 1105-174-0

PM 200 DOCTOR'S HANDPIECE CONTROL SYSTEM


SYMBOL	PM 200	IN	CM
F	•	5–1/4	13
G	•	6–1/2	16
Н	•	2	5
	•	6–9/16	17
J	•	2-1/2	6



STANDARD FEATURES INCLUDE

- 3-way air/water syringe
- Manual control for two handpieces one wet/one dry
- · Master on/off with indicator
- · Individual handpiece pressure adjustment
- · Handpiece water coolant adjustment
- Non-retraction water coolant on/off control
- · Panel mount holder bar assembly
- Handpiece tubing
 4-hole Midwest, asepsis style, straight sterling
- · Disc type foot control
- · Umbilical, 12'
- · Remote mount holder assembly





G -

PRICE DESCRIPTION
780.00 HANDPIECE CONTROL SYSTEM

PRICE DESCRIPTION(FACTORY OPTIONS)
360.00 STANDARD UTILITY CENTER
35.00 AIR COOLANT ADJUSTMENT

PART NO. 2000-061 PART NO. 1115-001 1100-223

HANDPIECE CONTROLS **F**



AIR AND WATER ACCESSORY PANEL



STANDARD FEATURES INCLUDE

- · 3/8" female air quick disconnect
- · Toggled shut off for air quick disconnect
- 1/4" female water quick disconnect
- · Flow adjustment for water quick disconnect
- 4" x 3" flush mount panel
- · White powder coated finish

PRICE DESCRIPTION PART NO. AIR AND WATER ACCESSORY PANEL 1300-030 145.00

LT 100 LAB HANDPIECE CONTROL

STANDARD FEATURES INCLUDE

- · Manual control for one handpiece one dry
- · Master on/off with indicator
- · Handpiece pressure adjustment
- Handpiece tubing 2-hole Midwest, straight sterling
- Disc type foot control
- · Umbilical, 12'
- · Air pressure regulator

PRICE DESCRIPTION 450.00 HANDPIECE CONTROL PART NO. 2000-011

LT 10 LAB HANDPIECE CONTROL KIT







INCLUDES

· Control for One Handpiece



STANDARD FEATURES

- · Air pressure regulator
- Handpiece tubing 2-hole Midwest, straight sterling
- · Handpiece holder
- · Mounting bracket
- · Disc type foot control

PRICE DESCRIPTION 250.00 HANDPIECE CONTROL KIT PART NO. 2000-010

5