

November 11, 2015

Re: Bravo Autoclave

Dear valued SciCan customer,

Thank you for taking the time to contact our company with your request for information regarding the BRAVO autoclave.

If the BRAVO is inactive for 30 minutes, the unit enters Stand-by Off mode, in which the heating elements are turned off, avoiding the useless consumption of electricity. The machine will be activated when the door is opened or a control panel button is pressed. The time after which the heating elements are turned off can be adjusted between 0 and 300 minutes (in 30 minute increments) in the Stand-by Options on the Advanced Menu.

When the sterilizer door is not opened at the end of the cycle (e.g. if a cycle completes when the office is closed), the vacuum pump and fan will run automatically for 5 mins in 30 minute intervals, until the chamber temperature drops below 50°C, to prevent condensation from accumulating in the sterilization chamber as it cools.

I trust this addresses your questions regarding the BRAVO. If you require anything further please do not hesitate to contact us.

Thank you,

Alana Ross
Product Manager, Sterilizers