
## Progeny®—the leader in intraoral x-ray systems in the U.S. and the preferred choice of dentists around the world.



Minimum Ceiling Requirement: 86.5"

Cabinet Pass Through 66" Arm: 31" ⊠ x 47.5" H 76" Arm: 41" ⊠ x 47.5" H

Cabinet Pass Through (Alternate) 66" Arm:26.75"⊠ x 47.5"H 76" Arm:36.75"⊠ x 47.5"H

Base: 26.9" ⋈ x 36" L Heiaht: 76.9" Arm reach: 44"

Line Voltage	110-230VAC +/- 10% 50-60 Hz	
Tube Potential	60, 65 or 70 kVp	
Tube Current	4, 5, 6, or 7 mA	
Exposure Time	20 ms to 2 seconds	
Cones Included	70 mm aperture with 20 cm SID 35 x 45 mm aperture with 20 cm SID	
Minimum Inherent Filtration	2 mm Al equivalent @ 70 kVp	
Focal Spot	0.4 mm (IEC 336)	
Duty Factor	1:15	
Shipping ⊠ eight	95 lbs.	
Applicable Documents Federal Regulations	21 CFR, UL, CE 0120	
⊠ arranty	2 Years, Limited	
☑ p tio ns Exposure hand switch 35 x 45 mm Long Cone	Remote Exposure Stations	
Two stud mounting plate	0 1	
Choice of 56", 66", 76", or 82	2" (total reach) and Mobile	
Sensor		0.1
Film Size equivalent	Size 1	Size 2
Dimensions	37 x 24 x 5.5 mm	43 x 32 x 5.5 mm
Active Area	30 x 20 mm	36 x 26 mm
Sensor Cable	90 cm	90 cm
Computer Interface	High Speed USB 2.0	High Speed USB 2.0

## Computer Requirements

Please see www.progenydental.com for operating system and hardware requirements

Contact your Midmark representative for additional system requirement information.

## midmark.com

675 Heathrow Drive Lincolnshire, II 60069 847-415-9800 Toll-free 888-924-3800

Midmark 2013 All Rights Reserved Printed in the USA ML-00041 Rev. D

