

DCI Edge

Back by popular demand, DCI is proud to launch our new expanded equipment line, DCI Edge.

The new line encompasses new and existing operatory equipment products such as delivery systems & compressors. Look for more new product launches in 2013 & 2014 as DCI innovates to give you the edge!



DCI Edge Operatory Equipment



15 – 17

Operatory Light



18 – 26

Series 4
Delivery Systems



27 – 28

Alternative
Delivery Systems



29 – 30

Specialty
Units



31 – 32

Compressors



33

Accessories



34

Surgical Vacuums



35 – 44

Dimensions & Illustrated
Parts Breakdown

DCI Catalog App

Looking for a digital version of the DCI Catalog

Visit the iTunes or Android app store and download it today.

View the app on your mobile device, iPAD or laptop!



ERGONOMIC
INTROD CING THE NEW
LED TECHNOLOGY
SERIES 5
ENERGY EFFICIENT
LED OPERATORY LIGHT
DEFINED LIGHT PATTERN



COMING SOON!

Post Mount Light
LP5000



Track Mount Light
LT5000



System Mount Light
LS5000



COMING
SOON!

System Mount Light
shown w/ optional
swing mount light pole
DS4002



Series 5 Order Guide

COMING SOON!

WHITE	GRAY	DESCRIPTION	PRICE
POST MOUNT SERIES 5 LIGHT			
LP5000	LP5050	Post Mount Series 5 Light, 120V	1, 0.00
LP5100	LP5150	Post Mount Series 5 Light, 220V	1, 0.00
21		Light Post Extension Kit	1 5.00
SYSTEM MOUNT SERIES 5 LIGHT			
LS5000	LS5050	System Mount Series 5 Light, 120V	1, 0.00
LS5100	LS5150	System Mount Series 5 Light, 220V	1, 0.00
TRACK MOUNT SERIES 5 LIGHT			
LT5000	LT5050	Track Mount Series 5 Light, 120V	2, 00.00
LT5100	LT5150	Track Mount Series 5 Light, 220V	2, 00.00
LT000	LT050	Track Light Suspension Tube for .5 Ceiling	NC
LT005	LT055	Track Light Suspension Tube for .5 Ceiling	NC
LT000	LT050	Track Light Suspension Tube for .5 Ceiling	NC
LT005	LT055	Track Light Suspension Tube for .5 Ceiling	NC
LT0010	LT0510	Track Light Suspension Tube for 10' Ceiling	NC

Specifications

LED Lamp 24Volt, 10 Watts

Input Voltage 110V/220V, 0/50 Hz

Color temperature 5000 Kelvin

Light Intensity Lux 30,000 24,000 15,000

Focal Length 2.5

Illumination Area 3 x

Series 4 Delivery Systems



(Delivery System shown with
A-dec® chair & light)

UPGRADE
your delivery system
WITHOUT
THE EXPENSE
of a new chair!

Compatible with New & Existing Chairs!

- A-dec® Cascade 1040
- A-dec® Priority 1005
- Belmont® -Calibur B-20, B-50
- Belmont® Acutrac B-, S-
- DCI Equipment 1235
- Healthco® Celebrity
- Marus® DC14 0, DC1 0, DC1 00, DC1 02
- Pelton & Crane® 1500, 1 00, 2000, 3000



(Delivery System shown with
Belmont® chair & light)

DCI's NEW Series 4 Delivery System
**offers a Left/Right &
Post Mount Configuration**


Series 4 Swing Mount Delivery System



Available w/
optional
monitor mount



New larger head
for high tech
integration
12 1/2" x 14"



(Delivery System shown with
Belmont® chair & light)

DIMENSIONS
PAGE 3

GRAY | WHITE

DS4050	DS4000	Swing Mount Pole w/ Auto Dental Unit & Light Pole	\$4,900.00
DS4051	DS4001	Swing Mount Pole w/ Auto Dental Unit	\$4,400.00
DS4052	DS4002	Swing Mount Light Pole Only	\$1,200.00

ORDER GUIDE ON PAGE 21

FEATURED NEW ITEM

Series 4 Post Mount Delivery System

FEATURED NEW ITEM



DIMENSIONS
PAGE 3

GRAY	WHITE
DO4050	DO4000
DO4051	DO4001
DO4052	DO4002
DW4050	DW4000

SIDE DELIVERY SYSTEM
 Fits into a 32" or larger opening
 Includes
 Utility Center & Stainless Steel Box
 Mounting Hardware for cabinet mount
 Optional wall board available



Over-the-Patient Auto Dental Unit w/ PMU/Asst. Pkg	\$4,200.00
Over-the-Patient Auto Dental Unit	\$3,500.00
Post Mount Utilities (PMU)	\$1,800.00
Side Delivery Auto Dental Unit	\$3,200.00

ORDER GUIDE ON PAGE 21

Series 4 Delivery System

	Over the Patient Auto Dental Unit with PMU & Assistant's Package	Over the Patient Auto Dental Unit	Post Mount Utilities (PMU) with Assistant's Package	Auto Dental Unit and light Pole	Swing Mount Pole with Auto Dental Unit	Swing Mount Pole with Auto Dental Unit and light Pole	DCI Swing Mount Light Pole Only	DCI Side Delivery Auto Dental Unit with Mount
White	DO4000	DO4001	DO4002	DS4000	DS4001	DS4002	DW4000	
Gray	DO4050	DO4051	DO4052	DS4050	DS4051	DS4052	DW4050	
Retail	\$4,200	\$3,500	\$1,800	\$4,900	\$4,400	\$1,200	\$3,200	

CHAIR ADAPTERS (For Traditional Chair Mounted Equipment)

A-dec® Cascade	8522	•	•	•				
A-dec® Priority	8515	•	•	•				
Belmont Bel 20/Healthco Celebrity	8493	•	•	•				
Belmont AcuTrac	8917	•	•	•				
Pelton Spirit/Schein Empire	8918	•	•	•				
Marus Pro/DCI Equipment	8561	•	•	•				

CHAIR ADAPTERS (For Swing Mount Equipment)

Pelton 1500, 1800, 2000, 3000, A-dec® Cascade 1040	4489				•	•	•	
Marus DC1490, DC1690, DC1700, DC1702, DCI Equipment 1235	4491				•	•	•	
Belmont X-Calibur B-20, B-50, Healthco Celebrity	4492				•	•	•	
AcuTrac B-99, S-99	4501				•	•	•	
A-dec® Priority Chair Adapter 1005	4498				•	•	•	

STANDARD FEATURES

Larger control head to accommodate accessories (Scaler, light source, etc.)		•	•		•	•		•
50" Flex arm		•	•		•	•		
3 Handpiece automatic delivery		•	•		•	•		•
Precision Comfort Air/Water syringe		•	•		•	•		•
Wet/Dry foot control		•	•		•	•		•
Oil collector for handpieces		•	•		•	•		•
Stainless steel tray		•	•		•	•		•
Asepsis design		•	•	•	•	•	•	•
Self contained water system		•	•	•	•	•		•
Mounts on 2" post (included)		•	•	•	•	•		
6' Umbilical with junction box and utilities		•	•	•	•	•		•
Air/Water quick disconnects		•		•	•			
Hub mount for light		•		•	•		•	
Assistants telescoping vacuum package (Syringe, HVE, SE) with canister		•		•				

OPTIONS

Add handpiece position	4467	\$110	•	•		•	•	
Delete handpiece position		-\$54.40	•	•		•	•	•
Delete syringe, tubing and holder		-\$105.80	•	•		•	•	•
Add circulating syringe and water heater	3575	\$449	•	•		•	•	•
Sub chassis for electric motor handpiece control	4230	\$45	•	•		•	•	•
ISO 5 hole light source system w/bulb, transformer and straight gray tubing	8744	\$665	•	•		•	•	•
ISO 6 pin light source system w/transformer and straight gray tubing	8749	\$600	•	•		•	•	•
Add chip air or scaler button on foot control		\$45	•	•		•	•	•
48" Post in lieu of 20" with Pelton Light Bushing		\$55	•	•	•			
Integrated scaler	8053	\$1,085	•	•		•	•	•
Delete utility center and junction box		-\$310	•	•	•	•	•	•
Delete assistant's telescoping vacuum package		-\$450	•		•			

Automatic Series 4 Arm Mounted Delivery System



COMPACT HEAD
11 3/4 W x 11 D

Shown w/ optional top post
mount adapter 8923



Thru Post Mount
4477 & 4478
(standard
configuration)

MANUAL UNIT AVAILABLE PON RE EST

DIMENSIONS
PAGE 3

GRAY | WHITE

4477	4478	Arm Mounted Automatic Control for 3 HP.....	\$1,955.00
8956	8923	Top Post Mount Adapter.....	\$150.00
	8035	Touch-Up Paint.....	\$20.00

PAGE 3



- Smooth lines for easy cleaning
- Automatic control for 3 handpieces
- Flex arm system mounts on 2" post
- Individual drive air & water coolant controls
- Stainless steel tray (PN 8013) & holder (PN 4280)
- Includes 8 Conductor Low Voltage Electrical Cable
- Precision Comfort™ Syringe

- Wet/Dry foot control
- Optional storage area for high-tech accessories
- Available in left or right hand configuration
- 7 ft. of 5/8" Umbilical
- Specify tubing color & type

Series 4 Cabinet or Wall Mounted Automatic Delivery System



COMPACT HEAD
11 3/4 W x 11 D



42
Wall mount
adapter included



PAGE 3
DIMENSIONS
PAGE 3

MANUAL UNIT AVAILABLE PON RE EST

PN 4470 Automatic Control for 3 HP	\$935.00
PN 8035 Touch-Up Paint (White)	\$20.00

- Smooth lines for easy cleaning
- Individual drive air & water coolant controls
- Precision Comfort™ Syringe
- Optional storage area for high-tech accessories
- Wet/Dry foot control
- 6 ft. of 5/8" Umbilical
- Specify tubing color & type

Series 4 Delivery Systems

	Arm Mounted	Cabinet Mounted
White	4478	4470
Gray	4477	
Price	\$1,955	\$935

STANDARD FEATURES

Automatic Controls for 3 Handpieces (expandable to 4)	●	●
Compact Size Delivery Head	●	●
Wet/Dry Foot Control	●	●
Precision Comfort® Syringe	●	●
Individual Non Retracting Water Coolant	●	●
Handpiece Purge Valve	●	●
50" Flex Arm for 2" Post w/ Air Brake	●	
Cabinet Mount		●
Stainless Steel Tray	●	
Oil Collector for Handpiece Exhaust	●	●

OPTIONS	PN	PRICE		
White Wall Mount Bracket	8647	\$150.00	●	
Gray Wall Mount Bracket	8349	\$150.00	●	
40" Chair Mount Post w/ P&C Bushing	8444	\$220.00	●	
White Top Post Mount Adapter	8923	\$150.00	●	
Gray Top Post Mount Adapter	8956	\$150.00	●	

CHAIR ADAPTERS *Refer to Mounts, Arms & Holders Section for complete listing*

Belmont Bel 20	8493	\$300.00	●	
A-dec® Priority	8515	\$300.00	●	
Pelton & Crane® 1500/Marus Empire	8918	\$300.00	●	

LIGHT SOURCE SYSTEMS *Refer to Handpiece Accessories Tab for complete listing*

ISO 5 hole light source system w/ bulb, transformer & 1 straight gray tubing	8744	\$665.00	●	●
ISO-C 6 pin light source system w/ transformer & 1 straight gray tubing	8749	\$600.00	●	●
Tubing (expandable to 3) 110V **Specify HP Position**				

CONTROL OPTIONS

Add Handpiece Position	4467	\$110.00	●	●
Delete Handpiece Position		-\$54.40	●	●
Delete Syringe, Tubing & Holder		-\$105.80	●	●
Circulating Syringe & Water Heater	3575	\$449.00	●	●
Add Chip Air or Scaler Button on Foot Control		\$45.00	●	●
Piezo Scaler Installed	8053	\$1,085.00	●	●
Sub Chassis for Electric Motor Handpiece Control - (Internal Size 6" x 1 3/4" x 2 1/2")	4230	\$45.00	●	●

ACCESSORIES *Refer to DCI Catalog for complete regulator & water system options*

Utility Center w/ White Box & Regulators	8329	\$398.00	●	●
2-Liter Asepsis Deluxe Water System-Cabinet Mount	8184	\$220.00		●
2-Liter Asepsis Deluxe Water System w/ 2" Post Clamp	8928	\$315.00	●	
White Telescoping Arm w/ Premium Vacuum Valves	5476	\$725.00	●	



Trouble Shooting a Unit? Find videos for diagnosing slow handpieces, combo valve & shut off block repair, needle valve & cartridge valve adjustment at the DCI Video Library!

Series 4 Cart Systems



PAGE 3

IPBDIMENSIONS
PAGE 3

PN 4511	Automatic Controls for 3 HP.....	\$2,450.00
PN 4540	Automatic Controls for 3 HP w/ Premium Vacuum.....	\$2,745.00
PN 4245	Manual Control for 2 HP.....	\$1,615.00
PN 4240	Manual Control for 2 HP w/ Premium Vacuum.....	\$2,050.00
PN 8035	Touch-Up Paint (White).....	\$20.00

- Sturdy all-metal control head w/ smooth lines for easy cleaning
- Automatic control for 3 handpieces or manual control for 2 handpieces
- Individual drive air & water coolant controls
- SS tray & holder
- Precision Comfort™ Syringe
- Wet/Dry foot control (Auto units only)
- Deluxe clean water system
- 1/4" Water Q.D. w/ flow control
- 7 ft. of 1 1/4" Umbilical
- Utility center w/ white box
- Specify tubing color & type



Operatory Support Cart

PN 4227	Operatory Support Cart w/ U-Frame.....	\$865.00
PN 8400	Storage Well w/ Organizer.....	\$270.00
PN 4220	Operatory Support Cart w/ Assistant's Package, U-Frame.....	\$1,425.00
PN 4228	Operatory Support Cart w/ H-Frame.....	\$865.00
PN 4221	Operatory Support Cart w/ Assistant's Package, H-Frame	\$1,425.00

- Solid 1" white laminate top w/ gray edge banding
- 19" x 22" work surface

Series 4 Cart Systems

	Automatic Carts	Manual Carts
	4511 4540	4240 4245
	\$2,450 \$2,745	\$2,050 \$1,615

STANDARD FEATURES

Automatic Controls for 3 Handpieces (expandable to 4)	●	●		
Manual Control for 2 Handpieces (1 wet & 1 dry)			●	●
Wet/Dry Foot Control	●	●		
Economy Foot Control			●	●
Precision Comfort® Syringe	●	●	●	●
Individual Non Retracting Water Coolant	●	●	●	●
Handpiece Purge Valve	●	●		
U-Frame	●	●	●	●
Stainless Steel Tray	●	●	●	●
Oil Collector for Handpiece Exhaust	●	●		
Deluxe Clean Water System	●	●	●	●
Utility Center w/ White Box & Regulators	●	●	●	●
1/4" Water Quick Disconnect	●	●	●	●
Precision Comfort® Premium HVE & SE Vacuum Package		●	●	

LIGHT SOURCE SYSTEMS *Refer to Handpiece Accessories Tab for complete listing*

ISO 5 hole light source system w/ bulb, transformer & 1 straight gray tubing	8744	\$665.00	●	●		
ISO-C 6 pin light source system w/ transformer & 1 straight gray tubing	8749	\$600.00	●	●		
Tubing (expandable to 3) 110V **Specify HP Position**						

CONTROL OPTIONS

Piezo Scaler Installed	8053	\$1,085.00	●	●	●	●
Add Handpiece Position	4467	\$110.00	●	●		
Delete Handpiece Position		-\$54.40	●	●		
Delete Syringe, Tubing & Holder		-\$105.80	●	●	●	●
Circulating Syringe & Water Heater	3575	\$449.00	●	●	●	●
Add Chip Air or Scaler Button on Foot Control		\$45.00	●	●		
Delete Utility Center & Box		-\$310.00	●	●	●	●
Delete Clean Water System		-\$90.00	●	●	●	●



Trouble Shooting a Unit? Find videos for diagnosing slow handpieces, combo valve & shut off block repair, needle valve & cartridge valve adjustment at the DCI Video Library!

ORDER GUIDE ON PAGE 30

Alternative Arm-Mounted Delivery System



43

PN 4376	Arm-Mounted Automatic Control for 3 Handpieces.....	\$1,795.00
PN 4128	Arm-Mounted Manual Control 4120 for 1 Wet & 1 Dry Handpiece	\$1,535.00
PN 4129	Arm-Mounted Manual Control 4121 for 2 Wet Handpieces	\$1,615.00
PN 8923	Top Post Mount Adapter (White).....	\$150.00
PN 8035	Touch-Up Paint (White).....	\$20.00

- Automatic control for 3 handpieces or manual control for 2 handpieces
- Flex arm system mounts on 2" post
- Individual drive air & water coolant controls
- Instrument tray holder
- Precision Comfort® Syringe
- Wet/Dry foot control (4376 only)
- Dimensions: Horizontal reach 62 3/4" Fully closed: 30 3/4" L x 21 1/4" W
- 7 ft. of 5/8" Umbilical
- Specify tubing color & type



Alternative Cabinet or Wall Mounted Automatic Delivery System



PN 4370	Alternative Automatic Control for 3 Handpieces.....	\$775.00
----------------	---	----------

- Individual drive air & water coolant controls
- Multiple mounting options (refer to pages 99–101)
- Precision Comfort® Syringe
- Wet/Dry foot control
- Dimensions: 12" x 8" x 5"
- 6 ft. of 1/2" Umbilical
- Specify tubing color & type



Alternative Cabinet or Wall Mounted Manual Delivery System



PN 4120	Manual Control for 1 Wet & 1 Dry Handpiece	\$525.00
PN 4121	Manual Control for 2 Wet Handpieces	\$570.00

- Manual control for 2 handpieces
- Precision Comfort® Syringe
- Disc type foot control
- Multiple mounting options (refer to pages 99–101)
- Dimensions: 10" x 3" x 4"
- 6 ft. of 1/2" Umbilical
- Specify tubing color & type

Panel & Horizontal Surface Mount Automatic Unit for 3 Handpieces



PN 4415	Panel-Mount Automatic Control	\$845.00
PN 4416	Horizontal-Mount Automatic Control	\$1,010.00
PN 4445	Optional Swing Arm Holder Bar	\$140.00



- Individual drive air & water coolant controls
- Precision Comfort® Syringe
- Wet/Dry foot control
- Control dimensions: 8 1/2" x 1 3/4" x 5"
- Cut-out dimensions: 7 1/2" x 1 1/2"
- Horizontal system can be mounted to top or bottom of surface
- Specify tubing color & type

Panel & Horizontal Surface Mount Manual Unit for 2 Handpieces



PN 4400	Panel-Mount Manual Control for 1 Wet & 1 Dry Handpiece	\$545.00
PN 4401	Horizontal-Mount Manual Control for 1 Wet & 1 Dry Handpiece	\$610.00
PN 4445	Optional Swing Arm Holder Bar	\$140.00



- Individual drive air & water coolant controls
- Precision Comfort® Syringe
- Wet/Dry foot control
- Control dimensions: 8 1/2" x 1 3/4" x 5"
- Cut-out dimensions: 7 1/2" x 1 1/2"
- Horizontal system can be mounted to top or bottom of surface
- Specify tubing color & type



DCI III Manual Unit

PN 4502 DCI III Manual Unit \$510.00

- Two HP manual selection for 1 Wet & 1 Dry Handpiece
- Precision Comfort® Syringe

- 6 ft. 1/2" Umbilical
- Specify tubing color & type



DCI II Manual Unit

PN 4102 DCI II Manual Unit \$440.00

- Two HP manual selection for 1 Wet & 1 Dry Handpiece
- Precision Comfort® Syringe

- Dimensions: 6" x 4 1/4" x 2 1/4"
- 6 ft. 1/2" Umbilical
- Specify tubing color & type



Air & Water Auxiliary Panel

PN 6559 Gray Auxiliary Panel \$175.00

PN 6569 White Auxiliary Panel \$175.00

- 1/4" Water quick connect w/ flow control
- Dimensions: 3" x 4"

- 3/8" Air quick connect w/ on/off switch
- Available in white or gray



Pelton & Crane Executive Conversion Kit

PN 4446 Executive Conversion Kit \$2,645.00

PN 4447 Executive Conversion Kit Less Utility Center \$2,435.00

- Precision Comfort™ Premium Autoclavable Vacuum Valves include two 5105 HVEs & one 5215 SE
- Includes 8323 utility center
- Automatic handpiece holders w/ manual lockout toggle

- Variable speed disc-type foot control w/ wet/dry toggle
- Precision Comfort® three-way syringe
- Dimensions: 12 1/2" x 18" x 4 1/2"
- 6 ft. Umbilical

FEATURED



Laboratory Handpiece Control Kit

PN 4198 Laboratory Handpiece Control Kit \$168.00

- Requires air only
- 4-hole handpiece tubing 432T
- Handpiece holder & bar 5910 & 5929

- Economy foot control 6301
- Air regulator w/ gauge 7220 & 7260

DCI Alternative Controls

	Arm Mounted			Cabinet Mounted			Specialty Units			
	4376	4128	4129	4370	4120	4121	4502	4102	4446	4447
	\$1,795	\$1,535	\$1,615	\$775	\$525	\$570	\$510	\$440	\$2,645	\$2,435

STANDARD FEATURES

Automatic Controls for 3 Handpieces	●			●					●	●
Manual Control for 2 Handpieces (1 wet & 1 dry)		●			●		●	●	●	●
Manual Control for 2 Handpieces (2 wet)			●			●				
Wet/Dry Foot Control	●			●		●			●	●
Economy Foot Control		●	●	●	●	●	●	●	●	●
Precision Comfort® Syringe	●	●	●	●	●	●	●	●	●	●
Individual Non Retracting Water Coolant	●	●	●	●	●	●	●	●	●	●
50" Flex Arm for 2" Post w/ Air Brake	●	●	●							
Cabinet Mount				●	●	●	●	●	●	●

OPTIONS

	PN	PRICE								
Wall Mount Bracket	8647	\$150.00	●	●	●					
40" Chair Mount Post w/ P&C Bushing	8444	\$220.00	●	●	●					
White Top Post Mount Adapter	8923	\$150.00	●	●	●					

CHAIR ADAPTERS *Refer to Mounts, Arms & Holders section for complete listing*

Belmont Bel 20	8493	\$300.00	●	●	●					
A-dec® Priority	8515	\$300.00	●	●	●					
Pelton & Crane® 1500/Marus Empire	8918	\$300.00	●	●	●					

LIGHT SOURCE SYSTEMS *Refer to Handpiece Accessories Tab for complete listing*

ISO 5-Hole Economy HP										
Illumination System	8801	\$252.00	●	●	●					
ISO 6-Pin Economy HP										
Illumination System	8791	\$190.00	●	●	●					

CONTROL OPTIONS

Piezo Scaler Installed	8053	\$1,085.00	●			●				
Delete Handpiece Position		-\$54.40	●			●				
Delete Syringe, Tubing & Holder		-\$105.80	●	●	●	●	●	●	●	●
Circulating Syringe & Water Heater	3575	\$449.00	●	●	●	●	●	●	●	●
Add Chip Air or Scaler Button on Foot Control		\$45.00	●			●				

ACCESSORIES *Refer to DCI Catalog for complete regulator & water system options*

Utility Center w/ Gray Box & Regulators	8321	\$375.00	●	●	●		●	●	●	
2-Liter Standard Water System Cabinet Mount	8144	\$155.00	●	●	●		●	●	●	
2-Liter Deluxe Water System Cabinet Mount	8143	\$190.00	●	●	●		●	●	●	
750ml Standard Water System Cabinet Mount	8145	\$150.00	●	●	●		●	●	●	
750ml Deluxe Water System Cabinet Mount	8139	\$185.00	●	●	●		●	●	●	
Adapter for 2" Post, White Powder Coat	8162	\$60.00	●	●	●					



Trouble Shooting a Unit? Find videos for diagnosing slow handpieces, combo valve & shut off block repair, needle valve & cartridge valve adjustment at the DCI Video Library!



SHOWN WITH
OPTIONAL
SO. ND COVER

DIMENSIONS
PAGE 44

DC-Series Deluxe Oil-less Air Compressors

PAGE 43

IPB



PN DC0103	3 ser / 1 HP / Single Head / 10 Gal Tank / 115 Vac	\$5,475.00
PN DC0203	3 ser / 1 HP / Single Head / 10 Gal Tank / 230 Vac	\$5,475.00
PN DC1203	3 ser / 1 HP / Single Head / 20 Gal Tank / 230 Vac	\$6,160.00
	Expandable to 2 Heads	
PN DC2206	ser / 2 HP / Dual Head / 20 Gal Tank / 230 Vac	\$8,585.00
PN DC3210	10 ser / 3 HP / Triple Head / 30 Gal Tank / 230 Vac	\$12,290.00

- Completely oil-less design
- Extremely quiet
- 115 / 230 dual voltage motors
- 3 stage filtration / drying system
- Emergency quick-connect
- 3 year warranty on compressor parts
- 5 year warranty on heads
- Includes remote air intake kit
- Includes 24 vac low voltage controls

NEW ITEM



DIMENSIONS
PAGE 44

C Series Oil-less Air Compressors

PAGE 43

IPB



PN C0103	3 ser / 1 HP / Single Head / 10 Gal Tank / 115 Vac	\$5,105.00
PN C0203	3 ser / 1 HP / Single Head / 10 Gal Tank / 230 Vac	\$5,105.00
PN C1203	3 ser / 1 HP / Single Head / 20 Gal Tank / 230 Vac	\$5,795.00
	Expandable to 2 Heads	
PN C2206	ser / 2 HP / Dual Head / 20 Gal Tank / 230 Vac	\$8,210.00
PN C3210	10 ser / 3 HP / Triple Head / 30 Gal Tank / 230 Vac	\$12,000.00

- Completely oil-less design
- Extremely quiet
- 115 / 230 dual voltage motors
- 3 stage filtration / drying system
- Emergency quick-connect
- 3 year warranty on compressor parts
- 5 year warranty on heads

NEW ITEM



Looking for a convenient, portable solution?

PN 2124 Portable Wheel Kit
for 10 Gallon Compressor.....\$325.00



P Series Portable Lubricated Air Compressor

PN **P050** 1 ser / 1/2 HP / Single Head / . 3 Gal Tank / 115 Vac. \$2,846.00

- 115 Volt 50- 0 Hertz
- quiet operation
- 2.1 CFM
- Only 44 lbs
- Compact size 11 x 11 x 11 /
- month warranty



203 DIMENSIONS PAGE 44



2030

Air Compressor Accessories

PN 2045	Compressor Head Upgrade Kit For Converting D C1203 to a D C220	\$2,500.00
PN 2068	Sound Cover for C0103 / DC0103, 115 Volt Fan, 10 Gal Tank	\$845.00
PN 2069	Sound Cover for C0203 / DC0203, 230 Volt Fan, 10 Gal Tank	\$845.00
PN 2043	Sound Cover for C1103 / DC1103, 115 Volt Fan, 20 Gal Tank	\$895.00
PN 2049	Sound Cover for C1203 / C220 / DC1203 / DC220 , 230 Volt Fan, 20 Gal Tank	\$895.00
PN 2059	Sound Cover for C3210, 230 Volt Only, 30 Gal Tank	\$950.00
PN 2036	Stacking Rack for Compressor or Vacuum	\$1,625.00
PN 2030	Low Voltage Control Box input 115 Vac or 230 Vac.	\$420.00
PN 2847	Buck / Boost Transformer 1.0 VA	\$950.00
PN 2033	Remote Air Intake Kit for C1103 C1203 / C220 , 20 Gal Tank Compressors	\$105.00
PN 2057	Remote Air Intake Kit for C3210, 30 Gal Tank Compressors	\$125.00

NOW THE OEM NEED THE DCI PN

use the DCI Cross Reference Database on www.dcionline.com to search from over 5,000 OEM PN's & find what you're looking for fast!

Alternating Control Panel

PN 2089 For Vacuum Pumps

\$1,900.00



- Control panel allows 2 systems to run in tandem
- Automatically shuts off one machine when demand is low
- Second machine turns back on when demand increases
- This panel alternates usage to provide equal run time

- L 50 A listed panel
- Automatic or manual mode
- One remote switch controls entire system
- 30 amp maximum per machine
- Requires 2 30 amp power sources
- Dimensions 1 H x 14 W x D



Refrigerated Air Drying System

PN 2797 10cfm 1 H x 1/4 D x 1 3/4 W

\$1,800.00

PN 2798 25cfm 20 H x 1/4 D x 1 3/4 W

\$2,110.00

- Can be used w/ any air compressor
- Available in 10 & 25 cfm models
- 1/2 NPT connections
- CFC-FREE R134A Refrigerant
- 115 volt AC

SHIPS VIA FREIGHT ONLY

Vacuum System Cleaner

PN 5835 Eco Vac 1-Pint Bottle

\$15.00

PN 5837 Eco Vac 1/2 Gallon Bottle

\$49.95

PN 5836 Eco Vac Case of 12 1-Pint Bottles

\$165.00



Eco Vac cleans vacuum systems and eliminates odors associated with organic waste by biologically degrading the waste into natural elements. Cleans and deodorizes without foaming or use of harsh chemicals. Eco Vac utilizes safe non-pathogenic, non-opportunistic, Class 1 bacteria.

One half ounce of Eco Vac mixed with one quart of warm water, drawn through the vacuum lines daily in each operatory, establishes a continuous biological action.

- Continuous cleaning action
- Environmentally safe
- Contains no hazardous materials
- Non-corrosive liquid safe to handle & store

- Easy to use, measured 1 pint dispenser holds 32 applications
- 1/2 gallon dispenser holds 12 applications
- Non-Toxic, non-pathogenic

DIMENSIONS
PAGE 44



Surgical Suction Systems

PN S1102	1-2 ser / Single 3/4 HP 20 Gal Tank / 115 Vac..	\$4,875.00
PN S1202	1-2 ser / Single 3/4 HP 20 Gal Tank / 230 Vac..	\$4,875.00
PN S2104	2-4 ser / Dual 3/4 HP 20 Gal Tank / 115 Vac..	\$6,995.00
PN S2204	2-4 ser / Dual 3/4 HP 20 Gal Tank / 230 Vac..	\$6,995.00
PN S3206	4- ser / Triple 3/4 HP 30 Gal Tank / 230 Vac..	\$8,990.00
PN 2507	Surgical Bottle Assembly..	\$500.00

14 Hg 22 Hg operating range
Oil-less pump design
For usage in all surgical applications

3 Year warranty on surgical vacuum parts
5 year warranty on pumps



Water Security Systems

PN W1100	3/4 System w/ Filter & 24 Volt Relay..	\$945.00
PN W2100	1 System w/ Filter & 24 Volt Relay..	\$1,100.00

Filters incoming water supply & shuts off water from control panel

See page 12 for Water Security System Replacement Parts

DIMENSIONS
PAGE 44

FEATURED

Control Panels

NEW & IMPROVED LED BACKLBS
NOW COME STANDARD




PN 2900	Single Switch Panel / Expandable to 4..	\$300.00
PN 2901	Dual Switch Panel / Expandable to 4..	\$340.00
PN 2902	Triple Switch Panel / Expandable to 4..	\$380.00
PN 2904	Quad Switch Panel..	\$420.00
PN 2723	Remote 24 Volt Red LED Light Switch..	\$40.00

Fits standard 2 gang electrical box
Includes custom labeling kit
24 Volt lighted LED rocker switches


TRYING TO
SOLVE A TOUGH
TECHNICAL
PROBLEM

Visit the DCI Service Forum at www.dcionline.com
or the DCI Video Library where you can collaborate
with fellow technicians and find the answers to the
most challenging dental repairs!


Series 5 LED Operatory System Mount & Post Mount Light




Series 5 LED Operatory Track Light





Series 4 Over the Patient Automatic Delivery System




Series 4 Swing Mount Pole
w/ Automatic Delivery System & Light Pole




Series 4 Side Delivery Automatic Delivery System w/ Mount

Series 4 Over the Patient Automatic Delivery System
w/ PM & Assistant's Package


**Series 4 Arm Mounted
Automatic Control
44 & 44**




**Series 4 Automatic Control
44 0**




Series 4 Automatic Cart




**DCI III Manual Unit
4502**




Series 4 Automatic Control




Series 4 Automatic Cart Unit 4511, 4540 & 4542




DCI III Manual nit 4502




**Alternative Manual Control
4120 & 4121**




Alternative Automatic Control 43 0




Alternative Automatic Cart 4411, 4440 & 4442




Alternative Manual Cart 4145, 4140 & 4142




Panel & Horizontal Surface Mount Control 4415 & 441




Panel & Horizontal Surface Mount Control 4400 & 4401



Oil-Less Compressor



Lubricated Compressor



MACHINE PHYSICAL DIMENSIONS & WEIGHT

Compressors

	LENGTH	DEPTH	HEIGHT	SHIPPING WEIGHT
C0103	32"	23"	27"	150 LBS
C0203	32	23	2	150 LBS
C1103	3	24	2 1/2	220 LBS
C1203	3	24	2 1/2	220 LBS
C2206	3	24	2 1/2	20 LBS
C3210	4	24	32	345 LBS
DC0103	32	23	2	10 LBS
DC0203	32	23	2	10 LBS
DC1103	3	24	31 1/2	250 LBS
DC1203	3	24	31 1/2	230 LBS
DC2206	3	24	31 1/2	20 LBS
DC3210	4	24	34	355 LBS

Surgical Vacuums

	LENGTH	DEPTH	HEIGHT	SHIPPING WEIGHT
S1102	37"	24"	29"	210 LBS
S1202	3	24	2	210 LBS
S2204	3	24	2	250 LBS
S3206	4	24	30	335 LBS

Stacking Rack

	LENGTH	DEPTH	HEIGHT
2036	39"	29"	34"

Water Security Systems

	LENGTH	DEPTH	HEIGHT
W1100	17"	24"	5"
W2100	1	24	5