



IMAGING | UTILITY ROOM | ACCESSORIES



PARTS
CATALOG
2012

AIR TECHNIQUES equipped for life
1962 50th Anniversary 2012

Table of Contents

Equipped for Life™ with Air Techniques

Air Techniques was founded in 1962 on the principal of providing top-of-the-line dental equipment and superb customer service. We are thrilled to be celebrating our 50th anniversary and appreciate the continued commitment from our network of dealers and customers.

Our mission at Air Techniques has stood true for 50 years: to manufacture products that are reliable and technologically advanced to support our customers' practices, while increasing their return on investment.

We understand that circumstances may arise, leading the way for parts to be replaced or exchanged and, therefore, causing your practice to experience downtime. To help prevent this situation from occurring, Air Techniques has established a Return Authorization Policy to ensure a quick and simple process. We also offer flexible options for the exchange of rebuilt parts.

We are proud to carry on our legacy as a U.S. manufacturer based in New York. Our team is committed to manufacturing products that will empower both our partners and customers.

We are truly grateful for the continued support from our dealer network and customers, and we thank you for choosing Air Techniques

Best Regards,

Richard S. Povak
President



Mission / Discontinued Parts	2
Return Authorization Policy / Contact Information	3
Replacement Parts Exchange Policy / Warranty	4
Air Compressors	6-41
AirStar with Membrane Dryer, AS10 to AS70	6-19
AirStar with Desiccant Dryer, AS10 to AS70	20-33
Buck/Boost Transformers, Contactors & Accessories for Air Systems	34
Vacuums	36-52
Mojave Dry Vacuum	36-37
STS Dry Vacuum	38-39
CAS, Acadia & Guardian Amalgam Separators	39-41
STS Accessories	42-43
VacStar VS20 to VS80H Wet Vacuum	44-51
Air/Water Separator, Hydromiser & Vacuum Equalizer	52
Accessories	54-55
Remote Water Control Valves	54
Remote Control Panels	54-55
Film Processing	56-72
A/T 2000XR	56-62
A/T 2000 Plus	62-66
Accessories for A/T 2000 Plus & 2000XR	66
Panoramic DLL for A/T 2000 Plus & 2000XR	66-67
Peri-Pro III	68-70
Peri-Pro & Accessories for Peri-Pro I & III	70-71
Daylight Loader & Duplicator / Daylight Loader -PPI,II,III	71-72
Main Transport Peri-Pro I & III	72
Intraoral X-Ray - Provecta 70	74-75
Digital Imaging	76-78
ScanX Digital Imaging System & Eraser	76-77
ScanX Flat Rate Repair/Upgrade Program	78
Video Imaging	80-81
Polaris Intraoral Camera	80
Acclaim Intraoral Camera	80
Spectra Caries Detection Aid	80
A/T SLC & VistaCam Omni & Omni IC4	80-81
Exchange Parts	82
STS, VacStar, AirStar, A/T 2000 Film Processors & ScanX	82

Air Techniques Company Mission

To develop and bring to our customers the products they value. Our products are superior by invention, are professionally supported and provide years of lasting service to their users.

Air Techniques endorses and supports high quality and customer service as a company-wide way of life, focusing on managing the prevention of problems and the continuous improvement of our products to assure customer requirements are met.

Our success is achieved by empowering employees on every level to prevent defects (even where this might include the temporary stoppage of work or shipments) and improve and refine our processes and products.

Discontinued Products and Parts

1. A/T 2000 Original Model 43000 manufactured from 1983-1994. Electronics Module P/N 43415X, Base PCB P/N 43420 and Electronics Module P/N 43400x are no longer available.
2. A-Series Compressor-Dri Air Compressors manufactured from 1972-1985. Replacement compressor heads are no longer available.
3. V-Series Compressor-Dri Air Compressors manufactured from 1985-1987. Replacement compressor heads are no longer available.
4. Turbo 2000 Dry Vacuum System manufactured from 1995-2001. Replacement pump assemblies are no longer available.
5. 70 Series Compressors manufactured from 1968-1998. Limited replacement parts availability.
6. 50/50,000 Series Compressors manufactured from 1969-1999. Limited replacement parts availability.
7. Dento Dri 3100R manufactured from 1962-1998. Limited replacement parts availability.
8. Dento Dri 3300G manufactured from 1963-1986. Limited replacement parts availability.
9. VistaCam Escort
10. Oil Type Compressor Parts
11. Original AirStar 1, 2-1, 2-2, 3, 5, 7 Heads
12. Vacstar 2, 4, 5, 8.
13. CAS 8-gallon PN 54300 Old Style.
14. 4-gallon CAS for Wet Suction Pumps.
15. 8-gallon CAS for Wet Suction Pumps.

Suggested Retail Prices are effective April 2, 2012.

This price list supersedes all other Parts Price Lists.

All prices are subject to change without notice.

Please see our website www.airtechniques.com for updates.

Return Authorization Policy

All part and component returns and replacement of equipment under warranty require an RMA.

Please follow this procedure so that we can process our returned equipment as soon as possible:

1. Before a product can be shipped back to Air Techniques, an RMA number must be obtained.
2. To obtain an RMA call our Technical Service department at (800) 247-8324 with:
 - a) a model number
 - b) a serial number for the defective product
 - c) detail on why the part is being returned
3. Only the Air Techniques product and quantity specified on the original RMA request can be returned with the RMA number issued.
 - a) If returning additional products to Air Techniques, a new RMA number will be required.
 - b) If we receive a shipment containing products not authorized for return on that RMA number, they will be returned as is.

Contact Information

For general and immediate questions please call:

Toll Free in the US and Canada: 1- 800-AIR TECH (247-8324)
Monday through Friday, 8:00am – 8:00pm (EST)

All shipments are FOB Melville, NY or Corona, CA

Replacement Parts Exchange Policy

Replacement parts are available on an "exchange basis" for some Air Techniques equipment. Parts available on this basis are listed in this price list. Replacement "exchange" parts may be ordered in either of two ways:

Advance Replacement

1. Order the replacement part needed.
2. The part is shipped and invoiced at the "outright sale price."
3. Return the old part (freight prepaid) to the factory, after obtaining an RMA number. See previous page.
4. Credit is issued for the difference between the "outright sale price" and the "exchange price."

Direct Exchange

1. Return the old part (freight prepaid) to the factory, after obtaining an RMA. See previous page.
2. Request that the part be replaced.
3. The replacement part is shipped (freight collect).
4. The "exchange price" is invoiced.

Warranty for Parts

Air Techniques parts are warranted to be free from defects in material and workmanship from the date of installation for a period of:


Replacement Compressor Heads, Oil-less	12 month limited warranty
Replacement Vacuum Pumps	12 month limited warranty
Rebuilt A/T 2000, vvPlus & XR Exchange Parts	6 month limited warranty
Peri-Pro Transports	6 month limited warranty
ScanX, ScanX Intraoral	6 month limited warranty
Provecta Intraoral X-Ray	6 month limited warranty
All Other Parts	3 month limited warranty

Any item returned to our factory in Melville, New York through an Air Techniques authorized Dealer will be repaired or replaced at our discretion at no charge provided that our inspection shall indicate it to have been defective. Dealer labor, shipping and handling charges are not covered by this warranty.

This warranty does not apply to damage due to shipping, misuse, careless handling or repairs by unauthorized service personnel. Warranty void if installed or serviced by anyone other than authorized Air Techniques Dealer service personnel.

Air Techniques, Inc. is not liable for indirect or consequential damage or loss of any nature in connection with this equipment. This warranty is in lieu of all other warranties expressed or implied. No representative or person is authorized to assume for us any liability in connection with the sale of our equipment.



Air Compressors with Membrane Dryer



Air Compressors with Membrane Dryer

Part Number	Description	Price
AirStar AS10 Parts Breakdown		
89948	1/4" OD, Black Nylon Tubing	\$ 2.00
86726	10mm Flamex Tubing, Cooler to Membrane Dryer	\$ 20.00
86727	12mm Flamex Tubing, Membrane Dryer to Main Tank Check Valve	\$ 24.00
85617	5 Micron Filter	\$ 251.00
87180	115V Cooler Assembly	\$ 384.00
85911	Drain Valve Assy	\$ 30.00
58038	Elbow, 3/8 MNPT x 3/8 FNPT, Tank Inlet	\$ 5.00
56121	Elbow, 1/4 MNPT x 1/4 Push, PConn	\$ 5.00
86729	Elbow, 10mm x 10mm, PConn	\$ 14.00
86747	Elbow, 10mm Push x 3/8", PConn	\$ 21.00
86749	Elbow, 12mm x 3/8", PConn	\$ 30.00
85603	Filter Element; Air: for Norgren	\$ 80.00
86193	Filter Element Replacement Kit, 5 Micron, Wilkerson	\$ 38.46
87067	Flex Braided Hose	\$ 81.00
S85482	Flex Pressure 6' Rubber Hook-Up Hose	\$ 38.00
89831	Intake Filter, Compressor Head	\$ 53.00
85471	Leveling Foot Kit	\$ 79.00
87045	Main Tank Check Valve	\$ 33.00
87175	Membrane Dryer Assembly	\$ 833.00
60013	Moisture Monitor Assy	\$ 23.00
85890	Motor Replacement (See Exchange Parts)	\$ 2,825.00
86019	Nipple to Tank, 1/2 MNPT x Brass, PConn	\$ 23.00
87036	Outlet Manifold	\$ 85.00
60011	Pressure Gauge	\$ 37.00
60017	Pressure Relief Valve, 125 PSI 60017 (Standard)	\$ 36.00
60510	Pressure Relief Valve, 140 PSI	\$ 54.00
87176	Pressure Switch, 85 to 115 PSI	\$ 53.00
85887	Tank Outlet Shutoff Valve	\$ 23.00
87108	Drain Tube and Drip Tray Kit	\$ 42.00
Electrical Breakdown		
85833	Capacitor, 65 uF, 370 VAC (2 per AS21)	\$ 62.00
117061	Circuit Breaker Low Voltage, 24 VAC, 1/2 Amp	\$ 16.00
85890	Compressor Head Motor Assembly, 115 VAC, 1 Cylinder	\$ 2,825.00
89345	Contactor, 40 Amp Line Voltage with 24 VAC Coil	\$ 159.00
87076-1	Cooler Assembly 115V Fan	\$ 28.00
87151	Line Cord Assembly, 115 VAC	\$ 61.00
87158	Motor Harness Replacement Kit	\$ 45.00
86065	Power On/Off DPST Switch and Circuit Breaker Motor	\$ 97.00
87176	Pressure Switch, 85 to 115 PSI	\$ 53.00
51952	Sound Cover Fan, 115 VAC	\$ 65.00
85961	Sound Cover Kit	\$ 950.00
87126	Sound Cover Receptacle	\$ 50.00
86637	Transformer, 115 VAC Step-down to 24 VAC	\$ 55.00
Additional Items		
85491	Remote Air Intake (Optional, Not Included with AirStar)	\$ 220.00
80183	Vibration Pad (Optional Not Included with AirStar)	\$ 41.00



Air Compressors with Membrane Dryer



Air Compressors with Membrane Dryer

Part Number	Description	Price
AirStar AS21 Parts Breakdown		
89948	1/4" OD, Black Nylon Tubing	\$ 2.00
86726	10mm Flamex Tubing, Cooler to Membrane Dryer	\$ 20.00
86727	12mm Flamex Tubing, Membrane Dryer to Main Tank Check Valve	\$ 24.00
85617	5 Micron Filter	\$ 251.00
87180	115V Cooler Assembly	\$ 384.00
85911	Drain Valve Assy	\$ 30.00
58038	Elbow, 3/8 MNPT x 3/8 FNPT, Tank Inlet	\$ 5.00
56121	Elbow, 1/4MNPT x 1/4 Push, PConn	\$ 5.00
86729	Elbow, 10mm x 10mm, PCConn	\$ 14.00
86747	Elbow, 10mm Push x 3/8", PCConn	\$ 21.00
86749	Elbow, 12mm x 3/8", PCConn	\$ 30.00
85603	Filter Element; Air: for Norgren	\$ 84.00
86193	Filter Element Replacement Kit, 5 Micron, Wilkerson	\$ 40.00
87049	Flex Braided Hose	\$ 87.00
S85482	Flex Pressure 6' Rubber Hook-Up Hose	\$ 38.00
89831	Intake Filter, Compressor Head	\$ 53.00
85471	Leveling Foot	\$ 79.00
87045	Main Tank Check Valve	\$ 33.00
87175	Membrane Dryer Assembly	\$ 833.00
60013	Moisture Monitor Assy	\$ 23.00
85960	Motor Replacement (See Exchange Parts)	\$ 3,638.00
86019	Nipple (to Tank), 1/2MNPT x Brass PConn	\$ 23.00
87036	Outlet Manifold	\$ 85.00
60011	Pressure Gauge	\$ 37.00
60017	Pressure Relief Valve, 125 PSI (Standard)	\$ 36.00
87176	Pressure Switch, 85 to 115 PSI	\$ 53.00
85887	Tank Outlet Shutoff Valve	\$ 23.00
87108	Drain Tube and Drip Tray Kit	\$ 42.00
Electrical Breakdown		
85833	Capacitor, 65 uF, 370 VAC (2 per AS21)	\$ 62.00
117061	Circuit Breaker, Low Voltage, 24 VAC, 1/2 Amp	\$ 16.00
85960	Compressor Head Replacement, 2 Cylinder, 115 VAC	\$ 3,638.00
89345	Contactor 40 Amp Line Voltage with 24 VAC Coil	\$ 159.00
87076-1	Cooler Assembly, 115V Fan	\$ 28.00
87151	Line Cord Assembly, 115 VAC	\$ 61.00
87156	Motor Harness Replacement Kit	\$ 45.00
56149	Power On/Off DPST Switch and Circuit Breaker Motor	\$ 89.00
87176	Pressure Switch, 85 to 115 PSI	\$ 53.00
51952	Sound Cover Fan, 115 VAC	\$ 65.00
85962-1	Sound Cover Kit	\$ 1,125.00
87126	Sound Cover Receptacle	\$ 50.00
86637	Transformer, 115 VAC Step-down to 24 VAC	\$ 55.00
Additional Items		
85492	Remote Air Intake (Optional, Not Included with AirStar)	\$ 240.00
80183	Vibration Pad (Optional Not Included with AirStar)	\$ 41.00



Air Compressors with Membrane Dryer



Air Compressors with Membrane Dryer



Part Number	Description	Price
AirStar AS22 Parts Breakdown		
89948	1/4" OD, Black Nylon Tubing	\$ 2.00
86726	10mm Flamex Tubing, Cooler to Membrane Dryer	\$ 20.00
86727	12mm Flamex Tubing, Membrane Dryer to Main Tank Check Valve	\$ 24.00
87160	230V Cooler Assembly	\$ 189.00
85617	5 Micron Filter	\$ 251.00
85911	Drain Valve Assy	\$ 30.00
85603	Filter Element; Air: for Norgren	\$ 80.00
86193	Filter Element Replacement Kit, 5 Micron, Wilkerson	\$ 38.46
87049	Flex Braided Hose	\$ 87.00
S85482	Flex Pressure 6' Rubber Hook-Up Hose	\$ 38.00
89831	Intake Filter, Compressor Head	\$ 50.00
85471	Leveling Foot	\$ 79.00
87045	Main Tank Check Valve	\$ 33.00
87175	Membrane Dryer Assembly	\$ 833.00
60013	Moisture Monitor Assy	\$ 23.00
85910	Motor Replacement (See Exchange Parts)	\$ 3,638.00
87036	Outlet Manifold	\$ 85.00
58038	Elbow, 3/8 MNPT x 3/8 FNPT, Tank Inlet	\$ 5.00
56121	Elbow, 1/4 MNPT x 1/4 Push, Pconn	\$ 5.00
86729	Elbow, 10mm x 10mm, Pconn	\$ 14.00
86747	Elbow, 10mm x 3/8", Pconn	\$ 21.00
86749	Elbow, 12mm x 3/8", Pconn	\$ 30.00
86019	Nipple (to Tank), Pconn	\$ 23.00
60011	Pressure Gauge	\$ 37.00
60017	Pressure Relief Valve, 125 PSI (Standard)	\$ 36.00
60510	Pressure Relief Valve, 140 PSI	\$ 54.00
87176	Pressure Switch, 85 to 115 PSI	\$ 53.00
85887	Tank Outlet Shutoff Valve	\$ 23.00
87108	Drain Tube and Drip Tray Kit	\$ 42.00
Electrical Breakdown		
85832	Capacitor, 55 uF, 370 VAC	\$ 53.00
117690	Circuit Breaker, Low Voltage, 24 VAC, 1/4 Amp	\$ 16.00
85910	Compressor Head Replacement, 2 Cylinder	\$ 3,638.00
89345	Contactor 40 Amp Line Voltage with 24 VAC Coil	\$ 159.00
87071-1	Cooler Assembly 230V Fan	\$ 28.00
87152	Line Cord Assembly, 208/230 VAC	\$ 66.00
87156	Motor Harness Replacement Kit	\$ 45.00
56148	Power On/Off DPST Switch and Circuit Breaker Motor	\$ 89.00
87176	Pressure Switch, 85 to 115 PSI	\$ 53.00
85132	Sound Cover Fan, 230 VAC	\$ 137.00
85962-2M	Sound Cover Kit	\$ 1,125.00
87125	Sound Cover Receptacle	\$ 50.00
86637	Transformer, 208/230 VAC Step-down to 24 VAC	\$ 55.00
Additional Items		
67002	Buck/Boost Transformer, 3.9 Kva	\$ 600.00
85492	Remote Air Intake (Optional, Not Included with AirStar)	\$ 240.00
80183	Vibration Pad (Optional Not Included with AirStar)	\$ 41.00

Air Compressors with Membrane Dryer



Part Number	Description	Price
AirStar AS30 Parts Breakdown		
89948	1/4" OD, Black Nylon Tubing	\$ 2.00
86726	10mm Flamex Tubing, Cooler to Membrane Dryer	\$ 20.00
86727	12mm Flamex Tubing, Membrane Dryer to Main Tank Check Valve	\$ 24.00
85617	5 Micron Filter	\$ 251.00
87160	230V Cooler Assembly	\$ 189.00
85911	Drain Valve Assy	\$ 30.00
58038	Elbow, 3/8 MNPT x 3/8 FNPT, Tank Inlet	\$ 5.00
56121	Elbow, 1/4 MNPT x 1/4 Push, Pconn	\$ 6.00
86729	Elbow, 10mm x 10mm, Pconn	\$ 14.00
86749	Elbow, 12mm x 3/8", Pconn	\$ 30.00
85603	Filter Element; Air: for Norgren	\$ 84.00
86193	Filter Element Replacement Kit, 5 Micron, Wilkerson	\$ 40.00
87041	Flex Braided Hose, 13.5" (Right Motor)	\$ 60.00
87042	Flex Braided Hose, 17.5" (Left Motor)	\$ 60.00
87047	Flex Braided Hose, 17" (To Cooler Assembly)	\$ 65.00
S85482	Flex Pressure 6' Rubber Hook-Up Hose	\$ 38.00
89831	Intake Filter, Compressor Head	\$ 53.00
85471	Leveling Foot	\$ 79.00
87045	Main Tank Check Valve	\$ 33.00
87175	Membrane Dryer Assembly	\$ 833.00
60013	Moisture Monitor Assy	\$ 23.00
85940	Motor Replacement (x2) (See Exchange Parts)	\$ 2,825.00
87036	Outlet Manifold	\$ 85.00
86019	Nipple (to Tank), Pconn	\$ 23.00
60011	Pressure Gauge	\$ 37.00
60017	Pressure Relief Valve, 125 PSI (Standard)	\$ 36.00
60510	Pressure Relief Valve, 140 PSI	\$ 54.00
87176	Pressure Switch, 85 to 115 PSI	\$ 53.00
85887	Tank Outlet Shutoff Valve	\$ 23.00
80048R	Tee, 1/4 FNPT x 1/4 FNPT, Pconn	\$ 6.00
87108	Drain Tube and Drip Tray Kit	\$ 42.00
Electrical Breakdown		
85831	Capacitor, 27.5 uF, 370 VAC (2 Per AS30)	\$ 27.00
117690	Circuit Breaker, Low Voltage, 24 VAC, 1/4 Amp (May 2006-Current)	\$ 16.00
85940	Compressor Head Motor Assembly, 1 Cylinder (2 Per AS30)	\$ 3,638.00
89345	Contactor 40 Amp Line Voltage with 24 VAC Coil	\$ 159.00
87071-1	Cooler Assembly 230V Fan	\$ 28.00
87152	Line Cord Assembly, 230 VAC	\$ 66.00
87157	Motor Harness (2 per AS30)	\$ 83.00
85838	Motor On/Off, 20 Amp DPST Switch (2 Per AS30)	\$ 26.00
86065	Power On/Off DPST Switch and Circuit Breaker Motor (May 2006-Current)	\$ 97.00
87176	Pressure Switch, 85 to 115 PSI	\$ 53.00
85132	Sound Cover Fan, 230 VAC	\$ 137.00
85963M	Sound Cover Kit	\$ 1,200.00
87125	Sound Cover Receptacle	\$ 50.00
86637	Transformer, 208/230 VAC Step-down to 24 VAC	\$ 55.00
Additional Items		
67002	Buck/Boost Transformer, 3.9 Kva	\$ 600.00
85492	Remote Air Intake (Optional, Not Included with AirStar)	\$ 240.00
80183	Vibration Pad (Optional Not Included with AirStar)	\$ 41.00



Air Compressors with Membrane Dryer



Air Compressors with Membrane Dryer

Part Number	Description	Price
AirStar AS50 Parts Breakdown		
89948	1/4" OD, Black Nylon Tubing	\$ 2.00
86726	10mm Flamex Tubing, Cooler to Membrane Dryer	\$ 20.00
86727	12mm Flamex Tubing, Membrane Dryer to Main Tank Check Valve	\$ 24.00
85617	5 Micron Filter	\$ 251.00
87160	230V Cooler Assembly	\$ 189.00
85911	Drain Valve Assy	\$ 30.00
56121	Elbow, 1/4 MNPT x 1/4 Push, Pconn	\$ 5.00
86729	Elbow, 10mm x 10mm, Pconn	\$ 14.00
86747	Elbow, 10mm x 3/8", Pconn	\$ 21.00
86749	Elbow, 12mm x 3/8", Pconn	\$ 30.00
58038	Elbow, 3/8 MNPT x 3/8 FNPT, Tank Inlet	\$ 5.00
85603	Filter Element; Air: for Norgren	\$ 84.00
86193	Filter Element Replacement Kit, 5 Micron, Wilkerson	\$ 40.00
87068	Flex Braided Hose, 19.5" (Right Motor)	\$ 59.00
87046	Flex Braided Hose, 20" (to Cooler Assembly)	\$ 69.00
87043	Flex Braided Hose, 23.5" (Left Motor)	\$ 74.00
S85482	Flex Pressure 6' Rubber Hook-Up Hose	\$ 38.00
89831	Intake Filter, Compressor Head	\$ 53.00
85471	Leveling Foot	\$ 79.00
87045	Main Tank Check Valve	\$ 33.00
87175	Membrane Dryer Assembly	\$ 833.00
60013	Moisture Monitor Assy	\$ 23.00
85910	Motor Replacement (x2) (See Exchange Parts)	\$ 3,638.00
86019	Nipple (to Tank), Pconn	\$ 23.00
87036	Outlet Manifold	\$ 85.00
60011	Pressure Gauge	\$ 37.00
60017	Pressure Relief Valve, 125 PSI (Standard)	\$ 36.00
60510	Pressure Relief Valve, 140 PSI	\$ 54.00
87176	Pressure Switch, 85 to 115 PSI	\$ 53.00
85887	Tank Outlet Shutoff Valve	\$ 23.00
80048R	Tee, 1/4 FNPT x 1/4 FNPT, Pconn	\$ 6.00
87108	Drain Tube and Drip Tray Kit	\$ 42.00
Electrical Breakdown		
85832	Capacitor, 55 uF, 370 VAC (2 Per AS50)	\$ 53.00
117690	Circuit Breaker, Low Voltage, 24 VAC, 1/4 Amp	\$ 16.00
85910	Compressor Head Motor Assembly, 2 Cylinder (2 Per AS50)	\$ 3,638.00
89345	Contactor 40 Amp Line Voltage with 24 VAC Coil	\$ 159.00
87160	230V Cooler Assembly	\$ 189.00
87155	Motor Harness (2 per AS50)	\$ 84.00
85838	Motor On/Off, 20 Amp DPST Switch (2 Per AS50)	\$ 26.00
86377	Power Connection Handy Box Assembly (4ft Long Pigtail)	\$ 63.00
86064	Power On/Off DPST Switch and Circuit Breaker Motor	\$ 119.00
87176	Pressure Switch, 85 to 115 PSI	\$ 53.00
85132	Sound Cover Fan, 230 VAC	\$ 137.00
89523M	Sound Cover Kit	\$ 1,500.00
87125	Sound Cover Receptacle	\$ 50.00
86637	Transformer, 208/230 VAC Step-down to 24 VAC	\$ 55.00
Additional Items		
67002	Buck/Boost Transformer, 3.9 Kva	\$ 600.00
85493	Remote Air Intake (Optional, Not Included with AirStar)	\$ 250.00
80183	Vibration Pad (Optional, Not Included with AirStar)	\$ 41.00

Air Compressors with Membrane Dryer




Air Compressors with Membrane Dryer

Part Number	Description	Price
AirStar AS70 Parts Breakdown		
89948	1/4" OD, Black Nylon Tubing	\$ 2.00
86726	10mm Flamex Tubing, Cooler to Membrane Dryer	\$ 20.00
87160	230V Cooler Assembly (x2)	\$ 189.00
85617	5 Micron Filter	\$ 251.00
87062	Braided Hose, 24"	\$ 66.00
85911	Drain Valve Assy	\$ 30.00
86747	Elbow (x5), 10mm x 3/8", Pconn	\$ 21.00
60083	Elbow, 1/4 FNPT x 1/4 FNPT, Pconn	\$ 6.00
56121	Elbow, 1/4 MNPT x 1/4 Push, Pconn	\$ 5.00
86729	Elbow, 10mm x 10mm, Pconn	\$ 14.00
86749	Elbow, 12mm x 3/8", Pconn	\$ 30.00
85603	Filter Element; Air: for Norgren	\$ 84.00
86193	Filter Element Replacement Kit, 5 Micron, Wilkerson	\$ 40.00
87042	Flex Braided Hose, 17.5" (x3)	\$ 60.00
87062	Flex Braided Hose, 24"	\$ 66.00
S85482	Flex Pressure 6' Rubber Hook-Up Hose	\$ 38.00
89831	Intake Filter, Compressor Head	\$ 53.00
85471	Leveling Foot	\$ 79.00
87045	Main Tank Check Valve	\$ 33.00
87175	Membrane Dryer Assembly (x2)	\$ 833.00
60013	Moisture Monitor Assy	\$ 23.00
85910	Motor Replacement (x3) (See Exchange Parts)	\$ 3,638.00
86019	Nipple (to Tank), Pconn	\$ 23.00
87036	Outlet Manifold	\$ 85.00
60011	Pressure Gauge	\$ 37.00
60017	Pressure Relief Valve, 125 PSI (Standard)	\$ 36.00
60510	Pressure Relief Valve, 140 PSI	\$ 54.00
87176	Pressure Switch, 85 to 115 PSI	\$ 53.00
85887	Tank Outlet Shutoff Valve	\$ 23.00
80048R	Tee (x3), 1/4 FNPT x 1/4 FNPT, Pconn	\$ 6.00
87108	Drain Tube and Drip Tray Kit	\$ 42.00
Electrical Breakdown		
85832	Capacitor, 55 uF, 370 VAC (3 Per AS70)	\$ 53.00
117690	Circuit Breaker, Low Voltage, 24 VAC, 1/4 Amp	\$ 16.00
85910	Compressor Head Replacement 2 Cylinder (3 Per AS70)	\$ 3,638.00
89345	Contactor 40 Amp Line Voltage with 24 VAC Coil	\$ 159.00
87071-1	Cooler Assembly 230V Fan	\$ 28.00
87155	Motor Harness (qty 1), Left	\$ 84.00
87157	Motor Harness (qty 2), Middle & Right	\$ 83.00
85838	Motor On/Off, 20 Amp DPST Switch (3 Per AS70)	\$ 26.00
86378	Power Connection Handy Box Assembly (4ft Long Pigtail)	\$ 63.00
86066	Power On/Off 22Amp DPST Switch and Circuit Breaker Motor	\$ 119.00
87176	Pressure Switch, 85 to 115 PSI	\$ 53.00
85132	Sound Cover Fan, 230 VAC	\$ 137.00
89574M	Sound Cover Kit	\$ 3,500.00
87125	Sound Cover Receptacle	\$ 50.00
86637	Transformer, 208/230 VAC Step-down to 24 VAC	\$ 55.00
Additional Items		
67000-1	Buck/Boost Transformer, 3.9 Kva	\$ 1,300.00
85494	Remote Air Intake (Optional, Not Included with AirStar)	\$ 330.00
80183	Vibration Pad (Optional, Not Included with AirStar)	\$ 41.00

Air Compressors with Membrane Dryer

AirStar Membrane Dryer Parts Breakdown



Air Compressors with Membrane Dryer

Part Number Description Price

Membrane Dryer Parts Breakdown

87174	Condensate Filter Bowl and O-Ring Kit	\$ 39.00
87367	Filter, Membrane Bottom	\$ 12.00
87366	Filter, Membrane Top	\$ 30.00
87364	Membrane, Top Cover	\$ 2.00
85529	O-Ring, Filter Bowl	\$ 2.00
87368	O-Ring, Top Cover	\$ 3.00

Maintenance Kits

87351	AS10 Membrane Maintenance Kit	\$ 75.00
87352	AS21 Membrane Maintenance Kit	\$ 120.00
87352	AS22 Membrane Maintenance Kit	\$ 120.00
87352	AS30 Membrane Maintenance Kit	\$ 120.00
87353	AS50 Membrane Maintenance Kit	\$ 220.00
87354	AS70 Membrane Maintenance Kit	\$ 325.00
87108	Drain Tube and Drip Tray Kit	\$ 42.00
85694	Tandem AirStar Kit	\$ 300.00

Sound Reducer Kits

86272	AS10 Intake and Sound Dampner Tube	\$ 182.00
86273	AS21/AS22/AS30 Intake and Sound Dampner Tube	\$ 286.00
86274	AS50 Intake and Sound Dampner Tube	\$ 390.00
86275	AS70 Intake and Sound Dampner Tube	\$ 494.00

Kits Include:


87367	Filter Membrane Bottom; Qty 1*	\$ 12.00
87366	Filter Membrane Top; Qty 1*	\$ 30.00
89831	Intake filters; Qty Model Specific (per number of filters)	\$ 53.00
87369	O-Ring Filter Bowl; Qty. 1*	\$ 2.00
87368	O-Ring Top Cover; Qty 1*	\$ 3.00

* AirStar 70 has a Qty of 2


Air Compressors with Desiccant Dryer

Air Compressors with Desiccant Dryer

AirStar 10 with Desiccant Dryer Parts Breakdown



AirStar 10 with Desiccant Dryer Electrical Breakdown



Part Number Description

AirStar AS10/AS10C Parts Breakdown

Part Number	Description	Price
89183	1/8" Tubing Replacement Kit (Includes new Regeneration Orifice and sleeves)	\$ 22.00
85617	5 Micron Outlet Air Filter	\$ 251.00
86193	5 Micron Filter Replacement	\$ 40.00
85890	Compressor Head Motor Assembly 1 Cylinder 115 VAC	\$ 2,825.00
85911	Drain Valve Assembly Main Tank	\$ 30.00
85373	Drying Chamber	\$ 462.00
85948	Flex Metal Motor Discharge Hose	\$ 210.00
S85482	Flex Pressure 6' Rubber Hook-up Hose	\$ 38.00
85408	Humidistat Assembly	\$ 515.00
89831	Intake Filter Replacement	\$ 53.00
85659	Main Tank Check Valve	\$ 53.00
85769	Main Tank Outlet Manifold Assembly (Does not include PN 85911+60013)	\$ 173.00
85887	Main Tank Outlet Shutoff Valve	\$ 23.00
60013	Moisture Monitor	\$ 23.00
60011	Pressure Gauge	\$ 37.00
85979	Pressure Switch	\$ 100.00
60017	Pressure Relief Valve 125 PSI **Standard	\$ 36.00
60510	Pressure Relief Valve 140 PSI **Standard	\$ 54.00
85957	Purge Control Valve	\$ 121.00
85252	Purge Tank Assembly	\$ 143.00
85411	Regeneration Orifice Kit	\$ 49.00
85834	Regeneration Solenoid Valve Assembly 24 VAC	\$ 173.00
85471	Rubber Leveling Feet Kit (x4)	\$ 79.00
85719	Unloader Purge Solenoid Valve 115 VAC	\$ 187.00
80055-1	Unloader Purge Solenoid Valve Muffler (replaces copper & foam tubing)	\$ 53.00



Electrical Breakdown

116088	24 VAC Regeneration Indicator LED (on models manufactured after May 2006)	\$ 6.00
85833	Capacitor 65 uF, 370 VAC	\$ 62.00
117061	Circuit Breaker Low Voltage 24 VAC 1/2 Amp (1996-April 2006)	\$ 16.00
117690	Circuit Breaker Low Voltage 24 VAC 1/4 Amp (May 2006-Current)	\$ 16.00
85890	Compressor Head Motor Assembly 115 VAC 1 Cylinder	\$ 2,825.00
89345	Contactor 40 Amp Line Voltage with 24 VAC Coil	\$ 159.00
85408	Humidistat Air Moisture Sensor 24 VAC Switch	\$ 515.00
86138	Line Cord Assembly 115 VAC	\$ 173.00
86372	Motor Harness Replacement Kit	\$ 100.00
85838	Motor On/Off 20 Amp DPST Switch (1996-April 2006)	\$ 26.00
56149	Power On/Off DPST Switch and Circuit Breaker Motor (May 2006-Current)	\$ 89.00
85979	Pressure Switch	\$ 100.00
85834	Regeneration Solenoid 24 VAC	\$ 173.00
51952	Sound Cover Fan 115 VAC	\$ 65.00
85961	Sound Cover Kit	\$ 950.00
89942	Sound Cover Receptacle	\$ 47.00
85275	Transformer 115 VAC Step-down to 24 VAC	\$ 147.00
85719	Unloader Purge Solenoid 115 VAC	\$ 187.00

Additional Items

85491	Remote Air Intake (Optional, Not Included with AirStar)	\$ 220.00
80183	Vibration Pad (Optional Not Included with AirStar)	\$ 41.00

Air Compressors with Desiccant Dryer



Air Compressors with Desiccant Dryer

Part Number	Description	Price
AirStar AS21/AS21C Parts Breakdown		
89183	1/8" Tubing Replacement Kit (Includes new Regeneration Orifice and sleeves)	\$ 22.00
85617	5 Micron Outlet Air Filter	\$ 251.00
85890	Compressor Head Motor Assembly 1 Cylinder 115 VAC	\$ 2,825.00
85911	Drain Valve Assembly Main Tank	\$ 30.00
85373	Drying Chamber	\$ 462.00
85948	Flex Metal Motor Discharge Hose	\$ 210.00
S85482	Flex Pressure 6' Rubber Hook-up Hose	\$ 38.00
85408	Humidistat Assembly	\$ 515.00
89831	Intake Filter Replacement	\$ 53.00
85659	Main Tank Check Valve	\$ 53.00
85769	Main Tank Outlet Manifold Assembly (Does not include PN 85911+60013)	\$ 173.00
85887	Main Tank Outlet Shutoff Valve	\$ 23.00
60013	Moisture Monitor	\$ 23.00
60011	Pressure Gauge	\$ 37.00
85979	Pressure Switch	\$ 100.00
60017	Pressure Relief Valve 125 PSI **Standard	\$ 36.00
60510	Pressure Relief Valve 140 PSI **Standard	\$ 54.00
85957	Purge Control Valve	\$ 121.00
85253	Purge Tank Assembly	\$ 415.00
85411	Regeneration Orifice Kit	\$ 49.00
85834	Regeneration Solenoid Valve Assembly 24 VAC	\$ 173.00
85471	Rubber Leveling Feet Kit (x4)	\$ 79.00
85719	Unloader Purge Solenoid Valve 115 VAC	\$ 187.00
80055-1	Unloader Purge Solenoid Valve Muffler (replaces copper & foam tubing)	\$ 53.00



Electrical Breakdown

89345	40 Amp Line Voltage Contactor with 24 VAC Coil	\$ 159.00
116088	24 VAC Regeneration Indicator LED (on models manufactured after May 2006)	\$ 6.00
85833	Capacitor 65 uF, 370 VAC	\$ 62.00
117061	Circuit Breaker Low Voltage 24 VAC 1/2 Amp (1996-April 2006)	\$ 16.00
117690	Circuit Breaker Low Voltage 24 VAC 1/4 Amp (May 2006-Current)	\$ 16.00
85960	Compressor Head Motor Assembly 2 Cylinder 115 VAC	\$ 3,638.00
85408	Humidistat Air Moisture Sensor 24 VAC Switch	\$ 515.00
86138	Line Cord Assembly 115 VAC	\$ 173.00
86372	Motor Harness Replacement Kit	\$ 100.00
85838	Motor On/Off 20 Amp DPST Switch (1996-April 2006)	\$ 26.00
56149	Power On/Off DPST Switch and Circuit Breaker Motor (May 2006-Current)	\$ 89.00
85979	Pressure Switch	\$ 100.00
85839	Regeneration Solenoid 24 VAC	\$ 173.00
51952	Sound Cover Fan 115 VAC	\$ 65.00
85962-1	Sound Cover Kit	\$ 1,050.00
89942	Sound Cover Receptacle	\$ 47.00
85275	Transformer 115 VAC Step-down to 24 VAC	\$ 147.00
85719	Unloader Purge Solenoid 115 VAC	\$ 187.00

Additional Items

85492	Remote Air Intake (Optional, Not Included with AirStar)	\$ 240.00
80183	Vibration Pad (Optional Not Included with AirStar)	\$ 41.00

Air Compressors with Desiccant Dryer



Air Compressors with Desiccant Dryer

Part Number	Description	Price
AirStar AS22/AS22C Parts Breakdown		
89183	1/8" Tubing Replacement Kit (Includes new Regeneration Orifice and sleeves)	\$ 22.00
85617	5 Micron Outlet Air Filter	\$ 251.00
85910	Compressor Head Motor Assembly 2 Cylinder 208/230 VAC	\$ 3,638.00
85911	Drain Valve Assembly Main Tank	\$ 30.00
85373	Drying Chamber	\$ 462.00
85948	Flex Metal Motor Discharge Hose	\$ 210.00
S85482	Flex Pressure 6' Rubber Hook-up Hose	\$ 38.00
85408	Humidistat Assembly	\$ 515.00
89831	Intake Filter Replacement	\$ 53.00
85659	Main Tank Check Valve	\$ 53.00
85769	Main Tank Outlet Manifold Assembly (Does not include PN 85911+60013)	\$ 173.00
85887	Main Tank Outlet Shutoff Valve	\$ 23.00
60013	Moisture Monitor	\$ 23.00
60011	Pressure Gauge	\$ 37.00
85979	Pressure Switch	\$ 100.00
60017	Pressure Relief Valve 125 PSI **Standard	\$ 36.00
60510	Pressure Relief Valve 140 PSI **Standard	\$ 54.00
85957	Purge Control Valve	\$ 121.00
85253	Purge Tank Assembly	\$ 415.00
85411	Regeneration Orifice Kit	\$ 49.00
85834	Regeneration Solenoid Valve Assembly 24 VAC	\$ 173.00
85471	Rubber Leveling Feet Kit (x4)	\$ 79.00
85944	Unloader Purge Solenoid Valve 115 VAC	\$ 207.00
80055-1	Unloader Purge Solenoid Valve Muffler (replaces copper & foam tubing)	\$ 53.00


Electrical Breakdown

116088	24 VAC Regeneration Indicator LED (on models manufactured after May 2006)	\$ 6.00
85832	Capacitor 65 uF, 370 VAC	\$ 62.00
117061	Circuit Breaker Low Voltage 24 VAC 1/2 Amp (1996-April 2006)	\$ 16.00
117690	Circuit Breaker Low Voltage 24 VAC 1/4 Amp (May 2006-Current)	\$ 16.00
85910	Compressor Head Motor Assembly 2 Cylinder 208/230 VAC	\$ 3,638.00
89345	Contactor 40 Amp Line Voltage with 24 VAC Coil	\$ 159.00
85408	Humidistat Air Moisture Sensor 24 VAC Switch	\$ 515.00
86177	Line Cord Assembly 208/230 VAC	\$ 103.00
86372	Motor Harness Replacement Kit	\$ 100.00
85838	Motor On/Off 20 Amp DPST Switch (1996-April 2006)	\$ 26.00
56148	Power On/Off DPST Switch and Circuit Breaker Motor (May 2006-Current)	\$ 89.00
85979	Pressure Switch	\$ 100.00
85839	Regeneration Solenoid 24 VAC	\$ 173.00
85132	Sound Cover Fan 230 VAC	\$ 137.00
85962-2	Sound Cover Kit	\$ 1,050.00
89942	Sound Cover Receptacle	\$ 47.00
85362	Transformer 208/230 VAC Step-down to 24 VAC	\$ 147.00
85944	Unloader Purge Solenoid 208/230 VAC	\$ 207.00

Additional Items

67002	Buck/Boost Transformer, 3.9 Kva	\$ 600.00
85492	Remote Air Intake (Optional, Not Included with AirStar)	\$ 240.00
80183	Vibration Pad (Optional Not Included with AirStar)	\$ 41.00


Air Compressors with Desiccant Dryer




Air Compressors with Desiccant Dryer

Part Number	Description	Price
AirStar AS30/AS30C Parts Breakdown		
89183	1/8" Tubing Replacement Kit (Includes new Regeneration Orifice and sleeves)	\$ 22.00
85617	5 Micron Outlet Air Filter	\$ 251.00
85940	Compressor Head Motor Assembly 1 Cylinder 208/230 VAC	\$ 2,825.00
85911	Drain Valve Assembly Main Tank	\$ 30.00
85373	Drying Chamber	\$ 462.00
85682	Flex Metal Motor Discharge Hose (one each hose)	\$ 72.00
s85482	Flex Pressure 6' Rubber Hook-up Hose	\$ 38.00
85408	Humidistat Assembly	\$ 515.00
89831	Intake Filter Replacement	\$ 53.00
85659	Main Tank Check Valve	\$ 53.00
85769	Main Tank Outlet Manifold Assembly (Does not include PN 85911+60013)	\$ 173.00
85887	Main Tank Outlet Shutoff Valve	\$ 23.00
60013	Moisture Monitor	\$ 23.00
60011	Pressure Gauge	\$ 37.00
85979	Pressure Switch	\$ 100.00
60017	Pressure Relief Valve 125 PSI **Standard	\$ 36.00
60510	Pressure Relief Valve 140 PSI **Standard	\$ 54.00
85957	Purge Control Valve	\$ 121.00
85253	Purge Tank Assembly	\$ 415.00
85411	Regeneration Orifice Kit	\$ 49.00
85834	Regeneration Solenoid Valve Assembly 24 VAC	\$ 173.00
85471	Rubber Leveling Feet Kit (x4)	\$ 79.00
85845	Selector / Shuttle Valve Assembly	\$ 142.00
64009-1	Selector / Shuttle Valve Ball	\$ 14.00
85944	Unloader Purge Solenoid Valve 208/230 VAC	\$ 207.00
80055-1	Unloader Purge Solenoid Valve Muffler (replaces copper & foam tubing)	\$ 53.00
Electrical Breakdown		
116088	24 VAC Regeneration Indicator LED (on models after May 2006)	\$ 6.00
85831	Capacitor 27.5 uF, 370 VAC (2 per AS 30)	\$ 27.00
117061	Circuit Breaker Low Voltage 24 VAC 1/2 Amp (1996-April 2006)	\$ 16.00
117690	Circuit Breaker Low Voltage 24 VAC 1/4 Amp (May 2006-Current)	\$ 16.00
85363	Circuit breaker Motor (1996-April 2006)	\$ 119.00
85940	Compressor Head Motor Assembly 1 Cylinder (2 per AS30)	\$ 2,825.00
89345	Contactor 40 Amp Line Voltage with 24 VAC Coil	\$ 159.00
85408	Humidistat Air Moisture Sensor 24 VAC Switch	\$ 515.00
86177	Line Cord Assembly 230 VAC	\$ 103.00
86263	Motor Harness Replacement Kit (Inc. both left & right motor harness)	\$ 56.00
85838	Motor On/Off 20 Amp DPST Switch (1996-April 2006)	\$ 26.00
86063	Power On/Off DPST Switch and Circuit Breaker Motor (May 2006-Current)	\$ 119.00
85979	Pressure Switch	\$ 100.00
85834	Regeneration Solenoid 24 VAC	\$ 173.00
85132	Sound Cover Fan 115 VAC	\$ 65.00
85963	Sound Cover Kit	\$ 1,166.00
89942	Sound Cover Receptacle	\$ 47.00
85362	Transformer 208/230 VAC Step-down to 24 VAC	\$ 147.00
85944	Unloader Purge Solenoid 208/230 VAC	\$ 207.00
Additional Items		
67002	Buck/Boost Transformer, 3.9 Kva	\$ 600.00
85492	Remote Air Intake (Optional, Not Included with AirStar)	\$ 240.00
80183	Vibration Pad (Optional Not Included with AirStar)	\$ 41.00

Air Compressors with Desiccant Dryer



AirStar 50 with Desiccant Dryer Electrical Breakdown



Air Compressors with Desiccant Dryer

Part Number	Description	Price
AirStar AS50/AS50C Parts Breakdown		
89185	1/8" Tubing Replacement Kit (Includes new Regeneration Orifice and sleeves)	\$ 42.00
85617	5 Micron Outlet Air Filter	\$ 251.00
85910	Compressor Head Motor Assembly 2 Cylinder 208/230 VAC	\$ 3,638.00
85911	Drain Valve Assembly Main Tank	\$ 30.00
85373	Drying Chamber	\$ 462.00
85949	Flex Metal Motor Discharge Hose (one each motor)	\$ 158.00
S85482	Flex Pressure 6' Rubber Hook-up Hose	\$ 38.00
85408	Humidistat Assembly	\$ 515.00
89831	Intake Filter Replacement	\$ 53.00
85659	Main Tank Check Valve	\$ 53.00
85769	Main Tank Outlet Manifold Assembly (Does not include PN 85911+60013)	\$ 173.00
85887	Main Tank Outlet Shutoff Valve	\$ 23.00
60013	Moisture Monitor	\$ 23.00
60011	Pressure Gauge	\$ 37.00
85979	Pressure Switch	\$ 100.00
60017	Pressure Relief Valve 125 PSI **Standard	\$ 36.00
60510	Pressure Relief Valve 140 PSI **Standard	\$ 54.00
85958	Purge Control Valve	\$ 131.00
85265	Purge Tank Assembly	\$ 276.00
85412	Regeneration Orifice Kit	\$ 44.00
85839	Regeneration Solenoid Valve Assembly 24 VAC	\$ 173.00
85471	Rubber Leveling Feet Kit (x4)	\$ 79.00
85845	Selector / Shuttle Valve Assembly	\$ 142.00
64009-1	Selector / Shuttle Valve Ball	\$ 14.00
85946	Unloader Purge Solenoid Valve 208/230 VAC	\$ 207.00
80055-1	Unloader Purge Solenoid Valve Muffler (replaces copper & foam tubing)	\$ 53.00


Electrical Breakdown

116088	24 VAC Regeneration Indicator LED (on models manufactured after May 2006)	\$ 6.00
85832	Capacitor 55 uF, 370 VAC (2 per AS50)	\$ 53.00
117061	Circuit Breaker Low Voltage 24 VAC 1/2 Amp (1996-April 2006)	\$ 16.00
117690	Circuit Breaker Low Voltage 24 VAC 1/4 Amp (May 2006-Current)	\$ 16.00
85910	Compressor Head Motor Assembly 2 Cylinder (2 per AS50)	\$ 3,638.00
89345	Contactor 40 Amp Line Voltage with 24 VAC Coil	\$ 159.00
86377	Handy Box Assembly (4ft pigtail)	\$ 63.00
85408	Humidistat Air Moisture Sensor 24 VAC Switch	\$ 515.00
86265	Motor Harness Replacement Kit	\$ 56.00
85838	Motor On/Off 20 Amp DPST Switch (1996-April 2006)	\$ 26.00
86064	Power On/Off DPST Switch and Circuit Breaker Motor (May 2006-Current)	\$ 119.00
85979	Pressure Switch	\$ 100.00
85839	Regeneration Solenoid 24 VAC	\$ 173.00
85132	Sound Cover Fan 230 VAC	\$ 137.00
85965	Sound Cover Kit	\$ 1,260.00
89942	Sound Cover Receptacle	\$ 47.00
85362	Transformer 208/230 VAC Step-down to 24 VAC	\$ 147.00
85946	Unloader Purge Solenoid 208/230 VAC	\$ 207.00


Additional Items

67002	Buck/Boost Transformer, 3.9 Kva	\$ 600.00
85493	Remote Air Intake (Optional, Not Included with AirStar)	\$ 250.00
80183	Vibration Pad (Optional Not Included with AirStar)	\$ 41.00

Air Compressors with Desiccant Dryer




AirStar 70 with Desiccant Dryer Electrical Breakdown



Air Compressors with Desiccant Dryer

Part Number	Description	Price
AirStar AS70 Parts Breakdown		
89185	1/8" Tubing Replacement Kit (Includes new Regeneration Orifice and sleeves)	\$ 42.00
85617	5 Micron Outlet Air Filter	\$ 251.00
85910	Compressor Head Motor Assembly 1 Cylinder 115 VAC	\$ 3,638.00
85911	Drain Valve Assembly Main Tank	\$ 30.00
85374	Drying Chamber	\$ 473.00
85948	Flex Metal Motor Discharge Hose	\$ 210.00
s85482	Flex Pressure 6' Rubber Hook-up Hose	\$ 38.00
85408	Humidistat Assembly	\$ 515.00
89831	Intake Filter Replacement	\$ 53.00
85692	Main Tank Check Valve	\$ 65.00
85769	Main Tank Outlet Manifold Assembly (Does not include PN 85911+60013)	\$ 173.00
85887	Main Tank Outlet Shutoff Valve	\$ 23.00
85698	Metal Braided Flex Right	\$ 81.00
85682	Metal Braided Flex Center	\$ 72.00
85696	Metal Braided Flex Left	\$ 90.00
60013	Moisture Monitor	\$ 23.00
60011	Pressure Gauge	\$ 37.00
85979	Pressure Switch	\$ 100.00
60017	Pressure Relief Valve 125 PSI **Standard	\$ 36.00
60510	Pressure Relief Valve 140 PSI **Standard	\$ 54.00
85958	Purge Control Valve	\$ 131.00
85257	Purge Tank Assembly	\$ 306.00
85412	Regeneration Orifice Kit	\$ 44.00
85839	Regeneration Solenoid Valve Assembly 24 VAC	\$ 173.00
85471	Rubber Leveling Feet Kit (x4)	\$ 79.00
85846	Selector / Shuttle Valve Assembly	\$ 210.00
64009-1	Selector / Shuttle Valve Ball	\$ 14.00
86383	Unloader Purge Solenoid Valve 115 VAC	\$ 207.00
80055-1	Unloader Purge Solenoid Valve Muffler (replaces copper & foam tubing)	\$ 53.00
Electrical Breakdown		
116088	24 VAC Regeneration Indicator LED (on models manufactured after May 2006)	\$ 6.00
85832	Capacitor 65 uF, 370 VAC	\$ 53.00
117061	Circuit Breaker Low Voltage 24 VAC 1/2 Amp (1996-April 2006)	\$ 16.00
117690	Circuit Breaker Low Voltage 24 VAC 1/4 Amp (May 2006-Current)	\$ 16.00
85367	Circuit Breaker Motor (1996 - April 2006)	\$ 243.00
85910	Compressor Head Motor Assembly 2 Cylinder (3 per AS70)	\$ 3,638.00
89345	Contactor 40 Amp Line Voltage with 24 VAC Coil	\$ 159.00
85408	Humidistat Air Moisture Sensor 24 VAC Switch	\$ 515.00
86267	Motor Harness Replacement Kit	\$ 83.00
85838	Motor On/Off 20 Amp DPST Switch (1996-April 2006)	\$ 26.00
86378	Power Connection Handy Box Assembly (4ft long pigtail)	\$ 63.00
86066	Power On/Off DPST Switch and Circuit Breaker Motor (May 2006-Current)	\$ 119.00
85979	Pressure Switch	\$ 100.00
85839	Regeneration Solenoid 24 VAC	\$ 173.00
85132	Sound Cover Fan 230 VAC	\$ 137.00
89574	Sound Cover Kit	\$ 3,255.00
89942	Sound Cover Receptacle	\$ 47.00
85362	Transformer 208/230 VAC Step-down to 24 VAC	\$ 147.00
86383	Unloader Purge Solenoid 208/230 VAC	\$ 207.00
Additional Items		
67000-1	Buck/Boost Transformer	\$ 1,300.00
85494	Remote Air Intake (Optional, Not Included with AirStar)	\$ 330.00
80183	Vibration Pad (Optional Not Included with AirStar)	\$ 41.00

Compressor Head Assembly 1 and 2 Cylinders Parts Breakdown



Part Number Description


Compressor Head Assembly 1 and 2 Cylinders

Part Number	Description	Price
85890	Compressor Head Assembly 1 Cylinder 115 VAC AS-10	\$2,825.00
85940	Compressor Head Assembly 1 Cylinder 208/230 VAC AS-30 (Special order)	\$2,825.00
85960	Compressor Head Assembly 2 Cylinders 115 VAC AS-21(Special order)	\$3,638.00
85910	Compressor Head Assembly 2 Cylinders 208/230 VAC AS-22	\$3,638.00
85610	Cylinder & Head Replacement Kit	\$ 277.00
85947	Cylinder & Head Seal and Valve Kit	\$ 95.00
89887	Eccentric & Piston Assembly 1 Cylinder AS-10, AS-30/	\$ 515.00
89893	Eccentric & Piston Assembly 2 Cylinders AS-21,22,50 & 70	\$1,066.00
89802	Fan Crankcase Cooling Replacement Kit	\$ 46.00
89841	Fan Crankcase Grill Guard Square	\$ 40.00
89549	Fan Motor Replacement Kit Side Hole	\$ 25.00
85601	Gray Valve Plate Assembly with Flapper	\$ 40.00
85525	Gray Valve Plate Gasket	\$ 6.00
89831	Intake Filter Replacement (1 ea cylinder)	\$ 53.00
89907	M6 Bolts (qty 3)	\$ 2.00
89839	Piston Compression Ring	\$ 38.00
89838	Piston Expander Ring	\$ 4.00
30897	Screw M6 x 12 Phillips Pan Head for Fan Crankcase Grill Guard Square (4 ea)	\$ 6.00
30010	Screw T-30 Torx Top Plate (3 ea)	\$ 6.00
86412	Shock Mount Rubber Front Crankcase End (1 ea)	\$ 11.00
86425	Shock Mount Rubber Rear Motor End (2 ea)	\$ 23.00

Note: Compressor Head Assembly 1 and 2 Cylinders are relevant to both Membrane and Desiccant Models of AirStar.

Air Compressors

Part Number	Description	Price
Contactors		
67000-1	Buck/Boost Kit, 7.8 KVA, 208/230	\$1,300.00
53180	Circuit Protector: 3 Amp, New Style - Blue Box	\$ 26.00
55407	Contactor, 20a For 57077, 57078 Turbo	\$ 72.00
53183	Contactor, 40 Amp New Style - Blue Box	\$ 272.00
53090	Contactor, 40a For #53082 Old Style-Brown Box	\$ 564.00
53080	Contactor;115v,Under 1-1/2 HP	\$ 729.00
53081	Contactor;230v,Under 2-1/2 HP	\$ 729.00
53082	Contactor;230v,Under 5 HP	\$ 729.00
89345	Contactor Kit,40 Amp,230 Volt	\$ 159.00
53082M	Contactor; Tandem Assy,40 Amp	\$ 709.00
90219	Fuse Holder Old Style - Brown Box	\$ 9.00
40719	Fuse, 2 Amp (Pkg Of 5) Old Style - Brown Box	\$ 11.00
53087	Transformer 115v/24v Old Style - Brown Box	\$ 294.00
53088	Transformer 230v/24v Old Style - Brown Box	\$ 291.00
53182	Transformer 24v,115/230 Primary New Style-Blue Bx	\$ 120.00
67002	Transformer Kit, 3.9 KVA, 230/208	\$ 600.00
Accessories For Air Systems		
85471	4 Leveling Feet For Compressor	\$ 79.00
85480	Flex Electrical Hook-Up Cable	\$ 64.00
S85482	Flex Pressure Hook-Up Hose	\$ 38.00
60018	Line Cord w/ Plug, 115v, 14/3, Sjo	\$ 34.00
64018	Line Cord w/ Plug, 230v, 14/3	\$ 61.00
85483	Nipple, Hex 1/2 Npt	\$ 13.00
80183	Pad,Vibration,4 X 4 X 1	\$ 41.00
56600	Utility Stack Rack Mini: (AS30,STS-5,VS50/50H & smaller - any config)	\$ 875.00
56650	Utility Stack Rack XL: (AS 50/70,STS-6/10,VS80/80H - any config)	\$1,095.00



Part Number	Description	Price
Mojave Accessories		
MMC	Mojave Master Controller	\$1,890.00
MT10	Mojave 10 Gallon Tank	\$2,490.00
CT20	Mojave Dual Continuum Tank	\$3,390.00
H5229	MT10 to CT20 Conversion Kit	\$2,490.00
H5231	Mojave Powered Auto Drain Pump	\$1,000.00
Replacement Parts		
H5213	V3 Replacement Pump Assembly	\$5,990.00
H5215	V5 Replacement Pump Assembly	\$7,990.00
H5223	V3 Replacement Power Module	\$ 690.00
H5225	V5 Replacement Power Module	\$ 690.00
H5216	V3/V5 Replacement Inlet Filter and Manifold	\$ 200.00
H5217	V3/V5 Replacement Bowl, Screen and Gasket for Inlet Filter	\$ 50.00
H5218	M10/CT20 Replacement Drain Check Valve	\$ 60.00
H5219	M10/CT20 Replacement Surge Suppressor	\$ 60.00
H5220	M10/CT20 Replacement Washout Solenoid Assembly	\$ 220.00
H5221	M10/CT20 Replacement Float Switch Assembly	\$ 60.00
H5222	M10/CT20 Replacement Tank outlet and Gasket	\$ 60.00
H5224	M10/CT20 Top Tank Replacement Washout Nozzle Assembly	\$ 60.00
H5226	CT20 Bottom Tank Replacement Washout Nozzle Assembly	\$ 70.00
H5227	MT10/CT20 Bottom Tank Drain Assembly	\$ 110.00
H5228	CT20 Replacement Clear Valve	\$ 70.00

Part Number	Description	Price
STS-3 Dry Vacuum		
54538	Adapter Plate, Vertical onto Horizontal Stacking	\$ 236.00
67002	Buck/Boost Transformer (+/- 10%)	\$ 600.00
54323	Capacitor 50mfd, 480vac 2.5" Dia x 3.8lg, ;.250 Tab	\$ 58.00
117690	Circuit Breaker 1/4 Amp (Low Voltage Circuit)	\$ 16.00
54247	Circuit Breaker, On/Off Switch, 15amp Dpst	\$ 144.00
80178-4P	Clamp, Capacitor 2-1/2" Dia	\$ 3.00
89324	Clamp, Hose 1-9/16" To 2-1/2" Dia	\$ 3.00
55407	Contactor 230 VAC, 24 VAC Coil STS-3 & STS-5	\$ 72.00
54445	Electrical Box Assy; STS-3	\$ 457.00
54440	Heat Exchanger	\$ 646.00
56057	Hose, Exhaust Vent 2-1/4" D X 17"L, Black Flex Rubber	\$ 40.00
54327	Inlet Plumbing Fittings STS-3 Pump	\$ 70.00
43106	Leveling Foot (4 Per Unit)	\$ 3.00
54383	Line Cord Assy 230v	\$ 65.00
54160	Pump Accessory Pack	\$ 58.00
30423	Screw;8-32 X 3/8,Phil Pan Hd,18-8 Ss (Motor To Shroud)	\$ 1.00
54442	Shipping Packaging Material	\$ 173.00
54435	Shroud;STS-3 & 5:Welded Assy	\$ 131.00
54064	STS-3 Chassis Assembly (Without Pump & Electric Box)	\$ 632.00
54400	Tandem Connection Accessory Kit	\$ 250.00
54623	Transformer 230vac To 24 VAC	\$ 89.00
54538	Vertical onto Horizontal Stacking Kit	\$ 236.00
55907	Vibration Shock Mount (4 Per Unit)	\$ 3.00
30052	Washer;#8,Flat,18-8 Ss (Motor To Shroud)	\$ 1.00
30263	Washer;#8,Lock,Split,18-8 Ss (Motor To Shroud)	\$ 1.00
30958	Washer;Flat,1/4",18-8 Ss (Motor To Chassis)	\$ 1.00
54068	Wiring Harness Assembly (Electrical Box To Pump)	\$ 23.00
STS-5 Dry Vacuum		
67002	Buck/Boost Transformer (+/- 10% Change)	\$ 600.00
54323	Capacitor 50mfd 440vac 2.5" Dia Sn259 & Above	\$ 58.00
54333	Capacitor 60 Mfd 440vac 2.5" Dia SN 258 & Below	\$ 58.00
54065	Chassis Assembly (Without Pump & Electrical Box)	\$ 632.00
117690	Circuit Breaker 1/4amp (Low Voltage Circuit) Changed, Was 1/2amp	\$ 16.00
54247	Circuit Breaker On/Off Switch 15amp Dpst	\$ 144.00
80178-4P	Clamp, Capacitor 2-1/2" Dia	\$ 3.00
89324	Clamp, Hose 1-9/16"- 2-1/2" Dia	\$ 3.00
55407	Contactor 230vac, 24vac Coil, STS-3 & STS-5	\$ 72.00
54250	Electrical Box Assy;STS-5 (Serial No. 259 & Above)	\$ 450.00
54310	Heat Exchanger	\$ 529.00
56057	Hose, Exhaust Vent 2 1/4" Dia X 17"L, Black Flex Rubber	\$ 40.00
54184	Hose; Vent 2-1/4"ID x 36"L,1"W Cuff: Rubber Compound	\$ 58.00
54239	Inlet Plumbing Fittings (Motor Side Of Flexible Connector)	\$ 84.00
54359	Inlet Plumbing Fittings, Tank Connect. Side Of Flexible Connector	\$ 84.00
43106	Leveling Foot (4 Per Unit)	\$ 3.00

Part Number	Description	Price
STS-15 Dry Vacuum		
54383	Line Cord Assy, 230v	\$ 65.00
49318	Pconn;Adapter:1/4 FNPT x 1/8 FNPT, Brass	\$ 3.00
57253	Pconn;Adapt,1-1/2" MNPT x 1-1/2" Barb: New Age #561-3944	\$ 15.00
55219	Pipe Support;1-1/2",Polypropylene, Black	\$ 13.00
54150	Pump Accessory Pack	\$ 56.00
30863t	Screw;6-32 X 3/8,Type F, Phil PN Hd, Electrical Box Cover	\$ 1.00
30453	Screw;8-32 X 1 1/4,Phil PN Hd,18-8 Ss (Elect. Box To Chassis)	\$ 1.00
30423	Screw;8-32 X 3/8,Phil Pan Hd,18-8 Ss (Elect. Box/Motor To Shroud)	\$ 1.00
54355	Shipping Packaging Material	\$ 194.00
54435	Shroud;STS-3 & 5:Welded Assy	\$ 131.00
54400	Tandem Connection Accessory Kit	\$ 250.00
54375	Transformer 230vac To 24vac (STS-5 & STS-15)	\$ 100.00
55907	Vibration Shock Mount (4 Per Unit)	\$ 3.00
30052	Washer;#8,Flat,18-8 Ss (Electrical Box/Motor To Shroud)	\$ 1.00
30263	Washer;#8,Lock,Split,18-8 Ss (Electrical Box/Motor To Shroud)	\$ 1.00
30958	Washer;Flat,1/4",18-8 Stainless Steel, Vacuum Connect Support	\$ 1.00
54071	Wiring Harness Assembly (Electrical Box To Pump)	\$ 17.00
CAS 4-Gallon		
57253	Adapter, 1-1/2" MNPT X 1-1/2" Barb Hose, PVC Black	\$ 15.00
54438	Bracket, Separator Tank Hold Down Support (2 Per Unit)	\$ 8.00
77231	Cap, Frame Tubing 1" Square (4 Per Unit)	\$ 1.00
55711	Cap, Vacuum Relief Valve	\$ 6.00
54076	CAS Leg Stand Replacement Kit (31"Long, Incl. Tube Cap Pn 77231)	\$ 26.00
54139	Check Valve Swing 1-1/2" PVC	\$ 55.00
57169	Clamp, Hose 1.31"-2.25" Dia (For 1-1/2" Black Flex Tubing)	\$ 4.00

Part Number	Description	Price
54075	Cover Assembly Replacement Kit For 4,8 & 20 Gal CAS (Stainless)	\$ 630.00
54074	Cover Assy Hardware Replacement Kit For 4/8/20 Gallon CAS	\$ 107.00
54660	Drain Assembly (Incl. 1-1/2" Check Valve PN 54139)	\$ 69.00
54148	Drip Shield Assembly Small	\$ 231.00
57234	Elbow 90 Degree, Street 1-1/2" MNPT X 1-1/2" FNPT PVC	\$ 6.00
54619	Gasket, CAS Cover	\$ 54.00
55616	Gauge, Vacuum 30"Hg 1/8"FNPT	\$ 72.00
54184	Hose; Vent 2-1/4"ID x 36"L,1"W Cuff: Rubber Compound	\$ 58.00
43106	Leveling Foot (4 Per Unit)	\$ 3.00
30989	Leveling Foot Standoff 5/16" Id X 5/8"L (4 Per Unit)	\$ 12.00
54038	Pad Neoprene (For Tank Hold Down Bracket)	\$ 8.00
54457	Shipping Packaging Material	\$ 70.00
54120	Tank Accessory Pack (4/8 Gallon CAS)	\$ 123.00
54281	Tank Support, Frame Front/Back (2 Per Unit)	\$ 56.00
54282	Tank Support, Frame Left/Right (2 Per Unit)	\$ 42.00
54061	Tank Washout Port Connection Kit	\$ 29.00
55865	Vacuum Relief Valve	\$ 158.00
CAS 8-Gallon New Style		
57253	Adapter 1-1/2" MNPT X 1-1/2" Barb Hose, PVC Black	\$ 15.00
54438	Bracket, Separator Tank Hold Down Support (2 Per Unit)	\$ 8.00
77231	Cap, Frame Tubing 1" Square (4 Per Unit)	\$ 1.00
55711	Cap, Vacuum Relief Valve	\$ 6.00
54076	CAS Leg Stand Replacement Kit (31" Long, Incl. Tube Cap Pn77231)	\$ 26.00
54139	Check Valve Swing 1-1/2" PVC	\$ 55.00
57169	Clamp, Hose 1.31" - 2.25" Dia (For 1-1/2"Black Flex Tubing)	\$ 4.00
54075	Cover Assembly Replacement Kit (For 4/8/20 Gallon CAS, Stainless)	\$ 630.00
54074	Cover Assy Hardware Repl Kit For 4/8/20 Gal CAS;3 Per Cover	\$ 107.00
54660	Drain Assembly (Incl. 1-1/2" Check Valve PN 54139)	\$ 69.00
54025	Drip Shield Assembly Large	\$ 96.00
57234	Elbow 90 Degree, Street, 1-1/2" MNPT X 1-1/2" FNPT PVC	\$ 6.00
54619	Gasket, CAS Cover	\$ 54.00
55616	Gauge, Vacuum 30"Hg 1/8" FNPT	\$ 72.00
54118	Hose; 1-1/2" Dia PVC, Black With Grey Spiral (Order Per Ft 10'inc)	\$ 17.00
54184	Hose; Vent 2-1/4"ID x 36"L, 1"W Cuff: Rubber Compound	\$ 58.00
54133	Large Drip Shield Kit, CAS Tank	\$ 170.00
43106	Leveling Foot (4 Per Unit)	\$ 3.00
30989	Leveling Foot Standoff 5/16"Id X 5/8"L (4 Per Unit)	\$ 12.00
54038	Pad Neoprene (For Tank Hold Down Bracket)	\$ 8.00
54428	Shipping Packaging Material	\$ 91.00
54120	Tank Accessory Pack, 4/8 Gallon CAS	\$ 123.00
54281	Tank Support, Frame Front/Back (2 Per Unit)	\$ 56.00
54282	Tank Support, Frame Left/Right (2 Per Unit)	\$ 42.00
54061	Tank Washout Port Connection Kit	\$ 29.00
55865	Vacuum Relief Valve	\$ 158.00


Part Number	Description	Price
CAS 20-Gallon		
57253	Adapter, 1-1/2" MNPT X 1-1/2"Barb Hose, PVC Black	\$ 15.00
54438	Bracket, Separator Tank Hold Down Support (4 Per Unit)	\$ 8.00
77231	Cap, Frame Tubing 1" Square (4 Per Unit)	\$ 1.00
55711	Cap, Vacuum Relief Valve	\$ 6.00
54077	CAS Leg Stand Replacement Kit 50"L Incl Tube Cap, Lev. Foot, Clip	\$ 43.00
54139	Check Valve Swing 1-1/2" PVC	\$ 55.00
57169	Clamp, Hose 1.31" - 2.25" Diameter (For 1-1/2"Black Flex Tubing)	\$ 4.00
54075	Cover Assembly Replacement Kit For 4/8/20 Gallon CAS (Stainless)	\$ 630.00
54074	Cover Assy Hardware Replacement Kit For 4/8/20 Gal CAS (Stainless)	\$ 107.00
54660	Drain Assembly (Includes 1-1/2" Check Valve PN 54139)	\$ 69.00
57234	Elbow 90 Degree, Street, 1-1/2" MNPT X 1-1/2" FNPT PVC	\$ 6.00
54619	Gasket, CAS Cover	\$ 54.00
55616	Gauge, Vacuum 30" Hg 1/8" FNPT	\$ 72.00
54118	Hose;1-1/2" Dia PVC, Black With Grey Spiral(Order Per Ft 15' Inc)	\$ 17.00
54184	Hose; Vent 2-1/4"ID x 36"L,1"W Cuff: Rubber Compound	\$ 58.00
43106	Leveling Foot (4 Per Unit)	\$ 3.00
54038	Pad, Neoprene (For Tank Hold Down Bracket)	\$ 8.00
54457	Shipping Packaging Material	\$ 70.00
54120	Tank Accessory Pack 4 & 8 Gallon CAS	\$ 123.00
54281	Tank Support, Frame Front/Back (2 Per Unit)	\$ 56.00
54282	Tank Support, Frame Left/Right (2 Per Unit)	\$ 42.00
54061	Tank Washout Port Connection Kit	\$ 29.00
55865	Vacuum Relief Valve	\$ 158.00
A1029	Wall Bracket Mount Left Side Tank	\$ 121.00
A1028	Wall Bracket Mount Right Side Tank	\$ 70.00
Acadia Amalgam Separator		
A1384	90-degree Elbow with O-ring	\$ 25.00
A1381	90-degree Non-Rotating Elbow with O-ring Spacer & Screw	\$ 25.00
A1223	1 1/2 Inch Hose Barb Adaptor (Qty. 2)	\$ 20.00
54521	1 1/2 inch ID Clear with White Spiral Hose	\$ 16.00
A1210	Amalgam Separator Assembly	\$ 320.00
A1297	Amalgam Separator Chassis/Mounting Bracket	\$ 65.00
A1250	Filter only w/o Recycler	\$ 950.00
A1390	Filter & Recycler Kit	\$ 415.00
A1282	Inlet Filter Housing	\$ 25.00
A1359	Installation Accessory Kit	\$ 80.00
A1283	Large Particle Inlet Filter Strainer	\$ 25.00
A1370	Recycler Kit w/o Filter	\$ 115.00
A1237	Retaining Clip (Qty. 5)	\$ 20.00
A1280	Solid Collection Filter	\$ 150.00

Part Number	Description	Price
Guardian Collector		
A1006	Adapter, 1"Msip X 25mm Fslip, Grey (Goes Inside Of PVC Inlet)	\$ 86.00
A1010	Amalgam Collection Container Replacement Kit	\$1,100.00
A1033	Bag, Ziplock 18" X 20" (To Seal Used Box Prior Shipping)	\$ 1.00
A1098	Cap, 1"Id X 7/8" Long, Inlet/Outlet Cover Assy (2 Per Unit)	\$ 2.00
A1093	Certificate Of Disinfection	\$ 2.00
A1092	O-Ring Lubricant & Container	\$ 2.00
54131	O-Ring, Large 1.112"Id X 0.103"Wide (On Adapter PN A1006/A1007)	\$ 1.00
51571	O-Ring, Small 0.864id X 0.70 Wide (Top Lid O-Rings, 2 Per Unit)	\$ 1.00
A1134	Shipping Carton With Label (To Ship Used Box To Recycler)	\$ 33.00
4 Or 8 Gal CAS With Guardian		
57253	Adapter 1-1/2" MNPT X 1-1/2" Barb Hose PVC Black (4 Per Unit)	\$ 15.00
54438	Bracket, Separator Tank Hold Down Support (2 Per Unit)	\$ 8.00
77231	Cap, Frame Tubing 1" Square (4 Per Unit)	\$ 1.00
55711	Cap, Vacuum Relief Valve	\$ 6.00
54074	CAS Cover Replacement Hard Kit	\$ 107.00
41211	Clamp Drain Hose 9/16"-1 1/4"Dia (For 3/4" Clear PVC, 4 Per Unit)	\$ 5.00
57169	Clamp, Hose 1.31"-2.25" Diameter (For 1-1/2" Black Flex Tubing)	\$ 4.00
54141	Drain Fitting Assembly w/ Inlet For Amalgam Collector	\$ 61.00
54437	Drip Shield Small 4 Gallon CAS	\$ 247.00
54133	Large Drip Shield Kit, CAS Tank	\$ 170.00
57234	Elbow 90 Degree, 1-1/2" MBPT X 1-1/2" FNPT PVC	\$ 6.00
54619	Gasket, CAS Cover	\$ 54.00
55616	Gauge, Vacuum 30"Hg 1/8" FNPT	\$ 72.00
55941	Hose, 3/4" Flex PVC Clear (Order Per Foot, 6' Incl. With Unit)	\$ 4.00
54184	Hose; Vent 2-1/4"ID x 36"L,1"W Cuff: Rubber Compound	\$ 58.00
43106	Leveling Foot (4 Per Unit)	\$ 3.00
30989	Leveling Foot Standoff 5/16" Id X 5/8" L (4 Per Unit)	\$ 12.00
A1035	Outlet Fitting Assembly For Guardian Amalgam Collector	\$ 131.00
54038	Pad, Neoprene (For Tank Hold Down Bracket)	\$ 8.00
54428	Shipping Packaging Material 8 Gallon CAS	\$ 91.00
54281	Tank Support, Frame Front/Back (4 Per Unit)	\$ 56.00
54282	Tank Support, Frame Left/Right (4 Per Unit)	\$ 42.00
54061	Tank Washout Port Connection Kit	\$ 29.00
55865	Vacuum Relief Valve	\$ 158.00
STS Accessories		
54400	Accessory Kit; Dual STS	\$ 250.00
57169	Clamp; Hose, 1/31" - 2.25" Dia.	\$ 4.00
89324	Clamp; Hose, 1-9/16" - 2 1/2" Max	\$ 3.00
54580	CleanShield Washout System	\$ 800.00
54521	Hose; 1 1/2" ID, Clear w/ White Spiral: Hi-Tech	\$ 17.00
54118	Hose; 1 1/2" Dia, PVC, Grey-Black (Order per Foot)	\$ 40.00
56057	Hose, 2 1/4" Dia x 17" L, 1" W Cuff	\$ 40.00

Part Number	Description	Price
54133	Large Drip Shield Kit, CAS Tank	\$ 170.00
30049	Nut; 1/4-20, Hex, Plated	\$ 1.00
57253	PConn Adapt, 1 1/2", MNPT 1 1/2" B	\$ 15.00
54150	Pump Accessory Pack; Dry Vacuum	\$ 56.00
56299	Switch,6VDC,With Markers	\$ 22.00
54350	Tank Accessory Pack; Dry Vacuum	\$ 173.00
54538	Mounting Kit (Old style to new style STS pump)	\$ 236.00
56600	Utility Stack Rack Mini: (AS30,STS-5,VS50/50H & smaller - any config)	\$ 875.00
56650	Utility Stack Rack XL: (AS 50/70,STS-6/10,VS80/80H - any config)	\$1,095.00

VacStar 20

BREAKDOWN




Part Number	Description	Price
VacStar VS20 Wet Vacuum (in chart)		
117061	1/2 amp Circuit Breaker	\$ 16.00
56361	Base Plate & O-Ring	\$ 64.00
55683	Capacitor Kit	\$ 62.00
56056	Check Valve	\$ 20.00
55857	Contactor Field Replacement Kit	\$ 98.00
55333	Drip Cover w/ Spacer Kit	\$ 80.00
90219	Fuse Holder	\$ 9.00
55869	Fuse	\$ 19.00
56149	Main Circuit Breaker	\$ 89.00
30851T	Screw	\$ 1.00
55863	Shock Mount Replacement Kit	\$ 29.00
55020	Solenoid Water Valve Kit	\$ 173.00
55880	Solids Collector Kit	\$ 40.00
55867	Transformer	\$ 131.00
55525	Vacuum Breaker	\$ 56.00
55517	Vacuum Breaker Seal Kit	\$ 37.00
55616	Vacuum Gauge	\$ 72.00
55865	Vacuum Relief Valve	\$ 158.00
55709	Vacuum Relief Valve Filter Cap	\$ 53.00
55909	Water Inlet Strainer	\$ 74.00
55505	Water System Replacement Kit	\$ 404.00

Additional Parts (not in chart)

55321	Base Plate O-Ring	\$ 3.00
56105	Check Valve & Nipples Kit	\$ 37.00
55292	Hose Assy Exhaust Manifold	\$ 47.00
55269	Hose, Barb Fitting, 1" Barb X 3/4" Fgh	\$ 22.00
55941	Hose, PVC 3/4" I.D. Per Feet, Reinforced, Suction/Discharge	\$ 4.00
H-2	Hydromiser Kit	\$ 995.00
55078	Inline Filter Kit 3/4"	\$ 102.00
55801	Installation Kit	\$ 243.00
55626	Manifold Plug	\$ 3.00
55802	Pconn; Adapter, 1/4 MNPT X 3/8 Poly	\$ 7.00
55110	Pump, 1 HP, Replacement, (See Exchange Parts)	\$ 2,100.00
55877	Ryton Base, 1HP Pump Housing Replacement	\$ 84.00
55519	Solenoid Valve Water (Square), 115 V	\$ 140.00
55387	Solenoid Valve, Water (Only), Round	\$ 99.00
55880C	Solids Collector, Case of 12	\$ 184.00
55862	Solid Collector Complete With Black Housing	\$ 50.00
55629	Solids Collector Bowl	\$ 26.00
55624	Solids Collector Filter Screen 40 Mesh	\$ 9.00
55625	Solids Collector Gasket	\$ 4.00
55861	Terminal Strip, Sub Assy Replacement	\$ 25.00

VacStar 40


BREAKDOWN




Part Number	Description	Price
VacStar VS40 Wet Vacuum (in chart)		
117690	1/4 amp Circuit Breaker	\$ 16.00
56361	Base Plate & O-Ring	\$ 64.00
55684	Capacitor Kit	\$ 65.00
56056	Check Valve	\$ 20.00
55857	Contactor Field Replacement Kit	\$ 98.00
55333	Drip Cover w/Spacer Kit	\$ 80.00
55524	Flow Regulator	\$ 78.00
55869	Fuse	\$ 19.00
56149	Main Circuit Breaker	\$ 89.00
30851T	Screw	\$ 1.00
55863	Shock Mount Replacement Kit	\$ 29.00
55040	Solenoid Water Valve Kit	\$ 117.00
55880	Solids Collector Replacement Kit	\$ 40.00
55866	Transformer	\$ 131.00
55525	Vacuum Breaker	\$ 56.00
55517	Vacuum Breaker Seal Kit	\$ 37.00
55616	Vacuum Gauge	\$ 72.00
55865	Vacuum Relief Valve	\$ 158.00
55709	Vacuum Relief Valve Filter Cap	\$ 53.00
55909	Water Inlet Strainer	\$ 74.00
55507	Water System Replacement Kit	\$ 404.00

Additional Items (not in chart)

55321	Base Plate O-Ring	\$ 3.00
56105	Check Valve & Nipples Kit	\$ 37.00
55292	Hose Assy Exhaust Manifold	\$ 47.00
55269	Hose, Barb Fitting, 1" Barb X 3/4" Fgh	\$ 22.00
55941	Hose, PVC 3/4" I.D. Per Feet, Reinforced, Suction/Discharge	\$ 4.00
H-4	Hydromiser Kit	\$ 995.00
55078	Inline Filter Kit 3/4"	\$ 102.00
55801	Installation Kit	\$ 243.00
55626	Manifold Plug	\$ 3.00
55802	Pconn; Adapter, 1/4 MNPT X 3/8 Poly	\$ 7.00
55111	Pump, 2 HP, Replacement, (See Exchange Parts)	\$2,625.00
55878	Ryton Base, 2HP Pump Housing Replacement	\$ 100.00
55518	Solenoid Valve Water (Square), 220 V	\$ 140.00
55107	Solenoid Valve, Water (Only), Round	\$ 99.00
55880C	Solids Collector, Case of 12	\$ 184.00
55862	Solid Collector Complete With Black Housing	\$ 50.00
55629	Solids Collector Bowl	\$ 26.00
55624	Solids Collector Filter Screen 40 Mesh	\$ 9.00
55625	Solids Collector Gasket	\$ 4.00



Part Number	Description	Price
VacStar VS50/50H Wet Vacuum (in chart)		
117690	1/4 amp Circuit Breaker	\$ 16.00
56361	Base Plate & O-Ring Kit	\$ 64.00
55683	Capacitor Kit	\$ 62.00
56105	Check Valve & Nipples Kit	\$ 37.00
56056	Check Valve	\$ 20.00
55857	Contactor Field Replacement Kit	\$ 98.00
55333	Drip Cover w/Spacer Kit	\$ 80.00
56043	Exhaust Manifold (VacStar 50 only)	\$ 210.00
55523	Flow Regulator (VacStar 50 only)	\$ 78.00
55115	Inlet Manifold	\$ 431.00
55092	J-Bolts	\$ 4.00
56148	Main Circuit Breaker	\$ 89.00
30851T	Screw	\$ 1.00
55612	Solids Collector Complete	\$ 67.00
55060	Solenoid Water Valve Kit	\$ 289.00
53182	Transformer (56325 for Canada)	\$ 120.00
55525	Vacuum Breaker	\$ 56.00
55517	Vacuum Breaker Seal Kit	\$ 37.00
55616	Vacuum Gauge	\$ 72.00
55865	Vacuum Relief Valve	\$ 158.00
55859	Vibration Mount	\$ 9.00
55909	Water Inlet Strainer	\$ 74.00
Additional Parts (not in chart)		
55178	1 HP replacement for VS50 or VS50H	\$1,995.00
55321	Base Plate O-Ring	\$ 3.00
53180	Circuit Protector 3 Amp	\$ 26.00
55142	Discharge Hose Assy	\$ 79.00
55117	Eductor Assy	\$ 378.00
56044	Exhaust Manifold (VacStar 50H only)	\$ 210.00
55292	Hose Assy Exhaust Manifold	\$ 47.00
56057	Hose, 2 1/4" Dia X 17'L X 1" Wide Cuff Suction/Discharge	\$ 40.00
55269	Hose, Barb Fitting, 1" Barb X 3/4" Fgh	\$ 22.00
55948	Hose; PVC 1" I.D. Per Feet, Reinforced	\$ 6.00
55087	Hydromiser Wallmount Kit	\$ 140.00
55149	Installation Kit	\$ 231.00
57278	Pconn, Coupling, 1" (Flex)	\$ 15.00
55802	Pconn; Adapter, 1/4" MNPT X 3/8" Poly	\$ 7.00
55877	Ryton Base, 1HP Pump Housing Replacement	\$ 84.00
55127	Separator Assy Hydromiser	\$ 225.00
55116	Solenoid Valve, Water (Only), Round	\$ 99.00
55114	Solids Collector Bowl Twin Units	\$ 45.00
55079	Solids Collector Complete, Remote Option For Twin Units	\$ 147.00
55112	Solids Collector Filter Screen 40 Mesh Twin Units	\$ 12.00
55113	Solids Collector Gasket Twin Units	\$ 2.00
55094C	Solids Collector ,Twin Units, Case of 8	\$ 210.00
85838	Switch Toggle 20 Amp Dpst	\$ 26.00
117403	Terminal Block	\$ 6.00
53182	Transformer 230v To 24 VAC	\$ 120.00
55093	Vacuum Relief Valve Filter	\$ 61.00
56010	Water Flow Regulator (1.0 GPM) (VacStar 50H only)	\$ 65.00
55096	Water System Replacement Kit	\$ 213.00



Part Number	Description	Price
VacStar VS80/80H Wet Vacuum (in chart)		
117690	1/4 amp Circuit Breaker	\$ 16.00
56361	Base Plate & O-Ring Kit	\$ 64.00
55684	Capacitor Kit	\$ 65.00
56105	Check Valve and Nipples Kit	\$ 37.00
56056	Check Valve	\$ 20.00
55857	Contactor Field Replacement Kit	\$ 98.00
55333	Drip Cover w/ Spacer Kit	\$ 80.00
56044	Exhaust Manifold (VacStar 80H only)	\$ 210.00
55524	Water Flow Regulator	\$ 78.00
55104	Inlet Manifold	\$ 446.00
55092	J-Bolts	\$ 4.00
56149	Main Circuit Breaker	\$ 89.00
53182	Transformer (56325 for Canada)	\$ 120.00
30851T	Screw	\$ 1.00
55612	Solids Collector Complete	\$ 67.00
55080	Solenoid Water Valve Kit	\$ 289.00
55525	Vacuum Breaker	\$ 56.00
55517	Vacuum Breaker Seal Kit	\$ 37.00
55616	Vacuum Gauge	\$ 72.00
55865	Vacuum Relief Valve	\$ 158.00
55093	Vacuum Relief Valve Filter	\$ 61.00
55859	Vibration Mount	\$ 9.00
55909	Water Inlet Strainer	\$ 74.00
Additional Parts (not in chart)		
55179	2 HP replacement for VS80 or VS80H - Cartridge Only	\$ 2,415.00
55321	Base Plate O-Ring	\$ 3.00
53180	Circuit Protector 3 Amp	\$ 26.00
55142	Discharge Hose Assy	\$ 79.00
55089	Eductor Assy	\$ 378.00
56043	Exhaust Manifold (VacStar 80 only)	\$ 210.00
55203	Filter, Inline Replacement	\$ 67.00
55292	Hose Assy Exhaust Manifold	\$ 47.00
56057	Hose, 2 1/4" Dia X 17" L, 1" Wide Cuff	\$ 40.00
55269	Hose, Barb Fitting, 1" Barb X 3/4" Fgh	\$ 22.00
55948	Hose; PVC 1" I.D. Reinforced Per Feet	\$ 6.00
55087	Hydromiser Wallmount Kit (Optional)	\$ 140.00
55149	Installation Kit	\$ 231.00
57278	Pconn, Coupling, 1" (Flex)	\$ 15.00
55802	Pconn; Adapter, 1/4" MNPT X 3/8" Poly	\$ 7.00
55179	Replacement Pump Assy 2 HP (See Exchange Parts)	\$ 2,415.00
55127	Separator Assy Hydromiser	\$ 225.00
55878	Ryton Base, 2HP Pump Housing Replacement	\$ 100.00
55107	Solenoid Valve Water, Round Only	\$ 99.00
55114	Solids Collector Bowl Twin Units	\$ 45.00
55079	Solids Collector Complete, Remote Option For Twin Units	\$ 147.00
55112	Solids Collector Filter Screen 40 Mesh Twin Units	\$ 12.00
55113	Solids Collector Gasket Twin Units	\$ 2.00
55094	Solids Collector, twin unit 1 filter	\$ 58.00
55094C	Solids Collector , Twin Units, Case of 8	\$ 210.00
85838	Switch Toggle 20 Amp	\$ 26.00
53182	Transformer 230v To 24 VAC	\$ 120.00
56011	Water Flow Regulator (1.5 GPM)	\$ 65.00
55099	Water System Replacement Kit (Incl. Vacbreaker & Flow Restrictor)	\$ 226.00

Vacuums

Part Number	Description	Price
Air/Water Separator		
55540	Air Water Separator (For Use w/All VacStar Models)	\$ 600.00
56057	Hose & Cuff Assy (2" Exhaust)	\$ 40.00
55235	Manifold Assy H-2 & H-4	\$ 116.00
55270	Separator Assy	\$ 273.00
Hydromiser		
49833	Clamp; Hose, S.S.	\$ 7.00
55245	Eductor Assy H-2 & H-5	\$ 431.00
55255	Eductor Assy H-4 & H-8	\$ 383.00
55764	Handle;3-Way Ball Valve All Models	\$ 11.00
56057	Hose Assy/ (2" Exhaust) - All Models	\$ 40.00
55269	Hose Barb Assy; 1barb X 3/4 FGH	\$ 22.00
H-2	Hydromiser For VacStar 20	\$ 995.00
H-4	Hydromiser For VacStar 40	\$ 995.00
56041	Hydromiser For VacStar 50	\$ 800.00
56042	Hydromiser For VacStar 80	\$ 800.00
55087	Hydromiser Wall Mount Kit For VacStar 50H & 80H	\$ 140.00
55235	Manifold Assy H-2 & H-4	\$ 116.00
55265	Manifold Assy H-5 & H-8	\$ 194.00
55780	Regulator Assy All Models, 60 PSI	\$ 269.00
55776-2	Regulator Gauge All Models	\$ 54.00
55776	Regulator;Pressure,60PSI All Models	\$ 156.00
55127	Separator Assy 50H, 80H	\$ 225.00
55225	Separator Assy H-2,H-4,H-5,H-8	\$ 278.00
55761	Valve, 3-Way Ball 1/8 FNPT All Models	\$ 158.00
Stack Racks		
56600	Utility Stack Rack Mini: (AS30,STS-5,VS50/50H & smaller - any config)	\$ 875.00
56650	Utility Stack Rack XL: (AS 50/70,STS-6/10,VS80/80H - any config)	\$1,095.00
Vacuum Equalizer		
56230	Remote Switch Panel	\$ 290.00
53202-1	Switch Assembly, Push Button, New Style	\$ 54.00
56233	Switch, Control Panel, Old Style Rocker	\$ 45.00
56200	Vacuum Equalizer	\$1,450.00
Accessories		
53240	Two Pump Water Lockout Relay Kit	\$ 95.00

Accessories

Part Number	Description	Price
3/4" Remote Water Control Valve		
53020-1	3/4" Remote Water Control Valve W/Filter;115v	\$ 1,000.00
53020	3/4" Remote Water Control Valve W/Filter;24v	\$ 1,071.00
53020T	3/4" Remote Water Control Valve W/Filter;24v (inc. Transform. Assy)	\$ 1,460.00
53021	3/4" Remote Water Control Valve W/O Filter;24v	\$ 860.00
53021-1	3/4"Remote Water Control Valve W/O Filter;115v	\$ 850.00
53021T	3/4" Remote Water Control Valve W/Filter;24v (inc. Transform Assy.)	\$ 1,120.00
53058	Housing, Filter Assy (Ametek); Upper Housing & Lower Bowl	\$ 106.00
53019	O-Ring For Filter Assy-Old Style (Clear Bowl), Cuno	\$ 16.00
53152	O-Ring Kit-For Filters w/On-Off Turn Valve (Cuno)	\$ 90.00
53106	O-Ring, New Style Ametek	\$ 9.00
53035	Replacement Filter Element (Pkg Of 2) 5 Micron	\$ 76.00
53031-2	Solenoid Coil,115v,Fits Metal Housing Only(Black Coil)(Old Style)	\$ 218.00
53031-5	Solenoid Coil,115v,Fits Plastic Housing Only(Green Coil)(New Style)	\$ 218.00
53031-3	Solenoid Coil,24v,Fits Metal Housing Only (Black Coil)(Old Style)	\$ 218.00
53031-4	Solenoid Coil,24v,Fits Plastic Housing Only(Green Coil)(New Style)	\$ 310.00
53031-1	Solenoid Valve, 3/4 Npt, 115v, Asco	\$ 461.00
53031	Solenoid Valve, 3/4 Npt, 24v, Asco	\$ 461.00
53150	Transformer Assembly for 53020 & 53021	\$ 350.00
53025	Valve, Gate Brass, 3/4"Npt	\$ 34.00
53097	Water Filter Replacement Kit(Housing,Bowl,Nipple) & Filter	\$ 201.00
1" Remote Water Control Valves		
53173	1" Remote Water Control Valve W/O Filter;24v	\$ 1,141.00
53173T	1" Remote Water Control Valve W/O Filter;24v (inc. Transformer Assy.)	\$ 1,384.00
53171	1" Remote Water Control Valve W/Filter;115v	\$ 1,260.00
53170	1" Remote Water Control Valve W/Filter;24v	\$ 1,260.00
53170T	1" Remote Water Control Valve W/Filter;24v (inc. Transformer Assy.)	\$ 1,600.00
53174	1" Remote Water Control Valve W/O Filter;115v	\$ 1,100.00
53153	Cartridge, Water Filter Replacement, Pkg Of 1	\$ 62.00
53172	Filter, 1", Bowl	\$ 234.00
53172-2	Housing; Water Filter, 1", Bowl	\$ 137.00
53172-1	O-Ring For Filter Bowl (1")	\$ 24.00
53195	Solenoid Coil 115v (1")	\$ 194.00
53194	Solenoid Coil 24v (1")	\$ 222.00
53155	Solenoid Valve: 1" Npt 115v	\$ 673.00
53154	Solenoid Valve: 1" Npt 24v	\$ 621.00
53150	Transformer Assembly for 53170 & 53173	\$ 350.00
53156	Valve, Gate 1" Npt	\$ 84.00
Remote Control Panels		
53225-1	Cover Plate, 2-Gang, Four Switch Hole, Etched	\$ 34.00
53222-1	Cover Plate, 2-Gang, One Switch Hole, Etched	\$ 34.00
23224-1	Cover Plate, 2-Gang, Three Switch Hole, Etched	\$ 34.00
53223-1	Cover Plate, 2-Gang, Two Switch Hole, Etched	\$ 34.00

Accessories

Part Number	Description	Price
53202-1	Equalizer Push Button Switch	\$ 51.00
53050	Label, Moisture Monitor	\$ 4.00
60013	Moisture Monitor	\$ 23.00
53251	Remote Control Panel;24v Switches; Air/Vac (2 switch)	\$ 245.00
53250	Remote Control Panel;24v Switches; Air/Vac/Water (3 switch)	\$ 300.00
53111	Remote Control Panel;24v Switches; Air/Vac/Water (1 switch)	\$ 210.00
53133	Remote Control Panel;24v Switches; Air/Water/Vac 1/Vac2 (4 switch)	\$ 325.00
53232	Ring; Designation, Air	\$ 2.00
53232-2	Ring; Designation, Air-2	\$ 2.00
53232-3	Ring; Designation, Air-3	\$ 2.00
53232-4	Ring; Designation, Air-4	\$ 2.00
53231	Ring; Designation, Vacuum	\$ 2.00
53231-2	Ring; Designation, Vacuum-2	\$ 2.00
53231-3	Ring; Designation, Vacuum-3	\$ 2.00
53231-4	Ring; Designation, Vacuum-4	\$ 2.00
53233	Ring; Designation, Water	\$ 2.00
30692	Screw:6-32 X 1/2 Oval Head - Phillips Ss	\$ 1.00
53202-2	Switch Assy, 6V, CAS Washout Timer Control	\$ 75.00
53201-1	Switch Assy; Push Btn, PNL MT, Grn 24V LED	\$ 56.00
53048	Switch, Control Panel, 24v (Rocker)	\$ 22.00
56233	Switch, Control Panel, 6v Rocker	\$ 45.00

Film Processing

Part Number	Description	Price
A/T 2000XR Major Assemblies		
45200	A/T 2000XR Base Assembly	\$5,039.00
45800	A/T 2000XR Cover Assembly	\$4,603.00
45500	Developer Rack, AT2000XR (See Exchange Parts)	\$ 898.00
45505	Fixer Rack, AT2000XR (See Exchange Parts)	\$ 872.00
45100	Leveling Base Assembly	\$ 194.00
43106	Leveling Feet	\$ 3.00
45510	Wash/Dry Rack, AT2000XR (See Exchange Parts)	\$1,628.00
A/T 2000XR Cover Assembly PN 45800		
43811	Air Heater Assy, 350w	\$ 154.00
43122	Bushing;Nylon,Split:Heyco#2872	\$ 2.00
43343-1	Cable Assy;Dental Fan Panel	\$ 47.00
45800M	Cover Assy, No Elec, w/Fan Panel	\$1,761.00
45800-1	Cover Assy. W/O Dryer & Module	\$1,090.00
45910	Cover Interlock Switch	\$ 81.00
45351	Cover/Tray Insert With Light Seal	\$ 68.00
45360	Cover; Top With Inserts, Gray Co	\$ 484.00
43530	Dryer Fan	\$ 395.00
43330	Dryer Heater Assembly	\$ 674.00
45415	Electronics Module (See Exchange Parts)	\$2,846.00
43309-1	Film Receptacle Divider	\$ 2.00
44482	Film Sensor Assembly (Shutter PCB)	\$ 105.00
43587	Hinge Block	\$ 8.00
45301	Molded Cover w/Hardware Inserts	\$ 455.00
43121	Plate; Vapor Barrier,304 S/S	\$ 8.00
44422	Power Switch	\$ 6.00
45982	Shipping Material For Cover	\$ 47.00
43584	Shutter Assembly	\$ 40.00
43658	Shutter Guard	\$ 56.00
A/T 2000XR Base Assembly PN 45200		
45200-1	Base Assy, W/O Drive Block & Circ Assy	\$3,122.00
45230	Base Assembly Complete with Water Recirculation	\$5,865.00
45201	Base Molded w/Hardware Inserts & Drain Elbows (No Tubing/Wiring)	\$ 736.00
45440	Base Wiring Harness, Complete	\$ 394.00
43232	Bearing, Circulator Magnet Plastic As Of S/N 464706	\$ 30.00
43885	Bottom Plate Assembly	\$ 71.00
43122	Bushing;Nylon,Split:Heyco#2872	\$ 2.00
29207	Cable Clamp	\$ 3.00
45262	Circulator Drawer Assembly W/O Tank Magnets	\$ 340.00
45691	Circulator Drawer Replacement Kit (Includes Tank Magnets)	\$ 420.00
43233	Circulator Magnet (Chem Tanks)	\$ 35.00
43267	Circulator Magnets (On Cir Motors)	\$ 45.00
43695	Circulator Motor Assy, w/ Fan	\$ 111.00

Film Processing

Part Number	Description	Price
43223	Clamp, Drive Shaft Black: Duerr	\$ 5.00
49963	Clamp, Tank Drain	\$ 2.00
41771	Clamp;Cable,3/8": Richco # E-6	\$ 1.00
45445	Connector Plate, Etched	\$ 40.00
44063	Dev Tube & Cap Assy	\$ 35.00
43384	Developer Circulator Cover - Black	\$ 9.00
43357	Drain Hose Strap	\$ 6.00
44361	Drain Hose, Developer (Black) (6ft 1/2" Id)	\$ 23.00
44362	Drain Hose, Fixer (Red) (6ft 1/2"Id)	\$ 32.00
44363	Drain Hose, Water (White) (6ft 1/2"Id)	\$ 23.00
44206	Drain Manifold, Molded	\$ 42.00
43286	Drain Tube Assy - Black	\$ 18.00
43287	Drain Tube Assy - Red	\$ 18.00
43270	Drain Tube Assy - Tan	\$ 16.00
43271	Drain Tube Washer	\$ 2.00
45240	Drive Block Assembly 2000 XR	\$1,134.00
43259	Electrical Compartment Cover	\$ 70.00
45975	EMI/RFI Filter Replacement, IEC Power Entry	\$ 42.00
43578	Fan Circulator Motor	\$ 5.00
43238-1	Film Feeder, Guide Block	\$ 6.00
45380	Film Inlet Assy (Includes P/N 45567 & P/N 43238-1)	\$ 32.00
45567	Film Inlet, Silkscreened	\$ 26.00
43385	Fixer Circulator Cover - Red	\$ 9.00
44064	Fixer Tube & Cap Assy	\$ 26.00
90219	Fuse Holder	\$ 9.00
41914	Fuse, Delayed Action, 8 Amp, (Pkg. Of 5)	\$ 24.00
43531	Garden Hose Washer & Mesh	\$ 2.00
43258	Gear Cover, Molded	\$ 51.00
45950	Heater Pad Assy Replacement	\$ 252.00
43859	Heater Pad Hold Down	\$ 9.00
43436	Line Cord And Plate Assy	\$ 33.00
45920	Line Cord Replacement, 115 V Plug IEC Socket	\$ 26.00
45043	Locking Ring For 45042	\$ 11.00
43220	Motor, Circulator (Without Motor Magnet)	\$ 57.00
30249	Nut;8-32,Hex,18-8 S.S.	\$ 1.00
45985	PCB Internal Recirculating Pump	\$ 72.00
45437	Plate; Adapter, IEC, Inside Mount	\$ 4.00
45436	Plate;Adapter,IEC,Outside Mount	\$ 3.00
45436P	Plate; Adapter, Powder Coated	\$ 40.00
45278	Plate; Dual Fuse Mount,Blk Anod	\$ 6.00
45439	Plate; Terminal Block Mount, A/T	\$ 17.00
43121	Plate; Vapor Barrier,304 S/S, AT	\$ 8.00
49793	Poppet Valve Kit (Includes 2 Valves & 2 O-Rings)	\$ 63.00
45550	Pressure Regulator Kit (Up To 40 PSI)	\$ 107.00
45070	Recirculator Pump	\$ 359.00


Film Processing

Part Number	Description	Price
44265	Replenisher Assembly	\$ 448.00
45935	Replenisher Assy w/ Circ. Drawer Assy	\$ 794.00
44044	Replenisher Developer Inlet Tubing (Blk) - 96"	\$ 11.00
44042	Replenisher Developer Outlet Tubing (Blk) - 26"	\$ 4.00
44043	Replenisher Fixer Inlet Tubing (Red) - 96"	\$ 14.00
44039	Replenisher Tubing Clamp	\$ 3.00
44041	Replenisher Fixer Outlet Tubing (Red) - 18"	\$ 3.00
30998	Rivet; Push: Heyco #9054	\$ 2.00
45042	Rt Angle Cap (White) Dev Replenisher	\$ 14.00
30568	Screw, #4 X 3/4"L Phillips Panhead Stainless (Gear Cover)	\$ 1.00
30550	Screw, #4 X 3/8"L Phillips Flathead Stainless (Bottom Plate)	\$ 3.00
30539	Screw, #4 X 3/8"L Phillips Panhead Stainless (Gear Cover)	\$ 1.00
30562	Screw, #7 X 1/2"L Phillips Panhead Stainless (Electrical Cover)	\$ 1.00
30569	Screw, 6-32 X 5/8" Phillips Panhead Stainless (Drive Block Mount)	\$ 1.00
30549	Screw;#4 X 5/8,Type Ab, Phil Pa	\$ 1.00
30535	Screw;#7 X 1/2,Hi-Lo Thrd Form	\$ 2.00
30541	Screw;#7 X 3/8,Type Ab, Phil Pa	\$ 2.00
30668	Screw;4-40 X 7/16,Phil Oval Hd	\$ 3.00
30431	Screw;6-32 X 3/8,Phil Pan Hd,1	\$ 1.00
30463	Screw;8-32 X 5/8,Phil Pan Hd,1	\$ 1.00
45952	Shipping Material All 3 Racks	\$ 71.00
45922	Shipping Material For Base, For Returns	\$ 34.00
45982	Shipping Material For Cover	\$ 47.00
41127	Snap Ring:M8 Truarc #5100 31-H	\$ 5.00
43980	Solenoid Assembly (No Flow Restrictor)	\$ 103.00
45980	Solenoid, Water	\$ 95.00
45323	Standoff, Interlock Actuator	\$ 29.00
41286	Teflon Washer (Between Circulator Magnet & Bearing)	\$ 1.00
43437	Terminal Block; Socket Type: Fem	\$ 6.00
45970	Thermistor Replacement Kit	\$ 56.00
30552	Thumb Screw, Circulator Cover	\$ 1.00
45048	Tubing cut to 9'	\$ 15.00
44210	Tubing Fixer Drain	\$ 4.00
44211	Tubing Water Drain	\$ 4.00
44209	Tubing, Developer Drain	\$ 16.00
45945	Varistor Replacement Kit,150v	\$ 55.00
43189	Wash Tank Drain Tube Plug Kit	\$ 10.00
30052	Washer;#8,Flat,18-8 Ss	\$ 1.00
30263	Washer;#8,Lock,Split,18-8 Ss	\$ 1.00
49165	Water Hose	\$ 34.00
45221	Water Inlet J Tube	\$ 9.00
45450	Water Level Sensor	\$ 51.00
45070	Water Recirculating Pump Assy	\$ 359.00

Film Processing

Part Number	Description	Price
	A/T 2000XR Drive Block Assembly PN 45240	
43244	Bearing, Shaft	\$ 3.00
44275	Capacitor, Drive Motor	\$ 14.00
45240	Drive Block Assy	\$ 1,134.00
45235	Drive Block Modified	\$ 230.00
43247	Drive Chain	\$ 39.00
45472	Drive Motor Assembly	\$ 721.00
43253	Drive Motor Mounting Plate	\$ 5.00
43223	Drive Shaft Clamp - Black	\$ 5.00
43246	Drive Shaft Clamp - Red	\$ 5.00
43243	Driver Gear (12 Tooth)	\$ 6.00
43242	Main Drive Shaft	\$ 91.00
30249	Nut;8-32,Hex,18-8 S.S.	\$ 1.00
43254	Rack Latch - Black	\$ 4.00
30539	Screw;#4 X 3/8:Type Ab, Phil Pa	\$ 1.00
30463	Screw;8-32 X 5/8,Phil Pan Hd,1	\$ 1.00
41127	Snap Ring	\$ 5.00
45970	Temp Sensor	\$ 56.00
30052	Washer;#8,Flat,18-8 Ss	\$ 1.00
30263	Washer;#8,Lock,Split,18-8 Ss	\$ 1.00
45450	Water Level Sensor	\$ 51.00
	A/T 2000XR Wash & Dry Assembly PN 45510	
43508	Axle, Shaft	\$ 14.00
43513	Drive Gear (14 Tooth)	\$ 6.00
43503	Gear (28 Tooth)	\$ 3.00
43504	Gear (36 Tooth)	\$ 3.00
43520	Gear (60 Tooth)	\$ 5.00
43736	Gear Replacement Kit (28,36,60 Teeth Gears Incl. 1/8" Drillbit)	\$ 32.00
43694	Gear Replacement Kit Idler (14 Tooth Gear Incl. #14 Drillbit)	\$ 42.00
43527	Guide Roller	\$ 34.00
41292	Idler Gear (14 Tooth) 3/16" Thick	\$ 2.00
41307	Inner Air Baffle Flat	\$ 15.00
43526	Outer Air Baffle U-Shaped	\$ 66.00
43518	Plastic Roller	\$ 34.00
43518W	Plastic Roller - Wash/Dry Rack Only	\$ 34.00
43512	Right End Plate, (Grey)	\$ 133.00
43519	Rubber Roller	\$ 65.00
41821	Rubber Roller Spacer Black	\$ 2.00
41284	Rubber Roller Spacer White	\$ 1.00
31006	Screw #7 X 1/4"L Phillips Panhead Stainless (For Spacer)	\$ 2.00
30567	Screw #7 X 5/8"L Phillips Pan Head Stainless	\$ 2.00
43509	Shaft Gear (24 Tooth)	\$ 5.00
45952	Shipping Material For All 3 Racks For Returns	\$ 71.00
41285	Spacer Bar	\$ 11.00

A/T 2000XR Parts Breakdown



Part Number	Description	Price
43517	Spacer Retainer	\$ 4.00
43524	Spacer Washer	\$ 4.00
45501	Spacer, Nylon	\$ 3.00
43514	Spacer; End Racks, Black	\$ 7.00

A/T 2000XR Fixer Rack PN 45505

43508	Axle Shaft	\$ 14.00
43513	Drive Gear (14 Tooth)	\$ 6.00
43507	End Plate (Red)	\$ 80.00
43503	Gear (28 Tooth)	\$ 3.00
43504	Gear (36 Tooth)	\$ 3.00
43736	Gear Replacement Kit (28,36,60 Teeth Gears, Incl. 1/8" Drillbit)	\$ 32.00
43694	Gear Replacement Kit, Idler 14 Tooth (Incl. #14 Drillbit)	\$ 42.00
43527	Guide Roller	\$ 34.00
43516	Handle	\$ 11.00
41292	Idler Gear (14 Tooth) 3/16" Thick	\$ 2.00
43518	Plastic Roller	\$ 34.00
41821	Rubber Roller Spacer Black	\$ 2.00
41284	Rubber Roller Spacer, White	\$ 1.00
31006	Screw #7 X 1 1/4" L Phillips Panhead Stainless (For Spacer)	\$ 2.00
30595	Screw #7 X 1" L Phillips Panhead Stainless	\$ 3.00
30567	Screw #7 X 5/8" L Phillips Panhead Stainless	\$ 2.00
43509	Shaft Gear (24 Tooth)	\$ 5.00
45952	Shipping Material For All 3 Racks For Returns	\$ 71.00
41285	Spacer Bar	\$ 11.00
43517	Spacer Retainer	\$ 4.00
43524	Spacer Washer	\$ 4.00
45501	Spacer, Nylon	\$ 3.00

A/T 2000XR Developer Rack PN 45500

43508	Axle Shaft	\$ 14.00
43513	Drive Gear (14 Tooth)	\$ 6.00
43502	End Plate (Black)	\$ 75.00
43503	Gear (28 Tooth)	\$ 3.00
43504	Gear (36 Tooth)	\$ 3.00
43736	Gear Replacement Kit (28,36,60 Teeth Gears, Incl. 1/8" Drillbit)	\$ 32.00
43694	Gear Replacement Kit, Idler (14 Tooth, Incl. #14 Drillbit)	\$ 42.00
43527	Guide Roller	\$ 34.00
43516	Handle	\$ 11.00
41292	Idler Gear (14 Tooth) 3/16" Thick	\$ 2.00
43518	Plastic Roller	\$ 34.00
43519	Rubber Roller	\$ 65.00
41284	Rubber Roller Spacer, White	\$ 1.00
31006	Screw #7 X 1 1/4" L Phillips Panhead Stainless (For Spacer)	\$ 2.00
30595	Screw #7 X 1" L Phillips Panhead Stainless	\$ 3.00
30567	Screw #7 X 5/8" L Phillips Panhead Stainless	\$ 2.00

Film Processing

Part Number	Description	Price
43509	Shaft Gear (24 Tooth)	\$ 5.00
45952	Shipping Material For All 3 Racks For Returns	\$ 71.00
41285	Spacer Bar	\$ 11.00
43517	Spacer Retainer	\$ 4.00
43524	Spacer Washer	\$ 4.00
45501	Spacer, Nylon	\$ 3.00
A/T 2000 Plus Major Assemblies		
43500	Developer Rack (See Exchange Parts)	\$ 898.00
43505	Fixer Rack (See Exchange Parts)	\$ 872.00
45100	Leveling Base Assembly	\$ 194.00
43106	Leveling Feet	\$ 3.00
43510	Wash/Dry Rack (See Exchange Parts)	\$1,628.00
A/T 2000 Plus Cover Assembly		
43811	Air Heater Assembly	\$ 154.00
43343-1	Cable Assy; Dental Fan Panel	\$ 47.00
43530	Dryer Fan	\$ 395.00
43330	Dryer Heater Assembly	\$ 674.00
44415	Electronics Module (See Exchange Parts)	\$2,846.00
44482	Film Sensor Assembly PCB	\$ 105.00
43587	Hinge Block	\$ 8.00
44422	Power Switch	\$ 6.00
45982	Shipping Material For Cover For Returns	\$ 47.00
43584	Shutter Assembly	\$ 40.00
43658	Shutter Guard	\$ 56.00
A/T 2000 Plus Base Assembly		
45201	Base Sub-Assy; A/T 2000 XR	\$ 736.00
44440	Base Wiring Harness, Complete	\$ 194.00
43232	Bearing, Ceramic	\$ 30.00
43885	Bottom Plate Assembly	\$ 71.00
43689	Bulkhead Fitting-Water Inlet	\$ 20.00
29207	Cable Clamp	\$ 3.00
45262	Circulator Draw Assy; Grey Color	\$ 340.00
45691	Circulator Drawer Replacement Kit (Includes Tank Magnets)	\$ 420.00
43233	Circulator Magnet	\$ 35.00
43220	Circulator Motor	\$ 57.00
43695	Circulator Motor Assy.	\$ 111.00
43578	Circulator Motor Fan	\$ 5.00
49963	Clamp, Tank Drain	\$ 2.00
43384	Developer Circulator Cover -Black	\$ 9.00
43357	Drain Hose Strap	\$ 6.00
44361	Drain Hose, Developer (Black)	\$ 23.00
44362	Drain Hose, Fixer (Red)	\$ 32.00

Film Processing

Part Number	Description	Price
44363	Drain Hose, Water (White)	\$ 23.00
44206	Drain Manifold, Molded	\$ 42.00
43286	Drain Tube Assy - Black	\$ 18.00
43287	Drain Tube Assy - Red	\$ 18.00
43270	Drain Tube Assy - Tan	\$ 16.00
43271	Drain Tube Washer	\$ 2.00
44240	Drive Block Assembly	\$1,050.00
43267	Drive Motor Magnet	\$ 45.00
43259	Electrical Compartment Cover	\$ 70.00
43567	Film Feeder, Molded	\$ 28.00
43380	Film Inlet Assy	\$ 131.00
43385	Fixer Circulator Cover - Red	\$ 9.00
90219	Fuse Holder	\$ 9.00
41914	Fuse, Delayed Action, 8 Amp. (Pkg.Of 5)	\$ 24.00
43258	Gear Cover, Molded	\$ 51.00
43290	Heater Bar Assembly	\$ 507.00
43436	Line Cord	\$ 33.00
43362	O-Ring For Heater Bar	\$ 2.00
49793	Poppet Valve Kit (Incl. 2 Valves & 2 O-Rings) Gorman Rupp Pumps	\$ 63.00
43437	Receptacle, 3-Pin Connector Plug	\$ 6.00
43224-1	Remote Water Solenoid Assembly	\$ 368.00
44265	Replenisher Assembly	\$ 448.00
44044	Replenisher Developer Inlet Tubing (Blk)	\$ 11.00
44042	Replenisher Developer Outlet Tubing (Blk)	\$ 4.00
44043	Replenisher Fixer Inlet Tubing (Red)	\$ 14.00
44041	Replenisher Fixer Outlet Tubing (Red)	\$ 3.00
44039	Replenisher Tubing Clamp	\$ 3.00
45935	Replenisher Assy w/ Circ. Drawer Assy	\$ 794.00
43531	Screen For Solenoid	\$ 2.00
30568	Screw, #4 X 3/4" L Phillips Panhead Stainless (Gear Cover)	\$ 1.00
30539	Screw, #4 X 3/8" L Phillips Panhead Stainless (Gear Cover)	\$ 1.00
30550	Screw, #4 X 3/8" L Phillips Panhead Stainless (Bottom Plate)	\$ 3.00
30562	Screw, #7 X 1/2" L Phillips Panhead Stainless (Electrical Cover)	\$ 1.00
30569	Screw, 6-32 X 5/8" L Phillips Panhead Stainless (Drive Block Mount)	\$ 1.00
45922	Shipping Material For Base For Returns	\$ 34.00
43931	Solenoid Repl Kit	\$ 100.00
43179	Solenoid, Elbow	\$ 5.00
41286	Teflon Washer (Between Circulator Magnet And Bearing)	\$ 1.00
30552	Thumb Screw, Circulator Cover	\$ 1.00
44210	Tubing Fixer Drain	\$ 4.00
44211	Tubing Water Drain	\$ 4.00
44209	Tubing, Developer Drain	\$ 16.00
43189	Wash Tank Drain Tube Plug Kit	\$ 10.00
41320	Washer Garden Hose	\$ 2.00
43397	Water Inlet Tube Assembly	\$ 12.00

Part Number	Description	Price
43447	Water Level Sensor	\$ 51.00
A/T 2000 Plus Base Assembly		
45922	Shipping Material; AT2000 Base	\$ 34.00
A/T 2000 Plus Drive Block Assembly		
43244	Bearing, Shaft	\$ 3.00
44275	Capacitor, Drive Motor (ECM Motor)	\$ 14.00
44240	Drive Block Assy	\$1,050.00
45235	Drive Block Modified	\$ 230.00
43247	Drive Chain	\$ 39.00
43243	Drive Gear (12 Tooth)	\$ 6.00
44472	Drive Motor Assembly	\$ 720.00
43253	Drive Motor Mounting Plate	\$ 5.00
43223	Drive Shaft Clamp - Black	\$ 5.00
43246	Drive Shaft Clamp - Red	\$ 5.00
43242	Main Drive Shaft	\$ 91.00
43254	Rack Latch - Black	\$ 4.00
41127	Snap Ring	\$ 5.00
43447	Water Level Sensor	\$ 51.00
A/T 2000 Plus Wash Rack		
43508	Axle, Shaft	\$ 14.00
43513	Drive Gear (14 Tooth)	\$ 6.00
43512	End Plate	\$ 133.00
43514	End Spacer	\$ 7.00
43503	Gear (28 Tooth)	\$ 3.00
43504	Gear (36 Tooth)	\$ 3.00
43520	Gear (60 Tooth)	\$ 5.00
43736	Gear Replacement Kit (28,36,60 Teeth Gears, Incl. 1/8" Drillbit)	\$ 32.00
43694	Gear Replacement Kit, Idler, 14 Tooth (Includes #14 Drillbit)	\$ 42.00
43527	Guide Roller	\$ 34.00
41292	Idler Gear (14 Tooth) 3/16" Thick	\$ 2.00
41307	Inner Air Baffle Flat	\$ 15.00
43526	Outer Air Baffle U-Shape	\$ 66.00
43518	Plastic Roller	\$ 34.00
43518W	Plastic Roller - Wash/Dry Rack Only	\$ 34.00
43519	Rubber Roller	\$ 65.00
41821	Rubber Roller Spacer Black	\$ 2.00
41284	Rubber Roller Spacer White	\$ 1.00
30550	Screw, #4 X 3/8" L Phillips Flathead Stainless	\$ 3.00
30567	Screw, #7 X 5/8" L Phillips Panhead Stainless	\$ 2.00
43509	Shaft Gear (24 Tooth)	\$ 5.00
45952	Shipping Material For All 3 Racks For Returns	\$ 71.00
41285	Spacer Bar	\$ 11.00
A/T 2000 Plus Fixer Rack		
43508	Axle Shaft	\$ 14.00
43513	Drive Gear (14 Tooth)	\$ 6.00
43507	End Plate	\$ 80.00
43514	End Spacer	\$ 7.00
43503	Gear (28 Tooth)	\$ 3.00
43504	Gear (36 Tooth)	\$ 3.00
43736	Gear Replacement Kit (28,36,60 Teeth Gears, Incl. 1/8" Drillbit)	\$ 32.00
43694	Gear Replacement Kit, Idler, 14 Tooth (Incl. #14 Drillbit)	\$ 42.00
43527	Guide Roller	\$ 34.00
43516	Handle	\$ 11.00
41292	Idler Gear (14 Tooth) 3/16" Thick	\$ 2.00
43518	Plastic Roller	\$ 34.00
43519	Rubber Roller	\$ 65.00
41821	Rubber Roller Spacer Black	\$ 2.00
41284	Rubber Roller Spacer, White	\$ 1.00
30550	Screw, #4 X 3/8" L Phillips Flathead Stainless	\$ 3.00
30595	Screw, #7 X 1" L Phillips Panhead Stainless	\$ 3.00
30567	Screw, #7 X 5/8" L Phillips Panhead Stainless	\$ 2.00
43509	Shaft Gear (24 Tooth)	\$ 5.00
45952	Shipping Material For All 3 Racks For Returns	\$ 71.00
41285	Spacer Bar	\$ 11.00
A/T 2000 Plus Developer Rack		
43508	Axle Shaft	\$ 14.00
43513	Drive Gear (14 Tooth)	\$ 6.00
43502	End Plate	\$ 75.00
43514	End Spacer	\$ 7.00
43503	Gear (28 Tooth)	\$ 3.00
43504	Gear (36 Tooth)	\$ 3.00
43736	Gear Replacement Kit (28,36,60 Teeth Gears, Incl. 1/8" Drillbit)	\$ 32.00
43694	Gear Replacement Kit Idler, 14 Tooth (Includes #14 Drillbit)	\$ 42.00
43527	Guide Roller	\$ 34.00
43516	Handle	\$ 11.00
41292	Idler Gear (14 Tooth) 3/16" Thick	\$ 2.00
43518	Plastic Roller	\$ 34.00
43519	Rubber Roller	\$ 65.00
41284	Rubber Roller Spacer, White	\$ 1.00
30550	Screw, #4 X 3/8" L Phillips Flathead Stainless	\$ 3.00
30595	Screw, #7 X 1" L Phillips Panhead Stainless	\$ 3.00



Film Processing

Part Number	Description	Price
30567	Screw, #7 X 5/6" L Phillips Panhead Stainless	\$ 2.00
43509	Shaft Gear (24 Tooth)	\$ 5.00
45952	Shipping Material For All 3 Racks For Returns	\$ 71.00
41285	Spacer Bar	\$ 11.00
43517	Spacer Retainer	\$ 4.00
43524	Spacer Washer	\$ 4.00
Accessories for A/T 2000 Plus, & 2000XR		
40270	Air Tech Developer & Fixer	\$ 87.00
40280	Air Tech Starter Kit; Case Of 6 12-Oz Bottles	\$ 66.00
44170	Chemistry Recovery Containers-Two 2 1/2 Gal.	\$ 89.00
44160	Chemistry Recovery Containers-Two 5 Gal.	\$ 109.00
45150	Chemistry Tray	\$ 27.00
40140	Cleaning Film - Pkg Of 12	\$ 60.00
45110	Daylight Loader; A/T 2000XR; Gray	\$1,379.00
43500	Developer Rack	\$ 898.00
43309-1	Film Divider	\$ 2.00
43505	Fixer Rack (See Exchange Parts)	\$ 872.00
43945	Formula 2000 Cleanser, Twin Pack	\$ 57.00
45100	Leveling Base Kit, A/T 2000XR (Grey)	\$ 194.00
43106	Leveling Feet	\$ 3.00
43115	Prop Stick	\$ 47.00
43699	Recovery Drain Hose Kit(Separate Hoses For Dev+Fix)-A/T 2000	\$ 155.00
40035	Sponge, Black	\$ 2.00
40034	Sponge, Orange	\$ 2.00
43965	Spray 2000	\$ 45.00
43999	Video A/T 2000	\$ 39.00
44999	Video A/T 2000 Plus	\$ 39.00
44054	Wall Chart For A/T 2000 Plus	\$ 5.00
45054	Wall Chart; At 2000 XR	\$ 5.00
43510	Wash/Dry Rack	\$1,628.00
45550	Water Pressure Regulator	\$ 107.00
45999	Video; A/T 2000 XR	\$ 39.00
45900-DVD	DVD; A/T 2000 XR	\$ 38.00
Panoramic DLL For 2000 Plus, & 2000XR		
43117-4	Adapter; Mount Used To Mount DLL To Base (Between Machine & DLL)	\$ 60.00
44637	Amber View Glass	\$ 23.00
40155	Cloth Film Bag	\$ 77.00
45630	Cover w/View Glass (Grey)	\$ 351.00
40153	Cuff Assembly (One Per Unit)	\$ 234.00
43117-3	Felt Gasket Used To Mount DLL To Base, (Between Machine & DLL)	\$ 8.00
43117-1	Foam Gasket (Between Machine & DLL)	\$ 4.00
44622	Gasket; Cover, DLL	\$ 6.00
40154	Gasket-Holds Cuff	\$ 74.00

Film Processing

Part Number	Description	Price
43114	Hardware, Mounting	\$ 13.00
45110	Panoramic Daylight Loader; Fits All A/T2000 Models	\$1,379.00
44647	Red View Glass S/N 6720 & Above	\$ 99.00
45118	Support Frame Assy Kit, Daylight Loader	\$ 209.00
43995	View Glass; Red, A/T 2000 DLL Old Style	\$ 75.00

Peri-Pro III Parts Breakdown



Part Number	Description	Price
Peri-Pro III		
94640	Access Door With Label	\$ 16.00
94993	Base w/Drain Port & Hose Kit	\$ 119.00
94310	Cover Drive Motor	\$ 79.00
90213	Cover Plate, Bent	\$ 24.00
90212	Cover Plate, Flat	\$ 15.00
90750	Daylight Loader Hole Plug Kit	\$ 26.00
94120	Drain Socket, Elbow	\$ 2.00
94994	Drain Valve Replacement Kit w/Panel Nut, With Hoses	\$ 53.00
94211	Dryer Fan Shroud	\$ 28.00
94630	Dryer Heater Replacement Kit, 370w (35 Ohms) Black Leads	\$ 84.00
90400	Dryer Transport	\$ 132.00
94996	Dryer Transport Gear Train Replacement Kit	\$ 40.00
90880	Fan & Motor Assy, Dryer, 115v	\$ 121.00
90555	Film Inlet Grill	\$ 7.00
90810	Film Receptacle W/3 Dividers	\$ 21.00
90219	Fuse Holder	\$ 9.00
90265	Fuse. 5 Amp (Pkg Of 5)	\$ 12.00
94525	Gear Idler, 18 Tooth, Dovetail Lock (Dryer Transfer)	\$ 3.00
94519	Gear Idler, 21 Tooth, Dovetail Lock (Dryer Transfer)	\$ 3.00
94545	Gear Train Cover	\$ 4.00
94995	Gear, Main Drive, 66 Tooth	\$ 79.00
90611	Heater Grill	\$ 7.00
94410	Heater Harness Assy	\$ 110.00
94970	Heater Plate Replacement Kit	\$ 252.00
94925	Light Suppressor Kit	\$ 2.00
94991	Main Cover	\$ 306.00
90147	Main Drive Motor Replacement Kit	\$ 285.00
94500	Main Transport, Peri Pro III	\$1,015.00
90226-1	Middle Plate	\$ 67.00
94450	Pc Board	\$ 228.00
94275	Plate Switches, Stamped Peri-Pro III	\$ 17.00
94134	Power Cord 115v	\$ 84.00
94269	Power Switch/On Light , Amber	\$ 16.00
94977	Process Switch/Ready Light, Green	\$ 21.00
90127	Retaining Ring For Dryer Fan	\$ 4.00
94111	Side Plate	\$ 69.00
117327	Socket, Thermal Fuse	\$ 2.00
94278	Thermal Fuse Dryer - 98 Deg.C	\$ 16.00
94976	Thermistor Arm Replacement Kit	\$ 26.00
94975	Thermistor Probe Replacement Kit	\$ 77.00
94461	Transformer, 40va, 120v-24v, Chemistry Heating System	\$ 72.00
94124	Tube, Drain, 48", White	\$ 7.00
94139	Valve Handle Only	\$ 2.00
94130	Wash Tank W/60604 O-Ring	\$ 37.00

Film Processing

Film Processing

Part Number	Description	Price
94039	Wash Tank Replacement Kit, w/o Drain System	\$ 32.00
94125	Water Bottle	\$ 11.00
94920	Water Bottle With Neck	\$ 26.00
Peri-Pro III Main Transport		
90522	Bearing,"U" Transfer Arm	\$ 4.00
90523	Bushing "L" Transfer Arm	\$ 3.00
90524	Bushing "U" Transfer Arm	\$ 2.00
94550	Gear & Shaft Assy, Large "L" Transfer Arm, Developer	\$ 24.00
90550	Gear & Shaft Assy;"L" Transfer Arm W/#90516 Gear	\$ 16.00
90519	Gear Idler, 24 Tooth	\$ 4.00
90516	Gear Transfer Arm, 24 Tooth	\$ 4.00
90518	Gear Transfer Arm, 24 Tooth	\$ 5.00
90517	Gear Transfer Arm, 60 Tooth	\$ 7.00
94501	Gear, 30 Tooth Speed Reducer, On Top Of 90517 Developer Gear	\$ 6.00
90521	Gear, Timer	\$ 11.00
94540	Hub & Pin Assy	\$ 37.00
94530	L Transfer Arm Assy (Lower Arm)	\$ 34.00
90520	L Transfer Arm Assy (Upper Arm)	\$ 35.00
94500	Main Transport, Peri-Pro III	\$1,015.00
60604	O-Ring, For Wash Tank Replacement, 5/16id X 1/16w:#011	\$ 1.00
90505	Pin, Top Plate	\$ 16.00
90543	Plastic Snap Ring	\$ 2.00
90546	Rubber Bumper - Red	\$ 2.00
90510	Shutter Plate Assy	\$ 90.00
90548	Shutter Release Bar	\$ 14.00
90512	Shutter Return Spring	\$ 5.00
30392	Shutter Spring Pin Rear	\$ 1.62
30391	Spring Pin, Shutter Release Bar	\$ 4.00
90545	Spring Shutter Release Bar	\$ 2.00
90590	Track Plate, Top Assy	\$ 144.00
94988	Transport Repair Kit	\$ 134.00
Peri-Pro		
90113	Bottom Plate	\$ 12.00
94310	Cover Modification; Drive Motor Cover, Peri-Pro III	\$ 79.00
90213	Cover Plate, Bent	\$ 24.00
90212	Cover Plate, Flat	\$ 15.00
90750	Daylight Loader Hole Plug Kit	\$ 26.00
91993	Dryer Drive Gear W/#30389 Set Screw	\$ 37.00
90223	Dryer Drive Motor Assy 115v	\$ 129.00
90610	Dryer Heater Assy	\$ 81.00
90400	Dryer Transport	\$ 132.00
90880	Fan & Motor Assy, 115v	\$ 121.00
90102	Fan Only For Dryer	\$ 7.00

Part Number	Description	Price
90555	Film Inlet Grille	\$ 7.00
90810	Film Receptacle w/Dividers	\$ 21.00
90219	Fuse Holder	\$ 9.00
90265	Fuse, 5 Amp (Pkg Of 5)	\$ 12.00
90611	Heater Grill	\$ 7.00
90997	Main Cover	\$ 235.00
90995	Main Drive Gear, 30 Tooth, #304, W/#30390 Set Screw	\$ 21.00
90147	Main Drive Motor Repl' Kit	\$ 285.00
90500	Main Transport	\$1,004.00
90221	Power Switch Kit Serial #14999 & Below (Push Button)	\$ 37.00
90262	Power Switch, Rocker Style: Serial #15000 & Above	\$ 29.00
90127	Retaining Ring For Dryer Fan	\$ 4.00
90111	Side Plate	\$ 63.00
117327	Socket, Thermal Fuse	\$ 2.00
90278	Thermal Fuse Assy Kit,(112 C)P	\$ 16.00
Accessories Peri-Pro I & III		
94640	Access Door; With Filling Label	\$ 16.00
90985	Adult Anterior Film #1 Size Carriers, Pkg Of 5	\$ 49.00
90990	Bite Wing Film #3 Size Carriers, Pkg Of 5	\$ 49.00
90117	Developer Tank, Black	\$ 29.00
90555	Film Inlet Grille	\$ 7.00
90810	Film Receptacle w/ Three Dividers	\$ 21.00
90118	Fixer Tank, Red	\$ 29.00
43945	Formula 2000 Cleanser, Twin Pack	\$ 57.00
90602	Inlet Cover	\$ 4.00
94500	Main Transport, Peri Pro III	\$1,015.00
90500	Main Transport, Peri-Pro & Peri-Pro II	\$1,004.00
90831	Occlusal #4 Size Carriers, Pkg Of 3	\$ 49.00
60604	O-Ring For Wash Tank	\$ 1.00
90980	Pedo Film #0 Size Carriers, Pkg Of 5	\$ 49.00
43975	Peri-Pro Transport Cleaning Kit	\$ 129.00
90121	Service Tray For Peri-Pro & Peri-Pro III	\$ 41.00
91992	Service Tray w/ Heater Bar Supports For Peri-Pro II	\$ 87.00
94130	Wash Tank Assy, White w/Drain Port & (60604) O-Ring For Pp III	\$ 37.00
90119	Wash Tank, White, Peri Pro & Peri Pro II	\$ 32.00
90703	Waste Deposit Box	\$ 11.00
94920	Water Bottle w/Spout	\$ 26.00
Daylight Loader For Peri-Pro		
90730	Cuff Assembly, 1 Each	\$ 79.00
90735	Cuff Assembly, 1 Pair	\$ 156.00
90090	Daylight Loader For Peri Pro	\$ 360.00
90750	Daylight Loader Hole Plug Kit	\$ 26.00
90713	Felt Gasket	\$ 5.00
90714	Handle	\$ 7.00

Film Processing

Part Number	Description	Price
30412	Handle Screw	\$ 1.00
90711	Hinge Pin, Left	\$ 5.00
90712	Hinge Pin, Right	\$ 5.00
90716	Hinge Socket	\$ 5.00
90720	View Glass, Amber - Special Order Only	\$ 120.00
90780	View Glass, Red	\$ 101.00
Duplicator/Daylight Loader For Peri-Pro		
49161	Bullet Catch With Strike Plate	\$ 3.00
92750	Duplicator Door Assy	\$ 349.00
92000	Film Duplicator/Daylight Loader	\$ 799.00
90714	Handle Dupl Assy	\$ 7.00
92002	Hardware Kit For Mounting	\$ 18.00
92720	Light Filter Kit	\$ 12.00
116574	Magnetic Switch	\$ 21.00
92746	P.C. Board & Reflector Assy	\$ 317.00
92786	Power Cord, Lugged	\$ 42.00
92765	Prism (Cover For Lamp Assy)	\$ 25.00
105A	Replacement Bulb	\$ 11.00
92755	Tray Base Assy (Of Dup Door Incl Foam Sheet #2 Film Holder Grid)	\$ 142.00
Main Transport Peri-Pro I & III		
90522	Bearing, "U" Transfer Arm	\$ 4.00
90546	Bumper, Rubber	\$ 2.00
90523	Bushing, "L" Transfer Arm	\$ 3.00
90524	Bushing, "U" Transfer Arm	\$ 2.00
90550	Gear & Shaft Assy, "U" Transfer Arm, W/#90516 Gear	\$ 16.00
90519	Gear Idler, 24 Tooth	\$ 4.00
90518	Gear, "U" Transfer Arm, 24 Tooth	\$ 5.00
90525	Gear, "L" Transfer Arm, 15 Tooth	\$ 4.00
90516	Gear, "U" Transfer Arm, 24 Tooth	\$ 4.00
90517	Gear, "U" Transfer Arm, 60 Tooth	\$ 7.00
90521	Gear, Timer	\$ 11.00
90540	Hub & Pin Assy	\$ 37.00
90520	L Transfer Arm, Serial #9520 & Above	\$ 35.00
90500	Main Transport	\$1,004.00
90505	Pin, Top Plate	\$ 16.00
90543	Plastic Snap Ring	\$ 2.00
90510	Shutter Plate Assy	\$ 90.00
90548	Shutter Release Bar	\$ 14.00
30391	Spring Pin, Shutter Release Bar	\$ 4.00
90545	Spring Shutter Release Bar	\$ 2.00
90512	Spring, Shutter Return	\$ 5.00
90988	Transport Repair Kit	\$ 129.00
90530	U Transfer Arm Assy, S/N 9520 & Above (See 90572 Sn9519 & Below)	\$ 39.00

Intraoral X-Ray

Part Number	Description	Price
Provecta 70 Intraoral X-Ray		
A2920	16" On Center Mounting Plate Standard, 12 3/4"H X 19"W	\$ 209.00
A2070	16" On Center Mounting Plate, Large (Optional) 16 3/4"X20"W	\$ 273.00
A2140-3	16" Straight Arm replacement Kit	\$ 683.00
A2934	2 Position Electric Mains Terminal Block	\$ 8.00
A2140-2	24" Straight Arm Replacement Kit	\$ 683.00
A2140-1	35" Straight Arm Replacement Kit	\$ 683.00
115952	630ma,250v,Fast-Acting Fuse	\$ 1.00
A2930	Bearing Block Replacement Kit	\$ 446.00
A2739	Bearing Block Replacement Kit for Mobile	\$ 281.00
A2970	Bearing Block Straight Arm Rotational Brake Assy Replacement Kit	\$ 27.00
A2705	Brake Kit for Mobile X-Ray	\$ 65.00
A2345	Cable Assy. Ribbon	\$ 52.00
A2921	Cap 1 3/4" Hole Plug For 16" Wallplate 16 3/4"Hx20"W & Upper Yok	\$ 47.00
A2985	Circuit Breaker Replacement Kit (Incl. 2 Breakers)	\$ 19.00
A2915	Cover Sleeve Tubehead Yoke With Screws Replacement Kit	\$ 39.00
A2808	End Cap Replacement Kit, Scissor Arm	\$ 63.00
115951	Fuse; 10a,250v,Fast-Acting,3ab	\$ 1.00
A2950	Handheld Exposure Holder Replacement Kit	\$ 37.00
A2992	Jacking Plate Adjustment Hex Wrench 3/16" T-Handle	\$ 17.00
A2983	Jacking Plate Screw	\$ 2.00
A2925	Keypad Membrane (Display) Replacement Kit	\$ 208.00
A2424	Line Cord Conversion Kit, 115 V	\$ 50.00
A2715	Lock Arm Kit	\$ 126.00
30387	Mobile Arm Clip 2"	\$ 9.00
A2733	Electronics Assembly for Mobile	\$2,415.00
A2986	Mounting Hardware & Tools Kit	\$ 152.00
A2020	Plastic Cap and Cover Kit, 1 1/8"	\$ 30.00
A2975	Power Switch Replacement Kit	\$ 17.00
A2055	Provecta 70 Packaging	\$ 308.00
A2290	Remote (Hall Mount) X-Ray Electronics Replacement Kit	\$1,141.00
A2991	Remote (Hall Mount) Control X-Ray Cover Replacement Kit	\$ 102.00
A2290-1	Remote (Wallmount) Conversion Kit (To Convert From Standard Unit	\$1,400.00
A2090	Remote Exposure Switch (Doorbell)	\$ 350.00
A2980R	Remote(Hall Mount) X-Ray Exposure Switch Handheld Replacement Kit	\$ 194.00
A2993	Scissor Arm Adjustment Hex Wrench 5/16" T-Handle	\$ 63.00
A2025	Scissor Arm Cap(S) Replacement Kit Includes Painted Gary Screws	\$ 131.00
A2997	Scissor Arm Cap(S) Replacement Screw Painted Grey	\$ 4.00
A2120	Scissor Arm Replacement Kit	\$1,118.00
A2960	Scissor Arm Rotational Brake Assembly Replacement Kit for Wall Mount	\$ 95.00
A2900	Standard Wall Mount Electronics Replacement Kit	\$1,628.00
A2980	Standard X-Ray Exposure Handheld Replacement Kit	\$ 194.00
A2911	Standard X-Ray Control Cover Replacement Kit	\$ 60.00
A2040	Straight Arm End Cap Replacement Kit	\$ 58.00
A2943	Straight Arm Rotation Stop Screw	\$ 63.00

Intraoral X-Ray

Part Number	Description	Price
A2969	Straight Arm Wire Access Cover With Screws Replacement Kit	\$ 16.00
A2999	Touch-up Paint Kit	\$ 29.00
A2908	Tubehead Lower Yoke Cap	\$ 58.00
A2905	Tubehead Retaining C-Clip & Grounding Clip Replacement Kit	\$ 90.00
A2982	Tubehead Yoke Rotation Adjustment Tool	\$ 32.00
A2254	X-Ray Tubehead Nose Cone Ring	\$ 58.46
A2250	X-Ray Tubehead Replacement Kit	\$2,924.00
A2990	X-Ray Tubehead Simulator	\$ 8.00

Part Number	Description	Price
A/T ScanX, ScanX Intraoral and ScanX Duo		
73248-0	Barrier Envelopes, Size 0 (Qty. 100)	\$ 19.00
73248-1	Barrier Envelopes, Size 1 (Qty. 100)	\$ 19.00
73248-2	Barrier Envelopes, Size 2 (Qty. 300)	\$ 50.00
73248-3	Barrier Envelopes, Size 3 (Qty. 100)	\$ 25.00
73248-4	Barrier Envelopes, Size 4 (Qty. 50)	\$ 25.00
73554	Cable, EPP, Packaged	\$ 76.00
B3554	Cable, USB, Packaged	\$ 15.00
B2010	Cleaning Sheets, ScanX, Pkg 12	\$ 140.00
B2020	Cleaning Sheets, ScanX, Pkg 25	\$ 245.00
B3010	Cleaning Sheets, ScanX I/O, Pkg 12	\$ 115.00
B3020	Cleaning Sheets, ScanX I/O, Pkg 25	\$ 195.00
D1010	Cleaning Sheets, ScanX Duo, Pkg 12	\$ 60.00
D1020	Cleaning Sheets, ScanX Duo, Pkg 25	\$ 115.00
74050	Dust Cover Kit	\$ 82.00
B7264	Electronic Assy; Main/Interface, ScanX ILE	\$1,806.00
B7464	Electronic Assy; Main/Interface, ScanX IO ILE	\$1,645.00
73596-10	EPP Communication Cable, 10ft	\$ 55.00
73596-20	EPP Communication Cable, 20ft	\$ 113.00
73596-35	EPP Communication Cable, 35ft	\$ 155.00
B7910	Fiber Op to Ribbon Kit	\$ 483.00
73431	Holder, Inner LED	\$ 1.62
73432	Holder, Outer LED	\$ 1.62
D1146	Light Seal Kit, Inlet, ScanX Duo	\$ 169.00
D1356	Light Seal Kit, Outlet, ScanX Duo	\$ 39.00
74610	Light Seal Kit, Inlet, ScanX ILE	\$ 168.00
73096	Line Cord	\$ 30.00
73767	Lower Cover Assy	\$ 129.00
74364	Main Assy PCB, Replacement Kit, ScanX Classic	\$1,646.00
B3464	Main Assy PCB, Replacement Kit, ScanX I/O	\$1,645.00
B7264	Main Assy PCB, Replacement Kit, ScanX ILE	\$1,806.00
B7464	Main Assy PCB, Replacement Kit, ScanX I/O ILE	\$1,645.00
D1364	Main Assy PCB, Replacement Kit, ScanX Duo	\$1,499.00
B7176	PCB, 20-24 LED Control, ScanX ILE	\$ 188.00
B7625	PCB, 20 LED / Roller Holder	\$ 558.00
73329	PCB, Stepper Motor Control Assy	\$ 184.00
74588	PCB, Start of Line, Replacement Kit	\$ 131.00
74365	PCB, Stepper Motor, Replacement Kit	\$ 210.00
73445-0	Phosphor Storage Plate (PSP), Size 0 (Qty.2)	\$ 62.00
73445-1	Phosphor Storage Plate (PSP), Size 1 (Qty.2)	\$ 62.00
73445-2	Phosphor Storage Plate (PSP), Size 2 (Qty.4)	\$ 126.00
73445-3	Phosphor Storage Plate (PSP), Size 3 (Qty.2)	\$ 74.00
73445-4	Phosphor Storage Plate (PSP), Size 4 (Qty. 1)	\$ 96.00
73578-57	Phosphor Storage Plate (PSP), TMJ 5"X7" (Qty.1)	\$ 578.00
73578-6	Phosphor Storage Plate (PSP), Panoramic, 15cm X 30cm (Qty.1)	\$ 980.00

Part Number	Description	Price
A/T ScanX Eraser		
73578-5	Phosphor Storage Plate (PSP), Panoramic, 5"X12" (Qty.1)	\$ 950.00
73578-8	Phosphor Storage Plate (PSP), Cephalometric, 8"X10" (Qty.1)	\$1,045.00
B8910	Phosphor Storage Plate Cleaning Wipes (Box of 50)	\$ 60.00
B8900	Phosphor Storage Plate Cleaning Wipes (10 Boxes of 50 count each)	\$ 600.00
73566-0	Plate Guide, Size 0 (Qty. 1)	\$ 30.00
73566-1	Plate Guide, Size 1 (Qty. 1)	\$ 30.00
73566-2	Plate Guide, Size 2 (Qty. 1)	\$ 30.00
73566-3	Plate Guide, Size 3 (Qty. 1)	\$ 30.00
73566-4	Plate Guide, Size 4 (Qty. 1)	\$ 30.00
73470	Plate Transfer Box	\$ 34.00
B7095	Power Supply, ScanX ILE (100w)	\$ 197.00
73097	Power Supply, ScanX Non-ILE (60w)	\$ 175.00
73471	Power Supply, ScanX Non-ILE (60w), with Line Cord	\$ 196.00
D1150	Transport Arch, ScanX Duo	\$2,300.00
73050	Transport Arch Replacement Kit, ScanX Non-ILE	\$4,013.00
B7760	Transport Arch Replacement Kit, ScanX ILE	\$4,200.00
D1920	Transport Belt Maintenance Kit, ScanX Duo, Set of 2	\$ 126.00
B7920	Transport Belt Maintenance Kit, ScanX, Set of 4	\$ 247.00
D1932	Transport Belt Drive Assy Kit, ScanX Duo, Set of 2	\$ 626.00
B7434	Transport Belt Drive Assy Kit, ScanX ILE, Set of 4	\$1,575.00
74634	Transport Belt Drive Assy Kit, ScanX Non-ILE, Set of 4	\$1,188.00
B7995	Upper Arch Cover Replacement Kit, ScanX ILE	\$ 289.00
74362	Upper Arch Cover Replacement Kit, ScanX Non-ILE	\$ 371.00
74078	Ballast Replacement Kit	\$ 116.00
74079	Bulb Replacement Kit 24 Watt/Ea	\$ 83.00
74080	Circuit Breaker Replacement Kit	\$ 21.00
73403-1	Magnet Only	\$ 45.00
73403-2	Magnetic Sensor Assy	\$ 47.00
74081	Pc Board Replacement Kit For Model 73800	\$ 189.00

Digital Imaging

Description	Price
ScanX Flat Rate Repair/ Upgrade Program	
Program Summary:	
1. All Dates referenced go by Date of Manufacture.	
2. ScanX Warranty periods:	
a. Upgrade warranty: 1-year	
b. Repairs & refurbished: 6 months (or remainder of original warranty, if greater)	
3. Dealer will be invoiced at full price as an advance replacement, upon receipt of the original unit by Air Techniques in New York, dealer will be credited accordingly.	
B7200U ScanX 73400 Non ILE Upgrade (New ILE unit not Refurbished)	\$12,995.00
B7100U ScanX 83100 Non ILE Upgrade (New ILE unit not Refurbished)	\$11,995.00
Repair B7200/B7100 or newer	\$4,995.00
Repair B7200/B7100 or older	\$6,995.00
Repair 73400/83100 USB Models	\$7,995.00
Hot Swap B7200 2008 or newer (Replace with an B7200R)	\$4,995.00
Hot Swap B7200 2007 or older (Replace with an B7200R)	\$7,995.00
Hot Swap B7100 2008 or newer (Replace with an B7100R)	\$4,995.00
Hot Swap B7100 2007 or older (Replace with an B7100R)	\$6,995.00
Hot Swap D1000 (Replace with a D1000R)	\$3,995.00

Part Number	Description	Price
Polaris Intraoral Camera System		
D7014	Blue Connection Grip	\$ 3.23
D7060	Polaris USB Handpiece	\$4,725.00
D7040	Polaris USB Umbilical (10')	\$ 85.00
D7520	Polaris Universal Camera Handpiece	\$4,725.00
D7560	Polaris Universal Control Module	\$1,299.00
D7540-6	Polaris Universal Umbilical (6')	\$ 307.00
D7540-9	Polaris Universal Umbilical (9')	\$ 340.00
D7540-15	Polaris Universal Umbilical (15')	\$ 390.00
D7540-30	Polaris Universal Umbilical (30')	\$ 515.00
D7540-50	Polaris Universal Umbilical (50')	\$ 650.00
D7170	Polaris Holder with Chairside Adapter	\$ 60.00
D7120	Polaris Sheaths - Box of 500	\$ 300.00
D7125	Polaris Sheaths - Box of 100	\$ 80.00
Acclaim Portable Intraoral Camera		
A5120	Acclaim Handpiece Holder With Chairside Adaptor	\$ 58.00
A5125	Adapter, USB To Gameport	\$ 54.00
A5115	Cable Assy, 10ft, Footswitch To Gameport	\$ 200.00
B3554	Cable USB, High Speed Cable, 6ft	\$ 15.00
77389	Cable, RCA Composite Video, 6ft	\$ 18.00
77379	Cable, S-Video, 6ft (For Acclaim With Capture Card)	\$ 38.00
A5105	Cable, Umbilical, 15ft, Handpiece To Control Module	\$ 339.00
A5105-9	Handpiece To Control Module Cable, 9 Ft, Acclaim Camera	\$ 323.00
A5105-30	Handpiece To Control Module Cable, 30 Ft, Acclaim Camera	\$ 541.00
A5165	Handpiece To Control Module Cable, 50 Ft, Acclaim Camera	\$ 683.00
A5105-6	Handpiece To Control Module Cable, 6 Ft, Acclaim Camera	\$ 323.00
A5105	Handpiece To Control Module Cable, 15 Ft, Acclaim Camera	\$ 323.00
A5145	Holder; Handpiece, Acclaim, Boxed	\$ 23.00
A5110	Infection Control Sheaths (Box Of 100)	\$ 61.00
A5100	Infection Control Sheaths (Box Of 500)	\$ 219.00
Spectra Caries Detection Aid		
D4060	Spectra Handpiece	\$4,725.00
D4020	Spectra Umbilical (10')	\$ 85.00
D4170	Spectra Holder with Chairside Adapter	\$ 60.00
D4120	Spectra Sheaths - Box of 500	\$ 300.00
D4125	Spectra Sheaths - Box of 100	\$ 80.00
D4140	Spectra Spacers (10mm) - Box of 25	\$ 125.00
A/T SLC		
77186	Infection Control Sheaths For Handpiece (Box Of 100)	\$ 75.00
77100	Infection Control Sheaths For Handpiece (Box Of 500)	\$ 258.00
A3112	Replacement Lamp, 100w, 12v	\$ 45.00
B3554	USB 2.0 Cable	\$ 15.00

Part Number	Description	Price
Vistacam Omni & Omni Ic4		
77243	Hospital Grade Line Cord, 10 Ft: For Omni And Omni Ic4	\$ 32.00
77186	Infection Control Sheaths For Handpiece (Box Of 100)	\$ 75.00
77100	Infection Control Sheaths For Handpiece (Box Of 500)	\$ 258.00
77121	Replacement Lamp, 150w	\$ 53.00

Exchange Parts

Part Number		Outright Sale Price	Exchange Sale Price
STS			
54646	Replacement Pump Assembly STS-3	\$ 4,996.00	\$ 2,520.00
54146	Replacement Pump Assembly STS-5	\$ 7,990.00	\$ 2,830.00
54110	Replacement Pump Assembly STS-15	\$17,730.00	\$ 7,098.00
Replacement VacStar Pumps			
55110	1 HP for VacStar 2, 5 or 20 with Electrical Box	\$ 2,100.00	\$ 1,676.00
55178	1 HP for VS50 or VS50H - Catridge Only	\$ 1,995.00	\$ 1,786.00
55111	2 HP for VacStar 4, 8 or 40 with Electrical Box	\$ 2,625.00	\$ 2,040.00
55179	2 HP for VS80 or VS80H - Catridge Only	\$ 2,415.00	\$ 2,166.00
AirStar Heads			
85890	AS10	\$ 2,825.00	\$ 2,359.00
85960	AS21	\$ 3,638.00	\$ 3,082.00
85910	AS22, AS50 & AS70	\$ 3,638.00	\$ 3,082.00
85940	AS30	\$ 2,825.00	\$ 2,359.00
A/T 2000, A/T 2000 Plus, A/T 2000XR			
45415	Electronics Assy; A/T 2000 XR	\$ 2,846.00	\$ 1,180.00
44415	Electronics Assy; A/T 2000 Plus	\$ 2,846.00	\$ 1,180.00
43500	Developer Rack - A/T 2000 & A/T 2000 Plus	\$ 898.00	\$ 733.00
45500	Rack Assy; Developer, AT 2000XR	\$ 898.00	\$ 733.00
43505	Fixer Rack - A/T 2000 & A/T 2000 Plus	\$ 872.00	\$ 665.00
45505	Fixer Rack - AT 2000XR	\$ 872.00	\$ 695.00
43510	Wash/Dry Rack - A/T 2000 & A/T 2000 Plus	\$ 1,628.00	\$ 1,208.00
45510	Wash/Dry Rack - A/T 2000XR	\$ 1,628.00	\$ 1,208.00
ScanX			
B7200U	ScanX Classic ILE Upgrade	\$21,995.00	\$11,020.00
B7100U	ScanX IO ILE Upgrade	\$13,695.00	\$11,020.00
73050	Transport Arch Assy Exchange	\$ 4,013.00	
	If Unit is 0 to 2 years old		N/C
	If Unit is 2 to 3 years old		\$ 1,902.00
	If Unit is 3 to 4 years old		\$ 2,481.00
	If Unit is over 4 years old		\$ 3,060.00
B7760	Transport Arch Replacement Kit ILE	\$ 4,200.00	
	If Unit is 0 to 2 years old		N/C
	If Unit is 2 to 3 years old		\$ 2,205.00
	If Unit is 3 to 4 years old		\$ 2,812.00
	If Unit is over 4 years old		\$ 3,368.00

We Serve You
NATIONWIDE

Professional Sales Associates
At Your Service




Mark Asplund




Victor Yoon



Rick Orme



Lloyd Scott



Dan Knox




Pat Moylan



Sam Paola



Paul Myers



Mel Seagle




Michael Park



Bill Reinhardt



Matt Highland



Aaron Renner



Michel Cabanie



Jim Domitrz



Don Rogers, III



Matt Langdon



Milt Baisley



Milton Lorring



Charlie Clarke



Bill Highland



Tom O'Berry



Rich Ellis



Keith O'Berry



Robert Rodriguez



Our 200,000 square foot manufacturing facility in Melville, New York. Innovation and reliability since 1962.

Air Techniques is a leading manufacturer of dental equipment for 5 decades. We have earned our reputation through our commitment to ensuring customer satisfaction, and we are proud to provide superior quality to dealers and dental practitioners worldwide.

Air Techniques' products are distributed through authorized dealers only.

Corporate Headquarters

1295 Walt Whitman Road | Melville, NY 11747-3062 | Phone 800.247.8324 | Fax 888.247.8481

Western Facility

291 Bonnie Lane, Suite 101 | Corona, CA 92880-2804 | Phone 800.247.8324 | Fax 951.898.7646

www.airtechniques.com

