

Parts Book





Table of Contents

Table of Contents	2
Table of Contents	3
Statim 2000 Parts Front View	4
Statim 2000 Parts Rear View	5
Statim 2000 Parts Bottom View	6
Statim 2000 Parts Top View No Cover	7
Statim 2000 Revision 5XX & Higher Software Alex Parts	8
Statim 2000 Individual Parts Pictures	9-16
Statim 5000 Parts Front View	17
Statim 5000 Parts Rear View	18
Statim 5000 Parts Top View No Cover	19
Statim 5000 Revision 5XX & Higher Software Alex Parts	20
Statim 5000 Individual Parts Pictures	21-27
Statim 7000 Front View	28
Statim 7000 Front Accessories & Rear View	29
Statim 7000 Front & Side View No Cover	30
Statim 7000 Rear View No Cover	31
Statim 7000 Cassette	32
Statim 7000 Condenser, Waste & Overflow Bottles	33
Statim 900 Parts Front View	34
Statim 900 Parts Rear View	35
Statim 900 Parts Left Side View	36
Statim 900 Parts Right Side View	37
Statim 900 Parts Top View	38
Statim 900 Individual Parts Pictures	39-42
Bravo Accessories	43
Bravo Parts Table 1	44
Bravo Parts List for Table 1	45
Bravo Parts Table 2	46
Bravo Parts List for Table 2	47
Bravo Parts Table 3	48
Bravo Parts List for Table 3 & Bravo Software	49
Bravo Parts Table 4	50
Bravo Parts List for Table 4	51
Bravo Parts Table 4A & 4B	52
Bravo Parts List for Tables 4A & 4B	53
Bravo Parts Table 5	54
Bravo Parts List for Table 5	55
Bravo Parts Table 6	56
Bravo Parts List for Table 6	57
Bravo Parts Table 7	58
Bravo Parts List for Table 7	59
Bravo Parts PCB Fuses	60
Bravo Individual Parts Pictures	61-65
Bravo Water Intake Pull Kit #01-111774S	66
Bravo Direct to Drain Kit #01-111775S	67
Quantum Parts Front View Door Closed	68
Quantum Parts Front View Door Open	69

Table of Contents


Quantim Parts Rear View with Cover	70
Quantim Parts Right Side No Cover	71
Quantim Parts Right Side Front View	72
Quantim Parts Rear View No Cover	73
Quantim Parts Left Side Front & Back No Cover	74
Quantim Parts Left Side Top Front No Cover	75
Quantim Parts Bottom View	76
Hydrim C51W Parts Front View Door Closed	77
Hydrim C51W Parts Front View Door Open	78
Hydrim C51W Parts Right Side	79
Hydrim C51W Right Side View & Inside Chamber with Dryer	80
Hydrim C51W Parts Bottom	81
Hydrim C51W Spare Parts Update 06-TSB-230	82
Hydrim C51W Bur Screen	83
Hydrim C51W w/Touchscreen Front View Door Closed	84
Hydrim C51W w/Touchscreen Front View Door Open w/4 Cassette Rack	85
Hydrim C51W w/Touchscreen Front View Door Open w/Baskets	86
Hydrim C51W w/Touchscreen Rear View	87
Hydrim C51W w/Touchscreen Front View Door Open No Trolley	88
Hydrim C51W w/Touchscreen Right Side View No Cover	89
Hydrim C51W w/Touchscreen Right Side View Behind Chemical Bracket	90
Hydrim C51W w/Touchscreen Left Side View No Cover	91
Hydrim C51W w/Touchscreen Top View No Cover	92
Hydrim C51W w/Touchscreen View Behind Controller Facia	93
Hydrim C51W w/Touchscreen Bottom View	94
Hydrim C51W w/Touchscreen Individual Parts Pictures	95-97
Hydrim L110W Parts Front View Door Closed	98
Hydrim L110W Parts Front View Door Open	99
Hydrim L110W Parts Left Side No Cover	100
Hydrim L110W Parts Right Side No Cover	101
Hydrim L110W Parts Inlet Solenoids & Controllers	102
Hydrim L110W Preventative Maintenance Kit	103
Hydrim L110W w/Glass Door Front View	104
Hydrim L110W w/Glass Door Chamber w/Trolleys	105
Hydrim L110W w/Glass Door Chamber	106
Hydrim L110W w/Glass Door Behind Front Kickplate	107
Hydrim L110W w/Glass Door Control Boards	108
Hydrim L110W w/Glass Door Right Side Top View	109
Hydrim L110W w/Glass Door Right Side Bottom View	110
Hydrim L110W w/Glass Door Rear View	111
Hydrim L110W w/Glass Door Top View	112
Statmatic Parts	113-117
Statmatic Couplers	118
Aquastat Boiler	119
Aquastat Chassis	120
Aquastat Reservoir	121
Aquastat Individual Parts Pictures	122-125
VistaBrite	126-128

Statim 2000 Front View




Item #	Part #	Description
1	01-101611S	Cover without Facia
2	01-100564S	Cover Complete for Model 101102
2	01-106121S	Cover Complete for Model 01-121101
3	01-100563S	Facia Complete for Model 101102
3	01-106122S	Facia Complete for Model 01-121101
3	01-103518S	Facia without LCD for Model 101102
3	01-106124S	Facia without LCD for Model 01-121101
4	01-100008A	Cassette Complete
5	01-104381S	LCD includes Decal (located inside Facia)
6	01-106950S	LCD Decal Only
7	01-101783S	Reservoir Cap & Filter
7	01-109300S	Filter Only, Reservoir (not shown under Reservoir Cap)
8	01-101707S	Keypad (old style not shown)
8	01-109680S	Keypad (pictured above)
9	01-100842S	Facia Gasket (not shown located inside Facia)
10	01-106675S	Label Kit Qty 10
11	01-103856S	Isoplate (not shown, not field replaceable requires alignment jig)

Statim 2000 Rear View




Item #	Part #	Description
1	01-100564S	Cover Complete (Model #101102)
1	01-106121S	Cover Complete (Model #01-121101)
2	01-102119S	Biological Filter (not on all units)
3	01-104284S	Biological Filter Bracket (behind Biological Filter)
4	01-100207S	Air Filter (not shown behind cover on Air Compressor)
5	01-107817S	Fitting 3/8" for Biological Filter Tubes (Pkg. of 2)
6	01-100782S	Push In Fitting
7	01-106087S	Line Filter
8	01-100573S	Power Switch
9	01-101783S	Reservoir Cap & Filter
9	01-109300S	Filter Only, Reservoir (not shown under Reservoir Cap)

Statim 2000 Bottom View



Item #	Part #	Description
1	01-100780S	Bumper
2	01-101553S	Leveler
2	01-104180S	Leveler Repair Kit (includes Leveler and Insert for Leveler)
3	01-108152S	Drain Tubing Kit
4	01-104343S	Plug Drain Tubing
5	74-100335A	Drain Tube 1"
6	01-112091S	Drain Tube Clip


Statim 2000 Revision 4XX & Lower Software Top View No Cover




Item #	Part #	Description
1	01-100684S	Bubble Level
2	01-100561S	Armature (not field replaceable requires alignment jig)
3	01-100574S	Check Valve
4	01-100560S	Compressor
5	01-106787S	Pressure Relief Valve
6	01-104461S	Thermal Fuse (not shown, under Boiler)
7	01-108979S	Boiler (software must be 5XX or higher to replace in field)
8	01-100571S	Microswitch (not shown, behind Boiler)
9	01-100570S	Probe Bracket (see note)
10	01-100782S	Push In Fitting
11	01-100557S	Solenoid
12	01-103571S	Water Quality Sensor
13	01-100568S	Pump
13	01-110445S	Pump Ulka (6XX> Software)
14	01-106087S	Line Filter
15	01-100573S	Power Switch
16	01-108541S	Power Cable Filter/Switch
17	01-106349S	PCB Power Cable
18	01-110868S	Reservoir with Water Quality Sensor
18	01-100565S	Reservoir with Float Switch (not shown)
18	01-106459S	Reservoir with Float Switch & Sensor (not shown)
19	01-109365S	PCB Rev. 7 (see note)

Note: This part is not field replaceable requires complete calibration

Statim 2000 Revision 5XX & Higher Software Alex Parts



Item #	Part #	Description
1	01-108979S	Boiler (software must be 5XX or higher to replace in field)
2	01-108983S	Validation Thermocouple 2000



Item #	Part #	Description
1	01-109365S	PCB Rev. 7 (see note)
2	01-108981S	Temperature Adapter Board for Alex
3	01-108985S	PCB Adapter Board Revision 5XX Software Alex 2000
4	01-108976S	Microprocessor Kit Revision 5XX
4	01-109385S	Microprocessor Kit Revision 6XX for Rev. 7 PCB

Note: This part is not field replaceable requires complete calibration

Statim 2000 Parts



Air Filter
#01-100207S



Alex Kit 2XX Software
#01-108996S



Alex Kit 4XX Software
#01-108993S



Alex Kit PCB Rev 3.0 &
4.2 #01-109716S




Allen Key 9/64
#01-103469S



Armature
#01-100561S



Biological Filter
#01-102119S



Biological Filter Bracket
#01-104284S



Boiler Alex
#01-108979S



Bubble Level
#01-100684S



Bumper
#01-100780S



Cassette Complete
#01-100008A

Statim 2000 Parts



Cassette Lid
#01-100834S



Cassette Seal & Lube Kit
#01-100028S



Cassette Tray Only
#01-103475S



Cassette Tray with Rack
#01-100271A



Chassis Kit (Frame)
#01-108388S



Check Valve
#01-100574S



Compressor
#01-100560S



Conductivity Meter
#01-103139S



Control Box
#103141S



Copper Tubing Validation Thermocouple
#01-109128S



Cover Complete Model 01-121101
#01-106121S



Cover Complete Model 101102
#01-100564S

Statim 2000 Parts



Cover without Facia
#01-101611S



Drain Tube 1"
#74-100335A



Drain Tubing Kit
#01-108152S




Exhaust Tubing
#01-100204S



Facia Complete
Model 01-121101
#01-106122S



Facia Complete
Model 101102
#01-100563S



Facia Gasket
#01-100842S



Facia without LCD
Model 01-121101
#01-106124S



Facia without LCD
Model 101102
#01-103518S



Filter Reservoir
#01-109300S



Fitting 3/8 for Biological
Filter Tubes
#01-107817S



Fuse 2A
#01-103498S


Statim 2000 Parts




Fuse 15A
#01-103472S




Fuse .25A
#01-103497S



Isoplate
#01-103856S



Keypad Model 01-121101
#01-109680S



Keypad Model 101102
#01-101707S



Label Kit
#01-106675S



LCD
#01-104381S



LCD Decal
#01-106950S



Leveler
#01-101553S



Leveler Repair Kit
#01-104180S



Line Filter
#01-106087S




Microprocessor 4XX
#01-108742S

Statim 2000 Parts



Microprocessor 5XX
#01-108976S



Microswitch
#01-100571S



Operator's Manual
#01-106938S



Packaging with Inserts
#01-108348S



PCB 4XX
#01-108741S



PCB 5XX
#01-109043S



PCB Adapter Board 4XX
#01-108775S



PCB Adapter Board 5XX
#01-108985S



PCB Power Cable
#01-106349S



Plug Drain Tubing
#01-104343S



Power Cable Filter/Switch
#01-108541S



Power Cord
#01-100572S

Statim 2000 Parts




Power Cord 2000
#01-101647S



Power Switch
#01-100573S



Pressure Relief Valve (PRV)
#01-106787S



Probe Bracket
#01-100570S



Probe Bracket Gasket
#01-100865S



Pump
#01-100568S



Pump Bracket
#01-108363S



Pump Tester Bottle
#01-100713S



Push In Fitting
#01-100782S



Rack Unwrapped
#01-103945S



Rack Wrapped (Mesh)
#01-106653




Reservoir Cap with Filter
#01-101783S

Statim 2000 Parts



Reservoir with Float
#01-100565S



**Reservoir with
Water Quality Sensor**
#01-110868S



Solenoid Coil
#01-104303S



Solenoid Plunger Wrench
#01-103471S




Solenoid Valve Complete
#01-100557S



Solenoid Valve Repair Kit
#01-100998S



Strain Relief
#53-100219A



**Temperature Adapter
Board**
#01-108981S



Test Strip Indicators
#SCI134



Thermal Fuse
#01-104461S



Validation Thermocouple
#01-108983S



Waste Bottle Complete
#01-100812S

Statim 2000 Parts



Waste Bottle Fitting
#01-100735S




Waste Bottle Only
#01-100724S




Water Quality Sensor
#01-103571S

Statim 5000 Front View




Item #	Part #	Description
1	01-101783S	Reservoir Cap & Filter
1	01-109300S	Filter Only, Reservoir (not shown under Reservoir Cap)
2	01-104300S	Cover Complete
3	01-104299S	Facia Complete
4	01-101650S	Facia Gasket (not shown located inside Facia)
5	01-101613S	Cassette Complete
6	01-103857S	Isoplate (not shown, not field replaceable requires alignment jig)
7	01-104382S	LCD (located inside Facia)
8	01-109671S	Keypad
9	01-101625S	Cover without Facia
10	01-101658S	Cassette Handles (Lid & Tray)

Statim 5000 Rear View




Item #	Part #	Description
1	01-101755S	Push In Fitting
2	01-101652S	Air Compressor Filter (not on all units)
3	01-104284S	Biological Filter Bracket
4	01-102119S	Biological Filter
5	01-103993S	Power Socket (older units)
5	01-106087S	Line Filter (newer units)
6	01-100573S	Power Switch

Statim 5000 Parts Top View No Cover




Item #	Part #	Description
1	01-100684S	Bubble Level
2	01-101623S	Armature (not field replaceable requires alignment jig)
3	01-106787S	Pressure Relief Valve
4	01-108980S	Boiler (software must be 5XX or higher to replace in field)
5	01-104462S	Thermal Fuse (not shown, under boiler)
6	01-101628S	Solenoid
7	01-101627S	Check Valve
8	01-101621S	Microswitch
9	01-101618S	Probe Bracket (see note)
10	01-101652S	Air Compressor Filter (not on all units)
11	01-102119S	Biological Filter
12	01-101761S	Pump
13	01-101619S	Compressor
13	01-112665S	Compressor Kit Thomas to Medo
14	01-103571S	Water Quality Sensor
15	01-110869S	Reservoir with Water Quality Sensor
15	01-101622S	Reservoir with Float Switch & Sensor (not shown)
16	01-109395S	PCB Rev. 7 for 5000 no printer USA (see note)
16	01-109823S	PCB Rev. 7 for 5000 with printer USA (see note)

Statim 5000 Revision 5XX & Higher Software Alex Parts



Item #	Part #	Description
1	01-108980S	Boiler (software must be 5XX or higher to replace in field)
2	01-108984S	Validation Thermocouple 5000



Item #	Part #	Description
1	01-109395S	PCB Rev. 7 for 5000 no printer USA (see note)
1	01-109823S	PCB Rev. 7 for 5000 with printer USA (see note)
2	01-108981S	Temperature Adapter Board for Alex
3	01-108987S	PCB Adapter Board Alex 5000 5XX & Higher Software
4	01-108978S	Microprocessor Kit Revision 5XX
4	01-109391S	Microprocessor Kit Revision 6XX for Rev. 7 PCB

Note: This part is not field replaceable requires complete calibration

Statim 5000 Parts



Air Compressor Filter
#01-101652S



Air Compressor
Filter Holder
#01-104329S



Alex Kit 2XX
#01-108998S



Alex Kit 4XX
#01-108995S



Armature
#01-101623S



Boiler Alex
#01-108980S



Bubble Level
#01-100684S




Bumper
#01-100780S




Cassette Complete
#01-101613S



Cassette Extended
#01-104104



Cassette Handles
#01-101658S



Cassette Handles Extended
#01-104472S

Statim 5000 Parts



Cassette Lid
#01-101757S



Cassette Seal & Lube Kit
#01-101649S



Cassette Tray with Rack
#01-101614S



Chassis Kit (Frame)
#01-108389S



Check Valve
#01-101627S



Compressor
#01-101619S



Conductivity Meter
#01-103139S



Control Box
#01-103141S



Copper Tubing
Validation Thermocouple
#01-109129S



Cover Complete
#01-104300S




Cover without Facia
#01-101625S



Drain Tube 1'
74-100335A

Statim 5000 Parts



Drying Plates
#01-103935



Exhaust Tubing
#01-100204S



Facia Complete
#01-104299S



Facia Gasket
#01-101650S



Filter, Reservoir
#01-109300S



Fuse 2A
#01-103498S



Fuse 15A
#01-103472S




Fuse .25A
#01-103497S



Isoplate
#01-103857S



Keypad
#01-109671S



LCD
#01-104382S



Leveler
#01-101553S

Statim 5000 Parts



**Leveler Repair Kit
#01-104180S**




**Line Filter
#01-106087S**



**Microprocessor 4XX
#01-108746S**



**Microprocessor 5XX
#01-108978S**



**Microswitch
#01-101621S**



**Operator's Manual
#01-106938S**



**Packaging with Inserts
#01-108327S**



**PCB 4XX
#01-108744S**



**PCB 5XX
#01-109045S**



**PCB Adapter Board 4XX
#01-108777S**



**PCB Adapter Board 5XX
#01-108987S**



**PCB Power Cable
#01-106350S**

Statim 5000 Parts



Plug Drain Tubing
#01-104343S



Power Cable Filter/Switch
#01-108541S



Power Cord
#01-101647S



Power Socket
#01-103993S



Power Switch
#01-100573S




Pressure Relief Valve
#01-106787S



Printer 5XX> Software
#01-210000



Printer Case with Damper
#01-103470S



Printer Door
#01-102159S



Printer Paper (Pkg. 10)
#01-101657S



Printer Paper Holder
#01-101973S



Printer Spring/Latch
#01-101971S


Statim 5000 Parts




Probe Bracket
#01-101618S



Probe Bracket Gasket
#01-100865S



Pump
#01-101761S



Pump Recovery Kit
#01-104313S




Pump Tester Bottle
#01-100713S



Push In Fitting
#01-101755S



Rack Unwrapped Mesh
#01-101709S



Rack Wrapped
#01-102054S



Reservoir w/Sensor
#01-110869S



Reservoir Cap with Filter
#01-101783S



Solenoid Coil
#01-104303S



Solenoid Plunger Wrench
01-103471S


Statim 5000 Parts



Solenoid Valve Complete
01-101628S



Solenoid Valve Repair Kit
01-100998S




**Temperature Adapter
Board**
01-108981S



Test Strip Indicators
SCI134



Thermal Fuse
01-104462S



Validation Thermocouple
#01-108984S



Waste Bottle Complete
01-100812S



Waste Bottle Fitting
01-100735S



Waste Bottle Only
01-100724S







Water Quality Sensor
01-103571S

Statim 7000 Front View







Item #	Part #	Description
1	01-110298S	Cover Complete
2	01-110523S	Cover Facia Front Cassette
3	01-110520S	Cover Reservoir Compartment
4	01-110522S	Cover Facia Front Left
5	01-110521S	Cover Logo SciCan
6	01-110526S	Cover Keypad Complete
7	01-110525S	Cover LCD Only
8	01-110524S	Cover Facia Gasket (not shown located inside cover)
9	01-110293S	Reservoir Complete
10	01-110328S	Reservoir Lid
11	01-110512S	Leveling Foot
12	01-110288S	Cassette Complete

Statim 7000 Front Accessories & Rear View





Item #	Part #	Description
1	SCWF1	Water Filter Cartridge
	SCWF6	Water Filter Cartridge (Package of 6)
	SCWFBP	Water Bypass Cartridge (not shown)
2	01-110293S	Reservoir Complete
3	01-102119S	Biological Filter
4	01-110503S	Push In Fitting
5	01-101652S	Air Compressor Filter
6	01-104329S	Air Compressor Filter Holder
7	01-104343S	Plug Drain Tubing
8	01-110505S	Line Filter (Plug Receptacle)
9	01-110504S	Power Switch
10	01-110327S	Reservoir Plunger Kit w/Large Outer O-ring 7000

Statim 7000 Parts









Item #	Part #	Description
1	01-102119S	Biological Filter
2	01-110472S	Bio Filter Holder
3	01-110496S	Armature
4	01-110498S	Armature Isoplate
5	01-110493S	Reservoir Inlet
6	01-110488S	PCB RFID Antenna Assembly
7	01-110502S	Armature Microswitch
8	01-110515S	Thermocouple Validation
9	01-110471S	Pump, Water
10	01-110497S	Armature Guide Rails (Pkg. Of 2)
11	01-110490S	Fan
12	01-110482S	PCB Power Board Auxiliary Heater
13	01-110716S	Fuse 10A 250V (Pkg. Of 2)
14	01-103498S	Fuse 2A
15	01-110477S	Auxiliary Heater Triac (Solid State Relay)
16	01-109816S	Fuse 0.8A
17	01-110483S	PCB RFID Board
18	01-110480S	PCB 7000 (not field replaceable)

Statim 7000 Rear View No Cover






Item #	Part #	Description
1	01-110474S	Auxiliary Heater
2	01-110527S	Solenoid Chamber (Valve 2)
3	01-110513S	Thermocouple Inlet/Outlet
4	01-110528S	Solenoid Validation (Valve 1)
5	01-110514S	Thermocouple Outlet (not field replaceable)
6	01-110472S	Boiler 7000
7	01-110510S	Check Valve
8	01-110471S	Pump, Water
9	01-110478S	Compressor

Statim 7000 Cassette




Item #	Part #	Description
1	01-110288S	Cassette Complete
2	01-110290S	Cassette Lid
3	01-110289S	Cassette Tray
4	01-110291S	Cassette Inlet Coupling
5	01-110292S	Cassette Outlet Coupling
6	01-110296S	Cassette Coupling Seal Kit (Seals for both Couplings included)
7	01-110294S	Cassette Perforated Rack
8	01-110330S	Cassette Tray Handle
9	01-110297S	Cassette Tray Exhaust Duct
10	01-110825S	Cassette Steam Deflector
11	01-106848S	Filter Exhaust Duct Insert
12	01-110824S	Cassette Steam Deflector Clip
13	01-110295S	Cassette Seal Kit
14	01-110329S	Cassette Lid Handle
Not Shown	01-111590S	Cassette O-ring Lubricator Kit

Statim 7000 Condenser, Waste & Overflow Bottles



Item #	Part #	Description
1	01-110506S	Bottle Waste Complete
2	01-100735S	Waste Bottle Fitting
3	01-111497S	Bottle Tubing Condenser to Waste (Includes Quick Disconnect)
4	01-112922S	Stem 1/4" X 5/16"
5	01-111660S	Bottle Condenser Fitting 3/8
6	01-111395S	Bottle Condenser
7	01-110508S	Bottle Waste Condenser Fitting Quick Disconnect
8	01-110300S	Bottle Waste Overflow Kit

Statim 900 Parts Front View




Item #	Part #	Description
1	01-101554S	Bubble Level
2	01-103527S	Keypad
3	01-103508S	Facia
4	01-103636S	Wire Basket W/Insert
5	01-103608S	Sterilization Chamber
6	01-101560S	Drawer Front
7	01-101562S	Drawer
8	01-103559S	Reservoir Cap
9	01-101558S	Cover

Statim 900 Parts Rear View




Item #	Part #	Description
1	01-104343S	Plug Drain Tubing
2	01-100573S	Power Switch
3	01-106087S	Line Filter
4	01-100782S	Push In Fitting

Statim 900 Parts Left Side View




Item #	Part #	Description
1	01-103508S	Facia
2	01-101542S	Reservoir
3	01-103608S	Sterilization Chamber
4	01-101562S	Drawer
5	01-110242S	Boiler (software must be 5XX or higher to replace in field)
6	01-102038S	Pump
7	01-101550S	Solenoid
8	01-100782S	Push In Fitting

Statim 900 Parts Right Side View



Item #	Part #	Description
1	01-103508S	Facia
2	01-101542S	Reservoir
3	01-103608S	Sterilization Chamber
4	01-101562S	Drawer
5	01-101552S	Microswitch
6	01-101552S	Microswitch

Statim 900 Parts Top View



Item #	Part #	Description
1	01-103508S	Facia
2	01-101542S	Reservoir
3	01-110242S	Boiler (software must be 5XX or higher to replace in field)
4	01-103085S	Pressure Relief Valve
5	01-101550S	Solenoid
6	01-100782S	Push In Fitting
7	01-110462S	PCB Rev. 7 (see note)
8	01-102038S	Pump
9	01-102061S	Motor Control PCB
10	01-104461S	Thermal Fuse (not shown, under boiler)
11	01-101545S	Motor

Note: This part is not field replaceable requires complete calibration

Statim 900 Parts



Basket
#01-103636S



Boiler Alex
#01-110242S



Bubble Level
#01-101554S



Cap Replacement
#01-101582S




Chamber Sterilization
#01-103608S



Cover
#01-101558S



Drawer
#01-101562S



Drawer Front
#01-101560S



Facia
#01-103508S



Filter Insert Exhaust
#01-106848S



Fuse .25A
#01-103497S



Fuses 15A
#01-103472S

Statim 900 Parts



Fuse 2A
#01-103498S



Keypad
#103527S



LCD
#01-103138S



Leveler
#01-101553S



Leveler Repair Kit
#01-104180S



Manual Op
#01-106847S



Microswitch
#01-101552S



Motor
#01-101545S



Motor Control PCB
#01-102061S



Plug Drain Tubing
#01-104343S



Power Cord
#01-101647S



Power Switch
#01-100573S

Statim 900 Parts



Pressure Relief Valve
#01-103085S



Pump
#01-102038S



Push-in Fitting
#01-100782S



Reservoir
#01-101542S



Reservoir Cap
#01-103559S



Seal (Exhaust Only)
#01-103144S



Seal (Pressure Only)
#01-103143S



Seal Assembly
(Complete w/Plate)
#01-106898S



Seal Frame
#01-102064S



Seal Kit
(Exhaust & Pressure)
#01-101543S



Seal Plate Flat
#01-102065S



Seal Screw
#01-102063S

Statim 900 Parts



Solenoid Coil
#01-104303S



Solenoid Valve Complete
#01-101550S



Solenoid Valve Repair Kit
#01-100998S



Thermal Fuse
#01-104461S




Waste Bottle Complete
#01-100812S




Water Filter In-Line
#01-106637S

Bravo Accessories



Item #	Part #	Description
1	C1XP076000Y	Bravo Tray for 17 & 17V
1	C1BG345000Y	Bravo Tray for 21V
2-1	STXX0080000	Bravo Tray Handle
2-2	49800030000	Bravo Wrench 12MM for Drain Filter
2-3	49800010000	Bravo Allen Wrench
3	A0XP0010000	Bravo Fill Bottle 2500 ml
4	110000003W0	Bravo Drain Tubing w/Fitting
5	47200010000	Bravo Biological Filter
6	47200030000	Bravo Chamber Filter Screen Metal
7	261000009MN	Bravo Elbow Barbed Fitting 1/8"
8	STXX0250000	Bravo Printer Paper Thermal Roll
Not Shown	SD-306	Bravo Operator's Manual
Not Shown	18300000400	Bravo Chamber Cleaner – Stahl Fix Classic
Not Shown	A7BP5920000	Bravo 17 Packaging & Foam
Not Shown	A7BG5380000	Bravo 21 Packaging & Foam


Bravo Parts Table 1



Parts List for Table 1

Item #	Table #	Part #	Description
D	1	A9BM5450000	Standard Kit Bravo 17 & 17V
D	1	A9BL0090000	Standard Kit Bravo 21V
1	1	C1XP076000Y	Tray for Bravo 17 & 17V
1	1	C1BG345000Y	Tray for Bravo 21V
2-1	1	STXX0080000	Tray Extractor Handle Bravo
2-2	1	49800030000	Wrench 12MM Bravo
2-3	1	49800010000	Allen Wrench Bravo
3	1	A0XP0010000	Reservoir Fill Bottle 2500 ml Bravo
4	1	110000003W0	Drain Tubing w/Fitting Bravo
5	1	47200010000	Biological Filter Bravo
6	1	47200030000	Water Filter Screen Bravo
7	1	STXX0250000	Thermal Paper Roll
8	1	C1BP583000Y	Rack for Chamber Bravo 17 & 17V
8	1	C1BG534000Y	Rack for Chamber Bravo 21V
9	1	C1BH0060001	Cover for Bravo 17 & 17V (1 Piece Cover)
9	1	C1BG4680001	Cover for Bravo 21V (1 Piece Cover)
9a	1	C1BP6230002	Cover Top/Right Side Bravo 17 & 17V (2 Piece Cover)
9a	1	C1BG5760002	Cover Top/Right Side Bravo 21V (2 Piece Cover)
9b	1	C1BP6240002	Cover Left Side Bravo 17 & 17V (2 Piece Cover)
9b	1	C1BG5770002	Cover Left Side Bravo 21V (2 Piece Cover)
10b	1	C3BP1850000	Door Cover Plastic Bravo (Original Style)
11	1	C5BP5500000	Printer Thermal Bravo (Printer only)
11	1	A9BP5430000	Printer Kit Bravo Gam (to upgrade Bravo with no printer)
11	1	A9JP0230000	Printer Kit Bravo Troll (to upgrade Bravo with no printer)
12a	1	C3JP0140000	Printer Door Bravo (New Style)
12b	1	C3BP1840000	Printer Door Bravo (Original Style)
13	1	230B04L20K0	Printer Door Pin 4 x 20 Bravo
14	1	42000070000	Power Switch Bravo
15b	1	C3BP1830000	Facia Front Plastic Bravo (Original Style)
16	1	C5BP1230000	LCD Bravo
17b	1	C6JP0010000	Keypad Bravo (Original Style)
18b	1	C6JM0010000	Label Bravo 17 (Original Style)
18b	1	C6JP0020000	Label Bravo 17V (Original Style)
18b	1	C6JG0010000	Label Bravo 21V (Original Style)
19	1	47200010000	Biological Filter Bravo
20	1	C1BH005000M	Left Crossbar Bravo 17 & 17V (1 Piece Cover)
20	1	C1BP655000M	Left Crossbar Bravo 17 & 17V (2 Piece Cover)
20	1	C1BG317000M	Left Crossbar Bravo 21V (1 Piece Cover)
20	1	C1BG575000M	Left Crossbar Bravo 21V (2 Piece Cover)
21	1	C1BP152000M	Right Crossbar Bravo 17 & 17V (1 Piece Cover)
21	1	C1BP654000M	Right Crossbar Bravo 17 & 17V (2 Piece Cover)
21	1	C1BG318000M	Right Crossbar Bravo 21V (1 Piece Cover)
21	1	C1BG574000M	Right Crossbar Bravo 21V (2 Piece Cover)


Bravo Parts Table 2



Parts List for Table 2

Item #	Table #	Part #	Description
22	2	C4BP2740000	Door Thermo insulating panel Bravo
23	2	C0BP260000P	Door Front Metal Bravo
24	2	C7BP2640000	Door Spring Bravo
25	2	C0BP060000Z	Door Fork Bravo
26	2	49000050000	Door Bushing Bravo
27	2	C0BP0500000	Door Fork Pin Bravo
28	2	203C06L30K0	Door Center Screw Bravo
29	2	C0BP079000N	Door Plate Adjustment Screw Bravo
30	2	C4BP5690000	Door Plate Thermo Insulating Bravo
31	2	A1BP027000Y	Door Seal Plate Bravo
32	2	48000050000	Door Gasket Bravo
33	2	C0BP2720000	Door Hinge Pin Bravo
34	2	A1BP5520000	Door Hinge Bravo
35	2	C2BP0010000	Chamber Bravo 17 & 17V
35	2	C2BG5300000	Chamber Bravo 21 & 21V
36	2	C0BP0430000	Door Locking System Pin Bravo
37	2	A1BP2710000	Door Locking System Complete Bravo
38	2	A0BP6250000	Door Locking System Motor Assy. Complete Bravo
39	2	232A03L25K0	Door Locking System Pin 3 x 25 Bravo
40	2	C0BP0450000	Door Locking System Motor Coupling Bravo
41	2	C0BP0460000	Door Locking System Intermediate Coupling Bravo
42	2	C0BP0590000	Door Locking System Closing Hook Bravo
43	2	43300010000	Door Locking System Hook Microswitch (open) Bravo
44	2	43300010000	Door Locking System Door Microswitch Bravo
45	2	43300080000	Door Locking System Hook Microswitch (close) Bravo
46	2	28600160000	Temperature Sensor Ferrule Bravo
47	2	2860014000N	Temperature Sensor Nut Bravo
48	2	43000190000	Temperature Sensor PT1000 – PT1 Bravo
49	2	41000180000	Chamber Heater 17 & 17V 1500W 120V Bravo
49	2	41000280000	Chamber Heater 17 & 17V 2000W 220V Bravo
49	2	41000190000	Chamber Heater 21V 1500W 120V Bravo
49	2	41000290000	Chamber Heater 21V 2000W 220V Bravo
50a	2	43000180000	Temperature Sensor PT1000 – PT3 Bravo
50b	2	43000180000	Temperature Sensor PT1000 – PT4 Bravo
50c	4	43000180000	Temperature Sensor PT1000 – PT2 Bravo
51	2	43000140000	Safety Thermostat TG400 Chamber Heater Bravo
52	2	C4BP2140100	Thermo Insulating Cover for Chamber Bravo 17 & 17V Bravo
52	2	C4BG3280000	Thermo Insulating Cover for Chamber Bravo 21 & 21V Bravo
53	2	C4BP5570000	Thermo Insulating Ring for Chamber Bravo
63	2	28600290000	Fitting for Chamber Filter Bravo
65	2	47200030000	Chamber Filter Screen Metal Bravo


Bravo Parts Table 3



Parts List for Table 3

Item #	Table #	Part #	Description
54	3	A0BP2900000	Reservoir Assembly Complete Bravo
54a	3	C3BP1700000	Reservoir and Waste Water Tank Bravo
54b	3	43200060000	Float Switch Waste Water Tank Bravo
54c	3	43200050000	Float Switch Clean Water Reservoir Bravo
55	3	C1BP175000M	Support Bar Reservoir Bravo 17 & 17V
55	3	C1BG320000M	Support Bar Reservoir Bravo 21 & 21V
56	3	A9BP5860000	Transformer Complete Bravo 120V (Gam Board)
56	3	A9BP5970000	Transformer Complete Bravo 220V (Gam Board)
56a	3	41200140000	Transformer Only Bravo 120V (Gam Board)
56a	3	41200080000	Transformer Only Bravo 220V (Gam Board)
56a	3	41200170000	Transformer Only Bravo 120V (Troll Board)
56a	3	41200110000	Transformer Only Bravo 220V (Troll Board)
56b	3	C1BP161000M	Transformer Support Plate Bravo
57	3	43100060000	Pressure Switch Bravo
58	3	C1BP538000M	Pressure Switch Support Bravo
59	3	260000009MN	Pressure Switch/Transducer Quick Connect Fitting Bravo
60	3	43100090000	Pressure Transducer Bravo
61a	3	C5BP5630000	Basic Board (Gam) Bravo 120V/220V
61a	3	C5BM1400200	CPU Board (Troll) Bravo 120V/220V
61b	3	C5BP4350000	PCB Filter Board (GAM) Bravo 120V
61b	3	C5BP1420000	PCB Filter Board (GAM) Bravo 220V 50/60Hz
61b	3	C5BM1420000	PCB Filter Board (Troll) Bravo 120V/220V
61c	3	C5BM1430000	PCB Printer Board Bravo
Not	Shown	SB17V120V	Bravo Software (Gam Board 17V/120V)
Not	Shown	SB17V220V	Bravo Software (Gam Board 17V/220V)
Not	Shown	SB17120V	Bravo Software (Gam Board 17/120V)
Not	Shown	SB17220V	Bravo Software (Gam Board 17/220V)
Not	Shown	SB21V120V	Bravo Software (Gam Board 21V/120V)
Not	Shown	SB21V220V	Bravo Software (Gam Board 21V/220V)
Not	Shown	A9BP5430000	Printer Upgrade Kit (Gam) Bravo
Not	Shown	A9JP0230000	Printer Upgrade Kit (Troll) Bravo
Not	Shown	BRAVO-GT2P120	Bravo Conversion Kit Gam to Troll2 w/Printer 120V
Not	Shown	BRAVO-GT2NP120	Bravo Conversion Kit Gam to Troll2 no Printer 120V
Not	Shown	BRAVO-GT2P220	Bravo Conversion Kit Gam to Troll2 w/Printer 220V
Not	Shown	BRAVO-GT2NP220	Bravo Conversion Kit Gam to Troll2 no Printer 220V

Bravo Parts Table 4




Parts List for Table 4


Item #	Table #	Part #	Description
63	2	28600290000	Filter Fitting Chamber Front Bravo
65	2	47200030000	Filter Screen Chamber Front LP1 Bravo
67	4	48300010000	Gasket Rubber Bravo
68	4	260000006V0	Fitting CPC Female Water Tap Chamber Front Bravo
69	4	A2BP2380000	Printer Plug RS232 9 pin w/Cable Bravo
70	4	C3BP5610000	Plate ABS Bravo
71	4	25600000400	Foot Front PVC Bravo
72	4	25600001400	Foot Rear Plastic Bravo
73	4	55V702016000	Plate Ground Bravo
74	4	C1BP147000M	Support Electronic Angle Bravo
75	4	C1BP269000M	Support Front Top Angle Bar Bravo
76	4	C1BP270000M	Support Front Side Angle Bar Bravo
77	4	C1BP150000M	Chassis Bravo 17 & 17V
77	4	A1BG350000M	Chassis Bravo 21 & 21V
78	4	A0BP2860000	Water Pump Complete w/Solenoid Bravo
78a	4	40000220000	Water Pump 24V Bravo
78b	4	40100370000	Valve EV6 for Water Pump Bravo
79	4	261000009MN	Fitting Elbow for Water Pump Solenoid Bravo
80	4	47200020000	Water Filter Plastic for Water Pump Bravo
81	4	111000006W0	Water Pump Adhesive Mousse Bravo
82	4	25600000800	Water Pump Mounting Bracket Foot Bravo
83	4	40000050000	Water Pump Mounting Bracket Bravo
84	4	A9BM5510000	Steam Generator Complete Bravo 120V
84	4	A9BP5250000	Steam Generator Complete Bravo 220V
84a	4	C1BP136000M	Steam Generator Support Bravo
84b	4	C0BP1310000	Steam Generator Support Spacer Insulator Bravo
84c	4	C4BP2360000	Steam Generator Thermo Insulating Panel End Bravo
84d	4	C4BP5700000	Steam Generator Thermo Insulating Cover Bravo
84e	4	C0BP548000P	Steam Generator Lower Part Bravo
84f	4	48100080000	Steam Generator O-ring Bravo
84g	4	C0BP547000P	Steam Generator Upper Part Bravo
85	4	41000170000	Steam Generator Heating Element Bravo 115V 850W
85	4	41000260000	Steam Generator Heating Element Bravo 220V 1000W
86	4	43000180000	Temperature Sensor PT1000 – PT2 Bravo
87	4	2860014000N	Steam Generator Temperature Sensor Nut Bravo
88	4	28600160000	Steam Generator Temperature Sensor Ferrule Bravo
89	4	43000150000	Steam Generator Safety Thermostat TG400 Bravo
90	4A	A9BM5520000	Vacuum Pump Small Complete Bravo 17/120V
90	4A	A9BM5270000	Vacuum Pump Small Complete Bravo 17/220V
90	4B	40000420000	Vacuum Pump EVO10 Complete Bravo 17V/120V
90	4B	40000360000	Vacuum Pump EVO10 Complete Bravo 17V/220V
90	4B	A9BP5870000	Vacuum Pump 8011 Large Complete Bravo 21V/120V
90	4B	A9BP5270000	Vacuum Pump 8011 Large Complete Bravo 21V/220V

Note: See Table 4a & 4b for Vacuum Pump parts.

Bravo Parts Table 4A



Bravo Parts Table 4B



Parts List for Table 4A

Item #	Table #	Part #	Description
90	4A	A9BM5520000	Vacuum Pump Small Complete Bravo 17/120V
90	4A	A9BM5270000	Vacuum Pump Small Complete Bravo 17/220V
90a	4A	40000500000	Vacuum Pump Small Membrane with Screw Bravo
90b	4A	40000450000	Vacuum Pump Single Head Repair Kit
90b	4A	A8BMR020000	Vacuum Pump Double Head Repair Kit
90c	4A	40000460000	Vacuum Pump Small Complete Head Bravo


Note: Single Head Repair Kit has 4 O-rings & 2 Shutters.

Double Head Repair Kit has 8 O-rings & 4 Shutters.

Parts List for Table 4B

Item #	Table #	Part #	Description
90	4B	40000420000	Vacuum Pump EVO10 Complete Bravo 17V/120V
90	4B	40000360000	Vacuum Pump EVO10 Complete Bravo 17V/220V
90	4B	A9BP5870000	Vacuum Pump 8011 Complete Bravo 21V/120V
90	4B	A9BP5270000	Vacuum Pump 8011 Complete Bravo 21V/220V
90a	4B	40000470000	Vacuum Pump EVO10 Membrane
90a	4B	40000260000	Vacuum Pump 8011 Membrane
90b	4B	40000480000	Vacuum Pump EVO10 Minimum Maintenance Kit
90b	4B	40000250000	Vacuum Pump 8011 Shutter
90c	4B	40000490000	Vacuum Pump EVO10 Full Maintenance Kit
90c	4B	48100150000	Vacuum Pump 8011 O-ring


Bravo Parts Table 5



Parts List for Table 5

Item #	Table #	Part #	Description
80	5	47200020000	Water Filter Plastic Bravo
91	5	C1BP151000M	Chassis Rear Bravo
92a	5	C1BS020000P	Pressure Relief Valve Support Bravo
92b	5	47000020000	Pressure Relief Valve Only Bravo
93	5	C1BP6220000	Heat Exchanger Complete Bravo
93a	5	40400040000	Heat Exchanger Protection Grid Bravo
93b	5	C1BP6220000	Heat Exchanger Only Bravo
93c	5	C1BP159000M	Heat Exchanger Support Bravo
93d	5	40400030000	Heat Exchanger Fan
94	5	40100240000	Valve EV5 3 Way 24V Bravo
95	5	55V432011000	Valve EV5 Support Bravo
96	5	A0BP2820000	Valve EV1 Complete Bravo
97	5	C1BP241000M	Valve EV1 Support Bracket Bravo
100	5	C1BP173000M	Valve EV2/EV3/EV4 Support Bracket Bravo
101	5	40000220000	Water Pump 24V Bravo
102	5	47200020000	Water Filter Plastic for Water Pump Bravo
103	5	01-112545S	Fuse Holder Bravo 120V/220V (Pkg of 2)
104	5	41700330000	Fuse 15A 250V Bravo
105	5	261000014V0	Tee Fitting Bravo
106	5	A2BP2630000	Level Sensor Cable External Bravo
107	5	27700003500	Strain Relief Bravo
108	5	A2BP5890000	Power Cord Bravo 120V
108	5	A2BP6140000	Power Cord Bravo 220V


Bravo Parts Table 6



Parts List for Table 6

Item #	Table #	Part #	Description
79	6	261000009MN	Fitting Elbow Barbed 1/8" Bravo
94a	6	40100240000	Valve EV5 3 ways 24V Bravo
96	6	A0BP2820000	Valve EV1 Complete Bravo
96a	6	40100240000	Valve EV1 3 ways 24V Bravo
98a	6	40100370000	Valve EV2 2 ways 24V Bravo
98b	6	40100230000	Valve EV3 2 ways 24V Bravo
99a	6	40100260000	Valve EV4 3 ways 24V Bravo
110	6	263000002MN	Fitting Tee Female/Female/Female Bravo
111	6	2860034000N	Fitting Pipe Barbed Bravo
113	6	263000003MN	Fitting Tee Female/Male/Female Bravo
114	6	265000009MN	Fitting Coupling Bravo
115	6	262000003MN	Fitting Elbow Female/Male Bravo
116	6	2860034000N	Fitting Elbow Barbed Bravo
117	6	261000004MN	Fitting Pipe Barbed Bravo
118	6	47000010000	Fitting One Way Valve Connector Bravo


Bravo Parts Table 7



Parts List for Table 7

Item #	Table #	Part #	Description
A1	7	A2BP2220000	Cable Microswitch Bravo
A2	7	A2BP4010000	Cable Safety Thermostat Bravo
A3	7	A2BG3940000	Cable Main Switch Bravo
A4	7	A2BM3390000	Cable Vacuum Pump (small) Bravo
A4	7	A2BP4030000	Cable Vacuum Pump (large) Bravo
A5	7	A2BG3520000	Cable Motor Reducer Bravo
A6	7	A2BP2270000	Cable Water Pump Bravo
A7	7	A2BP2280000	Cable Electro Valve Bravo
A8	7	A2BP4070000	Cable Filter Bravo
A9	7	A2BP4820000	Cable Ground Chamber Bravo
A10	7	A2BP4830000	Cable Ground Steam Generator Bravo
A11	7	A2BP4810000	Cable Ground Rear Frame Bravo
A12	7	A2BP4700000	Printer Cable Power Bravo
A13	7	A2BP4690000	Printer Cable Signal Bravo
P1	7	110000005W0	Silicone Tube 6 x 12 Black Bravo
P2	7	110000011W0	Silicone Tube 4 x 7 Transparent Bravo
P3	7	110000003W0	Silicone Tube 6 X 10 w/Quick Connect for Drain Bravo
P4	7	110000002V0	Teflon Tube 4 X 2.5 Bravo
P5	7	110000018W0	Silicone Tube 6 X 10 Bravo (Steam Generator to Chamber)

Bravo Fuses



Basic board: GAM version

F1: T 5A 250V (secondary transformer winding)
 F2: T 4A 250V (120V Bravo)(Primary transformer winding)
 F2: T 2A 250V (220V Bravo)(Primary transformer winding)
 F3: T 200mA 250V (door-lock accidental activation)
 F4: F 1.25A 250V (door-lock overload)

part # Fuse-T-5A
 part # Fuse-T-4A
 part # Fuse-TT-2A
 part # Fuse-F-200MA
 part # Fuse-F-1.25A

Basic board: TROLL version

F1: T 6.3A 250V (secondary transformer winding)
 F2: T 3.15A 250V (primary transformer winding)

part # Fuse-T-6.3A
 part # Fuse-T-3.15A

Printer PS board

F1 PTR: T 5A 250V (printer protection)

part # Fuse-T-5A

Line Fuses 15A 250V ABC (Fast Acting)(not shown)

part # 41700330000

Bravo Individual Parts Pictures



Bravo ABS Plate
#C3BP5610000



Bravo Basic Board (Gam)
#C5BP5630000
Bravo Basic Board (Troll)
C5BM1400200



Bravo Cable Flat 9 Pin
#A2BP2380000



Bravo Calibration Board
#STXX0120000



Bravo Chamber Heater
Safety Thermostat TG400
#43000140000



Bravo Chamber Thermo-
Insulating Cover 17
#C4BP2140100
Bravo Chamber Thermo-
Insulating Kit 17
#A8BPR140000



Bravo Chamber Thermo-
Insulating Cover 21
#C4BG3280000
Bravo Chamber Thermo-
Insulating Kit 21
#A8BGR010000




Bravo Chamber Thermo-
Insulating Panel Rear
#C4BP5570000



Bravo Door Fork
#C0BP060000Z



Bravo Door Fork Pin
#C0BP0500000



Bravo Door Locking
System Complete
#A1BP2710000



Bravo Door Locking
System Pin 3 x 25
#232A03L25K0



Bravo Door Motor
Complete
#A0BP6250000



Bravo Door Plate
Adjustment Screw
#C0BP079000N



Bravo Door Thermo-
Insulating Plate
#C4BP2740000



Bravo Drain Tubing
With Fitting 6 x 10
#110000003W0

Bravo Individual Parts Pictures



Bravo Facia Front Plastic
(Original Style)
#C3BP1830000



Bravo Fitting for
Chamber Filter
#28600290000



Bravo Fitting CPC
Female Water
#260000006V0



Bravo Float Switch Clean
Water Reservoir Min-Max
#43200050000



Bravo Float Switch Waste
Water Tank
#43200060000



Bravo Foot Front PVC
#25600000400



Bravo Foot Rear Plastic
#25600001400



Bravo Heat Exchanger
Fan
#40400030000



Bravo Keypad
#C6JP0010000



Bravo Label
Bravo 17
#C6JM0010000



Bravo Label
Bravo 17V
#C6JP0020000



Bravo Label
Bravo 21V
#C6JG0010000



Bravo PCB Filter
Board GAM 120V
#C5BP4350000



Bravo PCB Filter
Board GAM 220V
#C5BP1420000



Bravo Plate Ground
#55V702016000



Bravo Plug Plastic
#C3BP2790000

Bravo Individual Parts Pictures



Bravo Power Switch
#42000070000



Bravo Pressure Relief
Valve Group
#A0BS0520000



Bravo Pressure Relief
Valve Only
#47000020000



Bravo Pressure Switch
#43100060000



Bravo Pressure Switch
Transducer Fitting
#260000009MN



Bravo Printer Cable Power
#A2BP4700000



Bravo Printer Cable Signal
#A2BP4690000



Bravo Printer Door
(Original Style)
#C3BP1840000



Bravo Printer Door
Pin 4 x 20
#230B04L20K0



Bravo Printer Paper
Thermal 57mm
#STXX0250000



Bravo Printer Thermal
#C5BP5500000



Bravo Reservoir Fill
Bottle 2500ml
#A0XP0010000



Bravo Vacuum Pump 8011
Shutter
#40000250000



Bravo Steam Generator Heating
Element 115V 850W
#41000170000



Bravo Steam Generator
Heating Element 220V 1000W
#41000260000



Bravo Steam Generator
O-ring
#48100080000

Bravo Individual Parts Pictures



Bravo Steam Generator
Lower Part
#C0BP548000P



Bravo Steam Generator
Safety Thermostat TG400
#43000150000



Bravo Steam Generator
Temperature Sensor Nut
#2860014000N



Bravo Steam Generator
Temperature Sensor Ferrule
#28600160000



Bravo Steam Generator
Upper Part
#C0BP547000P



Bravo Temperature
Sensor PT1000-PT1
#43000190000



Bravo Temperature Sensor
PT1000-PT2/PT3/PT4
#43000180000



Bravo Transformer Only
GAM 120V
#41200140000



Bravo Transformer Only
GAM 220V
#41200080000



Bravo Tray Handle
#STXX0080000



Bravo Vacuum Pump 8011
Large 21V/120V
#A9BP5870000



Bravo Vacuum Pump 8011
Large 21V/220V
#A9BP5270000



Bravo Vacuum Pump
Small Membrane
w/Screw
#40000500000



Bravo Vacuum Pump
Small Shutter
#40000180000



Bravo Valve EV1 & EV5
Parker 3ways 24V
#40100240000

Bravo Individual Parts Pictures



Bravo Valve EV1
Elbow Fitting
#261000009MN



Bravo Valve EV1 Support
#C1BP241000M



Bravo Valve EV2
Parker 2way 24V
#4010037000



Bravo Valve EV3
Parker 2way 24V
#40100230000



Bravo Valve EV4
Parker 3ways 24V
#40100260000



Bravo Valve EV5 One Way
Valve Connector
#47000010000




Bravo Water Pump Group
With Solenoid
#A0BP2860000




Bravo Water Pump
Mounting Bracket Foot
#25600000800


**Bravo Water Intake Pull Kit
#01-111774S**



Bravo Direct to Drain Kit
#01-111775S




Quantim Front View Door Closed



Item #	Part #	Description
1	289124	Door Molding Kit
2	279159	Door Microswitch 2 Wire (Not Shown Under Door Molding)
2	279494	Door Microswitch 3 Wire (Not Shown Under Door Molding)
3	289162	Decal Quantim 16 (New Style no Liquids Cycle Not Pictured)
4	289125	Printer Door with Hinges


Note: Door Microswitch 2 Wire, wires are connected to switch terminals.
Door Microswitch 3 Wire, switch terminals are connected to door PCB.

Quantim Front View Door Open




Item #	Part #	Description
1	279011	Door Gasket (Includes Door Plate Seals)
2	279490	Acorn Nut (Pkg 2) (Pictured Above)
2	279491	Dome Nut (Pkg 2) (Not Shown in Picture)
3	279107	O-ring Kit Entry Port White (Pkg 2)
3	279467	Bolt & O-ring Entry Port
4	279527	Tubing & Fitting Kit
5	279102	Drain Tap
6	289132	Rack for Chamber 3 Tray (Wire)
7	309077	Air Filter
8	289006	Temperature Probe Chamber

Quantum Rear View with Cover




Item #	Part #	Description
1	319003	Cover w/Feet
2	319001	Power Cord
3	319014	Fan Assembly Rear
4	279469	Cover Circular
5	279503	Waste Bottle Complete Quantum 16 (Bottle Not Shown)

Quantum Right Side No Cover




Item #	Part #	Description
1	289125	Printer Door with Hinges
2	279537	Front Bezel
3	279170	Door Opening Solenoid
4	309039	Reservoir Single Use (Quantum with Waste Bottle)
4	279117	Reservoir Multi Use (Quantum no Waste Bottle)
5	309077	Air Filter & Clip
6	299051	Air Filter Nozzle
7	289134	Silencer Vacuum Pump
8	289118	Vacuum Pump

Quantum Right Side Front View




Item #	Part #	Description
1	279225	Display Module
2	279171	Door Solenoid Spring Kit (2 Springs)
3	279170	Door Opening Solenoid

Quantum Rear View No Cover




Item #	Part #	Description
1	309039	Reservoir Single Use (Quantum with Waste Bottle)
1	279117	Reservoir Multi Use (Quantum no Waste Bottle)
2	309077	Air Filter & Clip
3	319015	Fuse 15A
4	319001	Power Cord
5	319006	Control Module

Quantum Left Side Lower Front No Cover




Item #	Part #	Description
1	289011	Temperature Probe Boiler
2	279155	Insulator Heating Element
3	289136	Transformer
4	279544	Boiler Thermistor

Quantum Left Side Lower Back No Cover




Item #	Part #	Description
1	279157	ASME Safety Valve***
2	279155	Insulator Heating Element
***Note: Boiler must be removed to replace ASME valve. A new boiler to vessel seal (#319022) will need to be installed when reinstalling boiler.		

Quantim Left Side Top Front No Cover




Item #	Part #	Description
1	289006	Temperature Probe Chamber
2	289015	Door Opening Spring
3	289028	Door Spring Pillars (Qty 2)

Quantum Bottom View



Item #	Part #	Description
1	319019	Valve Block Kit
2	279176	Water Filter
3	319022	Seal Boiler to Vessel (Not Shown in Picture)
4	319024	Boiler w/PRV
5	279266	Thermal Cutout
6	279203	Thermal Cutout Mounting Kit
7	279155	Insulator Heating Element
8	289136	Transformer
9	289117	Heater Pad (Not Shown in Picture)

Hydrim C51W Parts Front View Door Closed







Item #	Part #	Description
1	01-107974S	Door Facia
2	01-108790S	Door Actuator Kit
3	01-108917S	Facia Control Panel
3	01-109147S	Facia Control Panel w/Dryer (HD01FK0001 & higher)
4	01-107935S	Chemical Door w/Hinge
5	01-107972S	Cover, Top (no vents)
5	01-110205S	Cover, Top (with vents)
5	01-112016S	Cover, Top (serial #210110C00001 & higher)
Not Pictured	01-108924S	Screw Kit (Cover Screws)

Hydrim C51W Parts Front View Door Open





Item #	Part #	Description
1	01-107805S	Wash Arm Upper
2	01-108154S	Plug Upper Arm
3	01-107786S	Seal Door
4	01-107804S	Wash Arm Lower
5	01-107806S	Drain Screen
6	01-107812S	Door Latch
7	01-107935S	Chemical Door New Style
Not Pictured	01-108090S	Seal Kit Door Bottom
Not Pictured	01-107807S	Drain Filter

Hydrim C51W Parts Right Side




Item #	Part #	Description
1	01-108030S	Quick Connect Female
2	01-108031S	Tube Chemical (HD01DF9999 & lower)
2	01-108699S	Tube Chemical (HD01DG0001 & higher)
3	01-108122S	Drip Tray
4	01-109058S	Dosing Reservoir Upgrade Kit (Old Style Reservoir not available must upgrade)
5	01-107790S	Dosing Pump
5	01-111848S	Dosing Pump (serial #210110C00001 & higher)
6	01-108700S	Fuse Holders (Pkg. 3)
7	01-108699S	Tube Chemical (HD01DG0001 & higher)
8	01-108253S	Dosing Reservoir New Style
9	01-109909S	Tube Dosing Pump Internal
10	01-108351S	Float Dosing Reservoir
11	01-107938S	Switch Pressure Assembly
12	01-107803S	Switch Overflow Chamber
13	01-107937S	Switch Overflow Bottom
Not Pictured	01-107800S	Fuse 10 Amp (Pkg. 10)
Not Pictured	01-107799S	Fuse 2 Amp (Pkg. 5)

**Hydrim C51W Right Side View & Inside Chamber with Dryer
HD01FK0001 & Newer**



Item #	Part #	Description
1	01-109143S	Dryer Tubing
2	01-109142S	Dryer Assembly
3	01-109145S	Vent Assembly Dryer
4	01-109144S	Dryer Fitting

Hydrim C51W Parts Bottom



Item #	Part #	Description
1	01-107794S	Pump Recirculation Lower Arm
2	01-108701S	Tubing (blue) Chamber
3	01-107797S	Pump Drain
4	01-107808S	Water Heater
5	01-107815S	Valve Hot Water
6	01-107801S	Valve Cold Water
7	01-107809S	Thermostat Heater NC Sensor
7	01-110304S	Thermostat Heater for Touchscreen
8	01-107795S	Pump Recirculation Upper Arm
9	01-107810S	Pressure Switch Heater

Note: See next page for upper and lower pump impellers and upper pump Y fittings with stainless steel barbs.

Hydrim C51w Spare Parts Update

06-TSB-230

Please be advised of the new Hydrim C51w spare parts.


The Impellers (items # 4&5) have been set-up separately due to cost concerns of changing circulation pumps. The Y-fittings have been upgraded and are now made of stainless steel. The part numbers for the items have remained the same.



1. Stainless Steel Fitting For Y's
p/n 01-109590S (Package of 2)

2. Y-Fitting, Upper Arm Pump
p/n 01-109303S

3. Y- Fitting, Upper Arm
p/n 01-108702S




4. Impeller - Upper
p/n 01-109447S

5. Impeller - lower
p/n 01-109448S

Hydrim C51W Bur Screen
Part #01-111828S



Hydrim C51W w/Touchscreen
Front View Door Closed



Item #	Part #	Description
1	01-110205S	Cover, Top (with vents)
1	01-112016S	Cover, Top (serial #210110C00001 & higher)
2	01-110201S	Door Facia w/Label
3	01-110202S	Facia Controller w/Label
4	01-110200S	Detergent Door w/Label
5	01-109785S	Touch Display
6	01-100573S	Power Switch
Not Visible	01-108924S	Screw Kit (Cover Screws)
Not Visible	01-107975S	Adjustable Feet (Pkg of 4)

**Hydrim C51W w/Touchscreen
Front View Door Open with 4 Cassette Rack**




Item #	Part #	Description
1	01-109894S	Trolley
2	01-108309S	Kickplate
Not Visible	01-109153	4 Cassette Rack

**Hydrim C51W w/Touchscreen
Front View Door Open with Baskets**



Item #	Part #	Description
1	01-109894S	Trolley
2	01-107240	Basket Statim 2000
Not Shown	01-107241	Basket Statim 5000
Not Shown	01-108241	Support Rack (hooks to bottom of trolley holds 2 additional Trays)

Hydrim C51W w/Touchscreen Rear View




Item #	Part #	Description
1	01-107801S	Inlet Valve Cold Water
2	01-107815S	Inlet Valve Hot Water
3	01-109305S	Power Cord
4	01-108700S	Fuse Holder (Pkg of 3)
5	01-107800S	Fuse 10A (Pkg of 10)
6	01-107799S	Fuse 2A (Pkg of 5)
Not Shown	01-107789S	Drain Hose
Not Shown	01-107788S	Inlet Hose (Hot or Cold)

Hydrim C51W w/Touchscreen Front View Door Open No Trolley




Item #	Part #	Description
1	01-111700S	Door Latch Retro-Fit Kit (Includes Door Latch Assy. in Unit & Latch Arm in Door)
2	01-109790S	Upper Spray Arm
3	01-109145S	Dryer Vent Assy. (Pkg of 2)
4	01-107786S	Door Seal
Not Visible	01-108090S	Door Seal Kit Bottom
5	01-107804S	Lower Wash Arm
6	01-107806S	Screen Drain
Not Visible	01-108310S	Mesh Drain Screen (small round mesh screen fits in center of Screen Drain)
Not Visible	01-107807S	Filter Drain (sits under Screen Drain)
Not Visible	01-111828S	Screen Assy. Sump
7	01-109144S	Dryer Fitting
8	01-108795S	Plastic Trim Edge

Hydrim C51W w/Touchscreen
Right Side View No Cover




Item #	Part #	Description
1	01-109789S	Rope Door Latch
2	01-108030S	Quick Connect Female
3	01-108699S	Tube Chemical
4	01-108122S	Drip Tray
5	01-108253S	Reservoir, Dosing
6	01-109143S	Dryer Tubing
7	01-109835S	Dryer Assy.
8	01-107790S	Dosing Pump
8	01-111848S	Dosing Pump (Serial #210110C00001 & higher)
Not Visible	01-109909S	Tube Dosing Pump Internal (Black)

Hydrim C51W w/Touchscreen
Right Side View Behind Chemical Bracket




Item #	Part #	Description
1	01-107791S	Valve Salt Regeneration
2	01-107938S	Full Switch Assy.
3	01-107803S	Overflow Switch
4	01-111619S	Styrofoam Float Overflow Bottom

**Hydrim C51W w/Touchscreen
Left Side View No Cover**




Item #	Part #	Description
1	16943.06	Joint Left Hand Side
Not Shown	16943.05	Joint Right Hand Side
2	16943.04	Frame Left Hand Side
Not Shown	16943.03	Frame Right Hand Side
3	01-107794S	Recirculation Pump Lower Arm
4	01-108797S	Support Bracket Lower Pump
5	01-107808S	Water Heater

**Hydrim C51W w/Touchscreen
Top View No Cover**




Item #	Part #	Description
1	01-109776S	Top Arm Fitting
2	01-110204S	Tubing Upper Arm Fitting
3	01-109143S	Dryer Tubing
4	01-109780S	Door Solenoid
5	01-111700S	Door Latch Retro-Fit Kit (Includes Door Latch Assy. in Unit & Latch Arm in Door)

Hydrim C51W w/Touchscreen View Behind Controller Facia



Item #	Part #	Description
1	01-109782S	Power Supply Board
2	01-110206S	Logic Board C51W
3	01-109783S	I/O Board
Not Visible	01-109785S	Touch Display

**Hydrim C51W w/Touchscreen
Bottom View**



Item #	Part #	Description
1	01-111619S	Styrofoam Float Overflow Bottom
2	01-110259S	Pump Upper Arm without Fitting
3	01-107801S	Inlet Valve Cold Water
4	01-107815S	Inlet Valve Hot Water
5	01-108700S	Fuse Holder (Pkg of 3)
6	01-107808S	Water Heater
7	01-107794S	Recirculation Pump Lower Arm
8	01-110230S	Single Temperature Sensor C51W
9	01-107797S	Waste Pump
10	01-107809S	NC Sensor (Water Heater)

Hydrim C51W w/Touchscreen Parts



01-107786S
Seal Door, J



01-107788S
Inlet Hose N.A., J



01-107789S
Drain Hose, J



01-107790S
Dosing Pump, J

01-111848S (S#210110C0001)



01-107791S
Valve Salt Regeneration, J



01-107794S
Pump Recirculation Lower Arm, J



01-110259S
Pump Upper Arm w/o
Fitting, J



01-107797S
Drain Pump N.A. 230V, 60Hz, J



01-107799S
2A Fuse, J



01-107800S
10A Fuse, J



01-107801S
Inlet Valve Cold Water, J



01-107803S
Switch Overflow Chamber, J



01-107804S
Lower Wash Arm, J



01-107806S
Screen Drain, J



01-107807S
Filter Drain, J



01-107808S
Water Heater, J



01-107810S
Switch Pressure Heater, J



01-107812S
Clip Door, J



01-107815S
Inlet Valve Hot Water, J



01-107938S
Switch, Pressure Assembly
(Chamber Full Switch Assembly)

Hydrim C51W w/Touchscreen Parts



01-107802S
Switch Full Chamber, J



01-107975S
Adjustable Feet, J



01-108030S
Quick Connect Female, J



01-108699S
Tube, chemical



01-108121S
Cap Quick Disconnect, J



01-108122S
Drip Tray, J



01-108253S
Reservoir



01-108309S
Kickplate, J



01-108310S
Mesh Drain, J



01-108700S
Fuse Holder, J



01-108795S
Plastic Trim Edge, J



01-108797S
Support Bracket 1-st Pump, J



01-109143S
Dryer Tubing, J



01-109144S
Dryer Fitting, J



01-109145S
Dryer Vent Assy, J



01-109776S
Top Arm Fitting Hydrim, S



01-110204S
Tubing, Upper Arm Fitting, C51W, J


Hydrim C51W w/Touchscreen Parts



01-109780S
Door Solenoid C51WD, S




01-110202S
Fascia Controller w/
Label, C51W, J



01-109782S
Power Supply Board, S



01-110206S
Logic Board Hydrim, C51W, J



01-109783S
I/O Board Hydrim, S



01-109785S
Touch Display Hydrim, S



01-110203S
Single Temperature
Sensor C51W, J



01-109789S
Rope, Door Latch, C51WD, S



01-109790S
Upper Spray Arm
Hydrim, S



01-110201S
Fascia Door w. Label, C51W, J



01-110200S
Detergent Door w. Label
C51W, J



01-107809S
NC Sensor



01-110205S
Cover Top, C51W, J



01-107973S
Packaging Hydrim
C51W, S




01-109834S
Seals, Door Bottom Hydrim, S



01-109835S
Dryer Assy Hydrim, S


01-112016S
(S#210110C00001)

Hydrim L110W Solid Door Parts Front View Door Closed




Item #	Part#	Description
1	50750.10	Cover Top
2	13790.06	Decal Front (Set of 2)
3	13459.00	Power Switch
3	13471.01	Button ON/OFF Silver
4	16945.00	Handle Door Lock
5	16979.01	Cover Front Door Complete
6	16648.01	Kick Plate Front
7	16695.00	Cover Right Side
8	16662.06	Door Cleaning Solution
9	16663.06	Cover Slave Controller (No Decal)
10	13471.02	Button Program Silver
11	13476.25	Front Facia Complete with Decal
Not Pictured	16696.00	Cover Left Side
Not Pictured	16656.00	Cover Back
Not Pictured	13633.00	Screw Kit (Cover Screws)

Hydrim L110W Solid Door Parts Front View Door Open




Item #	Part#	Description
1	13492.00	Wash Arm Middle
2	13456.00	Seal Peripheral Door Seal
3	13638.00	Lid for Salt Container
4	13481.00	Wash Arm Bottom
5	13611.06	Door Inner Stainless Steel
6	13466.00	Filter Bottom Flat
7	13465.00	Filter Bottom Round
8	13629.00	Pipe Water Feed to Wash Arm

Hydrim L110W Solid Door Parts Left Side No Cover





Item #	Part#	Description
1	13462.00	Side Bag
2	13455.00	Drain Tubing
3	13614.00	Solenoid Valve Drain
4	13615.00	Solenoid Valve Regeneration
5	13447.01	Emitter Housing Red
Not Pictured	13447.00	Emitter Housing Clear

Hydrim L110W Solid Door Parts Right Side No Cover




Item #	Part #	Description
1	16590.01	Dryer Unit Complete
2	13319.10	Dryer Air Filter
3	16590.10	Dryer Air Turbine Kit (to replace turbine #18466.00)
3	13220.00	Dryer Air Turbine (for Hydrims with kit #16590.10 installed)
4	13416.01	Check Valve Dryer 1"
5	01-111568S	Flow Switch w/Bracket
6	13599.00	Foot Front
7	21624.02	Dosing Pump Clamp External
8	13634.00	Dosing Pump
9	21970.00	Dosing Pump Hose
10	21027.00	Tubing Exhaust Air


Hydrim L110W Solid Door Parts Inlet Solenoids & Controllers



Item #	Part#	Description
1	N/A	Valve Solenoid Inlet***
***Note: If Solenoid Inlet Valve is defective install Preventative Maintenance Kit #01-109707S see next page.		




Item #	Part#	Description
1	16945.00	Handle Door Lock
2	13475.03	Controller Master



Item #	Part#	Description
1	13475.04	Controller Slave
2	23126.01	Relay Dryer

Hydrim L110W Solid Door Preventative Maintenance Kit

Part #01-109707S



Item #	Description
1	Inlet tubing's for hot and cold water (01-109861S)
2	Tubing assembly with sensor from inlet valve block to sidebag
3	Water inlet valve block
4	Vent tubing for top of sidebag
5	Tubing's from flow switch to dosing pump and dosing pump to chamber
6	Flow switch
7	Flow switch mounting bracket
8	Wiring harness for new valve block
9	Plate for water inlet valve block
10	Screws to mount valve block plate
11	Hardware to mount flow switch
12	Tie wraps for wiring harness

Hydrim L110W w/Glass Door Front View




Item #	Part #	Description
1	01-111469S	Cover Top
2	01-111484S	Screw 3mm Hex Head M5x16 Back Cover/Controller Door
3	01-110219S	Decal Kit
4	01-112024S	Rocker (Power) Switch
5	01-111470S	Door Cleaning Solution
6	01-111468S	Cover Right Side
7	01-111485S	Screw Kickplate M5x12 Slot Round Head
8	01-111476S	Kickplate Front
9	01-111466S	Cover Left Side
10	01-111483S	Screw 2.5mm Hex Head M4x8 Top Cover
Not Shown	01-111467S	Cover Rear
Not Shown	01-111477S	Kickplate Rear

Hydrim L110W w/Glass Door Chamber with Trolleys




Item #	Part #	Description
1	01-111491S	Trolley Upper
2	01-109963S	5 Cassette Rack Full Size
3	01-111490S	Trolley Lower

Hydrim L110W w/Glass Door Chamber




Item #	Part #	Description
1	01-111559S	Door Seal Kit
2	01-109790S	Wash Arm Upper & Lower Metal
3	01-111495S	Wash Arm Middle Complete
4	01-109144S	Dryer Fitting
5	13448.22	Salt Cap
6	01-111589S	Mesh Filter Wash Chamber
7	01-111705S	Plugs for Rinse Valve

Hydrim L110W w/Glass Door Control Boards




Item #	Part #	Description
1	01-109783S	I/O Board
2	01-110217S	Logic Board
3	01-109785S	Touch Display

**Hydrim L110W w/ Glass Door
Behind Front Kickplate Left Side**




Item #	Part #	Description
1	22815.03	Switch Safety (behind plate)
2	01-112080S	Door Latch Cord Assembly
3	01-111702S	Sensor, Chamber Temperature
4	01-111489S	Temperature Sensor PT-100

Behind Front Kickplate Right Side




Item #	Part #	Description
1	01-111471S	Dosing Pump
Not Shown	01-109909S	Tube Dosing Pump Internal

Hydrim L110W Right Side Top View




Item #	Part #	Description
1	01-109782S	Power Supply Board
2	01-111408S	Switch Full (Level Control)
3	01-109142S	Dryer Assembly
4	01-111474S	Flow Switch

Hydrim L110W w/Glass Door Right Side Bottom View




Item #	Part #	Description
1	01-111474S	Flow Switch
2	01-111580S	Tension Door Spring (light door no cutouts)
2	18156.00	Door Tension Spring (heavy door with half moon cutouts)
2	01-112729S	Dual Spring Retro-fit Kit
2	01-112730S	Dual Springs (replacement springs only)
3	01-111412S	Drain Pump
4	01-111494S	Valve Single Hot Water
5	01-111493S	Valve Double Cold Water
6	01-111473S	EMI Filter Kit
Not Shown	01-112026S	Check Valve Drain

Hydrim L110W w/Glass Door Rear View No Cover





Item #	Part #	Description
1	01-111409S	Switch Pressure
2	01-111585S	Breather
3	01-109142S	Dryer Assembly
4	01-111492S	Drain Tubing
5	01-107788S	Water Inlet Hose
6	01-111480S	Power Cord
7	01-111465S	Circulation Pump

Hydrim L110W w/Glass Door Top View No Cover




Item #	Part #	Description
1	01-111478S	Nut Top Cover for 2.5mm Screw
2	01-111783S	Door Latch Assembly
3	01-109782S	Power Supply Board
4	01-111408S	Switch Full (Level Control)
5	01-109142S	Dryer Assembly
6	01-111585S	Breather
7	01-111409S	Switch Pressure
8	01-112080S	Door Latch Cord Assembly

Statmatic Parts





Item #	Part #	Description
1	S30153	Cover, Top with keypad
2	S30123	Keypad
3	S30141	Spray Can Black Plastic Adapter (no brass sliding lock ring)
4	S30143	Door Collet





Item #	Part #	Description
5	S30113	Power Inlet Assembly
6	S30112	Regulator, Air
7	S30102	Absorption pad (6 large pads & 2 bottom pads under large pad)
8	S30120	Seal Kit Front Door

Statmatic Parts





Item #	Part #	Description
9	S30103	Hinge (for front door, no nuts or screws)
10	S30108	Rubber Foot





Item #	Part #	Description
11	S30101-CX	Coupling Multiflex/Connex (for SciCan & Kavo high speeds)
12	S30101-E	Coupling E-Type (for SciCan & Kavo low speed contra angles)
13	S30109	Chuck Adapter Complete
14	S30118	O-ring Intramatic (o-ring set for E-type coupler)
15	S30144	Clamping Nut
16	S30150	Recognition Flap
17	S30116	O-ring Set Multiflex/Connex Coupling

Statmatic Parts




Item #	Part #	Description
18	S30115	Fuse Drawer
19	S30117	Fuse 1A





Item #	Part #	Description
20	S30107	Comp. Air Cylinder (air solenoid for solution bottle)
21	S30105	Solenoid Valve Block
22	S30130	Power Cable (from power inlet to power supply board)
23	S30129	Cable PE (ground wire for power inlet assembly)
24	S30132	Tubing Blue (2" piece of blue air line tubing)
25	S30121	Control Board
26	S30119	Power Supply Board
27	S30151	Coupler Receptacle with 5 Hole Configuration
27	S30140	Connection Bush (Coupler receptacle with 4 hole configuration)
28	S30106	Microswitch Assembly (3 Switches with Wiring)
29	S30125	Tubing (clear 10" tubing)
30	S30110	Hose Set (complete tubing assembly including "Y" fittings)
31	S30139	Hex Nut 12 X 1 (hex-nut which holds the coupler receptacles)

Statmatic Parts



Item #	Part #	Description
32	S30107	Air Solenoid for Solution Bottle
33	S30127	Tubing Labeled (4" yellow tubing labeled "oil")
34	S30145	Sensor Collet Complete (with oil tubes, optical sensor & wiring)
34	S30148	Sensor Collet (with oil tubes & brass elbow fitting)
35	S30104	Fluid Sensor (optical fluid sensor)
36	S30137	Clamp Bugel (wire fastener the holds cable to fluid sensor)
37	S30138	Cable Photo Sensor (cable from fluid sensor to main PCB)
38	S30126	Tubing Labeled (10" yellow tubing labeled "oil")
39	S30133	Tubing Switch (door magnet sensor and cable)
40	S30131	Cable SNT/Elect. (cable from power supply board to main PCB)

Statmatic Parts



Item #	Part #	Description
41	S30124	Spring (for microswitch activating)
42	S30142	Guide Set (plunger assembly for solution bottle)

Statmatic Couplers



S30101-BA
BienAir



S30101-CX
Multiflex/Connex



S30101-E
E-type



S30101-KVH
Kavo Intra head



S30101-MW
Midwest Slow speed



S30101-NSK
Coupling NSK



S30101-ST1I
Sirona T1 Slow speed
Attachments




S30101-ST1T
Sirona T1 High speed




S30101-W&H
Coupling W&H

Aquastat Boiler




Item #	Part #	Description
1	P492-2A	Boiling Chamber Lid (No Gasket)
2	P167-72	Gasket for Boiling Chamber Lid
3	KP05EC001	Baffle Nut
4	KP03WS010	Baffle without O-ring
Not Shown	KP32GA162	Baffle Gasket
Not Shown	KP32GA163	Baffle O-ring
5	P186-1X	Boiler Vessel with Heater Only
6	KP28TS091	Thermostat
7	KP64CO001U	Control Board
8	P338-26	Strix Connector Boiler
9	KP44WP003	Boiler Base Bottom with Grey Feet
10	P128-134	Spring
11	P81-174	Button
12	P118-96	Retainer
13	P19-283	Grey Feet

Aquastat Chassis



Item #	Part #	Description
17	P8-34	Vent Piece for Top Cover
18	P72-8	Cover Top without Vent Piece
19	KP76PL001	Fan Blade Clear
20	KP63MT019	Motor 110V
21	P450-42	Clip Thermostat
22	P185-97	Thermostat Condensing Coil
23	P3-192	Elbow
24	P26-93	Hose without Vent Hole
25	P84-126	Condensing Coil
26	P26-93D	Hose with Vent Hole
27		Item #27, 29, 32, 34 See 35 Chassis Kit
28	KP03WS015	Spout
30	P338-26	Strix Connector Boiler
30	P314-26	Otter Connector Boiler
31	P16-78	Strix Connector Chassis
31	P13-78	Otter Connector Chassis
33	P118-14	Power Socket Main Frame
35	P35-197	Chassis Kit Strix (oval opening)

Aquastat Reservoir



Item #	Part #	Description
35	P160-89M	Reservoir Cap
36	W50211	Replacement Reservoir
37	P96-104	Faucet for Reservoir
38	W8010	Carbon Filters (6 Per Pack)

Note: Carbon Filters should be used for making drinking water only

Aquastat Parts



Baffle Gasket
#KP32GA162



Baffle Nut
#KP05EC001



Baffle O-ring
#KP32GA163



Baffle without O-ring
#KP03WS010



Boiler Base Bottom
KP44WP005



Boiling Chamber Lid
(no gasket)
#P492-2A



**Boiling Chamber New
Style
Oval Connector**
#P206-1H



**Boiling Vessel
with Heater Only**
#P186-1X



Button
#P81-174



Carbon Filters 6 per Pack
#W8010



Chassis Kit
#P35-197



Clip Thermostat
#P450-42

Aquastat Parts




Condensing Coil
#P84-126



Control Board
#KP64CO001U



**Cover Top without
Vent Piece**
#P72-8



**Distiller Cleaner &
Descaler**
#D1-500200S



Distiller Wiring Kit
#P376-90K



Elbow
#P3-192



Fan Blade Clear
#KP76PL001




Faucet for Reservoir
#P96-104



Feet Gray
#P19-283



**Gasket for Boiling
Chamber Lid**
#P167-72



Hose with Vent Hole
#P26-93D



Hose without Vent Hole
#P26-93

Aquastat Parts



Motor, 110V
#KP63MT019



Nut
#P340-81



Otter Connector Boiler
#P314-26



Otter Connector Chassis
#P13-78



Owners Manual
#L5446



Packaging
#J3366K



Power Cord
#P313-74



**Power Socket Main
Frame**
#P118-14



Replacement Reservoir
#W50211



Reservoir Cap
#P160-89M



Retainer
#P118-96



Screw
#KP13MS570

Aquastat Parts



Spacer
#KP12SN254



Spout
#KP03WS013



Spring
#P128-134



Strix Connector Boiler
#P338-26



Strix Connector Chassis
#P16-78



Thermostat Otter
#KP28TS091



Thermostat Strix
#P197-97




**Thermostat Condensing
Coil**
#P185-97



Vent Piece for Top Cover
#P8-34

VistaBrite Part Identification




Name	Part #	Description
Supply Hose Connectors	S4231	Garden hose adaptor – Female hose x 1/4" tube
VistaCheck™ Backflow Preventer	VC250	Backflow check valve - In-line 1/4" FNPT
System Check Valve	S6265P	In-line check valve – 1/4" tube
System Regulator	S9130	Automatic shut-off valve - Hydraulic
Temperature Gauge	S6435	Temp gauge – 1/4" MNPT back mount (0-250 F/C)
Filter Sump Brackets	S9170	Filter mounting bracket
Filter Sump Heads	S9171	Filter head
Inlet Pressure Gauge	S6431	Pressure gauge – Back mount 1/4" MNPT (0-100 PSI)
Storage Tank Pressure Gauge	S6431	Pressure gauge – Back mount 1/4" MNPT (0-100 PSI)
System Blending Manifold	S4310	System manifold ONLY (does not include fittings)
Carbon Prefilter	R4151	Replacement carbon prefilter (every 6 months)
Hyper Prefilter	R4155	Replacement hyperfiltration filter (approx 2 yrs)*
Drain Line Flow Restrictor	S9165	Waste water flow restrictor – 1/4" tube (1500 ml/min)
Sump Drain Elbow	S6121P	Male adaptor elbow – 1/8" NPT x 1/4" tube
Drain Saddle	S9175	Drain saddle – 1/4" tube
Valved Hose Wye	S9180	Garden hose valved wye
Male-Male Hose Adaptor	S9181	Garden hose male to male adaptor
Female Hose Adaptor	S4230	Garden hose adaptor – Female hose x 3/8" tube
Storage Tank (Standard 7 GAL)	S9111	Pressure storage tank w/ stand – 7 gallon nominal
Storage Tank (Optional 14 GAL)	S9113	Pressure storage tank w/ stand – 14 gallon nominal
Tank Valve Assembly	*****	(see top-down view on following page)

* Hyperfiltration filter should be changed every 2 years, or when TDS reduction is less than 70%

VistaBrite Part Identification




(Storage Tank Top-Down View)



Name	Part #	Description
Storage Tank (Standard 7 GAL)	S9111	Pressure storage tank w/ stand – 7 gallon nominal
Storage Tank (Optional 14 GAL)	S9113	Pressure storage tank w/ stand – 14 gallon nominal
Tank Valve	S6160P	Elbow shut-off valve – 1/4" tube x 1/4" FNPT
1/4" x 3/8" Tube Union	S6310P	Union – 1/4" tube x 3/8" tube
1/4" Cross	S6175P	Cross – 1/1" tube
Test Cock Valve	S6155P	Valve – In-line – 1/4" tube

VistaBrite Part Identification

(continued)



Item #	Part Number	Part Description
1	S6112P	Swivel Plug-In Elbow – 1/4" tube x 1/4" stem
2	S6110P	Elbow Union – 1/4" tube x 1/4" tube
3	S9130	Automatic Shut-Off Valve – Hydraulic / (Regulator)
4	S6145P	Valve – Straight Male – 1/4" MNPT x 1/4" tube
5	S6570	Wall Mount "C-Clip" – 2"
6	S6122	Male Adapter Straight – 1/4" NPT x 1/4" tube
7	S6120P	Male Adapter Elbow – 1/4" NPT x 1/4" tube
8	S9172	Filter Head Self-Tapping Screws
9	S6171P	Male Adapter Elbow – 3/8" NPT x 1/4" tube
10	S6123P	Male Adapter Straight – 3/8" NPT x 1/4" tube
11	S6170P	Male Adapter Straight – 3/8" NPT x 3/8" tube