




TECHWEST INC



Tech West Liquid Ring Pump
Pre-Installation Guide


All installations must conform to local codes.

Revised 1/16

Junction Box Installation

Box Mounted Control Panel Diagram

- Mount a three gang junction box to a solid wall stud. It should be placed so that it's front surface is set back 1/2" to 5/8".
- Run one 18-3 bell wire from each piece of equipment up to the installed junction box.
- Cut and install sheet rock.
- Connectors are provided to complete the wiring of the control panel by connecting wires to corresponding colored wires from the low voltage contactor box(es).
- Position the panel, then insert and tighten screws.
- 1 & 2 switch panels use 2 switch mud ring box.
3 or more switch panels use 3 switch mud ring box.



Switch Assembly


Part#
CPLV-100 (vacuum lense)
CPLA-100 (air lense)
CPLH-100 (water lense)

Part#
CPSV-100 (vacuum switch)
CPSA-100 (air switch)
CPSH-100 (water switch)

Part#
CPB-100



Wiring Schematic (4 Switch Panel)




B = Blue
R = Red
W = White

As shown above, control panel equipped with a water control switch are wired such that water switch must be in the "ON" position in order to start the vacuum pump(s). This feature is provided as a standard protective measure to ensure that the vacuum pump(s) are not operated without water and thereby damaged.

In panels without a water switch, the vacuum switch(es) are wired using the same configuration as shown on the air switch above.

Electrical Connections




Ventilation Requirements

Equipment Room Temperature

The Tech West Inc. equipment must be used in an environment that is temperature controlled. Room temperature must remain between 40 and 105 degrees Fahrenheit. Adequate forced ventilation must be provided across the equipment by placing an adequate exhaust fan opposite an equivalent air intake vent. The fan should be higher than the associated intake vent.

Exhaust Vent Protection


If the exhaust piping is venting to the outside of the building, precautions must be taken to protect the equipment room from weather elements and animal intrusion. This can be accomplished by using one of the three methods shown below.



Make sure to use the correct pipe type for associated system

Plumbing Installation

Sub Floor Installation: The sub-floor plumbing layout shown below is the recommended layout for Tech West Inc. system installations and should be used whenever possible.



Notes:

- 1 See optional drain connections shown below
- 2 8 foot maximum height from main line to pump
- 3 Consult dental unit manufacturer's guidelines for correct reduced size and height of termination of vacuum line inside junction box.
- 4 Limit branches. Orient main line under junction box or cabinet
- 5 When main line is 1 1/2" ID or larger, use 45 degree Y's and elbows only.
- 6 Long radius 90 degree elbows can be used as alternates to 45 degree elbows.
- 7 A total of 8 feet of 3/4 inch hose is supplied with Tech West Inc. units. This hose must be shared between inlet and drain.


All installations must conform to local codes

Plumbing Installation

Overhead Installation: The overhead plumbing layout shown below is the alternate layout for Tech West Inc.'s system installation and should be used only when unable to use the sub-floor plumbing layout.


Main Line

See plumbing requirement for main line diameters needed for specific units.
See chart A2



Connection Details All Installations

Making Turns




Main Line


- Use only 45° elbows or sweeping (sanitary) 90° elbows to make turns in main line.
- If piping is diverted to clear an obstruction, **DO NOT MAKE A TRAP.**

Drain Options

Direct connection to vented drain. No traps before vent.




Indirect connection (Air gap) with a P-trap.



Target Room Layout

Outside Air Pipe

2-inch pipe for air intake.
Must be protected from rain and animals.




Important Information

A. Safty First -

Always disconnect main power supply before installation.


B. Wire or Plug in -

Wire or plug in outlet to building power supply. (see product deimension chart)
NOTE: Buck boost transformers (BT) are only required when voltage drops below 208V or above 240V. On pump utilizing 115v BT is only required when voltage is below 110V and above 120V.



★ Warning:

Vent Outside



Electrical Information

C. All vacuum pumps need to have dedicated circuit breakers. For breaker size and electrical connection type please refer to the product specifications / Dimensions chart on last page.

Water Plumbing Connection

D. 1/4 FNPT Shut-off valve and 5ft. pressure hose (supplied)

E. 9ft of Flexible PVC hose 3/4" supplied.

F. Exhaust needs to be vented outside with 2" schedule 40 PVC pipe.

★ Warning: Avoid any bends resulting in a downward slope. Condensation could cause water to collect in vent pipe. (See diagram above)

G. Open floor sink - use 3/4" flexible PVC hose.

Service Clearances

H. Allow 12" on all sides for all units

Chart A2 Main Vacuum Line
Vacuum Line Pipe Diameter

Operatories	PVC sch 40	Copper
1	3/4"	3/4"
2	1"	1"
3	1 1/4"	1"
4	1 1/4"	1 1/4"
5	1 1/4"	1 1/2"
6	1 1/4"	1 1/2"
7	1 1/2"	1 1/2"
8	1 1/2"	1 1/2"
9	1 1/2"	2"
10	2"	2"
11	2"	2"
12	2"	2"

Ambient Tempertures

H. Must not exceed 105°F Must remain above 41°F

Product Dimensions

WHIRLWIND PUMP	MODEL	USERS	VOLTAGE	BREAKER	H.P.	HEIGHT	DEPTH	WIDTH	WEIGHT
	VPL2SS	2	115/230	20	1	16	13	13	59
	VPL2SSR	2	115/230	20	1	16	13	13	59
	VPL3SS	3	115/230	20	1.5	18	13	13	66
	VPL3SSR	3	115/230	20	1.5	18	13	13	66
	VPL4S2	4	208/230	20	2	18	13	13	66
	VPL4S2R	4	208/230	20	2	18	13	13	66
	VPL4D2	4	208/230	2 X 20	2 X 1	30	20	30	171
	VPL4D2R	4	208/230	2 X 20	2 X 1	30	20	30	171
	VPL6D2	6	208/230	2 X 20	2 X 1.5	30	20	30	180
	VPL6D2R	6	208/230	2 X 20	2 X 1.5	30	20	30	180
	VPL8D2	8	208/230	2 X 20	2 X 2	30	20	30	182
	VPL8D2R	8	208/230	2 X 20	2 X 2	30	20	30	182
	VPL6T2	6	208/230	3 X 20	3 X 1	30	20	41	255
	VPL6T2R	6	208/230	3 X 20	3 X 1	30	20	41	255
	VPL9T2	9	208/230	3 X 20	3 X 1.5	30	20	41	260
	VPL9T2R	9	208/230	3 X 20	3 X 1.5	30	20	41	260
	VPL12T2	12	208/230	3 X 20	3 X 2	30	20	41	264
	VPL12T2R	12	208/230	3 X 20	3 X 2	30	20	41	264

GOLDEN VAC PUMP	MODEL	USERS	VOLTAGE	BREAKER	H.P.	HEIGHT	DEPTH	WIDTH	WEIGHT
	VPLG3SS	3	115/230	20	1	18	13	13	69
	VPLG3SSR	3	115/230	20	1	18	13	13	69
	VPLG5SS	5	115/230	20	2	18	13	13	69
	VPLG5SSR	5	115/230	20	2	18	13	13	69
	VPLG6D2	6	208/230	2 X 20	2 X 1	30	20	28	188
	VPLG6D2R	6	208/230	2 X 20	2 X 1	30	20	28	188
	VPLG10D2	10	208/230	2 X 20	2 X 2	30	20	28	188
	VPLG10D2R	10	208/230	2 X 20	2 X 2	30	20	28	188
	VPLG15T2	15	208/230	3 X 20	3 X 2	30	20	41	273
	VPLG15T2R	15	208/230	3 X 20	3 X 2	30	20	41	273

ECOPUMP	MODEL	USERS	VOLTAGE	BREAKER	H.P.	HEIGHT	DEPTH	WIDTH	WEIGHT
	ESL2S	2-3	208/230	20	1	19	13	13	59
	ESL4D	4-6	208/230	2 X 20	2 X 1	30	20	30	166
	ESL5S	4-5	208/230	20	2	19	13	13	59
	ESL10D	8-10	208/230	2 X 20	2 X 2	30	20	30	166