


SNAP RING PLIERS

009.007.00

¥41.50




- Simplifies installation and removal of both internal and external snap rings
- Fits all snap ring sizes used in A-dec equipment

A-DEC SLEEVE TOOL

98.0072.00

¥7.25




- Use to secure the 1/4" and 1/8" tubing sleeves when installing tubing on barb fittings

HEMOSTAT

009.008.00

¥6.25



- Temporarily stops air or water flow through tubing while troubleshooting or repairing the unit

UMBILICAL TUBING STRINGER

009.015.00

¥90.00




- Use to string new or additional tubing and/or wiring into an existing umbilical assembly
- 12" of stringer with threading holes at each end

O-RING TOOL KIT

009.013.00

¥25.00




- For quick field repair
- Includes three tools, which fit the four smallest O-ring sizes used in A-dec equipment

SYRINGE WRENCH

98.0900.00


¥3.00



TUBING PLIERS

009.014.00

¥41.50




- Simplifies installation of 1/8" OD tubing to barbed fittings

TUBING CUTTERS

009.036.00


¥35.00



TUBING STRIPPERS

009.035.00


¥6.00



FIBER-OPTIC INSTALLATION TOOL KIT

90.0383.00

¥110.00



- For installing dual voltage intraoral light source and for adjusting voltage

PANEL MOUNT GAUGE

026.118.00

¥25.00




- 0 - 100 psi
- 1 1/2" diameter
- Used to check air/water pressure valves at various system locations
- 1/8" MPT x 10 - 32 threads

WATER RETRACTION TESTER

98.0051.00

N/C




- Use to determine whether or not a handpiece control system retracts water
- Includes instructions

HANDPIECE PRESSURE GAUGE

50.0271.00

¥180.00




- 0 - 60 psi
- Fits all handpiece couplers except Borden 3-hole

HEX WRENCH SET

009.018.00

¥10.00



- Includes hex wrenches required in servicing A-dec equipment
- Plastic holder

CABLE TIES

025.036.01 (10) Small ties

¥1.25

025.002.01 (10) Medium ties

¥2.80

025.003.01 (10) Large ties

¥7.00



- Three sizes:
Small: 3 3/4" (9.5 cm)
Medium: 6 3/4" (17 cm)
Large: 14" (35 cm)

BALL END, HEX DRIVER

009.037.00 1/8," Small handle


¥6.00

009.005.00 3/32," Std handle

¥6.00

008.023.00 5/64," Std handle

¥6.00




CLIPS AND CABLE TIE MOUNT

025.017.01 Cord clip (10)

¥12.75

025.019.01 Mount (10)

¥13.50



Adhesive-backed cord clip




Adhesive-backed cable tie mount

SYRINGE TOOL

98.0010.00

☒14.25



- Use to remove/replace valve cores on Uni-Flo and Tri-Flo Syringes
- Use 3/32" to adjust syringe flow on Century II
- Use 5/32" to remove/replace the syringe tip retaining nut

VALVE TEST SYRINGE

98.0050.01 (4)

☒13.75



- For making quick tests of pilot-operated valves
- Applies a static pressure of 5 to 75 psi

A-DEC SILICONE LUBRICANT

98.0090.01 (6)

☒4.00

**TOUCH-UP PAINT**

90.0351.00 Surf 4

☒10.50

90.0351.01 Surf 6

☒10.50

90.0351.02 White 2

☒10.50

90.0351.03 Gray 3

☒10.50

90.0350.09 Surf

☒10.50



- 1/2 fluid ounce
- Includes brush

LOCTITE®

060.001.00 Red ☒271

☒4.15

060.002.00 Blue ☒242

☒4.15

**SERVICE TIP**

A-dec's silicone lubricant is a high quality silicone-based grease. It is ideal for lubricating A-dec O-rings. Before installing O-rings always apply a light coating of silicone lubricant. This makes the installation of the O-ring easier and prevents damage to the O-ring. A-dec silicone lubricant is available in a package of 6 as Part Number 98.0090.01.

CAUTION

The use of petroleum-based lubricants (jellies, creams, or lotions) causes the O-rings in the syringe to swell and/or deteriorate, resulting in syringe malfunction or failure.