TROUBLE-SHOOTING

PROBLEM

POSSIBLE CAUSE

CORRECTIVE MEASURES

I KODLEMI	CAODE	
Compressor motor does not start.	1. No power	1. Check circuit breaker
2. Compressor motor runs, but, is noisy.	Plastic intake hose not installed properly or is cracked.	Replace; make sure plastic hose does not touch chrome tubing.
3. Compressor motor tries to start but Compresso-Dri circuit breaker trips.	3a. Voltage too low.	3a. A minimum of 208/105 volts is needed. If voltage is 208/105V or below, install a boost transformer. Consult your dealer.
	3b. Unloader valve not opening when compression cycle stops.	3b. Check unloader valve.
	3c. Defective Compresso- Dri circuit breaker.	3c. Replace
4. Compressor motor runs, but will not build up pressure to 100 PSI.	4a. Unloader valve not closing when compressor runs.	4a. Refer to service technician.
•	4b. Intake may be restricted.	4b. Remove intake filter. If this corrects problem, replace filter.
	4c. Leak in compressor.	4c. Close the main tank shut-off valve. Check compressor discharge tubing, cylinder heads, bolts, packing gland in tank shut-off valve, pressure relief valve and all fittings for leaks. Repair as necessary. (Verify pump-up time for specific unit, see p. 7)
	4d. Pressure switch needs adjustment.	4d. Turn off power to unit. Drain tank slowly until pressure switch clicks. Tank pressure gauge should read 80 PSI. Turn on power and verify pump-up time for specific unit on data sheet. (See p. 7) Call service technician if adjustment is needed.

PROBLEM

POSSIBLE CAUSE

CORRECTIVE MEASURES

5. Compressor cycles with no air being used.	5a. Total regeneration valve is in "ON" position.	5a. After insuring moisture monitor is blue, turn total regeneration valve "OFF"
	5b. Leak in office air system.	5b. Close main tank shut-off valve. Pump-up main tank to 100 PSI. If pressure is maintained at 100 PSI for 15-20 minutes, leak is in air system, not compressor. Contact plumber.
	5c. Leak in compressor.	5c. See 4d.
6. Moisture monitor is not blue.	6a. Compressor running too frequently or cycles without air being used.	 6a. 1. Compressor undersized for installation—check with dealer. 2. Leak in Compresso-Dri or air system—locate and repair.
	6b. Unloading system not functioning properly.	6b. Check to see that the unloader closes when the compressor runs and opens when the compressor stops. Purge tank pressure should decrease from 100 to 0 PSI within 40-45 seconds of shut-off. If time is not correct, check for clogged or sticking pressure switch unloader valve.