

# CHAPMAN-HUFFMAN, INC.

## ASEPSIS 3-WAY SYRINGE - PIN TYPE



*Syringe Only* 40-300-00

*Extra Asepsis 3-way Syringe* 40-312-00  
*Head Assembly - Pin Type, Less Handle*

**SYRINGE WITH COILED ASEPSIS TUBING**  
*Black* 40-301-00

*Gray* 40-302-00

*Surf* 40-303-00

*Sterling* 40-304-00

*Surf 6* 40-305-00

**SYRINGE W/ STRAIGHT ASEPSIS TUBING**  
*Black* 40-306-00

*Gray* 40-307-00

*Surf* 40-308-00

*Sterling* 40-309-00

*Surf 6* 40-310-00

- **Autoclavable** head assembly (after shutting off air & water supply to syringe and bleeding of pressure) head assembly can be removed for autoclaving at 275° F
- Autoclavable quick disconnect aerator tip can separately be removed for sterilization
- Locked-in cartridge style
- Replaceable air/water cartridges
- 360° tip rotation
- Precise mist, air and water control
- **Matte chrome finish**

# CHAPMAN-HUFFMAN, INC.

## AUTOCLAVABLE SYRINGE WITH QUICK DISCONNECT TIP



MANUFACTURED BY  
CHAPMAN-HUFFMAN, INC.  
IN THE USA

*Syringe Only* 40-700-00

*Extra Autoclavable Syringe  
Head Assembly With Quick  
Disconnect Tip, Less Handle* 40-712-00

### SYRINGE W/ COILED ASEPSIS TUBING

*Black* 40-701-00

*Gray* 40-702-00

*Surf* 40-703-00

*Sterling* 40-704-00

*Surf 6* 40-705-00

### SYRINGE W/ STRAIGHT ASEPSIS TUBING

*Black* 40-706-00

*Gray* 40-707-00

*Surf* 40-708-00

*Sterling* 40-709-00

*Surf 6* 40-710-00

- **Autoclavable** head assembly (after shutting off air & water supply to syringe and bleeding of pressure) head assembly can be removed for autoclaving at 275° F
- Quick disconnect tip from head by pushing in locking ring
- Locked in cartridge style air/water buttons
- 360° tip rotation
- Precise mist, air and water control
- **Matte chrome finish**
- **DCI cross ref. # 3430**

SYRINGES

**800-225-4621**

**184** FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: [ORDERS@CHAPMAN-HUFFMAN.COM](mailto:ORDERS@CHAPMAN-HUFFMAN.COM)

# CHAPMAN-HUFFMAN, INC.

---

## 3-WAY SYRINGE - PIN TYPE



**Syringe Only** 40-200-00

### SYRINGE W/ COILED ASEPSIS TUBING

**Black** 40-201-00

**Gray** 40-202-00

**Surf** 40-203-00

### SYRINGE W/ STRAIGHT ASEPSIS TUBING

**Black** 40-206-00

**Gray** 40-207-00

**Surf** 40-208-00

- Locked-in cartridge style
- Replaceable air/water cartridges
- 360° tip rotation
- Precise mist, air and water control
- Autoclavable quick disconnect aerator tip for easy sterilization

---

## 3-WAY CIRCULATING SYRINGE - PIN TYPE



**Syringe Only** 40-250-00

### SYRINGE W/ COILED ASEPSIS TUBING

**Black** 40-251-00

**Gray** 40-252-00

**Surf** 40-253-00

### SYRINGE W/ STRAIGHT ASEPSIS TUBING

**Black** 40-256-00

**Gray** 40-257-00

**Surf** 40-258-00

- Water circulating type, keeps warm water at syringe head when used with existing water warmer
- Locked-in cartridge style
- Replaceable air/water cartridges
- 360° tip rotation
- Precise mist, air and water control
- Autoclavable quick disconnect aerator tip for easy sterilization

## 3-WAY SYRINGE



*Syringe Only* 40-100-00

### SYRINGE W/ COILED ASEPSIS TUBING

*Black* 40-101-00

*Gray* 40-104-00

*Surf* 40-111-00

### SYRINGE W/ STRAIGHT ASEPSIS TUBING

*Black* 40-105-00

*Gray* 40-108-00

*Surf* 40-112-00

- **Conventional valve core design**
- 360° tip rotation
- Precise mist, air and water control
- Replaceable air/water buttons
- Autoclavable quick disconnect aerator tip for easy sterilization
- **CAUTION: Do not use sterile sleeves or disposable tips with this syringe.**

## SINGLE SHOT - AIR OR WATER SYRINGE



*Syringe Only* 40-400-00

### SYRINGE WITH COILED TUBING

*Black* 40-401-00

*Gray* 40-402-00

*Surf* 40-405-00

### SYRINGE WITH STRAIGHT TUBING

*Black* 40-403-00

*Gray* 40-404-00

*Surf* 40-406-00

- **Conventional valve core design**
- Single function for air or water
- 360° tip rotation
- Replaceable air/water button
- Autoclavable quick disconnect aerator tip for easy sterilization

# CHAPMAN-HUFFMAN, INC.

---

## RESTORATIVE ULTRA DRY AIR SYRINGE



**Syringe Only** 40-450-00

### SYRINGE WITH COILED TUBING

**Black** 40-451-00

**Gray** 40-452-00

**Surf** 40-455-00

### SYRINGE WITH STRAIGHT TUBING

**Black** 40-453-00

**Gray** 40-454-00

**Surf** 40-456-00

- For air only-single function
- **Conventional valve core design**
- For use in restorative procedures where ultra dry air is needed
- Built-in air filter (housed inside handle)
- 360° tip rotation
- Replaceable air button
- Autoclavable quick disconnect aerator tip for easy sterilization

---

## AUTOCLAVABLE CONTINENTAL 3-WAY SYRINGE



**Syringe Only** 40-600-00

### SYRINGE WITH 7' STRAIGHT SHEATHED ASEPSIS TUBING ASSEMBLY

**Gray** 40-607-00

**Surf** 40-608-00

**Sterling** 40-609-00

**Surf 6** 40-610-00

- **Autoclavable head assembly** (after shutting off air & water supply to syringe and bleeding off pressure) head assembly can be removed for autoclaving at 275° F
- Autoclavable quick disconnect aerator tip can separately be removed for sterilization
- 360° tip rotation
- Locked-in cartridge style
- Replaceable air/water cartridges
- Precise mist, air and water control
- **Matte chrome finish**

**800-225-4621**

# CHAPMAN-HUFFMAN, INC.

## QUICK DISCONNECT TIP ADAPTER KIT



**NEW!**

**40-730-00**

- Fits 40-700-00 Autoclavable syringe with quick disconnect tip
- Fits DCI Syringes  
DCI cross ref. # 3088

## SYRINGE BUTTON RETAINER SCREW



**40-616-00**

- Fits autoclavable continental 3-way syringes  
40-600-00, 40-607-00, 40-608-00, 40-609-00 and  
40-610-00
- Pkg. of 5

## NYLON WASHER



**95-139-00**

- Used with 8-32 fittings
- Pkg. of 25

## BARB FITTING



**95-119-00**

- Used on syringe head (in) and (out) ports
- 8-32 X 1/16 barb
- Plated
- Pkg. of 10

## CIRCULATING MANIFOLD ASSEMBLY



**40-551-00**

- Brass Construction
- (3) 1/16 Barbs

## SYRINGE TIP KIT AUTOCLAVABLE CHROME FINISH



**40-035-00**

- Tip adapter and O-ring
- Tip nut and O-ring
- Hex adapter wrench, 5/32
- Tip and cone
- Tip nut wrench
- Fits syringe 40-100-00, 40-200-00, 40-250-00,  
40-400-00 and 40-450-00

### Syringe Tip Kit

**40-040-00**

- **Note: Syringe Tip Kit adapts to Marco, Engle and Zirc style syringes**
- Same features as 40-035-00

# CHAPMAN-HUFFMAN, INC.

---

## SYRINGE TIP KIT AUTOCLAVABLE MATTE CHROME FINISH

40-328-00



- Fits 40-300-00 and 40-600-00 series syringe
- Matte chrome finish on tip and tip nut

---

## SYRINGE VALVE O-RING KIT

95-351-00



**Valve Core Wrench Only**

019-20-00

- This kit used for rebuilding 40-100-00 syringe buttons and valve cores

---

## SYRINGE TIP AUTOCLAVABLE CHROME FINISH

40-038-00



- Syringe tip only
- Interchangeable with MARUS, ENGLE and older ADEC syringes
- Fits syringe 40-100-00, 40-200-00, 40-250-00, 40-400-00 and 40-450-00

---

## SYRINGE TIP AUTOCLAVABLE MATTE CHROME FINISH

40-326-00



- Fits 40-300-00 and 40-600-00 series syringes
- **Matte chrome finish**

---

## AUTOCLAVABLE QUICK DISCONNECT SYRINGE TIP

40-740-00



**NEW!**

- Syringe tip only
- Fits 40-700-00 series autoclavable syringe with quick disconnect tip
- **Matte chrome finish**
- DCI cross ref. # 3060

**800-225-4621**

FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: [ORDERS@CHAPMAN-HUFFMAN.COM](mailto:ORDERS@CHAPMAN-HUFFMAN.COM)

SYRINGES

# CHAPMAN-HUFFMAN, INC.

## 3-WAY SYRINGE CARTRIDGE

40-213-00



- Fits (40-200-00) 3-way syringe - pin type and (40-250-00) 3-way circulating syringe - pin type
- Fits ENGLE DENTAL Syringe Button cross ref. # 20-380
- DCI cross ref. # 3085
- Requires 95-138-00 Conical Compression Spring (pkg. of 10) sold below
- Chrome plated button

## CONICAL COMPRESSION SPRING

95-138-00



- Fits under cartridge assembly in 40-200-00, 40-250-00, 40-300-00, 40-600-00 and 40-700-00 Series 3-way syringes
- Returns cartridge to normal position after having been depressed
- **Pkg. of 10**

## SYRINGE BUTTON FOR AIR-WATER

40-045-00



- Syringe button assembly fits syringe 40-100-00, 40-400-00 and 40-450-00

## SYRINGE VALVE CORE

40-032-00



- Fits syringe 40-100-00, 40-400-00, 40-450-00
- Interchangeable with old style ADEC short valve cores cross ref. # 026.033.00
- **Pkg. of 10**

## AUTOCLAVABLE 3-WAY SYRINGE CARTRIDGE - PIN TYPE

40-220-00



- ADEC style replacement soft touch syringe button assembly
- Chrome plated button
- Autoclavable 275°
- Requires 95-138-00 Conical Compression Spring (pkg. of 10) sold above
- ADEC Cross ref. # 23.1098.00 (Prior to Oct. 1995)
- MARUS cross ref. # 21R428
- DCI cross ref. # 3086

SYRINGES

**800-225-4621**



# CHAPMAN-HUFFMAN, INC.

---

## **AUTOCLAVABLE 3-WAY SYRINGE CARTRIDGE - PIN TYPE**

**40-230-00**



- Fits Chapman-Huffman 40-300-00, 40-600-00 and 40-700-00 series syringes and also fits ADEC syringe
- Autoclavable 275°
- Requires 95-138-00 Conical Compression Spring (pkg. of 10) sold on pg. 190
- **Matte chrome finish on button**
- ADEC **Old** cross ref. # 23.1098. 00 (**Beginning October 1995**)
- ADEC **New** cross ref. # 23.1232.01

Note: The new cross ref. # 23.1232.01 ADEC sells pkg. of 2. However Chapman-Huffman sells one at a time.

---

## **REPLACEMENT AIR-WATER SYRINGE CARTRIDGE**

**40-039-00**



- Syringe cartridge assembly for air or water
- Fits older style MARUS, MARCO, ENGLE and ZIRC syringes

---

## **FOREST STYLE SYRINGE CARTRIDGE**



**Non-Autoclaveable Black O-ring**  
**Autoclavable Brown O-rings**

**40-015-00**

**40-020-00**

- Fits FORSTAR-3 syringe
- Cross ref.  
# 1107-070 **non-autoclavable**  
# 1107-330 **autoclavable**
- Requires 40-010-00 service tool below for installation

---

## **FOREST STYLE SYRINGE SERVICE TOOL**

**40-010-00**



- Fits FORSTAR-3 syringe for removal of cartridge, tip adapter and barbs
- Cross ref. # 1107-057
- Stainless steel construction

**800-225-4621**

# CHAPMAN-HUFFMAN, INC.

## SYRINGE BUTTON REMOVAL TOOL

019-50-00



- Assists with the removal of syringe buttons that are stuck. The tapered design forces the button up and out as you push the tool under the button head
- Plastic construction

## MINIATURE FLOW RESTRICTOR

14-245-00



- For air or water flow control, ideal for handpiece drive air adjustment, or syringe flow adjustment
- Fits 1/4" or 1/8" OD tubing to restrict flow
- Brass construction, nickel plated

## SYRINGE FLOW RESTRICTOR

14-250-00



- Adapts to standard and asepsis syringe 2 hole tubing for flow adjustment
- Brass construction, nickel plated

## ADJUSTABLE SYRINGE CONNECTOR BLOCK


40-058-00



- Nickel plated
- Brass construction
- Metering screws allow independent adjustment of air and water
- In ports - Barbs fit 1/8" ID X 1/4" OD tubing
- Out ports - Barbs fit 1/16" ID X 1/8" OD tubing
- Mounting screws

## ADJUSTABLE SYRINGE CONNECTOR BLOCK

40-055-00



- Plastic construction
- Metering screws allow independent adjustment of air and water
- In ports - Barbs fit 1/16" ID X 1/8" OD tubing
- Out ports - Barbs fit 1/16" ID X 1/8" OD tubing
- Mounting screws

SYRINGES

**800-225-4621**

192 FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: [ORDERS@CHAPMAN-HUFFMAN.COM](mailto:ORDERS@CHAPMAN-HUFFMAN.COM)