SPARE-PARTS MANUAL

Publication part number 10006517 Version 8 Published 2009-09



PLANMECA INTRA

GENERATOR BOX	3
Generator assembly,	
220-240V~ code 06601005	
110-115V~ code 06601004	
100V~ code 06601003	4
CONTROL PANEL	5
EXTENSION ARM	6
EXTENSION ARM	
with new arm cover plugs, from ser. no. IXRF67494	7
Bracket arm assembly	8
Bracket arm assembly	
with new arm cover plugs, from ser. no. IXRF67494	,
code 06601100	9
Tube head support frame assembly	10
TUBE HEAD	П

OTHER INSTALLATION ALTERNATIVES

CEILING MOUNTING

Column	12
Adjustable ceiling adapter	13

ACCESSORIES	
for Delight operating light	4
ACCESSORIES	
for Delight S operating light	5
ACCESSORIES	
for LED operating light	6
SINGLE STUD ASSEMBLY	7
FIXED CONTROL PANEL	8
HOOK FOR REMOTE CONTROL/	
DATA WALL SOCKET	9
WALL EXPOSURE SWITCH2	0
MOUNTING ACCESSORIES 2	ı
PM INTRA MOBILE BASE2	2
MOBILE INTRA ACCESSORIES 2	3

PROSTYLE COMPACT DENTAL UNIT MOUNTING

New arm,	
from Prostyle Compact serial no. 338272	24
Prostyle Compact dental unit mounting,	
up to serial no. 338271	25
Prostyle Compact dental unit mounting,	
up to serial no. 332114	26
Light post, new arm,	
from Prostyle Compact serial no. 338272	27

Light post,

ıр	to	serial	no.	33827	I	28	
----	----	--------	-----	-------	---	----	--


PLANMECA DIXI

DIXI3 DIGITAL IMAGING SYSTEM	29
PM ETHERNET INTERFACE	30
PLANMECA DIXI SENSOR HOLDER KIT	3 I
SPARE-PART SETS FOR PLANMECA	
SENSOR HOLDERS	32

Planmeca pursues a policy of continual product development. Although every effort is made to produce up-to-date product documentation this publication should not be regarded as an infallible guide to current specifications. We reserve the right to make changes without prior notice.


COPYRIGHT PLANMECA HELSINKI FINLAND 2009 PUBLICATION PART NUMBER 10006517 VERSION 8 RELEASED: 17 SEPTEMBER 2009 Planmeca Intra GENERATOR BOX

Generator Box




Planmeca Intra GENERATOR BOX


Generator assembly




Control Panel




Extension Arm





Planmeca Intra Extension ARM

Extension Arm with new arm cover plugs, from serial no. IXRF67494



10015853 Extension arm assembly 385

10015855 Extension arm assembly 535


10015852 Extension arm assembly 760

10015847 Extension arm assembly 835

10015851 Extension arm assembly 935


10015854 Extension arm assembly 990

10015848 Extension arm assembly 1160




EXTENSION ARM

Bracket Arm Assembly




Planmeca Intra EXTENSION ARM

Bracket Arm Assembly with new arm cover plugs, from serial no. IXRF67494




10014716 Arm cover plug


Planmeca Intra Extension Ari




Planmeca Intra TUBE HEAD

Tube Head







Adjustable Ceiling Adapter




Planmeca Intra CEILING MOUNTING




Planmeca Intra CEILING MOUNTING




Planmeca Intra CEILING MOUNTING




Single Stud Assembly





Fixed Control Panel





Wall exposure switch





MOUNTING ACCESSORIES

Mounting Accessories

10006362 Mounting accessories set









ø8.4/17 DIN 125





00004806 Allen key 2.5mm



00327018 Expansion

00001721 Cable tie


10009515 Intra single stud accessory set




M8x25 DIN 912



M4x30 DIN 7981



10001262 Mounting accessory set, intra arm



10000250 Arc limiter




10000368 Arc limiter



(b)

10001342 Plug





10002238 Fastener tool

06320011 Extension cable 12m


06311308 Extension

cable PCB 113-10-08




02200030 Tape fastener


PM Intra mobile base



Mobile Intra accessories




New Arm, from Prostyle Compact serial no. 338272




Prostyle Compact Dental Unit Mounting,


up to Prostyle Compact serial no. 338271




Prostyle Compact Dental Unit Mounting, up to Prostyle Compact unit serial no. 332114




Light Post, New Arm, from Prostyle Compact serial no. 338272



Light Post



Dixi3 digital imaging system




PM Ethernet interface



Planmeca Dixi sensor holder kit

10009880 Planmeca Dixi sensor holder kit V4, size 1 10009881 Planmeca Dixi sensor holder kit V4, size 2



Spare-part sets for Planmeca sensor holders



10007914 Vertical sensor holder set V4, size 1 (10 pcs) 10007915 Vertical sensor holder set V4, size 2 (10 pcs)




10007906 Periapical bite piece set (10 pcs)




10007910 Horizontal sensor holder set, size 1 (10 pcs) 10007912 Horizontal sensor holder set, size 2 (10 pcs)



10007908 Endodontic bite piece set (5+5 pcs)



10007913 Aiming set (1 aiming ring, 2 aiming handles)



10007907 Bitewing bite piece set (10 pcs)

10007916 Sensor sheats (smooth, covers both sensor and holder, 300 pcs)

Spare parts for Universal sensor holder:

10006630 Bite tab holder-red (50 pcs/pkg)

10006631 Bitewing holder-blue (50pcs/pkg)

10006632 Endodontic holder-green (50pcs/pkg)

10006633 Anterior periapical holder-white (50pcs/pkg)

10006634 Posterior periapical holder-yellow (50pcs/pkg)

Alternative sensor sheats:

10006635 Sensor sheaths size 2, tight fit (300pcs/pkg)

10006636 Sensor sheaths size 1 and 0, tight fit (300pcs/pkg)



www.planmeca.com