

CHAPMAN-HUFFMAN, INC.

EXHAUST SEPARATOR ASSEMBLY

30-270-00



- For use with single vacuum pumps up to 1.5 HP
- Separates air/water discharged from a dental water injected vacuum pump
- Separates evacuated gases from liquids by means of gravity. Liquids exit separator in a controlled stream, eliminating hazardous waste splashing at drain location
- Enables contaminated air to be vented from vacuum pump location
- Elbow accepts 3/4 PVC pipe or 3/4 PVC **Flex Pipe** from Chapman-Huffman; part number 10-350-00
- Includes mounting hardware

AIR GAP / P-TRAP ASSEMBLY

30-260-00



- For use with water injected vacuum pumps, to be used in conjunction with an exhaust separator
- 1 1/2" ABS slip connection to P-trap
- Air gap bridge assembly with 90 degree 3/4" FPT PVC elbow with 360 degree rotation of complete elbow assembly
- Includes (1) 3/4" MPT X 5/8" barb adapter

VACUUM

UNIVERSAL CENTRAL VACUUM AUTOCLAVABLE LEVER VALVE KIT

30-006-00



- Two autoclavable central vacuum lever valves with 5' ASEPSI-FLEX tubing
- One autoclavable saliva ejector tip valve w/ 5' tubing
- Gray central vacuum canister with removable solids collector screen
- Three position holder bar w/ gray universal holders
- Mounting hardware

800-225-4621

FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

159

CHAPMAN-HUFFMAN, INC.

UNIVERSAL CENTRAL VACUUM AUTOCLAVABLE LEVER VALVE KIT

30-007-00



- Same as 30-006-00, except with (1) HVE, (1) SE, (2) Holders and bar

STANDARD CENTRAL VACUUM AUTOCLAVABLE LEVER VALVE KIT

30-001-00



- Two autoclavable central vacuum lever valves with 5' **Asepsi-Flex** tubing
- One autoclavable saliva ejector tip valve with 5' tubing
- Gray central vacuum canister with removable solids collector screen
- Three position holder bar with (1) standard holder and (2) HVE holders
- Mounting hardware

VACUUM

UNIVERSAL CENTRAL VACUUM LEVER VALVE KIT

Gray

30-008-00



- Central vacuum canister with removable solids collector screen
- Two universal central vacuum lever valves with 5' tubing
- One saliva ejector tip valve with 5' tubing
- Three position holder bar with universal holders
- Mounting hardware


800-225-4621

160 FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

CHAPMAN-HUFFMAN, INC.

MANUAL ON-OFF WATER POWERED SALIVA EJECTOR

30-330-00



- Air operated water On-Off valve with Venturi
- S.E. tip with tubing
- Check valve
- Toggle valve
- Drain Y-fitting
- Supply tubing
- Installation instructions

WATER VENTURI ASSEMBLY

30-334-00



- Water powered Venturi used in 30-330-00
- 1/8 FPT both ends
- Saliva ejector vacuum tube (on side of assembly) fits 3/16 ID or 7/32 ID saliva ejector tubing
- Brass construction, Nickel plated

Y-FITTING

30-340-00



- Y-fitting for drain lines
- Fits 5/8 ID vacuum tubing
- Smaller tube fits 3/8 OD Polyurethane tubing such as 10-300-00
- Brass construction, Nickel plated

AIR/WATER SALIVA EJECTOR

30-360-00



- Air/Water operated Saliva Ejector
- Vacuum Venturi Assembly
- Installation kit with instructions
- Recommended input pressure: Air - 40 to 70 psi
Water - 50 psi minimum

VACUUM

AIR/WATER SALIVA EJECTOR KIT

30-365-00



- Air/Water Operated Saliva Ejector
- Vacuum Venturi
- Saliva Ejector with 5' of tubing with check valve
- Holder for Saliva Ejector tip
- 5' of drain tubing
- On/Off toggle valve
- Installation kit with instructions
- Recommended input pressure: Air - 40 to 70 psi
Water - 50 psi minimum

800-225-4621

FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

161

CHAPMAN-HUFFMAN, INC.

MALE GARDEN HOSE TO MALE PIPE



3/4 x 1/2 Brass	14-011-00
3/4 x 3/4 Brass	14-013-00
3/4 x 3/4 White Nylon	14-014-00

- Garden hose thread X MPT adapter

VACUUM ADAPTER



14-166-00

- Commonly used on vacuum pumps
- 3/4 MPT x 1 dia. straight vacuum adapter with installation notch
- Brass construction
- Installation tool sold below

VACUUM ADAPTER INSTALLATION TOOL



019-03-00

- Used for the installation of 14-166-00 vacuum adapter
- Zinc plated steel

RUBBER GARDEN HOSE WASHER



95-126-00

- Fits 3/4" garden hose threaded nut
- Rubber
- Pkg. of 10

HOSE COUPLING ASSEMBLY 90°



95-001-00

- Commonly used in twin vacuum pumps
- Brass construction
- 3/4" barb X 3/4" female thread garden hose nut
- 90° bend
- 95-126-00 rubber garden hose washer not incl.

STRAIGHT HOSE COUPLING ASSEMBLY



017-03-00

- Commonly used in twin vacuum pumps
- Brass construction
- 3/4" barb X 3/4" female thread garden hose nut
- 95-126-00 rubber garden hose washer not incl.

800-225-4621

162 FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

CHAPMAN-HUFFMAN, INC.

SWIVEL GARDEN HOSE CONNECTOR TO FPT



3/4 Female hose to 1/2 FPT 017-06-00
3/4 Female hose to 3/4 FPT 017-07-00

- Swivel connector female hose to female pipe thread
- Brass construction
- 95-126-00 rubber garden hose washer not incl.

HOSE CLAMP ASSEMBLY



20-025-00

- Adjustable from 1/2" to 1 1/4"
- 1/2" wide
- Stainless steel construction

CUSPIDOR DRAIN VALVE



30-600-00

- Quiet operation
- Automatic drain for dental cuspidor
- Connects to central vacuum system
- Eliminates conventional plumbing drain
- Installation kit with instructions
- TRANNEL original

VACUUM

SINK DRAIN VALVE ASSEMBLY



30-605-00

- Quiet operation
- Automatic drain for standard or operatory sinks
- Connects to central vacuum system
- Eliminates conventional plumbing drain
- Installation kit with instructions
- TRANNEL original

DIAPHRAGM



30-614-00

- Fits **molded** vacuum cuspidor and sink valves 30-600-00 and 30-605-00
- Incl. on above assemblies beginning Sept. 1997

800-225-4621

FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

163

CHAPMAN-HUFFMAN, INC.

LONG SCREEN



30-612-00

- Fits 30-600-00 cuspidor drain valve assembly and 30-605-00 sink drain valve assembly
- Included on above assemblies beginning September 1997

GASKET



30-609-00

- Fits Clear Bowls 30-611-00 and 30-613-00
- Fits vacuum cuspidor and sink valves 30-600-00 and 30-605-00
- Also fits DENTAL-EZ SIGNATURE Unit

LONG CLEAR BOWL



30-613-00

- Fits 30-600-00 cuspidor drain valve assembly and 30-605-00 sink drain valve assembly
- Included on above assemblies beginning September 1997

VACUUM

DIAPHRAGM



30-608-00

- Fits vacuum cuspidor and sink valves 30-600-00 and 30-605-00 **manufactured prior to Sept. 1997**

SCREEN



30-610-00

- Fits vacuum cuspidor and sink valves 30-600-00 and 30-605-00 **manufactured prior to Sept. 1997**
- Also fits DENTAL-EZ SIGNATURE Unit

CLEAR BOWL



30-611-00

- Fits vacuum cuspidor and sink valves 30-600-00 and 30-605-00 **manufactured prior to Sept. 1997**
- Also fits DENTAL-EZ SIGNATURE Unit

800-225-4621

164 FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

CHAPMAN-HUFFMAN, INC.

STRAINER FILTER ASSEMBLY



30-125-00

- Commonly used on DENTAL-EZ Vacuum pumps
- Strains incoming water supply
- Bronze construction
- (2) 1/8 FPT ports
- DENTAL-EZ cross ref. # 6-4545-002

ADJUSTABLE VACUUM RELIEF VALVE ASSEMBLY



30-120-00

- Fits VACUDENT ESP Series vacuum pumps - cross ref. # 0951032
- Also fits other vacuum systems where an adjustable vacuum relief valve is needed
- Adjustable from 5" to 27" Hg.
- Steel and plastic construction
- 3/8 MPT

SWING CHECK VALVE ASSEMBLY



1/2 X 1/2 FPT

30-119-00

3/4 X 3/4 FPT

30-121-00

1 X 1 FPT

30-122-00

- Used in vacuum pump systems
- Flap with rubber seat
- Brass construction
- DENTAL-EZ cross ref. # 64622008
(For 3/4 valve only)

VACUUM

ANTI-SIPHON VACUUM BREAKER



1/4 FPT Ports

30-205-00

3/8 FPT Ports

30-200-00

1/2 FPT Ports

30-206-00

- Brass and plastic construction
- APOLLO cross ref. # PV V-50595
(for 30-205-00 only)
- DENTSPLY cross ref. # 1542601
(for 30-205-00 only)
- Fits DOUBLE, JR., 80 and 84 series vacuum pumps

800-225-4621

FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

165

CHAPMAN-HUFFMAN, INC.

ANTI-SIPHON VACUUM BREAKER REPAIR KIT



30-210-00

- All parts needed to rebuild 30-200-00 or 30-205-00 vacuum breakers

ANTI-SIPHON VACUUM BREAKER REPAIR KIT



30-211-00

- All parts need to rebuild 30-206-00 **1/2" FPT** Vacuum Breaker

ANTI-SIPHON VACUUM BREAKER



30-201-00

- Fits MATRIX vacuum system cross ref. # 7700-5139 (vacuum breaker only)
- Brass construction
- **Note: You must use flow restrictor and other fittings from existing vacuum breaker**
- 1/4 FPT ports
- Vacuum breaker can be installed horizontally or vertically

WATER FLOW CONTROL VALVE



0.13 GPM Flow

19-302-00

- DOLE type directional flow control
- 0.13 GPM flow
- This valve restricts the flow of water to a water injected vacuum pump impeller housing
- Commonly used on DENTAL-EZ 1 HP vacuum pumps when used with a water recirculator
- 3/8 FPT both ends
- DENTAL-EZ cross ref. # 6-4622-018

WATER FLOW CONTROL VALVE



0.19 GPM Flow

19-306-00

- DOLE type directional flow control
- 0.19 GPM flow
- This valve restricts the flow of water to a water injected vacuum pump impeller housing
- 3/8 FPT both ends
- APOLLO DENTAL cross ref. # PV V50505

800-225-4621

166 FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

CHAPMAN-HUFFMAN, INC.

WATER FLOW CONTROL VALVE



0.25 GPM Flow

19-308-00

- DOLE type directional flow control
- 0.25 GPM flow
- This valve restricts the flow of water to a water injected vacuum pump impeller housing
- Commonly used on DENTAL-EZ vacuum pumps when used with water recirculator, DENTAL-EZ cross ref. # 6-4622-019
- 3/8 FPT both ends

WATER FLOW CONTROL VALVE



1/2 GPM Flow

19-312-00

1 GPM Flow

19-313-00

- DOLE type directional flow control
- This valve restricts the flow of water to a water injected vacuum pump impeller housing
- Commonly used on DENTAL-EZ vacuum pumps and others
- 3/8 FPT both ends
- DENTAL-EZ cross ref. # 's
1/2 GPM - 6-4622-009
1 GPM - 6-4622-010

WATER FLOW CONTROL VALVE



3/8 MPT X 3/8 FPT 1/2 GPM Flow

19-309-00

3/8 MPT X 3/8 FPT 1 GPM Flow

19-314-00

- DOLE type directional flow control
- This valve restricts the flow of water to a water injected vacuum pump impeller housing
- Commonly used on DENTAL-EZ vacuum pumps
- **Cross ref. # same as earlier models except with 3/8 MPT on one end and 3/8 FPT on the other end 1/2 GPM # 6-4622-009 1 GPM # 6-4622-010**

VACSTAR TYPE WATER SYSTEM CONVERSION KIT

95-230-00



- Fits AIR-TECHNIQUES VAC STAR vacuum system
- Recommended to replace the N9 straight type vacuum breaker originally installed on system.
- Installation instructions
- Cross ref. # 55510

800-225-4621

FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

167

VACUUM

CHAPMAN-HUFFMAN, INC.

VACUUM CANISTER SOLIDS SEPARATOR ASSEMBLY



SIDE DRAIN PORT

<i>Black</i>	30-010-00
<i>Gray</i>	30-725-00
<i>Lid O-Ring Only</i>	007-35-00

- Used with central vacuum
- Removable screen
- Accepts (2) 1/2" ID vacuum valve tubing
- Accepts (1) saliva ejector tubing
- One 5/8" drain port

VACUUM CANISTER SOLIDS SEPARATOR ASSEMBLY



BOTTOM DRAIN PORT

<i>Black</i>	30-015-00
<i>Gray</i>	30-729-00
<i>Lid O-Ring Only</i>	007-35-00
<i>Mtg. Bracket - Black</i>	30-011-00
<i>Mtg. Bracket - Gray</i>	30-009-00

- Used with central vacuum
- Removable screen
- Accepts (2) 1/2" ID vacuum valve tubing
- Accepts (1) saliva ejector tubing
- One 5/8" drain port
- Mounting bracket

VACUUM

TWIST LOCK LID WITH O-RING



<i>Black</i>	95-153-00
<i>Gray</i>	95-154-00

- Fits Vacuum Canister Solids Separator Assemblies 30-010-00/30-015-00 black and 30-725-00/30-729-00 gray
- Includes lid and O-ring installed

2" POST MOUNT KIT



95-170-00

- Mounts 30-015-00 and 30-729-00 Vacuum Canister Solids Separator assembly to 2" diameter post

800-225-4621

168 FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

CHAPMAN-HUFFMAN, INC.

DISPOSABLE VACUUM CANISTER SCREEN



95-150-00

- Replacement screens for 30-010-00, 30-015-00, 30-725-00 and 30-729-00 plastic canisters
- 2 1/8" diameter screen
- ADEC cross ref. # 17.0593.00
- Pkg. of 100

TWIST LOCK LID CANISTER ASSEMBLY



VERTICAL MOUNTING BRACKET



HORIZONTAL MOUNTING BRACKET



Black 30-012-00

Gray 30-013-00

Vertical Mounting Bracket - Gray 30-035-00

Vertical Mounting Bracket - Black 30-036-00

Horizontal Mounting Bracket - Gray 30-038-00

Horizontal Mounting Bracket - Black 30-039-00

2" Post Mount for kit above 50-030-00

Lid O-Ring Only 007-35-00

- Twist lock removable lid
- Removable clean-out basket
- Will accept two HVE tubings and one SE tubing
- One 1/2" drain port
- Adapter for attaching 5/8" drain tubing
- Plastic construction
- Mounting bracket and screws

TWIST LOCK LID WITH O-RING



Black 95-151-00

Gray 95-152-00

- Fits Twist Lock Lid Vacuum Canister Assembly 30-012-00 black and 30-013-00 gray
- Includes lid and O-Ring installed

DISPOSABLE VACUUM CANISTER SCREEN

95-160-00



- Replacement screen for 30-012-00 and 30-013-00 plastic canisters
- 1 7/8" diameter screen
- Pkg. of 100

800-225-4621


FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

VACUUM

169

CHAPMAN-HUFFMAN, INC.


SOLID COLLECTOR COVER O-RING



007-40-00

- Fits only DENTECH solids collector cover # 01-04-176
- Cross ref. # for O-ring only DENTECH # 01-06-033

PELTON & CRANE DEBRIS CUP COVER O-RING



007-41-00

- Fits PELTON & CRANE SPIRIT debris cup cover
- PELTON & CRANE cross ref. # 42 05 337

DISPOSABLE COVERED VACUUM CANISTER SCREEN



**Each 30-169-00
Pkg. Of 100 95-169-00**

- Replacement screen for 30-012-00 and 30-013-00 plastic canisters when covered screen assembly is desired
- 1 7/8" diameter screen



**Each 30-159-00
Pkg. Of 100 95-159-00**

- Replacement screen for 30-010-00, 30-015-00, 30-725-00 and 30-729-00 plastic canisters when covered screen assembly is desired
- 2 1/8" diameter screen

CUSPIDOR BOWL O-RING KITS



**ADEC Cross ref. # 030.222.00 95-352-00
ADEC Cross ref. # 030.228.00 95-353-00**

- Fits ADEC cuspidor bowls
- Pkg. of 6

RUBBER PLUG



006-20-00

- Plugs vacuum at swivel adapter on autoclavable vacuum handpiece and saliva ejector
- Plugs HVE and saliva ejector ports on vacuum canisters

800-225-4621

170 FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

VACUUM

CHAPMAN-HUFFMAN, INC.

AUTOCLAVABLE DRY ORAL CUP



30-315-00

- Autoclavable
- Fits all vacuum handpiece valves that accept 7/16" OD tip
- 4 1/8" Dia. X 4" High
- White plastic construction
- Uses 30-316-00 liners sold below

LINER FOR DRY ORAL CUP



30-316-00

- Pre-cut
- Used in 30-315-00
- Pkg of 1000

INLINE VACUUM ON-OFF VALVE



1/2 Barbs 30-110-00

5/8 Barbs 30-115-00

- On-off rotary control valve
- Enables easy removal of vacuum cannister lid without shutting down building's vacuum system
- Fits 1/2" ID or 5/8" ID vacuum tubing
- Black plastic

STAINLESS STEEL ASPIRATOR TIP

NEW!

6" 30-460-00

4" 30-440-00



- The **Green Tip** helping save our landfills from plastic disposable aspirator tips
- Stainless steel construction
- Autoclavable
- 7/16" diameter OD
- 6" and 4" lengths available

VACUUM

UNIVERSAL CENTRAL VACUUM ROTARY VALVE



Black 30-025-00

Gray 30-026-00

- On/off and flow control with rotary type valve
- Swivel connector (7/8" OD) allows 360° rotation of entire valve and adapts to any 1/2" ID central vacuum tubing
- Universal valve has adapter added for use with universal type handpiece holders
- Valve accepts 7/16" OD tips

800-225-4621

FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

171

CHAPMAN-HUFFMAN, INC.

STANDARD CENTRAL VACUUM ROTARY VALVE



Black	30-027-00
Gray	30-028-00

- On/off and flow control with rotary type valve
- Swivel connector (7/8" OD) allows 360° rotation of entire valve and adapts to any 1/2" ID central vacuum tubing
- Valve accepts 7/16" OD tips

UNIVERSAL CENTRAL VACUUM LEVER VALVE



Gray	30-060-00
------	-----------

- Finger tip on/off and flow control with lever operated design
- Swivel quick disconnect (7/8" OD) allows 360° rotation of entire valve and adapts to any 1/2" ID central vacuum tubing
- Has adapter added for use with universal type handpiece holders
- Accepts 7/16" OD tips

STANDARD CENTRAL VACUUM LEVER VALVE




Gray	30-065-00
------	-----------

- Finger tip on/off and flow control with lever operated design
- Swivel quick disconnect (7/8" OD) allows 360° rotation of entire valve and adapts to any 1/2" ID central vacuum tubing
- Valve accepts 7/16" OD tips

AUTOCLAVABLE UNIVERSAL VACUUM LEVER VALVE



Valve With Standard Swivel	30-580-00
Valve With Tapered Swivel	30-535-00



- Autoclavable one piece valve body (275° F)
- Finger tip lever control for on/off and flow control
- Swivel quick disconnect allows 360 degrees of tubing rotation
- Quick disconnect anodized aluminum swivel remains in tubing while valve body is autoclaved
- Fits 1/2" ID vacuum tubing (works best with 1/2" asepsis vacuum tubing such as 10-239-00, 10-271-00, 10-371-00 and 10-471-00)
- Valve with Tapered swivel 30-535-00 fits ADEC newer style holder with taper
- Valve accepts 7/16" OD tips

800-225-4621

172 FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

CHAPMAN-HUFFMAN, INC.

UNIVERSAL AUTOCLAVABLE VACUUM LEVER VALVE



30-500-00

Same as 30-500-00, less swivel adapter 30-505-00


- Finger tip lever control for on/off and flow control
- Quick disconnect anodized aluminum swivel adapter remains in tubing while valve body is autoclaved
- 360° swivel at tubing
- Fits 1/2" ID vacuum tubing
- Anodized aluminum body designed for repeated autoclaving
- Easy, no tools needed for assembly
- 275° F
- Valve accepts 7/16 OD tips

STANDARD AUTOCLAVABLE VACUUM LEVER VALVE



Standard
Swivel

Valve With Standard Swivel 30-510-00



Tapered
Swivel

Valve with Tapered Swivel 30-533-00



Less
Swivel

Valve Less Swivel Adapter 30-513-00

- Autoclavable valve (275° F)
- Finger tip lever control for on/off and flow control
- Swivel quick disconnect allows 360 degrees of tubing rotation
- Quick disconnect anodized aluminum swivel remains in tubing while valve body is autoclaved
- Fits 1/2" ID vacuum tubing (works best with 1/2" asepsis vacuum tubing such as 10-239-00, 10-271-00, 10-371-00 and 10-471-00)
- Valve with Tapered swivel 30-535-00 fits ADEC newer style holder with taper
- Valve accepts 7/16" OD tips

RIGHT ANGLE AUTOCLAVABLE VACUUM LEVER VALVE (RA 360)



30-550-00

Same as 30-550-00, less swivel adapter 30-560-00

- Right angle (105°) design
- Autoclavable
- 360° tip rotation
- Valve accepts 7/16 OD tip
- Fits 1/2" ID vacuum tubing
- 360° anodized aluminum swivel at tubing (Quick disconnect swivel adapter remains in tubing while valve body is autoclaved)
- Anodized aluminum body is designed for repeated autoclaving
- Fits universal or vacuum holders
- Easy, no tools needed for assembly
- 275° F

800-225-4621

FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

VACUUM

173

CHAPMAN-HUFFMAN, INC.

LEVER



Saliva Ejector	30-075-00
Central Vac. (HVE)	30-504-00
• Autoclavable (275° F)	
• Gray plastic construction	

LEVER



Saliva Ejector	30-094-00
Central Vac. (HVE)	30-032-00
• Black plastic construction	

PLASTIC SWIVEL



Tapered Saliva Ejector	30-078-00
Saliva Ejector	30-073-00
Central Vac. (HVE 1/2)	30-503-00
Central Vac (HVE 5/8)	30-507-00
• White Plastic construction	
• Saliva Ejector swivels are 5/8" OD and adapt to 3/16" ID or 7/32" ID tubing	
• The Central Vacuum swivels are 7/8" OD and adapt to any 1/2" ID vacuum tubing (or using 30-507-00 swivel adapt to 5/8" vacuum tubing)	
• Fits saliva ejector valves 30-070-00, 30-076-00, 30-080-00 and 30-082-00	
• Fits central vacuum valves 30-500-00, 30-505-00, 30-510-00, 30-513-00, 30-533-00, 30-535-00 and 30-580-00	

VACUUM

METAL SWIVEL



Tapered Saliva Ejector	30-068-00
Saliva Ejector	30-083-00
Central Vacuum (HVE 1/2)	30-523-00
Tapered Central Vacuum (HVE 1/2)	30-530-00
Central Vacuum (HVE 5/8)	30-517-00
• Anodized aluminum construction	
• Saliva ejector swivels are 5/8" OD and adapt to 3/16" ID or 7/32" ID tubing	
• The central vacuum swivels are 7/8" OD and adapt to any 1/2" ID vacuum tubing (or using 30-517-00 swivel adapt to 5/8" vacuum tubing)	
• Fits saliva ejector valves 30-070-00, 30-076-00, 30-080-00 and 30-082-00	
• Fits central vacuum valves 30-500-00, 30-505-00, 30-510-00, 30-513-00, 30-533-00, 30-535-00 and 30-580-00	

800-225-4621

174 FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

CHAPMAN-HUFFMAN, INC.


UNIVERSAL TIP ADAPTER



30-040-00

- Adapts to existing **plastic** Central Vacuum Valves to provide use with Universal type handpiece holder
- Accepts 7/16" OD tips


AUTOCLAVABLE SALIVA EJECTOR LEVER VALVE W/ PUSH-ON TIP



Valve With Standard Swivel	30-070-00
Valve With Tapered Swivel	30-085-00
Valve Less Swivel	30-076-00

- **Push-on type tips**
- Quick disconnect valve from swivel for autoclaving (275°F)
- Variable lever control for on/off and flow control
- 360° valve rotation
- Plastic screen
- Anodized aluminum body is designed for repeated autoclaving
- Valve with tapered swivel 30-085-00 fits ADEC newer style holder with taper
- Connects to 3/16" ID or 7/32" ID tubing

AUTOCLAVABLE SALIVA EJECTOR LEVER VALVE WITH SCREW - ON TIP



Valve With Standard Swivel	30-170-00
Valve With Tapered Swivel	30-175-00
Valve Less Swivel	30-176-00

- Screw-on gray autoclavable saliva ejector tip
- Quick disconnect valve from swivel for autoclaving
- Variable lever control for on/off and flow control
- 360° valve rotation
- Plastic screen
- Anodized aluminum swivel
- Anodized aluminum body is designed for repeated autoclaving (275°F)
- Valve with tapered swivel 30-175-00 fits ADEC newer style holder with taper
- Connects to 3/16" ID or 7/32" ID tubing

VACUUM


800-225-4621

FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

175

CHAPMAN-HUFFMAN, INC.

AUTOCLAVABLE SALIVA EJECTOR LEVER VALVE



<i>Valve With Standard Swivel</i>	30-080-00
<i>Valve With Tapered Swivel</i>	30-095-00
<i>Valve Less Swivel Adapter</i>	30-082-00

- Anodized aluminum saliva ejector valve and swivel
- No rubber tip or screen to replace
- Disposable saliva ejector tip is retained by O-ring friction
- Quick disconnect valve from swivel for autoclaving (275°F)
- Variable lever control for on/off and flow control
- 360° valve rotation
- Anodized aluminum body is designed for repeated autoclaving
- Valve with tapered swivel 30-095-00 fits ADEC newer style holders with taper
- Connects to 3/16" ID or 7/32" ID tubing

UNIVERSAL SALIVA EJECTOR LEVER VALVE



Gray **30-150-00**

- Quick disconnect valve from swivel for cleaning
- Variable lever on/off control
- Plastic screen
- 360° valve rotation
- Connects to 3/16" ID or 7/32" ID tubing

UNIVERSAL SALIVA EJECTOR ROTARY VALVE



Gray **30-021-00**

- Variable on/off control
- Larger knob for easier control when using rubber gloves
- Plastic screen
- Connects to 3/16" ID or 7/32" ID tubing
- ADEC style taper

SALIVA EJECTOR TIP ASSEMBLY



NEW!

30-118-00

- Use with water powered type saliva ejector system
- Black **screw-on** style tip
- Plastic screen

800-225-4621

176 FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

CHAPMAN-HUFFMAN, INC.

SALIVA EJECTOR TIP ASSEMBLY



30-018-00

- Use with water powered type saliva ejector systems
- Black **push-on** style tip
- Plastic screen

SALIVA EJECTOR CHECK VALVE



19-311-00

- Commonly used with water powered saliva ejectors to **control backflow** in tubing
- It is recommended that 7/32 ID (10-235-00) tubing be used with this valve
- Size 2 3/8" long x 5/8" OD
- Plastic construction

SALIVA EJECTOR TIP ADAPTER ASSEMBLY



Black Tip 30-087-00
Gray Tip 30-088-00

- Adapts to HVE vacuum valve
- Anodized aluminum body

AUTOCLAVABLE SALIVA EJECTOR TIPS ONLY



Push-On - Black 30-023-00
Push-On - Gray 30-074-00

Screw-On - Black 30-105-00
Screw-On - Gray 30-107-00

- Replacement Autoclavable saliva ejector tips
- **Push-on** and **screw-on** type

VACUUM

SALIVA EJECTOR SCREEN



95-121-00

- Plastic screen fits push-on and screw-on style saliva ejector assemblies
- **Pkg. of 10**

800-225-4621

FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

177

CHAPMAN-HUFFMAN, INC.

SALIVA EJECTOR SCREEN METAL



30-017-00

- Metal construction
- Plated

VACUUM DRAIN ADAPTER



14-003-00

- Adapts 1/2 ID vacuum tubing to 5/8 OD vacuum canister drain port

SALIVA EJECTOR ADAPTER



30-056-00

- Adapts 3/16" ID or 7/32" ID saliva ejector tubing to central vacuum handpiece using 7/16" OD tip

SALIVA EJECTOR TUBING ADAPTER



14-001-00

- Saliva ejector tubing barb X 1/2 ID
- Adapts 10-234-00 or 10-235-00 saliva ejector tubing to vacuum canisters, 1/2 HVE port
- White plastic construction

SALIVA EJECTOR TUBING ADAPTER



14-002-00

- Saliva ejector tubing barb X 1/4 ID
- Adapts 10-234-00 or 10-235-00 saliva ejector tubing to saliva ejector tip
- White plastic construction

VACUUM RELIEF MUFFLER



1/2 FPT

30-215-00

3/4 FPT

30-220-00

- Use with wet type vacuum pump
- Filters and reduces noise from pump air inlet

800-225-4621

178 FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

CHAPMAN-HUFFMAN, INC.

VACUUM INLINE FILTER ASSEMBLY



30-700-00

- 3/4" FPT ports
- Removable stainless steel screen (30-301-00) and gasket (30-302-00)
- Clear see-through bowl
- Plastic construction

CLEAR BOWL REPLACEMENT



30-701-00

- Fits 30-700-00 Vacuum Inline Filter Assembly

FILTER SCREEN



30-301-00

- Fits 30-700-00 vacuum inline filter assembly
- Fits MDT cross ref. # 3-20-2092-10
- Exact replacement for APOLLO DENTAL vacuum pump inline 3/4 FPT filter assembly
- Fits AIR-TECHNIQUES VACSTAR cross ref. # 55624
- Stainless steel construction
- Fits wet type vacuum pump filters w/ 3/4" FPT ports
- Commonly used on single vacuum pumps
- 1 1/2" Dia. X 3 11/16" Long

VACUUM

FILTER SCREEN



30-303-00

- **Coarse mesh** filter screen has a more coarse mesh than 30-301-00 screen
- Also fits 30-700-00 vacuum in-line filter assembly if a coarse mesh filter screen is desired
- Fits MATRX cross ref. # 77005077
- Stainless steel construction
- Fits wet type vacuum pump filters with 3/4 MPT ports commonly used on single vacuum pumps
- 1 1/2" dia. X 3 11/16" long

800-225-4621

FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

179

CHAPMAN-HUFFMAN, INC.

RUBBER GASKET



30-302-00

- Fits 30-700-00 vacuum inline filter assembly
- Fits wet type vacuum pump filters with 3/4" FPT ports

WALL MOUNTING BRACKET



30-702-00

- Fits vacuum inline filter assembly 30-700-00, enabling it to be wall mounted
- Includes Zip-It™ wall anchors and screws

VACUUM INLINE FILTER ASSEMBLY



30-710-00

- 1" FPT ports
- Removable stainless steel screen (30-401-00) and gasket (30-402-00)
- Clear see-through bowl
- Plastic Construction

VACUUM

CLEAR BOWL REPLACEMENT



30-712-00

- Clear see-through bowl
- Plastic construction
- Fits 30-710-00 vacuum inline filter assembly

FILTER SCREEN



30-401-00

- Fits 30-710-00 vacuum inline filter assembly
- Fits DENTAL-EZ cross ref. # 6-4609-013
- Also fits MDT cross ref. # 3-20-0457-10
- Stainless steel construction
- Fits wet type vacuum pump filters with 1" FPT ports
- Commonly used on double vacuum pumps
- 2 1/4" Dia. X 4 7/32" Long

800-225-4621

180 FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

CHAPMAN-HUFFMAN, INC.

RUBBER GASKET



30-402-00

- Fits 30-710-00 vacuum inline filter assembly
- Fits wet type vacuum pump filters with 1" FPT

WALL MOUNTING BRACKET



30-711-00

- Fits Vacuum Inline Filter Assembly 30-710-00, enabling it to be wall-mounted
- Includes Zip-It™ wall anchors and screws

VACUUM ADAPTER



14-166-00

- Commonly used on vacuum pumps
- 3/4 MPT X 1 dia. straight vacuum adapter
- Brass construction with installation notches
- Vacuum Adapter Installation Tool 019-03-00 found on pg. 109

MALE GARDEN HOSE TO MALE PIPE



3/4 X 1/2 Brass 14-011-00

3/4 X 3/4 Brass 14-013-00

3/4 X 3/4 White Nylon 14-014-00

- Garden hose thread X MPT

VACUUM

FILTER SCREEN



30-425-00

- Fits DENTSPLY MVS series 80 and 84 vacuum pumps with white plastic solids collector
- Stainless steel construction
- 2 3/4" Dia. X 5 1/8" Long


800-225-4621

FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM

181

CHAPMAN-HUFFMAN, INC.

RUBBER GASKET



30-426-00

- Fits DENTSPLY MVS series 80 and 84 vacuum pumps with white plastic solids collector

FILTER SCREEN



30-427-00

- Fits & is exact replacement for APOLLO DENTAL vacuum pump inline filter assembly w/ 1" FPT ports
- Stainless steel construction
- 2 1/4" dia. X 4 7/32" long

FILTER SCREEN BASKET



30-305-00

- Fits DENTSPLY MVS, MVS JUNIOR, SPECTRUM II & III vacuum pump solids collectors with glass jar
- One end closed
- Stainless steel construction
- 1 3/8" ID X 5 7/16" Long

VACUUM

MVS TYPE GASKET



30-411-00

- Fits DENTSPLY MVS and MVS JUNIOR vacuum pumps solids collector with glass jar
- Adhesive with protective paper covering on one side of gasket, which is removed at installation
- Installation instructions

SPECTRUM TYPE GASKET



30-412-00

- Fits DENTSPLY SPECTRUM II and III vacuum pumps solids collector with glass jar
- Adhesive with protective paper covering on one side of gasket, which is removed at installation
- Installation instructions

800-225-4621

182 FAX 800-836-2765 • WWW.CHAPMAN-HUFFMAN.COM • EMAIL: ORDERS@CHAPMAN-HUFFMAN.COM