



---

**Forest Dental Equipment Catalog**

---

DentalEZ® Integrated Solutions



# **Got Color? We Do.**

**9 POWDER COAT PAINT COLORS:** Forest Grey, Shadow, White, Pewter, Midnight, Copper, Charcoal, Bronze, Arctic Silver

**17 ULTRALEATHER® UPHOLSTERY COLORS:** Raven Wing, Cave, Stone, Doe, Papyrus, Charcoal, Diplomat Blue, Dove Grey, Eucalyptus, Riviera, Pearl Midnight, Baltic, Parrot, Red, Apricot, Pearl Pewter, Pearl Abalone

**3 CARBON FIBERS:** Laser Leather, Lemans, Greywood

**2 PLASTICS:** White, Forest Grey

**11 SOLID SURFACES:** Rocky Road Granite, Silver Mist, Platinum Granite, Moonstone, Sanibel Granite, Goldrush Gemstone, Cinder Granite, Starry Starry Night, Arctic Melange, Baja Melange, Antarctica

**4 LAMINATES:** Surf, Folkstone, Oyster Grey, White Blanc

# GET HYGGE!

HYGGE is a Danish word that defines a quality of cosiness and comfortable conviviality. It engenders a feeling of contentment or well-being. Forest Dental recognizes that all treatment rooms need to GET HYGGE—for both doctor and patient! When a patient is reclining in the operatory chair, neck and knee support lend to added comfort and peace of mind. Doctors and assistants have recognized this for years and improvised support by rolling towels to place under the neck and knees. You no longer need to roll the towel! Forest Dental has designed a set of neck and knee pillows that will not only make your patient rest at ease but will be easy on the eyes.






3900 Patient Chair Powder Coat Paint Color: Copper Dual Plush Upholstery Pearl Abalone and Copper Ultraleather® w/ Brown Contrast Deviated Stitching



# THE 7020 PRO SERIES

ANY WAY YOU SPIN IT, IT'S ALL ABOUT  
REAR MOUNT SIDE DELIVERY.




**3900 Patient Chair** Powder Coat Paint Color: Shadow  
**Upholstery** Standard Plush: Riviera Ultraleather®  
**7020 PRO** Powder Coat Paint Color: Bronze - Solid Surface: Baja Melange

# Table of Contents

<b>Pivot Chair Mount Packages</b>	<b>12</b>
Standard IC Control Head	12
Euro Control Head	14
Rear Standard IC Control Head	16
<b>Fixed Chair Mount with Sidebox Packages</b>	<b>18</b>
Standard IC Control Head	18
Euro Control Head	20
<b>Chairs and Stools</b>	<b>25</b>
3900 Chair Accessories	25
3900 Hydraulic Patient Chair	26
3900 Chair Upholstery Options	27
Economy and Deluxe Stools	28
Comfort Stools	29
<b>Lights Monitor Mounts</b>	<b>33</b>
Post Mount Lights	33
Track Mount Lights	34
Ceiling Mount Lights	35
Cabinet and Wall Mount Lights	36
Monitor Mount Assemblies	37
<b>Rear Mount Side Delivery Systems</b>	<b>41</b>
Cabinetry Mounting Configurations	41
7020PRO Duo Swing Delivery	42
7021PRO Duo Swing Delivery	44
7022PRO Assistant's Delivery	46
7023PRO Assistant's Delivery	48
7010 Duo Swing Delivery	50
3486 Assistants Delivery	52
Storage Organizer, Trays and Monitor Mount	54
<b>Ancillary Ordering and Options</b>	<b>57</b>
Standard IC Delivery Unit	57
Euro Delivery Unit	58
<b>Control Head Accessories</b>	<b>59</b>
Ancillary Products	59-61
<b>The Fusion Package</b>	<b>65</b>
The Fusion Pivot Package	66
The Fusion Fixed with Sidebox Package	68
The Fusion Delivery System	70
<b>Simplified</b>	<b>74</b>
Fixed Mount with Sidebox, External Umbilical	74
Fixed Mount without Sidebox, External Umbilical	76
Package FE Standard IC Control Head	78
Fixed Chair Post Mount, External Umbilical	80
Rear Beneath Shelf or Work Surface Mount	81
Side Cabinet and Wall Mount	82
Rear or Side Beneath Counter Mount	83
Rear Beneath Counter Mount	84
Panel Mounts	86
Rear or Side Mount Compact Controls	87
Mobile Cart	88
<b>Factory Options and Parts</b>	<b>90</b>
<b>General Information</b>	<b>96</b>
Dimensions	96-111
Equipment Maintenance	112
Terms and Warranty	113
Professional Sales Associates (PSA)	115
Designer Friendly® Colors and Surfaces	116

## Chair Mount Packages

---



## Five Packages Configured by Forest

To create a Forest package, first start with the delivery system mounting method & control head of your choice. The delivery system is what differentiates one package from another. Purchase the exact package shown or add/delete items to make it your own.



**WORD**  
OF MOUTH

We feel good in a warm and friendly environment. We tell our friends about a restaurant with great food at a reasonable price in the perfect setting. The same is true with dentistry. We've seen that within a year of upgrading, practices experience up to 40% increase in new patients.

(Shown on Front)

**3900 Patient Chair** Powder Coat Paint Color: Bronze

**Upholstery** Duo Color Standard Plush: Doe and Lagoon Ultraleather®


**9078 Light** Powder Coat Paint Color: Bronze

**Standard IC** Powder Coat Paint Color: Bronze - Carbon Fiber: Laser Leather

**1095 Assistant's** Powder Coat Paint Color: Bronze - Carbon Fiber: Laser Leather


## Pivot Chair Mount

Pivot Chair Mounts easily move left to right on large roller bearings to accommodate left and right-handed practitioners. Unlike a Fixed Chair Mount delivery system, the mounting arm comes from underneath the patient rather than over the top, making it less obtrusive and more patient-friendly. Pivot Assistant's Instrumentation, Operatory Lights and Monitor Mounts are popular options to support your style of practice.



## Fixed Chair Mount With Sidebox

Fixed Chair Mounts With Sidebox allow for a Cuspidor, Assistant's Instrumentation, a Monitor and Light to mount on a sidebox, which is large enough to accommodate a host of ancillary support products, finish the aluminum side covers in one of our nine powder coat paint colors and choose a complementary carbon fiber trim to finish the top cap. If you use a cuspidor, select our beautiful high-gloss ceramic black for a customized and beautiful patient-friendly delivery system.



# Pivot Chair Mount Packages

## Standard IC Control Head

### Standard Features:

#### **Delivery System 4485PI**

Automatic control for three handpieces  
asepsis 3-way air/water syringe w/tubing, 5'  
handpiece tubing, 4-hole Midwest, asepsis style, 5'  
pneumatic master control block w/diaphragm  
master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air/water coolant adjustment  
handpiece flush  
handpiece exhaust canister/oil collector  
dual handle brake release function  
self-contained clean water system w/quick disconnect,  
single bottle, 1 1/2 liter  
wet/dry foot control  
Flexible control arm w/pneumatic brake  
vertical travel 28" w/horizontal travel 62"  
Cable for unit-mounted touchpad, 3900 chair  
Internal umbilical, 5' w/3' external from chair to utility center  
Contoured utility center  
air pressure regulator w/filter  
master air shut-off valve  
Chair adapter, pivot



#### **Assistant's Instrumentation, Rear Pivot 1095 (option not shown)**

Central vacuum accessory package  
central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'  
Internal umbilical, 5' w/ 3' external  
Rear chair pivot mount

#### **Cuspidor, Rear Pivot 1580 (option not shown)**

One-piece ceramic bowl w/touch lever water activator  
Automatic cup fill  
Timed bowl rinse  
Swing away or over patient's lap  
Water pressure regulator w/filter  
Master water shut-off valve  
Internal umbilical, 5' w/ 3' external  
Rear chair pivot mount

#### **Cuspidor, Rear Pivot 1590 (option not shown)**

Same as 1580 above, w/additional features:  
Swivel assistant's quad holder bar

Assistant's instrumentation  
central vacuum canister w/ disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'

#### **Light, Halogen 9083**

3-Axis light head movement  
Automatic on/off activation w/override  
Color shade matching temperature setting 4200 to 5500 Kelvin  
Shadow free illumination

#### **Chair, Hydraulic 3900**

Ultraleather® upholstery w/wide back for elbow support  
Double articulating pinch-style headrest  
Thin back w/full lumbar support  
Right/left patient access  
Swing-down locking armrests  
Hydraulic cantilever mechanism  
Emergency shut-off  
30 degree base rotation w/lock  
Multifunction touchpad w/4 programmable preset positions  
Seat height: minimum 17" maximum 29"

### Popular Equipment Choices



Assistant's  
Instrumentation 1095



Cuspidor 1580



Cuspidor 1590



LED Light  
9078



Monitor Mount  
1300RM



Office Comfort  
Stool 6198



Doctor's Comfort  
Stool 6176



Assistant's Comfort  
Stool 6187

	Factory Options	Part No	Retail	Price
<b>Base Package</b> 4485PI Delivery System 3900 Chair 9083 Light 1112-644 Light Post	<b>\$ 17,390</b>			
<b>As Shown Package</b> Additionally equipped with: 1110-468 Swivel Tray, Arm-Mounted 3914-998-CS Chair Control, Touchpad, Unit-Mounted	<b>\$ 17,735</b>			
<b>Designer Friendly® Surface Treatments</b>				
<b>STEP 1 SELECT UNIT PLASTIC FINISHES:</b>				
<input type="radio"/> Forest Grey <input type="radio"/> Greywood <input type="radio"/> Laser Leather				
<input type="radio"/> Lemans <input type="radio"/> White				
<b>STEP 2 SELECT POWDER COAT PAINT FOR UNIT &amp; LIGHT:</b>				
<input type="radio"/> Arctic Silver <input type="radio"/> Bronze <input type="radio"/> Charcoal				
<input type="radio"/> Copper <input type="radio"/> Forest Grey <input type="radio"/> Midnight				
<input type="radio"/> Pewter <input type="radio"/> Shadow <input type="radio"/> White				
<b>STEP 3 SELECT STANDARD POWDER COAT PAINT (FOR CHAIR):</b>				
<input type="radio"/> Forest Grey <input type="radio"/> Shadow				
<b>UPGRADE:</b> Chair Base/Back Powder Coat Paint (+ \$ 1,500):				
<b>UPGRADE:</b> PRO Chair (+ \$ 2,570) (Less Utility Center Housing and Cover)				
<input type="radio"/> Arctic Silver <input type="radio"/> Bronze <input type="radio"/> Charcoal				
<input type="radio"/> Copper <input type="radio"/> Midnight <input type="radio"/> Pewter				
<input type="radio"/> White				
<b>STEP 4 SELECT UTRALEATHER® UPHOLSTERY COLOR CHOICES FOR CHAIR AND DELUXE/ECONOMY STOOLS:</b>				
<input type="radio"/> Apricot <input type="radio"/> Baltic <input type="radio"/> Cave <input type="radio"/> Charcoal				
<input type="radio"/> Diplomat Blue <input type="radio"/> Doe <input type="radio"/> Dove Grey				
<input type="radio"/> Eucalyptus <input type="radio"/> Papyrus <input type="radio"/> Parrot				
<input type="radio"/> Raven Wing <input type="radio"/> Red <input type="radio"/> Riviera <input type="radio"/> Stone				
Pearlized Colors: <input type="radio"/> Abalone <input type="radio"/> Midnight <input type="radio"/> Pewter				
<b>OR SELECT</b> Special Color Ultraleather® (Call for leadtime)				
<input type="radio"/> Chair + \$ 290 3907-999-SP				
<input type="radio"/> Deluxe/Economy Stool + \$ 105 6007-999-SP				
<b>STEP 5 UPHOLSTERY SELECTIONS:</b>				
Plush Selections:				
<input type="radio"/> Standard Back, Wide/Thin N/C 3907-PL-WB				
<input type="radio"/> Narrow Back N/C 3907-PL-NB				
<input type="radio"/> Narrow Back w/Slings + \$ 420 3907-PL-NS				
Plush Additional Options:				
<input type="radio"/> Contrast/Deviated Stitching + \$ 290 3907-999-DS				
<input type="radio"/> Dual Color + \$ 290 3907-999-DC				
<input type="radio"/> Dual Color /Contrast/Deviated Stitching				
Call Customer Service for information (1.800.423.3555)				
<input type="radio"/> Massage Option + \$ 555 3900-100				
Asepsis Selections:				
<input type="radio"/> Standard Back, Wide/Thin N/C 3907-AS-WB				
<input type="radio"/> Narrow Back N/C 3907-AS-NB				
<input type="radio"/> Narrow Back w/Slings + \$ 415 3907-AS-NS				
Asepsis Additional Options:				
<input type="radio"/> Massage Option + \$ 555 3900-100				
Memory Comfort Selections:				
<input type="radio"/> Wide/Thin + \$ 825 3907-MC-WB				
<input type="radio"/> Dual Color + \$ 1115 3907-MC-DC				
<b>DELIVERY SYSTEM WITH STANDARD IC.</b> ..... 4485PI* \$ 6605				
Handpiece illumination system				
One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS	715			
Add'l handpiece illumination tubing (5-hole) ..... 0017-615-CS	275			
One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS	715			
Add'l handpiece illumination tubing (6-pin) ..... 0017-603-CS	195			
<b>Ancillaries (use existing handpiece position) (see page 57)</b>				
Write in selection here .....				
Write in selection here .....				
Additional handpiece position, wet ..... 1100-189-CS	425			
Additional handpiece position, air drive only ..... 1100-190-CS	215			
Assistant's hygiene instrumentation (control head mount) ... 1408-CS	385			
Quick disconnect with flow control 1/2" (for water) .... 0016-454-CS	145			
Quick disconnect with flow control 3/8" (for air) .... 0016-456-CS	145			
Swivel tray, arm-mounted (over control head) ..... 1110-468*	345			
Chrome accents ..... 9000-998-CS	140			
Additional umbilical (maximum 9 feet) (per foot) ..... 1115-999-CS	30			
Routing valve kit with water shutoff & regulator/factory installation ..... 1115-315-CS	305			
Routing valve kit for use with cuspidor/factory installation 1115-316-CS	145			
<b>ASSISTANT'S INSTRUMENTATION, REAR PIVOT</b> ..... ○ 1095	1605			
<b>OR CUSPIDOR, REAR PIVOT</b> ..... ○ 1580	3160			
<b>OR CUSPIDOR WITH ASSISTANT'S INSTRUMENTATION,</b>				
<b>REAR PIVOT</b> ..... ○ 1590	3595			
Quick disconnect with flow control 1/4" (specify air or water) (1095) 0016-455-CS	145			
Additional autoclavable standard HVE				
with tubing, 5' & autoclavable oral dry cup ..... 1400-113-CS	90			
Self-contained clean water system with quick disconnect, single bottle, 1 1/2 liter (1095) ..... 1700-319	275			
Touchpad, unit-mounted (1095) for 3900 chair ..... 1700-088-CS	515			
Tray assembly (1095) (see page 94) ..... 1110-680	575			
Dual articulating quad holder assembly, standard HVE (1095) 0024-227-CS	260			
Deluxe cuspidor unit-mounted activation (1580/1090) .. 1100-480-CS	270			
<b>LIGHT, HALOGEN</b> (mounting options, see pages 33) ..... ○ 9083 *	2195			
<b>OR LED</b> ..... ○ 9078	2995			
• Light post, pivot (does not include chair adaptor)..... 1112-644 *	600			
• Monitor post (for top post mount) (does not include chair adaptor). 1331-644	600			
Monitor mount assembly (monitor screen not included) (see page 37)				
Ring mount (used with light post) ..... 1300RM	1065			
Top post mount (used with monitor post) ..... 1300TPM	1065			
<b>CHAIR, HYDRAULIC</b> ..... ○ 3900 *	7990			
• One chair control included (specify type below) ..... 3914-998-CS *	N/C			
○ Footswitch control				
○ Touchpad, chair-mounted				
○ Touchpad, unit-mounted				
• Additional chair control (limit 2 on chair, limit 2 on unit) 3914-999-CS	515			
○ Footswitch control				
○ Touchpad, chair-mounted				
○ Touchpad, unit-mounted				
• Additional chair options (See page 25 for pricing, write selection here)				
<b>STOOL, DR'S, ECONOMY/DELUXE/COMFORT</b> 6156C/6159C/6176	530/895/895			
<b>STOOL, ASST'S, ECONOMY/DELUXE</b> ..... 6157C/6160C	680/1010			
<b>STOOL, ASST'S/BACK, DELUXE/COMFORT</b> ..... 6161C/6187	1195/1195			
<b>STOOL, OFFICE, COMFORT</b> ..... 6198	1060			
Fully assembled, blanket wrapped (chair, delivery, light).... BLANKET	450			
* ITEM included in AS SHOWN PACKAGE price				
• Choose only one post when sold with unit				
Often overlooked				
		<b>TOTAL</b>		

Equipment dimensions, see page 96

Shipping weight estimate:

Delivery unit & chair 450 lb / 204 kg

Base package 580 lb / 264 kg

As shown package 630 lb / 286 kg

# Pivot Chair Mount Packages

## Euro Control Head

Standard Features:

### **Delivery System 4987PI**

Automatic control for four handpieces  
asepsis 3-way air/water syringe w/tubing, 3'  
handpiece tubing, 4-hole Midwest, asepsis style, 3'  
pneumatic master control block w/diaphragm  
master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air/water coolant adjustment  
handpiece flush  
handpiece exhaust canister/oil collector  
dual handle brake release function  
swivel instrument tray  
self-contained clean water system w/ quick disconnect,  
single bottle, 1 1/2 liter  
wet/dry foot control  
Flexible control arm w/pneumatic brake  
vertical travel 28" w/horizontal travel 62"  
Cable for unit-mounted touchpad, 3900 chair  
Internal umbilical, 5' w/3' external from chair to utility center  
Contoured utility center  
air pressure regulator w/filter  
master air shut-off valve  
Chair adapter, pivot

### **Cuspidor, Rear Pivot 1590**

Same as 1580 below, w/additional features:  
Swivel assistant's quad holder bar  
Assistant's instrumentation  
central vacuum canister w/disposable  
solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'

### **Cuspidor, Rear Pivot 1580 (Option Not Shown)**

One-piece ceramic bowl w/touch lever water activator  
Automatic cup fill  
Timed bowl rinse  
Swing away or over patient's lap  
Water pressure regulator w/filter  
Master water shut-off valve  
Internal umbilical, 5' w/ 3' external  
Rear chair pivot mount

### **Assistant's Instrumentation,**

### **Rear Pivot 1095 (Option Not Shown)**

Central vacuum accessory package  
central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'  
Internal umbilical, 5' w/3' external  
Rear chair pivot mount

### **Light, Halogen 9083**

3-Axis light head movement  
Automatic on/off activation w/override  
Color shade matching temperature setting 4200 to 5500 Kelvin  
Shadow free illumination

### **Chair, Hydraulic 3900**

Ultraleather® upholstery w/wide back for elbow support  
Double articulating pinch-style headrest  
Thin back w/full lumbar support  
Right/left patient access  
Swing-down locking armrests  
Hydraulic cantilever mechanism  
Emergency shut-off  
30 degree base rotation w/lock  
Multifunction touchpad w/4 programmable preset positions  
Seat height: minimum 17" maximum 29"



## Popular Equipment Choices



Assistant's  
Instrumentation 1095



Cuspidor 1580



Add'l HVE Position  
With Oral Dry Cup  
1400-113-CS



LED Light  
9078



Doctor's Comfort  
Stool 6176



Assistant's Comfort  
Stool 6187



Doctor's Deluxe  
Stool 6159C



Assistant's Deluxe  
Stool 6161C

	Factory Options	Part No	Retail	Price
<b>Base Package</b> 4987PI Delivery System 3900 Chair 9083 Light 1112-644 Light Post	<b>\$ 18,460</b>	<b>DELIVERY SYSTEM WITH EURO.....</b>	<b>4987PI*</b>	<b>\$ 7675</b>
	Handpiece illumination system			
	One handpiece (5-hole), accommodates 1 add'l tubing	1700-021-CS	715	
	Add'l handpiece illumination tubing (5-hole).....	0017-615-CS	275	
	One handpiece (6-pin), accommodates 1 add'l tubing .	1700-019-CS	715	
	Add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	195	
<b>As Shown Package</b> Additionally equipped with: 1590 Deluxe Cuspidor 1300RM Monitor Assembly 3914-998-CS Chair Control, Chair Mounted	<b>\$ 23,120</b>	Ancillaries (use existing handpiece position) (see page 57)		
	Write in selection here .....			
	Write in selection here .....			
	Quick disconnect with flow control 1/4" (for water) ...	0016-454-CS	145	
	Quick disconnect with flow control 3/8" (for air) .....	0016-456-CS	145	
	Chrome accents .....	9000-998-CS	140	
	Additional umbilical (maximum 9 feet) (per foot) ....	1115-999-CS	30	
	Routing valve kit with water shutoff & regulator/factory installation .....	1115-315-CS	305	
	Routing valve kit for use with cuspidor/factory installation	1115-316-CS	145	
	<b>ASSISTANT'S INSTRUMENTATION, REAR PIVOT .....</b>	<b>○ 1095</b>	<b>1605</b>	
	<b>OR CUSPIDOR, REAR PIVOT .....</b>	<b>○ 1580</b>	<b>3160</b>	
	<b>OR CUSPIDOR WITH ASSISTANT'S INSTRUMENTATION,</b>			
	<b>REAR PIVOT .....</b>	<b>○ 1590*</b>	<b>3595</b>	
	Quick disconnect with flow control 1/4"			
	(specify air or water) (1095) .....	0016-455-CS	145	
	Additional autoclavable standard HVE			
	with tubing, 5' & autoclavable oral dry cup.....	1400-113-CS	90	
	Self-contained clean water system with quick disconnect,			
	single bottle, 1 1/2 liter (1095). ....	1700-319	275	
	Touchpad, unit-mounted (1095) for 3900 chair .....	1700-088-CS	515	
	Tray assembly (1095) (see page 94) .....	1110-680	575	
	Dual articulating quad holder assembly,			
	standard HVE (1095) .....	0024-227-CS	260	
	Deluxe cuspidor unit-mounted activation (1580/1590)	1100-480-CS	270	
	<b>LIGHT, HALOGEN</b> (mounting options, see pages 33) .....	<b>○ 9083*</b>	<b>2195</b>	
	<b>OR LED.....</b>	<b>○ 9078</b>	<b>2995</b>	
	● Light post, pivot (does not include chair adaptor) .....	1112-644*	600	
	• Monitor post (for top post mount)			
	(does not include chair adaptor) .....	1331-644	600	
	Monitor mount assembly (monitor screen not included) (see page 37)			
	Ring mount (used with light post).....	1300RM*	1065	
	Top post mount (used with monitor post).....	1300TPM	1065	
	<b>CHAIR, HYDRAULIC.....</b>	<b>○ 3900*</b>	<b>7990</b>	
	● One chair control included (specify type below) .....	3914-998-CS*	N/C	
	○ Footswitch control			
	○ Touchpad, chair-mounted			
	○ Touchpad, unit-mounted			
	● Additional chair control (limit 2 on chair, limit 2 on unit) ..	3914-999-CS	515	
	○ Footswitch control			
	○ Touchpad, chair-mounted			
	○ Touchpad, unit-mounted			
	Additional chair options (See page 25 for pricing, write selection here)			
	<b>STOOL, DR'S, ECONOMY/DELUXE/COMFORT</b>	<b>6156C/6159C/6176</b>	<b>530/895/895</b>	
	<b>STOOL, ASST'S, ECONOMY/DELUXE .....</b>	<b>6157C/6160C</b>	<b>680/1010</b>	
	<b>STOOL, ASST'S WITH BACK, DELUXE/COMFORT .....</b>	<b>.6161C/6187</b>	<b>1195/1195</b>	
	<b>STOOL, OFFICE, COMFORT.....</b>	<b>6198</b>	<b>1060</b>	
	● Fully assembled, blanket wrapped (chair, delivery, light)....	BLANKET	450	
	* ITEM included in AS SHOWN PACKAGE price			
	• Choose only one post when sold with unit			
	● Often overlooked			
		<b>TOTAL</b>		

Equipment dimensions, see page 96

Shipping weight estimate:

Delivery unit & chair 445 lb / 203 kg

Base package 533 lb / 264 kg

As shown package 630 lb / 286 kg

# Pivot Chair Mount Packages

## Rear Standard IC Control Head



### Standard Features:

#### **Delivery System 1500PI**

Automatic control for three handpieces

- central vacuum canister w/disposable solids collector
- asepsis 3-way air/water syringe w/tubing, 5'
- autoclavable standard HVE lever valve w/tubing, 5'
- autoclavable saliva ejector lever valve w/tubing, 5'
- handpiece tubing, 4-hole Midwest, asepsis style, 5'
- pneumatic master control block w/diaphragm
- master on/off
- individual handpiece pressure adjustment
- individual water coolant adjustment
- non-retraction water coolant
- air coolant adjustment
- handpiece flush
- handpiece exhaust canister/oil collector
- self-contained clean water system w/quick disconnect, single bottle, 1 1/2 liter
- wet/dry foot control

Cable for unit-mounted touchpad, 3900 chair

Internal umbilical, 5' w/3' external from chair to utility cent

Contoured utility center

- air pressure regulator w/filter
- master air shut-off valve

Rear chair pivot mount

#### **Light, Halogen 9083**

3-Axis light head movement

Automatic on/off activation w/override

Color shade matching temperature setting

- 4200 to 5500 Kelvin

Shadow free illumination

#### **Chair, Hydraulic 3900**

Ultraleather® upholstery w/wide back for elbow support

Double articulating pinch-style headrest

Thin back w/full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation w/lock

Multifunction touchpad w/ 4 programmable preset positions

Seat height: minimum 17" maximum 29"

### Popular Equipment Choices



Swivel Tray,  
Arm Mounted  
1110-731



Monitor Mount  
1300RM



Add'l HVE Position  
With Oral Dry Cup  
1400-113-CS



LED Light  
9078



Office Comfort  
Stool 6198



Doctor's Comfort  
Stool 6176



Doctor's Deluxe  
Stool 6159C



Doctor's Economy  
Stool 6156C

## Factory Options

## Part No

## Retail

## Price

**Base Package** \$ 16,830

1500PI Delivery System  
3900 Chair  
9083 Light  
1112-644-3 Light Post

**As Shown Package** \$ 16,995

Additionally equipped with:  
0024-054 Assistant's Quad Holder Bar  
3914-998-CS Chair Control, Chair Mounted

## Designer Friendly® Surface Treatments

**STEP 1 SELECT UNIT PLASTIC FINISHES:**

- Forest Grey  Greywood  Laser Leather
- Lemans  White

**STEP 2 SELECT POWDER COAT PAINT FOR UNIT & LIGHT:**

- |                                     |                                   |                                |
|-------------------------------------|-----------------------------------|--------------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze      | <input type="radio"/> Charcoal |
| <input type="radio"/> Copper        | <input type="radio"/> Forest Grey | <input type="radio"/> Midnight |
| <input type="radio"/> Pewter        | <input type="radio"/> Shadow      | <input type="radio"/> White    |

**STEP 3 SELECT STANDARD POWDER COAT PAINT (FOR CHAIR):**

- Forest Grey  Shadow

**UPGRADE: Chair Base/Back Powder Coat Paint (+ \$ 1,500):****UPGRADE: PRO Chair (+ \$ 2,570) (Less Utility Center Housing and Cover)**

- |                                     |                                |                                |
|-------------------------------------|--------------------------------|--------------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze   | <input type="radio"/> Charcoal |
| <input type="radio"/> Copper        | <input type="radio"/> Midnight | <input type="radio"/> Pewter   |
| <input type="radio"/> White         |                                |                                |

**STEP 4 SELECT UTRALEATHER® UPHOLSTERY COLOR CHOICES FOR CHAIR AND DELUXE/ECONOMY STOOLS:**

- |                                     |                               |                                 |                                |
|-------------------------------------|-------------------------------|---------------------------------|--------------------------------|
| <input type="radio"/> Apricot       | <input type="radio"/> Baltic  | <input type="radio"/> Cave      | <input type="radio"/> Charcoal |
| <input type="radio"/> Diplomat Blue | <input type="radio"/> Doe     | <input type="radio"/> Dove Grey |                                |
| <input type="radio"/> Eucalyptus    | <input type="radio"/> Papyrus | <input type="radio"/> Parrot    |                                |
| <input type="radio"/> Raven Wing    | <input type="radio"/> Red     | <input type="radio"/> Riviera   | <input type="radio"/> Stone    |

Pearlized Colors:  Abalone  Midnight  Pewter

**OR SELECT Special Color UltraLeather® (Call for leadtime)**

- |   |             |
|---|-------------|
| <input type="radio"/> Chair + \$ 290                | 3907-999-SP |
| <input type="radio"/> Deluxe/Economy Stool + \$ 105 | 6007-999-SP |

**STEP 5 UPHOLSTERY SELECTIONS:**

Plush Selections:

- |   |            |
|---|------------|
| <input type="radio"/> Standard Back, Wide/Thin N/C  | 3907-PL-WB |
| <input type="radio"/> Narrow Back N/C               | 3907-PL-NB |
| <input type="radio"/> Narrow Back w/Slings + \$ 420 | 3907-PL-NS |

Plush Additional Options:

- |   |             |
|---|-------------|
| <input type="radio"/> Contrast/Deviated Stitching + \$ 290    | 3907-999-DS |
| <input type="radio"/> Dual Color + \$ 290                     | 3907-999-DC |
| <input type="radio"/> Dual Color /Contrast/Deviated Stitching |             |

Call Customer Service for information (1.800.423.3555)

- |   |          |
|---|----------|
| <input type="radio"/> Massage Option + \$ 555 | 3900-100 |
|---|----------|

Asepsis Selections:

- |   |            |
|---|------------|
| <input type="radio"/> Standard Back, Wide/Thin N/C  | 3907-AS-WB |
| <input type="radio"/> Narrow Back N/C               | 3907-AS-NB |
| <input type="radio"/> Narrow Back w/Slings + \$ 415 | 3907-AS-NS |

Asepsis Additional Options:

- |   |          |
|---|----------|
| <input type="radio"/> Massage Option + \$ 555 | 3900-100 |
|---|----------|

Memory Comfort Selections:

- |  |            |
|--|------------|
| <input type="radio"/> Wide/Thin + \$ 825   | 3907-MC-WB |
| <input type="radio"/> Dual Color + \$ 1115 | 3907-MC-DC |

**DELIVERY SYSTEM REAR PIVOT/STANDARD IC..... 1500PI\*** \$ 5460

Handpiece illumination system

- |   |             |     |
|---|-------------|-----|
| One handpiece (5-hole), accommodates 1 add'l tubing | 1700-021-CS | 715 |
| Add'l handpiece illumination tubing (5-hole).....   | 0017-615-CS | 275 |
| One handpiece (6-pin), accommodates 1 add'l tubing  | 1700-019-CS | 715 |
| Add'l handpiece illumination tubing (6-pin).....    | 0017-603-CS | 195 |

Ancillaries use existing handpiece position (see page 57)

Write in selection here.....

Write in selection here.....

- |  |             |     |
|--|-------------|-----|
| Additional handpiece position, wet ..... | 1100-189-CS | 425 |
|--|-------------|-----|

- |  |             |     |
|--|-------------|-----|
| Additional handpiece position, air drive only..... | 1100-190-CS | 215 |
|--|-------------|-----|

Additional autoclavable standard HVE

- |   |             |    |
|---|-------------|----|
| with tubing, 5' & autoclavable oral dry cup ..... | 1400-113-CS | 90 |
|---|-------------|----|

- |                                   |          |     |
|-----------------------------------|----------|-----|
| Assistant's quad holder bar ..... | 0024-054 | 165 |
|-----------------------------------|----------|-----|

- |   |             |     |
|---|-------------|-----|
| Quick disconnect with flow control, 1/4" (specify air or water) | 0016-455-CS | 145 |
|---|-------------|-----|

- |                                |          |     |
|--------------------------------|----------|-----|
| Swivel tray, arm-mounted ..... | 1110-731 | 345 |
|--------------------------------|----------|-----|

- |                     |             |     |
|---------------------|-------------|-----|
| Chrome accents..... | 9000-998-CS | 140 |
|---------------------|-------------|-----|

- |  |             |    |
|--|-------------|----|
| Additional umbilical (maximum 9 feet) (per foot) ..... | 1115-999-CS | 30 |
|--|-------------|----|

- |   |             |     |
|---|-------------|-----|
| Routing valve kit w/ water shutoff & regulator/factory installation | 1115-315-CS | 305 |
|---|-------------|-----|

- |  |          |      |
|--|----------|------|
| <b>LIGHT, HALOGEN</b> (mounting options, see pages 33) ..... | ○ 9083 * | 2195 |
|--|----------|------|

- |                     |        |      |
|---------------------|--------|------|
| <b>OR LED</b> ..... | ○ 9078 | 2995 |
|---------------------|--------|------|

• Light post, pivot (includes chair adaptor) .....

• Light post, pivot if sold with monitor post (does not include chair adaptor) .....

- |       |          |     |
|-------|----------|-----|
| ..... | 1112-644 | 600 |
|-------|----------|-----|

- |  |            |      |
|--|------------|------|
| Monitor post (for top post mount) (includes chair adaptor) | 1331-644-3 | 1185 |
|--|------------|------|

- |  |  |  |
|--|--|--|
| Monitor mount assembly (monitor screen not included) (see page 37) |  |  |
|--|--|--|

- |  |        |      |
|--|--------|------|
| Ring mount (used with light post)..... | 1300RM | 1065 |
|--|--------|------|

- |  |         |      |
|--|---------|------|
| Top post mount (used with monitor post)..... | 1300TPM | 1065 |
|--|---------|------|

- |                               |          |      |
|-------------------------------|----------|------|
| <b>CHAIR, HYDRAULIC</b> ..... | ○ 3900 * | 7990 |
|-------------------------------|----------|------|

• One chair control included (specify type below) .....

• Footswitch control .....

• Touchpad, chair-mounted .....

• Touchpad, unit-mounted .....

• Additional chair control (limit 2 on chair, limit 2 on unit) .....

- |             |   |     |
|-------------|---|-----|
| 3914-999-CS | * | 515 |
|-------------|---|-----|

• Footswitch control .....

• Touchpad, chair-mounted .....

• Touchpad, unit-mounted .....

Additional chair options (See page 25 for pricing, write selection here)

- |  |                  |             |
|--|------------------|-------------|
| <b>STOOL, DR'S, ECONOMY/DELUXE/COMFORT</b> | 6156C/6159C/6176 | 530/895/895 |
|--|------------------|-------------|

- |  |             |          |
|--|-------------|----------|
| <b>STOOL, ASST'S, ECONOMY/DELUXE</b> ..... | 6157C/6160C | 680/1010 |
|--|-------------|----------|

- |  |            |           |
|--|------------|-----------|
| <b>STOOL, ASST'S WITH BACK, DELUXE/COMFORT</b> ..... | 6161C/6187 | 1195/1195 |
|--|------------|-----------|

• Fully assembled, blanket wrapped (chair, delivery, light).....

BLANKET .....

450 .....

\* ITEM included in AS SHOWN PACKAGE price

• Choose only one post when sold with unit

• Often overlooked

**TOTAL**

**Equipment dimensions, see page 96**

Shipping weight estimate:

Delivery unit & chair 445 lb / 203 kg

Base package 533 lb / 242 kg

As shown package 538 lb / 245 kg

# Pivot Chair Mount Packages

## Fixed Chair Mount with Sidebox

### Standard Features:

#### **Delivery System 4485SI**

Automatic control for three handpiece  
asepsis 3-way air/water syringe w/tubing, 5'  
handpiece tubing, 4-hole Midwest, asepsis style, 5'  
pneumatic master control block w/diaphragm  
master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air coolant adjustment  
handpiece flush  
handpiece exhaust canister/oil collector  
single handle brake release function  
wet/dry foot control  
Flexible control arm w/pneumatic brake  
vertical travel 28" w/horizontal travel 62"  
Sidebox w/self-contained clean water system  
w/quick disconnect, single bottle, 1 1/2 liter  
Water quick disconnect 1/4" w/flow control  
Air quick disconnect 3/8" w/flow control  
Toggle valve for city/bottle water routing  
Telescoping arm w/swivel quad holder bar  
Central vacuum accessory package  
central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'  
Cable for unit-mounted touchpad, 3900 chair  
Internal umbilical, 5' w/3' external from chair to utility center  
Contoured utility center  
air/water pressure regulators w/filter  
master air/water shut-off valves  
Chair adapter, fixed

#### **Cuspidor, Deluxe 1577** (option not shown)

One-piece ceramic bowl w/touch lever water activator  
Automatic cup fill  
Timed bowl rinse

#### **Cuspidor, Basic 1679** (option not shown)

Removable ceramic bowl  
Manual cup fill  
Timed bowl rinse

#### **Light, Halogen 9083**

3-Axis light head movement  
Automatic on/off activation w/override  
Color shade matching temperature setting 4200 to 5500 Kelvin  
Shadow free illumination

#### **Chair, Hydraulic 3900**

Ultraleather® upholstery w/wide back for elbow support  
Double articulating pinch-style headrest  
Thin back w/full lumbar support  
Right/left patient access  
Swing-down locking armrests  
Hydraulic cantilever mechanism  
Emergency shut-off  
30 degree base rotation w/lock  
Multifunction touchpad w/4 programmable preset positions  
Seat height: minimum 17" maximum 29"



### Popular Equipment Choices



Cuspidor 1577  
(Sidebox Comes  
With Delivery)



Cuspidor 1679  
(Sidebox Comes  
With Delivery)



Monitor Mount  
1300RM



LED Light  
9078



Doctor's Comfort  
Stool 6176



Assistant's Comfort  
Stool 6187



Doctor's Deluxe  
Stool 6159C



Assistant's Deluxe  
Stool 6161C

## Factory Options

## Part No

## Retail

## Price

**Base Package** \$ 18,380

4485SI Delivery System  
3900 Chair  
9083 Light  
1110-900 Light Post

**As Shown Package** \$ 18,725

Additionally equipped with:  
1110-468 Swivel Tray, Arm-Mounted  
3914-998-CS Chair Control, Touchpad, Unit-Mounted

## Designer Friendly® Surface Treatments

**STEP 1 SELECT UNIT PLASTIC FINISHES:**

- Forest Grey  Greywood  Laser Leather
- Lemans  White

**STEP 2 SELECT POWDER COAT PAINT FOR UNIT & LIGHT:**

- |                                     |                                   |                                |
|-------------------------------------|-----------------------------------|--------------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze      | <input type="radio"/> Charcoal |
| <input type="radio"/> Copper        | <input type="radio"/> Forest Grey | <input type="radio"/> Midnight |
| <input type="radio"/> Pewter        | <input type="radio"/> Shadow      | <input type="radio"/> White    |

**STEP 3 SELECT STANDARD POWDER COAT PAINT (FOR CHAIR):**

- Forest Grey  Shadow

**UPGRADE:** Chair Base/Back Powder Coat Paint (+ \$ 1,500):**UPGRADE:** PRO Chair (+ \$ 2,570) (Less Utility Center Housing and Cover)

- |                                     |                                |                                |
|-------------------------------------|--------------------------------|--------------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze   | <input type="radio"/> Charcoal |
| <input type="radio"/> Copper        | <input type="radio"/> Midnight | <input type="radio"/> Pewter   |
| <input type="radio"/> White         |                                |                                |

**STEP 4 SELECT UTRALEATHER® UPHOLSTERY COLOR CHOICES FOR CHAIR AND DELUXE/ECONOMY STOOLS:**

- |                                     |                               |                                 |                                |
|-------------------------------------|-------------------------------|---------------------------------|--------------------------------|
| <input type="radio"/> Apricot       | <input type="radio"/> Baltic  | <input type="radio"/> Cave      | <input type="radio"/> Charcoal |
| <input type="radio"/> Diplomat Blue | <input type="radio"/> Doe     | <input type="radio"/> Dove Grey |                                |
| <input type="radio"/> Eucalyptus    | <input type="radio"/> Papyrus | <input type="radio"/> Parrot    |                                |
| <input type="radio"/> Raven Wing    | <input type="radio"/> Red     | <input type="radio"/> Riviera   | <input type="radio"/> Stone    |

Pearlized Colors:  Abalone  Midnight  Pewter**OR SELECT** Special Color UltraLeather® (Call for leadtime)

- |   |             |
|---|-------------|
| <input type="radio"/> Chair + \$ 290                | 3907-999-SP |
| <input type="radio"/> Deluxe/Economy Stool + \$ 105 | 6007-999-SP |

**STEP 5 UPHOLSTERY SELECTIONS:**

Plush Selections:

- |   |            |
|---|------------|
| <input type="radio"/> Standard Back, Wide/Thin N/C  | 3907-PL-WB |
| <input type="radio"/> Narrow Back N/C               | 3907-PL-NB |
| <input type="radio"/> Narrow Back w/Slings + \$ 420 | 3907-PL-NS |

Plush Additional Options:

- |   |             |
|---|-------------|
| <input type="radio"/> Contrast/Deviated Stitching + \$ 290    | 3907-999-DS |
| <input type="radio"/> Dual Color + \$ 290                     | 3907-999-DC |
| <input type="radio"/> Dual Color /Contrast/Deviated Stitching |             |

Call Customer Service for information (1.800.423.3555)

- |   |          |
|---|----------|
| <input type="radio"/> Massage Option + \$ 555 | 3900-100 |
|---|----------|

Asepsis Selections:

- |   |            |
|---|------------|
| <input type="radio"/> Standard Back, Wide/Thin N/C  | 3907-AS-WB |
| <input type="radio"/> Narrow Back N/C               | 3907-AS-NB |
| <input type="radio"/> Narrow Back w/Slings + \$ 415 | 3907-AS-NS |

Asepsis Additional Options:

- |   |          |
|---|----------|
| <input type="radio"/> Massage Option + \$ 555 | 3900-100 |
|---|----------|

Memory Comfort Selections:

- |  |            |
|--|------------|
| <input type="radio"/> Wide/Thin + \$ 825   | 3907-MC-WB |
| <input type="radio"/> Dual Color + \$ 1115 | 3907-MC-DC |

**DELIVERY SYSTEM WITH STANDARD IC** ..... 4485SI\* \$ 7975

Handpiece illumination system

- |   |             |     |
|---|-------------|-----|
| One handpiece (5-hole), accommodates 1 add'l tubing | 1700-021-CS | 715 |
| Add'l handpiece illumination tubing (5-hole).....   | 0017-615-CS | 275 |
| One handpiece (6-pin), accommodates 1 add'l tubing  | 1700-019-CS | 715 |
| Add'l handpiece illumination tubing (6-pin).....    | 0017-603-CS | 195 |

Ancillaries (use existing handpiece position) (see page 57)

Write in selection here .....

Write in selection here .....

- |  |             |     |
|--|-------------|-----|
| Additional handpiece position, wet ..... | 1100-189-CS | 425 |
| Additional handpiece position .....      | 1100-190-CS | 215 |

Swivel tray, arm-mounted option (over control head) .... 1110-468\*

Less sidebox telescoping arm and assistant's instrumentation 1400-000-CS -350

Chrome accents .....
 9000-998-CS | 140 |

Additional autoclavable standard HVE

with tubing, 5' & autoclavable oral dry cup .....
 1400-113-CS | 90 |

Additional umbilical (maximum 9 feet) (per foot) .... 1115-999-CS 30

**CUSPIDOR** (sold only with delivery)**DELUXE** ..... ○ 1577 1780**OR BASIC** ..... ○ 1679 1070

Deluxe cuspidor unit-mounted activation (1577) .... 1100-480-CS 270

**LIGHT, HALOGEN** (mounting options, see pages 33) ..... ○ 9083\* 2195**OR LED** ..... ○ 9078 2995

• Light post, 44" ..... 1110-900\* 220

• Monitor post (for top post mount)..... 1331-998 220

Monitor mount assembly (monitor screen not included) (see page 37)

Ring mount (used with light post)..... 1300RM 1065

Top post mount (used with monitor post)..... 1300TPM 1065

**CHAIR, HYDRAULIC** ..... ○ 3900\* 7990

One chair control included (specify type below).... 3914-998-CS \* N/C

○ Footswitch control

○ Touchpad, chair-mounted

○ Touchpad, unit-mounted

Additional chair control (limit 2 on chair, limit 2 on unit) 3914-999-CS 515

○ Footswitch control

○ Touchpad, chair-mounted

○ Touchpad, unit-mounted

Additional chair options (See page 25 for pricing, write selection here)

**STOOL, DR'S, ECONOMY/DELUXE/COMFORT** 6156C/6159C/6176 530/895/895**STOOL, ASST'S, ECONOMY/DELUXE** ..... 6157C/6160C 680/1010**STOOL, ASST'S WITH BACK, DELUXE/COMFORT** ..... 6161C/6187 1195/1195**STOOL, OFFICE, COMFORT** ..... 6198 1060

Fully assembled, blanket wrapped (chair, delivery, light) .. BLANKET 450

\* ITEM included in AS SHOWN PACKAGE price

• Choose only one post when sold with unit

● Often overlooked

**TOTAL****Equipment dimensions, see page 97**

Shipping weight estimate:

Delivery unit &amp; chair 440 lb / 202 kg

Base package 556 lb / 253 kg

As shown package 561 lb / 255 kg

# Fixed Chair Mount With Sidebox Packages

## Euro Control Head

### Standard Features:

#### **Delivery System 4987SI**

Automatic control for four handpieces

asepsis 3-way air/water syringe w/tubing, 3' handpiece tubing, 4-hole Midwest, asepsis style, 3' pneumatic master control block w/diaphragm master on/off

individual handpiece pressure adjustment

individual water coolant adjustment

non-retraction water coolant

air coolant adjustment

handpiece flush

handpiece exhaust canister/oil collector

dual handle brake release function

swivel instrument tray

wet/dry foot control

Flexible control arm w/pneumatic brake

vertical travel 28"/horizontal travel 62"

Sidebox w/self-contained clean water system

w/quick disconnect, single bottle, 1 1/2 liter

Water quick disconnect 1/4" w/flow control

Air quick disconnect 3/8" w/flow control

Toggle valve for city/bottle water routing

Telescoping arm w/swivel quad holder bar

Central vacuum accessory package

central vacuum canister w/disposable solids collector

asepsis 3-way air/water syringe w/tubing, 5'

autoclavable standard HVE lever valve w/tubing, 5'

autoclavable saliva ejector lever valve w/tubing, 5'

Cable for unit-mounted touchpad, 3900 chair

Internal umbilical, 5' w/3' external from chair to utility center

Contoured utility center

air/water pressure regulators w/filter

master air/water shut-off valves

Chair adapter, fixed

#### **Cuspidor, Deluxe 1577** (option not shown)

One-piece ceramic bowl w/touch lever water activator

Automatic cup fill

Timed bowl rinse

#### **Cuspidor, Basic 1679** (option not shown)

Removable ceramic bowl

Manual cup fill

Timed bowl rinse

#### **Light, Halogen 9083**

3-Axis light head movement

Automatic on/off activation w/override

Color shade matching temperature setting 4200 to 5500 Kelvin

Shadow free illumination

#### **Chair, Hydraulic 3900**

Ultraleather® upholstery w/wide back for elbow support

Double articulating pinch-style headrest

Thin back w/full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation w/lock

Multifunction touchpad w/4 programmable preset positions

Seat height: minimum 17" maximum 29"



### Popular Equipment Choices



Cuspidor 1577  
(Sidebox Comes  
With Delivery)



Cuspidor 1679  
(Sidebox Comes  
With Delivery)



Monitor Mount  
1300RM



LED Light  
9078



Doctor's Comfort  
Stool 6176



Assistant's Comfort  
Stool 6187



Doctor's Deluxe  
Stool 6159C



Assistant's Deluxe  
Stool 6161C

## Factory Options

## Part No

## Retail

## Price

**Base Package** \$ 19,455

4987SI Delivery System  
3900 Chair  
9083 Light  
1110-900 Light Post

**As Shown Package** \$ 19,455

3914-998-CS Chair Control, Chair Mounted

**DELIVERY SYSTEM WITH EURO** ..... 4987SI\* \$ 9050

Handpiece illumination system

One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS 715  
Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS 275  
One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS 715  
Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS 195

**Designer Friendly® Surface Treatments****STEP 1 SELECT UNIT PLASTIC FINISHES:**

- Forest Grey  Greywood  Laser Leather
- Lemans  White

**STEP 2 SELECT POWDER COAT PAINT FOR UNIT & LIGHT:**

- |                                     |                                   |                                |
|-------------------------------------|-----------------------------------|--------------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze      | <input type="radio"/> Charcoal |
| <input type="radio"/> Copper        | <input type="radio"/> Forest Grey | <input type="radio"/> Midnight |
| <input type="radio"/> Pewter        | <input type="radio"/> Shadow      | <input type="radio"/> White    |

**STEP 3 SELECT STANDARD POWDER COAT PAINT (FOR CHAIR):**

- Forest Grey  Shadow

**UPGRADE:** Chair Base/Back Powder Coat Paint (+ \$ 1,500):

**UPGRADE:** PRO Chair (+ \$ 2,570) (Less Utility Center Housing and Cover)

- |                                     |                                |                                |
|-------------------------------------|--------------------------------|--------------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze   | <input type="radio"/> Charcoal |
| <input type="radio"/> Copper        | <input type="radio"/> Midnight | <input type="radio"/> Pewter   |
| <input type="radio"/> White         |                                |                                |

**STEP 4 SELECT UTRALEATHER® UPHOLSTERY COLOR CHOICES FOR CHAIR AND DELUXE/ECONOMY STOOLS:**

- |                                     |                               |                                 |                                |
|-------------------------------------|-------------------------------|---------------------------------|--------------------------------|
| <input type="radio"/> Apricot       | <input type="radio"/> Baltic  | <input type="radio"/> Cave      | <input type="radio"/> Charcoal |
| <input type="radio"/> Diplomat Blue | <input type="radio"/> Doe     | <input type="radio"/> Dove Grey |                                |
| <input type="radio"/> Eucalyptus    | <input type="radio"/> Papyrus | <input type="radio"/> Parrot    |                                |
| <input type="radio"/> Raven Wing    | <input type="radio"/> Red     | <input type="radio"/> Riviera   | <input type="radio"/> Stone    |

Pearlized Colors:  Abalone  Midnight  Pewter

**OR SELECT** Special Color Ultraleather® (Call for leadtime)

- |   |             |
|---|-------------|
| <input type="radio"/> Chair + \$ 290                | 3907-999-SP |
| <input type="radio"/> Deluxe/Economy Stool + \$ 105 | 6007-999-SP |

**STEP 5 UPHOLSTERY SELECTIONS:**

Plush Selections:

- |   |            |
|---|------------|
| <input type="radio"/> Standard Back, Wide/Thin N/C  | 3907-PL-WB |
| <input type="radio"/> Narrow Back N/C               | 3907-PL-NB |
| <input type="radio"/> Narrow Back w/Slings + \$ 420 | 3907-PL-NS |

Plush Additional Options:

- |   |             |
|---|-------------|
| <input type="radio"/> Contrast/Deviated Stitching + \$ 290    | 3907-999-DS |
| <input type="radio"/> Dual Color + \$ 290                     | 3907-999-DC |
| <input type="radio"/> Dual Color /Contrast/Deviated Stitching |             |

Call Customer Service for information (1.800.423.3555)

- |   |          |
|---|----------|
| <input type="radio"/> Massage Option + \$ 555 | 3900-100 |
|---|----------|

Asepsis Selections:

- |   |            |
|---|------------|
| <input type="radio"/> Standard Back, Wide/Thin N/C  | 3907-AS-WB |
| <input type="radio"/> Narrow Back N/C               | 3907-AS-NB |
| <input type="radio"/> Narrow Back w/Slings + \$ 415 | 3907-AS-NS |

Asepsis Additional Options:

- |   |          |
|---|----------|
| <input type="radio"/> Massage Option + \$ 555 | 3900-100 |
|---|----------|

Memory Comfort Selections:

- |  |            |
|--|------------|
| <input type="radio"/> Wide/Thin + \$ 825   | 3907-MC-WB |
| <input type="radio"/> Dual Color + \$ 1115 | 3907-MC-DC |

**DELIVERY SYSTEM WITH EURO** ..... 4987SI\* \$ 9050

Handpiece illumination system

One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS 715  
Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS 275  
One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS 715  
Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS 195

**(i) Ancillaries (use existing handpiece position) (see page 57)**

Write in selection here .....

Write in selection here .....

Less sidebox telescoping arm and assistant's instrumentation 1400-000-CS -350

Chrome accents ..... 9000-998-CS 140

Additional autoclavable standard HVE

with tubing, 5' & autoclavable oral dry cup..... 1400-113-CS 90

Additional umbilical (maximum 9 feet) (per foot) .... 1115-999-CS 30

**CUSPIDOR (sold only with delivery)**

**DELUXE** ..... ○1577 1780

**OR BASIC** ..... ○1679 1070

Deluxe cuspidor unit-mounted activation (1577) .... 1100-480-CS \* 270

**LIGHT, HALOGEN** (mounting options, see pages 33) ..... ○ 9083 2195

**OR LED** ..... ○ 9078\* 2995

**(i) Light post, 44"** .....

• Monitor post (for top post mount)..... 1331-998 220

Monitor mount assembly (monitor screen not included) (see page 37)

Ring mount (used with light post)..... 1300RM 1065

Top post mount (used with monitor post)..... 1300TPM\* 1065

**CHAIR, HYDRAULIC** ..... 3900\* 7990

One chair control included (specify type below).... 3914-998-CS N/C

Footswitch control

Touchpad, chair-mounted

Touchpad, unit-mounted

**(i) Additional chair control (limit 2 on chair, limit 2 on unit) .** 3914-999-CS 515

Footswitch control

Touchpad, chair-mounted

Touchpad, unit-mounted

Additional chair options (See page 25 for pricing, write selection here)

**STOOL, DR'S, ECONOMY/DELUXE/COMFORT** ..... 6156C/6159C/6176 530/895/895

**STOOL, ASST'S, ECONOMY/DELUXE** ..... 6157C/6160C 680/1010

**STOOL, ASST'S WITH BACK, DELUXE/COMFORT** ..... 6161C/6187 1195/1195

**STOOL, OFFICE, COMFORT** ..... 6198 1060

Fully assembled, blanket wrapped (chair, delivery, light) .. BLANKET 450

\* ITEM included in AS SHOWN PACKAGE price

• Choose only one post when sold with unit

**(i) Often overlooked**

**TOTAL**

**Equipment dimensions, see page 97**

Shipping weight estimate:

Delivery unit & chair 445 lb / 203 kg

Base package 551 lb / 251 kg

As shown package 551 lb / 251 kg

## Surface Treatments

Forest is the only dental equipment with a WOW factor. It really is a kick watching a dentist react to Forest equipment when they first see it. A common reaction is "wow, that is beautiful."

### Carbon Fiber

#### Control Heads

- Standard IC - Chassis
- Euro - Lid

Sidebox - Cover

Basic Cuspidor - Cover

Assistant's Instrumentation - Arm Cover

### Powder Coat Paint

Dentists and interior designers love being able to put something other than white into their beautiful patient friendly decors.

#### Control Heads

- Standard IC - Lid, Brake Handles and Holder Bar
- Euro - Chassis, Whip Arms, Brake Handles and Swivel Tray
- 7000 Series Rear Mount, Side Delivery and Mobile Cart
- 7020 PRO Series Rear Mount, Side Delivery and Mobile Cart

Sidebox - Housing

Basic Cuspidor - Housing and Holder Bar

Assistant's Instrumentation - Arm

### Ultraleather® (standard)

Once you see a Forest chair, delivery system and light package with color coordinated paint and upholstery in a well designed treatment room, you wonder why it's taken so long for a manufacturer to offer dental equipment with a full range of surface treatments.

#### 3900 Chair

- Plush Upholstery
- Asepsis Upholstery
- Knee and Neck Pillow
- Deluxe and Economy Stools

### Solid Surface

#### Control Heads

- 7000 Series Rear Mount Side Delivery and Mobile Cart (standard)
- 7020 PRO Series Rear Mount Side Delivery and Mobile Cart (standard)
- 3486 Rear Mount Side Delivery (standard)
- 5386/5886 Mobile Cart (optional)
- 2184BC/2186BC Rear and Side Delivery (optional)

### Laminate

#### Control Heads

- 3486 Rear Mount Side Delivery (optional)
- 5386/5886 Mobile Cart (standard)
- 2184BC/2186BC Rear and Side Delivery (standard)



---

Chairs and Stools



## BUILT TO LAST

New doesn't have to mean plastic. While our competitors continue to move toward more and more plastic parts, Forest uses metal whenever possible so you get greater value from your investment.

### Patient Chair

Material: cast aluminum chair back, chair base & utility center base

Benefit: strong, impervious to fading or cracking

Material: pressure-formed plastic chair & utility center covers

Benefit: warp, impact & stain resistant

Material: aluminum hydraulic manifold & speed adjustment screws

Benefit: independent speed adjustment of base & back both up & down

Material: optional soft-touch back massage with integrated compressor

Benefit: saves central compressor capacity for delivery systems

Material: Ultraleather® upholstery standard

Benefit: withstands heavy use & cleaning, with glove-soft leather feel

Material: forest grey or shadow powder coat standard on 3900 chair / seven optional Designer Friendly® colors

Benefit: complements treatment room interior design color palette

### Comfort Stool

Material: tension woven material with rocking back

Benefit: eliminates pressure points, massages & supports the seat & back

Material: all-surface 22mm polyurethane casters

Benefit: smooth rolling eliminates squeaks / long-life

(Shown on Front)

**3900 Patient Chair** Powder Coat Paint Color: Shadow

**Upholstery** Plush: Ultraleather® Pro Blueberry and Green Apple

**Neck Pillow**

Provides ergonomic support for patients with neck discomfort and/or during long procedures.  
Neck support offers patients a comfortable relaxed experience.  
Specify upholstery color

3907-318 \$ 125

**Knee Pillow**

Knee support ideal for patients with joint discomfort and/or long procedures.  
Zero gravity posture, eliminates pressure on knees.  
Specify upholstery color

3907-317 \$ 390

Dual color upholstery option (upcharge)  
Contrast/Deviated stitching option (upcharge)

3907-317-DC 80  
3907-317-DS 80**Pedo Booster Seat** (to accommodate smaller patients)

A cushion that takes the lumbar out of the chair, essentially turning it into a bed, ideal for children.  
Specify upholstery color

3907-180 \$ 365

**Chair Control, Footswitch**

Additional Chair Control, Footswitch



Factory Installed	3914-998-CS	N/C
Factory Installed	3914-999-CS	\$ 515
After Market	3914-010-SH	515

4 programmable preset positions

**3900 Remote Mount Touchpad with 25' Cable****FIXED HORIZONTAL MOUNT**

White	1700-079	\$ 645
Forest Grey	1700-079-G	645
Designer Friendly® (specify color)	1700-079-C	645

**SWING HORIZONTAL MOUNT**

White	1700-089-HM	645
Forest Grey	1700-089-HM-G	645
Designer Friendly® (specify color)	1700-089-HM-C	645

**SWING VERTICAL MOUNT**

White	1700-089-VM	645
Forest Grey	1700-089-VM-G	645
Designer Friendly® (specify color)	1700-089-VM-C	645

- 4 programmable preset positions
- Included w/4485PI, 4987PI, 1500PI, 4485SI and 4987SI

**Contoured Utility Center Housing & Cover**

Factory option:	3900-101	\$ 435
Forest Grey	3900-101-G	435
Shadow	3900-101-SH	435
Designer Friendly® (specify color)	3900-101-C	535

Included w/PI and SI chair mounted packages

# Chairs and Stools

## 3900 Hydraulic Patient Chair




**Plush**



**Model 3900 \$ 7990**

**Asepsis**



**Part No**   **Retail**   **Price**

### Standard Features

#### Hydraulic Chair 3900

Ultraleather® upholstery w/  
wide/thin back for elbow support  
Plastic toe guard, 25"  
Double articulating pinch-style headrest  
Thin back with full lumbar support  
Right/left patient access  
Swing-down locking armrests  
Hydraulic cantilever mechanism  
Emergency shut-off  
30 degree base rotation with lock  
Multifunction touchpad  
w/4 programmable preset positions,  
accommodates 3 additional controls  
Seat height: minimum 17" maximum 29"

**CHAIR** (Forest Grey or Shadow)..... 3900 \$ 7990

**CHAIR** Designer Friendly® ..... 3900C 9490

**CHAIR PRO** Designer Friendly® Includes: Plush dual color and deviated/  
contrast stitching upholstery, massage, contoured utility center housing  
and cover, knee pillow and neck pillow ..... 3900C-PRO  
**CHAIR PRO** Designer Friendly® less utility center housing and cover

..... 3900C-PRO-LUC

1170 10560

One chair control included (specify type below)..... 3914-998-CS

N/C

Footswitch control

Touchpad, chair-mounted

Additional chair control (limit 2 on chair, limit 2 on unit) 3914-999-CS

515

Footswitch control

Touchpad, chair-mounted

Chair adapter, fixed (for external umbilical) ..... 1115-195

445

Chair adapter, fixed (for internal umbilical) ..... 1115-530

585

Chair adapter, pivot (for internal umbilical) ..... 1115-528

585

Contoured utility center housing & cover ..... 3900-101

435

Contoured utility center housing & cover

Designer Friendly® ..... 3900-101-C

535

### Designer Friendly® Surface Treatments

#### STEP 1 SELECT STANDARD POWDER COAT PAINT (CHAIR):

Forest Grey    Shadow

**UPGRADE:** Chair Base/Back Powder Coat Paint (+ \$ 1,500):

**UPGRADE:** PRO Chair (+ \$ 2,570)

(Less Utility Center Housing and Cover)

<input type="radio"/> Arctic Silver	<input type="radio"/> Bronze	<input type="radio"/> Charcoal
<input type="radio"/> Copper	<input type="radio"/> Midnight	<input type="radio"/> Pewter
<input type="radio"/> White		

#### UPHOLSTERY SELECTIONS:

Write in selection here (see page 27) .....

.....

Write in selection here (see page 27) .....

.....

Pedo booster seat (see page 25) ..... 3907-180

365

Knee pillow (see page 25) ..... 3907-317

390

Neck pillow (see page 25) ..... 3907-318

125

Soft-touch back massage option

(not available with Memory Comfort) ..... 3900-100

555

Plastic toe guard 33" in lieu of 25" ..... 3907-028-CS

25

#### STEP 2 SELECT UTRALEATHER® UPHOLSTERY COLOR CHOICES

FOR CHAIR AND DELUXE /ECONOMY STOOLS:

<input type="radio"/> Apricot	<input type="radio"/> Baltic	<input type="radio"/> Cave	<input type="radio"/> Charcoal
<input type="radio"/> Diplomat Blue	<input type="radio"/> Doe	<input type="radio"/> Dove Grey	
<input type="radio"/> Eucalyptus	<input type="radio"/> Papyrus	<input type="radio"/> Parrot	
<input type="radio"/> Raven Wing	<input type="radio"/> Red	<input type="radio"/> Riviera	<input type="radio"/> Stone

Pearlized Colors:  Abalone  Midnight  Pewter

**OR SELECT** Special Color Ultraleather® (Call for leadtime)

Chair 3907-999-SP + \$ 290

#### LIGHT POSTS (for unit mount lights/post mount monitor assemblies see pages 33)

Fixed (chair adaptor sold separately)

44" light post without slot, chair mount ..... 1110-794

220

44" light post with slot for monitor cord ..... 1110-050

220

Pivot (includes chair adaptor)

Pivot light post, 3900 (light not included) ..... 1112-644-3

1185

Pivot monitor post, 3900

(monitor mount assembly not included) ..... 1331-644-3

1185

Monitor mount assembly, top post mount (see page 37) .. 1300TPM

1065

Monitor mount assembly, ring mount (see page 37) ..... 1300RM

1065

**TOTAL** \_\_\_\_\_

## 3900 Chair Upholstery Options

**Memory Comfort****Narrow Back or Narrow Back with Slings**

Standard Features	Factory Options	Part No	Retail	Price
<b>Memory Comfort Upholstery</b> Designed for all day/ultra-long procedures or greater general patient comfort 3 layers of foam topped with 1.5 inches of memory foam for extra comfort Side-bolsters on seat cushion support legs Memory foam headrest	<b>UPHOLSTERY SELECTIONS:</b> Plush selections: Standard back, wide/thin ..... 3907-PL-WB N/C Narrow back ..... 3907-PL-NB N/C Narrow back with slings ..... 3907-PL-NS \$ 420  Plush additional options: Contrast/deviated stitching ..... 3907-999-DS 290 _____ Dual color ..... 3907-999-DC 290 _____  Dual color with contrast/deviated stitching Call customer service for information (1.800.423.3555) Massage option ..... 3900-100 555 _____			
<b>Narrow Back</b> Perfect for Pediatric or Orthodontic dentistry Unrestricted patient access	Asepsis selections: Standard back, wide/thin ..... 3907-AS-WB N/C Narrow Back ..... 3907-AS-NB N/C Narrow Back with Slings ..... 3907-AS-NS 415 _____			
<b>Narrow Back with Slings</b> Same as above, with additional features: Slings provide alternative for additional patient elbow support Slings create a barrier between patient and practitioner	Asepsis additional options: Massage option ..... 3900-100 555 _____  Memory comfort selections: Wide/thin ..... 3907-MC-WB 825 _____ Dual color ..... 3907-MC-DC 1115 _____			
			<b>TOTAL</b>	_____



3900 Chair Equipment dimensions, see page 98

Shipping weight estimate:

335 lb / 152 kg

# Chairs + Stools

## Economy and Deluxe Stools

**Model 6156C \$ 530****Model 6157C \$ 680****Model 6159C \$ 895****Model 6160C \$ 1010****Model 6161C \$ 1195****Standard Features****Factory Options****Part No****Retail****Price****Economy and Deluxe Stools**

Ultraleather® Upholstery  
Injection molded foam  
5-leg fully polished base  
Twin wheel hard floor casters  
Accommodates maximum of 300 lbs.

**Doctor's Economy 6156C**

Backrest, non-adjusting  
Round seat with single lever height adjustment  
Pneumatic cylinder height adjustment, 18" - 23"

**Assistant's Economy 6157C**

Round seat with single lever height adjustment  
Torso support, non-ratcheting, with height adjustment  
Pneumatic cylinder height adjustment, 21 1/2" - 28"  
Foot ring affixed to base

**Doctor's Deluxe 6159C**

Backrest with height adjustment  
Contoured seat with triple lever adjustment  
Pneumatic cylinder height adjustment, 18 1/2" - 23"

**Assistant's Deluxe 6160C**

Contoured seat with single lever adjustment  
Torso support with ratcheting and height adjustment  
Pneumatic cylinder height adjustment, 21 1/2" - 28"  
Foot ring with adjustment

**Assistant's Deluxe with Back 6161C**

Backrest with height adjustment  
Contoured seat with triple lever adjustment  
Torso support with ratcheting and height adjustment  
Pneumatic cylinder height adjustment, 21 1/2" - 28"  
Foot ring with adjustment

**STOOL ..... MODEL CHOICE ..... PRICE:**

Vinyl in lieu of Ultraleather® ..... 6007-999-CS \* \$ - 100  
Special color upholstery (upcharge) ..... 6007-999-SP 105

**► CYLINDER CHOICES:**

Doctor's short cylinder, 16 " - 21" (6156C/6159C) .... 6000-628-CS N/C  
Doctor's tall cylinder, 21 1/2" - 28" (6156C/6159C) foot ring  
with support highly recommended ..... 6000-604-CS N/C  
Assistant's short cylinder, 18" - 23"  
(6157C/6160C/6161C) ..... 6000-605-CS N/C  
Foot ring with support (6156C/6159C) ..... 6000-602 120

**TOTAL****► Often overlooked**

\* Specify vinyl manufacturer and color  
(not all colors may be available.)  
(Call Forest for ordering details 1-800-423-3555.)

**Equipment dimensions, see page 98**

Shipping weight estimate:

6156C/6159C 35 lb / 15 kg

6161C 45 lb / 21 kg

6157C/6160C 40 lb / 18 kg

**Designer Friendly® Surface Treatments****STEP 1 SELECT ULTRALEATHER® UPHOLSTERY COLOR CHOICES FOR CHAIR AND DELUXE /ECONOMY STOOLS:**

Apricot  Baltic  Cave  Charcoal  Diplomat Blue  Doe  Dove Grey  
 Eucalyptus  Papyrus  Parrot  Raven Wing  Red  Riviera  Stone

Pearlized Colors:  Abalone  Midnight  Pewter**OR SELECT Special Color Ultraleather® (Call for leadtime)** Deluxe/Economy Stool 6007-999-SP + \$ 105

**Model 6176 \$ 895****Model 6187 \$ 1195****Model 6198 \$ 1060**

Standard Features	Factory Options	Part No	Retail	Price
<b>Comfort Stools</b> Air circulating weave construction Active tilt® seat (actively moves with user) Tilt seat lever function Contoured seat 5-leg base Twin wheel hard floor casters Accommodates maximum of 300 lbs.	<b>STOOL .....</b> MODEL CHOICE _____ Cylinder choices: _____ Doctor's/Office short cylinder, 16 " – 21" (6176/6198) ... 6000-628-CS N/C Doctor's/Office tall cylinder, 21 1/2" – 28" (6176/6198) foot ring with support highly recommended ..... 6000-604-CS N/C Assistant's short cylinder, 18" – 23" (6187) ..... 6000-605-CS N/C Foot ring with support (6176/6198) ..... 6000-602 \$ 120			
<b>Doctor's and Office Comfort 6176 / 6198</b> Triple lever adjustment - independent back and seat Backrest with angle & height adjustment Fully adjustable office arms (6198) Pneumatic cylinder height adjustment, 19 1/4"– 24 1/2"	Often overlooked		TOTAL	_____
<b>Assistant's Comfort 6187</b> Triple lever adjustment - independent back and seat Backrest with angle & height adjustment Fully adjustable ratcheting assistant arm Foot ring with height adjustment Pneumatic cylinder height adjustment, 23"– 30 1/2"				

**Equipment dimensions, see page 98**

Shipping weight estimate:

6176 48 lb / 23 kg

6187 60 lb / 28 kg

6198 54 lb / 24 kg

## Notes

Lights and Monitor Mounts




## Powder Coat



Forest Grey

Shadow

White



Midnight

Arctic Silver

Pewter



Copper

Charcoal

Bronze

## Halogen

Material: reflector material embedded into glass / lexan light shield

Benefit: reflector will not peel or crack / shield is scratch resistant

Material: oil impregnated bronze bushing / 3-axis light head joint

Benefit: no light drift / 3-axis light keeps light pattern parallel with the long axis of mouth while freeing up space above patient's head

## LED

Material: ten 1-watt LED lamps

Benefit: creates a light pattern with feathered edges reducing eye fatigue /

consumes approximately 75 percent less power than halogen lights / no heat / 22 -year projected LED life span

Material: cast aluminum light head

Benefit: allows for nine Designer Friendly® powder coat paint colors

(Shown on Front)

**LED** Powder Coat Paint Color: Bronze

**Halogen** Powder Coat Paint Color: Shadow



Model 9083 \$ 2195



Model 9078 \$ 2995

Standard Features	Factory Options	Part No	Retail	Price
<b>Halogen Light 9083</b> Automatic on/off activation w/override 3-Axis light head movement Additional bulb Color shade matching temperature setting 4200 to 5500 Kelvin Designer Friendly® powder coat paint light head	<b>LIGHT</b> ..... <b>MODEL CHOICE</b> ..... <b>PRICE:</b> ..... Replacement bulb (9083/9093)..... 0014-300 \$ 40 Short rigid arm 20" in lieu of standard 32"..... 1110-227-CS N/C <b>LIGHT POSTS:</b> 44" with slot for sidebox access ..... 1110-900 220 44" without slot, chair mount ..... 1110-794 220 35" with post clamp ..... 1110-820 290 Pivot (mounts on Forest 3900 chair only)..... 1112-644 600 Pivot (includes 3900 cast chair adapter)..... 1112-644-3 1185 <b>LIGHT POSTS, WITH MONITOR CORD SLOT:</b> 44" with dual slot for sidebox access ..... 1110-170 220 44" with single slot, chair mount ..... 1110-050 220 35" with single slot & post clamp ..... 1110-045 290			
<b>LED Light 9078</b> Automatic on/off activation w/override 3-Axis light head movement Easy release autoclavable handles Light source 10 LED lamps Color temperature, 4,200 to 6,000 Kelvin Designer Friendly® powder coat paint light head				
			<b>TOTAL</b>	



Designer Friendly® Surface Treatments



Equipment dimensions, see page 99

Shipping weight estimate:

35 lb / 16 kg

## STEP 1 SELECT POWDER COAT PAINT FOR UNIT &amp; LIGHT

- |                                     |                                |                                |                              |
|-------------------------------------|--------------------------------|--------------------------------|------------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze   | <input type="radio"/> Charcoal | <input type="radio"/> Copper |
| <input type="radio"/> Forest Grey   | <input type="radio"/> Midnight | <input type="radio"/> Pewter   | <input type="radio"/> Shadow |

- White





**Model 9083-C \$ 2830**  
Halogen Light



**Model 9078-C \$ 3995**  
LED Light

### Lights

Ceiling mounting plate and shroud w/drop post  
Specify mounting height

#### Halogen Light 9083-C

Automatic on/off activation w/override  
3-Axis light head movement  
Additional bulb  
Color temperature, 5,500 Kelvin  
Designer Friendly® powder coat paint light head

#### LED Light 9078-C

Automatic on/off activation w/override  
3-Axis light head movement  
Easy release autoclavable handles  
Light source 10 LED lamps  
Color temperature, 4,200 to 6,000 Kelvin  
Designer Friendly® powder coat paint light head

**LIGHT . . . . . MODEL CHOICE . . . . . PRICE:**  
Replacement bulb (9083/9093) . . . . . 0014-300 \$ 40

**TOTAL** \_\_\_\_\_

### Standard Mounting Options

#### STEP 1 CHOOSE YOUR STANDARD MOUNTING OPTION:

##### MOUNTING HEIGHT:

- |  |                            |
|--|----------------------------|
| <input type="radio"/> 8'                           | 12"                        |
| <input type="radio"/> 8' 6"                        | 18"                        |
| <input type="radio"/> 9'                           | 24"                        |
| <input type="radio"/> 9' 6"                        | 30"                        |
| <input type="radio"/> 10'                          | 36"                        |
| <input type="radio"/> 10' 6"                       | 42"                        |
| <input type="radio"/> 11'                          | 48"                        |
| <input type="radio"/> _____ Custom mounting height | _____ * Custom post length |
- \* Cannot exceed 48"

##### POST LENGTH:

### Designer Friendly® Surface Treatments

#### STEP 2 SELECT POWDER COAT PAINT FOR UNIT & LIGHT:

- |                                     |                                |                                |                              |
|-------------------------------------|--------------------------------|--------------------------------|------------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze   | <input type="radio"/> Charcoal | <input type="radio"/> Copper |
| <input type="radio"/> Forest Grey   | <input type="radio"/> Midnight | <input type="radio"/> Pewter   | <input type="radio"/> Shadow |
| <input type="radio"/> White         |                                |                                |                              |



Equipment dimensions, see page 99

Shipping weight estimate:

50 lb / 23 kg

# Lights and Monitor Mounts

## Cabinet and Wall Mount Lights



**Model 9079C \$ 4730**  
(Standard) 9079C-T (Tall)  
LED Cabinet

**Model 9093C \$ 3355**  
(Standard) 9093C-T (Tall)  
Halogen Cabinet

**Model 9093W \$ 3355**  
(Standard) 9093W-T (Tall) Halogen Wall

**Model 9079W \$ 4730**  
(Standard) 9079W-T (Tall) LED Wall

### Standard Features

### Factory Options

### Part No

### Retail

### Price

#### **Lights**

Cabinet mount, internal vertical mount w/post  
Wall mount, vertical mounting plate  
Standard rigid arm  
Tall rigid arm (additional 1' drop)

#### **Halogen Light 9093**

Automatic on/off activation w/override  
3-Axis light head movement  
Additional bulb  
Color shade matching temperature  
setting 4200 to 5500 Kelvin  
Designer Friendly® powder coat paint light head

#### **LED Light 9079**

Automatic on/off activation w/override  
3-Axis light head movement  
Easy release autoclavable handles  
Light source 10 LED lamps  
Color temperature, 4,200 to 6,000 Kelvin  
Designer Friendly® powder coat paint light head

**LIGHT** ..... **MODEL CHOICE** \_\_\_\_\_ **PRICE:** \_\_\_\_\_  
Replacement bulb (9083/9093) ..... 0014-300 \$ 40

Cabinet mounting options: (for dimensions see page 100)  
Horizontal top cabinet mount in lieu of internal mount 1110-881-CS N/C  
Vertical exterior cabinet mount in lieu of internal mount 1110-886-CS N/C

**TOTAL** \_\_\_\_\_

► \*\* SPECIFY CABINET/WALL MOUNTING HEIGHT \*\* \_\_\_\_\_

► Often overlooked



Equipment dimensions, see page 100

Shipping weight estimate:

Cabinet 60 lb / 28 kg

Wall 65 lb / 30 kg



Designer Friendly® Surface Treatments

#### STEP 1 SELECT POWDER COAT PAINT FOR UNIT & LIGHT:

Arctic Silver

Bronze

Charcoal

Copper

Forest Grey

Midnight

Pewter

Shadow

White

**Model 1300RM Ring Post Mount \$ 1065****Model 1300TPM Top Post Mount \$ 1065**

Monitor Screens And Chair-Mounted Posts Not Included

Standard Features	Factory Options	Part No	Retail	Price
<b>Mount Assembly,</b> <b>Ring Post Mount 1300RM</b> ICW Quicklink monitor mount Rex positioning handles 22" ring mount rigid arm	<b>MONITOR MOUNT ASSEMBLY .....</b> MODEL CHOICE _____ PRICE : _____			
<b>Monitor Mount Assembly,</b> <b>Top Post Mount 1300TPM</b> ICW Quicklink monitor mount Rex positioning handles 22" top post mount rigid arm	◀ Light Post (1300RM): 44" with dual slot, fixed sidebox mount, with Forest light bushing . 1110-170 \$ 220 35" with slot & post clamp, fixed mount, with Forest light bushing 1110-045 290 Pivot light posts & requirements ..... (call to order)			
Monitor Requirements: 12 volt with power supply Weight: 8 1/2 to 13 lbs max, width of 19"	◀ MONITOR POST (1300TPM): 30" with dual slot, fixed sidebox mount, with Oilite bushing 1331-998 220 30" with slot & post clamp, fixed mount, with Oilite bushing 1331-997 290 Pivot monitor posts & requirements ..... (call to order)			
			<b>TOTAL</b>	_____
	◀ Often overlooked			



Equipment dimensions, see page 100

Shipping weight estimate:

1300RM, 1300TPM 15 lb / 7 kg

## Designer Friendly® Surface Treatments

## STEP 1 SELECT POWDER COAT PAINT FOR UNIT &amp; LIGHT:

- |                                     |                                |                                |                              |
|-------------------------------------|--------------------------------|--------------------------------|------------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze   | <input type="radio"/> Charcoal | <input type="radio"/> Copper |
| <input type="radio"/> Forest Grey   | <input type="radio"/> Midnight | <input type="radio"/> Pewter   | <input type="radio"/> Shadow |
| <input type="radio"/> White         |                                |                                |                              |

## Notes



Rear Mount Side Delivery Systems



## Powder Coat



Forest Grey



Shadow



White



Midnight



Arctic Silver



Pewter



Copper



Charcoal



Bronze

## 7000 Series Rear and Side Delivery

Material: North America's only all metal control head chassis

Benefit: resistant to UV & cleaners

Material: hinged solid surface lid

Benefit: long-lasting surface gives easy access to ancillary products

Material: optional exclusive brass handpiece lockout valves

Benefit: allows assistant to remove a handpiece from its holder to change a burr while dentist continues working

Material: four steel mounting configurations

Benefit: wall, under cabinet, beneath counter & cart

Material: secondary aluminum work surface

Benefit: dedicated work surface for burr block, electric handpiece attachment, mirror, explorer

Material: steel height adjust mechanism

Benefit: adjustable work surface & handpiece delivery arm height

Material: nine Designer Friendly® powder coat paint colors including white & forest grey


Benefit: dental equipment with a WOW factor

(Shown on Front)

**7000** Powder Coat Paint Color: Charcoal - Solid Surface: Starry Starry Night

**Made for Each Other**

From the ground up, Forest delivery systems are designed for easy cabinet integration.

**3486 Beneath Counter Mount (BC)**

Start with a strong foundation

- Your choice of 3/4" hardened steel 11" x 14-1/2" or optional 10" x 24" mounting plate that bolts to the bottom of your counter top.
- Followed by a 6" x 6" leveling plate with friction adjustment.

**3486 Under Cabinet Floor Mount (UC)**

- Low profile 1/2" steel floor mount plate allows room for your stool casters.
- Needle roller bearings with friction adjustment eliminates drift.
- Durable 2-1/2" steel tube post and arm assembly with needle roller bearings installed between pivoting sections.

**7000 Series Beneath Counter Mount (BC)**


Start with a strong foundation

- Your choice of 3/4" hardened steel 11" x 14-1/2" or optional 10" x 24" mounting plate that bolts to the bottom of your counter top
- Followed by a 6" x 6" leveling plate and a set of needle roller bearings with friction adjustment.
- Durable 2-1/2" steel tube dual arm assembly with needle roller bearings installed between pivoting arm sections.

**7000 Series Under Cabinet Floor Mount (UC)**

Start with a strong foundation

- Low profile 1/2" steel floor mount plate allows room for your stool casters.
- Needle roller bearings with friction adjustment eliminates drift.
- Durable 2-1/2" steel tube post and arm assembly with needle roller bearings installed between pivoting sections.




# Rear Mount Side Delivery Systems

## 7020PRO Duo Swing Delivery System



7020PRO



**Model 7020PRO-BC \$ 8800**  
**Beneath Counter Mount**

Universal heavy-duty mounting plate 11" x 14 1/2" included (mounts to under counter)

Specify Cabinet Manufacturer

### Standard Features

#### Duo Delivery System 7020PRO Base Model

##### Duo control

- automatic control for three handpieces
- handpiece lockouts
- asepsis 3-way air water syringe w/tubing, 5'
- handpiece tubing, 4-hole Midwest, asepsis style, 5'
- pneumatic master control block w/diaphragm
- master on/off
- individual handpiece pressure adjustment
- individual water coolant adjustment
- non-retraction water coolant
- air coolant adjustment
- handpiece flush
- handpiece exhaust canister/oil collector
- wet/dry foot control

Self-contained clean water system w/quick disconnect,  
single bottle, 1 1/2 liter

"Burr Buddy" tray holder and covered Zirc tray

Chrome accents

Central vacuum accessory package

- central vacuum canister with disposable solids collector
- asepsis 3-way air/water syringe w/tubing, 5'
- autoclavable standard HVE lever valve w/tubing, 5'
- autoclavable saliva ejector lever valve w/tubing, 5'
- dual assistant's quad holder assembly

Solid work surface, 17" x 30"

Pneumatic work surface easy lift/close mechanism

Height adjustable work surface, 5"

Large and small organizer bins w/frame, dividers  
and sharps container

Procedure shelf

Keyboard tray

Stainless trash/CPU tray (MC)

Universal heavy duty mounting plate 11" x 14 1/2",  
beneath counter (BC)

Wall mount panel 16" x 18" (WM)

Standard utility center housing and cover (MC)

Air pressure regulator w/filter

Master air shut-off valve

Umbilical, 6'



**Beneath Counter Mount**



**Mobile Cart**



**Under Cabinet floor Mount**




**Wall Mount**


Mounting plate dimensions, see pages 101

Shipping weight estimate:


7020PRO BC, MC, UC, WM 100 lb / 46 kg



**Model 7020PRO-MC \$ 8800**  
Mobile Cart Mount



**Model 7020PRO-UC \$ 8800**  
Under Cabinet Floor Mount



**Model 7020PRO-WM \$ 8800**  
Wall Mount

Factory Options	Part No	Retail	Price
Designer Friendly® Surface Treatments			
<b>STEP 1 SELECT POWDER COAT PAINT:</b>			
<input type="radio"/> Arctic Silver <input type="radio"/> Bronze <input type="radio"/> Charcoal			
<input type="radio"/> Copper <input type="radio"/> Forest Grey <input type="radio"/> Midnight			
<input type="radio"/> Pewter <input type="radio"/> Shadow <input type="radio"/> White			
<b>STEP 2 SOLID WORK SURFACE COLORS:</b>			
Corian® <input type="radio"/> Antarctica			
Gibraltar®: <input type="radio"/> Arctic Melange <input type="radio"/> Baja Melange			
Meganite® <input type="radio"/> Cinder Granite <input type="radio"/> Goldrush Gemstone			
<input type="radio"/> Moonstone <input type="radio"/> Platinum Granite			
<input type="radio"/> Rocky Road Granite <input type="radio"/> Sanibel Granite			
<input type="radio"/> Silver Mist <input type="radio"/> Starry Starry Night			
<b>DELIVERY SYSTEM .....</b> <b>MODEL CHOICE.....</b>			<b>PRICE:</b>
Handpiece Illumination System			
One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS	\$ 715		
Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS	275		
One handpiece (6-pin), accommodates 1 add'l tubing . 1700-019-CS	715		
Add'l handpiece illumination tubing (6-pin) ..... 0017-603-CS	195		
<b>■ Ancillaries use existing handpiece position (tubing lengths may vary) (see page 57)</b>			
Write in selection here .....			
Write in selection here .....			
Additional handpiece position, wet .....	1100-189-CS	425	
Additional handpiece position, air drive only .....	1100-190-CS	215	
Additional standard autoclavable HVE			
with tubing, 5' & autoclavable oral dry cup .....	1400-113-CS	90	
Quick disconnect with flow control 1/4" (for water).... 0016-454-CS	145		
Quick disconnect with flow control 3/8" (for air) .... 0016-456-CS	145		
<b>■ Monitor mount assembly (monitor screen not included) 1300-025-CS</b>	450		
<b>■ Swivel tray .....</b> 1110-909-CS	500		
Large bin divider kit (white or black)....0010-860 / 0010-860-BLK	175		
Small bin divider kit (white or black)....0010-861 / 0010-861-BLK	95		
Sharps container kit (white or black)....0010-870 / 0010-870-BLK	95		
<b>■ Touchpad for 3900 chair (25 ft. Cable) .....</b> 1700-081-1-CS	515		
Standard utility center housing & cover (BC/UC/WM).... 1115-020	145		
In floor recessed utility center housing & cover..... 1115-305	280		
Additional umbilical (per foot)..... 1115-999-CS	30		
Routing valve kit with water shutoff &			
regulator/factory installation .....	1115-315-CS	305	
Universal mounting plate, 10" x 24"			
in lieu of 11" x 141/2" (BC)..... 0010-084-CS	130		
Special color solid work surface .....	0023-856-SSV		
(call for price and leadtime)			
<b>TOTAL</b>			

<b>Designer Friendly® Surface Treatments</b>			
<b>STEP 1 SELECT POWDER COAT PAINT:</b>			
<input type="radio"/> Arctic Silver <input type="radio"/> Bronze <input type="radio"/> Charcoal			
<input type="radio"/> Copper <input type="radio"/> Forest Grey <input type="radio"/> Midnight			
<input type="radio"/> Pewter <input type="radio"/> Shadow <input type="radio"/> White			
<b>STEP 2 SOLID WORK SURFACE COLORS:</b>			
Corian® <input type="radio"/> Antarctica			
Gibraltar®: <input type="radio"/> Arctic Melange <input type="radio"/> Baja Melange			
Meganite® <input type="radio"/> Cinder Granite <input type="radio"/> Goldrush Gemstone			
<input type="radio"/> Moonstone <input type="radio"/> Platinum Granite			
<input type="radio"/> Rocky Road Granite <input type="radio"/> Sanibel Granite			
<input type="radio"/> Silver Mist <input type="radio"/> Starry Starry Night			
<b>DELIVERY SYSTEM .....</b> <b>MODEL CHOICE.....</b>			<b>PRICE:</b>
Handpiece Illumination System			
One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS	\$ 715		
Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS	275		
One handpiece (6-pin), accommodates 1 add'l tubing . 1700-019-CS	715		
Add'l handpiece illumination tubing (6-pin) ..... 0017-603-CS	195		
<b>■ Ancillaries use existing handpiece position (tubing lengths may vary) (see page 57)</b>			
Write in selection here .....			
Write in selection here .....			
Additional handpiece position, wet .....	1100-189-CS	425	
Additional handpiece position, air drive only .....	1100-190-CS	215	
Additional standard autoclavable HVE			
with tubing, 5' & autoclavable oral dry cup .....	1400-113-CS	90	
Quick disconnect with flow control 1/4" (for water).... 0016-454-CS	145		
Quick disconnect with flow control 3/8" (for air) .... 0016-456-CS	145		
<b>■ Monitor mount assembly (monitor screen not included) 1300-025-CS</b>	450		
<b>■ Swivel tray .....</b> 1110-909-CS	500		
Large bin divider kit (white or black)....0010-860 / 0010-860-BLK	175		
Small bin divider kit (white or black)....0010-861 / 0010-861-BLK	95		
Sharps container kit (white or black)....0010-870 / 0010-870-BLK	95		
<b>■ Touchpad for 3900 chair (25 ft. Cable) .....</b> 1700-081-1-CS	515		
Standard utility center housing & cover (BC/UC/WM).... 1115-020	145		
In floor recessed utility center housing & cover..... 1115-305	280		
Additional umbilical (per foot)..... 1115-999-CS	30		
Routing valve kit with water shutoff &			
regulator/factory installation .....	1115-315-CS	305	
Universal mounting plate, 10" x 24"			
in lieu of 11" x 141/2" (BC)..... 0010-084-CS	130		
Special color solid work surface .....	0023-856-SSV		
(call for price and leadtime)			
<b>TOTAL</b>			

**■ Often overlooked**


Cabinet mounting configurations, see page 41

# Rear Mount Side Delivery Systems

## 7021PRO Duo Swing Delivery System



**7021PRO**



**Model 7021PRO-BC \$ 7575**  
**Beneath Counter Mount**

Universal heavy-duty mounting plate 11" x 14 1/2" included (mounts to under counter)

Specify Cabinet Manufacturer

### Standard Features

#### Duo Delivery System 7021PRO

Duo control

- automatic control for three handpieces
- handpiece lockouts
- asepsis 3-way air water syringe w/tubing, 5'
- handpiece tubing, 4-hole Midwest, asepsis style, 5'
- pneumatic master control block w/diaphragm
- master on/off
- individual handpiece pressure adjustment
- individual water coolant adjustment
- non-retraction water coolant
- air coolant adjustment
- handpiece flush
- handpiece exhaust canister/oil collector
- wet/dry foot control
- Self-contained clean water system w/quick disconnect, single bottle, 1 1/2 liter
- "Burr Buddy" tray holder and covered Zirc tray
- Chrome accents
- Central vacuum accessory package
  - central vacuum canister with disposable solids collector
  - asepsis 3-way air/water syringe w/tubing, 5'
  - autoclavable standard HVE lever valve w/tubing, 5'
  - autoclavable saliva ejector lever valve w/tubing, 5'
  - dual assistant's quad holder assembly
- Solid work surface, 17" x 30"
- Pneumatic work surface easy lift/close mechanism
- Height adjustable work surface, 5"
- Large and small organizer bins w/frame, dividers and sharps container
- Universal heavy duty mounting plate 11" x 14 1/2", beneath counter (BC)
- Wall mount panel 16" x 18" (WM)
- Standard utility center housing and cover (MC)
- Air pressure regulator w/filter
- Master air shut-off valve
- Umbilical, 6'



**Beneath Counter Mount**



**Mobile Cart**



**Under Cabinet floor Mount**



**Wall Mount**


Mounting plate dimensions, see pages 101

Shipping weight estimate:


7021PRO BC, MC, UC, WM 100 lb / 46 kg

## Rear Mount Side Delivery Systems


## 7021PRO Duo Swing Delivery System



## **Model 7021PRO-MC \$ 7575 Mobile Cart Mount**



## **Model 7021PRO-UC \$ 7575 Under Cabinet Floor Mount**



**Model 7021PRO-WM \$ 7575**  
**Wall Mount**

Designer Friendly® Surface Treatments

## **STEP 1 SELECT POWDER COAT PAINT:**

- Arctic Silver
  - Bronze
  - Charcoal
  - Copper
  - Forest Grey
  - Midnight
  - Pewter
  - Shadow
  - White

## **STEP 2 SOLID WORK SURFACE COLORS:**

- Corian®  Antarctica

Gibraltar®: Arctic Melange

- Baja Melange

Meganite®

- Cinder Granite
  - Goldrush Gemstone
  - Moonstone
  - Platinum Granite
  - Rocky Road Granite
  - Sanibel Granite
  - Silver Mist
  - Starry Starry Night

Factory Options	Part No	Retail	Price
<b>DELIVERY SYSTEM . . . . . MODEL CHOICE _____</b>		<b>PRICE :</b>	
Handpiece Illumination System			
One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS	\$715		
Add'l handpiece illumination tubing (5-hole) . . . . . 0017-615-CS	275		
One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS	715		
Add'l handpiece illumination tubing (6-pin) . . . . . 0017-603-CS	195		
<b> Ancillaries use existing handpiece position (tubing lengths may vary) (see page 57)</b>			
Write in selection here . . . . .			
Write in selection here . . . . .			
Additional handpiece position, wet . . . . . 1100-189-CS	425		
Additional handpiece position, air drive only. . . . . 1100-190-CS	215		
Additional standard autoclavable HVE			
with tubing, 5' & autoclavable oral dry cup . . . . . 1400-113-CS	90		
Quick disconnect with flow control 1/4" (for water) . . 0016-454-CS	145		
Quick disconnect with flow control 3/8" (for air) . . . . 0016-456-CS	145		
<b> Procedure shelf. . . . . 0010-864-CS</b>	775		
<b> Keyboard tray . . . . . 0010-862-CS</b>	475		
<b> Stainless trash/CPU tray, H-frame (MC) . . . . . 0010-863</b>	150		
<b> Monitor mount assembly (monitor screen not included)</b>			
base model requires purchase of . .0010-864-CS/1300-025-CS	450		
<b> Swivel tray. . . . . 1110-909-CS</b>	500		
Large bin divider kit (white or black) . . . . . 0010-860 / 0010-860-BLK	175		
Small bin divider kit (white or black) . . . . . 0010-861 / 0010-861-BLK	95		
Sharps container kit (white or black) . . . . . 0010-870 / 0010-870-BLK	95		
<b> Touchpad for 3900 chair (25 ft. Cable) . . . . . 1700-081-1-CS</b>	515		
Standard utility center housing & cover (BC/UC/WM) . . 1115-020	145		
In floor recessed utility center housing & cover . . . . 1115-305	280		
Additional umbilical (per foot). . . . . 1115-999-CS	30		
Routing valve kit with water shutoff &			
regulator/factory installation. . . . . 1115-315-CS	305		
Universal mounting plate, 10" x 24" in lieu of 11" x 14 1/2" (BC)			
. . . . . 0010-084-CS	130		
Special color solid work surface . . . . . 0023-856-SSV			
(call for price and leadtime)			


**Cabinet mounting configurations, see page 41**

### Often overlooked

**TOTAL**

# Rear Mount Side Delivery Systems

## 7022PRO Assistant's Delivery System



**Model 7022PRO-BC \$ 7300**  
**Beneath Counter Mount**

Universal Heavy-Duty Mounting Plate 11" x 14 1/2" included (mounts to under counter)  
Specify Cabinet Manufacturer

### Standard Features

#### **Assistant's Delivery System 7022PRO**

Central vacuum accessory package  
central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'  
dual assistant's quad holder assembly  
Self-contained clean water system w/quick disconnect,  
single bottle, 1 1/2 liter  
Solid work surface, 17" x 30"  
Pneumatic work surface easy lift/close mechanism  
Height adjustable work surface, 5"  
Large and small organizer bins w/frame, dividers  
and sharps container  
Procedure shelf  
Keyboard tray  
Stainless trash/CPU tray (MC)  
Universal heavy duty mounting plate 11" x 14 1/2",  
beneath counter (BC)  
Wall mount panel 16" x 18" (WM)  
Standard utility center housing and cover (MC)  
Air pressure regulator w/filter  
Master air shut-off valve  
Umbilical, 6'



**Beneath Counter Mount**



**Mobile Cart**



**Under Cabinet floor Mount**




**Wall Mount**

**Mounting plate dimensions, see pages 101**

Shipping weight estimate:

7022PRO BC, MC, UC, WM 90 lb / 42 kg



**Model 7022PRO-MC \$ 7300**  
Mobile cart mount

**Model 7022PRO-UC \$ 7300**  
Under Cabinet Floor Mount

**Model 7022PRO-WM \$ 7300**  
Wall Mount



### Designer Friendly® Surface Treatments

#### STEP 1 SELECT POWDER COAT PAINT:

- |                                     |                                   |                                |
|-------------------------------------|-----------------------------------|--------------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze      | <input type="radio"/> Charcoal |
| <input type="radio"/> Copper        | <input type="radio"/> Forest Grey | <input type="radio"/> Midnight |
| <input type="radio"/> Pewter        | <input type="radio"/> Shadow      | <input type="radio"/> White    |

#### STEP 2 SOLID WORK SURFACE COLORS:

Corian®  Antarctica

Gibraltar®:  Arctic Melange  
 Baja Melange

Meganite®  Cinder Granite  
 Goldrush Gemstone  
 Moonstone  
 Platinum Granite  
 Rocky Road Granite  
 Sanibel Granite  
 Silver Mist  
 Starry Starry Night

### Factory Options

### Part No      Retail      Price

DELIVERY SYSTEM .....	MODEL CHOICE.....	PRICE : .....
Additional standard autoclavable HVE		
with tubing, 5' & autoclavable oral dry cup .....	1400-113-CS	\$ 90
Quick disconnect with flow control 1/4" (for water)....	0016-454-CS	145
Quick disconnect with flow control 3/8" (for air) .....	0016-456-CS	145
Monitor mount assembly (monitor screen not included)	1300-025-CS	450
Swivel tray .....	1110-909-CS	500
Large bin divider kit (white or black).....	0010-860 / 0010-860-BLK	175
Small bin divider kit (white or black).....	0010-861 / 0010-861-BLK	95
Sharps container kit (white or black).....	0010-870 / 0010-870-BLK	95
Touchpad for 3900 chair (25 ft. Cable) .....	1700-081-1-CS	515
Standard utility center housing & cover (BC/UC/WM)....	1115-020	145
In floor recessed utility center housing & cover.....	1115-305	280
Additional umbilical (per foot) .....	1115-999-CS	30
Routing valve kit with water shutoff & regulator/factory		
installation .....	1115-315-CS	305
Universal mounting plate, 10" x 24" in lieu of 11" x 14 1/2" (BC) ...	0010-084-CS	130
Special color solid work surface .....	0023-856-SSV	
(call for price and leadtime)		
<b>TOTAL</b>		


► Often overlooked

# Rear Mount Side Delivery Systems

## 7023PRO Assistant's Delivery System



7023PRO



**Model 7023PRO-BC \$ 6075**  
**Beneath Counter Mount**

Universal heavy-duty mounting plate 11" x 14 1/2" included (mounts to under counter)  
Specify Cabinet Manufacturer

### Standard Features

#### Assistant's delivery system 7023PRO

Central vacuum accessory package

central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'  
dual assistant's quad holder assembly

Self-contained clean water system w/quick disconnect,  
single bottle, 1 1/2 liter

Solid work surface, 17" x 30"

Pneumatic work surface easy lift/close mechanism

Height adjustable work surface, 5"

Large and small organizer bins w/frame,  
dividers and sharps container

Universal heavy duty mounting plate 11" x 14 1/2",  
beneath counter (BC)

Wall mount panel 16" x 18" (WM)

Standard utility center housing and cover (MC)

Air pressure regulator with filter

Master air shut-off valve

Umbilical, 6'



**Beneath Counter Mount**



**Mobile Cart**



**Under Cabinet floor Mount**




**Wall Mount**

**Mounting plate dimensions, see pages 101**

Shipping weight estimate:

7023PRO BC, MC, UC, WM 90 lb / 42 kg



**Model 7023PRO-MC \$ 6075**  
Mobile Cart Mount

**Model 7023PRO-UC \$ 6075**  
Under Cabinet Floor Mount

**Model 7023PRO-WM \$ 6075**  
Wall Mount

Factory Options	Part No	Retail	Price
<b>DELIVERY SYSTEM .....</b> <b>MODEL CHOICE</b> _____ <b>PRICE:</b> _____			
Additional standard autoclavable HVE			
with tubing, 5' & autoclavable oral dry cup .....	1400-113-CS	\$ 90	_____
Quick disconnect with flow control 1/4" (for water)....	0016-454-CS	145	_____
Quick disconnect with flow control 3/8" (for air) .....	0016-456-CS	145	_____
<b>Procedure shelf.</b> ....	0010-864-CS	775	_____
<b>Keyboard tray</b> .....	0010-862-CS	475	_____
<b>Stainless trash/CPU tray, H-frame (MC)</b> .....	0010-863	150	_____
<b>Monitor mount assembly (monitor screen not included)</b> base model requires purchase of 0010-864-CS .....	1300-025-CS	450	_____
<b>Swivel tray</b> .....	1110-909-CS	500	_____
Large bin divider kit (white or black).....	0010-860 / 0010-860-BLK	175	_____
Small bin divider kit (white or black).....	0010-861 / 0010-861-BLK	95	_____
Sharps container kit (white or black).....	0010-870 / 0010-870-BLK	95	_____
<b>Touchpad for 3900 chair (25 ft. Cable)</b> .....	1700-081-1-CS	515	_____
<b>Standard utility center housing &amp; cover (BC/UC/WM)</b> ....	1115-020	145	_____
In floor recessed utility center housing & cover .....	1115-305	280	_____
Additional umbilical (per foot) .....	1115-999-CS	30	_____
Routing valve kit with water shutoff & regulator/factory installation .....	1115-315-CS	305	_____
Universal mounting plate, 10" x 24" in lieu of 11" x 14 1/2" (BC) .....	0010-084-CS	130	_____
Special color solid work surface .....	0023-856-SSV (call for price and leadtime)		_____
<b>TOTAL</b> _____			

**Often overlooked**

# Rear Mount Side Delivery Systems

## 7010 Duo Swing Delivery System




### 1110-185-CS and 1110-908-CS

Shown with optional storage organizer and swivel tray



Shown with optional handpiece lockouts 0010-748  
and additional assistant's quad holder assembly 0024-156



### Model 7010BC \$ 6630 Beneath Counter Mount

Universal heavy-duty mounting plate 11" x 14 1/2" included (mounts to under counter)  
Specify Cabinet Manufacturer

## Standard Features

### Duo Swing Delivery System 7010

#### Duo control

- automatic control for three handpieces
- asepsis 3-way air water syringe w/tubing, 4'
- handpiece tubing, 4-hole Midwest, asepsis style, 4'
- pneumatic master control block w/diaphragm
- master on/off
- individual handpiece pressure adjustment
- individual water coolant adjustment
- non-retraction water coolant
- air coolant adjustment
- handpiece flush
- handpiece exhaust canister/oil collector
- self-contained clean water system w/quick disconnect, single bottle, 1 1/2 liter
- wet/dry foot control

#### Height adjustable work surface, 5" (BC/UC/WM)

#### "Burr Buddy" tray holder and covered Zirc tray

#### Pneumatic work surface easy lift/close mechanism

#### Central vacuum accessory package

- central vacuum canister w/disposable
- solids collector
- asepsis 3-way air/water syringe w/tubing, 5'
- autoclavable standard HVE lever valve w/tubing, 5'
- autoclavable saliva ejector lever valve w/tubing, 5'

#### Solid work surface, 19" x 24"

#### Universal heavy duty mounting plate

11" x 14 1/2", beneath counter (BC)

#### Wall mount panel 16" x 18" (WM)

#### Adjustable height frame (MC)

#### Twin wheel hard floor casters (MC)

#### Standard utility center housing and cover (MC)

#### Air pressure regulator w/filter

#### Master air shut-off valve

#### Umbilical, 7'

#### External umbilical, 6' (MC)



Beneath Counter Mount



Mobile Cart



Under Cabinet floor Mount




Wall Mount

Mounting plate dimensions, see pages 105


Shipping weight estimate:

7010BC, MC, UC, WM 130 lb / 59 kg


## 7010 Duo Swing Delivery System



**Model 7010MC \$ 5990**  
Mobile Cart Mount



**Model 7010UC \$ 6630**  
Under Cabinet Floor Mount



**Model 7010WM \$ 6630**  
Wall Mount



### Designer Friendly® Surface Treatments

#### STEP 1 SELECT POWDER COAT PAINT:

- |                                     |                                   |                                |
|-------------------------------------|-----------------------------------|--------------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze      | <input type="radio"/> Charcoal |
| <input type="radio"/> Copper        | <input type="radio"/> Forest Grey | <input type="radio"/> Midnight |
| <input type="radio"/> Pewter        | <input type="radio"/> Shadow      | <input type="radio"/> White    |

#### STEP 2 SOLID WORK SURFACE COLORS:

Corian®  Antarctica

Gibraltar®:  Arctic Melange  
 Baja Melange

Meganite®  Cinder Granite  
 Goldrush Gemstone  
 Moonstone  
 Platinum Granite  
 Rocky Road Granite  
 Sanibel Granite  
 Silver Mist  
 Starry Starry Night

#### Factory Options

#### Part No

#### Retail

#### Price

#### DELIVERY SYSTEM ..... MODEL CHOICE ..... PRICE:

##### Handpiece Illumination System

- |   |        |
|---|--------|
| One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS | \$ 715 |
| Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS   | 275    |
| One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS  | 715    |
| Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS    | 195    |

► Ancillaries use existing handpiece position (tubing lengths may vary) (see page 57)

Write in selection here .....

Write in selection here .....

- |   |             |     |
|---|-------------|-----|
| Additional handpiece position, wet .....            | 1100-189-CS | 425 |
| Additional handpiece position, air drive only ..... | 1100-190-CS | 215 |
| Lockout feature for 3 handpieces .....              | 0010-748    | 215 |

Additional standard autoclavable HVE  
with tubing, 5' & autoclavable oral dry cup .....

1400-113-CS 90

Additional assistant's quad holder assembly .....

0024-156 205

Quick disconnect with flow control 1/4" (for water) ...

0016-454-CS 145

Quick disconnect with flow control 3/8" (for air) .....

0016-456-CS 145

► Storage organizer with bins & frame..... 1110-185-CS

625

► Swivel tray (requires purchase of 1110-185-CS)..... 1110-908-CS

350

► Swivel tray (mounts without storage organizer)..... 1110-079-CS

475

Monitor mount assembly (monitor screen not included)

(requires purchase of 1110-185-CS) .....

1300-025-CS 450

Touchpad for 3900 chair (25 ft. Cable) .....

1700-081-1-CS 515

Chrome accents .....

9000-998-CS 140

Standard utility center housing & cover (BC/UC/WM)....

1115-020 145

In floor recessed utility center housing & cover.....

1115-305 280

Additional umbilical (per foot) .....

1115-999-CS 30

Routing valve kit with water shutoff & regulator/factory installation .....

1115-315-CS 305

Universal mounting plate, 10" x 24" in lieu of 11" x 14 1/2" (BC) .....

0010-084-CS 130

Special color solid work surface .....

0023-776-SSV

(call for price and leadtime)

**TOTAL**

# Rear Mount Side Delivery Systems

## 3486 Swing Delivery System



**Model 3486**

Shown with optional Storage Organizer and Swivel Tray  
1110-179-CS and 1110-908-CS

### Standard Features

#### **Assistant's Swing Delivery System 3486 Series**

Central vacuum accessory package

arm-mounted central vacuum canister

w/disposable solids collector

asepsis 3-way air/water syringe w/tubing, 5'

autoclavable standard HVE lever valve w/tubing, 5'

autoclavable saliva ejector lever valve w/tubing, 5'

self-contained clean water system w/quick disconnect,  
single bottle, 1 1/2 liter

Assistant's quad holder assembly

Solid work surface, 19" x 24"

Beneath counter rigid arm assembly

w/heavy duty mounting plate, 11" x 14 1/2" (BC)

Under cabinet floor mount delivery post (UC)

Wall mount panel, 16" x 18" (WM)

Air pressure regulator w/filter

Master shut-off valve

Umbilical, 7'



**Beneath Counter Mount**






**Under Cabinet floor Mount**

**Wall Mount**

Mounting plate dimensions, see pages 107-108

Shipping weight estimate:

3486 115 lb / 53 kg



**Model 3486BC \$ 3225**  
Beneath Counter Mount  
Specify Cabinet Manufacturer

**Model 3486UC \$ 3690**  
Under Cabinet Floor Mount

**Model 3486WM \$ 3225**  
Wall Mount

Factory Options	Part No	Retail	Price
DELIVERY SYSTEM ..... MODEL CHOICE _____			PRICE : _____
Additional standard autoclavable HVE			
with tubing, 5' & autoclavable oral dry cup.....	1400-113-CS	\$ 90	_____
Additional assistant's quad holder assembly .....	0024-149	205	_____
Quick disconnect with flow control 1/4" (for water) .....	0016-454-CS	145	_____
Quick disconnect with flow control 3/8" (for air) .....	0016-456-CS	145	_____
Storage organizer with bins & frame.....	1110-179-CS	525	_____
Swivel tray (requires purchase of 1110-179-CS).....	1110-908-CS	350	_____
Swivel tray (mounts without storage organizer).....	1110-482-CS	425	_____
Monitor mount assembly (monitor screen not included)			
(requires purchase of 1110-179-CS) .....	1300-025-CS	450	_____
Touchpad for 3900 chair (25 ft. Cable) .....	1700-081-1-CS	515	_____
Standard utility center housing & cover .....	1115-020	145	_____
In floor recessed utility center housing & cover.....	1115-305	280	_____
Additional umbilical (per foot) .....	1115-999-CS	30	_____
Routing valve kit with water shutoff & regulator/factory			
installation.....	1115-315-CS	305	_____
Standard color laminate work surface in lieu of solid surface 0023-003-LAM		-325	_____
Special color laminate or solid work surface upgrade .....	0023-003-X/0023-003-SSV		_____
(call for price and leadtime)			
<b>TOTAL</b> _____			

Additional standard autoclavable HVE

with tubing, 5' & autoclavable oral dry cup..... 1400-113-CS \$ 90

Additional assistant's quad holder assembly .....

0024-149 205

Quick disconnect with flow control 1/4" (for water) .....

0016-454-CS 145

Quick disconnect with flow control 3/8" (for air) .....

0016-456-CS 145

Storage organizer with bins & frame.....

1110-179-CS 525

Swivel tray (requires purchase of 1110-179-CS).....

1110-908-CS 350

Swivel tray (mounts without storage organizer).....

1110-482-CS 425

Monitor mount assembly (monitor screen not included)

(requires purchase of 1110-179-CS) .....

1300-025-CS 450

Touchpad for 3900 chair (25 ft. Cable) .....

1700-081-1-CS 515

Standard utility center housing & cover .....

1115-020 145

In floor recessed utility center housing & cover.....

1115-305 280

Additional umbilical (per foot) .....

1115-999-CS 30

Routing valve kit with water shutoff & regulator/factory

    installation.....

1115-315-CS 305

Standard color laminate work surface in lieu of solid surface 0023-003-LAM

Special color laminate or solid work surface upgrade .....

0023-003-X/0023-003-SSV


(call for price and leadtime)

Often overlooked

# Rear Mount Side Delivery Systems

## Storage Organizer, Trays and Monitor Mount

### Storage Organizer with Bins and Frame



Factory option model 7010                            1110-185-CS                            \$ 625  
Factory option models: 2184BC/2186BC/3486/5386/5886                    1110-179-CS                            525

- Metal powder coated organizer frame
- Three storage bins with removable covers
- Storage dividers

### Swivel Tray



Factory option models:

1110-908-CS                            \$ 350

2184BC/2186BC/3486/5386/5886/7010 requires purchase of 1110-179-CS or 1110-185-CS

- Inserts inside storage organizer tower
- Metal powder coated tray holder
- Plastic instrument tray

### Swivel Tray with Left or Right Mounting Bracket



Factory option model 7010                            1110-079-CS                            \$ 475  
Factory option models: 2184BC/2186BC/3486/5386/5886                    1110-482-CS                            425

- Mounts without storage organizer
- Metal powder coated tray holder with tower and bracket
- Plastic instrument tray

### Swivel Tray with Left and Right Mounting Brackets



Factory option models:                                    1110-909-CS                            \$ 500

- Metal powder coated tray holder, towers and brackets
- Plastic instrument tray

PRO Series 2184BC/286BC/3486 mounts without storage organizer 2184BC/2186BC/3486 and 7010 third tray addition mounts with organizer

### Monitor Mount



Factory option models:                                    1300-025-CS                            \$ 450

- Inserts inside storage organizer or PRO Series procedure shelf tower
- Forest's unique locking tilt mechanism with knob 1/2" pin


PRO Series 7021 and 7023 requires purchase of 0010-864-CS/2184BC/2186BC/3486/5386/5886 and 7010 requires purchase of 1110-179-CS or 1110-185-CS

#### Monitor requirements (not included):

12 volt with power supply

Weight 8-1/4 to 13 lbs. max, width of 19"

## Delivery Units and Ancillary Products



## Carbon Fiber




Laser Leather




Lemans



Greywood



White



Forest Grey

## Plastics

### Standard IC

Material: aluminum hinged lid / 1/8" ABS plastic chassis with durable DR acrylic protective coating

Benefit: hinge allows quick access to controls & ancillary devices / lid & chassis are resistant to UV, impact & compatible surface cleaners / aluminum lid is strong & lightweight

Material: electro-less nickel plated brass air & water Ts, routing valves, pneumatic control & distribution blocks

Benefit: lifetime warranty on all parts "under the hood" (on all forest control heads) / trouble-free handpiece operation

Material: aluminum handpiece holder bar with glass-filled nylon injection molded handpiece holders

Benefit: allows for up to four handpiece positions, a syringe & ancillary touchpad at your desired angle & spacing

Material: hinged aluminum control head lid

Benefit: available in nine two-coat powder coat paint colors

Material: plastic chassis

Benefit: available in three carbon fiber look & two solid color finishes

### Euro

Material: whip arms rotate on nylon wheels with aluminum & brass components

Benefit: feather-touch whip movement / ultralight handpiece feel / long, trouble free life

Material: control head swivels at its center point

Benefit: 240-degrees of control head rotation allows for beside-or over-the-patient placement

Material: forward mounted aluminum whip arm system

Benefit: allows for greater reach of handpiece & syringe expanding control head placement options.

Material: plastic control head lid

Benefit: available in three carbon fiber looks & two solid color finishes

Feature: autoclavable silicone handpiece pad available in black & white

Benefit: blends with Designer Friendly® choices

(Shown on Front)

**Forest Euro** Powder Coat Paint Color: Charcoal - Carbon Fiber: Greywood


**Forest Standard IC** Powder Coat Paint Color: Bronze - Carbon Fiber: Laser Leather

**Standard IC**

Please use graph below to determine number of ancillaries and options possible for the Standard IC.

**Standard Features**

Automatic control for three handpieces  
 asepsis 3-way air/water syringe with tubing, 5' handpiece tubing, 4-hole Midwest, asepsis style, 5' pneumatic master control block with diaphragm master on/off  
 individual handpiece pressure adjustment  
 individual water coolant adjustment  
 non-retraction water coolant  
 air/water coolant adjustment  
 handpiece flush  
 handpiece exhaust canister/oil collector  
 handle brake release function

**% Standard IC Control Head Space Absorbed per Option****Please note:**

If your ancillary and option choices exceed 100%, please call customer service for additional information.

\* Uses existing handpiece position

★ Standard on SE and SI deliveries

● Standard on 1500PI delivery

+ 3900 Chair touchpad cable standard on PI and SI deliveries

◆ Two or more require purchase of 0010-310 ancillary shelf (see page 61)

# Ancillary Ordering and Options

## Euro Delivery Unit

### Euro


Please use graph below to determine number of ancillaries and options possible for the Euro.



### Standard Features

Automatic control for four handpieces  
asepsis 3-way air/water syringe with tubing, 3' handpiece tubing, 4-hole Midwest, asepsis style, 3' pneumatic master control block with diaphragm master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air/water coolant adjustment  
handpiece flush  
handpiece exhaust canister/oil collector  
dual handle brake release function  
swivel instrument tray

% Euro Control Head Space Absorbed per Option



#### Please note:

If your ancillary and option choices exceed 100%, please call customer service for additional information.

\* Uses existing handpiece position

★ Standard on SE and SI deliveries

+ 3900 Chair touchpad cable standard on PI and SI deliveries

**Satelec SP**

Factory option	SATELEC-S	\$ 2310
Factory option	SATELEC-S-LED	2900



- Please call 800-289-63367 for warranty and technical support of the Scaler
- Uses existing handpiece position
- Satelec SP Newtron LED lit Piezo Ultrasonic Scaler w/power supply
- Color coding system ensures proper tip and power setting
- Three tips: #1, #10Z, H3 (Intro Scaling Kit)
- Tips available for periodontics, endodontics, scaling, restorative, and cavity preparation
- Cruise Control® automatically adjusts frequency amplitude
- On/off controlled by rheostat
- First ever autoclavable LED lit handpiece – 100,000 LUX
- Newtron technology means power, control and reliability
- Satelec warranty 1 year / tips 30 days from date of purchase

**Scaler**

Factory option	CAVITRON	\$ 4065
----------------	----------	---------



- Please call 800-800-7202 for warranty and technical support of the Scaler
- Uses existing handpiece position Dentsply Cavitron Magnetostrictive Ultrasonic scaler
- Blue Zone® for patient comfort
- SPS® (Sustained Performance System)  
constant balance between scaling efficiency and patient comfort
- Handpiece w/330 degree swivel and lavage control (reduces cord drag)
- Autoclavable Steri-mate® (infection control)
- PowerLINE 1000 30K insert
- Frequency 30 kHz
- Cavitron warranty 2 years from date of purchase

# Control Head Accessories

## Ancillary Products

### HP Illumination System, 5 Hole Adapter with Krypton Lamp

#### HP Illumination System, 6 Pin



Factory option

1700-021-CS \$ 715

Factory option

1700-019-CS 715

- Uses existing handpiece position
- Dual channel outputs individually adjustable from 2.5VDC to 4.0VDC
- Accommodates one additional handpiece tubing
- Illumination delay (adjustable between 0-15 seconds) can be set individually for 2 HP
- Automatic activation of HP illumination w/the dental unit foot control
- Simple installation into any new or existing dental unit
- Indicator LEDs confirms correct operation of the control module
- Powers all brands of optic handpieces and air-motors
- ISO standard handpiece terminal with stainless steel nut
- Asepsis style handpiece tubing, (4' to 5' installed)
- Power transformer, output 24VAC, 1.05 Amps
- Installation hardware and instructions

### HP Illumination Tubing, 5 Hole Adapter with Krypton Lamp



Factory option

0017-615-CS \$ 275

Used to replace power optic handpiece tubing  
ISO standard power optic handpiece terminal with stainless steel nut  
5 hole adapter with krypton lamp  
Asepsis style handpiece tubing, (4' to 5' factory installed)

### HP Illumination Tubing, 6 Pin



Factory option

0017-603-CS \$ 195

- Used to replace power optic handpiece tubing
- ISO standard power optic handpiece terminal with stainless steel adapter and nut
- Asepsis style handpiece tubing, (4' to 5' factory installed)

**Chair Control, Touchpad, Unit-Mounted**

Factory option

3914-998-CS

N/C

- 4 programmable preset positions

**Additional Chair Control, Touchpad, Unit-Mounted**

Factory option

3914-999-CS

\$ 515

**Assistant's Hygiene Instrumentation**

Factory option

1408-CS

\$ 385

- Canister with 1/2" drain
- Mounting bracket
- Drain tubing
- Autoclavable standard HVE w/tubing, 5'
- Autoclavable saliva ejector w/tubing, 5'
- Standard vacuum holder
- Universal holder
- Installed on Standard IC and Value control heads (not available on Euro)

**Ancillary Shelf**

Factory option

0010-310

\$135

- Accommodates additional ancillary products
- Designed for easy access to installed devices
- Powder coated aluminum
- Installed on Standard IC control head only

## Notes

Coming Soon



## THE FUSION PACKAGE

Combining elements of strength, design and affordability.

JUNE 2019




## Powder Coat



Greige

## NaugaHyde®



Black




Cocoa



Deep Sapphire



Saddle



Sand



Taupe

(Shown on Front)

The Fusion Package Powder Coat Paint Color: Greige  
Upholstery NaugaHyde®: Black

## Patient Chair

Material: hydraulic cantilever mechanism

Benefit: allows for smooth and steady heavy lifting

Material: cast aluminum chair back, base, armrest and headrest

Benefit: strong, impervious to fading or cracking

Material: flip down armrests

Benefit: allows for easy patient entry and exit

## LED Light

Material: hands free motion activation

Benefit: easy on/off while maintaining asepsis

Material: cast aluminum light head

Benefit: allows for better aesthetics and durability

## Delivery System

Material: ABS plastic control cover w/tough DR acrylic protective coating

Benefit: resistant to UV, impact & compatible surface cleaners/easy removal for quick access to controls & ancillary devices

Material: aluminum chassis base plate

Benefit: strong support for handpiece holder bar, positioning handles & control lid

Material: chrome-plated or powder coated handpiece holder bar

Benefit: rich look enhancing a value-based control head

# The Fusion Package

## Pivot Package

### Popular Equipment Choices



Cuspidor 1590

Ceiling Light



Shipping weight estimate:

Delivery Unit & Chair 445 lb / 203 kg

Base package 445 lb / 203 kg

As shown package 585 lbs / 266 kg

## Standard Features

## Factory Options

## Part No

## Retail Price

**Delivery System 4295PI**

Automatic control for three handpieces  
 asepsis 3-way air/water syringe w/tubing, 5'  
 handpiece tubing, 4-hole Midwest, asepsis style, 5'  
 pneumatic master control block w/diaphragm  
 master on/off  
 individual handpiece pressure adjustment  
 individual water coolant adjustment  
 non-retraction water coolant  
 air coolant adjustment  
 handpiece flush  
 handpiece exhaust canister/oil collector  
 dual handle brake release function  
 wet/dry foot control  
 self-contained clean water system w/quick disconnect,  
 single bottle, 1 1/2 liter w/on/off  
 toggle valve for city/bottle water routing  
 Flexible control arm w/pneumatic brake  
 vertical travel 28" / horizontal travel 62"  
 Internal umbilical, 5' w/3' external/chair to utility center  
 Utility center air/water pressure regulators w/filter  
 master air/water shut-off valves  
 Chair adapter, pivot

**Assistant's Instrumentation,  
Rear Pivot 1295**

Central vacuum accessory package  
 central vacuum canister w/disposable solids collector  
 asepsis 3-way air/water syringe w/tubing, 5'  
 autoclavable standard HVE lever valve w/tubing, 5'  
 autoclavable saliva ejector lever valve w/tubing, 5'  
 Internal umbilical, 5' w/3' external  
 Rear chair pivot mount

**Cuspidor, Rear Pivot 1582** (Option Not Shown)  
 One-piece ceramic bowl w/touch lever water activator  
 Automatic cup fill  
 Timed bowl rinse  
 Swing away or over patient's lap  
 Water pressure regulator w/filter  
 Master water shut-off valve  
 Internal umbilical, 5' w/3' external  
 Rear chair pivot mount

**Cuspidor, Rear Pivot 1592** (Option Not Shown)  
 Same as 1580 above, w/additional features:  
 Swivel assistant's quad holder bar  
 Assistant's instrumentation  
 central vacuum canister w/disposable solids collector  
 asepsis 3-way air/water syringe w/tubing, 5'  
 autoclavable standard HVE lever valve w/tubing, 5'  
 autoclavable saliva ejector lever valve w/tubing, 5'

**Simplicity™ LED**

Smooth rotation and 3 axis adjustment allow for a wide range of light positioning to better view the oral cavity  
 Control manually or w/hands free sensor  
 Ergonomic handles for comfort of operator  
 Cast aluminum light head for better aesthetics/durability  
 Low, medium and high settings for flexibility from 18,000, 27,000 and 36,000 lux  
 Rectangular light pattern (6 1/2" x 3 1/2") lights the oral cavity w/minimal shadows  
 Color temperature ranges from 5,000 to 5,700 Kelvin

**Patient Chair**

Hydraulic cantilever mechanism  
 Heavy duty cast iron chair frame  
 Aluminum casting chair back, base, armrest and headrest  
 Emergency shut-down control  
 Double articulating headrest that adjusts easily to fit the patient  
 Flip down armrests for easy patient access  
 Standard NaugaHyde® upholstery w/wide back  
 30-degree swivel rotation w/lock

**Base Package**

4295PI Delivery System  
 3800 Chair (nontubed)

**\$ 13,320****As Shown Package \$ 17,517**

Additionally equipped with:  
 9700 Light  
 1295 Assistant's Instrumentation

**FUSION DELIVERY SYSTEM . . . . . 4295PI \* \$ 6605**

Handpiece illumination system

One handpiece (5-hole), accommodates 1 add'l tubing	1700-021-CS	715
Add'l handpiece illumination tubing (5-hole) . . . . .	0017-615-CS	275
One handpiece (6-pin), accommodates 1 add'l tubing	1700-019-CS	715
Add'l handpiece illumination tubing (6-pin) . . . . .	0017-603-CS	195

► Ancillaries use existing handpiece position (see page 59)

Write in selection here . . . . .

Write in selection here . . . . .

Additional handpiece position, wet . . . . .	1100-189-CS	425
--	-------------	-----

Additional handpiece position, air drive only . . . . .	1100-190-CS	215
---	-------------	-----

Assistant's hygiene instrumentation (control head mount) . . . . .	1408-CS	385
--	---------	-----

Quick disconnect with flow control 1/4"		
---	--	--

(specify air or water) . . . . .	0016-445-CS	145
----------------------------------	-------------	-----

Swivel tray, arm-mounted (over control head) . . . . .	1110-468	345
--	----------	-----

Chrome accents . . . . .	9000-998-CS	140
--------------------------	-------------	-----

Additional umbilical (maximum 9 feet) (per foot) . . . . .	1115-999-CS	30
--	-------------	----

<b>ASSISTANT'S INSTRUMENTATION, REAR PIVOT . . . . .</b>	<input type="radio"/> 1295 *	1605
--	------------------------------	------

<b>OR CUSPIDOR, REAR PIVOT . . . . .</b>	<input type="radio"/> 1582	3160
--	----------------------------	------

<b>OR CUSPIDOR WITH ASSISTANT'S INSTRUMENTATION,</b>		
--	--	--

<b>REAR PIVOT . . . . .</b>	<input type="radio"/> 1592	3595
-----------------------------	----------------------------	------

Quick disconnect with flow control 1/4"		
---	--	--

(specify air or water) (1095) . . . . .	0016-455-CS	145
---	-------------	-----

Additional autoclavable standard HVE		
--------------------------------------	--	--

with tubing, 5' & autoclavable oral dry cup . . . . .	1400-113-CS	90
---	-------------	----

Self-contained clean water system with quick disconnect,		
--	--	--

single bottle, 1 1/2 liter (1095) . . . . .	1700-319	275
---	----------	-----

Tray assembly (1095) (see page 74) . . . . .	1110-680	575
--	----------	-----

Dual articulating quad holder assembly, standard HVE (1095) 0024-227-CS		260
---	--	-----

Deluxe cuspidor unit-mounted activation (1580/1590) . . . . .	1100-480-CS	270
---	-------------	-----

<b>SIMPLICITY™ LED</b>		
------------------------	--	--

Chair Mounted Light (includes pivot light post) . . . . .	<input type="radio"/> 9700 *	2592
---	------------------------------	------

Ceiling Mounted Light, 8' . . . . .	<input type="radio"/> 9700-C-S	2840
-------------------------------------	--------------------------------	------

Ceiling Mounted Light, 8' 6" - 10' . . . . .	<input type="radio"/> 9700-C-L	2840
--	--------------------------------	------

<b>PATIENT CHAIR</b>		
----------------------	--	--

Patient Chair (pretubed) . . . . .	<input type="radio"/> 3800PT	7195
------------------------------------	------------------------------	------

Patient Chair (nontubed) . . . . .	<input type="radio"/> 3800 *	6715
------------------------------------	------------------------------	------

Wired Foot Control . . . . .	<input type="radio"/> 3810-010	515
------------------------------	--------------------------------	-----

Additional Touchpad . . . . .	<input type="radio"/> 3814-102	515
-------------------------------	--------------------------------	-----

\* ITEM included in AS SHOWN PACKAGE price

**TOTAL**

► Often overlooked

Designer Friendly® Surface Treatments

UNIT PLASTIC FINISHES:  GreigePOWDER COAT PAINT FINISHES:  GreigeUPHOLSTERY COLORS:  Black  Cocoa  Deep Sapphire  Saddle  Sand  Taupe

# The Fusion Package

Fixed with Sidebox Package



## Popular Equipment Choices



Cuspidor 1577  
(Sidebox Comes  
with Delivery)



Ceiling Light



Shipping weight estimate:

Delivery unit & chair 440 lb / 202 kg

Base package 440 lb / 202 kg

As shown package 561 lbs / 255 kg

Standard Features	Factory Options	Part No	Retail	Price
<b>Delivery System 4295SI</b> Automatic control for three handpieces asepsis 3-way air/water syringe w/tubing, 5' handpiece tubing, 4-hole Midwest, asepsis style, 5' pneumatic master control block w/diaphragm master on/off individual handpiece pressure adjustment individual water coolant adjustment non-retraction water coolant air/water coolant adjustment handpiece flush handpiece exhaust canister/oil collector handle brake release function wet/dry foot control Flexible control arm w/pneumatic brake vertical travel 28" w/horizontal travel 62" Sidebox w/self-contained clean water system w/quick disconnect, single bottle, 1 1/2 liter Water quick disconnect 1/4" w/flow control Air quick disconnect 3/8" w/flow control Toggle valve for city/bottle water routing Telescoping arm w/swivel quad holder bar Central vacuum accessory package central vacuum canister w/disposable solids collector asepsis 3-way air/water syringe w/tubing, 5' autoclavable standard HVE lever valve w/tubing, 5' autoclavable saliva ejector lever valve w/tubing, 5' Internal umbilical, 5' w/3' external from chair to utility center Utility center air/water pressure regulators w/filter master air/water shut-off valves Chair adapter	<b>Base Package</b> 4295SI Delivery System 3800 Chair (nontubed)	<b>\$ 14,690</b>		
<b>Cuspidor, Deluxe 1577</b> (option not shown) One-piece ceramic bowl w/touch lever water activator Automatic cup fill Timed bowl rinse	<b>As Shown Package</b> additionally equipped with: 9700 Light 1113-900 Light Post	<b>\$ 17,479</b>		
<b>Cuspidor, Basic 1679</b> (option not shown) Removable ceramic bowl Manual cup fill Timed bowl rinse	<b>FUSION DELIVERY SYSTEM</b> ..... 4295SI * \$ 7975 Handpiece illumination system One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS 715 Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS 275 One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS 715 Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS 195			
<b>Simplicity™ LED</b> Smooth rotation and 3 axis adjustment allow for a wide range of light positioning to better view the oral cavity Control manually or w/hands free sensor Ergonomic handles for comfort of operator Cast aluminum light head for better aesthetics/durability Low, medium and high settings for flexibility from 18,000, 27,000 and 36,000 lux Rectangular light pattern (6 1/2" x 3 1/2") lights the oral cavity w/minimal shadows Color temperature ranges from 5,000 to 5,700 Kelvin	<b>Ancillaries (use existing handpiece position) (see page 59)</b> Write in selection here ..... Write in selection here ..... Additional handpiece position, wet ..... 1100-189-CS 425 Additional handpiece position ..... 1100-190-CS 215 Swivel tray, arm-mounted option (over control head) ..... 1110-468 345 Less sidebox telescoping arm and assistant's instrumentation 1400-000-CS -350 Additional autoclavable standard HVE with tubing, 5' & autoclavable oral dry cup ..... 1400-113-CS 90 Additional umbilical (maximum 9 feet) (per foot) ..... 1115-999-CS 30 <b>CUSPIDOR</b> , (sold only with delivery) <b>DELUXE</b> ..... ○ 1577 <b>1780</b> <b>OR BASIC</b> ..... ○ 1679 <b>1070</b> Deluxe cuspidor unit mount activation (1577) ..... 1100-480-CS 270 <b>SIMPPLICITY™ LED</b> Chair Mounted Light ..... ○ 9700 * 2592 Ceiling Mounted Light, 8' ..... ○ 9700-C-S 2840 Ceiling Mounted Light, 8' 6" - 10' ..... ○ 9700-C-L 2840 Light Post (fixed) ..... ○ 1113-900 * 197 <b>PATIENT CHAIR</b> Patient Chair (pretubed) ..... ○ 3800PT 7195 Patient Chair (nontubed) ..... ○ 3800 * 6715 Wired Foot Control ..... ○ 3810-010 515 Additional Touchpad ..... ○ 3814-102 515			
<b>Patient Chair</b> Hydraulic cantilever mechanism Heavy duty cast iron chair frame Aluminum casting chair back, base, armrest and headrest Emergency shut-down control Double articulating headrest that adjusts easily to fit the patient Flip down armrests for easy patient access Standard NaugaHyde® upholstery w/wide back 30-degree swivel rotation w/locks	* ITEM included in AS SHOWN PACKAGE price Often overlooked		<b>TOTAL</b>	



Designer Friendly® Surface Treatments

UNIT PLASTIC FINISHES:  Greige

POWDER COAT PAINT FINISHES:  Greige

CHAIR UPHOLSTERY COLORS:  Black  Cocoa  Deep Sapphire  Saddle  Sand  Taupe

# The Fusion Package

## Delivery System

### Delivery System


Please use graph below to determine number of ancillaries and options possible for the delivery system.



### Standard Features

Automatic control for three handpieces  
asepsis 3-way air/water syringe w/tubing, 5' handpiece tubing, 4-hole Midwest, asepsis style, 5' pneumatic master control block w/diaphragm master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air coolant adjustment  
handpiece flush  
handpiece exhaust canister/oil collector  
handle brake release function

### % Control Head Space Absorbed per Option



#### Please note:

If your ancillary and option choices exceed 100%, please call customer service for additional information.

\*: Uses existing handpiece position

★ Standard on SI deliveries



---

Simplified Equipment



Forest Grey



White

(Shown on Front)

**Standard IC** Powder Coat Paint Color: Forest Grey and White

# Ordering Forest Equipment Made Easy **Simplified**

Too many options can sometimes be overwhelming. That's why this select product offering has just two color choices. You can make upgrading a treatment room or matching existing equipment easy.

# Simplified Delivery Systems, Chair Mounted

Fixed Mount with Sidebox, External Umbilical



## Choose Your Control Head

### Model Price

Package offerings see pages 18, 20 and 68

**4485SE Standard IC \$ 7975**

**4987SE Euro \$ 9050**

### Standard Features

Pneumatic  
3 HP  
Syringe

Pneumatic  
4 HP  
Syringe

Delivery Systems can be mounted on any 2" fixed post mount brand of chair

### Optional Equipment Choose One:



- 1577 + \$ 1780**  
**Cuspidor, Deluxe**  
(Sidebox comes with delivery)



- 1679 + \$ 1070**  
**Cuspidor, Basic**  
(Sidebox comes with delivery)



- 9083 + \$ 2195**  
**Light, Halogen**



- 1110-900 + \$ 220**  
**Light Post**



- 9078 + \$ 2995**  
**Light, LED**

- 1110-900 + \$ 200**  
**Light Post**

## Standard Features

## Factory Options

## Part No

## Retail Price

### **Delivery System**

#### **Choice of Control Head**

Automatic control for handpieces

asepsis 3-way air/water syringe with tubing, 5' (Euro tubing 3')

handpiece tubing, 4-hole Midwest, asepsis style, 5' (Euro tubing 3')

pneumatic master control block w/diaphragm master on/off

individual handpiece pressure adjustment

individual water coolant adjustment

non-retraction water coolant

air coolant adjustment

handpiece flush

handpiece exhaust canister/oil collector

single handle brake release function

dual handle brake release function (Euro)

swivel instrument tray (Euro)

wet/dry foot control

Flexible control arm with pneumatic brake

vertical travel 28" / horizontal travel 62"

Sidebox with self-contained clean water system w/ quick disconnect, single bottle, 1 1/2 liter electrical outlet/UL

Water quick disconnect 1/4" w/flow control

Air quick disconnect 3/8" w/flow control

Toggle valve for city/bottle water routing

Telescoping arm with swivel quad holder bar

Central vacuum accessory package

central vacuum canister with disposable solids collector

asepsis 3-way air/water syringe w/tubing, 5'

autoclavable standard HVE lever valve w/tubing, 5'

autoclavable saliva ejector lever valve w/tubing, 5'

External umbilical, 5'

Deluxe utility center

air pressure regulator w/filter

master air shut-off valve

Chair adapter, fixed (specify chair make/model)

### **Deluxe Cuspidor 1577**

One-piece ceramic bowl with touch

lever water activator

Automatic cup fill

Timed bowl rinse

### **Basic Cuspidor 1679**

Removable ceramic bowl

Manual cup fill

Timed bowl rinse

### **DELIVERY SYSTEM ..... MODEL CHOICE**

**PRICE :** \_\_\_\_\_

Handpiece illumination system

One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS \$ 715

Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS 275

One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS 715

Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS 195

► Ancillaries use existing handpiece position (see page 57)

Write in selection here .....

Write in selection here .....

Additional handpiece position, wet (4485)..... 1100-189-CS 425

Additional handpiece position, air drive only(4485) .... 1100-190-CS 215

Swivel tray, arm-mounted (Euro Standard)..... 1110-468 345

Less telescoping arm and assistant's instrumentation .. 1400-000-CS - 350

Chrome accents..... 9000-998-CS 140

Additional autoclavable standard HVE

with tubing, 5' & autoclavable oral dry cup ..... 1400-113-CS 90

Additional umbilical (per foot)..... 1115-999-CS 30

Without chair adapter ..... 1115-099-CS -150

### **CUSPIDOR, (sold only with delivery)**

**DELUXE** ..... ○ 1577 **1780**

**OR BASIC** ..... ○ 1679 **1070**

Deluxe cuspidor unit mount activation (1577) ..... 1100-480-CS 270

**LIGHT, HALOGEN** (mounting options, see pages 33) ..... ○ 9083 **2195**

**OR LED** ..... ○ 9078 **2995**

►● Light post (specify light make/model)

● Light post, 44" (Forest) ..... 1110-900 220

Light post, 44" (Pelton & Crane) ..... 1110-890 220

**MONITOR MOUNT ASSEMBLY** (monitor screen not included) (see page 37)

Ring mount (used with light post)..... 1300RM 1065

Top post mount (used with monitor post)..... 1300TPM 1065

● Monitor post (for top post mount)..... 1331-998 220

**TOTAL** \_\_\_\_\_

● Choose only one post when sold with unit

► Often overlooked

**Equipment dimensions, see page 97**

Shipping weight estimate:

4485SE Delivery Unit 105 lb / 48 kg

4987SE Delivery Unit 110 lb / 50 kg



### Simplified Surface Treatments

#### UNIT PLASTIC FINISHES:

○ Forest Grey ○ White

#### POWDER COAT PAINT FINISHES:

○ Forest Grey ○ White

# Simplified Delivery Systems, Chair Mounted

## Fixed Mount without Sidebox, External Umbilical



**4485FE Standard IC \$ 5455**

**4987FE Euro \$ 6615**

### Model Price

To order packages see pages 78

#### Standard Features

Pneumatic  
3 HP  
Syringe

Pneumatic  
4 HP  
Syringe

Delivery Systems can be mounted on any  
2" fixed post mount brand of chair



### Optional Equipment Choose One:

- 1075FE + \$ 930**  
**Assistant's Telescoping**  
**Arm with Instrumentation**



- 1636FE + \$ 1320**  
**Cuspidor, Basic**



- 1637FE + \$ 1815**  
**Cuspidor, Basic**  
**with Assistant's**  
**Instrumentation**



- 9083 + \$ 2195**  
**Light, Halogen**



- 1110-820 + \$ 290**  
**Light Post**



- 9078 + \$ 2995**  
**Light, LED**

- 1110-820 + \$ 290**  
**Light Post**

Standard Features	Factory Options	Part No	Retail	Price
<b>Delivery System</b>	<b>DELIVERY SYSTEM</b> ..... <b>MODEL CHOICE</b> _____		<b>PRICE :</b>	
Choice of Control Head	Handpiece illumination system			
Automatic control for handpieces	One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS	\$ 715		
asepsis 3-way air/water syringe w/tubing, 5' (Euro tubing 3')	Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS	275		
handpiece tubing, 4-hole Midwest, asepsis style, 5' (Euro tubing 3')	One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS	715		
master control block	Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS	195		
master on/off	<b>► Ancillaries use existing handpiece position (see page 57)</b>			
individual handpiece pressure adjustment	Write in selection here .....			
individual water coolant adjustment	Write in selection here .....			
non-retraction water coolant	Additional handpiece position, wet (4485) ..... 1100-189-CS	425		
air coolant adjustment	Additional handpiece position, air drive only (4485) .. 1100-190-CS	215		
handpiece flush	<b>Assistant's Hygiene Instrumentation</b>			
handpiece exhaust canister/oil collector	(control head mount) ..... 1408-CS	385		
single handle brake release function	Swivel tray, arm-mounted (standard for Euro) ..... 1110-468	345		
self-contained clean water system w/ quick disconnect, single bottle, 1 1/2 liter	Chrome accents ..... 9000-998-CS	140		
wet/dry foot control	Quick disconnect with flow control, 1/4" (specify air or water)..... 0016-455-CS	145		
Flexible control arm w/pneumatic brake	Additional umbilical (per foot)(9' max) ..... 1115-999-CS	30		
vertical travel 28" / horizontal travel 62"	<b>► Deluxe utility center</b>			
External umbilical, 5'	Routing valve kit with water shutoff & regulator/factory installation ..... 1115-315-CS	305		
<b>► Chair adapter, fixed (specify chair make/model)</b>	Routing valve kit for use with cuspidor/factory installation ..... 1115-316-CS	145		
<b>Assistant's Instrumentation 1075FE</b>	Without chair adapter ..... 1115-099-CS	- 150		
Central vacuum accessory package	<b>ASSISTANT'S INSTRUMENTATION</b> ..... <b>1075FE</b>	<b>930</b>		
central vacuum canister w/disposable solids collector	Additional autoclavable standard HVE with tubing, 5' and autoclavable oral dry cup ..... 1400-113-CS	90		
asepsis 3-way air/water syringe w/tubing, 5'	<b>► CUSPIDOR, BASIC</b> ..... <input type="radio"/> <b>1636FE</b>	<b>1320</b>		
autoclavable standard HVE lever valve w/tubing, 5'	with assistant's instrumentation ..... <input type="radio"/> <b>1637FE</b>	<b>1815</b>		
autoclavable saliva ejector lever valve w/tubing, 5'	Additional autoclavable standard HVE with tubing, 5' & autoclavable oral dry cup (1637FE) 1400-113-CS	90		
external umbilical, 5'	<b>LIGHT, HALOGEN</b> (mounting options, see pages 33) ..... <input type="radio"/> <b>9083</b>	<b>2195</b>		
Post mount, 2"	<b>OR LED</b> ..... <input type="radio"/> <b>9078</b>	<b>2995</b>		
<b>Cuspidor, Basic 1636FE</b>	<b>► LIGHT POST</b> with 2" post clamp (specify light make/model)			
Gravity drain cuspidor w/swing away arm	Light post, 35" (Forest) ..... 1110-820	290		
timed bowl rinse	Light post, 35" (Pelton & Crane) ..... 1110-822	290		
manual cup fill	<b>► MONITOR MOUNT ASSEMBLY</b>			
removable ceramic bowl	(monitor screen not included) (see page 37)			
Post mount adapter, 2"	Ring mount (used with light post) ..... 1300RM	1065		
External 5' umbilical	Top post mount (used with monitor post) ..... 1300TPM	1065		
Water pressure regulator w/filter	<b>• Monitor post (for top post mount)</b> ..... 1331-997	290		
Master water shut-off valve	<b>► Choose only one post when sold with unit</b>			
<b>Cuspidor, Basic 1637FE</b>	<b>► When optional cuspidor added, standard utility center housing and cover is used to accommodate two umbilicals.</b>			
Same as 1636FE above, w/additional features:	<b>► Often overlooked</b>			
Individual assistant's instrumentation holders				
Central vacuum accessory package				
central vacuum canister w/disposable solids collector				
asepsis 3-way air/water syringe w/tubing, 5'				
autoclavable standard HVE lever valve w/tubing, 5'				
autoclavable saliva ejector lever valve w/tubing, 5'				



#### Simplified Surface Treatments

**UNIT PLASTIC FINISHES:**

Forest Grey  White

**POWDER COAT PAINT FINISHES:**

Forest Grey  White

**Equipment dimensions, see page 97**

Shipping weight estimate:

4485FE Delivery Unit 95 lb / 44 kg

4987FE Delivery Unit 100 lb / 46 kg

# Simplified Package

## Package FE Standard IC Control Head

Standard Features:

### 4485FE Delivery System (option shown)

Choice of Control Head

Automatic control for handpieces

asepsis 3-way air/water syringe w/tubing, 5'  
(Euro tubing 3')

handpiece tubing, 4-hole Midwest, asepsis style, 5'  
(Euro tubing 3')

master control block

master on/off

individual handpiece pressure adjustment

individual water coolant adjustment

non-retraction water coolant

air coolant adjustment

handpiece flush

handpiece exhaust canister/oil collector

single handle brake release function Euro-dual

self-contained clean water system w/quick disconnect,  
single bottle, 1 1/2 liter

wet/dry foot control

Flexible control arm w/pneumatic brake

vertical travel 28" / horizontal travel 62"

External umbilical, 5'

► Deluxe utility center

air pressure regulator w/filter

master air shut-off valve

Chair adapter, fixed

### Assistant's Instrumentation 1075FE

(option shown)

Central vacuum accessory package

central vacuum canister w/disposable solids collector

asepsis 3-way air/water syringe w/tubing, 5'

autoclavable standard HVE lever valve w/tubing, 5'

autoclavable saliva ejector lever valve w/tubing, 5'

External umbilical, 5'

Post mount, 2"

### Cuspidor, Basic 1636FE/1637FE (option not shown)

Removable ceramic bowl

Manual cup fill

Timed bowl rinse

Water pressure regulator w/filter

Master water shut-off valve

Post mount, 2"

### 1637FE includes:

Individual assistant's instrumentation holders

Assistant's instrumentation

central vacuum canister w/disposable  
solids collector

asepsis 3-way air/water syringe w/tubing, 5'

autoclavable standard HVE lever valve w/tubing, 5'

autoclavable saliva ejector lever valve w/tubing, 5'

### Light, Halogen 9083

3-Axis light head movement

Automatic on/off activation w/override

Color shade matching temperature

setting 4200 to 5500 Kelvin

Shadow free illumination

### Chair, Hydraulic 3900

Ultraleather® upholstery w/wide back for elbow support

Double articulating pinch-style headrest

Thin back w/full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation w/lock

Multifunction touchpad w/4 programmable preset positions

Seat height: minimum 17" maximum 29"



### Control Head Choice



4987FE Euro (option not shown)

### Popular Equipment Choices



1636FE Cuspidor, Basic



1637FE Cuspidor, Basic  
with Assistant's  
Instrumentation



Doctor's Economy  
Stool 6156C



Led Light  
9078

		Factory Options	Part No	Retail	Price
<b>Base Package</b> 4485FE Delivery 3900 Chair 9083 Light 1110-820 Light Post	<b>\$ 15,930</b>	<b>DELIVERY SYSTEM (CHOOSE ONE)</b> .....			
		Standard IC Control .....	4485FE *	\$5455	
		Euro Control .....	4987FE	6615	
		Handpiece illumination system			
		One handpiece (5-hole), accommodates 1 add'l tubing	1700-021-CS	715	
		Add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	275	
		One handpiece (6-pin), accommodates 1 add'l tubing	1700-019-CS	715	
		Add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	195	
		Ancillaries use existing handpiece position (see page 57)			
		Write in selection here .....			
		Write in selection here .....			
		Additional handpiece position, wet (4485) .....	1100-189-CS	425	
		Additional handpiece position, air drive only (4485) ..	1100-190-CS	215	
		Assistant's Hygiene Instrumentation			
		(control head mount) (4485) .....	1408-CS	385	
		Swivel tray, arm-mounted (standard for Euro) .....	1110-468	345	
		Chrome accents .....	9000-998-CS	140	
		Quick disconnect with flow control, 1/4" (specify air or water)	0016-455-CS	145	
		Additional umbilical (maximum 9 feet) (per foot) ....	1115-999-CS	30	
		Routing valve kit with water shutoff & regulator/factory installation	1115-315-CS	305	
		Routing valve kit for use with cuspidor/factory installation	1115-316-CS	145	
		<b>ASSISTANT'S INSTRUMENTATION</b> .....	<b>1075FE *</b>	<b>930</b>	
		Additional autoclavable standard HVE			
		with tubing, 5' & autoclavable oral dry cup.....	1400-113-CS	90	
		► <b>CUSPIDOR, BASIC</b> .....	<b>1636FE</b>	<b>1320</b>	
		with assistant's instrumentation .....	1637FE	1815	
		Additional autoclavable standard HVE			
		with tubing, 5' & autoclavable oral dry cup (1637FE) 1400-113-CS	90		
		<b>LIGHT, HALOGEN</b> (mounting options, see pages 34-35) .....	<b>9083 *</b>	<b>2195</b>	
		<b>OR LED</b> .....	<b>9078</b>	<b>2995</b>	
		► • <b>LIGHT POST WITH 2" POST CLAMP</b> (specify light make/model)			
		Light post, 35" (Forest) .....	1110-820 *	290	
		Light post, 35" (Pelton & Crane) .....	1110-822	290	
		<b>MONITOR MOUNT ASSEMBLY</b> (monitor screen not included) (see page 37)			
		Ring mount (used with light post) .....	1300RM	1065	
		Top post mount (used with monitor post) .....	1300TPM	1065	
		• Monitor post (for top post mount) .....	1331-997	290	
		<b>CHAIR, HYDRAULIC</b> .....	<b>3900 *</b>	<b>7990</b>	
		► One chair control included (specify type below) .....	3914-998-CS *	N/C	
		<input type="radio"/> Footswitch control			
		<input type="radio"/> Touchpad, chair-mounted			
		<input type="radio"/> Touchpad, unit-mounted			
		Additional chair control (limit 2 on chair, limit 2 on unit)	3914-999-CS	515	
		<input type="radio"/> Footswitch control			
		<input type="radio"/> Touchpad, chair-mounted			
		<input type="radio"/> Touchpad, unit-mounted			
		► Additional chair options (See page 25 for pricing, write selection here)			
		<b>STOOL, DR'S, ECONOMY/DELUXE/COMFORT</b> .....	<b>6156C/6159C/6176</b>	<b>530/895/895</b>	
		<b>STOOL, ASST'S, ECONOMY/DELUXE</b> .....	<b>6157C/6160C</b>	<b>680/1010</b>	
		<b>STOOL, ASST'S WITH BACK, DELUXE/COMFORT</b> .....	<b>6161C/6187</b>	<b>1195/1195</b>	
		<b>STOOL, OFFICE, COMFORT</b> .....	<b>6198</b>	<b>1060</b>	
		► Fully assembled, blanket wrapped (chair, delivery, light) ..	BLANKET	450	
		<b>TOTAL</b>			

\* ITEM included in AS SHOWN PACKAGE price

• Choose only one post when sold with unit

► When optional cuspidor added, standard utility center housing and cover is used to accommodate two umbilicals.

► Often overlooked

#### Equipment dimensions, see page 97

Shipping weight estimate:

4485FE Delivery Unit 95 lb / 44 kg

4987FE Delivery Unit 100 lb / 46 kg

# Simplified Cuspidors and Assistant's Instrumentation

## Fixed Chair Post Mount, External Umbilical



**Model 1636 \$ 1320**



**Model 1637 \$ 1815**



**Model 1075 \$ 930**

Standard Features	Factory Options	Part No	Retail	Price
<b>Cuspidor, Basic 1636</b> Gravity drain cuspidor w/swing away arm timed bowl rinse manual cup fill removable ceramic bowl Post mount adapter, 2" External umbilical, 5'  ► Does not include post or chair adapter	<b>CUSPIDOR, BASIC .....</b> MODEL CHOICE _____ Additional autoclavable standard HVE with tubing, 5' & autoclavable oral dry cup (1637). 1400-113-CS \$ 90 Quick disconnect with flow control, 1/4" (specify air or water)..... 0016-455-CS..... 145 ► Self-contained clean water system with quick disconnect 1700-310-CS 275 Standard utility center housing & cover ..... 1115-020 145 ► Standard utility center ..... 1115-001 500 Additional umbilical (per foot)..... 1115-999-CS 30 ► MOUNTING POSTS 2" diameter: Control post with cap, 18" ..... 1110-487 115 ► LIGHT POSTS: 44" 2" diameter Forest model ..... 1110-794 220 Pelton & Crane ..... 1110-785 220 ASSISTANT'S INSTRUMENTATION ..... 1075 930 Additional standard autoclavable HVE with tubing, 5' & autoclavable oral dry cup..... 1400-113-CS 90 ► Self-contained clean water system with quick disconnect 1700-323-CS 275 Assistant's air shut-off & regulator utility center w/o housing & cover..... 1115-008 235 Standard utility center housing & cover ..... 1115-020 145 Additional umbilical (per foot)..... 1115-998-CS 20			
<b>Cuspidor, Basic 1637</b> Same as 1636 above, w/additional features: Individual assistant's instrumentation holders Central vacuum accessory package central vacuum canister w/disposable solids collector asepsis 3-way air/water syringe w/tubing, 5' autoclavable standard HVE lever valve w/tubing, 5' autoclavable saliva ejector lever valve w/tubing, 5'  ► Does not include post or chair adapter	TOTAL _____			
<b>Post Mount Vacuum Telescoping Arm 1075</b> Post mount adapter, 2" Central vacuum accessory package post mounted central vacuum canister w/disposable solids collector asepsis 3-way air/water syringe w/tubing, 5' autoclavable standard HVE lever valve w/tubing, 5' autoclavable saliva ejector lever valve w/tubing, 5' External umbilical, 8'  ► Does not include post or chair adapter	► Standard utility center includes: housing, cover, air & water pressure regulator with filter, master air and water shut-off valves. ► Often overlooked			

**WARNING:** Isolated water sources (as opposed to city water hook-up) are highly recommended to allow full control of the source and quality of water circulating through the syringe and cuspidor water lines. If you do not have a water quality protocol, the optional self-contained water system allows you to isolate this product from the municipal water supply, however, it requires monitoring, maintenance and disinfection of the water lines, reservoir water bottle and source water.

Equipment dimensions, see page 97

Shipping weight estimate:

1636 17 lb / 8 kg

1637 20 lb / 9 kg

1075 15 lb / 7 kg



### Simplified Surface Treatments

UNIT PLASTIC FINISHES:

Forest Grey

White

POWDER COAT PAINT FINISHES:

Forest Grey

White

# Simplified Cuspidors and Assistant's Instrumentation

## Rear Beneath Shelf or Work Surface Mount

Mounts under floating work surface



Mounts beneath shelf



### **Model 1284 \$ 1130**

Shown with optional Clean Water System 1700-312 not included in price  
Assistants Utilities 1115-008 sold separately

### **Model 1286 \$ 1220**

Optional Clean Water System 1700-312 and  
Assistants Utilities 1115-008 sold separately

Standard Features	Factory Options	Part No	Retail	Price
<b>Horizontal Assistant's Delivery System 1284/1286</b>	<b>DELIVERY SYSTEM .....</b> MODEL CHOICE _____		<b>PRICE:</b> _____	
Central vacuum accessory package, swing mount arm-mounted central vacuum canister w/disposable solids collector	Additional standard autoclavable HVE position			
asepsis 3-way air/water syringe w/tubing, 5' autoclavable standard HVE lever valve w/ tubing, 5'	with tubing, 5' & autoclavable oral dry cup..... 1400-113-CS	\$ 90		
autoclavable saliva ejector lever valve w/ tubing, 5'	Slide mount, 7 5/8" extension (1284) ..... 1700-142	290		
Horizontal surface mount	█ Standard utility center housing & cover ..... 1115-020	145		
Link assembly, arm to plate (1286)	█ Self-contained clean water system with quick disconnect, single bottle..... 1700-312	275		
External umbilical, 7'	█ Assistant's air shut-off & regulator utility center without housing & cover ..... 1115-008	235		
			<b>TOTAL</b>	_____

█ Often overlooked

**⚠ WARNING:** Isolated water sources (as opposed to city water hook-up) are highly recommended to allow full control of the source and quality of water circulating through the syringe and cuspidor water lines. If you do not have a water quality protocol, the optional self-contained water system allows you to isolate this product from the municipal water supply, however, it requires monitoring, maintenance and disinfection of the water lines, reservoir water bottle and source water.



### Simplified Surface Treatments

Mounting plate dimensions, see page 111

Shipping weight estimate:

1284 12 lb / 5 kg

1286 14 lb / 6 kg

UNIT PLASTIC FINISHES:

Forest Grey

POWDER COAT PAINT FINISHES:

Forest Grey

# Simplified Rear and Side Delivery Systems

## Side Cabinet & Wall Mount



**Model 3235 \$ 4435**  
12" Rigid Arm Cabinet Mount



**Model 3230 \$ 4580**  
22" Rigid Arm Wall Mount

Standard Features	Factory Options	Part No	Retail	Price
<p><b>Doctor Standard IC Delivery System 3235/3230</b></p> <p>Standard IC control head automatic control for three handpieces asepsis 3-way air/water syringe w/tubing, 5' handpiece tubing, 4-hole Midwest, asepsis style, 5' pneumatic master control block w/diaphragm master on/off individual handpiece pressure adjustment individual water coolant adjustment non-retraction water coolant air coolant adjustment handpiece flush handpiece exhaust canister/oil collector single handle brake release function self-contained clean water system w/quick disconnect, single bottle, 1 1/2 liter wet/dry foot control Flexible control arm w/pneumatic brake (3235) cabinet vertical mount 12" rigid arm &amp; 18" flex arm Flexible control arm w/pneumatic brake (3230) wall vertical mount 22" rigid arm &amp; 18" flex arm External umbilical, 5' Standard utility center housing and cover (3230) Air pressure regulator w/filter Master air shut-off valve</p>	<p><b>DELIVERY SYSTEM . . . . . MODEL CHOICE</b> _____</p> <p>Handpiece Illumination System</p> <p>One handpiece (5-hole), accommodates</p> <p>1 add'l tubing ..... 1700-021-CS \$ 715</p> <p>Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS 275</p> <p>One handpiece (6-pin), accommodates</p> <p>1 add'l tubing ..... 1700-019-CS 715</p> <p>Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS 195</p> <p>Assistant's hygiene instrumentation (control head mount) . 1408-CS 385</p> <p>Quick disconnect with flow control, 1/4"</p> <p>(specify air or water) ..... 0016-455-CS 145</p> <p>Swivel tray, arm-mounted (over control head) ..... 1110-468 345</p> <p>Standard utility center housing &amp; cover (3235) ..... 1115-020 145</p> <p>Additional umbilical (per foot)..... 1115-999-CS 30</p> <p>Routing valve kit with water shutoff &amp; regulator/factory installation</p> <p>..... 1115-315-CS 305</p> <p>Wall mount panel, 16" x 18"..... 1110-586-CS 285</p> <p>Ring mount rigid arm in lieu of top post mount</p> <p>22" without vertical mount..... 1110-536-RM 100</p> <p>26" without vertical mount..... 1110-537-RM 100</p> <p>26" rigid with 24" flex without vertical mount.... 1110-771-RM 100</p>			

**TOTAL**

► Often overlooked

Mounting plate dimensions, see page 110

Shipping weight estimate:  
3235, 3230 50 lb / 23 kg



Simplified Surface Treatments

UNIT PLASTIC FINISHES:

Forest Grey

POWDER COAT PAINT FINISHES:

Forest Grey

# Simplified Rear and Side Delivery Systems

Rear or Side Beneath Counter Mount



## **Model 2400 \$ 3055**

Beneath Counter Mount with swivel included  
Countertop/work surface not included

### Standard Features

**Doctor Standard IC Delivery System 2400**  
Beneath counter swivel mount  
(work surface not included)  
Standard IC control head  
automatic control for three handpieces  
asepsis 3-way air/water syringe w/tubing, 5'  
handpiece tubing, 4-hole Midwest, asepsis style, 5'  
pneumatic master control block w/diaphragm  
master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air coolant adjustment  
handpiece flush  
handpiece exhaust canister/oil collector  
self-contained clean water system w/quick  
disconnect, single bottle, 1 1/2 liter  
wet/dry foot control  
External braided umbilical, 7'  
Air pressure regulator w/ filter  
Master air shut-off valve

### Optional Beneath Counter Slide 1700-140

### Optional Folding Arm Assembly 1110-115-CS

Standard Features	Factory Options	Part No	Retail	Price
	<b>DELIVERY SYSTEM .....</b> <b>MODEL CHOICE: 2400</b> <b>PRICE:</b> _____			
	Handpiece Illumination System			
	One handpiece (5-hole), accommodates			
	1 add'l tubing ..... 1700-021-CS     \$ 715			
	Add'l handpiece illumination tubing (5-hole).....0017-615-CS     275			
	One handpiece (6-pin), accommodates			
	1 add'l tubing ..... 1700-019-CS     715			
	Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS     195			
	Assistant's hygiene instrumentation (control head mount) . 1408-CS     385			
	Quick disconnect with flow control, 1/4"			
	(specify air or water) ..... 0016-455-CS     145			
	Slide mount with swivel, 7 5/8" extension ..... 1700-140     290			
	Folding arm assembly, 8" x 12" with wall mount ..... 1110-115-CS     345			
	Standard utility center housing & cover ..... 1115-020     145			
	Routing valve kit with water shutoff &			
	regulator/factory installation ..... 1115-315-CS     305			
	<b>TOTAL</b> _____			

Often overlooked



### Simplified Surface Treatments

Mounting plate dimensions, see page 110

Shipping weight estimate:  
2400 28 lb / 13kg

UNIT PLASTIC FINISHES:

Forest Grey

POWDER COAT PAINT FINISHES:

Forest Grey

# Simplified Rear and Side Delivery Systems

## Rear Beneath Counter Mount



Panel Mount is recessed in cabinet separate from duo articulating instrument delivery 2184.



**Model 2184 \$ 4430**

Mounts under floating work surface.



**Model 2186 \$ 3895**

Mounts under floating work surface.

Standard Features	Factory Options	Part No	Retail	Price
<b>Duo Delivery System 2184/2186</b> Automatic panel control for three handpieces asepsis 3-way air/water syringe w/tubing, 5' handpiece tubing, 4-hole Midwest, asepsis style, 5' pneumatic master control block w/diaphragm master on/off individual handpiece pressure adjustment individual water coolant adjustment non-retraction water coolant air coolant adjustment handpiece flush wet/dry foot control Central vacuum accessory package, swing mount arm-mounted central vacuum canister w/disposable solids collector asepsis 3-way air/water syringe w/tubing, 5' autoclavable standard HVE lever valve w/tubing, 5' autoclavable saliva ejector lever valve w/tubing, 5' Articulating left/right instrument positioning Self-contained clean water system w/quick disconnect, single bottle, 1 1/2 liter External umbilical, 7' Air pressure regulator w/filter Master air shut-off valve	<b>DELIVERY SYSTEM .....</b> MODEL CHOICE ..... Handpiece Illumination System One handpiece (5-hole), accommodates 1 add'l tubing ..... 1700-021-CS \$ 715 Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS 275 One handpiece (6-pin), accommodates 1 add'l tubing ..... 1700-019-CS 715 Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS 195 Additional standard autoclavable HVE position with tubing, 5' & autoclavable oral dry cup..... 1400-113-CS 90 Slide mount, 7 5/8" extension (2184) ..... 1700-142 290 Touchpad for 3900 chair (25 ft. Cable) ..... 1700-081-1-CS 515 Standard utility center housing & cover ..... 1115-020 145 ► Quick disconnect for air or water..... dealer installed (see page 95) ► Routing valve kit for city water ..... dealer installed (see page 95)			
				<b>TOTAL</b>
	► Often overlooked			

Mounting plate dimensions, see pages 108

Panel cut-out, see page 108

Shipping weight estimate:

2184 30 lb / 14 kg

2186 30 lb / 14 kg



### Simplified Surface Treatments

#### UNIT PLASTIC FINISHES:


Forest Grey

#### POWDER COAT PAINT FINISHES:

Forest Grey

# Simplified Rear and Side Delivery Systems

Rear Beneath Counter Mount



Panel Mount is recessed in cabinet  
separate from duo articulating  
instrument delivery 2184BC.



## Model 2184BC \$ 5430

Universal Heavy-Duty Mounting Plate 11" x 14 1/2"  
included (mounts to under counter)  
Specify cabinet manufacturer

## Model 2186BC \$ 4895

Universal Heavy-Duty Mounting Plate 11" x 14 1/2"  
included (mounts to under counter)  
Specify cabinet manufacturer

Standard Features	Factory Options	Part No	Retail	Price
<p><b>Duo Delivery System 2184BC/2186BC</b></p> <p>Automatic panel control for three handpieces asepsis 3-way air/water syringe w/tubing, 5' handpiece tubing, 4-hole Midwest, asepsis style, 5' pneumatic master control block w/diaphragm master on/off individual handpiece pressure adjustment individual water coolant adjustment non-retraction water coolant air coolant adjustment handpiece flush wet/dry foot control</p> <p>Central vacuum accessory package, swing mount arm-mounted central vacuum canister w/disposable solids collector asepsis 3-way air/water syringe w/tubing, 5' autoclavable standard HVE lever valve w/tubing, 5' autoclavable saliva ejector lever valve w/tubing, 5'</p> <p>Articulating left/right instrument positioning</p> <p>Self-contained clean water system w/quick disconnect, single bottle, 1 1/2 liter</p> <p>Laminate work surface, 19" x 24"</p> <p>► Beneath counter rigid arm assembly w/heavy duty mounting plate, 11" x 14 1/2"</p> <p>External umbilical, 7'</p> <p>Air pressure regulator w/filter</p> <p>Master air shut-off valve</p>	<p><b>DELIVERY SYSTEM</b> ..... <b>MODEL CHOICE</b> ..... <b>PRICE:</b></p> <p>Handpiece Illumination System</p> <ul style="list-style-type: none"> <li>One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS \$ 715</li> <li>Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS 275</li> <li>One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS 715</li> <li>Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS 195</li> </ul> <p>Additional standard autoclavable HVE position</p> <ul style="list-style-type: none"> <li>with tubing, 5' &amp; autoclavable oral dry cup..... 1400-113-CS 90</li> <li>Quick disconnect with flow control 1/4" (for water) ... 0016-454-CS 145</li> <li>Quick disconnect with flow control 3/8" (for air) .... 0016-456-CS 145</li> </ul> <p>■ Storage organizer with bins &amp; frame..... 1110-179-CS 525</p> <p>■ Swivel tray (requires purchase of 1110-179-CS).... 1110-908-CS 350</p> <p>■ Swivel tray (mounts without storage organizer).... 1110-482-CS 425</p> <p>■ Monitor mount assembly (monitor screen not included)</p> <ul style="list-style-type: none"> <li>(requires purchase of 1110-179-CS) ..... 1300-025-CS 450</li> </ul> <p>Touchpad for 3900 chair (25 ft. Cable) ..... 1700-081-1-CS 515</p> <p>Standard utility center housing &amp; cover ..... 1115-020 145</p> <p>Routing valve kit with water shutoff &amp; regulator/factory installation ..... 1115-315-CS 305</p> <p>Standard color solid work surface upgrade ..... 0023-003-SS 325</p> <p>Special color laminate or solid work surface upgrade ..... 0023-003-X / 0023-003-SSV (call for price and leadtime)</p>			

**TOTAL**

■ Often overlooked



### ► Beneath Counter Mount

Mounting plate dimensions, see pages 109

Panel cut-out, see page 108

Shipping weight estimate:

2184BC 80 lb / 36 kg

2186BC 80 lb / 36 kg

### Simplified Surface Treatments

UNIT PLASTIC FINISHES:

Forest Grey

POWDER COAT PAINT FINISHES:

Forest Grey

LAMINATE WORK SURFACE :

FORMICA®:  Surf  Folkstone  Oyster Gray  White Blanc

SOLID WORK SURFACE UPGRADE COLORS 0023-003-SS \$325

CORIAN®:  Antarctica

GIBRALTAR®:  Arctic Melange  Baja Melange

MEGANITE®:  Cinder Granite  Goldrush Gemstone  Moonstone  Platinum Granite

Rocky Road Granite  Sanibel Granite  Silver Mist  Starry Starry Night

# Simplified Rear and Side Delivery Systems

## Panel Mounts



**Model 2500 \$ 2525**

Optional Assistant's Hygiene  
Instrumentation 1408-CS

**Model 2520 \$ 2095**

Standard Features	Factory Options	Part No	Retail	Price
<b>Doctor Control 2500</b> Automatic control for three handpieces asepsis 3-way air/water syringe w/tubing, 7' handpiece tubing, 4-hole Midwest, asepsis style, 7' pneumatic master control block w/diaphragm master on/off individual handpiece pressure adjustment individual water coolant adjustment non-retraction water coolant air coolant adjustment self-contained clean water system w/quick disconnect, single bottle, 1 1/2 liter wet/dry foot control Remote mount holder assembly w/control tubing, 6' External umbilical, 6' Air pressure regulator w/filter Master air shut-off valve	<b>DELIVERY SYSTEM . . . . . MODEL CHOICE</b> Handpiece Illumination System One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS Add'l handpiece illumination tubing (5-hole).....0017-615-CS 275 One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS 715 Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS 195 Standard utility center housing & cover ..... 1115-020 145 Assistant's hygiene instrumentation with extended 10 1/2" holder bar (Drain Tubing 12') standard autoclavable HVE, saliva ejector with 2 holders and canister ..... 1408-CS 385 ► Quick disconnect for air or water..... dealer installed (see page 95) ► Routing valve kit for city water ..... dealer installed (see page 95)			
<b>Doctor Control 2520</b> Automatic control for two handpieces asepsis 3-way air/water syringe w/tubing, 7' handpiece tubing, 4-hole Midwest, asepsis style, 7' pneumatic master control block w/diaphragm master on/off individual handpiece pressure adjustment non-retraction water coolant air coolant adjustment self-contained clean water system w/quick disconnect, single bottle, 1 1/2 liter wet/dry foot control Remote mount holder assembly w/control tubing, 6' External umbilical, 6' Air pressure regulator w/filter Master air shut-off valve	► Often overlooked			
			<b>TOTAL</b>	

Wall/cabinet cut-out dimensions, see pages 108, 109  
Shipping weight estimate:

2500 18 lb / 8 kg  
2520 15 lb / 7 kg



Simplified Surface Treatments

UNIT PLASTIC FINISHES:

Forest Grey

POWDER COAT PAINT FINISHES:

Forest Grey

# Simplified Rear and Side Delivery Systems

Rear or Side Mount Compact Controls



**Model MT600 \$ 2205**

Work surface not included



**Model MT610 \$ 2105**

Work surface not included

Standard Features	Factory Options	Part No	Retail	Price
<b>Doctor Control with 3 HP</b> Top Shelf Mount Beneath Shelf Mount w/Swivel	<b>DELIVERY SYSTEM .....</b> MODEL CHOICE _____ Handpiece Illumination System One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS Slide mount with swivel, 7 5/8" extension (order BS) .... 1700-140 Folding arm assembly, 8" x 12" with wall mount (order TS) 1110-115-CS Standard utility center housing & cover ..... 1115-020 Assistant's hygiene instrumentation with extended 10 1/2" holder bar (Drain Tubing 12') standard autoclavable HVE, saliva ejector with 2 holders & canister ..... 1408-CS ■► Quick disconnect for air or water..... dealer installed (see page 95) ■► Routing valve kit for city water ..... dealer installed (see page 95)	MT600TS MT600BS	PRICE : \$ 715 275 715 195 290 345 145 385 _____	_____
<b>Doctor Control with 2 HP</b> Top Shelf Mount Beneath Shelf Mount w/Swivel		MT610TS MT610BS		
<b>All Units Include the Following Features:</b> Automatic control for three handpieces (MT600), or Automatic control for two handpieces (MT610) asepsis 3-way air/water syringe w/tubing, 4' handpiece tubing, 4-hole Midwest, asepsis style, 4' pneumatic master control block w/diaphragm master on/off individual handpiece pressure adjustment non-retraction water coolant air coolant adjustment self-contained clean water system w/quick disconnect, single bottle, 1 1/2 liter wet/dry type foot control External umbilical, 12' Air pressure regulator w/filter Master air shut-off valve	■► Often overlooked		TOTAL _____	



## Simplified Surface Treatments

Mounting plate dimensions, see page 110  
Shipping weight estimate:  
MT600, MT610      20 lb / 10 kg

UNIT PLASTIC FINISHES:

Forest Grey

POWDER COAT PAINT FINISHES:

Forest Grey

# Simplified Rear and Side Delivery Systems

## Mobile Cart



Shown with optional storage organizer and swivel tray  
1110-179-CS and 1110-908-CS



**Model 5386 \$ 3355**

**Model 5886 \$ 5455**

Standard Features	Factory Options	Part No	Retail	Price
<b>Assistant's Cart with Instrumentation 5386</b>	<b>CART .....</b> <b>MODEL CHOICE_____</b> <b>PRICE:</b>			
Central vacuum accessory package, swing mount arm-mounted central vacuum canister w/disposable solids collector asepsis 3-way air/water syringe w/tubing, 5' autoclavable standard HVE lever valve w/tubing, 5' autoclavable saliva ejector lever valve w/tubing, 5' self-contained clean water system w/quick disconnect, single bottle 1 1/2 liter	Handpiece Illumination System (5886) One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS \$ 715 Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS 275 One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS 715 Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS 195			
Assistant's holder quad assembly Laminate work surface, 19" x 24" Adjustable height frame Twin wheel hard floor casters Electrical outlet/UL External umbilical, 7' Assistant's utility center w/housing and cover air pressure regulator w/filter master shut-off valve	Additional standard autoclavable HVE with tubing, 5' & autoclavable oral dry cup..... 1400-113-CS 90 Additional assistant's quad holder bar ..... 0024-149 205 Quick disconnect with flow control 1/4" (for water).... 0016-454-CS 145 Quick disconnect with flow control 3/8" (for air) .... 0016-456-CS 145 ➡ Storage organizer with bins & frame..... 1110-179-CS 525 ➡ Swivel tray (requires purchase of 1110-179-CS)..... 1110-908-CS 350 ➡ Swivel tray (mounts without storage organizer)..... 1110-482-CS 425 ➡ Monitor mount assembly (monitor screen not included) (requires purchase of 1110-179-CS) ..... 1300-025-CS 450 Additional umbilical (per foot) ..... 1115-999-CS 30 Routing valve kit with water shutoff & regulator/factory installation ..... 1115-315-CS 305 Standard color solid work surface upgrade ..... 0023-003-SS 325 Special color laminate or solid work surface upgrade ..... 0023-003-X / 0023-003-SSV			
<b>Duo Swing Cart 5886</b> Same as above, w/additional features: Standard IC Control Head automatic control for three handpieces asepsis 3-way air/water syringe w/tubing, 4' handpiece tubing, 4-hole Midwest, asepsis style, 4' pneumatic master control block w/diaphragm master on/off individual handpiece pressure adjustment individual water coolant adjustment non-retraction water coolant air coolant adjustment handpiece flush handpiece exhaust canister/oil collector wet/dry foot control Standard utility center housing and cover Air pressure regulator w/filter Master air shut-off valve	Often overlooked			<b>TOTAL</b>

➡ Often overlooked

## Simplified Surface Treatments

**UNIT PLASTIC FINISHES:**  Forest Grey

**POWDER COAT PAINT FINISHES:**  Forest Grey

**LAMINATE WORK SURFACE :**

**FORMICA®:**  Surf  Folkstone  Oyster Gray  White Blanc

**SOLID WORK SURFACE UPGRADE COLORS 0023-003-SS \$325**

**CORIAN®:**  Antarctica

**GIBRALTAR®:**  Arctic Melange  Baja Melange

**MEGANITE®:**  Cinder Granite  Goldrush Gemstone  Moonstone  Platinum Granite

Rocky Road Granite  Sanibel Granite  Silver Mist  Starry Starry Night

**Equipment dimensions, see page 111**


Shipping weight:

5386 115 lb / 53 kg

5886 120 lb / 55 kg



Model 5350 \$ 3775

Shown with optional swivel tray,  
arm mounted 1110-731

Model 5355 \$ 4115

Standard Features	Factory Options	Part No	Retail	Price
<b>Doctor Standard IC Cart 5350</b>	<b>CART .....</b> MODEL CHOICE _____		<b>PRICE:</b>	
Standard IC Control Head	Handpiece Illumination System			
automatic control for three handpieces	One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS	\$ 715		
asepsis 3-way air/water syringe w/tubing, 4'	Add'l handpiece illumination tubing (5-hole)..... 0017-615-CS	275		
handpiece tubing, 4-hole Midwest, asepsis style, 4'	One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS	715		
pneumatic master control block w/diaphragm	Add'l handpiece illumination tubing (6-pin)..... 0017-603-CS	195		
master on/off				
individual handpiece pressure adjustment				
individual water coolant adjustment				
non-retraction water coolant				
air coolant adjustment				
handpiece flush				
handpiece exhaust canister/oil collector				
self-contained clean water system w/quick				
disconnect, single bottle 1 1/2 liter				
wet/dry foot control				
Adjustable height frame				
Twin wheel hard floor casters				
External 7' umbilical				
Standard utility center housing and cover				
Air pressure regulator w/filter				
Master air shut-off valve				
			<b>TOTAL</b>	

► Often overlooked

**Duo Standard IC Cart 5355**

Same as above, w/additional features:

Central vacuum canister w/disposable solids collector  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'

## Simplified Surface Treatments

## Equipment dimensions, see page 111

Shipping weight estimate:

5350 50 lb / 23 kg

5355 55 lb / 25 kg

## UNIT PLASTIC FINISHES:

 Forest Grey

## POWDER COAT PAINT FINISHES:

 Forest Grey

# Factory Options and Parts

	Part No	Retail
--	---------	--------

## Standard Chair Adapters - for 2" Post

3900 fixed accommodates external umbilical	1115-195	\$ 445
3900 fixed sidebox accommodates internal umbilical	1115-530	585
3900 pivot accommodates internal umbilical	1115-528	585

---

## Monitor Posts - 2" Diameter with Oilite Bushing (used with 1300TPM)

30" with dual slot, fixed sidebox mount	1331-998	\$ 220
30" with single slot, fixed (includes post clamp)	1331-997	290

---

## Light Posts - 2" with Light Bushing

35" fixed (includes post clamp)		
Forest	1110-820	\$ 290
Pelton & Crane	1110-822	290
35" fixed with slot (for monitor cord) (includes post clamp)		
Forest	1110-045	290
44" fixed		
Forest	1110-794	220
Pelton & Crane	1110-785	220
44" fixed with slot (for monitor cord)		
Forest	1110-050	220
44" fixed with slot (sidebox access)		
Forest	1110-900	220
Pelton & Crane	1110-890	220
44" fixed with dual slot (for monitor cord & sidebox access)		
Forest	1110-170	220

	Part No	Retail
--	---------	--------



### 250 Watt (Universal Power Supply)

Used when two or more ancillary devices are ordered w/a Forest Operatory Light (excluding the 3900 soft-touch massager power source)

- Provides power for MULTIPLE electrical accessories such as: fiber-optics operatory lights, scalers, curing lights, and electric micro motors
- Size: 7.5" - 5.5" - 2.5". Easily fits in all Forest floor boxes
- Quick disconnect connectors for all of Forest ancillary products
- Available in both 230 volt and 110 volt
- Auto-resettable Klixon switch: inbuilt over current and over temperature protection, eliminating the need for manually changing fuses or circuit breakers

Please note: built in air-electric switch, power source turns on/off w/pilot air only



### Handpiece Illumination System, 5 hole Adapter w/ Krypton Lamp

Factory option


1700-021-CS \$ 715



### Handpiece Illumination System, 6 Pin

Factory option

1700-019-CS 715



- Dual channel outputs individually adjustable from 2.5VDC to 4.0VDC
- Accommodates one additional handpiece tubing
- Illumination delay (adjustable between 0-15 seconds) can be set individually for two handpieces
- Automatic activation of handpiece illumination w/the dental unit foot control
- Simple installation into any new or existing dental unit
- Indicator LED's confirms correct operation of the control module
- Powers all brands of optic handpieces and air-motors
- ISO standard handpiece terminal w/stainless steel nut
- Asepsis style handpiece tubing, 7' long (4' to 5' installed)
- Power transformer, output 24VAC, 1.05 Amps
- Installation hardware and instructions



### Handpiece Illumination Tubing, 5 hole Adapter w/ Krypton Lamp

Factory option

0017-615-CS \$ 275

Used to replace power optic handpiece tubing

- ISO standard power optic handpiece terminal w/stainless steel nut
- 5 hole adapter w/krypton lamp
- Asepsis style handpiece tubing, 7' (4' to 7' factory installed)



### Handpiece Illumination Tubing, 6 pin

Factory option

0017-603-CS \$ 195

Used to replace power optic handpiece tubing

- ISO standard power optic handpiece terminal w/stainless steel adapter and nut
- Asepsis style handpiece tubing, 7' (4' to 7' factory installed)

	Part No	Retail
<b>Standard Deluxe Air/Water Syringe with Quick-release Tip</b>		
Factory option	1109-010	\$ 90
<b>with Sheathed Tubing Assembly, 5'</b>		
Sterling	1107-755	140
Shadow	1107-755-SH	140
Black	1107-755-BLK	140
<b>Straight deluxe Air/Water Syringe with Quick-release Tip</b>		
Factory option	1109-025	110
<b>with Sheathed Tubing Assembly, 4'</b>		
Sterling	1109-585-STG	160
Shadow	1109-585-SHW	160
Black	1109-585-BLK	160
<b>Standard Deluxe One Button Syringe with Quick-release Tip</b>		
Factory option	1109-035	105
<b>with Sheathed Tubing Assembly, 5'</b>		
Sterling	1109-485-STG	180
Shadow	1109-485-SHW	180
Black	1109-485-BLK	180

<b>Additional HVE with Tubing 5' and Autoclavable Oral Dry Cup</b>		
Factory option	1400-113-CS	\$ 90
• Standard vacuum valve		

<b>Dual Articulating Quad Holder Assembly</b>		
Factory option	0024-227-CS	\$ 260
4 standard, 4 universal, White	0024-227	320
4 standard, 4 universal, Forest Grey	0024-227-G	320
4 standard, 4 universal, Designer Friendly® (specify color)	0024-227-C	320
8 universal, White	0024-224	320
8 universal, Forest Grey	0024-224-G	320
8 universal, Designer Friendly® (specify color)	0024-224-C	320

<b>Support Shelf</b>		
Factory option	0023-665	\$ 185
Forest Grey	0023-665-G	185
• Provides adequate support for use of foot control while allowing cart mobility.		

	Part No	Retail
		
<b>Contoured Utility Center Housing and Cover</b>		
Factory option	3900-101	\$ 435
Forest Grey	3900-101-G	435
Shadow	3900-101-SH	435
Designer Friendly® (specify color)	3900-101-C	535
<b>Complete Contoured Utility Center</b>		
Factory option	3900-102	790
Forest Grey	3900-102-G	790
Shadow	3900-102-SH	790
Designer Friendly® (specify color)	3900-102-C	890
<ul style="list-style-type: none"> <li>• Air and water pressure regulator w/filter</li> <li>• Master air and water shut-off valves</li> </ul>		
		
<b>Deluxe Utility Center Housing and Cover</b>		
Forest Grey	1115-007-G	\$ 235
Shadow	1115-007-SH	235
Black	1115-007-BLK	235
White	1115-007-WH	235
<ul style="list-style-type: none"> <li>• Right umbilical access</li> </ul>		
Forest Grey	1115-009-G	235
Shadow	1115-009-SH	235
Black	1115-009-BLK	235
White	1115-009-WH	235
<ul style="list-style-type: none"> <li>• Rear umbilical access</li> </ul>		
		
<b>Complete Deluxe Utility Center</b>		
Forest Grey	1115-005-G	675
Shadow	1115-005-SH	675
Black	1115-005-BLK	675
White	1115-005-WH	675
<ul style="list-style-type: none"> <li>• Right umbilical access</li> <li>• Air and water pressure regulator w/filter</li> <li>• Master air and water shut-off valves</li> </ul>		
		
<b>Standard Utility Center Housing and Cover</b>		
Factory option	1115-020	\$ 145
Forest Grey	1115-020-G	145
Designer Friendly® (specify color)	1115-020-C	145
<b>Complete Standard Utility Center</b>		
Factory option	1115-001	500
Forest Grey	1115-001-G	500
Designer Friendly® (specify color)	1115-001-C	500
<ul style="list-style-type: none"> <li>• Air and water pressure regulator w/ filter</li> <li>• Master air and water shut-off valves</li> </ul>		
<b>Recessed utility center housing and cover</b>	1115-305	\$ 280

Part No	Retail
---------	--------



**Tray Holder Assembly, Use with 1095 Rear Pivot  
Assistant's Instrumentation**

Factory option	1110-680	\$ 575
Forest Grey	1110-680-G	575
Designer Friendly® (specify color)	1110-680-C	575

Includes:

- Pivoting powder coated tray holder, 9 3/4" x 13 1/2"
- Clamp assembly adapts to 1 1/2" horizontal or vertical post



Swivel Bracket

**Slide with Swivel 7 5/8" Extension, use w/ 1284, 2184, 2400, MT's**

Factory option (2400/MT's)	1700-140	\$ 290
Forest Grey (2400/MT's)	1700-140-G	290
Factory option (1284, 2184) (Does not include swivel bracket)	1700-142	290



**Folding Arm Assembly, use with 2400, MT's**

Factory option	1110-115-CS	\$ 345
Forest Grey	1110-115-G	345

Includes:

- 8" folding arm with 1/2" pin and umbilical retainer
- 4" x 4" wall mounting plate which accepts 1/2" pin
- 12" folding arm with 1/2" pin and umbilical retainer
- Control mount with 1/2" pin



**Universal Heavy Duty Mounting Plate**

10" x 24" beneath counter	0010-084-BLK	\$ 330
11" x 14 1/2" beneath counter	0010-086-BLK	200
10" x 24" in lieu of 11" x 14 1/2"	0010-084-CS	130



**26" Rigid with 24" Flexible Arm In Lieu of Top Post Mount**

Factory Option	1110-771-RM	\$ 100
----------------	-------------	--------

**22" Rigid with 18" Flexible Arm In Lieu of Top Post Mount,**

Factory Option	1110-536-RM	100
----------------	-------------	-----

**26" Rigid With 18" Flexible Arm In Lieu Of Top Post Mount,**

Factory Option	1110-537-RM	100
----------------	-------------	-----

- Designed for use with 2" diameter post.



**Floating Work Surface Arm Assembly, 20"**

Forest Grey	0023-815-G	\$ 755
Designer Friendly® (specify color)	0023-815-C	755





	Part No	Retail
<b>Quick Disconnects</b>		
<b>1/2" Female and Male QD with Flow Control</b>		
Factory Installed	0016-454-CS	\$ 145
<b>1/4" Female QD with flow control, black/grey/shadow/white</b>		
Factory Installed	0016-454-BLK /-G /-SH /-WH	115
<b>1/4" Female and Male QD With Flow Control, Arm and Control Head Mount,</b>		
Factory Installed	0016-455-CS	145
<b>1/4" Female QD with Flow Control, Arm and Control Head Mount, Black /Grey</b>		
Factory Installed	0016-455-BLK /-G	115
<b>1/4" Female and Male QD With Flow Control, Factory Installed</b>		
Factory Installed	0016-456-CS	145
<b>1/4" Female QD with Flow Control, Black /Grey/Shadow/White</b>		
Factory Installed	0016-456-BLK /-G /-SH /-WH	115
<b>Routing Valve Kits</b>		
— Check with Local Authorities About Installation Requirements		
Factory installed with water shut-off and regulator	1115-315-CS	\$ 305
Factory installed without water shut-off and regulator	1115-316-CS	145
<b>Dealer Installed Includes;</b>		
<b>Toggle Valve/Bracket/Tubing/Water Shut-Off/Regulator</b>		
black / grey / shadow / white	1115-310-BLK / -G / -SH / -WH	250
<b>Dealer Installed Includes; Toggle Valve, Bracket, And Tubing</b>		
black / grey / shadow / white	1115-312-BLK / -G / -SH / -WH	70

	<b>High Density Bottle</b>		
1100-431	1 1/2" liter with quick disconnect adapter	1100-431	\$ 70
1100-423	1 1/2" liter with quick disconnect upgrade kit	1100-423	105
	For use with clean water system		
	<ul style="list-style-type: none"> <li>• Opaque bottle</li> <li>• High Density Polyethylene (HDPE)</li> </ul>		

# General Information





## Dimensions




### Pivot Delivery Systems



Text below the drawing: "From tip of assistant's holder to end of pivot light post is 111"

### 1500PI Rear Pivot Delivery




**Fixed with Sidebox Delivery Systems****Fixed without Sidebox Delivery Systems****1636/1637 Basic**


# General Information

## Dimensions


### 3900 Chair




### 6176/6198 Doctor's Office Comfort Stool




### 6159C Doctor's Deluxe Stool




### 6156C Doctor's Economy Stool




Deluxe      Standard      Recessed




### 6187 Assistant's Comfort Stool






### 6160C/6161C Assistant's Deluxe Stool



### 6157C/6158C Assistant's Economy Stool




**Chair-Mounted Light****Ceiling-Mounted Light****Track Mount Light/Monitor**


# General Information

## Dimensions


### Standard Wall & Cabinet Mounted Light




### Tall Wall & Cabinet Mounted Light




### Cabinet Mount





### Optional Cabinet Mount



### 1300RM





### 1300TPM





**7020/7022PRO-BC Beneath Counter**




\*7022 does not include doctor's delivery.

**7021/7023PRO-BC Beneath Counter**

\*7023 does not include doctor's delivery.



- Large contoured work surface standard for 7020/7021/7022/7023
- Work surface 9/16" thick
- Shown with organization features
- Height adjustable work surface, 5"





**7020/7021/7022/7023PRO-BC Mounting Plates**

# General Information

## Dimensions



### 7020/7022PRO-WM Wall Mount

\*7022 does not include doctor's delivery.




### 7021/7023PRO-WM Wall Mount


\*7023 does not include doctor's delivery.




- Large contoured work surface standard for 7020/7021/7022/7023
- Work surface 9/16" thick
- Shown with organization features
- Height adjustable work surface, 5"



**7020/7021/7022/7023 WM  
Mounting Plate**





**7020/7021/7022/7023 WM  
Leveling Mounting Board**





**7020/7022PRO-UC Under Cabinet Floor Mount**




\*7022 does not include doctor's delivery.

**7021/7023PRO-UC Under Cabinet Floor Mount**

\*7023 does not include doctor's delivery.



- Large contoured work surface standard for 7020/7021/7022/7023
- Work surface 9/16" thick
- Shown with organization features
- Height adjustable work surface, 5"



**7020/7021/7022/7023PRO-UC Under Cabinet Floor Mount Leveling/Mounting Plate****1300-025-CS Monitor Assembly Option**

# General Information

## Dimensions



### 7020/7022PRO-MC Duo Cart

\*7023 does not include doctor's delivery.




### 7021/7023PRO-MC Duo Cart


\*7023 does not include doctor's delivery.




- Large contoured work surface standard for 7020/7021/7022/7023
- Work surface 9/16" thick
- Shown with organization features
- Height adjustable work surface, 5"












### 1300-025-CS Monitor Assembly Option



**7010BC Beneath Counter**


- Work surface 9/16" thick (semi-round), 7010
- Shown with optional storage organizer and swivel tray
- Shown with additional optional asst's quad holder assembly
- Height adjustable work surface, 5"

**7010BC Mounting Plates****0010-086****7010WM Mounting Plate****7010WM Wall Mount****0010-084-CS Option****7010WM Leveling/Mounting Board****1300-025-CS Monitor Assembly Option**


# General Information

## Dimensions


- Work surface 9/16" thick (semi-round), 7010
- Shown with optional storage organizer and swivel tray
- Shown with additional optional asst's quad holder assembly
- Height adjustable work surface, 5"




## 1300-025-CS Monitor Assembly Option




## 7010UC Under Cabinet Floor Mount











## 7010UC Under Cabinet Floor Mount Leveling/Mounting Plate




## 7010MC Duo Cart



**3486BC and WM Rear Beneath Counter Mount****3486BC Heavy Duty Mounting Plate****3486WM Wall Mount Hub****3486WM Leveling/Mounting Board****3486BC and WM Mounting Arm**


- Semi-round work surface standard for 2184BC/2186BC/3486/5386/5886
- Work surface 9/16" thick (semi-round and rectangular)
- Shown with optional storage organizer and swivel tray
- Shown with additional optional asst's quad holder assembly




# General Information

## Dimensions


### 3486UC Under Cabinet Floor Mount Assistant's




- Semi-round work surface standard for 2184BC/2186BC/3486/5386/5886
- Work surface 9/16" thick (semi-round and rectangular)
- Shown with optional storage organizer and swivel tray
- Shown with additional optional asst's quad holder assembly




### 2184 Rear Beneath Work Surface Duo Unit




#### 2184 MOUNTING PLATE




### 3486UC Under Cabinet Floor Mount Leveling/Mounting Plate




### 1300-025-CS Monitor Assembly Option











### 2184/2500 Wall Cut-Out



### 1700-142 SLIDE MOUNT OPTION, 1284/2184




**2186 Rear Beneath Work Surface Duo Unit****2186 Control Mounting****2184BC Beneath Counter Duo Unit****2184BC/2186BC Beneath Counter Duo Unit****2186BC Rear Beneath Counter Duo Unit****2500 PANEL MOUNT****2520 PANEL MOUNT****2520 Wall Cut-Out**

# General Information


## Dimensions

### 3235 Cabinet Side Delivery




Minimum cabinet opening 18" deep x 30" wide


### 3230 Wall Side Delivery




### 2400 Rear Mount Beneath Counter




### 2400 Mounting Plate




### 3230/3235 Leveling Mounting Plate




### MT600/MT610 Compact Controls











### MT600/610BS Leveling/Mounting Plate



### 1700-140 Slide Mount With Swivel Option



**1284 Rear, Beneath Work Surface****1286 Rear, Beneath Shelf****1700-142 Slide Mount Option, 1284/2184****1284/1286 Mounting Plate****5386 Assistant's Cart****5886 Duo Cart****5350/5355 Standard IC Cart****5350/5355 Work Surface**

## General Information

### Equipment Maintenance

Proper maintenance will prolong the life of the dental equipment. It is important to adhere to routine practices for patient safety and avoid premature aging of the dental equipment. Following the step by step details outlined in the Forest Operator Guides included with the product will provide valuable information to achieve optimum success maintaining patient safety and the lifespan of the treatment room equipment.

We recommend barrier protection be used for surfaces that may be contacted by dental practitioners or instruments during a procedure. In addition to infection control, barrier protectors also preserve the finish, appearance and life of the dental equipment.

The water reservoir system installed on the Forest dental unit does not control bacteria levels in dental unit water line "DUWL" or the reservoir bottle therefore bacteria levels must be maintained and monitored by the dental team to prevent biofilm build-up in procedural water quality. Isolated water sources (as opposed to city water hook-up) are highly recommended to allow full control of the source and quality of water circulating through the dental unit water line.

The Centers for Disease Control and Prevention's (CDC) recommendations for the cleaning, disinfection, and sterilization of dental equipment can be found in:

Guidelines for Infection Control in Dental HealthCare Settings — 2003  
[www.cdc.gov/mmwr/PDF/rr/rr5217.pdf](http://www.cdc.gov/mmwr/PDF/rr/rr5217.pdf)

Guidelines for Disinfection and Sterilization in Healthcare Facilities – 2008  
[www.cdc.gov/infectioncontrol/pdf/guidelines/disinfection-guidelines.pdf](http://www.cdc.gov/infectioncontrol/pdf/guidelines/disinfection-guidelines.pdf)

Summary of Infection Prevention Practices in Dental Settings – 2016  
<https://www.cdc.gov/oralhealth/infectioncontrol/pdf/safe-care2.pdf>

The American Dental Association's recommendations can be found at: [www.ada.org](http://www.ada.org)

The Organization for Safety and Asepsis Procedures recommendations can be found at: [www.osap.org](http://www.osap.org)

FDA Guidelines for Dental Unit Water Line Maintenance can be found at: [www.fda.gov](http://www.fda.gov)

## Terms

All prices unless otherwise listed are F.O.B. Hillsboro, Oregon U.S.A. Payments are net 30 days of invoice date. A service charge of 1.5% per month (annual rate 18%) will be charged on past due accounts, \$1.00 minimum. All prices are subject to change without notice.

## Return Merchandise

Authorized Forest dealers wishing to return merchandise for credit consideration must:

Obtain prior factory authorization.

Provide basic reason for return.

Provide serial number, purchase order or invoice number of original order.

Provide a purchase order for warranty items when necessary. Credit will be issued by Forest upon determination of warranty coverage in accordance with the stated Warranty Policy.

Return product within a period of 90 days—if product is not returned within time period, authorization becomes invalid.

Forest will approve return and credit amount after evaluation of product. If product is used or modified, the reason for return will be used to determine if replacement product or credit will be allowed.

### Handling Charges:

25% or \$15.00 minimum restock fee will be charged on returned equipment approved for credit.

10% or \$15.00 minimum restock fee will be charged on new and unopened parts/components.

Returned items must be shipped freight prepaid.

## Warranty

Forest warrants its products to be free from defects in material and workmanship only. No other warranties are expressed or implied. This warranty shall extend for five years unless otherwise stated. Written notice of breach must be given to Forest within this period. Buyer's remedy for breach of this warranty is limited to repair parts or replacement of items by Forest. Factory repairs are covered for a period of 90 days. No claim for labor or consequential damages or freight will be allowed. Forest is not responsible for ensuring equipment or installation conformity to local code or regulatory requirements. Check with local code or regulatory requirements before purchase and installation of all Forest product. Field modifications that alter the electrical and /or mechanical safety of Forest products conflict with agency construction file requirements and are not sanctioned by Forest.

### Warranty is void if:

Items are carelessly used or improperly maintained or installed.

Damage results from using cleaning, disinfecting or sterilization chemicals and processes.

Damage is the result of freight damage and as such claims must be filed with the carrier by dealer.

### Warranty term exceptions cover items that are normally expected service items and items not covered by Forest's warranty:

Handpiece illumination tubings, heated syringe tubings & assemblies and light bulbs are warranted for six months.

Power supplies, syringe buttons, HVE & SE valves are warranted for one year.

Water bottles, gas springs and handpiece tubings are warranted for two years.

Air and Water tubing connectors (T's), Distribution Block, Master Control Block, Flow Valves, 7000 Chassis are warranted for a maximum of 20 years.

Upholstery warranty, excluding wear and tear, for two years.

## Cancellation Policy

### Standard orders:

- Fax cancellation request required 10 working days prior to promised ship date, otherwise a 25% restock fee will be charged.
- Expedited ship date requests are non-cancellable.

## Custom Orders

- Custom configurations are those that deviate from Forest's published product offering:

Orders are non-cancellable and non-returnable.


May be subject to an additional "Deviation" charge and/or a non-refundable engineering charge.

- Custom powder coat, custom plastic surface treatments and/or custom treatments such as multiple color combinations or custom colors may result in an upcharge of 10-15% per model and are non-cancellable and non-returnable.

## Notes

## General Information

Professional Sales Associates (PSA)



**Steve Carlin**  
952.491.1223  
[scarlin@profsales.com](mailto:scarlin@profsales.com)

**John Miller**  
503.758.9895  
[jmiller@profsales.com](mailto:jmiller@profsales.com)

 **Skip Ellender**  
617.894.1224  
[sellender@profsales.com](mailto:sellender@profsales.com)

**Sam Paola**  
310.569.1178  
[spaola@profsales.com](mailto:spaola@profsales.com)

**Jim Schlosser**  
909.215.3277  
[ischlosser@profsales.com](mailto:ischlosser@profsales.com)

**Rick Orme**  
916.837.3839  
[rorme@profsales.com](mailto:rorme@profsales.com)

**Victor Yoon**  
778.688.8429  
[vyoon@profsales.com](mailto:vyoon@profsales.com)

**Scott Rumble**  
403.801.5390  
[srumble@profsales.com](mailto:srumble@profsales.com)

 **Bill Highland**  
810.623.7297  
[bhighland@profsales.com](mailto:bhighland@profsales.com)

**Charlie Clarke**  
859.509.6298  
[cclarke@profsales.com](mailto:cclarke@profsales.com)

 **Matt Highland**  
901.606.3633  
[mhighland@profsales.com](mailto:mhighland@profsales.com)

 Michael Park  
281.702.3992  
mpark@profsales.com

 Stephan Nugent  
214-394-0074  
[snugent@profsales.com](mailto:snugent@profsales.com)

 **Pat Moylan**  
816.729.2900  
[pmoylan@profsales.com](mailto:pmoylan@profsales.com)

**Brett Hobbs**  
630.772.4321  
[bhobbs@profsales.com](mailto:bhobbs@profsales.com)

**Aaron Renner**  
519.319.7457  
[arenner@profsales.com](mailto:arenner@profsales.com)

 **Michel Cabanie**  
514.793.3368  
[mcabanie@profsales.com](mailto:mcabanie@profsales.com)

 Jim Domitrz  
617.921.4580  
[jdomitrz@profsales.com](mailto:jdomitrz@profsales.com)

 Matt Langdon  
516.316.5913  
[mlangdon@profsales.com](mailto:mlangdon@profsales.com)

**Milton Lorring**  
216.396.2341  
[mlorring@profsales.com](mailto:mlorring@profsales.com)

**Don Rogers, III**  
717.201.7332  
[drogers@profsales.com](mailto:drogers@profsales.com)

 **Rich Ellis**  
301.502.7424  
[rellis@profsales.com](mailto:rellis@profsales.com)

336.210.1429  
koberry@profsales.com

205.529.0243  
mgullage@profsales.com



Raven Wing



Charcoal



Cave



Stone



Pearl Pewter



Papyrus



Doe



Pearl Abalone



Baltic



Dove Grey



Riviera



Pearl Midnight



Diplomat Blue



Parrot



Eucalyptus



Red



Apricot


Ultraleather®



Laser Leather




Lemans




Greywood

Carbon Fiber



White



Forest Grey

Plastics




Forest Grey



Shadow



White



Arctic Silver



Pewter



Midnight



Copper




Charcoal




Bronze

Powder Coat

Formica®



Surf



Folkstone




Rocky Road Granite




Silver Mist



Platinum Granite



Oyster Grey



White Blanc



Moonstone



Sanibel Granite




Goldrush Gemstone

Gibraltar®



Arctic Melange



Baja Melange

Corian®



Antarctica



Cinder Granite



Starry Starry Night



---

Phone

- 800 . 423 . 3555
  - 503 . 640 . 3012
- 

Fax

- 503 . 693 . 9715
- 

Digital

- [sales@forestdental.com](mailto:sales@forestdental.com)
  - [forestdental.com](http://forestdental.com)
- 

Address

- 6200 ne cherry drive hillsboro oregon 97124

Ultraleather® is a registered trademark of Ultrafabrics, Inc.®

Naugahyde® is a registered trademark of Uniroyal Engineered Products, LLC®

Simplicity™ is a trademark of ©DentalEZ Inc.



©DTE 2019 Oregon, Inc.

98-2019-EDC .Rev-