



assistina



Weekly checklist




**Check liquid level!
Refill if necessary!**

Service oil:
500 ml – 02675930

Maintenance solution:
1000 ml – 02680230





Check function!



Check o-rings!


**Look for right position
and / or
damage!**





General hints



**When maintaining after every patient – Press program button once!
When maintaining twice a day – Press program button twice!**




Check connection is tight!



**After consumption of one
500 ml – bottle of Service oil
[approx. 7000 cycles]**

Changing of filters!



Once a year

Manufacturer:
W&H Dentalwerk Bürmoos Ges.m.b.H.
A-5111 Bürmoos, Postfach 1
Tel.: +43/6274/6236-0, Fax +43/6274/6236-55

Distributor:
A-dec Inc.
2601 Crestview Drive, Newberg, Oregon 97132
Phone: (800) 547-1883, Fax (503) 537-2752