



COMPLETE NITROUS OXIDE/OXYGEN SEDATION SOLUTIONS

2019 EQUIPMENT AND ACCESSORIES CATALOG



CROSSTEX
A CANTEL MEDICAL COMPANY



All Accutron products are made in the U.S.A. We take pride in providing innovative, high-quality products for over 45 years. Safety for patients, clinicians, dental office staff, and installation experts is our first concern.

The complete Accutron product line is available through authorized dental distributors.



Products

Flowmeter Systems & Flowmeter Components

Digital Ultra™ Flushmount Flowmeter & Accessories	4-5
Digital Ultra™ Flowmeter	6-7
Digital Newport™ Flowmeter System	7
Ultra PC™ % Cabinet Mount Flowmeter & Accessories	8-9
Newport™ Flowmeter System	9
Ultra PC™ % Flowmeter	10-11
Ultra DC™ Flowmeter	12-13
Flowmeters without Pre-Install Kits & Pre-Install Kit Components	14
Oral Surgery Flowmeter Systems	14-15
Oral Surgery Accessories	16
Oral Surgery Suction Accessories	16

Vacuum Controllers & Mounts

Vacuum Controllers	17
Universal Cabinet Slide Mounts	17
Flowmeter Mounts & Flowmeter Mounting Adapters	18
RFS™ Remote Flow System Chairmount Kits	19-21

Central System Schematic & Manifold Systems

Central System Schematic (Gas Cylinder Room)	22-23
Digi-Flo™ Digital Automatic Switching Manifold Systems	24-25
Guardian Monitor™ Analog Manifold Systems	25-26
Manifold Accessories & Components	26
Portable System Manifolds & Components	27

Hoses & Outlets

Hoses & Vacuum Connectors	28
Outlets	29

Oxygen Equipment, Components & Accessories

Emergency Oxygen Equipment	30-31
----------------------------------	-------

Nasal Masks, Scavenging Systems, Rubber Goods & Scavenging Circuit Components

Axess™ Low Profile Nasal Masks	32
Axess™ Capnography Adapters	32
Axess™ Scavenging Circuits	33
ClearView™ Single-Use Nasal Masks	34
ClearView™ Capnography Single-Use Nasal Masks	34
ClearView™ Scavenging Circuits	35
ClearView™ Adaptor Easy Order Guide	35
PIP+™ Single-Use Nasal Masks	36
PIP+™ Multi-Use (steam-sterilizable) Masks	36
PIP+™ Scavenging Circuits	37
PIP+™ Adaptor Easy Order Guide	37
Rubber Goods	38
Scavenging Circuit Components	38

Diagrams, Templates & Service/Warranty>Returns

Scavenging Circuit Diagrams	39
Flowmeter Diagrams & Templates	40-42
Outlet Templates	43
Service/Warranty>Returns	44
Notes	45

Toll-free: 800.531.2221

Local Phone: 623.780.2020 • Fax: 623.780.0444

accutron-inc.com



1 Digital Ultra™ Flushmount Flowmeter

- Sealed flush surface makes it easy to disinfect or barrier-protect
- N₂O Ratio Control with 70% Max
- Dual Readouts for Each Gas
- Audible & Visual Low O₂ Alarms
- Oxygen Fail-safe; Oxygen Flush
- Oxygen Resuscitation Connector
- Directional Check Valve
- Emergency Air Valve
- Space-saving Low Profile Unit
- Mounting Options

Bag Tee &
Remote BlockRFS™ Remote
Flow System DRFS™ Remote
Flow System AAccu-Vac™
Vacuum ControllerScavenging Circuit II
(Standard)RFS™ System D
& Accu-Vac™
Vacuum Controller

**RFS™ Chairmount Kits
for 56800-CHAIR
on pages 19-21**

Product Specifications

Dimensions (W x H x D) Flushmount Control Module:	9 1/4" x 4" x 1 1/2" 23,5 cm x 10,8 cm x 3,8 cm
Dimensions (W x H x D) Remote Main Module:	8" x 4" x 7 1/2" 18,4 cm x 8,9 cm x 19,1 cm
Mounted Dimensions (W x H x D):	9 1/2" x 4 1/4" x 1 1/4"
Weight Flushmount Control Module:	1.3 lbs (0.6 kg)
Weight Remote Main Module:	6.1 lbs (2.8 kg)
Mixed Gas Outlet Diameter:	7/8" Nominal (22mm)
Oxygen Inlet:	Male DISS (CGA 1240)
Nitrous Oxide Inlet:	Male DISS (CGA 1040)
Nitrous Oxide Mixed Ratio Range:	0-70%
Oxygen Flow Range:	1.0-9.9 LPM
Nitrous Oxide Flow Range:	0-6.9 LPM
Inlet Pressure (each gas):	50-55 psi
Resuscitator Connector:	Quick Connect (1/4" Imperial) at Remote Outlet
Resuscitator Connector Flow:	100 LPM Oxygen
Power Requirements:	100-240v AC, 50-60 Hz, 30 Watts

Digital % Flushmount

Central Systems

Order No.	Product Description	Retail List
56600	Digital Ultra™ Flushmount Package G Digital Ultra™ Flushmount Flowmeter, Mounting Hardware, Check Valve Risers, Bag Tee, Scavenging Circuit II (standard) (Includes 1 2 6 8)	\$6,152
56800	Digital Ultra™ Flushmount Package H Digital Ultra™ Flushmount Flowmeter, Mounting Hardware, Check Valve Risers, RFS™ Remote Flow System D (with O ₂ port), RFS™ Flush Panel Mount Faceplate, Scavenging Circuit II (for RFS™-D) (Includes 1 3 7 8)	\$6,447
56800-CHAIR	Digital Ultra™ Flushmount Package H CHAIR Digital Ultra™ Flushmount Flowmeter, Mounting Hardware, RFS™ Remote Flow System A, Oxygen Outlet, O ₂ and N ₂ O Check Valve Risers, 25' of Vacuum Hose, 25' of Mixed Gas Hose, One 6' N ₂ O Hose, Two 6' O ₂ Hoses, and Scavenging Circuit II (for RFS™-A) (Includes 1 4 7 8 9) <i>Requires separate purchase of appropriate chair brand Chairmount Kit – see pages 19-21 for Chairmount kits</i>	\$6,732
56950	Digital Ultra™ Flushmount Package I Digital Ultra™ Flushmount Flowmeter, Mounting Hardware, Check Valve Risers, Accu-Vac™ Vacuum Controller, Accu-Vac™ Vacuum Controller Flush Panel Mount Faceplate, Scavenging Circuit II (for Accu-Vac™ Vacuum Controller) (Includes 1 5 7 8)	\$6,830

To order Flushmount packages without pre-install kits and flowmeter pre-install kit components see page 14. Scavenging Circuit II includes a complimentary sampling of single-use nasal masks.



Flowmeter not included.

Digital Ultra™ Flushmount Flowmeter Accessories

Order No.	Product Description	Retail List
23173-FRU	O ₂ Check Valve Riser - DISS (for operatory); includes 28085	\$148
23174-FRU	N ₂ O Check Valve Riser - DISS (for operatory); includes 28086	\$148
25043	Universal Light Pole Mount Kit for RFS™ - D	\$184
27193	Bracket for Digital Ultra™ Flushmount Flowmeter (use for mounting flowmeter in wall)	\$134
27793-FRU	Oxygen Outlet (for use with Digital Ultra™ Flushmount, RFS™ System A & RFS™ Chairmount Kits)	\$274
27328	Main Module Mounting Bracket	\$72
28085-FRU	O ₂ DISS x DISS Demand Check Valve Assembly	\$74
28086-FRU	N ₂ O DISS x DISS Demand Check Valve Assembly	\$74
32015	Universal Cabinet Slide, Accu-Vac™ Vacuum Controller (slide mounts to left, right, or under cabinet)	\$544
32018	Universal Cabinet Slide, Digital RFS™ Remote Flow System (slide mounts to left, right, or under cabinet)	\$544
32019	Universal Cabinet Slide, Digital Bag Tee (slide mounts to left, right, or under cabinet)	\$544
32055	Surface Mount for Digital Ultra™ Flushmount (flowmeter not included)	\$322
54500	Remote Control for Digital Ultra™ Flowmeter	\$492
54600	Printer for Digital Ultra™ Flowmeter	\$658
54650	Printer Paper for Digital Ultra Printer (10 rolls)	\$64





1 Digital Ultra™ Flowmeter

- Sealed flush surface makes it easy to disinfect or barrier-protect
- N₂O Ratio Control with 70% Max
- Dual Readouts for Each Gas
- Audible & Visual Low O₂ Alarms
- Oxygen Fail-safe; Oxygen Flush
- Oxygen Resuscitation Connector
- Directional Check Valve
- Emergency Air Valve

Product Specifications

Dimensions (W x H x D):	7 1/2" x 7" x 8 3/4" 19,1 cm x 17,8 cm x 22,2 cm
Weight:	6.3 lbs
Mixed Gas Outlet Diameter:	7/8" Nominal (22mm)
Oxygen Inlet:	Male DISS (CGA 1240)
Nitrous Oxide Inlet:	Male DISS (CGA 1040)
Nitrous Oxide Mixed Ratio Range:	0-70%
Oxygen Flow Range:	1.0-9.9 LPM
Nitrous Oxide Flow Range:	0-6.9 LPM
Inlet Pressure (each gas):	50-55 psi
Resuscitator Connector:	Quick Connect (1/4" Imperial)
Resuscitator Connector Flow:	100 LPM Oxygen
Power Requirements:	100-240v AC, 50-60 Hz, 30 Watts

Digital % Flowmeter



4-Cylinder Portable System

Height with Flowmeter:
36 3/8"
W x D
(with cylinders):
22" x 22"



Scavenging Circuit II (Standard)


Flowmeter Only

Order No.	Product Description	Retail List
51975	Digital Ultra™ Flowmeter (Includes 1)	\$3,902

Portable Systems

Order No.	Product Description	Retail List
56100	Digital Ultra™ Flowmeter Package B Digital Ultra™ Flowmeter 4-Cylinder Portable System Scavenging Circuit II (Includes 1 2 3)	\$7,487

See Digital Newport™ Flowmeter System (4-cylinder portable) on page 7.
Scavenging Circuit II includes a complimentary sampling of single-use nasal masks.



Telescoping Wall Arm
(Extends 12"-17")

Central Systems

Order No.	Product Description	Retail List
56200	Digital Ultra™ Flowmeter Package C Digital Ultra™ Flowmeter Mobile Stand & Scavenging Circuit II (Includes 1 3 4)	\$4,940
56200-ADJS	Digital Ultra™ Flowmeter Package C ADJS Digital Ultra™ Flowmeter Adjustable Mobile Stand – Short Scavenger Circuit II (Includes 1 3 5)	\$ 4,940
56300	Digital Ultra™ Flowmeter Package D Digital Ultra™ Flowmeter 9-inch Wall Arm & Mounting Adapter 24552 Scavenging Circuit II (Includes 1 3 6)	\$4,940
56400	Digital Ultra™ Flowmeter Package E Digital Ultra™ Flowmeter 16-inch Bi-fold Wall Arm & Mounting Adapter 24552 Scavenging Circuit II (Includes 1 3 7)	\$4,940
56500	Digital Ultra™ Flowmeter Package F Digital Ultra™ Flowmeter Telescoping Wall Arm & Mounting Adapter 24552 Scavenging Circuit II (Includes 1 3 8)	\$4,940

See page 28 for a complete listing of Quick Connect Hoses and Vacuum Male Connectors.

Digital Portable

Digital Newport™ Flowmeter System

The **Digital Newport™ Flowmeter System** has all of the operating features of the Digital Ultra™ Flowmeter, plus its enclosed 4-cylinder design makes it a user-friendly portable that can easily be moved from operatory to operatory. The unit's height allows it to be stored under counter space that is open below; if not stored its flat surface can function as additional counter space.



Order No.	Product Description	Retail List
51000	Newport™ Digital Flowmeter System 4-Cylinder Enclosed Portable Scavenging Circuit II (Includes 1 2)	\$9,103



Scavenging Circuit II (Standard)

Accu-Vac™ Vacuum Controller and RFS™ Remote Flow System not available with Newport™ Digital Flowmeter System

Ultra PC™ % Cabinet Mount Flowmeter & Accessories

Ultra PC™ % Cabinet Mount Flowmeter & Accessories



1 Ultra PC™ % Cabinet Mount Flowmeter

- Single-knob Total Flow Control
- N₂O Ratio Control with 70% Max
- Oxygen Fail-safe; Oxygen Flush
- Oxygen Resuscitation Connector
- Emergency Air Valve
- Directional Check Valve
- Space-saving Low Profile Unit
- Continuous O₂ & N₂O Flow Tubes
- Neutral Color

Product Specifications

Dimensions (W x H x D):	9 3/4" x 4 3/4" x 3 1/2" 24,8 cm x 12,1 cm x 8,9 cm
Mounted Dimensions (W x H x D):	9 3/4" x 4 3/4" x 1 1/4" 24,8 cm x 12,1 cm x 3,2 cm (to end of knobs) 9 3/4" x 4 3/4" x 1/4" 24,8 cm x 12,1 cm x 0,6 cm (front plate only)
Weight:	6 lb (2.7 kg)
Mixed Gas Outlet Diameter:	7/8" nominal (22 mm)
Oxygen Inlet:	Male DISS (CGA 1240)
Nitrous Oxide Inlet:	Male DISS (CGA 1040)
Nitrous Oxide Mix Ratio Range:	0-70%
Oxygen Flow Range:	1-10 LPM
Nitrous Oxide Flow Range:	0-7 LPM
Inlet Pressure (each gas):	50-55 PSI
Resuscitator Connector:	Quick Connect (1/4" Imperial)
Resuscitator Connector Flow:	100 LPM Oxygen

RFS™
Chairmount Kits
for 36800-CHAIR
on pages 19-21

Analog % Flushmount

Central Systems



26679

Bag Tee



36855

RFS™ Remote Flow System A



43003

Scavenging Circuit II (Standard)



43004

Scavenging Circuit II
(RFS™ System A)

Order No.

Product Description

Retail List

36600

Ultra PC™ % Cabinet Mount Flowmeter Package G

\$5,300

36800

Ultra PC™ % Cabinet Mount Flowmeter Package H

\$5,724

36800-CHAIR

Ultra PC™ % Cabinet Mount Flowmeter Pkg. H CHAIR

\$6,012

To order above Ultra PC™ % Cabinet Mount packages without pre-install
kits and flowmeter pre-install kit components see page 14.

Scavenging Circuit II includes a complimentary sampling of single-use nasal masks.



23173-FRU
23174-FRU


O₂ & N₂O
Check Valve Risers

28085
28086

O₂ DISS & N₂O DISS x DISS Demand
Check Valve Assembly

Universal Cabinet Slides

3 Cabinet Mounting Positions



32008
Universal Cabinet
Slide, Analog Bag Tee



32014
Universal Cabinet Slide,
Remote Flow System (RFS™ System A)

Analog Portable


Newport™ Flowmeter System

The **Newport™ Flowmeter System** has all of the operating features of the Ultra PC™ % Flowmeter, plus its enclosed 4-cylinder design makes it a user-friendly portable that can easily be moved from operatory to operatory. The unit's height allows it to be stored under counter space that is open below; if not stored its flat surface can function as additional counter space.

1 Newport™ Flowmeter System



43003
Scavenging Circuit II (Standard)



50000
Newport™ Flowmeter System
4-Cylinder Enclosed Portable
Scavenging Circuit II
(Includes 1 2)

Ultra PC™ % Cabinet Mount Accessories

Order No.	Product Description	Retail List
23173-FRU	O ₂ Check Valve Riser - DISS includes 28085	\$148
23174-FRU	N ₂ O Check Valve Riser - DISS includes 28086	\$148
28085-FRU	O ₂ DISS x DISS Demand Check Valve Assembly	\$74
28086-FRU	N ₂ O DISS x DISS Demand Check Valve Assembly	\$74
25043	Universal Light Pole Mount Kit for Remote Flow System (RFS™ System A)	\$184
32008	Universal Cabinet Slide, Analog Bag Tee (slide mounts to left, right, or under cabinet)	\$544
32014	Universal Cabinet Slide, Analog RFS™ System (slide mounts to left, right, or under cabinet)	\$544

Central Systems


Order No.	Product Description	Retail List
36200	Ultra PC™ % Flowmeter Package C Ultra PC™ % Flowmeter Mobile Stand Scavenging Circuit II (Includes 1 4 5)	\$4,216
36200-ADJS	Ultra PC™ % Flowmeter Package C ADJS Ultra PC™ % Flowmeter Adjustable Mobile Stand – Short Scavenger Circuit II (Includes 1 4 6)	\$4,216
36300	Ultra PC™ % Flowmeter Package D Ultra PC™ % Flowmeter 9-inch Wall Arm & Mounting Adapter 24552 Scavenging Circuit II (Includes 1 4 7)	\$4,216
36400	Ultra PC™ % Flowmeter Package E Ultra PC™ % Flowmeter 16-inch Bi-fold Wall Arm & Mounting Adapter 24552 Scavenging Circuit II (Includes 1 4 8)	\$4,216
36500	Ultra PC™ % Flowmeter Package F Ultra PC™ % Flowmeter Telescoping Wall Arm & Mounting Adapter 24552 Scavenging Circuit II (Includes 1 4 9)	\$4,216
36500-US	Ultra PC™ % Flowmeter Package U Ultra PC™ % Flowmeter Universal Cabinet Slide, PC/DC (slide mounts to left, right, or under cabinet), and Mounting Adapter 24552 Scavenging Circuit II (Includes 1 4 10)	\$4,216

Order quick connect hoses for Packages C - U separately (page 28)



Universal Cabinet Slides, Ultra PC™ % Flowmeters & Ultra DC™ Flowmeters

3 Cabinet Mounting Positions



Analog % Flowmeter

Flowmeter Only

Order No.	Product Description	Retail List
31975	Ultra PC™ % Flowmeter (Includes 1)	\$3,048

Portable Systems

Order No.	Product Description	Retail List
36000	Ultra PC™ % Flowmeter Package A Ultra PC™ % Flowmeter 2-Cylinder Portable System Scavenging Circuit II (Includes 1 2 4)	\$6,034
36100	Ultra PC™ % Flowmeter Package B Ultra PC™ % Flowmeter 4-Cylinder Portable System Scavenging Circuit II (Includes 1 3 4)	\$6,646

See Newport™ Flowmeter System (4-cylinder portable) on page 9.

Scavenging Circuit II includes a complimentary sampling of single-use nasal masks.



- 1 Ultra PC™ % Flowmeter**
- Single-knob Total Flow Control
 - N₂O Ratio Control with 70% Max
 - Oxygen Fail-safe; Oxygen Flush
 - Oxygen Resuscitation Connector
 - Emergency Air Valve
 - Directional Check Valve
 - DISS Inlet Fittings
 - Neutral Color

Product Specifications

Dimensions (W x H x D):	4 1/4" x 10 1/2" x 8" 11 cm x 27 cm x 20 cm
Weight:	4.5 lb (2.0 kg)
Mixed Gas Outlet Diameter:	7/8" nominal (22 mm)
Oxygen Inlet:	Male DISS (CGA-1240)
Nitrous Oxide Inlet:	Male DISS (CGA-1040)
Nitrous Oxide Mix Ratio Range:	0-70%
Oxygen Flow Range:	1-10 LPM
Nitrous Oxide Flow Range:	0-7 LPM
Inlet Pressure (each gas):	50-55 psi
Resuscitator Connector:	Quick Connect (1/4" Imperial)
Resuscitator Connector Flow:	100 LPM Oxygen



2

Height with Flowmeter:
39 3/8"
W x D
(with cylinders):
22" x 22"

2-Cylinder Portable System



3

Height with Flowmeter:
min: 31 3/8"
max: 41 3/8"
W x D
(with cylinders):
22" x 22"

4-Cylinder Portable System



4

43003

Scavenging Circuit II



1 Ultra DC™ Flowmeter

- O₂ Control with On/Off Switch
- Independent O₂ & N₂O Controls
- Oxygen Fail-Safe; Oxygen Flush
- Oxygen Resuscitation Connector
- Emergency Air Valve
- Directional Check Valve
- DISS Inlet Fittings
- Neutral Color

Product Specifications

Dimensions (W x H x D):	4 1/4" x 10 1/2" x 8" 11 cm x 27 cm x 20 cm
Weight:	4.5 lb (2.0 kg)
Mixed Gas Outlet Diameter:	7/8" nominal (22 mm)
Oxygen Inlet:	Male DISS (CGA-1240)
Nitrous Oxide Inlet:	Male DISS (CGA-1040)
Oxygen Flow Range:	3-10 LPM
Nitrous Oxide Flow Range:	0-7 LPM
Inlet Pressure (each gas):	50-55 psi
Resuscitator Connector:	Quick Connect (1/4" Imperial)
Resuscitator Connector Flow:	100 LPM Oxygen

Independent O₂ & N₂O Flow Control



2-Cylinder Portable System

Height with Flowmeter:
39 3/8"
W x D
(with cylinders):
22" x 22"



4-Cylinder Portable System

Height with Flowmeter:
min: 31 3/8"
max: 41 3/8"
W x D
(with cylinders):
22" x 22"



Scavenging Circuit II

Scavenging Circuit II includes a complimentary sampling of single-use nasal masks.

Flowmeter Only

Order No.	Product Description	Retail List
31980	Ultra DC™ Flowmeter (Includes 1)	\$2,162

Portable Systems

Order No.	Product Description	Retail List
46000	Ultra DC™ Flowmeter Package A Ultra DC™ Flowmeter 2-Cylinder Portable System Scavenging Circuit II (Includes 1 2 4)	\$5,038
46100	Ultra DC™ Flowmeter Package B Ultra DC™ Flowmeter 4-Cylinder Portable System Scavenging Circuit II (Includes 1 3 4)	\$5,698



Mobile Stand



Adjustable Mobile Stand



9" Wall Arm



16" Bi-fold Wall Arm



Telescoping Wall Arm
(Extends 12"-17")



Universal Cabinet Slide,
Ultra PC™ % Flowmeters
& Ultra DC™ Flowmeters




Central Systems


Order No.	Product Description	Retail List
46200	Ultra DC™ Flowmeter Package C Ultra DC™ Flowmeter Mobile Stand Scavenging Circuit II (Includes 1 4 5)	\$3,316
46200-ADJS	Ultra DC™ Flowmeter Package C ADJS Ultra DC™ Flowmeter Adjustable Mobile Stand – Short Scavenger Circuit II (Includes 1 4 6)	\$ 3,316
46300	Ultra DC™ Flowmeter Package D Ultra DC™ Flowmeter 9-inch Wall Arm & Mounting Adapter 24552 Scavenging Circuit II (Includes 1 4 7)	\$3,316
46400	Ultra DC™ Flowmeter Package E Ultra DC™ Flowmeter 16-inch Bi-fold Wall Arm & Mounting Adapter 24552 Scavenging Circuit II (Includes 1 4 8)	\$3,316
46500	Ultra DC™ Flowmeter Package F Ultra DC™ Flowmeter Telescoping Wall Arm & Mounting Adapter 24552 Scavenging Circuit II (Includes 1 4 9)	\$3,316
46500-US	Ultra DC™ Flowmeter Package U Ultra DC™ Flowmeter Universal Cabinet Slide, PC/DC (slide mounts to left, right, or under cabinet), and Mounting Adapter 24552 Scavenging Circuit II (Includes 1 4 10)	\$3,316

Additional mounts on pages 17 & 18.
See page 28 for a complete listing of Quick Connect Hoses and Vacuum Male Connectors.

Universal Cabinet Slides, Ultra PC™ % Flowmeters & Ultra DC™ Flowmeters

3 Cabinet Mounting Positions





Flowmeters without Pre-installation Kits


Order No.	Product Description	Retail List
Cabinet Mount Flowmeters (without pre-installation kits):		
36600-NS	Ultra PC™ % Cabinet Mount Flowmeter Pkg. G (includes mounting hardware)	\$5,052
Digital Flushmount Flowmeters (without pre-installation kits):		
56600-NS	Digital Ultra™ Flushmount Flowmeter Pkg. G (includes mounting hardware)	\$5,864
56800-NS	Digital Ultra™ Flushmount Flowmeter Pkg. H (includes mounting hardware, Remote Flow System & mounting hardware)	\$5,484
56950-NS	Digital Ultra™ Flushmount Flowmeter Pkg. I (includes mounting hardware, Accu-Vac™ Vacuum Controller & Accu-Vac™ Flush Mount Faceplate)	\$6,540

Flowmeter Pre-installation Kit Components

Order No.	Product Description	Retail List
Pre-Installation Components:		
23173-FRU	O ₂ Check Valve Riser - DISS (for operatory); includes 28085	\$148
23174-FRU	N ₂ O Check Valve Riser - DISS (for operatory); includes 28086	\$148
28085-FRU	O ₂ DISS x DISS Demand Check Valve Assembly	\$74
28086-FRU	N ₂ O DISS x DISS Demand Check Valve Assembly	\$74

Digital Ultra™ Oral Surgery Flowmeter

Order No.	Product Description	Retail List
Portable Packages		
56100-OS	4-Cylinder Portable (includes portable manifold) Height with Flowmeter: 36 3/8" W x D: 22" x 22"	\$9,100
Central Packages		
56200-OS	Mobile Stand Height with Flowmeter: 36 3/8" W x D: 22" x 22"	\$6,210
56200-OS-ADJT	Adjustable Mobile Stand – Tall Height with Flowmeter: 36 3/8" - 47 3/8" W x D: 22" x 22"	\$6,210
56300-OS	9" Wall Arm	\$6,210
56400-OS	16" Bi-fold Wall Arm	\$6,210
56500-OS	Telescoping Wall Arm	\$6,210
56600-OS	Flushmount	\$7,436



Ultra PC™ % Oral Surgery Flowmeter

Order No.	Product Description	Retail List
Portable Packages		
36000-OS	2-Cylinder Portable (includes portable manifold) Height with Flowmeter: 39 3/8" W x D: 22" x 22"	\$7,266
Central Packages		
36100-OS	4-Cylinder Portable (includes portable manifold) Height with Flowmeter: 39 3/8" W x D: 22" x 22"	\$7,994
36200-OS	Mobile Stand Height with Flowmeter: 39 3/8" W x D: 22" x 22"	\$5,450
36200-OS-ADJT	Adjustable Mobile Stand – Tall Height with Flowmeter: 39 3/8" - x 50 3/8" W x D: 22" x 22"	\$5,450
36300-OS	9" Wall Arm	\$5,450
36400-OS	16" Bi-fold Wall Arm	\$5,450
36500-OS	Telescoping Wall Arm	\$5,450
36600-OS	Cabinet Mount	\$6,582



Ultra DC™ Oral Surgery Flowmeter

Order No.	Product Description	Retail List
Portable Packages		
46000-OS	2-Cylinder Portable (includes portable manifold) Height with Flowmeter: 39 3/8" W x D: 22" x 22"	\$6,270
Central Packages		
46100-OS	4-Cylinder Portable (includes portable manifold) Height with Flowmeter: 39 3/8" W x D: 22" x 22"	\$6,932
46200-OS	Mobile Stand Height with Flowmeter: 39 3/8" W x D: 22" x 22"	\$4,548
46200-OS-ADJT	Adjustable Mobile Stand – Tall Height with Flowmeter: 39 3/8" - 50 3/8" W x D: 22" x 22"	\$4,548
46300-OS	9" Wall Arm	\$4,548
46400-OS	16" Bi-fold Wall Arm	\$4,548
46500-OS	Telescoping Wall Arm	\$4,548

Oral Surgery Accessories



Order No.	Product Description	Retail List
47000	Oral Surgery Accessories Package	\$1,222
23854-FRU	Universal Female Adapter (22 mm male x 15 mm female)	\$23
47050	Directional Y Valve	\$745
47075	Hanks Elbow with Pop-off	\$335
33014	Corrugated Tube – 3-foot Grey	\$52
22137-FRU	Corrugated Tube – 6-Foot Grey	\$79
33029-10	Pediatric Cannula with 7' Tubing (pkg. of 10)	\$96
33030-10	Adult Cannula with 7' Tubing (pkg. of 10)	\$96
38008	Full Facemask – Adult	\$90
38009	Full Facemask – Child	\$90
42006	Adjustable Mobile Stand - Tall H: 28 1/4" - 40" W x D: 22" x 22"	\$514

Oral Surgery Suction Accessories




Order No.	Product Description	Retail List
34014	Surgical E-Vac Kit Includes: Five 1500 ml single-use canisters, wall mount plate and mounting hardware (requires vacuum outlet 34011)	\$1,290
34016	Replacement Canisters for Surgical E-Vac Kit Includes: Five 1500 ml single-use canisters	\$110
26986	6' Tubing – 3/8" ID, Clear PVC	\$46
26987	Canister Holder	\$58
26988	Vacuum Regulator	\$580
34017	Vacuum Regulator Mounting Bracket (for 26988)	\$102
26989	Vacuum Slide	\$82
34011	Recessed Q/C Vacuum Outlet	\$286



Vacuum Controllers

Order No.	Product Description	Retail List
36855	RFS™ Remote Flow System – Analog (for use with Ultra PC™ % Cabinet Mount Flowmeters, and Digital Ultra™ Flushmount Flowmeters with RFS™ Chairmount Kits)	\$1,044
56855	RFS™ Remote Flow System – Digital (for use with Digital Ultra™ Flushmount without RFS™ Chairmount Kits)	\$1,044
56955	Accu-Vac™ Vacuum Controller (for use with Digital Ultra™ Flushmount Flowmeter)	\$1,314



Universal Cabinet Slide Mounts

Order No.	Product Description	Retail List
	Accu-Vac™ Mount: (for Digital Ultra™ Flushmount Flowmeter)	
32015	Universal Cabinet Slide, Accu-Vac™ Vacuum Controller (slide mounts to left, right, or under cabinet)	\$544
32018	Universal Cabinet Slide, Digital RFS™ System (slide mounts to left, right, or under cabinet)	\$544
	RFS™ Remote Flow System A Mounts: (for Ultra PC™ % Cabinet Mount Flowmeter)	
25043	Universal Light Pole Mount Kit (not shown)	\$184
32014	Universal Cabinet Slide, Analog RFS™ System (slide mounts to left, right, or under cabinet)	\$544
	Universal Slide Mount (for Ultra PC™ % Flowmeters & Ultra DC™ Flowmeters)	
32004	Universal Cabinet Slide for Ultra™ PC/DC with Mounting Adapter #24552 (slide mounts to left, right, or under cabinet)	\$544
	Analog Bag Tee Mount: (for Ultra PC™ % Cabinet Mount Flowmeter)	
32008	Universal Cabinet Slide, Analog Bag Tee (slide mounts to left, right, or under cabinet)	\$544
	Digital Bag Tee Mount: (for Digital Ultra™ Flushmount Flowmeter)	
32019	Universal Cabinet Slide, Digital Bag Tee (slide mounts left, right, or under cabinet)	\$544

Mounts continues on page 18.



Flowmeter Mounts		
Order No.	Product Description	Retail List
32001	16" Bi-fold Wall Arm & Mounting Adapter 24552	\$462
32002	9" Wall Arm & Mounting Adapter 24552	\$430
32003	Telescoping Wall Arm (extends 12" - 17") & Mounting Adapter 24552	\$456
32004	Universal Cabinet Slide for Ultra PC/DC (slide mounts to left, right, or under cabinet), includes Mounting Adapter 24552	\$544
32005	Mobile Stand (W x D: 22" x 22" H: 31 1/2")	\$462
32021	Wall bracket with mounting post for the Ultra PC™ Flowmeters and Ultra DC™ Flowmeters	\$238
32055	Surface Mount for Digital Ultra™ Flushmount (flowmeter not included)	\$322
42005-ADJS	Adjustable Mobile Stand – Short (W x D: 22" x 22" H: 19 1/2 - 30")	\$462
42006-ADJT	Adjustable Mobile Stand - Tall (32 1/4" - 43 1/4")	\$514




Flowmeter Mounting Adapters		
Order No.	Product Description	Retail List
24552	Mounting Adapter for Ultra PC and Ultra DC Flowmeters	\$90
32010	Portable and Mobile Stand Adapter – Digital Ultra™ Flowmeters, Ultra PC™ % Flowmeters & Ultra DC™ Flowmeters Flowmeters to Porter™ Flowmeter	\$84
32011	Wall Arm Adapter – Ultra Ultra PC & Ultra DC to Matrix® Wall Arms	\$84

RFS™ Remote Flow System Chairmount Kits



RFS™ Remote Flow System – Not So Remote!

This Illustration depicts the combination of an RFS™ Chairmount Kit and a Flowmeter Package (RFS™ Chairmount Kits are to be purchased separately).



Digital Ultra™ Flushmount Flowmeter

56800-CHAIR
56800-CHAIR-NS

Also available:
Ultra PC™ % Cabinet Mount Flowmeter

36800-CHAIR
36800-CHAIR-NS

See page 21 for illustration of Ultra PC™ Cabinet Mount package with RFS™ mounted.

Chairside Delivery Systems

RFS™ Remote Flow System Mounted at Patient Chair

- Eliminate long mixed gas and vacuum hoses draping from adjacent cabinets
- Create uncluttered, safer work zone for doctor and staff
- Minimize displacement of patient mask caused by drag/pull and weight of long hoses
- Store circuit and rubber goods when not in use



RFS™ Remote Flow System mounts to patient chairs without drilling new holes and can be installed in the field – retrofit where possible or new construction. Kits are available for a wide variety of dental chair brands.

Beginning on page 40 see installation instructions and templates.



A-dec® Chairs



Belmont® Chairs



DCI Series 5 Chairs




DentalEZ® Chairs



Boyd Oral Surgery Chairs



Forest™ Chairs



Marus® Chairs

RFS™ Remote Flow System Chairmount Kits for Chair Brands


Order No.	Product Description	Retail List
31005-1	A-dec® Cascade (1040) RFS™ Chairmount Kit	\$350
31005-2	A-dec® Decade (1021) & Decade + (1221) RFS™ Chairmount	\$350
31005-3	A-dec® 511 RFS™ Chairmount Kit	\$350
31005-4	A-dec® 311 RFS™ Chairmount Kit	\$350
31005-5	A-dec® 400 Series RFS™ Chairmount Kit	\$350
31005-6	A-dec® Performer 8000 Chairmont Kit	\$350
31006-1	Belmont® X-Calibur B-20N RFS™ Chairmount Kit	\$350
31006-2	Belmont® X-Calibur V B-50 RFS™ Chairmount Kit	\$350
31006-3	Belmont® Quolis™ RFS™ Chairmount Kit	\$350
31014-1	Boyd™ Oral Surgery S2600 Series RFS Chairmount Kit	\$350
31016-1	Brumaba™ Oral Surgery OMS Series RFS Chairmount Kit	\$350
31013-1	DCI™ Series 5 RFS™ Chairmount Kit	\$350
31009-1	DentalEZ® Simplicity RFS™ Chairmount Kit	\$350
31009-2	DentalEZ® NuSimplicity RFS™ Chairmount Kit	\$350
31009-3	DentalEZ® J/V Generation RFS™ Chairmount Kit	\$350
31009-4	DentalEZ® Core RFS™ Chairmount Kit	\$350
31018-1	Dexta™ Oral Surgery (20X-1, 20X-3, 25X, 9S, 52X02, 52X-3) RFS™ Chairmount Kit	\$350
31008-1	Engle™ 300 RFS™ Chairmount Kit	\$350
31001-1M	Forest™ 3900 (no extras, metric thread) RFS™ Chairmount Kit	\$350
31001-1S	Forest™ 3900 (touch pad, standard thread), RFS™ Chairmount Kit	\$350
31001-2M	Forest™ Forest 3900 (touch pad, metric thread) RFS™ Chairmount Kit	\$350
31001-2S	Forest™ 3900 (touch pad, standard thread) RFS™ Chairmount Kit	\$350
31001-3	Forest™ 3900 (assist pkg.) RFS™ Chairmount Kit	\$350
31001-4	Forest™ 3900 (touch pad & assist. Pkg.) RFS™ Chairmount Kit	\$350
31004-1	Marus MaxStar™ (no extras) RFS™ Chairmount Kit	\$350
31004-3	Marus NuStar™ (no extras) RFS™ Chairmount Kit	\$350



Note: Chairmount Kits do not include a Remote Flow System. See RFS™ Vacuum Controllers on page 17.

RFS™ Remote Flow System Chairmount Kits for Chair Brands (continued)

Order No.	Product Description	Retail List
31002-1	Midmark® Ultra RFS™ Chairmount Kit	\$350
31002-2	Midmark® Elevance® RFS™ Chairmount Kit	\$350
31002-3	Midmark® 641 Surgery RFS™ Chairmount Kit	\$350
31017-1	MTI™ Oral Surgery 721 RFS™ Chairmount Kit	\$350
31010-1	OSSI Dino™ RFS™ Chairmount Kit	\$350
31012-1	Pediatric Bench™ (standard) RFS™ Chairmount Kit	\$350
31003-1	Pelton & Crane™ 1800 Series RFS™ Chairmount Kit	\$350
31003-2	Pelton & Crane™ 3000 Series RFS™ Chairmount Kit	\$350
31003-3	Pelton & Crane™ 1700 RFS™ Chairmount Kit	\$350
31007-1	Royal® Alliant RFS™ Chairmount Kit	\$350
31007-2	Royal® Signet R16 RFS™ Chairmount Kit	\$350
31007-3	Royal® PD2/OR2/OR3/GP2 (w/o light) RFS™ Chairmount Kit	\$350
31007-4	Royal® PD2/OR2/OR3/GP2 (with light) RFS™ Chairmount Kit	\$350
31007-5	Royal® Domain RFS™ Chairmount Kit	\$350
31007-6	Royal® OS1 Oral Surgery RFS™ Chairmount Kit	\$350
31015-1	Westar™ Oral Surgery (OS VIII & OS III) RFS™ Chairmount Kit	\$350

This Illustration depicts the combination of an RFS™ Chairmount Kit and a Flowmeter Package (RFS™ Chairmount Kits are to be purchased separately).





Cylinder Information: (Approx)	"E" Cylinder: - 4.5" DIA, x 29.5" HT. - O ₂ Volume: 24.9 ft ³ , 690 L - N ₂ O Volume: 56.2 ft ³ , 1,560 L	"G" Cylinder: - 8.5" DIA, x 55" HT. - O ₂ Volume: 222 ft ³ , 6,150 L - N ₂ O Volume: 489 ft ³ , 13,550 L	"H" Cylinder: - 9" DIA, x 56" HT. - O ₂ Volume: 249 ft ³ , 6,900 L - N ₂ O Volume: 558 ft ³ , 15,490 L
--------------------------------	--	---	---

Medical Gas systems must be installed by a licensed medical gas installer (i.e. ASSE 6010) and be verified by a licensed medical gas verifier (i.e. ASSE 6030) and must be installed in compliance with NFPA 99 requirements.

Cylinder Storage Room Requirements:


No more than 3000 cubic feet of gas connected and in storage at one time (excluding nitrogen). Cylinders should be stored in a locked room or cabinet.

- A separate regulator for each gas cylinder
- A pressure relief valve set at 75 psi
- A check valve down stream of each regulator
- System line pressure set at 50 psi
- When the supply system is remote, the switchover shall be automatic

Flexible tubing must not go through wall(s) or partition(s) and have a minimum burst pressure of 1000 psi. When nonflammable gas supply (N₂O-O₂) cylinders are located inside building, they should be in a separate room or enclosure, and separated from the rest of the building by not less than one-hour fire-resistant construction. Two vents of not less than 72 square inches in area should be included. One vent within 6 inches of the floor and one vent within 6 inches of the ceiling. Exception: When an exterior wall cannot be provided see code for requirements. Openings between room or enclosure and interior shall be self-closing smoke and draft control, having a fire protection rating of not less than one hour.

All information provided is intended to assist in the installation of the medical gas system. It is not intended to be complete or eliminate the need to reference NFPA 99 regulations.

Outlet Stations



General Piping Requirements:

Tubes, valves, fittings, outlet stations and other piping components in medical gas systems shall have been cleaned for oxygen service prior to installation. Piping shall be hard-drawn seamless medical gas tube, type K or L (ASTM B819), and bear one of the following markings: OXY, MED, OXY/MED, ACR/OXY or ACR/MED. Piping for non-flammable medical gas systems shall be suitable for oxygen service. Each length of tube shall be permanently labeled and delivered plugged and capped. Fittings, valves and other devices shall be sealed and marked. The installer shall furnish documentation certifying that all installed piping materials comply with the requirements and that the plumber has the appropriate medical gas license.

Testing Requirements:

Piping Integrity and System Verification Final Testing include:

- Initial Pressure Test
- Blow Down
- Piping Purge
- Static Pressure Test
- Pressure Relief Test
- Cross Connection Test
- Alarm Testing for Gas System

Refer to NFPA 99. The responsible facility authority will review all inspection and testing prior to the use of all systems covered by this requirement.

Note: At the time of this printing, NFPA 2005, 2012, 2015, and 2018 are in use in North America. It is important to consult with your local municipality to confirm which codes(s) they have currently adapted.

Digital Automatic Manifold/Alarm System

Pre-install mounting bracket allows effortless installation!



- Automatic cylinder changeover
- Alarm with audio and visual warnings
- Single or multiple zone valve and/or alarm panel capacity
- Clear and precise LED numeric line pressure readings
- Easily accessible components; field upgradeable software
- Optional conduit adapter box available if local codes require wiring inside conduit

Product Specifications

Dimensions with Cover:	18" x 9" x 3" 46 cm x 23 cm x 8 cm
Working Pressure:	50-55 psi
Hose Length:	4 ft.
N ₂ O Outlet:	N ₂ O DISS
O ₂ Outlet:	O ₂ DISS
Primary Voltage:	120 volts
Cable:	Category 6, Shielded, Ethernet 4 twisted pr.

Manifold Package Contents:

Manifold with pressure gauges & high/low pressure alarm switches, system indicator panel, 4 cylinder regulators (2 O₂ & 2 N₂O), four 4 ft. DISS hoses (2 O₂ & 2 N₂O), wall or desk office control panel, cylinder restraints and optional security system.

Digi-Flo™ Manifolds and Valves with Pre-Install Kits

Order No.	Product Description	Retail List
49150	Digi-Flo™ Automatic Switching Manifold/Alarm Pkg. A - Desk Alarm Panel (includes pre-install kit)	\$6,320
49250	Digi-Flo™ Automatic Switching Manifold/Alarm Pkg. B - Wall Alarm Panel (includes pre-install kit)	\$6,320
49350	Digi-Flo™ Automatic Switching Manifold/Alarm Pkg. C - Desk Alarm Panel & Security System (includes pre-install kit)	\$7,054
49450	Digi-Flo™ Automatic Switching Manifold/Alarm Pkg. D - Wall Alarm Panel & Security System (includes pre-install kit)	\$7,054
49550	Digi-Flo™ Zone Valve	\$1,550


Zone Valve required when in remote location.

Digi-Flo™ Manifolds and Alarms without Pre-Install Kits

Order No.	Product Description	Retail List
49150-NP	Manifold with Desk Alarm Panel Pkg. A	\$6,090
49250-NP	Manifold with Wall Alarm Panel Pkg. B	\$6,090
49350-NP	Manifold with Desk Alarm Panel & Security System Pkg. C	\$6,822
49450-NP	Manifold with Wall Alarm Panel & Security System Pkg. D	\$6,822

Digi-Flo™ Manifold Pre-Install Kits

Order No.	Product Description	Retail List
49051	Pre-installation Kit for Digi-Flo™ Automatic Manifold Pkgs. A & C Contents: Manifold template, manifold mounting plate assembly with brass risers, 100 ft. cable, wall outlet box for desk alarm panel connection	\$234
49052	Pre-installation Kit for Digi-Flo™ Automatic Manifold Pkgs. B & D Contents: Manifold template, manifold Mounting plate assembly with brass risers, 100 ft. cable, outlet box and mounting bracket for wall alarm panel	\$234
49053	Pre-Install Kit for Digi-Flo™ 2nd Wall Alarm	\$116
49054	Digi-Flo™ 2nd Wall Alarm (includes splitter for 2nd alarm)	\$810
49054-S	Digi-Flo™ 2nd Wall Alarm with Security System (includes splitter for 2nd alarm)	\$870
27985-FRU	Digi-Flo™ Optional Conduit Kit	\$290
27954-FRU	Matrix® Manifold Adapter Kit	\$158



Digi-Flo™ Manifold Accessories & Adapters

Order No.	Product Description	Retail List
27985-FRU	Digi-Flo™ Optional Conduit Kit	\$290
28191-FRU	Oxygen Add-on Dual Regulator Kit	\$1,534
27954-FRU	Matrix® Manifold Adapter Kit	\$158
28406-FRU	Porter™ Manifold Adapter Kit	\$158
28408-FRU	Belmed Manifold Adapter Kit	\$158
28401-FRU	Accutron™ (Old Style) Manifold Adapter Kit	\$158
49400	Zone Valve Adapter Box – Converts Normally Open (NO) switches on Zone Valves to Digi-Flo™ Manifold	\$520
49400-NC	Zone Valve Adapter Box – Converts Normally Closed (NC) switches on Zone Valves to Digi-Flo™ Manifold	\$520

Analog Manifold/Alarm System



- Manual cylinder changeover
- Pressure relief valve
- Alarm with audio and visual warnings
- Accommodates G & H cylinders
- Zone valve pkg. available




Guardian Monitor™ Manifold and Alarms with Pre-Install Kits

Order No.	Product Description	Retail List
38002	Manifold with Wall Alarm System A – 2 O ₂ + 2 N ₂ O (includes pre-install kit)	\$4,377
38003	Manifold with Desk Alarm System B – 2 O ₂ + 2 N ₂ O (includes pre-install kit)	\$4,377
38004	Manifold with Wall Alarm System C – 2 O ₂ + 1 N ₂ O (includes pre-install kit)	\$3,939
38005	Manifold with Desk Alarm System D – 2 O ₂ + 1 N ₂ O (includes pre-install kit)	\$3,939
39550	Zone Valve (for Guardian Monitor™ Conventional Manifold)	\$1,391

Note: Zone valve required when manifold is in a remote location.

See additional Guardian Monitor™ Manifold products on the following page.



Guardian Monitor™ Manifold and Alarms without Pre-Install Kits

Order No.	Product Description	Retail List
38002-NP	Conventional Manifold/Wall Alarm System A (2+2)	\$4,245
38003-NP	Conventional Manifold/Desktop Alarm System B (2+2)	\$4,245
38004-NP	Conventional Manifold/Wall Alarm System C (2+1)	\$3,805
38005-NP	Conventional Manifold/Desktop Alarm System D (2+1)	\$3,805
39550	Zone Valve (for Guardian Monitor™ Conventional)	\$1,391



Note: Zone valve required when manifold is in a remote location.

Guardian Monitor™ Manifold and Alarms Pre-Install Kits

Order No.	Product Description	Retail List
38051	Pre-Installation Kit for Guardian Monitor™ Conventional Manifold System B & D Contents: Manifold template, copper pipe risers, 100 ft. cable, outlet box for desk control panel	\$132
38052	Pre-Installation Kit for Guardian Monitor™ Conventional Manifold System A & C Contents: Manifold template, copper pipe risers, 100 ft. cable, outlet box for wall control panel	\$132

Manifold Accessories & Components

Order No.	Product Description	Retail List
23123-FRU	Zone Valve Box Trim/Frame Assembly (for zone valve #39550)	\$276
23189	Zone Valve Box Door/Cover (for #39550)	\$138
27945	Zone Valve Box Trim/Frame Assembly (for Zone Valve #49550)	\$850
27703	Zone Valve Box Door/Cover (for #49550)	\$150
23126	Gas Identification Labels (6 O ₂ and 6 N ₂ O)	\$40
23148	100 ft. CAT6 Black Cable and Connector	\$118
27224-FRU	Cylinder Room Caution Sticker	\$6
27319	Purge Fitting Set	\$86
34057	Oxygen Regulator	\$456
34061	Oxygen Hose (DISS 3 feet)	\$130
34058	Nitrous Oxide Regulator	\$456
34062	Nitrous Oxide Hose (DISS 3 feet)	\$130
48075	Nitrogen Regulator	\$506
48078	Nitrogen Hose DISS (DISS 3 feet)	\$148
34075	Cylinder Restraints (accommodates 4 cylinders)	\$260
48025	Cylinder Restraints (accommodates 2 cylinders)	\$154



Portable System Manifolds & Components

Order No.	Product Description	Retail List
20107-FRU	E Cylinder Wrench	\$29
24494	Replacement E Washers (pkg. of 10)	\$32
24549	E Yoke Handle	\$46
32010	Portable and Mobile Stand Adaptor – Accutron™ Digital Ultra™ Flowmeters, Ultra PC™ % Flowmeters & Ultra DC™ Flowmeters to Porter™ Flowmeter	\$84
42000	2-Cylinder Non-Digital Portable Manifold	\$2,454
44000	4-Cylinder Non-Digital Portable Manifold	\$3,218



Hoses & Vacuum Connectors


HOSE ENDS



Hoses & Vacuum Connectors

Order No.	Product Description	Retail List
Compatible with Accutron™ White Slim-line Outlets & Porter™ type Outlets:		
33968	Accutron™ E-Z Q/C X DISS Hose Set*	\$268
33968-90	Accutron™ E-Z Q/C X DISS 90° Hose Set*	\$290
33972	Accutron™ Vacuum Male Connector	\$138
33973	Scavenger Control, Outlet Mount	\$350
Compatible with Accutron™ (silver outlets), Puritan Bennett®, McKesson® & Fraser™ Outlets:		
34070	Puritan Bennett® O ₂ Q/C X DISS Hose*	\$138
34071	Puritan Bennett® N ₂ O Q/C X DISS Hose*	\$138
34072	Puritan Bennett® Vacuum Male Connector	\$138
Compatible with Qxyquip™ Outlets:		
34072	Oxyquip™ O ₂ Q/C X DISS Hose*	\$138
34073	Oxyquip™ N ₂ O Q/C X DISS Hose*	\$138
34084	Oxyquip™ Vacuum Male Connector	\$138
Compatible with Chemtron™, NCG & Coastal Outlets:		
34085	Chemtron™ O ₂ Q/C X DISS Hose*	\$138
34086	Chemtron™ N ₂ O Q/C X DISS Hose*	\$138
34088	Chemtron™ Vacuum Male Connector	\$138
Compatible with Matrx®, Ohio & Centurion Outlets:		
34074	Matrx® O ₂ Q/C X DISS Hose*	\$138
34087	Matrx® N ₂ O Q/C X DISS Hose*	\$138
34080	Matrx® Vacuum Male Connector	\$138
Compatible with Schrader® Inline Outlets:		
34090	Schrader® Inline O ₂ Q/C X DISS Hose	\$132
34091	Schrader® Inline N ₂ O Q/C X DISS Hose	\$132
Compatible with Schrader® Flush Outlets:		
34093	Schrader® Flush O ₂ Q/C X DISS Hose	\$132
34094	Schrader® Flush N ₂ O Q/C X DISS Hose	\$132
DISS Hoses		
34061	O ₂ DISS x DISS Hose	\$138
34061-HT	O ₂ DISS x DISS Hand-Tight	\$160
34061-90	O ₂ DISS x DISS 90°	\$150
34061-90HT	O ₂ DISS Hand-Tight x DISS 90°	\$170
34062	N ₂ O DISS x DISS Hose	\$138
34062-HT	N ₂ O DISS x DISS Hand-Tight	\$160
34062-90	N ₂ O DISS x DISS 90°	\$150
34062-90HT	N ₂ O DISS Hand-Tight x DISS 90°	\$170
Vacuum Connectors		
33972	Accutron™ Vacuum Male Connector	\$138
33974	Vacuum DISS Hand-Twist x 90° 1/2" Hose	\$138
33975	W.A.G.D. DISS Hand-Twist x 90° 1/2" Hose	\$138
Female Q/C Adapter		
35005	N ₂ Hand-Tight 90° x Schrader® Female Q/C	\$232


*Specify lengths up to 6 feet.



Recessed Outlet Stations

Order No.	Product Description	Retail List
34004	Recessed Q/C Triple Outlet (N ₂ O-O ₂ -Vac)*	\$682
34005	Recessed Q/C Dual Outlet (N ₂ O-O ₂)	\$436
34009	Recessed Q/C O ₂ Outlet	\$286
34011	Recessed Q/C Vacuum Outlet*	\$286
34010-I	Recessed Imperial Q/C N ₂ Outlet	\$616
34010-S	Recessed Schrader Q/C N ₂ Outlet	\$616

*Requires use of Vacuum Male Connector #33972 (page 28).



Surface Outlet Stations

Order No.	Product Description	Retail List
34003	Floor Q/C Dual Outlet (N ₂ O-O ₂)	\$602
34060	Surface Q/C Triple Outlet (N ₂ O-O ₂ -Vac)*	\$752
34076	Surface Q/C Dual Outlet (N ₂ O-O ₂)	\$480
34077	Surface Q/C O ₂ Outlet	\$336
34079	Surface Q/C Vacuum Outlet*	\$336
34078-I	Surface Imperial Q/C N ₂ Outlet	\$694
34078-S	Surface Schrader Q/C N ₂ Outlet	\$694

*Requires use of Vacuum Male Connector #33972 (page 28).



Micro Outlet Stations

Order No.	Product Description	Retail List
34012-3	Micro In-line Triple Outlet Kit – DISS* (includes Check Valve Risers, hoses and outlet station)	\$1,102
34012-3-NS	Micro In-line Triple Outlet – DISS* (for mounting within cabinetry)**	\$826
34013-3	Micro In-line Dual Outlet Kit – DISS (includes Check Valve Risers, hoses and outlet station)	\$848
34013-3-NS	Micro In-line Dual Outlet Kit – DISS** (for mounting within cabinetry)	\$554
28000	Micro Outlet Spacer (for Dual & Triple Micro Outlets)	\$26
27793	Oxygen Outlet – DISS with Female Q/C** (for use with Digital Ultra RFS™ Chairmount Kits)	\$276

*Does not include risers and hoses.

**Requires use of Vacuum Male Connector #33972 (page 28).

Emergency Oxygen Equipment



Order No.	Product Description	Retail List
38013	Deluxe Portable Oxygen System Includes: Demand Valve with O2 Q/C Fitting, Adult Facemask, Adjustable Liter Flow Regulator, In-line O2 Adapter, Cylinder Wrench, & Cart (gas cylinder not included)	\$1,680
38010	Portable Oxygen System Includes: Demand Valve with O2 Q/C Fitting, Adult Facemask, Cylinder Wrench, & Cart (gas cylinder not included)	\$1,352
38006	Portable Oxygen System Includes: Bag, Valve, Mask (Manual Resuscitator), Oxygen Tube, Cylinder Wrench, & Cart (gas cylinder not included)	\$886
38011	Demand Valve Resuscitator Includes: Demand Valve, 6' O2 Hose, O2 Male Q/C Fitting	\$808
38007	Manual Resuscitator with Bag, Valve, Mask & Oxygen Tube	\$330
26849-FRU	Oxygen Fitting - DISS with Female Q/C	\$96
38008	Full Facemask - Adult Replacement part for Demand Valve, Manual Resuscitator & Oral Surgery System	\$90
38009	Full Facemask - Child Replacement part for Demand Valve, Manual Resuscitator & Oral Surgery System	\$90
38001	E Cylinder Mobile Cart	\$246
33029-10	Pediatric Nasal Cannula with 7' Tubing (curved tips; *Non-Latex; 10 per pkg.)	\$96
33030-10	Adult Nasal Cannula with 7' Tubing (curved tips; *Non-Latex; 10 per pkg.)	\$96
35018	Deluxe "E" Yoke Regulator with Content Pressure Gauge and Features Two 50 psi DISS Outlets & an Internal Liter-Flow Control with O2 Barb	\$614
35013	"E" Yoke Regulator with Content Pressure Gauge and Features a 50 psi DISS DISS Outlet	\$494
35012	"E" Yoke Regulator, Liter-Flow Control Knob, and O2 DISS with Barb Adapter included (p/n 24589)	\$494

Emergency Oxygen Equipment (continued)

Order No.	Product Description	Retail List
35014-U	Oxygen Liter-Flow (for use with oxygen outlets that do not have a hex nut)	\$160
35014-HL	Oxygen Liter-Flow Meter (for use with oxygen outlets that have a hex nut, <i>including stand-alone outlet 27793</i>)	\$160
34018	Oxygen Liter-Flow Mounting Bracket (for 35014-U)	\$102
35015	Oxygen Liter-flow Selector Valve with Bracket (0-15 lpm) with O2 Barb & O2 Q/C Fitting (p/n 27775)	\$662
35016	Oxygen Liter-flow Selector Valve without Bracket (0-15 lpm) with O2 Barb & O2 Q/C Fitting (p/n 27958)	\$662
35017	Oxygen Liter-flow Selector Valve with 90° Ell (0-15 lpm) with O2 Barb, Set-up for RFS™ Remote Flow System D	\$662
27253-FRU	Oxygen Fitting - DISS with Male Q/C with Bracket	\$96
27254-FRU	Oxygen Fitting - DISS with Male Q/C without Bracket	\$96
27775-FRU	Oxygen Fitting - 1/8" Male Pipe Thread with Male Q/C with Bracket	\$96
27958-FRU	Oxygen Fitting - 1/8" Male Pipe Thread with Male Q/C without Bracket	\$96
27793-FRU	Oxygen Outlet - DISS with Female Q/C*	\$274
24589-FRU	Oxygen DISS with O2 Barb Fitting (green nut & nipple)	\$8

*Does not include riser and hoses.

Oxygen Male Q/C (MQC) Chart

Connectors / Outlets	O ₂ Outlet	Micro O ₂	Micro Double	Micro Triple	RFS™-A	RFS™-D	Accu-Vac™	Ultra PC™ CM	Ultra PC™/DC™	Newport™
O ₂ DISS x MQC w/ Bracket p/n 27253		X			X	X	X	X	X	
1/8" MPT x MQC w/ Bracket p/n 27775		X			X	X *	X	X	X	
O ₂ DISS x MQC w/o Bracket p/n 27254	X	X	X	X	X	X	X	X	X	X
1/8" MPT x MQC w/o Bracket p/n 27958	X	X	X	X	X	X	X	X	X	X

* Note: 35017 O₂ Selector Valve can be attached to system if bag tee is removed, or adding a 90° "ell" and nipple will angle Selector Valve up and out of the way of the bag tee.



Axess™ Low Profile Nasal Masks


Not made with natural rubber latex



"It really does smell like bubblegum!"



FREE SAMPLES!
Call us at **800-531-2221** or
accutron-inc.com



Axess™ Nasal Mask



PATIENT COMFORT AND CLINICAL EFFICIENCY A PERFECT COMBINATION!

- Low profile design provides unencumbered access to the oral cavity
- Fits snugly without tape or uncomfortable tubes protruding into the patient's nostril
- Lightweight, autoclavable scavenging circuit minimizes drag and pull, enhancing patient comfort
- Available in scented and unscented
- Translucent design allows visual monitoring of patient's breathing
- Individually wrapped for greater infection protection
- Lower cost than competing single-use disposable nasal masks
- Single-use minimizes cross-contamination; saves staff time

Axess™ Single-Use Nasal Masks (24 per pkg.)

Order No.	Product Description	Retail List
53037-9	Small Unscented Clear	\$140
53037-12	Small Birthday Bubblegum	\$140
53037-16	Small Fresh Mint	\$140
53035-9	Medium Unscented Clear	\$140
53035-12	Medium Birthday Bubblegum	\$140
53035-16	Medium Fresh Mint	\$140
53034-9	Large Unscented Clear	\$140
53034-12	Large Birthday Bubblegum	\$140
53034-16	Large Fresh Mint	\$140

Axess™ Capnography Adapters (24 per pkg.)

Order No.	Product Description	Retail List
53055	Axess™ Capnography Adapters	\$60

- Attaches via a Luer fitting to a capnography machine's sample line for the purpose of monitoring patient end-tidal CO₂

THE CIRCUIT CAN BE **STEAM-STERILIZED**, WHICH MEANS IT'S ENVIRONMENTALLY FRIENDLY AND COST EFFECTIVE!

The slender and lightweight tubing of the Axess Scavenging Circuit goes hand-in-hand with the design of the Axess Low Profile Mask:

- More comfortable for patients
- Increases the clinician's field of view
- Allows clinician unobstructed access to the oral cavity
- Minimizes drag and pull, so nasal mask stays in place
- Not made from natural rubber latex



Order No. 53006

Scavenging Circuit for RFS™ Remote Flow System and Accu-Vac™ Vacuum Controller



Order No. 28657-FRU

Axess™ Accessory Kit



Order No. 53005

Scavenging Circuit for In-line Vacuum Controller

Axess™ Scavenging Circuits

Order No.	Product Description	Retail List
53005	Axess™ Scavenging Circuit for In-line Vacuum Controller Includes 6-foot long low profile scavenging circuit with In-line vacuum controller, complimentary low profile nasal mask assortment, adapters for connecting most non-Accutron vacuum controllers, instructions for use and reprocessing.	\$300
53006	Axess™ Scavenging Circuit for RFS™ Remote Flow System and Accu-Vac™ Vacuum Controller Includes 6-foot long low profile scavenging circuit (to be used with RFS and Accu-Vac™ Vacuum Controllers – not provided), complimentary low profile nasal mask assortment, adapters for connecting most non-Accutron vacuum controllers, instructions for use and reprocessing.	\$250
53004	Axess™ Extension Tube 32" autoclavable circuit that allows extension of Axess Scavenging Circuits [i.e., Axess™ Scavenging Circuit for In-line Vacuum Controller and Axess™ Scavenging Circuit for Remote Flow System (RFS™) and Accu-Vac™].	\$32
28657-FRU	Axess™ Accessory Kit Includes white cap, 90° outspout connector, universal vacuum adapter.	\$50

See scavenging circuit parts diagram on page 39.



ClearView™ Nasal Mask



ClearView™ Nasal Masks

Not made with natural rubber latex



FREE SAMPLES!
Call us at **800-531-2221** or
accutron-inc.com

Secure Female Luer-Lok® Sample Line Connection



Connect to capnography monitor sample line and verify patient end-tidal CO₂

Luer-Lok® is a registered trademark of Becton Dickinson and Company.



ClearView™ Scavenging Circuits - for Standard Bag Tee

Order No.	Product Description	Retail List
43003	ClearView™ Scavenging Circuit II (includes Scavenging Circuit with Vacuum Control Valve, Corrugated Tube, Reservoir Bag and complimentary ClearView™ Single-Use Nasal Masks)	\$822
43005	ClearView™ Scavenging Circuit I (includes Scavenging Circuit with Vacuum Control Valve and complimentary ClearView™ Single-Use Nasal Masks)	\$715

ClearView™ Scavenging Circuits - for RFS™ Remote Flow System & Accu-Vac™ Vacuum Controller

Order No.	Product Description	Retail List
43004	ClearView™ Scavenging Circuit II (includes Scavenging Circuit Corrugated Tube, Reservoir Bag and complimentary ClearView™ Single-Use Nasal Masks)	\$776
43006	ClearView™ Scavenging Circuit I (includes Scavenging Circuit and complimentary ClearView™ Single-Use Nasal Masks)	\$675

See scavenging circuit parts diagram on page 39.

ClearView™ Adaptor Easy Order Guide

Accutron's ClearView™ Nasal Masks connect effortlessly to:

- ClearView™ Scavenging Circuits I (43005 & 43006)
- ClearView™ Scavenging Circuit II (43003 & 43004)
- Any existing PIP+® Scavenging Circuit
- Any Porter™ Brown Scavenging Circuit (no adaptor required)
- Matrix® Scavenging System (has one smooth tube on each side & comes with two connections) requires adaptor part number 35360 \$57



35360



PIP+™ Nasal Masks

- Single-use
- 8 soothing scents
- Plus unscented
- 3 comfortable sizes
- Not made with natural rubber latex
- Made in the USA



FREE SAMPLES!

Call us at **800-531-2221** or
accutron-inc.com



PIP+™ Nasal Mask

PIP+™ Single-Use Nasal Masks (24 per pkg.)

Order No.	Product Description	Retail List
33015-3	Large Variety Pack 3 (Mint, Peach & Vanilla)	\$103
33015-9	Large Unscented Grey	\$103
33015-14	Large French Vanilla	\$103
33015-15	Large Perfect Peach	\$103
33015-16	Large Fresh Mint	\$103
33016-1	Medium Variety Pack 1 (Orange, Strawberry & Bubblegum)	\$103
33016-2	Medium Variety Pack 2 (Piña Colada, Vanilla & Orange)	\$103
33016-3	Medium Variety Pack 3 (Mint, Peach & Vanilla)	\$103
33016-9	Medium Unscented Grey	\$103
33016-10	Medium Outlaw Orange	\$103
33016-11	Medium Sassy Strawberry	\$103
33016-12	Medium Birthday Bubblegum	\$103
33016-13	Medium Piña Colada	\$103
33016-14	Medium French Vanilla	\$103
33016-15	Medium Perfect Peach	\$103
33016-16	Medium Fresh Mint	\$103
33017-1	Small Variety Pack 1 (Orange, Strawberry & Bubblegum)	\$103
33017-9	Small Unscented Grey	\$103
33017-10	Small Outlaw Orange	\$103
33017-11	Small Sassy Strawberry	\$103
33017-12	Small Birthday Bubblegum	\$103
33017-17	Small Groovy Grape	\$103

Multi-Use Masks (steam sterilize)

Order No.	Product Description	Retail List
33018	Multi-Use Large Mask (1/ea.)	\$140
33019	Multi-Use Medium Mask (1/ea.)	\$131
33020	Multi-Use Small Mask (1/ea.)	\$118



Order No. 32007

PIP+™ and Multi-Use Scavenging Circuits - for Standard Bag Tee

Order No.	Product Description	Retail List
32007	PIP+™ Scavenging Circuit II (includes Scavenging Circuit with Vacuum Control Valve, Corrugated Tube, Reservoir Bag and Complimentary PIP+® Single-Use Nasal Masks)	\$822
33010	PIP+™ Scavenging Circuit I (includes Scavenging Circuit with Vacuum Control Valve and complimentary PIP+™ Single-Use Nasal Masks)	\$715
33011	Scavenger with Large Multi-Use Mask	\$741
33012	Scavenger with Medium Multi-Use Mask	\$741
33013	Scavenger with Small Multi-Use Mask	\$741

PIP+™ Scavenging Circuits - for RFS™ Remote Flow Systems & Accu-Vac™ Vacuum Controller

Order No.	Product Description	Retail List
32009	PIP+™ Scavenging Circuit II (includes Scavenging Circuit, Corrugated Tube, Reservoir Bag and complimentary PIP+™ Single-Use Nasal Masks)	\$776
33008	PIP+™ Scavenging Circuit I (includes Scavenging Circuit and complimentary PIP+™ Single-Use Nasal Masks)	\$675

See scavenging circuit parts diagram on page 39.

PIP+™ Adaptor Easy Order Guide

Accutron's Personal Inhaler Plus+™ Nasal Masks connect effortlessly to Accutron™ Scavenging Circuit I (33008 & 33010) and Scavenging Circuit II (32007 & 32009). They also connect easily to any brand of scavenging system. To order an adaptor, please look for these identifying features and note the Accutron™ order number.

35000	Porter™ Brown Scavenging System Has 4 tubes – 2 ports on each side, mask within mask bearing the name "Brown."	\$59	35350	Matrix® Scavenging System Has one smooth, white tube on each side. \$59
35100	Porter™ Clean Air Scavenging System Has tube inside tube and wings or vacuum sweeps to draw off excess gas.	\$23	35400	McKesson® Scavenging System Has clear mask and 1 tube on each side. \$23
35325	Matrix® Scavenging System Has spiral, white tubing.	\$59	35500	Accutron™ ClearView™ Scavenging System Has 4 tubes - 2 ports on each side, colored & scented mask within clear outer mask. \$59

Axess™, ClearView™ & PIP+™ Scavenging Circuit Diagrams

Rubber Goods & Scavenging Circuit Components



Rubber Goods


Order No.	Product Description	Retail List
22137-FRU	6' (72") Corrugated Tube	\$78
33014	3' (36") Corrugated Tube	\$52
33027	1.5' (18") Corrugated Tube	\$50
33025	Reservoir Bag – 3 Liter	\$47
33026	Reservoir Bag – 2 Liter	\$47

Scavenging Circuit Components

Order No.	Product Description	Retail List
21178	Double Connectors for PIP+™ Nasal Mask (2/pkg.)	\$33
21179	Single Connectors for PIP+™ Nasal Mask (2/pkg.)	\$33
21245	Large Y Connector	\$33
21246	Small Y Connector	\$33
22099	Scavenging Hub for PIP+™ Nasal Masks (clear disk)	\$56
22122	Vacuum Controller w/ Gauge (for standard bag tee) For ClearView™ and PIP+™ Nasal Masks	\$174
28530-FRU	Vacuum Controller w/ Gauge for Axess™ In-line Scavenging Circuit	\$167
28657-FRU	Axess™ Accessory Kit (Includes: 90° Outspot Adapter, White Cap, Matrix® Scavenger Adapter, Porter® In-line Scavenger Adapter & Porter® AVS Scavenger Adapter)	\$50
28463-FRU	Vinyl Cap (7/8" x 1")	\$6
28624-FRU	Axess™ RFS Hose Assembly	\$200
28625-FRU	Axess™ In-Line Hose Assembly	\$200
25590	Double Tubing for PIP+™ Scavenging Circuit (order 2)	\$33
53004	Axess™ Extension Tube (tubing that extends the Axess™ Scavenging Circuit by 32")	\$32
28559-FRU	Axess™ Scavenging Circuit to Accu-Vac™ Adapter	\$36
25591	Spiral Vacuum Tube – 3 1/2'	\$18
25592	Spiral Vacuum Tube – 1 1/2'	\$18
25593	5" Large Tube for PIP+™ Scavenging Circuit (order 2)	\$18
25594	7 1/4" Small Tube for PIP+™ Scavenging Circuit (order 2)	\$18
26252	PIP+™ Scavenging Plug	\$23
26253	Tube Slide Double	\$33
26254	Single Vacuum Connector	\$18
26255-FRU	90° Elbow (creates 90° angle off flowmeter for attachment of scavenging circuit)	\$26
27132	6' Spiral Vacuum Tube (for RFS™ Scavenging Circuit)	\$33
27693-FRU	Double Tubing for ClearView™ Scavenging Circuit (order 2)	\$33


See page 39 for scavenging circuit parts diagrams.

Axess™ Scavenging Circuit Components Diagram (standard in-line vacuum)




Note: For Axess™ Scavenging Circuits for Remote Flow System (RFS™) and Accu-Vac™ Vacuum Controller, part 28517-FRU is replaced with 28516-FRU, and parts 25591 and 28530-FRU are not included.

ClearView™ Scavenging Circuit Components Diagram (standard bag tee)



Note: For ClearView™ Scavenging Circuits for RFS™ Remote Flow System and Accu-Vac™ Vacuum Controller, parts 25591, 22122 and 25592 are replaced with a 6' Spiral Vacuum Tube (part 27132).


PIP+™ Scavenging Circuit Components Diagram (standard bag tee)



Note: For PIP+™ Scavenging Circuits for RFS™ Remote Flow System and Accu-Vac™ Vacuum Controller, parts 25591, 22122 and 25592 are replaced with a 6' Spiral Vacuum Tube (part 27132).



Digital Ultra™ Flushmount Flowmeter 56800-CHAIR



Template not to scale.

1) In advance of construction:

- Provide for tubing chase from flowmeter location to chair J-Box.
- Install Vacuum, Oxygen, and Nitrous Oxide risers.


2) Installation:

- Ensure cut-outs for the flowmeter and oxygen outlet are located correctly.
- Mount Control Module Bracket, with long slot down. Feed cable through the slot. Connect cable to controller and press controller into bracket. Make sure cable is not pinched. To remove controller from the bracket, insert screwdriver into slot at bottom corner of controller and pry outward.
- Pull the Mixed Gas and Vacuum lines up to the Remote Flow System mounting location on the chair.
- Install Remote Flow System to chair.
(Requires purchase of brand-specific Chairmount Kit - see page 19-20).
- Complete connections of all hoses.

OXYGEN OUTLET CUTOOT


<

Ultra PC™ % Cabinet Mount Flowmeter Diagram




Ultra PC™ % Cabinet Mount Flowmeter

Diagram is specific to 36800 Ultra PC™ % Cabinet Mount Flowmeter
with RFS™ Vacuum Controller.




Ultra PC™ % Cabinet Mount Flowmeter Template

Ultra PC™ %
Cabinet Mount
Flowmeter

- Vertical dimension of cut is 4 1/2" MAX
- Horizontal dimension of cut is 6 1/4" MAX
- Minimum Depth is 3 1/2"

Options for vacuum system controllers include: Standard Bag Tee and RFS™ Remote Flow System.


Micro In-Line Triple Outlet Template



Order Nos.
34012-3
34012-3-NS

Cut-out dimensions not to scale.


Micro In-Line Dual Outlet Template



Order Nos.
34013-3
34013-3-NS

Cut-out dimensions not to scale.

Recessed Outlet Station Template



Face Plate Outline

Order Nos.
34004
34005
34009
34010-I
34010-S
34011

Cut-out dimensions not to scale.

ACCUTRON 2-YEAR FLOWMETER LIMITED WARRANTY

IF AN ACCUTRON FLOWMETER NEEDS TO HAVE REPAIR WORK OR REPLACEMENT PARTS DURING THE 2-YEAR WARRANTY PERIOD DUE TO MANUFACTURING DEFECTS, ACCUTRON WILL PROVIDE THE PARTS AND LABOR AT NO CHARGE. THE FLOWMETER OWNER IS RESPONSIBLE ONLY FOR A \$35.00 SHIPPING AND HANDLING FEE, WHICH WILL BE ASSESSED EACH TIME A FLOWMETER IS RETURNED TO ACCUTRON FOR WARRANTY WORK.

WARRANTY TERMS

Limited Warranty and Disclaimer: Accutron ("Seller") warrants that its product will be free from manufacturing defects subject to the terms, conditions, and limitation set forth hereinafter, for a period of 2 years for flowmeters and 1 year for other equipment products. Seller's obligations under this limited warranty are contingent on Buyer's full payment of the product purchase price. Except as specifically set forth above, Seller and its affiliates make no warranties, expressed or implied, and specifically disclaim any warranties of merchantability or fitness for a particular purpose.

The liability of Seller and its affiliates for any claims, losses, damages, or expenses from any cause whatsoever (including acts or omissions of third parties) regardless of the form of the action, whether in tort, contract, or otherwise, shall not exceed the repair cost, replacement cost, or purchase price of the product that directly gives rise to the claim. Seller and its affiliates shall not be liable for any incidental, special, reliance, consequential, or indirect loss or damage rising out of this agreement or the products. As used in this paragraph, consequential damages include, but are not limited to, lost profits, lost revenues, property damage, personal injury damage to the Purchaser or third parties, loss of business or profits, and/or loss of business reputation. It is the sole responsibility of Purchaser to determine the suitability of the products for the Purchaser's intended use. Seller's obligation to repair, replace or refund, as set forth above shall be Buyer's exclusive remedy.

This warranty constitutes the entire warranty. This warranty and Seller's liability hereunder shall be construed according to the laws of the State of Arizona without regard to conflict of law principles.

To activate the Flowmeter Warranty, complete and mail the warranty registration card that accompanies Flowmeter.

Accutron warranties are subject to the following conditions: Accutron products and equipment are warranted to be free from defects in material and workmanship under normal use and service, including all component parts. This warranty shall not apply to defects resulting from accidents, alterations, or misuse. If modifications have affected the operation of the product to render it faulty, this warranty shall be void. This warranty shall be void if any part not of Accutron's manufacture or supply has been incorporated into the product.

THIS WARRANTY IS GIVEN IN PLACE OF ALL OTHER WARRANTIES EXPRESSED OR IMPLIED, OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR OTHERWISE.

No statement or claim about the product by any employee, agent, representative or dealer of Accutron shall constitute a warranty by Accutron or give rise to any liability or obligation of Accutron.

Warranty and Returned Goods Policy

All warranty resolution issues and merchandise returns will be handled through the local authorized Accutron distributor. Contact distributor where unit was purchased. Custom hoses cannot be returned.

Repair Service Policy

All service issues will be handled through the local authorized Accutron distributor. Contact distributor where unit was purchased.

Assistance

For Assistance, contact your local dental distributor or call Accutron Customer Service at:

Toll-free: (800) 531-2221

Local: (623) 780-2020

Fax: (623) 780-0444

Hours of operation: 7:00 AM – 4:30 PM MST

Service ship-to address:

Accutron


1733 W. Parkside Lane

Phoenix, AZ 85027

Visit our website:

www.accutron-inc.com

**NOTES:**



45⁺

YEARS OF INNOVATION



1733 W. Parkside Lane
Phoenix, AZ 85027

800-531-2221 | accutron-inc.com

CROSSTEX
A CANTEL MEDICAL COMPANY

© 2018 Crosstex International, Inc. Part No. 25969 - Rev. 21 - 12/2018
Acutron™, Axess™, ClearView™, PIP+™, Accu-Vac™, Guardian Monitor™, Digi-Flo™, Digital Ultra™, Ultra DC™, Ultra PC™, Newport™ and RFS™ are trademarks of Crosstex International, its affiliates or related companies. All other product names are trademarks of their respective companies.