

### SELF-CONTAINED WATER CAP ASSEMBLY

14.0299.00 Pre-Cascade systems/universal	☒320.00
14.0314.01 3072 Wall mounts	☒320.00
14.0406.00 For Cascade chair-mount	☒320.00
14.0409.00 With QD for cart and wall mounts	☒380.00
14.0424.00 For Radius chair-mount	☒320.00
14.0451.00 For Performer	☒320.00



- Provides water for handpieces and syringe
- Can be used with ICX or an approved maintenance protocol to maintain dental water lines
- Includes (2) 750 ml water bottles
- Includes mounting bracket


### SELF-CONTAINED WATER CAP ASSEMBLY

14.0408.00 Cap assembly	☒90.00
004.182.00 Washer for barb	☒0.25
004.137.00 Washer for cap	☒1.05




### SELF-CONTAINED WATER CAP ASSEMBLY

14.0308.02 Cap assembly	☒90.00
004.182.00 Washer for barb	☒0.25
004.137.00 Washer for cap	☒1.05



### SELF-CONTAINED WATER SYSTEM KIT


14.0416.00 System kit	☒30.00
14.0332.01 Tube (6)	☒25.00



- Includes bottle, bottle cap, washer, supply tube, and barb fitting with washer

### WATER BOTTLE KIT


90.0460.03	☒48.00
------------	--------



- (2) water bottles
- (2) caps
- Extra bottles for the system, they do not replace the system

### WASHER AND GASKET FOR SELF-CONTAINED WATER BOTTLE LIDS

004.137.00 Washer, .750 ID	☒1.05
17.0210.00 Gasket	☒34.50



- Washer for polyurethane water bottles (004.137.00)
- Gasket for portable stainless steel water bottles (17.0210.00)

### TAMPER PROOF BOTTLE CAP

052.015.00	☒0.10
------------	-------



### A-DEC ICX® WATERLINE TREATMENT TABLET


90.1064.00 Case, 0.7-liter tablet	☒1,470.00
-----------------------------------	-----------



- For 0.7-liter bottles
- Case contains (36) boxes with (50) tablets each

### SELF-CONTAINED WATER QD WITH FLOW CONTROL KIT


90.0370.00	☒73.00
------------	--------



- Includes male and female QDs

### MOUNTING BRACKET FOR SELF-CONTAINED WATER

14.0327.01	☒23.00
------------	--------



### SERVICE TIP

When installing a water bottle on the unit, turn the bottle until finger tight, then turn ☒ turn more. Do not overtighten or the threads on the bottle top strip.