



forest dental

Affordable Excellence + Designer Friendly™

## Affordable Excellence

---

Forest is committed to matching the feature set of the big name brands but selling for less. This exclusive offering is the key reason why Forest is today's growth leader through the down economy.

## Designer Friendly™

---

Forest has the industry exclusive feature of dental equipment with a WOW factor. Now more than fifty percent of Forest product is ordered in a color other than the industry standard white.



## Forest is the undisputed growth leader over the last five years of the down economy...

Forest dares to be different with our two brand strategies of Affordable Excellence & Designer Friendly™.

I've always appreciated how cool and classy a car can be. Why not dental equipment? Taking the elements of dental equipment to the level of a beautiful car seemed natural. Dental equipment should be a little more exciting to look at and a lot more patient friendly.

After I purchased Forest in 2004, I began implementing these ideas and before long some of the special elements that make cars so alluring found their way into our dental equipment: many small details done well, rich paint and upholstery color, chrome accents and carbon fiber trim. Thus was born Affordable Excellence and Designer Friendly™ dental equipment.

Affordable Excellence reflects our focus on mirroring many of the expensive features of big name brands with even more precision of fit and finish and selling it for less. Designer Friendly™ signifies that we cater to the taste and style of our customers. Everyone knows color is important in any interior design. We are the only company in the industry that gives you the choice of nine paint colors enhanced with chrome accents and carbon fiber surface treatments. Simply stated, it's beautiful.

Before you build your new clinic or modernize your existing office, please insist that one of our authorized dealers show you the Forest product line. Our customers are quick to say the money they save on Forest, especially in this economy, is a very pleasant surprise. I tell them spending their savings on practice-enhancing technology versus essential equipment is simply smart business.



Hank Barton  
Hank Barton, Owner

# Affordable Excellence

A lot of little things done well is our primary focus at Forest. Before you make your next equipment purchase, take some time to compare Forest to any brand, regardless of price. Our quality is unmatched. And the money you save on Forest chairs, lights, delivery systems and stools goes a long way toward other practice enhancing investments. **Simply smart business.**



1 Headrest



2 Chair Arm Rest



5 Vertical Post with Monitor Arm



6 Chair Mount Adapter & Pivot Delivery Arm

## 1 Headrest

- Stainless Steel Headrest Shaft
- Powder Coated Aluminum Headrest Housing
- Chromed Aluminum Release Button
- Double Articulation

Metal VS Plastic Means Long Life & Easy Function

## 2 Chair Arm Rest

- One-Piece Machined & Polished Stainless Steel Arm Support
  - Chromed Aluminum Release Button
- Unique One-Piece Construction Withstands Extreme Weight  
Button Release Drops Arm Down VS Out To Save Room

## 5 Vertical Post With Monitor Arm

- Aluminum Horizontal Monitor Arm
- Spring Steel Internal Lock Ring
- Steel Vertical Delivery Post

Tight Tolerances Because of Precise Machining  
Perfect Fit Eliminates Sag & Gaps for Great Function

## 6 Chair Mount Adapter & Pivot Delivery Arm

- Cast Iron Monocoque Chair Adapter
  - Needle Roller Bearings
- Chair Adapter Has Superior Rigidity VS Flat Steel  
Easy To Position Delivery & Light Arms



3 Pivot Points



4 Delivery Arms



7 Metal Control Head Chassis on 7000 Series



8 Under the Hood Parts

## A Lot of Little Things Done Well

### 3 Pivot Points

- Hardened Steel Set Screws
  - Self-Locking Nuts & Washers At Pivot Points
- Connected Parts Remain Tight

### 4 Delivery Arms

- Chromed Billet Aluminum End Caps
  - Oil Impregnated Bronze Bushings
  - Friction Adjustment Set Screws on Delivery Arms
- Metal End Caps Won't Break or Discolor Over Time  
Appropriate Bearings With Friction Adjustment Assure  
Adjustable Feather Touch Movement of Arms & No Drift

### 7 Metal Control Head Chassis On 7000 Series

- All Aluminum Control Chassis
- North America's Only All Aluminum Control Head Chassis  
Results In A Lifetime Warranty

### 8 Under The Hood Parts

- Brass Parts w. Electroless Nickel Plating: Air & Water T's, Distribution Block, Master Control Block, Flow Valves
- No Cracking or Warping of Metal Results in Lifetime Warranty

# Designer Friendly™

**What Exactly is Designer Friendly?**

**Designer Friendly is Patient Friendly.**

It is being able to choose from a palette of nine beautiful paint colors including the industry standard white, three carbon fiber look trim pieces and chrome accents.

Your chair, light and delivery systems can now become one with your interior design. Our exclusive Designer Friendly™ makes dental equipment fun to look at.

**Forest is the Only North American Dental Equipment Manufacturer with a WOW Factor.**





# Forest Dental Control Heads

## Affordable Excellence

### Standard IC Control Head

- Material: aluminum hinged lid / 1/8" ABS plastic chassis with super tough DR acrylic protective coating
- Benefit: hinge allows quick access to controls & ancillary products / lid & chassis are resistant to UV, impact & compatible surface cleaners / aluminum lid is strong & lightweight
- Material: electroless nickel plated brass air & water T's, routing valves, pneumatic control & distribution blocks
- Benefit: lifetime warranty on all parts "under the hood" (on all forest control heads) / trouble free handpiece operation
- Material: aluminum handpiece holder bar with glass filled nylon injection molded handpiece holders
- Benefit: allows for up to four handpiece positions, a syringe & ancillary touchpad at your desired angle & spacing

### Euro Control Head

- Material: whip arms rotate on nylon wheels with aluminum & brass components
- Benefit: feather-touch whip movement / ultra light handpiece feel / long, trouble free life
- Material: control head swivels at its center point
- Benefit: 240-degrees of control head rotation allows for beside or over-the-patient placement
- Material: forward mounted aluminum whip arm system
- Benefit: allows for greater reach of handpiece & syringe expanding control head placement options

### Value Control Head

- Material: PVC pinch valve control blocks
- Benefit: ultra-reliable handpiece functions (drive air & air/water coolant) with a long trouble-free life
- Material: ABS plastic control cover w/ tough DR acrylic protective coating
- Benefit: resistant to UV, impact & compatible surface cleaners / easy removal for quick access to controls & ancillary products
- Material: aluminum chassis base plate
- Benefit: strong support for handpiece holder bar, positioning handles & control lid


## Designer Friendly™

- Material: hinged aluminum control head lid
- Benefit: available in nine two-coat powder coat paint colors
- Material: plastic chassis
- Benefit: available in three carbon fiber look & two solid color finishes


- Material: plastic control head lid
- Benefit: available in three carbon fiber looks & two solid color finishes
- Feature: autoclavable silicone handpiece pad available in black & white
- Benefit: blends with Designer Friendly™ choices

- Material: plastic control head lid
- Benefit: available in three carbon fiber looks & two solid color finishes
- Material: chrome plated or powder coated handpiece holder bar
- Benefit: rich look enhancing a value based control head

Affordable Excellence



Designer Friendly™





Pivot Mount Package: Value Control Head, 3900 Chair (Designer Friendly™),  
9080 Light, 1095 Assistant's Instrumentation, 1300RM Monitor Mount

## Chair Mount Packages

### Create a Package

To create a Forest package, first start with the delivery system mounting method & control head of your choice. This is what differentiates one package from another.

---

#### Create a Custom Package

---

use pages 2-3 to build your own package

---

#### Nine Packages Configured By Forest

---

purchase the exact package shown in the picture  
or add/delete items to make it your own.

four pivot chair mount packages

- pivot mounted deliveries with four control head choices.

three fixed chair mount with sidebox packages

- sidebox mounted deliveries with three control head choices.

two fixed chair mount without sidebox packages

- post mounted deliveries with two control head choices.



forest dental

[forestdental.com](http://forestdental.com)

[sales@forestdental.com](mailto:sales@forestdental.com)



chair mount packages  
build a custom package

**CONFIGURE A CUSTOM PACKAGE USING PAGES 2-3**

**STEP 1. CHOOSE YOUR CHAIR MOUNT DELIVERY SYSTEM**

ADDITIONAL EQUIPMENT

LIGHT

LIGHT POST

STOOLS

**STEP 2. ADD UP YOUR TOTAL**

**PACKAGE CHOICE:**

DELIVERY SYSTEM

3900 HYDRAULIC PATIENT CHAIR

ADDITIONAL EQUIPMENT

OPERATORY LIGHT

LIGHT POST

STOOL

POPULAR ADD ONS AND FACTORY OPTIONS

**MODEL #**

**PRICE**

\_\_\_\_\_ \$ \_\_\_\_\_

**3900** \$ **7,390**

\_\_\_\_\_ \$ \_\_\_\_\_

\_\_\_\_\_ \$ \_\_\_\_\_

\_\_\_\_\_ \$ \_\_\_\_\_

\_\_\_\_\_ \$ \_\_\_\_\_

\_\_\_\_\_ \$ \_\_\_\_\_

\_\_\_\_\_ \$ \_\_\_\_\_

\_\_\_\_\_ \$ \_\_\_\_\_

\_\_\_\_\_ \$ \_\_\_\_\_

\_\_\_\_\_ \$ \_\_\_\_\_

\_\_\_\_\_ \$ \_\_\_\_\_

\_\_\_\_\_ \$ \_\_\_\_\_

PACKAGE TOTAL

\$ \_\_\_\_\_ x QTY: \_\_\_\_\_

GRAND TOTAL

\$ \_\_\_\_\_

**EXAMPLE PACKAGE:**

DELIVERY SYSTEM

3900 HYDRAULIC PATIENT CHAIR

ADDITIONAL EQUIPMENT

OPERATORY LIGHT

LIGHT POST

STOOL

POPULAR ADD ONS AND FACTORY OPTIONS

**MODEL #**

**PRICE**

**4987SI** \$ **8,145**

**3900** \$ **7,390**

**1678** \$ **810**

**9045** \$ **2,390**

**1110-900** \$ **200**

**6126C** \$ **790**

**6147C** \$ **1,150**

**1107-450-CS** \$ **640**

**1300RM** \$ **950**

**1700-019-CS** \$ **595**

**1700-081-CS** \$ **500**

**3914-081-CS** \$ **-500**

PACKAGE TOTAL

\$ **23,060** x QTY: **2**

GRAND TOTAL

\$ **46,120**

# chair mount package

## build a custom package

### PIVOT chair mount delivery system



See pages 4-11 for additional information



#### Choose one:

4485PI Standard IC  
Unit Price \$ 5,925  
Base Package Price \$ 15,705  
package example page 4

4987PI Euro  
Unit Price \$ 6,925  
Base Package Price \$ 16,705  
package example page 6

4195PI Value  
Unit Price \$ 4,725  
Base Package Price \$ 14,505  
package example page 8

Option:  
 1495PI Rear Pivot Large IC  
Unit Price \$ 4,995  
Base Package Price \$ 15,340  
package example page 10

#### ADDITIONAL EQUIPMENT



1095 \$ 1,425  
Assistant's Instrumentation  
rear pivot  
(not available with 1495PI)



1590 \$ 3,350  
Cuspidor, Deluxe  
rear pivot with assistant's  
instrumentation  
(not available with 1495PI)



1580 \$ 2,950  
Cuspidor, Deluxe  
rear pivot  
(not available with 1495PI)

### FIXED chair mount with sidebox delivery system



See pages 12-17 for additional information



4485SI Standard IC  
Unit Price \$ 7,145  
Base Package Price \$ 16,625  
package example page 12

4987SI Euro  
Unit Price \$ 8,145  
Base Package Price \$ 17,625  
package example page 14

4195SI Value  
Unit Price \$ 5,945  
Base Package Price \$ 15,425  
package example page 16



1576 \$ 1,450  
Cuspidor, Deluxe  
w/swivel assistant's quad  
holder bar in lieu of  
telescoping arm



1577 \$ 1,450  
Cuspidor, Deluxe



1678 \$ 810  
Cuspidor, Basic  
w/individual assistant's  
instrumentation  
holders in lieu of  
telescoping arm



1679 \$ 810  
Cuspidor, Basic

### FIXED chair mount without sidebox delivery system



See pages 18-21 for additional information



4485FI Standard IC  
Unit Price \$ 5,090  
Base Package Price \$ 14,630  
package example page 18

4195FI Value  
Unit Price \$ 3,890  
Base Package Price \$ 13,430  
package example page 20



1075FI \$ 840  
Telescoping Arm  
w/assistant's instrumentation



1636FI \$ 1,220  
Cuspidor, Basic



1637FI \$ 1,670  
Cuspidor, Basic  
w/assistant's  
instrumentation

# chair mount package

## build a custom package

### ADD TO ANY PACKAGE CONFIGURATION

#### OPERATORY LIGHT

Choose one:



9080 \$ 1,890  
Halogen unit mount light



9045 \$ 2,390  
Halogen unit mount light



9073 \$ 3,490  
LED unit mount light

COMING SUMMER 2013

See page 27 for additional information

#### LIGHT POST

Choose one:



1112-644 \$ 500  
1112-644-3 \$ 1,065 (1495PI)  
Pivot chair mount



1110-900 \$ 200  
Fixed chair mount with sidebox, 44"



1110-820 \$ 260  
Fixed chair mount without sidebox, 35" w/clamp

See pages 4-21 for additional information

#### STOOLS:



6176 \$ 840  
Doctor's comfort  
Air-circulating weave



6126C \$ 790  
Doctor's deluxe  
Ultraleather™ upholstery



6187 \$ 1,200  
Assistant's comfort  
Air-circulating weave



6137C \$ 905  
Assistant's value  
Ultraleather™ upholstery



6147C \$ 1,150  
Assistant's deluxe  
Ultraleather™ upholstery

See page 23 for  
additional information

#### POPULAR ADD ON SUGGESTIONS:

For additional options see pages 4-21



1700-081-CS \$ 500  
Touchpad, unit mount



1110-468 \$ 335  
Swivel tray (standard on Euro  
control — not available on 1495PI)



1700-019-CS \$ 595  
Handpiece illumination system (6 pin)



1700-021-CS \$ 595  
Handpiece illumination system (5 hole)



1400-113-CS \$ 85  
additional HVE with oral dry cup



1107-450-CS \$ 320  
Heated syringe conversion kit  
(1) per syringe (\$320 X 2 = \$640)



1300RM \$ 950  
Monitor mount  
(monitor not included)



SATELEC-S \$ 1845  
 SATELEC-S-LED \$ 2320  
SATELEC piezo scaler uses existing HP position

# pivot chair mount packages

## package P1 | standard IC control head

### STANDARD FEATURES:

#### DELIVERY SYSTEM 4485PI

Automatic control for three handpieces

asepsis 3-way air/water syringe w/tubing, 5'

handpiece tubing, 4-hole Midwest, asepsis style, 5'

pneumatic master control block w/diaphragm

master on/off

individual handpiece pressure adjustment

individual water coolant adjustment

non-retraction water coolant

air/water coolant adjustment

handpiece flush

handpiece exhaust filter/oil collector

dual handle brake release function

self-contained clean water system w/ quick disconnect,

single bottle 1½ liter

wet/dry foot control

Flexible control arm with pneumatic brake

vertical travel 28" / horizontal travel 62"

Internal umbilical, 5' with 3' external from chair to utility center

Contoured utility center

air pressure regulator with filter

master air shut-off valve

Chair adapter, pivot

### ASSISTANT'S INSTRUMENTATION,

#### REAR PIVOT 1095 (option not shown)

Central vacuum accessory package

central vacuum canister w/disposable solids collector

asepsis 3-way air/water syringe w/tubing, 5'

autoclavable standard HVE lever valve w/tubing, 5'

autoclavable saliva ejector lever valve w/tubing, 5'

Internal umbilical, 5' with 3' external

Rear chair pivot mount

#### CUSPIDOR, REAR PIVOT 1580 (option not shown)

One-piece ceramic bowl w/touch lever water activator

Automatic cup fill

Timed bowl rinse

Swing away or over patient's lap

Water pressure regulator with filter

Master water shut-off valve

Internal umbilical, 5' with 3' external

Rear chair pivot mount

#### CUSPIDOR, REAR PIVOT 1590 (option not shown)

*Same as 1580 above, with additional features:*

Swivel assistant's quad holder bar

Assistant's instrumentation

central vacuum canister w/disposable solids collector

asepsis 3-way air/water syringe w/tubing, 5'

autoclavable standard HVE lever valve w/tubing, 5'

autoclavable saliva ejector lever valve w/tubing, 5'

#### LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching temperature setting 4200K to 5500K

Shadow free illumination

#### CHAIR, HYDRAULIC 3900

Ultraleather™ upholstery with wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad w/4 programmable preset positions

Seat height: minimum 17" maximum 29"



## POPULAR EQUIPMENT CHOICES



ASSISTANT'S  
INSTRUMENTATION 1095



CUSPIDOR 1580



CUSPIDOR 1590



MONITOR MOUNT  
1300RM



DOCTOR'S COMFORT  
STOOL 6176



ASSISTANT'S COMFORT  
STOOL 6187

		FACTORY OPTIONS	PART NO	RETAIL	PRICE
<b>BASE PACKAGE</b>	<b>\$ 15,705</b>				
4485PI Delivery System		DELIVERY SYSTEM W/STANDARD IC *	4485PI	\$ 5925	
3900 Chair		Handpiece illumination system			
9080 Light		one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595	
1112-644 Light Post		add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265	
		one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595	
		add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190	
		Additional handpiece position, wet .....	1100-189-CS	405	
		Additional handpiece position, air drive only .....	1100-190-CS	195	
		Less handpiece position .....	1700-099-CS	- 100	
		Heated syringe conversion kit (1) per syringe .....	1107-450-CS	320	
		Assistant's Instrumentation (control head mount) .....	1409-CS	400	
		Quick disconnect w/flow control 1/4" (for water) .....	0016-454-CS	100	
		Quick disconnect w/flow control 3/8" (for air) .....	0016-456-CS	100	
		<b>Swivel tray</b> , * arm mounted (over control head) .....	<b>1110-468</b>	<b>335</b>	
		Touchpad * unit mount for 3900 chair .....	<b>1700-081-CS</b>	<b>500</b>	
		Chrome accents .....	9000-998-CS	135	
		Additional umbilical (per foot) .....	1115-999-CS	30	
		Less self-contained clean water system, w/quick disconnect,			
		single bottle, add water shutoff and regulator .....	1700-355-CS	N/C	
		Routing valve kit w/water shutoff & regulator/factory installation .....	1115-315-CS	235	
		Routing valve kit for use with cuspidor/factory installation .....	1115-316-CS	60	
		Water shut-off valve & regulator kit .....	1115-298	160	
		Ancillaries, uses existing handpiece position (see pages 62-63)			
		Write in selection here .....			
		Write in selection here .....			
		ASSISTANT'S INSTRUMENTATION, REAR PIVOT .....	○ 1095	1425	
		or CUSPIDOR, REAR PIVOT .....	○ 1580	2950	
		or CUSPIDOR W/ASSISTANT'S INSTRUMENTATION, REAR PIVOT .....	○ 1590	3350	
		Quick disconnect w/flow control 1/4" (specify air or water) (1095) .....	0016-455-CS	100	
		Additional autoclavable standard HVE			
		w/tubing, 5' & autoclavable oral dry cup .....	1400-113-CS	85	
		Self-contained clean water system w/quick disconnect,			
		single bottle 1½ liter (1095) .....	1700-319-CS	265	
		Touchpad, unit mount (1095) for 3900 chair .....	1700-088-CS	500	
		Tray assembly (1095) see page 68 .....	1110-680	560	
		Dual articulating quad holder assembly, standard HVE (1095) .....	0024-227-CS	235	
		Wet handpiece position,			
		includes self-contained clean water & foot control (1095) .....	1100-490-CS	895	
		Dry handpiece position includes foot control (1095) .....	1100-493-CS	640	
		Deluxe cuspidor unit mount activation (1580/1590) .....	1100-480-CS	200	
		<b>LIGHT</b> * Halogen (mounting options, see pages 27-30) .....	○ 9080	<b>1890</b>	
		or Halogen .....	○ 9045	2390	
		or LED (coming summer 2013) .....	○ 9073	3490	
		● • <b>Light post, pivot</b> * .....	<b>1112-644</b>	<b>500</b>	
		Monitor mount assembly (monitor screen not included) (see page 32)			
		Ring mount (used with light post) .....	1300RM	950	
		Top post mount (used with monitor post) .....	1300TPM	950	
		• Monitor post (for top post mount) .....	1331-644	500	
		<b>CHAIR, HYDRAULIC</b> * .....	○ 3900	<b>7390</b>	
		Additional touchpad .....	3914-102	500	
		Additional foot control .....	3914-010	500	
		Foot control in lieu of touchpad .....	3914-010-CS	N/C	
		● • <b>Touchpad</b> * deduct chair mount (option on unit) .....	<b>3914-081-CS</b>	<b>- 500</b>	
		Additional chair options (see pages 24-25 for pricing, write selection here) .....			
		Stool, doctor's deluxe / doctor's comfort .....	6126C/6176	790/840	
		Stool, assistant's value w/Ultraleather™ upholstery .....	6137C	905	
		Stool, assistant's deluxe / assistant's comfort .....	6147C/6187	1150/1200	
		* ITEM included in AS SHOWN PACKAGE price .....			
		• Choose only one post when sold with unit .....			
		● Often Overlooked .....			
				<b>TOTAL</b>	

1.800.423.3555

www.forestdental.com

prices subject to change without notice 5

# pivot chair mount packages

## package P2 | euro control head

### STANDARD FEATURES:

#### DELIVERY SYSTEM 4987PI

Automatic control for four handpieces

asepsis 3-way air/water syringe w/tubing, 3'  
handpiece tubing, 4-hole Midwest, asepsis style, 3'  
pneumatic master control block w/diaphragm  
master on/off

individual handpiece pressure adjustment  
individual water coolant adjustment

non-retraction water coolant  
air/water coolant adjustment

handpiece flush

handpiece exhaust filter/oil collector  
dual handle brake release function

swivel instrument tray  
self-contained clean water system, single bottle 1½ liter  
wet/dry foot control

Flexible control arm with pneumatic brake  
vertical travel 28" / horizontal travel 62"

Internal umbilical, 5' with 3' external from chair to utility center  
Contoured utility center

air pressure regulator with filter  
master air shut-off valve

Chair adapter, pivot

### CUSPIDOR, REAR PIVOT 1590

*Same as 1580 below, with additional features:*

Swivel assistant's quad holder bar

Assistant's instrumentation

central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'

### CUSPIDOR, REAR PIVOT 1580 (option not shown)

One-piece ceramic bowl w/touch lever water activator

Automatic cup fill

Timed bowl rinse

Swing away or over patient's lap

Water pressure regulator with filter

Master water shut-off valve

Internal umbilical, 5' with 3' external

Rear chair pivot mount

### ASSISTANT'S INSTRUMENTATION,

#### REAR PIVOT 1095 (option not shown)

Central vacuum accessory package

central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'

Internal umbilical, 5' with 3' external

Rear chair pivot mount

### LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching temperature setting 4200K to 5500K

Shadow free illumination

### CHAIR, HYDRAULIC 3900

Ultraleather™ upholstery w/wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad w/4 programmable preset positions

Seat height: minimum 17" maximum 29"



### POPULAR EQUIPMENT CHOICES



ASSISTANT'S  
INSTRUMENTATION 1095



CUSPIDOR 1580



ADD'L HVE POSITION  
W/ORAL DRY CUP  
1400-113-CS



DOCTOR'S DELUXE  
STOOL 6126C



ASSISTANT'S VALUE  
STOOL 6137C



ASSISTANT'S DELUXE  
STOOL 6147C

		FACTORY OPTIONS	PART NO	RETAIL	PRICE
<b>BASE PACKAGE</b>	<b>\$ 16,705</b>	<b>DELIVERY SYSTEM W/EURO *</b>	<b>4987PI</b>	<b>\$ 6925</b>	
4987PI Delivery System		Handpiece illumination system			
3900 Chair		one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595	
9080 Light		add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265	
1112-644 Light Post		one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595	
		add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190	
<b>AS SHOWN PACKAGE</b>	<b>\$ 21,005</b>	Additional handpiece position, wet .....	1100-189-CS	405	
Additionally equipped with:		Additional handpiece position, air drive only .....	1100-190-CS	195	
1590 Deluxe Cuspidor		Heated syringe conversion kit (1) per syringe .....	1107-450-CS	320	
1300RM Monitor Assembly		Quick disconnect w/flow control 1/4" (for water) .....	0016-454-CS	100	
		Quick disconnect w/flow control 3/8" (for air) .....	0016-456-CS	100	
		Touchpad unit mount .....	1700-081-CS	500	
		Chrome accents .....	9000-998-CS	135	
		Additional umbilical (per foot) .....	1115-999-CS	30	
		Less self-contained clean water system, w/quick disconnect,			
		single bottle, add water shutoff and regulator .....	1700-355-CS	N/C	
		Routing valve kit w/water shutoff & regulator/factory installation .....	1115-315-CS	235	
		Routing valve kit for use with cuspidor/factory installation .....	1115-316-CS	60	
		Water shut-off valve & regulator kit .....	1115-298	160	
<b>DESIGNER FRIENDLY™ SURFACE TREATMENTS</b>		<b>■ Ancillaries, uses existing handpiece postion (see pages 62-63)</b>			
standard unless noted		Write in selection here .....			
Unit plastic finishes:		Write in selection here .....			
<input type="radio"/> Forest Grey	<input type="radio"/> White	<input type="radio"/> Lemans			
<input type="radio"/> Greywood	<input type="radio"/> Laser Leather				
Powder coat paint for unit & light:					
<input type="radio"/> Forest Grey	<input type="radio"/> Shadow	<input type="radio"/> White			
<input type="radio"/> Pearl	<input type="radio"/> Charcoal	<input type="radio"/> Pewter			
<input type="radio"/> Sage	<input type="radio"/> Bronze	<input type="radio"/> Arctic Silver			
Powder coat paint for chair:					
<input type="radio"/> Forest Grey	<input type="radio"/> Shadow				
Add \$1,500:					
<input type="radio"/> White	<input type="radio"/> Pearl	<input type="radio"/> Charcoal			
<input type="radio"/> Pewter	<input type="radio"/> Sage	<input type="radio"/> Bronze			
<input type="radio"/> Arctic Silver					
Ultraleather™ upholstery color choices for chair & deluxe stools:					
<input type="radio"/> Baltic	<input type="radio"/> Dove Grey	<input type="radio"/> Charcoal			
<input type="radio"/> Diplomat Blue	<input type="radio"/> Doe	<input type="radio"/> Grenadine			
<input type="radio"/> Papyrus	<input type="radio"/> Parrot	<input type="radio"/> Raven Wing			
<input type="radio"/> Sassy	<input type="radio"/> Stone	<input type="radio"/> Walnut			
<input type="radio"/> Pearlized Mica	<input type="radio"/> Pearlized Green Tea				
<input type="radio"/> Pearlized Midnight	<input type="radio"/> Pearlized Pewter				
Chair, raffia or special color Ultraleather™	3907-999-SP add \$280				
<input type="radio"/> Durango Dust raffia	<input type="radio"/> Noche raffia	<input type="radio"/> Special color			
<input type="radio"/> Deluxe stool, special color Ultraleather™	6007-999-SP add \$100				
Upholstery choices (must choose one):					
<input type="radio"/> Standard Back, wide/thin	Standard / no charge				
<input type="radio"/> Memory Comfort (plush only)	3907-999-MC add \$800				
<input type="radio"/> Narrow Back	3907-999-NB add \$100				
<input type="radio"/> Narrow Back with slings	3907-999-NS add \$400				
<input type="radio"/> Dual color (plush only)	3907-999-DC add \$280				
<input type="radio"/> Raffia dual color (plush only)	3907-999-RD add \$560				
Upholstery styles (must choose one):					
<input type="radio"/> Plush	Standard / no charge				
<input type="radio"/> Asepsis	3907-999-AS no charge				
Equipment dimensions, see page 72					
Shipping weight estimate:					
delivery unit & chair	450 lb / 204 kg				
base package	580 lb / 264 kg				
as shown package	630 lb / 286 kg				
<b>FACTORY OPTIONS</b>					
<b>DELIVERY SYSTEM W/EURO *</b>					
Handpiece illumination system					
one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595			
add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265			
one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595			
add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190			
Additional handpiece position, wet .....	1100-189-CS	405			
Additional handpiece position, air drive only .....	1100-190-CS	195			
Heated syringe conversion kit (1) per syringe .....	1107-450-CS	320			
Quick disconnect w/flow control 1/4" (for water) .....	0016-454-CS	100			
Quick disconnect w/flow control 3/8" (for air) .....	0016-456-CS	100			
Touchpad unit mount .....	1700-081-CS	500			
Chrome accents .....	9000-998-CS	135			
Additional umbilical (per foot) .....	1115-999-CS	30			
Less self-contained clean water system, w/quick disconnect,					
single bottle, add water shutoff and regulator .....	1700-355-CS	N/C			
Routing valve kit w/water shutoff & regulator/factory installation .....	1115-315-CS	235			
Routing valve kit for use with cuspidor/factory installation .....	1115-316-CS	60			
Water shut-off valve & regulator kit .....	1115-298	160			
<b>■ Ancillaries, uses existing handpiece postion (see pages 62-63)</b>					
Write in selection here .....					
Write in selection here .....					
<b>ASSISTANT'S INSTRUMENTATION, REAR PIVOT</b>		<b>○ 1095</b>	<b>1425</b>		
or CUSPIDOR, REAR PIVOT		<b>○ 1580</b>	<b>2950</b>		
<b>or CUSPIDOR * W/ASSISTANT'S INSTRUMENTATION, REAR PIVOT</b>	<b>○ 1590</b>	<b>3350</b>			
Quick disconnect w/flow control 1/4" (specify air or water) (1095) . 0016-455-CS		100			
Additional autoclavable standard HVE					
w/tubing, 5' & autoclavable oral dry cup .....	1400-113-CS	85			
Self-contained clean water system w/quick disconnect,					
single bottle 1 1/2 liter (1095) .....	1700-319-CS	265			
Touchpad, unit mount (1095) for 3900 chair .....	1700-088-CS	500			
Tray assembly (1095) see page 68 .....	1110-680	560			
Dual articulating quad holder assembly, standard HVE (1095) .... 0024-227-CS		235			
Wet handpiece position,					
includes self-contained clean water & foot control (1095) .... 1100-490-CS		895			
Dry handpiece position includes foot control (1095) .....	1100-493-CS	640			
Deluxe cuspidor unit mount activation (1580/1590) .....	1100-480-CS	200			
<b>LIGHT * Halogen (mounting options, see pages 27-30)</b>	<b>○ 9080</b>	<b>1890</b>			
or Halogen .....	<b>○ 9045</b>	<b>2390</b>			
or LED (coming summer 2013) .....	<b>○ 9073</b>	<b>3490</b>			
<b>• Light post, pivot *</b>	<b>1112-644</b>	<b>500</b>			
Monitor mount assembly (monitor screen not included) (see page 32)					
<b>Ring mount *</b> (used with light post) .....	<b>1300RM</b>	<b>950</b>			
Top post mount (used with monitor post) .....	1300TPM	950			
• Monitor post (for top post mount) .....	1331-644	500			
<b>CHAIR, HYDRAULIC *</b>	<b>○ 3900</b>	<b>7390</b>			
Additional touchpad .....	3914-102	500			
Additional foot control .....	3914-010	500			
Foot control in lieu of touchpad .....	3914-010-CS	N/C			
<b>■ Touchpad deduct chair mount (option on unit) .....</b>	<b>3914-081-CS</b>	<b>- 500</b>			
Additional chair options (see pages 24-25 for pricing, write selection here) ...					
Stool, doctor's deluxe / doctor's comfort .....	6126C/6176	790/840			
Stool, assistant's value w/Ultraleather™ upholstery .....	6137C	905			
Stool, assistant's deluxe / assistant's comfort .....	6147C/6187	1150/1200			
<b>* ITEM included in AS SHOWN PACKAGE price</b>					
• Choose only one post when sold with unit					
<b>■ Often Overlooked</b>					
		<b>TOTAL</b>			

# pivot chair mount packages

## package P3 | value control head

### STANDARD FEATURES:

#### DELIVERY SYSTEM 4195PI

Automatic control for three handpieces

asepsis 3-way air/water syringe w/tubing, 5'  
handpiece tubing, 4-hole Midwest, asepsis style, 5'  
individual pinch valve control blocks  
master on/off

individual handpiece pressure adjustment  
individual water coolant adjustment

non-retraction water coolant

air coolant adjustment

handpiece flush

handpiece exhaust filter/oil collector

dual handle brake release function

self-contained clean water system w/quick disconnect,  
single bottle 1½ liter

wet/dry foot control

Flexible control arm with pneumatic brake

vertical travel 28" / horizontal travel 62"

Internal umbilical, 5' with 3' external from chair to utility center

Contoured utility center

air pressure regulator with filter

master air shut-off valve

Chair adapter, pivot

### ASSISTANT'S INSTRUMENTATION,

#### REAR PIVOT 1095

Central vacuum accessory package

central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'

Internal umbilical, 5' with 3' external

Rear chair pivot mount

#### CUSPIDOR, REAR PIVOT 1580 (option not shown)

One-piece ceramic bowl w/touch lever water activator

Automatic cup fill

Timed bowl rinse

Swing away or over patient's lap

Water pressure regulator with filter

Master water shut-off valve

Internal umbilical, 5' with 3' external

Rear chair pivot mount

#### CUSPIDOR, REAR PIVOT 1590 (option not shown)

*Same as 1580 above, with additional features:*

Swivel assistant's quad holder bar

Assistant's instrumentation

central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'

#### LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching temperature setting 4200K to 5500K

Shadow free illumination

#### CHAIR, HYDRAULIC 3900

Ultraleather™ upholstery w/wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad w/4 programmable preset positions

Seat height: minimum 17" maximum 29"



### POPULAR EQUIPMENT CHOICES



CUSPIDOR 1580



CUSPIDOR 1590



MONITOR MOUNT  
1300RM



DOCTOR'S DELUXE  
STOOL 6126C



ASSISTANT'S VALUE  
STOOL 6137C



ASSISTANT'S DELUXE  
STOOL 6147C

		FACTORY OPTIONS	PART NO	RETAIL	PRICE
<b>BASE PACKAGE</b>	<b>\$ 14,505</b>				
4195PI Delivery System		Handpiece illumination system	4195PI	\$ 4725	
3900 Chair		one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595	
9080 Light		add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265	
1112-644 Light Post		one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595	
		add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190	
<b>AS SHOWN PACKAGE</b>	<b>\$ 16,265</b>				
Additionally equipped with:		Additional handpiece position, wet .....	1100-189-CS	405	
1110-468 Arm Mounted Swivel Tray		Additional handpiece position, air drive only .....	1100-190-CS	195	
1095 Assistant's Instrumentation		Less handpiece position .....	1700-099-CS	-100	
		Heated syringe conversion kit (1) per syringe .....	1107-450-CS	320	
		Assistant's Instrumentation (control head mount) .....	1409-CS	400	
		Quick disconnect w/flow control 1/4" (for water) .....	0016-454-CS	100	
		Quick disconnect w/flow control 3/8" (for air) .....	0016-456-CS	100	
		<b>Swivel tray, * arm mounted (over control head)</b> .....	<b>1110-468</b>	<b>335</b>	
		Touchpad, unit mount for 3900 chair .....	1700-081-CS	500	
		Chrome accents .....	9000-998-CS	135	
		Additional umbilical (per foot) .....	1115-999-CS	30	
		Less self-contained clean water system, w/quick disconnect, single bottle, add water shutoff and regulator .....	1700-355-CS	N/C	
		Routing valve kit w/water shutoff & regulator/factory installation .....	1115-315-CS	235	
		Routing valve kit for use with cuspidor/factory installation .....	1115-316-CS	60	
		Water shut-off valve & regulator kit .....	1115-298	160	
		Ancillaries, uses existing handpiece position (see pages 62-63)			
		Write in selection here .....			
		Write in selection here .....			
		<b>ASSISTANT'S INSTRUMENTATION, * REAR PIVOT</b> .....	<b>1095</b>	<b>1425</b>	
		or CUSPIDOR, REAR PIVOT .....	1580	2950	
		or CUSPIDOR W/ASSISTANT'S INSTRUMENTATION, REAR PIVOT .....	1590	3350	
		Quick disconnect w/flow control 1/4" (specify air or water) (1095) . 0016-455-CS		100	
		Additional autoclavable standard HVE w/tubing, 5' & autoclavable oral dry cup .....	1400-113-CS	85	
		Self-contained clean water system w/quick disconnect, single bottle 1-1/2 liter (1095) .....	1700-319-CS	265	
		Touchpad, unit mount (1095) for 3900 chair .....	1700-088-CS	500	
		Tray assembly (1095) see page 68 .....	1110-680	560	
		Dual articulating quad holder assembly, standard HVE (1095) .....	0024-227-CS	235	
		Wet handpiece position, includes self-contained clean water & foot control (1095) .....	1100-490-CS	895	
		Dry handpiece position includes foot control (1095) .....	1100-493-CS	640	
		Deluxe cuspidor unit mount activation (1580/1590) .....	1100-480-CS	200	
		<b>LIGHT * Halogen (mounting options, see pages 27-30)</b> .....	<b>9080</b>	<b>1890</b>	
		or Halogen .....	9045	2390	
		or LED (coming summer 2013) .....	9073	3490	
		<b>• Light post, pivot *</b> .....	<b>1112-644</b>	<b>500</b>	
		Monitor mount assembly (monitor screen not included) (see page 32)			
		Ring mount (used with light post) .....	1300RM	950	
		Top post mount (used with monitor post) .....	1300TPM	950	
		• Monitor post (for top post mount) .....	1331-644	500	
		<b>CHAIR, HYDRAULIC *</b> .....	<b>3900</b>	<b>7390</b>	
		Additional touchpad .....	3914-102	500	
		Additional foot control .....	3914-010	500	
		Foot control in lieu of touchpad .....	3914-010-CS	N/C	
		<b>• Touchpad deduct chair mount (option on unit)</b> .....	3914-081-CS	-500	
		Additional chair options (see pages 24-25 for pricing, write selection here) .....			
		Stool, doctor's deluxe / doctor's comfort .....	6126C/6176	790/840	
		Stool, assistant's value w/Ultraleather™ upholstery .....	6137C	905	
		Stool, assistant's deluxe / assistant's comfort .....	6147C/6187	1150/1200	
		<b>* ITEM included in AS SHOWN PACKAGE price</b>			
		• Choose only one post when sold with unit			
		<b>• Often Overlooked</b>			
				<b>TOTAL</b>	

Equipment dimensions, see page 72

Shipping weight estimate:

delivery unit & chair	445 lb / 203 kg
base package	556 lb / 257 kg
as shown package	585 lb / 266 kg

1.800.423.3555

www.forestdental.com

prices subject to change without notice 9

# pivot chair mount packages

## package P4 | rear large IC control head

### STANDARD FEATURES:

#### DELIVERY SYSTEM 1495PI

Automatic control for three handpieces

central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'  
handpiece tubing, 4-hole Midwest, asepsis style, 5'  
pneumatic master control block w/diaphragm  
master on/off

individual handpiece pressure adjustment

individual water coolant adjustment

non-retraction water coolant

air coolant adjustment

handpiece flush

handpiece exhaust filter/oil collector

self-contained clean water system w/quick disconnect,  
single bottle 1½ liter

wet/dry foot control

Internal umbilical, 5' with 3' external from chair to utility center

Contoured utility center

air pressure regulator with filter

master air shut-off valve

Rear chair pivot mount

#### FLEX TRAY ASSEMBLY 1800RM

Flexible counterbalanced arm system 26" rigid x 24" flex  
with pneumatic brake; accommodates up to 13½ lbs

Stainless steel tray holder, 9¾" x 13 ½"

wraparound holder bar

dual handle brake release function  
on/off toggle, 3-way with exhaust

Air pressure regulator with filter

Master air shut-off valve

Ring post mount

#### LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching temperature setting 4200K to 5500K

Shadow free illumination

#### CHAIR, HYDRAULIC 3900

Ultraleather™ upholstery with wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad w/4 programmable preset positions

Seat height: minimum 17" maximum 29"



### POPULAR EQUIPMENT CHOICES



REAR CONTROL HEAD MTD  
ASSISTANT'S QUAD  
HOLDER 0024-054



MONITOR MOUNT  
1300RM



ADD'L HVE POSITION  
W/ORAL DRY CUP  
1400-113-CS



DOCTOR'S DELUXE  
STOOL 6126C



ASSISTANT'S VALUE  
STOOL 6137C



ASSISTANT'S DELUXE  
STOOL 6147C

		FACTORY OPTIONS	PART NO	RETAIL	PRICE
<b>BASE PACKAGE</b>	<b>\$ 15,340</b>	<b>DELIVERY SYSTEM REAR PIVOT LARGE IC *</b> ..... 1495PI \$ 4995			
1495PI Delivery System		Handpiece illumination system			
3900 Chair		one handpiece (5-hole), accommodates 1 add'l tubing ..... 1700-021-CS	595		
9080 Light		add'l handpiece illumination tubing (5-hole) ..... 0017-615-CS	265		
1112-644-3 Light Post		one handpiece (6-pin), accommodates 1 add'l tubing ..... 1700-019-CS	595		
		add'l handpiece illumination tubing (6-pin) ..... 0017-603-CS	190		
		Additional handpiece position, wet ..... 1100-189-CS	405		
		Additional handpiece position, air drive only ..... 1100-190-CS	195		
		Less handpiece position ..... 1700-099-CS	-100		
		Additional autoclavable standard HVE w/tubing, 5' & autoclavable oral dry cup ..... 1400-113-CS	85		
		Heated syringe conversion kit ..... 1107-450-CS	320		
		Assistant's quad holder bar ..... 0024-054	150		
		Quick disconnect w/flow control 1/4" (for water) ..... 0016-454-CS	100		
		Quick disconnect w/flow control 3/8" (for air) ..... 0016-456-CS	100		
		Touchpad, unit mount for 3900 chair ..... 1700-081-CS	500		
		Chrome accents ..... 9000-998-CS	135		
		Additional umbilical (per foot) ..... 1115-999-CS	30		
		Less self-contained clean water system, w/quick disconnect, single bottle, add water shutoff and regulator ..... 1700-355-CS	N/C		
		Routing valve kit w/water shutoff & regulator/factory installation ..... 1115-315-CS	235		
		Water shut-off valve & regulator kit ..... 1115-298	160		
		Ancillaries, uses existing handpiece postion (see pages 62-63)			
		Write in selection here			
		Write in selection here			
		<b>LIGHT *</b> Halogen (mounting options, see pages 27-30) ..... 9080	<b>1890</b>		
		or Halogen ..... 9045	2390		
		or LED (coming summer 2013) ..... 9073	3490		
		<b>Light post, pivot *</b> ..... 1112-644-3	<b>1065</b>		
		Light post, pivot if sold with monitor post ..... 1112-644	500		
		Monitor mount assembly (monitor screen not included) (see page 32)			
		Ring mount (used with light post) ..... 1300RM	950		
		Top post mount (used with monitor post) ..... 1300TPM	950		
		Monitor post (for top post mount) ..... 1331-644-3	1065		
		<b>Flex tray assembly,* ring mount</b> ..... 1800RM	<b>1890</b>		
		<b>CHAIR, HYDRAULIC *</b> ..... 3900	<b>7390</b>		
		Additional touchpad ..... 3914-102	500		
		Additional foot control ..... 3914-010	500		
		Foot control in lieu of touchpad ..... 3914-010-CS	N/C		
		Touchpad deduct chair mount (option on unit) ..... 3914-081-CS	-500		
		Additional chair options (see pages 24-25 for pricing, write selection here) ...			
		Stool, doctor's deluxe / doctor's comfort ..... 6126C/6176	790/840		
		Stool, assistant's value w/Ultraleather™ upholstery ..... 6137C	905		
		Stool, assistant's deluxe / assistant's comfort ..... 6147C/6187	1150/1200		
		<b>* ITEM included in AS SHOWN PACKAGE price</b>		<b>TOTAL</b>	
		<b>Often Overlooked</b>			
Upholstery choices (must choose one):					
<input type="radio"/> Standard Back, wide/thin		Standard / no charge			
<input type="radio"/> Memory Comfort (plush only)	3907-999-MC	add \$800			
<input type="radio"/> Narrow Back	3907-999-NB	add \$100			
<input type="radio"/> Narrow Back with slings	3907-999-NS	add \$400			
<input type="radio"/> Dual color (plush only)	3907-999-DC	add \$280			
<input type="radio"/> Raffia dual color (plush only)	3907-999-RD	add \$560			
Upholstery styles (must choose one):					
<input type="radio"/> Plush		Standard / no charge			
<input type="radio"/> Asepsis	3907-999-AS	no charge			
Equipment dimensions, see page 72					
Shipping weight estimate:					
delivery unit & chair	445 lb / 203 kg				
base package	533 lb / 242 kg				
as shown package	538 lb / 245 kg				

# fixed chair mount with sidebox packages

## package S1 | standard IC control head

### STANDARD FEATURES:

#### DELIVERY SYSTEM 4485SI

Automatic control for three handpieces

asepsis 3-way air/water syringe w/tubing, 5'  
handpiece tubing, 4-hole Midwest, asepsis style, 5'  
pneumatic master control block w/diaphragm  
master on/off

individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air coolant adjustment  
handpiece flush  
handpiece exhaust filter/oil collector  
single handle brake release function  
wet/dry foot control

Flexible control arm with pneumatic brake

vertical travel 28" / horizontal travel 62"

Sidebox with self-contained clean water system w/quick  
disconnect, single bottle, 1½ liter

Telescoping arm with swivel quad holder bar

Central vacuum accessory package

central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'

Internal umbilical, 5' with 3' external from chair to utility center

Utility center

contoured housing and cover  
air pressure regulator with filter  
master air shut-off valve

Chair adapter, fixed

### CUSPIDOR, DELUXE 1577/1576 (option not shown)

One-piece ceramic bowl w/touch lever water activator

Automatic cup fill

Timed bowl rinse

Water pressure regulator with filter

Master water shut-off valve

1576 includes:

Swivel assistant's quad holder bar in lieu of telescoping arm

### CUSPIDOR, BASIC 1678/1679 (option not shown)

Removable ceramic bowl

Manual cup fill

Timed bowl rinse

Water pressure regulator with filter

Master water shut-off valve

1678 includes:

Individual assistant's instrumentation holders in lieu  
of telescoping arm

### LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching temperature setting 4200K to 5500K

Shadow free illumination

### CHAIR, HYDRAULIC 3900

Ultraleather™ upholstery with wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad w/4 programmable preset positions

Seat height: minimum 17" maximum 29"



### POPULAR EQUIPMENT CHOICES



CUSPIDOR 1577



CUSPIDOR 1679



MONITOR MOUNT  
1300RM  
(monitor not included)



DOCTOR'S DELUXE  
STOOL 6126C



ASSISTANT'S VALUE  
STOOL 6137C



ASSISTANT'S DELUXE  
STOOL 6147C

	FACTORY OPTIONS	PART NO	RETAIL	PRICE
<b>BASE PACKAGE</b>	<b>\$ 16,625</b>			
4485SI Delivery System				
3900 Chair				
9080 Light				
1110-900 Light Post				
<b>AS SHOWN PACKAGE</b>	<b>\$ 16,960</b>			
Additionally equipped with:				
1110-468 Arm Mounted Swivel Tray				
1700-081-CS Unit Mount Touchpad				
3914-081-CS Deduct Chair Mount Touchpad				
<b>DESIGNER FRIENDLY™ SURFACE TREATMENTS</b> standard unless noted				
Unit plastic finishes:				
<input type="radio"/> Forest Grey	<input type="radio"/> White	<input type="radio"/> Lemans		
<input type="radio"/> Greywood	<input type="radio"/> Laser Leather			
Powder coat paint for unit & light:				
<input type="radio"/> Forest Grey	<input type="radio"/> Shadow	<input type="radio"/> White		
<input type="radio"/> Pearl	<input type="radio"/> Charcoal	<input type="radio"/> Pewter		
<input type="radio"/> Sage	<input type="radio"/> Bronze	<input type="radio"/> Arctic Silver		
Powder coat paint for chair:				
<input type="radio"/> Forest Grey	<input type="radio"/> Shadow			
Add \$1,500:				
<input type="radio"/> White	<input type="radio"/> Pearl	<input type="radio"/> Charcoal		
<input type="radio"/> Pewter	<input type="radio"/> Sage	<input type="radio"/> Bronze		
<input type="radio"/> Arctic Silver				
Ultraleather® upholstery color choices for chair & deluxe stools:				
<input type="radio"/> Baltic	<input type="radio"/> Dove Grey	<input type="radio"/> Charcoal		
<input type="radio"/> Diplomat Blue	<input type="radio"/> Doe	<input type="radio"/> Grenadine		
<input type="radio"/> Papyrus	<input type="radio"/> Parrot	<input type="radio"/> Raven Wing		
<input type="radio"/> Sassy	<input type="radio"/> Stone	<input type="radio"/> Walnut		
<input type="radio"/> Pearlized Mica	<input type="radio"/> Pearlized Green Tea			
<input type="radio"/> Pearlized Midnight	<input type="radio"/> Pearlized Pewter			
Chair, raffia or special color Ultraleather™	3907-999-SP add \$280			
<input type="radio"/> Durango Dust raffia	<input type="radio"/> Noche raffia	<input type="radio"/> Special color		
<input type="radio"/> Deluxe stool, special color Ultraleather™	6007-999-SP add \$100			
Upholstery choices (must choose one):				
<input type="radio"/> Standard Back, wide/thin	Standard / no charge			
<input type="radio"/> Memory Comfort (plush only)	3907-999-MC add \$800			
<input type="radio"/> Narrow Back	3907-999-NB add \$100			
<input type="radio"/> Narrow Back with slings	3907-999-NS add \$400			
<input type="radio"/> Dual color (plush only)	3907-999-DC add \$280			
<input type="radio"/> Raffia dual color (plush only)	3907-999-RD add \$560			
Upholstery styles (must choose one):				
<input type="radio"/> Plush	Standard / no charge			
<input type="radio"/> Asepsis	3907-999-AS no charge			
Equipment dimensions, see page 73				
Shipping weight estimate:				
delivery unit & chair	105 lb / 48 kg			
base package	556 lb / 253 kg			
as shown package	561 lb / 255 kg			
<b>DELIVERY SYSTEM W/STANDARD IC *</b>	<b>4485SI</b>	<b>\$ 7145</b>		
Handpiece illumination system				
one handpiece (5-hole), accommodates 1 add'l tubing	1700-021-CS	595		
add'l handpiece illumination tubing (5-hole)	0017-615-CS	265		
one handpiece (6-pin), accommodates 1 add'l tubing	1700-019-CS	595		
add'l handpiece illumination tubing (6-pin)	0017-603-CS	190		
Additional handpiece position, wet	1100-189-CS	405		
Additional handpiece position	1100-190-CS	195		
Less handpiece position	1700-099-CS	-100		
Heated syringe conversion kit (1) per syringe	1107-450-CS	320		
Additional brake handle activation (left)	0023-642-CS	140		
<b>Swivel tray,* arm mounted option (over control head)</b>	<b>1110-468</b>	<b>335</b>		
Additional autoclavable standard HVE				
w/tubing, 5' & autoclavable oral dry cup	1400-113-CS	85		
Quick disconnect w/flow control 1/4" (for water)	0016-454-CS	100		
Quick disconnect w/flow control 3/8" (for air)	0016-456-CS	100		
Less sidebox telescoping arm and assistant's instrumentation	1400-000-CS	-350		
<b>Touchpad *</b> unit mount for 3900 chair	<b>1700-081-CS</b>	<b>500</b>		
Chrome accents	9000-998-CS	135		
Additional umbilical (per foot)	1115-999-CS	30		
Less self-contained clean water system, w/quick disconnect,				
single bottle, add water shutoff and regulator	1700-355-CS	N/C		
Routing valve kit w/water shutoff & regulator/factory installation	1115-315-CS	235		
Routing valve kit for use with cuspidor/factory installation	1115-316-CS	60		
Water shut-off valve & regulator kit	1115-298	160		
<b>Ancillaries, uses existing handpiece postion (see pages 62-63)</b>				
Write in selection here				
Write in selection here				
CUSPIDOR, DELUXE (in lieu of telescoping arm - 1576)	○ 1576 ○ 1577	1450		
or BASIC (in lieu of telescoping arm - 1678)	○ 1678 ○ 1679	810		
Basic cuspidor automatic cup fill	1700-055-CS	195		
Deluxe cuspidor unit mount activation (1576/1577)	1100-480-CS	200		
<b>LIGHT *</b> Halogen (mounting options, see pages 27-30)	○ 9080	1890		
or Halogen	○ 9045	2390		
or LED (coming summer 2013)	○ 9073	3490		
<b>Light post, 44" *</b>	<b>1110-900</b>	<b>200</b>		
Monitor mount assembly (monitor screen not included) (see page 32)				
Ring mount (used with light post)	1300RM	950		
Top post mount (used with monitor post)	1300TPM	950		
Monitor post (for top post mount)	1331-998	200		
<b>CHAIR, HYDRAULIC *</b>	<b>3900</b>	<b>7390</b>		
Additional touchpad	3914-102	500		
Additional foot control	3914-010	500		
Foot control in lieu of touchpad	3914-010-CS	N/C		
<b>Touchpad *</b> deduct chair mount (option on unit)	<b>3914-081-CS</b>	<b>- 500</b>		
Additional chair options (see pages 24-25 for pricing, write selection here)...				
Stool, doctor's deluxe / doctor's comfort	6126C/6176	790/840		
Stool, assistant's value w/Ultraleather™ upholstery	6137C	905		
Stool, assistant's deluxe / assistant's comfort	6147C/6187	1150/1200		
<b>* ITEM included in AS SHOWN PACKAGE price</b>				
<b>Often Overlooked</b>				
		<b>TOTAL</b>		

# fixed chair mount with sidebox packages

## package S2 | euro control head

### STANDARD FEATURES:

#### DELIVERY SYSTEM 4987SI

Automatic control for four handpieces

asepsis 3-way air/water syringe w/tubing, 3'  
handpiece tubing, 4-hole Midwest, asepsis style, 3'  
pneumatic master control block w/diaphragm  
master on/off

individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air coolant adjustment  
handpiece flush  
handpiece exhaust filter/oil collector  
dual handle brake release function  
swivel instrument tray  
wet/dry foot control

Flexible control arm with pneumatic brake  
vertical travel 28" / horizontal travel 62"

Sidebox with self-contained clean water system, w/quick  
disconnect single bottle, 1½ liter

Telescoping arm with swivel quad holder bar

Central vacuum accessory package

central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'

Internal umbilical, 5' with 3' external from chair to utility center

Utility center

contoured housing and cover  
air pressure regulator with filter  
master air shut-off valve

Chair adapter, fixed

### CUSPIDOR, DELUXE 1577/1576 (1576 option shown)

One-piece ceramic bowl w/touch lever water activator

Automatic cup fill

Timed bowl rinse

Water pressure regulator with filter

Master water shut-off valve

1576 includes:

Swivel assistant's quad holder bar in lieu of telescoping arm

### CUSPIDOR, BASIC 1678/1679 (option not shown)

Removable ceramic bowl

Manual cup fill

Timed bowl rinse

Water pressure regulator with filter

Master water shut-off valve

1678 includes:

Individual assistant's instrumentation holders in lieu  
of telescoping arm

### LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching

temperature setting 4200K to 5500K

Shadow free illumination

### CHAIR, HYDRAULIC 3900

Ultraleather™ upholstery with wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad w/4 programmable preset positions

Seat height: minimum 17" maximum 29"



### POPULAR EQUIPMENT CHOICES



CUSPIDOR 1577



CUSPIDOR 1678



CUSPIDOR 1679



MONITOR MOUNT  
1300RM



DOCTOR'S COMFORT  
STOOL 6176



ASSISTANT'S COMFORT  
STOOL 6187

		FACTORY OPTIONS	PART NO	RETAIL	PRICE
<b>BASE PACKAGE</b>	<b>\$ 17,625</b>				
4987SI Delivery System		DELIVERY SYSTEM W/EURO * ..... 4987SI	\$ 8145		
3900 Chair		Handpiece illumination system			
9080 Light		one handpiece (5-hole), accommodates 1 add'l tubing ..... 1700-021-CS	595		
1110-900 Light Post		add'l handpiece illumination tubing (5-hole) ..... 0017-615-CS	265		
		one handpiece (6-pin), accommodates 1 add'l tubing ..... 1700-019-CS	595		
		add'l handpiece illumination tubing (6-pin) ..... 0017-603-CS	190		
<b>AS SHOWN PACKAGE</b>	<b>\$ 19,075</b>				
Additionally equipped with:		Additional handpiece position, wet ..... 1100-189-CS	405		
1576 Deluxe Cuspidor		Additional handpiece position, air drive only ..... 1100-190-CS	195		
		Heated syringe conversion kit (1) per syringe ..... 1107-450-CS	320		
		Additional autoclavable standard HVE			
		w/tubing, 5' & autoclavable oral dry cup ..... 1400-113-CS	85		
		Quick disconnect w/flow control 1/4" (for water) ..... 0016-454-CS	100		
		Quick disconnect w/flow control 3/8" (for air) ..... 0016-456-CS	100		
		Less sidebox telescoping arm and assistant's instrumentation ..... 1400-000-CS	-350		
		Touchpad unit mount ..... 1700-081-CS	500		
		Chrome accents ..... 9000-998-CS	135		
		Additional umbilical (per foot) ..... 1115-999-CS	30		
		Less self-contained clean water system, w/quick disconnect,			
		single bottle, add water shutoff and regulator ..... 1700-355-CS	N/C		
		Routing valve kit w/water shutoff & regulator/factory installation ..... 1115-315-CS	235		
		Routing valve kit for use with cuspidor/factory installation ..... 1115-316-CS	60		
		Water shut-off valve & regulator kit ..... 1115-298	160		
		Ancillaries, uses existing handpiece postion (see pages 62-63)			
		Write in selection here .....			
		Write in selection here .....			
		<b>CUSPIDOR, DELUXE *</b> (in lieu of telescoping arm - 1576) <input checked="" type="radio"/> 1576 <input type="radio"/> 1577	<b>1450</b>		
		or BASIC (in lieu of telescoping arm - 1678) ..... <input type="radio"/> 1678 <input checked="" type="radio"/> 1679	810		
		Basic cuspidor automatic cup fill ..... 1700-055-CS	195		
		Deluxe cuspidor unit mount activation (1576/1577) ..... 1100-480-CS	200		
		<b>LIGHT *</b> Halogen (mounting options, see pages 27-30) ..... <input type="radio"/> 9080	<b>1890</b>		
		or Halogen ..... <input type="radio"/> 9045	2390		
		or LED (coming summer 2013) ..... <input type="radio"/> 9073	3490		
		<b>Light post, 44" *</b> ..... 1110-900	<b>200</b>		
		Monitor mount assembly (monitor screen not included) (see page 32)			
		Ring mount (used with light post) ..... 1300RM	950		
		Top post mount (used with monitor post) ..... 1300TPM	950		
		Monitor post (for top post mount) ..... 1331-998	200		
		<b>CHAIR, HYDRAULIC *</b> ..... 3900	<b>7390</b>		
		Additional touchpad ..... 3914-102	500		
		Additional foot control ..... 3914-010	500		
		Foot control in lieu of touchpad ..... 3914-010-CS	N/C		
		Touchpad deduct chair mount (option on unit) ..... 3914-081-CS	-500		
		Additional chair options (see pages 24-25 for pricing, write selection here) ...			
		Stool, doctor's deluxe / doctor's comfort ..... 6126C/6176	790/840		
		Stool, assistant's value w/Ultraleather™ upholstery ..... 6137C	905		
		Stool, assistant's deluxe / assistant's comfort ..... 6147C/6187	1150/1200		
		<b>* ITEM</b> included in AS SHOWN PACKAGE price	<b>TOTAL</b>		
		Often Overlooked			
Equipment dimensions, see page 73					
Shipping weight estimate:					
delivery unit & chair		110 lb / 50 kg			
base package		551 lb / 251 kg			
as shown package		571 lb / 260 kg			

# fixed chair mount with sidebox packages

## package S3 | value control head

### STANDARD FEATURES:

#### DELIVERY SYSTEM 4195SI

Automatic control for three handpieces

asepsis 3-way air/water syringe w/tubing, 5'  
handpiece tubing, 4-hole Midwest, asepsis style, 5'

individual pinch valve control blocks

master on/off

individual handpiece pressure adjustment

individual water coolant adjustment

non-retraction water coolant adjustment

air coolant adjustment

handpiece flush

handpiece exhaust filter/oil collector

single handle brake release function

wet/dry foot control

Flexible control arm with pneumatic brake

vertical travel 28" / horizontal travel 62"

Sidebox with self-contained clean water system, w/quick disconnect single bottle, 1½ liter

Telescoping arm with swivel quad holder bar

Central vacuum accessory package

central vacuum canister w/disposable solids collector

asepsis 3-way air/water syringe w/tubing, 5'

autoclavable standard HVE lever valve w/tubing, 5'

autoclavable saliva ejector lever valve w/tubing, 5'

Internal umbilical, 5' with 3' external from chair to utility center

Utility center

contoured housing and cover

air pressure regulator with filter

master air shut-off valve

Chair adapter, fixed

### CUSPIDOR, BASIC 1678/1679 (1678 option shown)

Removable ceramic bowl

Manual cup fill

Timed bowl rinse

Water pressure regulator with filter

Master water shut-off valve

1678 includes:

Individual assistant's instrumentation holders in lieu of telescoping arm

### CUSPIDOR, DELUXE 1577/1576 (option not shown)

One-piece ceramic bowl w/touch lever water activator

Automatic cup fill

Timed bowl rinse

Water pressure regulator with filter

Master water shut-off valve

1576 includes:

Swivel assistant's quad holder bar in lieu of telescoping arm

### LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching temperature setting 4200K to 5500K

Shadow free illumination

### CHAIR, HYDRAULIC 3900

Ultraleather™ upholstery with wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad w/4 programmable preset positions

Seat height: minimum 17" maximum 29"



### POPULAR EQUIPMENT CHOICES



CUSPIDOR 1576



CUSPIDOR 1577



CUSPIDOR 1679



DOCTOR'S DELUXE  
STOOL 6126C



ASSISTANT'S VALUE  
STOOL 6137C



ASSISTANT'S DELUXE  
STOOL 6147C

		FACTORY OPTIONS	PART NO	RETAIL	PRICE
<b>BASE PACKAGE</b>	<b>\$ 15,425</b>				
4195SI Delivery System		DELIVERY SYSTEM W/VALUE * ..... 4195SI	\$ 5945		
3900 Chair		Handpiece illumination system			
9080 Light		one handpiece (5-hole), accommodates 1 add'l tubing ..... 1700-021-CS	595		
1110-900 Light Post		add'l handpiece illumination tubing (5-hole) ..... 0017-615-CS	265		
		one handpiece (6-pin), accommodates 1 add'l tubing ..... 1700-019-CS	595		
		add'l handpiece illumination tubing (6-pin) ..... 0017-603-CS	190		
<b>AS SHOWN PACKAGE</b>	<b>\$ 16,235</b>				
Additionally equipped with:		Additional handpiece position, wet ..... 1100-189-CS	405		
1700-081-CS Unit Mount Touchpad		Additional handpiece position, air drive only ..... 1100-190-CS	195		
1678 Basic Cuspidor		Less handpiece position ..... 1700-099-CS	-100		
3914-081-CS Deduct Chair Mount Touchpad		Heated syringe conversion kit (1) per syringe ..... 1107-450-CS	320		
		Additional brake handle activation (left) ..... 0023-642-CS	140		
		Swivel tray, arm mounted option (over control head) ..... 1110-468	335		
		Additional autoclavable standard HVE w/tubing, 5' & autoclavable oral dry cup ..... 1400-113-CS	85		
		Quick disconnect w/flow control 1/4" (for water) ..... 0016-454-CS	100		
		Quick disconnect w/flow control 3/8" (for air) ..... 0016-456-CS	100		
		Less sidebox telescoping arm and assistant's instrumentation ..... 1400-000-CS	-350		
		Touchpad * unit mount for 3900 chair ..... 1700-081-CS	500		
		Chrome accents ..... 9000-998-CS	135		
		Additional umbilical (per foot) ..... 1115-999-CS	30		
		Less self-contained clean water system, w/quick disconnect, single bottle, add water shutoff and regulator ..... 1700-355-CS	N/C		
		Routing valve kit w/water shutoff & regulator/factory installation ..... 1115-315-CS	235		
		Routing valve kit for use with cuspidor/factory installation ..... 1115-316-CS	60		
		Water shut-off valve & regulator kit ..... 1115-298	160		
		Ancillaries, uses existing handpiece postion (see pages 62-63)			
		Write in selection here .....			
		Write in selection here .....			
		<b>CUSPIDOR, BASIC *</b> (in lieu of telescoping arm - 1678) ○ 1678 ○ 1679	810		
		or Deluxe (in lieu of telescoping arm - 1576) ..... ○ 1576 ○ 1577	1450		
		Basic cuspidor automatic cup fill ..... 1700-055-CS	195		
		Deluxe cuspidor unit mount activation (1576/1577) ..... 1100-480-CS	200		
		<b>LIGHT *</b> Halogen (mounting options, see pages 27-30) ..... ○ 9080	1890		
		or Halogen ..... ○ 9045	2390		
		or LED (coming summer 2013) ..... ○ 9073	3490		
		Light post, 44" * ..... 1110-900	200		
		Monitor mount assembly (monitor screen not included) (see page 32)			
		Ring mount (used with light post) ..... 1300RM	950		
		Top post mount (used with monitor post) ..... 1300TPM	950		
		Monitor post (for top post mount) ..... 1331-998	200		
		<b>CHAIR, HYDRAULIC *</b> ..... ○ 3900	7390		
		Additional touchpad ..... 3914-102	500		
		Additional foot control ..... 3914-010	500		
		Foot control in lieu of touchpad ..... 3914-010-CS	N/C		
		Touchpad * deduct chair mount (option on unit) ..... 3914-081-CS	- 500		
		Additional chair options (see pages 24-25 for pricing, write selection here) ...			
		Stool, doctor's deluxe / doctor's comfort ..... 6126C/6176	790/840		
		Stool, assistant's value w/Ultraleather™ upholstery ..... 6137C	905		
		Stool, assistant's deluxe / assistant's comfort ..... 6147C/6187	1150/1200		
		* ITEM included in AS SHOWN PACKAGE price		TOTAL	
		Often Overlooked			
Equipment dimensions, see page 73					
Shipping weight estimate:					
delivery unit & chair	104 lb / 48 kg				
base package	547 lb / 249 kg				
as shown package	567 lb / 258 kg				

# fixed chair mount without sidebox packages

## package F1 | standard IC control head

### STANDARD FEATURES:

#### DELIVERY SYSTEM 4485FI

Automatic control for three handpieces

- asepsis 3-way air/water syringe w/tubing, 5'
- handpiece tubing, 4-hole Midwest, asepsis style, 5'
- pneumatic master control block w/diaphragm
- master on/off
- individual handpiece pressure adjustment
- individual water coolant adjustment
- non-retraction water coolant
- air coolant adjustment
- handpiece flush
- handpiece exhaust filter/oil collector
- single handle brake release function
- self-contained clean water system, w/quick disconnect single bottle, 1½ liter
- wet/dry foot control

Flexible control arm with pneumatic brake  
vertical travel 28" / horizontal travel 62"

Internal umbilical, 5' with 3' external from chair to utility center  
Utility center

Contoured housing and cover  
air pressure regulator with filter  
master air shut-off valve

Chair adapter, fixed

### ASSISTANT'S INSTRUMENTATION 1075FI

Central vacuum accessory package

- central vacuum canister w/disposable solids collector
- asepsis 3-way air/water syringe w/tubing, 5'
- autoclavable standard HVE lever valve w/tubing, 5'
- autoclavable saliva ejector lever valve w/tubing, 5'

Internal umbilical, 5' with 3' external from chair to utility center  
Post mount, 2"

### CUSPIDOR, BASIC 1636FI/1637FI (option not shown)

Removable ceramic bowl

Manual cup fill

Timed bowl rinse

Internal umbilical, 5' with 3' external from chair to utility center  
Water pressure regulator with filter

Master water shut-off valve

Post mount, 2"

1637FI includes:

- Individual assistant's instrumentation holders
- Assistant's instrumentation
- central vacuum canister w/disposable solids collector
- asepsis 3-way air/water syringe w/tubing, 5'
- autoclavable standard HVE lever valve w/tubing, 5'
- autoclavable saliva ejector lever valve w/tubing, 5'

### LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching temperature setting 4200K to 5500K

Shadow free illumination

### CHAIR, HYDRAULIC 3900

Ultraleather™ upholstery with wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad w/4 programmable preset positions

Seat height: minimum 17" maximum 29"



### POPULAR EQUIPMENT CHOICES



CUSPIDOR 1636FI



CUSPIDOR 1637FI



MONITOR MOUNT  
1300RM  
(monitor not included)



DOCTOR'S DELUXE  
STOOL 6126C



ASSISTANT'S VALUE  
STOOL 6137C



ASSISTANT'S DELUXE  
STOOL 6147C

		FACTORY OPTIONS	PART NO	RETAIL	PRICE
<b>BASE PACKAGE</b>	<b>\$ 14,630</b>				
4485FI Delivery System		<b>DELIVERY SYSTEM W/STANDARD IC *</b> ..... 4485FI	<b>\$ 5090</b>		
3900 Chair		Handpiece illumination system			
9080 Light		one handpiece (5-hole), accommodates 1 add'l tubing ..... 1700-021-CS	595		
1110-820 Light Post		add'l handpiece illumination tubing (5-hole) ..... 0017-615-CS	265		
		one handpiece (6-pin), accommodates 1 add'l tubing ..... 1700-019-CS	595		
		add'l handpiece illumination tubing (6-pin) ..... 0017-603-CS	190		
<b>AS SHOWN PACKAGE</b>	<b>\$ 15,805</b>				
Additionally equipped with:		Additional handpiece position, wet ..... 1100-189-CS	405		
1110-468 Arm Mounted Swivel Tray		Additional handpiece position, air drive only ..... 1100-190-CS	195		
1075FI Assistant's Instrumentation		Less handpiece position ..... 1700-099-CS	-100		
		Heated syringe conversion kit (1) per syringe ..... 1107-450-CS	320		
		Assistant's Instrumentation (control head mount) ..... 1409-CS	400		
		Quick disconnect w/flow control 1/4" (for water) ..... 0016-454-CS	100		
		Quick disconnect w/flow control 3/8" (for air) ..... 0016-456-CS	100		
		Additional brake handle activation (left) ..... 0023-642-CS	140		
		<b>Swivel tray,* arm mounted (over control head)</b> ..... 1110-468	<b>335</b>		
		Touchpad, unit mount for 3900 chair ..... 1700-081-CS	500		
		Chrome accents ..... 9000-998-CS	135		
		Deluxe utility center housing and cover in lieu of contoured ..... 1115-007-CS	-200		
		Additional umbilical (per foot) ..... 1115-999-CS	30		
		Less self-contained clean water system, w/quick disconnect, single bottle, add water shutoff and regulator ..... 1700-355-CS	N/C		
		Routing valve kit w/water shutoff & regulator/factory installation ..... 1115-315-CS	235		
		Routing valve kit for use with cuspidor/factory installation ..... 1115-316-CS	60		
		Water shut-off valve & regulator kit ..... 1115-298	160		
		Ancillaries, uses existing handpiece postion (see pages 62-63)			
		Write in selection here .....			
		Write in selection here .....			
		<b>ASSISTANT'S INSTRUMENTATION *</b> ..... 1075FI	<b>840</b>		
		or BASIC CUSPIDOR ..... 1636FI	1220		
		or BASIC CUSPIDOR w/assistant's instrumentation ..... 1637FI	1670		
		Additional autoclavable standard HVE w/tubing, 5' & autoclavable oral dry cup (1075FI/1637FI) ..... 1400-113-CS	85		
		Basic cuspidor automatic cup fill ..... 1700-055-CS	195		
		<b>LIGHT *</b> Halogen (mounting options, see pages 27-30) ..... 9080	<b>1890</b>		
		or Halogen ..... 9045	2390		
		or LED (coming summer 2013) ..... 9073	3490		
		<b>Light post, 35"</b> * ..... 1110-820	<b>260</b>		
		Monitor mount assembly (monitor screen not included) (see page 32)			
		Ring mount (used with light post) ..... 1300RM	950		
		Top post mount (used with monitor post) ..... 1300TPM	950		
		Monitor post (for top post mount) ..... 1331-997	260		
		<b>CHAIR, HYDRAULIC *</b> ..... 3900	<b>7390</b>		
		Additional touchpad ..... 3914-102	500		
		Additional foot control ..... 3914-010	500		
		Foot control in lieu of touchpad ..... 3914-010-CS	N/C		
		<b>Touchpad deduct chair mount (option on unit)</b> ..... 3914-081-CS	-500		
		Additional chair options (see pages 24-25 for pricing, write selection here) ...			
		Stool, doctor's deluxe / doctor's comfort ..... 6126C/6176	790/840		
		Stool, assistant's value w/Ultraleather™ upholstery ..... 6137C	905		
		Stool, assistant's deluxe / assistant's comfort ..... 6147C/6187	1150/1200		
		<b>* ITEM</b> included in AS SHOWN PACKAGE price		<b>TOTAL</b>	
		<b>Often Overlooked</b>			
Equipment dimensions, see page 73					
Shipping weight estimate:					
delivery unit & chair		95 lb / 44 kg			
base package		544 lb / 247 kg			
as shown package		562 lb / 255 kg			

# fixed chair mount without sidebox packages

## package F2 | value control head

### STANDARD FEATURES:

#### DELIVERY SYSTEM 4195FI

Automatic control for three handpieces

asepsis 3-way air/water syringe w/tubing, 5'  
handpiece tubing, 4-hole Midwest, asepsis style, 5'  
individual pinch valve control block

master on/off

individual handpiece pressure adjustment

individual water coolant adjustment

non-retraction water coolant

air coolant adjustment

handpiece flush

handpiece exhaust filter/oil collector

single handle brake release function

self-contained clean water system, w/quick  
disconnect single bottle, 1½ liter

wet/dry foot control

Flexible control arm with pneumatic brake

vertical travel 28" / horizontal travel 62"

Internal umbilical, 5' with 3' external from chair to utility center

Utility center

Contoured housing and cover

air pressure regulator with filter

master air shut-off valve

Chair adapter, fixed

#### CUSPIDOR, BASIC 1636FI/1637FI (1637FI option shown)

Removable ceramic bowl

Manual cup fill

Timed bowl rinse

Internal umbilical, 5' with 3' external from chair to utility center

Water pressure regulator with filter

Master water shut-off valve

Post mount, 2"

1637FI includes:

Individual assistant's instrumentation holders

Assistant's instrumentation

central vacuum canister w/disposable solids collector

asepsis 3-way air/water syringe w/tubing, 5'

autoclavable standard HVE lever valve w/tubing, 5'

autoclavable saliva ejector lever valve w/tubing, 5'

#### ASSISTANT'S INSTRUMENTATION 1075FI

(option not shown)

Central vacuum accessory package

central vacuum canister w/disposable solids collector

asepsis 3-way air/water syringe w/tubing, 5'

autoclavable standard HVE lever valve w/tubing, 5'

autoclavable saliva ejector lever valve w/tubing, 5'

Internal umbilical, 5' with 3' external from chair to utility center

Post mount, 2"

#### LIGHT, HALOGEN 9080

3-Axis light head movement

Automatic on/off activation with override

Color shade matching temperature setting 4200K to 5500K

Shadow free illumination

#### CHAIR, HYDRAULIC 3900

Ultraleather™ upholstery with wide back for elbow support

Double articulating pinch-style headrest

Thin back with full lumbar support

Right/left patient access

Swing-down locking armrests

Hydraulic cantilever mechanism

Emergency shut-off

30 degree base rotation with lock

Multifunction touchpad w/4 programmable preset positions

Seat height: minimum 17" maximum 29"



### POPULAR EQUIPMENT CHOICES



ASSISTANT'S  
INSTRUMENTATION 1075FI



CUSPIDOR 1636FI



MONITOR MOUNT  
1300RM



DOCTOR'S DELUXE  
STOOL 6126C



ASSISTANT'S VALUE  
STOOL 6137C



ASSISTANT'S DELUXE  
STOOL 6147C

		FACTORY OPTIONS	PART NO	RETAIL	PRICE
<b>BASE PACKAGE</b>	<b>\$ 13,430</b>				
4195FI Delivery System		<b>DELIVERY SYSTEM W/VALUE *</b> ..... 4195FI	<b>\$ 3890</b>		
3900 Chair		Handpiece illumination system			
9080 Light		one handpiece (5-hole), accommodates 1 add'l tubing ..... 1700-021-CS	595		
1110-820 Light Post		add'l handpiece illumination tubing (5-hole) ..... 0017-615-CS	265		
<b>AS SHOWN PACKAGE</b>	<b>\$ 15,100</b>	one handpiece (6-pin), accommodates 1 add'l tubing ..... 1700-019-CS	595		
Additionally equipped with:		add'l handpiece illumination tubing (6-pin) ..... 0017-603-CS	190		
1700-081-CS Unit Mount Touchpad		Additional handpiece position, wet ..... 1100-189-CS	405		
1637FI Basic Cuspidor		Additional handpiece position, air drive only ..... 1100-190-CS	195		
3914-081-CS Deduct Chair Mount Touchpad		Less handpiece position ..... 1700-099-CS	-100		
<b>DESIGNER FRIENDLY™ SURFACE TREATMENTS</b> <i>standard unless noted</i>		Heated syringe conversion kit (1) per syringe ..... 1107-450-CS	320		
Unit plastic finishes:		Assistant's Instrumentation (control head mount) ..... 1409-CS	400		
<input type="radio"/> Forest Grey	<input type="radio"/> White	Quick disconnect w/flow control 1/4" (for water) ..... 0016-454-CS	100		
<input type="radio"/> Greywood	<input type="radio"/> Laser Leather	Quick disconnect w/flow control 3/8" (for air) ..... 0016-456-CS	100		
Powder coat paint for unit & light:		Additional brake handle activation (left) ..... 0023-642-CS	140		
<input type="radio"/> Forest Grey	<input type="radio"/> Shadow	Swivel tray, arm mounted (over control head) ..... 1110-468	335		
<input type="radio"/> Pearl	<input type="radio"/> Charcoal	<b>Touchpad *</b> unit mount for 3900 chair ..... 1700-081-CS	<b>500</b>		
<input type="radio"/> Sage	<input type="radio"/> Bronze	Chrome accents ..... 9000-998-CS	135		
Powder coat paint for chair:		Deluxe utility center housing and cover in lieu of contoured ..... 1115-007-CS	-200		
<input type="radio"/> Forest Grey	<input type="radio"/> Shadow	Additional umbilical (per foot) ..... 1115-999-CS	30		
Add \$1,500:		Less self-contained clean water system, w/quick disconnect,			
<input type="radio"/> White	<input type="radio"/> Pearl	single bottle, add water shutoff and regulator ..... 1700-355-CS	N/C		
<input type="radio"/> Pewter	<input type="radio"/> Sage	Routing valve kit w/water shutoff & regulator/factory installation ..... 1115-315-CS	235		
<input type="radio"/> Arctic Silver		Routing valve kit for use with cuspidor/factory installation ..... 1115-316-CS	60		
<input type="radio"/> Arctic Silver		Water shut-off valve & regulator kit ..... 1115-298	160		
Ultraleather™ upholstery color choices for chair & deluxe stools:		<b>Ancillaries</b> , uses existing handpiece position (see pages 62-63)			
<input type="radio"/> Baltic	<input type="radio"/> Dove Grey	Write in selection here .....			
<input type="radio"/> Diplomat Blue	<input type="radio"/> Doe	Write in selection here .....			
<input type="radio"/> Papyrus	<input type="radio"/> Parrot	<b>CUSPIDOR, BASIC *</b> w/assistant's instrumentation ..... 1637FI	<b>1670</b>		
<input type="radio"/> Sassy	<input type="radio"/> Stone	or BASIC CUSPIDOR ..... 1636FI	1220		
<input type="radio"/> Pearlized Mica	<input type="radio"/> Pearlized Green Tea	or ASSISTANT'S INSTRUMENTATION ..... 1075FI	840		
<input type="radio"/> Pearlized Midnight	<input type="radio"/> Pearlized Pewter	Additional autoclavable standard HVE w/tubing,			
Chair, raffia or special color Ultraleather™	3907-999-SP add \$280	5" & autoclavable oral dry cup (1075FI/1637FI) ..... 1400-113-CS	85		
<input type="radio"/> Durango Dust raffia	<input type="radio"/> Noche raffia	Basic cuspidor automatic cup fill ..... 1700-055-CS	195		
<input type="radio"/> Deluxe stool, special color Ultraleather™	6007-999-SP add \$100	<b>LIGHT *</b> Halogen (mounting options, see pages 27-30) ..... 9080	<b>1890</b>		
Upholstery choices (must choose one):		or Halogen ..... 9045	2390		
<input type="radio"/> Standard Back, wide/thin	Standard / no charge	or LED (coming summer 2013) ..... 9073	3490		
<input type="radio"/> Memory Comfort (plush only)	3907-999-MC add \$800	<b>Light post, 35" *</b> ..... 1110-820	<b>260</b>		
<input type="radio"/> Narrow Back	3907-999-NB add \$100	Monitor mount assembly (monitor screen not included) (see page 32)			
<input type="radio"/> Narrow Back with slings	3907-999-NS add \$400	Ring mount (used with light post) ..... 1300RM	950		
<input type="radio"/> Dual color (plush only)	3907-999-DC add \$280	Top post mount (used with monitor post) ..... 1300TPM	950		
<input type="radio"/> Raffia dual color (plush only)	3907-999-RD add \$560	Monitor post (for top post mount) ..... 1331-997	260		
Upholstery styles (must choose one):		<b>CHAIR, HYDRAULIC *</b> ..... 3900	<b>7390</b>		
<input type="radio"/> Plush	Standard / no charge	Additional touchpad ..... 3914-102	500		
<input type="radio"/> Asepsis	3907-999-AS no charge	Additional foot control ..... 3914-010	500		
Equipment dimensions, see page 73		Foot control in lieu of touchpad ..... 3914-010-CS	N/C		
Shipping weight estimate:		<b>Touchpad *</b> deduct chair mount (option on unit) ..... 3914-081-CS	<b>- 500</b>		
delivery unit & chair	95 lb / 44 kg	Additional chair options (see pages 24-25 for pricing, write selection here) ...			
base package	554 lb / 247 kg	Stool, doctor's deluxe / doctor's comfort ..... 6126C/6176	790/840		
as shown package	564 lb / 256 kg	Stool, assistant's value w/Ultraleather™ upholstery ..... 6137C	905		
		Stool, assistant's deluxe / assistant's comfort ..... 6147C/6187	1150/1200		
<b>* ITEM included in AS SHOWN PACKAGE price</b>					<b>TOTAL</b>
<b>Often Overlooked</b>					





3900 Patient Chair



6176 Doctor's Comfort Stool

## Chairs & Stools

### Affordable Excellence

#### 3900 Patient Chair

- Material: cast aluminum chair back, chair base & utility center base
- Benefit: strong, impervious to fading or cracking
- Material: pressure-formed plastic chair & utility center covers
- Benefit: warp, impact & stain resistant
- Material: aluminum hydraulic manifold & speed adjustment screws
- Benefit: independent speed adjustment of base & back both up & down
- Material: optional soft-touch back massage with integrated compressor
- Benefit: saves central compressor capacity for delivery systems

#### 6176 Doctor's Comfort Stool

- Material: tension woven material with rocking back
- Benefit: eliminates pressure points, massages & supports the seat & back
- Material: all-surface 22mm polyurethane casters
- Benefit: smooth rolling eliminates squeaks / long-life




forest dental

[forestdental.com](http://forestdental.com)

[sales@forestdental.com](mailto:sales@forestdental.com)



## Ultraleather™



Memory Comfort Option / Pearlized Mica Ultraleather™



Standard Plush / Grenadine Ultraleather™

## Designer Friendly™

- Material: memory foam in Memory Comfort chair upholstery option
- Benefit: designed for all-day procedures keeping patient relaxed
- Material: ultraleather™ upholstery standard
- Benefit: withstands heavy use & cleaning with glove-soft leather feel
- Material: forest grey or shadow powder coat standard on 3900 chair / seven optional Designer Friendly™ colors
- Benefit: complements treatment room interior design color palette



MODEL 6126C  
\$ 790



MODEL 6137C  
\$ 905



MODEL 6147C  
\$ 1150



MODEL 6176  
\$ 840



MODEL 6187  
\$ 1200

**STANDARD FEATURES**

**DOCTOR'S DELUXE 6126C**

Ultraleather™ upholstery  
5-leg base  
Twin wheel hard floor casters  
Pneumatic cylinder height adjustment, 17 1/4" – 23 3/4"  
Backrest w/height adjustment  
Tilt seat lever function  
Contoured seat  
Accommodates maximum of 250 lbs.

**ASSISTANT'S STANDARD 6137C**

Ultraleather™ upholstery  
5-leg base  
Twin wheel hard floor casters  
Pneumatic cylinder height adjustment, 22 1/2" – 31 1/4"  
Contoured seat  
Torso support w/height adjustment & release lever  
Foot ring w/height adjustment  
Accommodates maximum of 250 lbs.

**ASSISTANT'S DELUXE 6147C**

*Same as 6137C stool with additional features:*  
Backrest w/height adjustment

**DOCTOR'S COMFORT 6176**

Air circulating weave construction  
5-leg base  
Twin wheel hard floor casters  
Pneumatic cylinder height adjustment, 17 1/4" – 23 3/4"  
Triple lever adjustment - independent back and seat  
Backrest w/angle & height adjustment  
Active tilt™ seat (actively moves with user)  
Tilt seat lever function  
Contoured seat  
Accommodates maximum of 300 lbs.

**ASSISTANT'S COMFORT 6187**

*Same as 6176 stool with additional features:*  
Pneumatic cylinder height adjustment, 22 1/2" – 31 1/4"  
Fully adjustable ratcheting assistant arm  
Foot ring w/height adjustment

Shipping weight:

6126C 34 lb / 16 kg

6137C 33 lb / 15 kg

6147C 42 lb / 19 kg

6176 48 lb / 15 kg

6187 60 lb / 19 kg

**FACTORY OPTIONS**

STOOL .....	MODEL CHOICE _____	PRICE: _____
Vinyl in lieu of Ultraleather™ 6126C/6137C/6147C* .....	6007-999-CS	\$ – 100 _____
Special color upholstery 6126C/6137C/6147C (upcharge) .....	6007-999-SP	100 _____
Twin wheel carpet casters .....	0010-401-CS	N/C _____
Cylinder choices (upcharge)		
Doctor's standard cylinder, 16 1/4" – 20 1/2" .....	6000-299-CS	N/C _____
Doctor's X-long cylinder, 19 1/2" – 27" (foot ring w/support highly recommended) .....	6000-297-CS	N/C _____
Assistant's short cylinder, 21 1/4" – 27" .....	6000-303-CS	N/C _____
Foot ring w/support .....	6000-575	115 _____

**TOTAL** \_\_\_\_\_

\*Specify vinyl manufacturer and color (Not all colors may be available. Call Forest for ordering details.)

**ULTRALEATHER™ UPHOLSTERY standard unless noted**


Ultraleather™:	<input type="radio"/> Baltic	<input type="radio"/> Papyrus	<input type="radio"/> Charcoal	<input type="radio"/> Parrot
	<input type="radio"/> Diplomat Blue	<input type="radio"/> Raven Wing	<input type="radio"/> Doe	<input type="radio"/> Sassy
	<input type="radio"/> Stone	<input type="radio"/> Grenadine	<input type="radio"/> Walnut	<input type="radio"/> Dove Grey
	<input type="radio"/> Pearlized Midnight	<input type="radio"/> Pearlized Pewter	<input type="radio"/> Pearlized Green Tea	<input type="radio"/> Pearlized Mica

# chairs + stools

## 3900 hydraulic patient chair



MODEL 3900 (110V)  
\$ 7390



### STANDARD FEATURES

#### HYDRAULIC CHAIR 3900

Ultraleather™ upholstery with wide back for elbow support  
Double articulating pinch-style headrest  
Thin back with full lumbar support  
Right/left patient access  
Swing-down locking armrests  
Hydraulic cantilever mechanism  
Emergency shut-off  
30 degree base rotation with lock  
Multifunction touchpad  
w/4 programmable preset positions, accommodates 3 additional controls  
Seat height: minimum 17" maximum 29"

### FACTORY OPTIONS

PART NO	RETAIL	PRICE
CHAIR (Forest Grey or Shadow) .....	3900	\$ 7390
CHAIR (Designer Friendly) .....	3900C	8890
Additional touchpad .....	3914-102	500
Foot control in lieu of touchpad .....	3914-010-CS	N/C
Additional foot control .....	3914-010	500
Touchpad deduct chair mount (option on unit) .....	3914-081-CS	- 500
Chair adapter, fixed (for external umbilical) .....	1115-195	430
Chair adapter, fixed (for internal umbilical) .....	1115-530	565
Chair adapter, fixed with sidebox (for internal umbilical) .....	1115-531	565
Chair adapter, pivot (for internal umbilical) .....	1115-528	565
Contoured utility center housing & cover .....	3900-101	420
UPHOLSTERY OPTIONS		
Write in selection here (see page 25)		
Write in selection here (see page 25)		
Pedo booster seat .....	3907-180	280
Soft-touch back massage option (only available on wide/thin back) .....	3900-100	540
IV Board .....	3900-200	640
LIGHT POSTS (for unit mount lights/post mount monitor assemblies see page 60)		
Fixed (chair adaptor sold separately)		
44" light post without slot, chair mount .....	1110-794	200
44" light post w/slot for monitor cord .....	1110-050	200
Pivot (includes chair adaptor)		
Pivot light post, 3900 (light not included) .....	1112-644-3	1065
Pivot monitor post, 3900 (monitor mount assembly not included) ..	1331-644-3	1065
Monitor mount assembly, top post mount (see page 32) .....	1300TPM	950
Monitor mount assembly, ring mount (see page 32) .....	1300RM	950
<b>TOTAL</b>		

### ULTRALEATHER™ UPHOLSTERY standard unless noted

Ultraleather™ upholstery color choices for chair:

- |  |   |                                  |
|--|---|----------------------------------|
| <input type="radio"/> Baltic             | <input type="radio"/> Dove Grey           | <input type="radio"/> Charcoal   |
| <input type="radio"/> Diplomat Blue      | <input type="radio"/> Doe                 | <input type="radio"/> Grenadine  |
| <input type="radio"/> Papyrus            | <input type="radio"/> Parrot              | <input type="radio"/> Raven Wing |
| <input type="radio"/> Sassy              | <input type="radio"/> Stone               | <input type="radio"/> Walnut     |
| <input type="radio"/> Pearlized Mica     | <input type="radio"/> Pearlized Green Tea |                                  |
| <input type="radio"/> Pearlized Midnight | <input type="radio"/> Pearlized Pewter    |                                  |

Raffia or special color Ultraleather™ 3907-999-SP add \$280:  
 Durango Dust     Noche     Special color

Upholstery choices (must choose one):

- |   |                           |
|---|---------------------------|
| <input type="radio"/> Standard wide/thin back   | Standard / no charge      |
| <input type="radio"/> Upgrade upholstery option | (see page 25 for pricing) |

Upholstery styles (must choose one):

- |                               |                         |
|-------------------------------|-------------------------|
| <input type="radio"/> Plush   | Standard / no charge    |
| <input type="radio"/> Asepsis | 3907-999-AS (no charge) |

Shipping weight: 335 lb / 152 kg



### DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Powder coat paint:	<input type="radio"/> Forest Grey	<input type="radio"/> Shadow			
Add \$1,500:	<input type="radio"/> White	<input type="radio"/> Pearl	<input type="radio"/> Charcoal	<input type="radio"/> Pewter	<input type="radio"/> Sage



Narrow Back or  
Narrow Back with Slings



#### STANDARD FEATURES

##### **MEMORY COMFORT UPHOLSTERY**

Designed for all day/ultra-long procedures or greater general patient comfort  
3 layers of foam topped with 1.5 inches of memory foam for extra comfort  
Side-bolsters on seat cushion support legs  
Memory foam headrest

##### **NARROW BACK**

Perfect for Pediatric or Orthodontic dentistry  
Unrestricted patient access

##### **NARROW BACK WITH SLINGS**

*Same as above, with additional features:*  
Slings provide alternative for additional patient elbow support  
Slings create a barrier between patient and practitioner

#### FACTORY OPTIONS

##### UPHOLSTERY UPGRADE OPTION (add pricing to page 24)

Memory Comfort (plush only, not available with back massage) .....	3907-999-MC	\$800
Narrow Back (not available with back massage) .....	3907-999-NB	100
Narrow Back with Slings (not available with back massage) .....	3907-999-NS	400

##### UPHOLSTERY COLOR OPTIONS (add pricing to page 24)

Raffia or special color Ultraleather™ .....	3907-999-SP	280
Dual color (plush only) .....	3907-999-DC	280
Raffia dual color (plush only) .....	3907-999-RD	560

#### PART NO

#### RETAIL

Shipping weight: 335 lb / 152 kg



# chairs + stools

## office stool



MODEL 6198

\$ 940

### STANDARD FEATURES

#### OFFICE COMFORT

- Air circulating weave construction
- 5-leg base
- Twin wheel hard floor casters
- Pneumatic cylinder height adjustment, 17 1/4" – 23 3/4"
- Triple lever adjustment - independent back and seat
- Backrest w/angle & height adjustment
- Active tilt™ seat (actively moves with user)
- Tilt seat lever function
- Contoured seat
- Fully adjustable office arms
- Accommodates maximum of 300 lbs.

### FACTORY OPTIONS

	PART NO	RETAIL	PRICE
STOOL .....	6198	940	_____
Twin wheel carpet casters .....	0010-401-CS	N/C	_____
► Cylinder choices (upcharge)			
Standard cylinder, 16 1/4" – 20 1/2" .....	6000-299-CS	N/C	_____
X-long cylinder, 19 1/2" – 27" .....	6000-297-CS	N/C	_____
		<b>TOTAL</b>	_____

Shipping weight:  
6198 54 lb / 18 kg

9080/9090 Halogen  
9045/9065 Halogen  
9073/9074 LED



9080/9090 Halogen



9045/9065 Halogen



9073/9074 LED

## Lights & Monitor Mounts

### Affordable Excellence

#### 9080/9090 Halogen

- Material: reflector material imbedded into glass / lexan light shield
- Benefit: reflector will not peel or crack / shield is scratch resistant
- Material: oil impregnated bronze bushing / 3-axis light head joint
- Benefit: no light drift / 3-axis light keeps light pattern parallel with the long axis of mouth while freeing up space above patient's head

#### 9045/9065 Halogen

- Material: aluminum light-defeat shield
- Benefit: foot candles (brightness) can be adjusted while maintaining 5500 kelvin / ideal for shade matching & color selection
- Material: quiet electric fan
- Benefit: keeps halogen bulb, wiring & reflector cool for long life

#### 9073/9074 LED

- Material: six-mode, polycarbonate-button keyboard
- Benefit: easily controlled with soft-touch buttons / resistant to cleaners
- Material: eight 1-watt led lamps
- Benefit: creates a light pattern with feathered edges reducing eye fatigue / consumes approximately 75 percent less power than halogen lights / no heat / 22 year projected LED life span



forest dental

[forestdental.com](http://forestdental.com)

[sales@forestdental.com](mailto:sales@forestdental.com)



## Powder Coat

arctic silver

bronze

charcoal

forest grey

pearl

pewter

sage

shadow

white



9080/9090 Halogen (Bronze Powder Coat)



9073/9074 LED (Sage Powder Coat)

## Designer Friendly™

- Material: cast aluminum light head (9080/9090/9073/9074)
- Benefit: allows for nine Designer Friendly™ powder coat paint colors

# lights + monitor mount assemblies

## post mount lights



MODEL 9080 (110V)  
HALOGEN LIGHT  
**\$ 1890**



MODEL 9045 (110V)  
HALOGEN LIGHT  
**\$ 2390**



MODEL 9073  
LED LIGHT  
**\$ 3490**  
COMING SUMMER 2013

### STANDARD FEATURES

#### LIGHTS

##### MODEL 9080

Automatic on/off activation with override  
3-Axis light head movement  
Additional bulb  
Color temperature, 4,200 to 5,500 Kelvin  
Designer Friendly™ powder coat paint light head

##### MODEL 9045

Automatic on/off activation with override  
Cooling fan  
Easy release autoclavable handles  
Additional bulb  
Color temperature, 5,000 Kelvin, constant  
Two-tone plastic light head

##### MODEL 9073

Automatic on/off activation with override  
Easy release autoclavable handles  
Light source 8 LED lamps  
Color temperature, 4,200 to 6,000 Kelvin  
Designer Friendly™ powder coat paint light head

### FACTORY OPTIONS

LIGHT	MODEL CHOICE	PRICE:	
Replacement bulb (9045/9080/9065/9090)	0014-300	\$ 35	
Short rigid arm 20" in lieu of standard 32"	1110-227-CS	N/C	
Light posts:			
44" w/slot for sidebox access	1110-900	200	
44" without slot, chair mount	1110-794	200	
35" w/post clamp	1110-820	260	
Light posts, with monitor cord slot:			
44" w/dual slot for sidebox access	1110-170	200	
44" w/single slot, chair mount	1110-050	200	
35" w/single slot & post clamp	1110-045	260	

**TOTAL**

Shipping weight: 35 lb / 16 kg

### DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Powder coat paint finishes:  Forest Grey  Shadow  White  Pearl  Charcoal  
 Pewter  Sage  Bronze  Arctic Silver

# lights + monitor mount assemblies

## track mount lights

Shown with optional 1300T monitor mount assembly (monitor screen not included)

MODEL 1300T  
MONITOR MOUNT ASSEMBLY  
(track not included)  
**\$ 2290**

MODEL 9080T (110V)  
HALOGEN LIGHT/TRACK  
**\$ 3690**

MODEL 9045T (110V)  
HALOGEN LIGHT/TRACK  
**\$ 4190**

MODEL 9073T  
LED LIGHT/TRACK  
**\$ 5290**  
COMING SUMMER 2013

STANDARD FEATURES	FACTORY OPTIONS	PART NO	RETAIL	PRICE
<b>LIGHT + TRACK ASSEMBLY</b> Pre-wired through flex arm Track/Trolley assembly Vertical light tube with 6 inch horizontal side offset for clearance of optional monitor Flex arm & light head assembly Trolley mounting plate and shroud Drop post (specify mounting height)	LIGHT / TRACK ..... MODEL CHOICE _____ Replacement bulb (9045/9080/9065/9090) ..... 0014-300 Monitor mount assembly (track not included) ..... 1300T S Video, cable pre-installed ..... 0014-282 RCA cables (3 conductor), pre-installed ..... 0014-283		PRICE: \$ 35 2290 40 40	_____ _____ _____ _____
			<b>TOTAL</b>	_____
<b>MONITOR MOUNT ASSEMBLY</b> Pre-wired from track through flex arm Coax cable Low voltage wiring (2 conductor 18 awg wire) VGA shielded cable ICW Quicklink monitor mount Rex positioning handles  Monitor Requirements: 12 volt with power supply weight: up to 18 lbs max, width of 19"	► **SPECIFY MOUNTING HEIGHT** _____  ► Often Overlooked			
Shipping weight: Light/track 115 lb / 53 kg 1300T 140 lb / 64 kg (includes mount/track light weight)	► DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted  Powder coat paint finishes: <input type="radio"/> Forest Grey <input type="radio"/> Shadow <input type="radio"/> White <input type="radio"/> Pearl <input type="radio"/> Charcoal <input type="radio"/> Pewter <input type="radio"/> Sage <input type="radio"/> Bronze <input type="radio"/> Arctic Silver			

# lights + monitor mount assemblies

## ceiling mount lights



MODEL 9080-C  
HALOGEN LIGHT  
\$ 2490



MODEL 9045-C  
HALOGEN LIGHT  
\$ 2990



MODEL 9073-C  
LED LIGHT  
\$ 4090  
COMING SUMMER 2013

STANDARD FEATURES	FACTORY OPTIONS	PART NO	RETAIL	PRICE
<b>LIGHTS</b> Ceiling mounting plate and shroud with drop post <b>►</b> specify mounting height  <b>MODEL 9080-C</b> Automatic on/off activation with override 3-Axis light head movement Additional bulb Color temperature, 4,200 to 5,500 Kelvin Designer Friendly™ powder coat paint light head  <b>MODEL 9045-C</b> Automatic on/off activation with override Cooling fan Easy release autoclavable handles Additional bulb Color temperature, 5,000 Kelvin, constant Two-tone plastic light head  <b>MODEL 9073-C</b> Automatic on/off activation with override Easy release autoclavable handles Light source 8 LEDs Color temperature, 4,200 to 6,000 Kelvin Designer Friendly™ powder coat paint light head	LIGHT ..... MODEL CHOICE _____ Replacement bulb (9045/9080/9065/9090) ..... 0014-300  <b>► **SPECIFY MOUNTING HEIGHT**</b> _____  <b>► Often Overlooked</b>			PRICE: \$ 35  <b>TOTAL</b> _____
Shipping weight: 50 lb / 23 kg	<b>► DESIGNER FRIENDLY™ SURFACE TREATMENTS</b> standard unless noted  Powder coat paint finishes: <input type="radio"/> Forest Grey <input type="radio"/> Shadow <input type="radio"/> White <input type="radio"/> Pearl <input type="radio"/> Charcoal <input type="radio"/> Pewter <input type="radio"/> Sage <input type="radio"/> Bronze <input type="radio"/> Arctic Silver			

# lights + monitor mount assemblies

## cabinet + wall mount lights



HALOGEN CABINET  
MODEL 9065C (standard)  
9065C-T (tall)  
**\$ 3490**



HALOGEN WALL  
MODEL 9065W (standard)  
9065W-T (tall)  
**\$ 3490**

COMING SUMMER 2013



LED CABINET  
MODEL 9074C (standard)  
9074C-T (tall)  
**\$ 4590**

HALOGEN CABINET  
MODEL 9090C (standard)  
9090C-T (tall)  
**\$ 2990**

HALOGEN WALL  
MODEL 9090W (standard)  
9090W-T (tall)  
**\$ 2990**

COMING SUMMER 2013



LED WALL  
MODEL 9074W (standard)  
9074W-T (tall)  
**\$ 4590**

### STANDARD FEATURES

### FACTORY OPTIONS

### PART NO

### PRICE

#### LIGHTS

Cabinet mount, internal vertical mount with post  
Wall mount, vertical mounting plate  
Standard rigid arm  
Tall rigid arm (additional 1' drop)

#### MODEL 9090

Automatic on/off activation with override  
3-Axis light head movement  
Additional bulb  
Color temperature, 4,200 to 5,500 Kelvin  
Designer Friendly™ powder coat paint light head

#### MODEL 9065

Automatic on/off activation with override  
Cooling fan  
Easy release autoclavable handles  
Additional bulb  
Color temperature, 5,000 Kelvin, constant  
Two tone plastic light head

#### MODEL 9074

Automatic on/off activation with override  
Easy release autoclavable handles  
Light source 8 LED lamps  
Color temperature, 4,200 to 6,000 Kelvin  
Designer Friendly™ powder coat paint light head

#### LIGHT .....

MODEL CHOICE .....  
Replacement bulb (9045/9080/9065/9090) .....

#### PRICE:

\$ 35

Cabinet mounting options: (for dimensions see page 76)

Horizontal top cabinet mount in lieu of internal mount .....

N/C

Vertical exterior cabinet mount in lieu of internal mount .....

N/C

#### TOTAL

#### ► \*\*SPECIFY CABINET/WALL MOUNTING HEIGHT\*\*

► Often Overlooked

#### ► DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Shipping weight:  
Cabinet 60 lb / 28 kg Wall 65 lb / 30 kg

Powder coat paint finishes:  Forest Grey  Shadow  White  Pearl  Charcoal  
 Pewter  Sage  Bronze  Arctic Silver

# lights + monitor mount assemblies

## cabinet + wall monitor mount assemblies

(monitor screens not included)



### CABINET MOUNT

MODEL 1310C (standard)

1310C-T (tall)

**\$ 2595**

### WALL MOUNT

MODEL 1310W (standard)

1310W-T (tall)

**\$ 2595**

#### STANDARD FEATURES

##### MONITOR MOUNT ASSEMBLY

Cabinet mount, internal vertical mount with post  
Wall mount, vertical mounting plate  
Standard rigid arm  
Tall rigid arm (additional 1' drop)  
Coax cable  
Low voltage wiring (2 conductor 18 awg wire)  
VGA shielded cable  
ICW Quicklink monitor mount  
Rex positioning handles

##### Monitor Requirements:

12 volt with power supply  
weight: 8 1/4 to 13 lbs max, width of 19"

#### FACTORY OPTIONS

MONITOR MOUNT ASSEMBLY ..... MODEL CHOICE \_\_\_\_\_  
S Video, cable pre-installed ..... 0014-282  
RCA cables (3 conductor), cable pre-installed ..... 0014-283

#### TOTAL

#### PRICE

#### ► \*\*SPECIFY CABINET/WALL MOUNTING HEIGHT\*\*

► Often Overlooked

#### ► DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Shipping weight:  
25 lb / 12 kg

Powder coat paint finishes:  Forest Grey  Shadow  White  Pearl  Charcoal  
 Pewter  Sage  Bronze  Arctic Silver

# lights + monitor mount assemblies

## ceiling + post monitor mount assemblies



MODEL 1300-C  
CEILING MOUNT  
\$ 2290

(monitor screens and chair mount posts not included)



MODEL 1300RM  
RING POST MOUNT  
\$ 950



MODEL 1300TPM  
TOP POST MOUNT  
\$ 950



ORGANIZER BIN  
MOUNT 1300-025-CS  
\$ 395

used with rear delivery systems and carts with storage organizer bin

### STANDARD FEATURES

### FACTORY OPTIONS

### PART NO

### PRICE

#### **MONITOR MOUNT ASSEMBLY**

##### MONITOR MOUNT ASSEMBLY, CEILING MOUNT 1300-C

Ceiling mounting plate and shroud  
Drop post (specify mounting height)  
Coax cable  
Low voltage wiring (2 conductor 18 awg wire)  
VGA shielded cable  
ICW Quicklink monitor mount  
Rex positioning handles

##### MONITOR MOUNT ASSEMBLY, RING POST MOUNT 1300RM

ICW Quicklink monitor mount  
Rex positioning handles  
22" ring mount rigid arm

##### MONITOR MOUNT ASSEMBLY, TOP POST MOUNT 1300TPM

ICW Quicklink monitor mount  
Rex positioning handles  
22" top post mount rigid arm

##### ORGANIZER BIN MOUNT 1300-025-CS

Forest's unique locking tilt mechanism with knob  
1/2" pin

#### Monitor Requirements:

12 volt with power supply  
weight: 8 1/4 to 13 lbs max, width of 19"

#### Shipping weight:

1300-C	25 lb / 12 kg
1300RM, 1300TPM	15 lb / 7 kg
1300-025-CS	5 lb / 2 kg

MONITOR MOUNT ASSEMBLY ..... MODEL CHOICE \_\_\_\_\_

PRICE :

S Video, cable pre-installed (1300-C) ..... 0014-282

\$ 40

RCA cables (3 conductor), cable pre-installed (1300-C) ..... 0014-283

40

#### ► LIGHT POST (1300RM)

44" w/dual slot, fixed sidebox mount, w/Forest light bushing ..... 1110-170

200

35" w/slot & post clamp, fixed mount, w/Forest light bushing ..... 1110-045

260

Pivot light posts & requirements ..... (call to order)

#### ► MONITOR POST (1300TPM)

30" w/dual slot, fixed sidebox mount, w/oilite bushing ..... 1331-998

200

30" w/slot & post clamp, fixed mount, w/oilite bushing ..... 1331-997

260

Pivot monitor posts & requirements ..... (call to order)


#### TOTAL

#### ► \*\*SPECIFY MOUNTING HEIGHT for ceiling mount 1300-C\*\*

#### ► Often Overlooked

#### ► DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Powder coat paint finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> Shadow	<input type="radio"/> White	<input type="radio"/> Pearl	<input type="radio"/> Charcoal
	<input type="radio"/> Pewter	<input type="radio"/> Sage	<input type="radio"/> Bronze	<input type="radio"/> Arctic Silver	



Rear Mount, Side Delivery 7000 Series

## Delivery Systems

### Affordable Excellence

#### Rear Mount, Side Delivery 7000 Series

- Material: north america's only all metal control head chassis
- Benefit: resistant to UV & cleaners
- Material: hinged solid surface lid
- Benefit: long lasting surface gives easy access to ancillary products
- Material: exclusive brass handpiece lockout valves
- Benefits: allows assistant to remove a handpiece from its holder to change a burr while dentist continues working
- Material: five steel mounting configurations
- Benefit: wall, under cabinet, under counter, floor mount & cart
- Material: secondary aluminum work surface
- Benefit: dedicated work surface for burr block, electric handpiece attachment, mirror, explorer
- Material: steel height adjust mechanism
- Benefit: adjustable work surface & handpiece delivery arm height



forest dental

[forestdental.com](http://forestdental.com)

[sales@forestdental.com](mailto:sales@forestdental.com)



## Powder Coat

arctic silver

bronze

charcoal

forest grey

pearl

pewter

sage

shadow

white

laser leather

lemons

greywood



Rear Mount, Side Delivery 7000 Series (Pearl Powder Coat)



Deluxe Cuspidor (Black)

## Designer Friendly™

- Material: white & forest grey, three carbon fiber-look plastic finishes & optional chrome accents
- Benefit: rich design enhances overall practice image & patient experience
- Material: nine designer friendly powder coat paint colors including white & forest grey
- Benefit: dental equipment with a WOW factor
- Material: black & white porcelain deluxe cuspidor bowl with chrome accents
- Benefit: shape & color enhances aesthetics of a typically utilitarian device

**delivery systems**page 34

---

**4485SE Standard IC / 4987SE Euro / 4195SE Value  
fixed chair mount with sidebox, external umbilical****4485FE Standard IC / 4195FE Value  
fixed chair mount without sidebox, external umbilical****cuspidors / assistant's instrumentation**page 38

---

**1636 / 1637 Basic Cuspidor / 1075 Vacuum Telescoping Arm  
fixed chair post mount, external umbilical****support systems**page 39

---

**1800RM / 1800TPM  
flex tray****rear + side delivery systems**page 40

---

**M&CC PFR8144 / M&CC DR Cabinet  
cabinetry mounting configurations****7000BC / 7000FM / 7000UC / 7000WM  
7000 Series: rear mount, side delivery****3386WM / 3486WM / 3686WM  
rear wall mount, side delivery****3386BC / 3486BC / 3686BC  
rear beneath counter mount, side delivery****3230 / 3235  
side cabinet & wall mount - standard IC control****3730 / 3735  
side cabinet & wall mount - value control****2400 / 2700  
rear or side beneath counter mount****1280 / 2180  
rear beneath counter mount****optional delivery systems**page 50

---

**2500 / 2505 / 2510 / 2515  
panel mounts****MT600 / MT610  
compact controls****1055  
beneath counter assistant's instrumentation**

# delivery systems without chair


## fixed chair mount with sidebox, external umbilical



MODEL	4485SE Standard IC	4987SE Euro	4195SE Value
PRICE	\$ 7145	\$ 8145	\$ 5945
STANDARD FEATURES	pneumatic, 3 HP, syringe	pneumatic, 4 HP, syringe	pinch, 3 HP, syringe

**DELIVERY SYSTEMS** can be mounted on any brand of chair

Choose your control head



### ADDITIONAL EQUIPMENT:

Choose one:

- 1576 Add \$ 1450  
Cuspidor, Deluxe  
w/swivel assistant's quad  
holder bar in lieu of  
telescoping arm



- 1577 Add \$ 1450  
Cuspidor, Deluxe



- 1678 Add \$ 810  
Cuspidor, Basic  
w/individual assistant's  
instrumentation  
holders in lieu  
of telescoping arm



- 1679 Add \$ 810  
Cuspidor, Basic



## STANDARD FEATURES

## FACTORY OPTIONS

## PART NO

## RETAIL

## PRICE

**DELIVERY SYSTEM**

Choice of Control Head  
 Automatic control for handpieces  
 asepsis 3-way air/water syringe w/tubing, 5'  
 (Euro tubing 3')  
 handpiece tubing, 4-hole Midwest, asepsis style, 5'  
 (Euro tubing 3')  
 master control block (with pneumatic or pinch valves)  
 master on/off  
 individual handpiece pressure adjustment  
 individual water coolant adjustment  
 non-retraction water coolant  
 air coolant adjustment  
 handpiece flush  
 handpiece exhaust filter/oil collector  
 single handle brake release function  
 dual handle brake release function (Euro)  
 swivel instrument tray (Euro)  
 wet/dry foot control  
 Flexible control arm with pneumatic brake  
 vertical travel 28" / horizontal travel 62"  
 Sidebox with self-contained clean water system, w/quick disconnect single bottle, 1½ liter  
 electrical outlet/UL  
 Telescoping arm with swivel quad holder bar  
 Central vacuum accessory package  
 central vacuum canister w/disposable solids collector  
 asepsis 3-way air/water syringe w/tubing, 5'  
 autoclavable standard HVE lever valve w/tubing, 5'  
 autoclavable saliva ejector lever valve w/tubing, 5'  
 External umbilical, 5'  
 Deluxe utility center  
 air pressure regulator with filter  
 master air shut-off valve  
 Chair adapter, fixed (specify chair make/model)

**CUSPIDOR, DELUXE 1577/1576**

One-piece ceramic bowl w/touch lever water activator  
 Automatic cup fill  
 Timed bowl rinse  
 Water pressure regulator with filter  
 Master water shut-off valve  
 1576 includes:  
 Swivel assistant's quad holder bar in lieu of telescoping arm

**CUSPIDOR, BASIC 1678/1679**

Removable ceramic bowl  
 Manual cup fill  
 Timed bowl rinse  
 Water pressure regulator with filter  
 Master water shut-off valve  
 1678 includes:  
 Individual assistant's instrumentation holders in lieu of telescoping arm

**DELIVERY SYSTEM** ..... MODEL CHOICE \_\_\_\_\_

Handpiece illumination system  
 one handpiece (5-hole), accommodates 1 add'l tubing ..... 1700-021-CS  
 add'l handpiece illumination tubing (5-hole) ..... 0017-615-CS  
 one handpiece (6-pin), accommodates 1 add'l tubing ..... 1700-019-CS  
 add'l handpiece illumination tubing (6-pin) ..... 0017-603-CS  
 Additional handpiece position, wet ..... 1100-189-CS  
 Additional handpiece position, air drive only ..... 1100-190-CS  
 Less handpiece position ..... 1700-099-CS  
 Heated syringe conversion kit (1) per syringe ..... 1107-450-CS  
 Additional brake handle activation (left) (euro standard) ..... 0023-642-CS  
 Swivel tray, arm mounted (euro standard) ..... 1110-468  
 Additional autoclavable standard HVE  
 w/tubing, 5' & autoclavable oral dry cup ..... 1400-113-CS  
 Quick disconnect w/flow control 1/4" (for water) ..... 0016-454-CS  
 Quick disconnect w/flow control 3/8" (for air) ..... 0016-456-CS  
 Less telescoping arm and assistant's instrumentation ..... 1400-000-CS  
 Chrome accents ..... 9000-998-CS  
 Standard utility center in lieu of deluxe ..... 1115-001-CS  
 Additional umbilical (per foot) ..... 1115-999-CS  
 Less self-contained clean water system, w/quick disconnect,  
 single bottle, add water shutoff and regulator ..... 1700-355-CS  
 Routing valve kit w/water shutoff & regulator/factory installation ..... 1115-315-CS  
 Routing valve kit for use with cuspidor/factory installation ..... 1115-316-CS  
 Water shut-off valve & regulator kit ..... 1115-298  
 Without chair adapter ..... 1115-099-CS

Ancillaries, uses existing handpiece postion (see pages 62-63)

Write in selection here .....

Write in selection here .....

**CUSPIDOR**

Deluxe w/swivel assistant's quad holder bar in lieu of  
 telescoping arm ..... ○ 1576 ..... 1450  
 Deluxe ..... ○ 1577 ..... 1450  
 Deluxe cuspidor unit mount activation (1576/1577) ..... 1100-480-CS ..... 200  
 Basic w/individual assistant's instrumentation holders in lieu of  
 telescoping arm ..... ○ 1678 ..... 810  
 Basic ..... ○ 1679 ..... 810  
 Basic cuspidor automatic cup fill ..... 1700-055-CS ..... 195  
 Vacuum drain cuspidor ..... 1400-115 ..... 420

LIGHT POST (specify light make/model)

Light post, 44" (Forest) ..... 1110-900 ..... 200  
 Light post, 44" (Pelton & Crane) ..... 1110-890 ..... 200

MONITOR MOUNT ASSEMBLY (monitor screen not included) (see page 32)

Ring mount (used with light post) ..... 1300RM ..... 950  
 Top post mount (used with monitor post) ..... 1300TPM ..... 950  
 Monitor post (for top post mount) ..... 1331-998 ..... 200

**TOTAL**

\_\_\_\_\_

Often Overlooked

**DESIGNER FRIENDLY™ SURFACE TREATMENTS** standard unless noted

Unit plastic finishes:  Forest Grey  White  Lemans  Laser Leather  Greywood

Powder coat paint for finishes:  Forest Grey  Shadow  White  Pearl  
 Pewter  Sage  Bronze  Arctic Silver  Charcoal

# delivery systems without chair


## fixed chair mount without sidebox, external umbilical



MODEL	4485FE Standard IC	4195FE Value
PRICE	\$ 5090	\$ 3890
STANDARD FEATURES	pneumatic, 3 HP, syringe	pinch, 3 HP, syringe

**DELIVERY SYSTEMS** can be mounted on any brand of chair

Choose your control head



### ADDITIONAL EQUIPMENT:

- 1075 Add \$ 840  
Assistant's telescoping arm w/instrumentation



- 1636FE Add \$ 1220  
Cuspidor, Basic



- 1637FE Add \$ 1670  
Cuspidor, Basic  
w/assistant's instrumentation



## STANDARD FEATURES

### DELIVERY SYSTEM

Choice of Control Head  
Automatic control for handpieces  
asepsis 3-way air/water syringe w/tubing, 5'  
(Euro tubing 3')  
handpiece tubing, 4-hole Midwest, asepsis style, 5'  
(Euro tubing 3')  
master control block (with pneumatic or pinch valves)  
master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air coolant adjustment  
handpiece flush  
handpiece exhaust filter/oil collector  
single handle brake release function  
self-contained clean water system, w/quick disconnect  
single bottle, 1½ liter  
wet/dry foot control  
Flexible control arm with pneumatic brake  
vertical travel 28" / horizontal travel 62"  
External umbilical, 5'  
Deluxe utility center  
air pressure regulator with filter  
master air shut-off valve  
 Chair adapter, fixed (specify chair make/model)

### ASSISTANT'S INSTRUMENTATION 1075

Central vacuum accessory package  
central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'  
external umbilical, 5'  
Post mount, 2"

### CUSPIDOR, BASIC 1636FE

Gravity drain cuspidor with swing away arm  
timed bowl rinse  
manual cup fill  
removable ceramic bowl  
Post mount adapter, 2"  
External 5' umbilical  
Water pressure regulator with filter  
Master water shut-off valve

### CUSPIDOR, BASIC 1637FE

*Same as 1636FE above, with additional features:*

Individual assistant's instrumentation holders  
Central vacuum accessory package  
central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'

## FACTORY OPTIONS

### PART NO

### RETAIL

### PRICE

DELIVERY SYSTEM	MODEL CHOICE	PRICE:
Handpiece illumination system		
one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595
add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265
one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595
add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190
Additional handpiece position, wet .....	1100-189-CS	405
Additional handpiece position, air drive only .....	1100-190-CS	195
Less handpiece position .....	1700-099-CS	-100
Heated syringe conversion kit (1) per syringe .....	1107-450-CS	320
Assistant's Instrumentation		
(control head mount) .....	1409-CS	400
Quick disconnect with flow control, 1/4" (specify air or water) .....	0016-455-CS	100
Additional brake handle activation (left) .....	0023-642-CS	140
Swivel tray, arm mounted .....	1110-468	335
Chrome accents .....	9000-998-CS	135
Standard utility center, housing and cover, in lieu of deluxe .....	1115-001-CS	-100
Additional umbilical (per foot) .....	1115-999-CS	30
Less self-contained clean water system, w/quick disconnect, single bottle, add water shutoff and regulator .....	1700-355-CS	N/C
Routing valve kit w/water shutoff & regulator/factory installation .....	1115-315-CS	235
Routing valve kit for use with cuspidor/factory installation .....	1115-316-CS	60
Water shut-off valve & regulator kit .....	1115-298	160
Without chair adapter .....	1115-099-CS	-150
 Ancillaries, uses existing handpiece position (see pages 62-63)		
Write in selection here .....		
Write in selection here .....		
ASSISTANT'S INSTRUMENTATION .....	1075	840
Additional autoclavable standard HVE w/tubing, 5' & autoclavable oral dry cup .....	1400-113-CS	85
CUSPIDOR, BASIC .....	<input type="radio"/> 1636FE <input type="radio"/> with assistant's instrumentation .....	1220 1670
Additional autoclavable standard HVE w/tubing, 5' & autoclavable oral dry cup (1635) .....	1400-113-CS	85
Cuspidor automatic cup fill .....	1700-055-CS	195
Vacuum drain cuspidor .....	1400-115	420
 LIGHT POST w/2" post clamp (specify light make/model)		
Light post, 35" (Forest) .....	1110-820	260
Light post, 44" (Pelton & Crane) .....	1110-822	260
 MONITOR MOUNT ASSEMBLY (monitor screen not included) (see page 32)		
Ring mount (used with light post) .....	1300RM	950
Top post mount (used with monitor post) .....	1300TPM	950
Monitor post (for top post mount) .....	1331-997	260

### TOTAL

 Often Overlooked

### DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Unit plastic finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> White	<input type="radio"/> Lemans	<input type="radio"/> Laser Leather	<input type="radio"/> Greywood
Powder coat paint for finishes:	<input type="radio"/> Forest Grey <input type="radio"/> Pewter	<input type="radio"/> Shadow <input type="radio"/> Sage	<input type="radio"/> White <input type="radio"/> Bronze	<input type="radio"/> Pearl <input type="radio"/> Arctic Silver	<input type="radio"/> Charcoal

# cuspidors / assistant's instrumentation

## fixed chair post mount, external umbilical



BASIC MODEL 1636  
\$ 1095



BASIC MODEL 1637  
\$ 1495



MODEL 1075  
\$ 840

### STANDARD FEATURES

#### CUSPIDOR, BASIC 1636

Gravity drain cuspidor with swing away arm  
timed bowl rinse  
manual cup fill  
removable ceramic bowl  
Post mount adapter, 2"  
External umbilical, 5'

(does not include post or chair adapter)

#### CUSPIDOR, BASIC 1637

Same as 1636 above, with additional features:

Individual assistant's instrumentation holders  
Central vacuum accessory package  
central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'

(does not include post or chair adapter)

#### POST MOUNT VACUUM TELESCOPING ARM 1075

Post mount adapter, 2"  
Central vacuum accessory package  
post mounted central vacuum canister  
w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'  
External umbilical, 8'

(does not include post or chair adapter)

Shipping weight:  
1636 17 lb / 8 kg  
1637 20 lb / 9 kg  
1075 15 lb / 7 kg

\* All optional utility centers include:  
housing and cover, air & water  
pressure regulator with filter master  
air and water shut-off valve, unless  
otherwise stated

### FACTORY OPTIONS

#### CUSPIDOR, BASIC .....

MODEL CHOICE ..... PRICE :  
Vacuum drain cuspidor ..... 1400-115 \$ 420

Basic cuspidor automatic cup fill ..... 1700-055-CS 195

Self-contained clean water system, w/quick disconnect

single bottle, 1½ liter ..... 1700-310-CS 265

Heated syringe conversion kit (1637) ..... 1107-450-CS 320

Additional autoclavable standard HVE

w/tubing, 5' & autoclavable oral dry cup (1637) ..... 1400-113-CS 85

Quick disconnect with flow control, 1/4" (specify air or water) ..... 0016-455-CS 100

Standard utility center housing & cover ..... 1115-020 130

Standard utility center \* ..... 1115-001 475

Additional umbilical (per foot) ..... 1115-999-CS 30

#### MOUNTING POSTS 2" diameter:

Control post with cap, 18" ..... 1110-487 110

#### LIGHT POSTS: 44", 2" diameter

Forest model ..... 1110-794 200

Pelton & Crane ..... 1110-785 200

#### ASSISTANT'S INSTRUMENTATION .....

Additional standard autoclavable HVE

w/tubing, 5' & autoclavable oral dry cup ..... 1400-113-CS 85

Heated syringe conversion kit ..... 1107-450-CS 320

Self-contained clean water system, w/quick disconnect

single bottle, 1½ liter ..... 1700-323-CS 265

#### Assistant's air shut-off & regulator utility center w/o housing & cover ..... 1115-008 225

Standard utility center housing & cover ..... 1115-020 130

Additional umbilical (per foot) ..... 1115-998-CS 20

#### TOTAL

#### Often Overlooked

#### DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Unit plastic finishes:  Forest Grey  White  Lemans  Laser Leather  Greywood

Powder coat paint for finishes:  Forest Grey  Shadow  White  Pearl  Charcoal  
 Pewter  Sage  Bronze  Arctic Silver



MODEL 1800RM  
\$ 1890



MODEL 1800TPM  
\$ 1890

#### STANDARD FEATURES

Flexible Tray Assembly  
Flexible counterbalanced arm system 26" rigid x 24" flex  
with pneumatic brake; accommodates up to 13½ lbs  
Stainless steel tray holder, 9¾" x 13½"  
wraparound holder bar  
dual handle brake release function  
on/off toggle, 3-way with exhaust  
Air pressure regulator with filter  
Master air shut-off valve  
RM: ring post mount (slides over 2" post)  
TPM: top post mount (inserts into 1½" ID bearing)  
Optional wall mount bracket available

#### FACTORY OPTIONS

FLEX TRAY .....	MODEL CHOICE _____	PRICE : _____
44" light post w/slot (w/Forest light bushing for RM) .....	1110-170	200
30" control post, with 1½" ID bronze bearing (for TPM) .....	1110-480	185
Wall mount bracket with 6" x 6" mounting plate (for TPM) .....	1110-386	185
Laminate wall mount panel 18 x 24" .....	1110-498	285
(wall mount bracket sold separately for TPM)		

**TOTAL** \_\_\_\_\_

Shipping weight:  
15 lb / 7 kg


#### DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Powder coat paint finishes:  Forest Grey  Shadow  White  Pearl  Charcoal  
 Pewter  Sage  Bronze  Arctic Silver

# rear + side delivery systems

## cabinetry mounting configurations

M&CC PFR8144




### good thinking

(low profile 1/2" hardened steel mount  
leaves room for your stool casters)

### made for each other


from the ground up, Forest delivery  
systems and M&CC cabinets are  
designed for easy integration




### start with a strong foundation

- your choice of 3/4" hardened steel 11"x 14 $\frac{1}{2}$ " (pre-installed on the PFR8144) or 10"x 24" mounting plate that bolts to the bottom of your counter top
- the 6"x 6" leveling plate is next, then a set of needle roller bearings with friction adjustment
- a super strong 2 $\frac{1}{2}$ " steel tube scissor arm assembly
- needle roller bearings are used where the two arm sections pivot

M&CC PFR8144



see page 45 for  
ordering information




\* shown with optional  
storage organizer,  
tray and additional  
assistant's holder bar

# rear + side delivery systems

## cabinetry mounting configurations


M&CC DR CABINET



3235

3735

see page 46-47 for  
ordering information



OPTIONAL

STANDARD

2400

2700


MT600BS/MT610BS

see page 48 for  
ordering information

see page 52 for  
ordering information

## why so many choices?

choose exactly what works best for YOU –  
Forest has the delivery system and mounting  
solutions for any cabinetry system




2180

1280

see page 49 for  
ordering information

# rear + side delivery systems

## 7000 Series: rear mount, side delivery



### STANDARD FEATURES

#### DUO SWING DELIVERY SYSTEM 7000 SERIES

##### Duo control

automatic control w/lockout feature for three handpieces  
asepsis 3-way air water syringe w/tubing, 4' handpiece tubing, 4-hole Midwest, asepsis style, 4' pneumatic master control block w/diaphragm master on/off individual handpiece pressure adjustment individual water coolant adjustment non-retraction water coolant air coolant adjustment handpiece flush handpiece exhaust filter/oil collector self-contained clean water system, w/quick disconnect single bottle, 1½ liter wet/dry foot control

##### Height adjustable work surface, 5"

##### "Burr Buddy" tray holder and covered Zirc tray

##### Central vacuum accessory package

central vacuum canister w/disposable solids collector asepsis 3-way air/water syringe w/tubing, 5' autoclavable standard HVE lever valve w/tubing, 5' autoclavable saliva ejector lever valve w/tubing, 5'

##### Assistant's quad holder assembly

##### Solid work surface, 19" x 24"

##### Universal heavy duty mounting plate

11" x 14 1/2", beneath counter (7000BC)

##### Laminate wall mount panel 18" x 24" (7000WM)

##### Air pressure regulator with filter

##### Master air shut-off valve

##### Umbilical, 7'

#### BENEATH COUNTER MOUNT

#### MODEL 7000BC

**\$ 5950**

Mounting plate dimensions, see page 77

Shipping weight:

7000BC, UC, FM, WM 130 lb / 59 kg

# rear + side delivery systems

**7000 Series: rear mount, side delivery**

CABINET MOUNTING CONFIGURATION  
(see page 40)



UNDER CABINET MOUNT  
MODEL 7000UC  
**\$ 5950**



FLOOR MOUNT  
MODEL 7000FM  
**\$ 5950**



WALL MOUNT  
MODEL 7000WM  
**\$ 5950**

## DESIGNER FRIENDLY™ SURFACE TREATMENTS

*standard unless noted*

Powder coat paint finishes:  Forest Grey  Shadow  White  
 Pearl  Charcoal  Pewter  Sage  
 Bronze  Arctic Silver

Solid work surface colors:

Corian:  Antartica  Blue Pebble  Doeskin

Fossil  Savannah

Gibraltar:  Artic Melange  Baja Melange  Frosty White Mirage

Light Beige Melange  Midnight Melange  Northern Melange

Meganite:  Goldrush Gemstone  Moonstone

Laminate work surface colors in lieu of solid surface 0023-776-LAM  
deduct —\$525

Formica:  Almond  Folkstone  Graphite

Oyster Gray  Surf  White Blanc

Wilsonart:  Shadow

Laminate wall mount panel (Model 7000WM)

Formica:  Almond  Folkstone  Graphite

Oyster Gray  Surf  White Blanc

Wilsonart:  Shadow

FACTORY OPTIONS	PART NO	RETAIL	PRICE
DELIVERY SYSTEM .....	MODEL CHOICE _____	PRICE: _____	
Handpiece Illumination System			
one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595	_____
add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265	_____
one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595	_____
add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190	_____
Additional handpiece position, wet .....	1100-189-CS	405	_____
Additional handpiece position, air drive only .....	1100-190-CS	195	_____
Less handpiece position .....	1700-099-CS	-100	_____
Additional standard autoclavable HVE			
w/tubing, 5' & autoclavable oral dry cup .....	1400-113-CS	85	_____
Heated syringe conversion kit (1) per syringe .....	1107-450-CS	320	_____
Additional assistant's quad holder bar .....	0024-156	150	_____
Quick disconnect w/flow control 1/4" (for water) .....	0016-454-CS	100	_____
Quick disconnect w/flow control 3/8" (for air) .....	0016-456-CS	100	_____
Storage organizer w/bins & frame .....	1110-179-CS	420	_____
Swivel tray .....	1110-909-CS	335	_____
Touchpad for 3900 chair (25 ft. Cable) .....	1700-081-1-CS	500	_____
Chrome accents .....	9000-998-CS	135	_____
Monitor mount assembly (monitor screen not included) .....	1300-025-CS	395	_____
Standard utility center housing & cover .....	1115-020	130	_____
In floor recessed utility center housing & cover .....	1115-305	270	_____
Additional umbilical (per foot) .....	1115-999-CS	30	_____
Less self-contained clean water system, w/quick disconnect,			
single bottle, add water shutoff and regulator .....	1700-355-CS	N/C	
Routing valve kit w/water shutoff & regulator/factory installation .....	1115-315-CS	235	_____
Water shut-off valve & regulator kit .....	1115-298	160	_____
Standard color laminate work surface in lieu of solid surface .....	0023-776-LAM	-525	_____
Special color solid or laminate work surface .....	0023-776-SSV / 0023-776-X	(call for price)	
Universal mounting plate, 10"x 24" in lieu of 11"x 14 1/2" .....	0010-084-CS	125	_____
Ancillaries, uses existing handpiece position (see pages 62-63)			
Write in selection here .....			_____
Write in selection here .....			_____
Often Overlooked			
		TOTAL	_____

# rear + side delivery systems

## rear wall mount, side delivery



**MODEL 3486WM**  
\$ 2750



**MODEL 3386WM**  
Work Surface without controls  
\$ 2250



**MODEL 3686WM**  
\$ 4750

### STANDARD FEATURES

#### ASSISTANT'S SWING DELIVERY SYSTEM 3486WM

Central vacuum accessory package  
arm mounted central vacuum canister  
w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'  
self-contained clean water system w/quick disconnect,  
single bottle 1½ liter

Assistant's quad holder assembly  
Wall mount rigid arm assembly

Laminate wall mount panel, 18" x 24"

Laminate work surface, 19" x 24"

Air pressure regulator with filter

Master shut-off valve

Umbilical, 7'

#### DUO SWING DELIVERY SYSTEM 3686WM

*Same as Unit above, with additional features:*

Standard IC control head  
automatic control for three handpieces  
asepsis 3-way air/water syringe w/tubing, 4'  
handpiece tubing, 4-hole Midwest, asepsis style, 4'  
pneumatic master control block w/diaphragm  
master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air coolant adjustment  
handpiece flush  
handpiece exhaust filter/oil collector  
wet/dry foot control

#### WORK SURFACE SYSTEM 3386WM

Wall mount rigid arm assembly  
Laminate wall mount panel, 18" x 24"  
Laminate work surface, 19" x 24"

### FACTORY OPTIONS

### PART NO RETAIL PRICE

DELIVERY SYSTEM .....	MODEL CHOICE .....	PRICE:
Handpiece Illumination System (3686)		
one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595
add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265
one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595
add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190
Additional handpiece position, wet (3686) .....	1100-189-CS	405
Additional handpiece position, air drive only (3686) .....	1100-190-CS	195
Less handpiece position (3686) .....	1700-099-CS	-100
Additional standard autoclavable HVE		
w/tubing, 5' & autoclavable oral dry cup .....	1400-113-CS	85
Heated syringe conversion kit (1) per syringe .....	1107-450-CS	320
Additional assistant's quad holder bar .....	0024-149	150
Quick disconnect w/flow control 1/4" (for water) .....	0016-454-CS	100
Quick disconnect w/flow control 3/8" (for air) .....	0016-456-CS	100
Storage organizer w/bins & frame .....	1110-179-CS	420
Swivel tray .....	1110-909-CS	335
Touchpad for 3900 chair (25 ft. Cable) .....	1700-081-I-CS	500
Chrome accents .....	9000-998-CS	135
Monitor mount assembly (monitor screen not included) .....	1300-025-CS	395
Less standard utility center housing & cover .....	1115-020-CS	-50
Additional umbilical (per foot) (3686) .....	1115-999-CS	30
Additional umbilical (per foot) (3486) .....	1115-998-CS	20
Less self-contained clean water system, w/quick disconnect, single bottle, add water shutoff and regulator .....	1700-355-CS	N/C
Routing valve kit w/water shutoff & regulator/factory installation .....	1115-315-CS	235
Water shut-off valve & regulator kit .....	1115-298	160
Standard color solid work surface upgrade .....	0023-003-SS	525
Special color laminate or solid work surface upgrade .....	0023-003-X / 0023-003-SSV	(call for price)
Ancillaries, uses existing handpiece position (see pages 62-63)		
Write in selection here .....		
Write in selection here .....		
Often Overlooked		
	<b>TOTAL</b>	



#### DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Unit plastic finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> White	<input type="radio"/> Lemans	<input type="radio"/> Greywood	<input type="radio"/> Laser Leather
Powder coat paint finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> Shadow	<input type="radio"/> White	<input type="radio"/> Pearl	<input type="radio"/> Charcoal
	<input type="radio"/> Pewter	<input type="radio"/> Sage	<input type="radio"/> Bronze	<input type="radio"/> Arctic Silver	<input type="radio"/> Oyster Gray
Laminate work surface and wall mount colors:	<input type="radio"/> Almond	<input type="radio"/> Graphite	<input type="radio"/> Surf	<input type="radio"/> Folkstone	<input type="radio"/> White Blanc
	<input type="radio"/> Oyster Gray	<input type="radio"/> White Blanc			<input type="radio"/> Shadow
Wilsonart:	<input type="radio"/> Shadow				

Mounting plate dimensions, see page 79

Shipping weight:

3686 120 lb / 55 kg      3486 115 lb / 53 kg

3386 100 lb / 46 kg

# rear + side delivery systems

rear beneath counter mount, side delivery



**MODEL 3486BC**  
\$ 2750



CABINET MOUNTING CONFIGURATION (see page 40)



**MODEL 3686BC**  
\$ 4750

## STANDARD FEATURES

### ASSISTANT'S SWING DELIVERY SYSTEM 3486BC

Central vacuum accessory package  
arm mounted central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'  
Self contained clean water system, single bottle 1½ liter  
Assistant's quad holder assembly  
Beneath counter rigid arm assembly  
Beneath counter heavy duty mounting plate, 11" x 14-1/2"  
Laminate work surface, 19" x 24"  
Air pressure regulator with filter  
Master shut-off valve  
Umbilical, 7'

### DUO SWING DELIVERY SYSTEM 3686BC

*Same as Unit above, with additional features:*

Standard IC control head  
automatic control for three handpieces  
asepsis 3-way air/water syringe w/tubing, 4'  
handpiece tubing, 4-hole Midwest, asepsis style, 4'  
pneumatic master control block w/diaphragm  
master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air coolant adjustment  
handpiece flush  
handpiece exhaust filter/oil collector  
wet/dry foot control

### WORK SURFACE SYSTEM 3386BC

Beneath counter rigid arm assembly  
Beneath counter heavy duty mounting plate, 11" x 14-1/2"  
Laminate work surface, 19" x 24"

## FACTORY OPTIONS

DELIVERY SYSTEM .....	MODEL CHOICE .....	PRICE :
Handpiece Illumination System (3686)		
one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595
add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265
one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595
add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190
Additional handpiece position, wet (3686) .....	1100-189-CS	405
Additional handpiece position, air drive only (3686) .....	1100-190-CS	195
Less handpiece position (3686) .....	1700-099-CS	-100
Additional standard autoclavable HVE		
w/tubing, 5' & autoclavable oral dry cup .....	1400-113-CS	85
Heated syringe conversion kit (1) per syringe .....	1107-450-CS	320
Additional assistant's quad holder bar .....	0024-149	150
Quick disconnect w/flow control 1/4" (for water) .....	0016-454-CS	100
Quick disconnect w/flow control 3/8" (for air) .....	0016-456-CS	100
► Storage organizer w/bins & frame .....	1110-179-CS	420
► Swivel tray .....	1110-909-CS	335
Touchpad for 3900 chair (25 ft. Cable) .....	1700-081-1-CS	500
Chrome accents .....	9000-998-CS	135
Monitor mount assembly (monitor screen not included) .....	1300-025-CS	395
► Universal heavy duty mounting plate, 10" x 24" in lieu of 11" x 14 1/2" .....	0010-084-CS	125
Standard utility center housing & cover .....	1115-020	130
Additional umbilical (per foot) (3686) .....	1115-999-CS	30
Additional umbilical (per foot) (3486) .....	1115-998-CS	20
Less self-contained clean water system, w/quick disconnect,		
single bottle, add water shutoff and regulator .....	1700-355-CS	N/C
Routing valve kit w/water shutoff & regulator/factory installation .....	1115-315-CS	235
Water shut-off valve & regulator kit .....	1115-298	160
Standard color solid work surface upgrade .....	0023-003-SS	525
Special color laminate or solid work surface upgrade .....	0023-003-X / 0023-003-SSV	(call for price)
► Ancillaries, uses existing handpiece position (see pages 62-63)		
Write in selection here .....		
Write in selection here .....		
► Often Overlooked		

**TOTAL**

## ► DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Unit plastic finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> White	<input type="radio"/> Lemans	<input type="radio"/> Greywood	<input type="radio"/> Laser Leather
Powder coat paint finishes:	<input type="radio"/> Charcoal	<input type="radio"/> Forest Grey	<input type="radio"/> Shadow	<input type="radio"/> White	<input type="radio"/> Pearl
		<input type="radio"/> Pewter	<input type="radio"/> Sage	<input type="radio"/> Bronze	<input type="radio"/> Arctic Silver
Laminate work surface colors:	Formica:	<input type="radio"/> Almond	<input type="radio"/> Graphite	<input type="radio"/> Surf	<input type="radio"/> Folkstone
		<input type="radio"/> Oyster Gray	<input type="radio"/> White Blanc		
	Wilsonart:	<input type="radio"/> Shadow			

Mounting plate dimensions, see page 79

Shipping weight:

3486 120 lb / 55 kg      3486 115 lb / 53 kg  
3386 100 lb / 46 kg



# rear + side delivery systems

## side cabinet & wall mount - value control



**12" rigid arm  
cabinet mount**



**22" rigid arm  
wall mount**

CABINET MOUNTING CONFIGURATIONS  
(see page 41)

MODEL 3735  
VALUE CONTROL HEAD  
**\$ 3420**

MODEL 3730  
VALUE CONTROL HEAD  
**\$ 3560**

### STANDARD FEATURES

#### DOCTOR VALUE DELIVERY SYSTEM 3730/3735

Value control head  
automatic control for three handpieces  
asepsis 3-way air/water syringe w/tubing, 5' handpiece tubing, 4-hole Midwest, asepsis style, 5' individual pinch valve control blocks  
master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air coolant adjustment  
handpiece flush  
handpiece exhaust filter/oil collector  
single handle brake release function  
self-contained clean water system w/quick disconnect, single bottle 1½ liter  
wet/dry foot control  
Flexible control arm with pneumatic brake (3735)  
cabinet vertical mount 12" rigid arm & 18" flex arm  
Flexible control arm with pneumatic brake (3730)  
wall vertical mount 22" rigid arm & 18" flex arm  
External umbilical, 5'  
Standard utility center housing and cover (3730)  
Air pressure regulator with filter  
Master air shut-off valve

### FACTORY OPTIONS

PART NO	RETAIL	PRICE
DELIVERY SYSTEM .....	MODEL CHOICE _____	PRICE: _____
Handpiece Illumination System		
one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595
add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265
one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595
add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190
Additional handpiece position, wet .....	1100-189-CS	405
Additional handpiece position, air drive only .....	1100-190-CS	195
Less handpiece position .....	1700-099-CS	-100
Heated syringe conversion kit .....	1107-450-CS	320
Assistant's Instrumentation (control head mount) .....	1408-CS	320
Quick disconnect with flow control, 1/4" (specify air or water) .....	0016-455-CS	100
Additional brake release handle activation (left) .....	0023-642-CS	140
Swivel tray, arm mounted (over control head) .....	1110-468	335
Touchpad for 3900 chair (25 ft. Cable) .....	1700-081-1-CS	500
Chrome accents .....	9000-998-CS	135
Standard utility center housing & cover (3735) .....	1115-020	130
Additional umbilical (per foot) .....	1115-999-CS	30
Less self-contained clean water system, w/quick disconnect, single bottle, add water shutoff and regulator .....	1700-355-CS	N/C
Routing valve kit w/water shutoff & regulator/factory installation .....	1115-315-CS	235
Water shut-off valve & regulator kit .....	1115-298	160
► Laminate wall mount panel, 18" x 24" .....	1110-498	285
Ring mount rigid arm in lieu of top post mount (see page 69)		
22" w/o vertical mount .....	1110-536-RM	N/C
26" w/o vertical mount .....	1110-537-RM	N/C
26" rigid w/24" flex w/o vertical mount .....	1110-771-RM	N/C
► Ancillaries, uses existing handpiece position (see pages 62-63)		
Write in selection here .....		_____
Write in selection here .....		_____
		<b>TOTAL</b> _____

► Often Overlooked

#### ► DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Unit plastic finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> White	<input type="radio"/> Lemans	<input type="radio"/> Laser Leather	<input type="radio"/> Greywood
Powder coat paint finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> Shadow	<input type="radio"/> White	<input type="radio"/> Pearl	<input type="radio"/> Charcoal
	<input type="radio"/> Pewter	<input type="radio"/> Sage	<input type="radio"/> Bronze	<input type="radio"/> Arctic Silver	

Mounting plate dimensions, see page 80

Shipping weight:  
3730, 3735    50 lb / 23 kg

# rear + side delivery systems

## rear or side beneath counter mount



**MODEL 2400**  
**\$ 2900**



(beneath counter slide  
optional see page 68)



**MODEL 2700**  
**\$ 2100**

### CABINET MOUNTING CONFIGURATION

(see page 41)

#### STANDARD FEATURES

##### DOCTOR IC DELIVERY SYSTEM 2400

Beneath counter swivel mount (work surface not included)  
Standard IC control head

automatic control for three handpieces  
asepsis 3-way air/water syringe w/tubing, 5'  
handpiece tubing, 4-hole Midwest, asepsis style, 5'  
pneumatic master control block w/diaphragm  
master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air coolant adjustment  
handpiece flush  
handpiece exhaust filter/oil collector  
self-contained clean water system w/quick disconnect,  
single bottle 1½ liter  
wet/dry foot control  
External braided umbilical, 7'

Air pressure regulator with filter  
Master air shut-off valve

##### DOCTOR VALUE DELIVERY SYSTEM 2700

Same as Unit above, with alternative control:

Value control head  
individual pinch valve control blocks

#### FACTORY OPTIONS

#### PART NO

#### RETAIL

#### PRICE

DELIVERY SYSTEM .....	MODEL CHOICE _____	PRICE : _____
Handpiece Illumination System		
one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595
add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265
one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595
add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190
Additional handpiece position, wet .....	1100-189-CS	405
Additional handpiece position, air drive only .....	1100-190-CS	195
Less handpiece position .....	1700-099-CS	-100
Heated syringe conversion kit .....	1107-450-CS	320
Assistant's Instrumentation (control head mount) .....	1408-CS	320
Quick disconnect with flow control, 1/4" (specify air or water) .....	0016-455-CS	100
Slide mount with swivel, 9" extension .....	1700-140	280
Folding arm assembly, 8" x 12" with wall mount .....	1110-115-CS	335
Touchpad for 3900 chair (25 ft. Cable) .....	1700-081-1-CS	500
Chrome accents .....	9000-998-CS	135
Standard utility center housing & cover .....	1115-020	130
Additional umbilical (per foot) .....	1115-998-CS	20
Less self-contained clean water system, w/quick disconnect, single bottle, add water shutoff and regulator .....	1700-355-CS	N/C
Routing valve kit w/water shutoff & regulator/factory installation .....	1115-315-CS	235
Water shut-off valve & regulator kit .....	1115-298	160
► Ancillaries, uses existing handpiece position (see pages 62-63)		
Write in selection here .....		
Write in selection here .....		
► Often Overlooked		
		<b>TOTAL</b> _____

DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted					
Unit plastic finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> White	<input type="radio"/> Lemans	<input type="radio"/> Laser Leather	<input type="radio"/> Greywood
Powder coat paint finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> Shadow	<input type="radio"/> White	<input type="radio"/> Pearl	<input type="radio"/> Charcoal
	<input type="radio"/> Pewter	<input type="radio"/> Sage	<input type="radio"/> Bronze	<input type="radio"/> Arctic Silver	

Mounting plate dimensions, see page 80

Shipping weight:

2400, 2700 28 lb / 13 kg


# rear + side delivery systems

## rear beneath counter mount

mounts to counter or floating work surface



**MODEL 1280**  
\$ 1590



CABINET MOUNTING CONFIGURATION  
(see page 41)



**MODEL 2180**  
\$ 3610

### STANDARD FEATURES

#### HORIZONTAL ASSISTANT'S DELIVERY SYSTEM 1280

Central vacuum accessory package, swing mount arm mounted central vacuum canister w/disposable solids collector asepsis 3-way air/water syringe w/tubing, 5' autoclavable universal HVE lever valve w/tubing, 5' autoclavable saliva ejector lever valve w/tubing, 5' self-contained clean water system w/quick disconnect, single bottle 1½ liter Horizontal surface mount External umbilical, 7' Air pressure regulator with filter Master air shut-off valve

#### HORIZONTAL DUO DELIVERY SYSTEM 2180

*Same as Unit above, with additional features:*

Automatic panel control for three handpieces Instrument delivery asepsis 3-way air/water syringe w/tubing, 5' handpiece tubing, 4-hole Midwest, asepsis style, 5' pneumatic master control block w/diaphragm master on/off individual handpiece pressure adjustment individual water coolant adjustment non-retraction water coolant air coolant adjustment wet/dry foot control Internal umbilical, 7'

Mounting plate dimensions, see page 81  
Panel cut-out, see page 83

Shipping weight:

1280 12 lb / 5 kg  
2180 30 lb / 14 kg

### FACTORY OPTIONS

### PART NO RETAIL PRICE

DELIVERY SYSTEM .....	MODEL CHOICE .....	PRICE : .....
Handpiece Illumination System (2180)		
one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595
add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265
one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595
add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190
Additional handpiece position, wet (2180) .....	1100-189-CS	405
Additional handpiece position, air drive only (2180) .....	1100-190-CS	195
Less handpiece position (2180) .....	1700-099-CS	-100
Additional universal autoclavable HVE position		
w/tubing, 5' & autoclavable oral dry cup .....	1400-119-CS	85
Dual articulating quad holder assembly, universal HVE .....	0024-224-CS	235
Heated syringe conversion kit (1) per syringe .....	1107-450-CS	320
Slide mount, 7 5/8" extension .....	1700-141	280
Touchpad for 3900 chair (25 ft. Cable) .....	1700-081-I-CS	500
Chrome accents (2180) .....	9000-998-CS	135
Floating work surface to under counter arm assembly, 20" .....	0023-806	730
Universal heavy-duty mounting plate, 11" x 14 1/2" .....	0010-086-BLK	195
Universal heavy-duty mounting plate 10" x 24" .....	0010-084-BLK	320
Standard color laminate work surface, 19" x 24" .....	0010-652-LAM	325
Standard color solid work surface, 19" x 24" .....	0010-652-SS	850
Special color laminate or solid work surface 19" x 24" .....	0010-652-X / 0010-652-SSV	(call for price)
Standard utility center housing & cover .....	1115-020	130
Additional umbilical (per foot) .....	1115-998-CS	20
Less self-contained clean water system, w/quick disconnect,		
single bottle, add water shutoff and regulator .....	1700-355-CS	N/C
Routing valve kit w/water shutoff & regulator/factory installation .....	1115-315-CS	235
Water shut-off valve & regulator kit .....	1115-298	160

► Ancillaries, uses existing handpiece position (see pages 62-63)

Write in selection here .....

Write in selection here .....

### TOTAL

► Often Overlooked

### DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Unit plastic finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> White	<input type="radio"/> Lemans	<input type="radio"/> Laser Leather	<input type="radio"/> Greywood
Powder coat paint finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> Shadow	<input type="radio"/> White	<input type="radio"/> Pearl	<input type="radio"/> Charcoal
	<input type="radio"/> Pewter	<input type="radio"/> Sage	<input type="radio"/> Bronze	<input type="radio"/> Arctic Silver	

# optional delivery systems

## panel mounts



**MODEL 2500**  
\$ 2295



**MODEL 2510**  
\$ 1650

### STANDARD FEATURES

### FACTORY OPTIONS

### PART NO RETAIL

### PRICE

#### DOCTOR CONTROL 2500

Automatic control for three handpieces  
 asepsis 3-way air/water syringe w/tubing, 7'  
 handpiece tubing, 4-hole Midwest, asepsis style, 7'  
 pneumatic master control block w/diaphragm  
 master on/off  
 individual handpiece pressure adjustment  
 individual water coolant adjustment  
 non-retraction water coolant  
 air coolant adjustment  
 self-contained clean water system w/quick disconnect,  
     single bottle 1½ liter  
     wet/dry foot control  
 Remote mount holder assembly with control tubing, 6'  
 External umbilical, 6'  
 Air pressure regulator with filter  
 Master air shut-off valve

#### DOCTOR CONTROL 2510

Manual control for two handpieces  
 asepsis 3-way air/water syringe w/tubing, 7'  
 handpiece tubing, 4-hole Midwest, asepsis style, 7'  
 individual module control block with manual selector,  
     1 wet and 1 dry handpiece position  
 individual handpiece pressure adjustment  
 non-retraction water coolant  
 self contained clean water system,  
     single bottle 1½ liter  
     wet/dry foot control  
 Remote mount holder assembly with control tubing, 6'  
 External umbilical, 6'  
 Air pressure regulator with filter  
 Master air shut-off valve

#### DELIVERY SYSTEM .....MODEL CHOICE

#### PRICE:

#### Handpiece Illumination System

one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595
add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265
one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595
add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190
Additional handpiece position, wet (2500) .....	1100-189-CS	405
Additional handpiece position, air drive only (2500) .....	1100-190-CS	195
Heated syringe conversion kit (1) per syringe .....	1107-450-CS	320
Standard utility center housing & cover .....	1115-020	130
Additional umbilical (per foot) .....	0017-819-CS	10
Less self-contained clean water system, w/quick disconnect, single bottle, add water shutoff and regulator .....	1700-355-CS	N/C
Water shut-off valve & regulator kit .....	1115-298	160
► Quick disconnect for air or water .....	dealer installed (see page 69)	
► Routing valve kit for city water .....	dealer installed (see page 69)	

#### TOTAL

#### ► Often Overlooked

Wall/cabinet cut-out dimensions, see page 83



#### DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Powder coat paint finishes:  Forest Grey  Shadow  White  Pearl  Charcoal  
 Pewter  Sage  Bronze  Arctic Silver

# optional delivery systems

panel mounts



**MODEL 2515**  
\$ 2995



**MODEL 2505**  
\$ 3945

## STANDARD FEATURES

### DUO CONTROL with ASSISTANT'S INSTRUMENTATION 2515

Automatic control for two handpieces  
asepsis 3-way air/water syringe w/tubing, 7'  
handpiece tubing, 4-hole Midwest, asepsis style, 7'  
handpiece illumination system (5-hole)  
pneumatic master control block w/diaphragm  
master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air coolant adjustment  
self contained clean water system,  
single bottle 1½ liter  
wet/dry foot control w/chip blower and 12' tubing  
Single assistant's instrumentation  
central vacuum canister w/disposable solids collector  
autoclavable universal HVE lever valve with tubing, 7'  
autoclavable saliva ejector lever valve with tubing, 7'  
Remote mount 5 position holder assembly  
External umbilical, 6'  
Air pressure regulator with filter  
Master air shut-off valve

### DUO CONTROL with DUAL ASSISTANT'S INSTRUMENTATION 2505

*Same as Unit above, with additional features:*

Automatic control for four handpieces  
additional handpiece illumination tubing (5-hole)  
Dual, remote mount, 2 position holder assemblies  
Dual assistant's instrumentation for left/right mounting  
central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe  
autoclavable universal HVE lever valve with tubing, 7'  
autoclavable saliva ejector lever valve with tubing, 7'

## FACTORY OPTIONS

DELIVERY SYSTEM .....	MODEL CHOICE _____	PRICE: _____
Handpiece Illumination System		
add'l handpiece illumination tubing (5-hole) (2515) .....	0017-615-CS	265
one handpiece (6-pin), in lieu of 5-hole .....	0017-603-1-CS	N/C
add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190
Heated syringe conversion kit (1) per syringe .....	1107-450-CS	320
Standard utility center housing & cover .....	1115-020	130
Additional umbilical (per foot) .....	0017-819-CS	10
Less self-contained clean water system, w/quick disconnect, single bottle, add water shutoff and regulator .....	1700-355-CS	N/C
Water shut-off valve & regulator kit .....	1115-298	160
► Quick disconnect for air or water .....	dealer installed (see page 69)	
► Routing valve kit for city water .....	dealer installed (see page 69)	

**TOTAL** \_\_\_\_\_

► Often Overlooked

Wall/cabinet cut-out dimensions, see page 83


Shipping weight:  
2515 17 lb / 8 kg      2505 20 lb / 10 kg

## DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Powder coat paint finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> Shadow	<input type="radio"/> White	<input type="radio"/> Pearl	<input type="radio"/> Charcoal
	<input type="radio"/> Pewter	<input type="radio"/> Sage	<input type="radio"/> Bronze	<input type="radio"/> Arctic Silver	

# optional delivery systems

## compact controls



STANDARD FEATURES	FACTORY OPTIONS	PART NO	RETAIL	PRICE
<b>DOCTOR CONTROL WITH 3 HP's</b>	DELIVERY SYSTEM ..... MODEL CHOICE _____ Handpiece Illumination System one handpiece (5-hole), accommodates 1 add'l tubing ..... 1700-021-CS add'l handpiece illumination tubing (5-hole) ..... 0017-615-CS one handpiece (6-pin), accommodates 1 add'l tubing ..... 1700-019-CS add'l handpiece illumination tubing (6-pin) ..... 0017-603-CS Heated syringe conversion kit (1) per syringe ..... 1107-450-CS Slide mount with swivel, 7 5/8" extension (order BS) ..... 1700-140 Folding arm assembly, 8" x 12" with wall mount (order TS) ..... 1110-115-CS Standard utility center housing & cover ..... 1115-020 Additional umbilical (per foot) ..... 0017-819-CS Less self-contained clean water system, w/quick disconnect, single bottle, add water shutoff and regulator ..... 1700-355-CS Water shut-off valve & regulator kit ..... 1115-298 Assistant Instrumentation Option (Drain Tubing 12') These options add extended holder bar to MT's with holders mounted Universal autoclavable with (2 holders) (HVE, saliva ejector & canister) ..... 1406-CS Standard autoclavable with (2 holders) (HVE, saliva ejector & canister) ..... 1405-CS Additional drain tubing (per foot) ..... 0017-933-CS Quick disconnect for air or water ..... dealer installed (see page 69) Routing valve kit for city water ..... dealer installed (see page 69)			
Top Shelf Mount Beneath Shelf Mount with Swivel	MT600TS MT600BS			
<b>DOCTOR CONTROL WITH 2 HP's</b>				
Top Shelf Mount Beneath Shelf Mount with Swivel	MT610TS MT610BS			
<i>All units include the following features:</i>				
Automatic control for three handpieces, or Automatic control for two handpieces asepsis 3-way air/water syringe w/tubing, 4' handpiece tubing, 4-hole Midwest, asepsis style, 4' pneumatic master control block w/diaphragm master on/off individual handpiece pressure adjustment non-retraction water coolant air coolant adjustment self-contained clean water system w/quick disconnect, single bottle 1½ liter wet/dry type foot control External umbilical, 12' Air pressure regulator with filter Master air shut-off valve				
Mounting plate dimensions, see page 83	DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted		<b>TOTAL</b>	
Shipping weight: MT600, MT610 20 lb / 10 kg	Plastic finishes: <input type="radio"/> Forest Grey <input type="radio"/> White Powder coat paint finishes: <input type="radio"/> Forest Grey <input type="radio"/> White			

# optional delivery systems

## beneath counter assistant's instrumentation



**MODEL 1055**


**\$ 1160**

(work surface not included)

STANDARD FEATURES	FACTORY OPTIONS	PART NO	RETAIL	PRICE
<b>BENEATH SHELF SWING MOUNT 1055</b>	DELIVERY SYSTEM ..... 1055 Heated syringe conversion kit ..... 1107-450-CS Additional standard autoclavable HVE w/tubing, 5' & autoclavable oral dry cup ..... 1400-113-CS • Quick disconnect with flow control, 1/4" (specify air or water) .... 0016-455-CS Dual articulating quad holder assembly, standard HVE ..... 0024-227-CS • Wet handpiece position, includes self-contained clean water & foot control ..... 1100-492-CS • Dry handpiece position includes foot control ..... 1100-491-CS Wall mount conversion ..... 1110-349-CS ⇒ Assistant's air shut-off & regulator utility center w/o housing & cover .... 1115-008 Standard utility center without housing & cover ..... 1115-004 Standard utility center housing & cover ..... 1115-020 Additional umbilical (per foot) ..... 1115-998-CS Self-contained clean water system w/quick disconnect, single bottle 1½ liter ..... 1700-301	1055 1107-450-CS 1400-113-CS 0016-455-CS 0024-227-CS 1100-492-CS 1100-491-CS 1110-349-CS 1115-008 1115-004 1115-020 1115-998-CS 1700-301	\$1160 320 85 100 235 640 310 115 225 365 130 20 265	_____ _____ _____ _____ _____ _____ _____ _____ _____ _____ _____ _____ _____ _____
			<b>TOTAL</b>	_____
	<ul style="list-style-type: none"> <li>• Cannot be purchased in combination. Choice of one ONLY.</li> </ul>			
	⇒ Often Overlooked			
Mounting plate dimensions, see page 83 Shipping weight: 12 lb / 6 kg	<b>DESIGNER FRIENDLY™ SURFACE TREATMENTS</b> standard unless noted Plastic finishes: <input type="radio"/> Forest Grey <input type="radio"/> White Powder coat paint finishes: <input type="radio"/> Forest Grey <input type="radio"/> White			

Mounting plate dimensions, see page 83  
Shipping weight:  
12 lb / 6 kg





7000MC Cart



5886 Cart

## Carts

### Affordable Excellence

#### Mobile Carts

- Material: north america's only all metal control head (7000) metal control lid (5200/5210/5400/5410/5886)
- Benefit: resistant to UV & cleaners finish in one of nine designer friendly™ paint color choices
- Material: stainless steel lid hinges (7000/5200/5210/5400/5410/5886)
- Benefit: quick access to controls & ancillary products
- Material: exclusive brass handpiece lockout valves (7000)
- Benefits: allows assistant to remove a handpiece from its holder & change a burr while dentist continues working
- Material: secondary aluminum work surface on 7000 cart (burr buddy)
- Benefit: dedicated work surface for burr block, electric handpiece attachment, mirror, explorer
- Material: powder coated aluminum tub holders, cassette trays & monitor mounts (optional on 5386/5886/7000)
- Benefit: expandable work surface accommodates additional armamentarium
- Material: thick steel-walled cart frames with high-quality all surface polyurethane casters on all mobile carts
- Benefit: extra stable mobile delivery moves to accommodate your position



forest dental

[forestdental.com](http://forestdental.com)

[sales@forestdental.com](mailto:sales@forestdental.com)



## Powder Coat

arctic silver

bronze

charcoal

forest grey

pearl

pewter

sage

shadow

white

laser leather

lemans

greywood

### Designer Friendly™

- Material: nine powder coat paint colors including white & forest grey
- Benefit: dental equipment with a **WOW** factor
- Material: white, forest grey & three carbon fiber-look plastic lids (5910/5915)
- Benefit: rich design improves overall image of practice & patient experience
- Material: 13 solid surface material colors (standard 7000 & optional 5386/5886)
- Benefit: tough, repairable surface enhancing operatory interior design



5200/5400 Cart (Arctic Silver Powder Coat, Chrome Option)



5910/5915 Cart (Greywood Carbon Fiber, Chrome Option)

## 7000 Series: mobile carts



Optional storage organizer and  
swivel tray  
1110-179-CS and 1110-909-CS



MODEL 7000MC  
\$ 5950

## STANDARD FEATURES

**DUO SWING CART**

Automatic control w/lockout feature for three handpieces  
asepsis 3-way air water syringe w/tubing, 4' handpiece tubing, 4-hole Midwest, asepsis style, 4' pneumatic master control block w/diaphragm master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air coolant adjustment  
handpiece flush  
handpiece exhaust filter/oil collector  
wet/dry foot control  
Assistant's quad holder assembly  
"Burr Buddy" sliding tray holder and covered Zirc tray  
Pneumatic work surface easy lift/close mechanism  
Central vacuum accessory package  
canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5'  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'  
Self contained clean water system, single bottle  
1½ liter w/quick disconnect bottle assembly  
Solid work surface, 19" x 24"  
Adjustable height frame, 24" x 46"  
Twin wheel hard floor casters  
External umbilical, 6'  
Standard utility center housing and cover  
Air pressure regulator with filter  
Master air shut-off valve

**DESIGNER FRIENDLY™ SURFACE TREATMENTS**  
*standard unless noted*

## ► Solid work surface colors:

- |            |   |  |   |
|------------|---|--|---|
| Corian:    | <input type="radio"/> Antarctica          | <input type="radio"/> Blue Pebble      | <input type="radio"/> Doeskin             |
|            | <input type="radio"/> Fossil              | <input type="radio"/> Savannah         |   |
| Gibraltar: | <input type="radio"/> Arctic Melange      | <input type="radio"/> Baja Melange     | <input type="radio"/> Frosty White Mirage |
|            | <input type="radio"/> Light Beige Melange | <input type="radio"/> Midnight Melange | <input type="radio"/> Northern Melange    |
| Meganite:  | <input type="radio"/> Goldrush Gemstone   | <input type="radio"/> Moonstone        |   |

Shipping weight: 119 lb / 54 kg

## FACTORY OPTIONS

PART NO	RETAIL	PRICE
CART .....	7000MC	\$ 5950
Handpiece Illumination System		
one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595
add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265
one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595
add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190
Additional handpiece position, wet .....	1100-189-CS	405
Additional handpiece position, air drive only .....	1100-190-CS	195
Less handpiece position .....	1700-099-CS	-100
Additional standard autoclavable HVE w/tubing, 5' & autoclavable oral dry cup .....	1400-113-CS	85
Heated syringe conversion kit (1) per syringe .....	1107-450-CS	320
Additional assistant's quad holder bar .....	0024-156	150
Quick disconnect w/flow control 1/4" (for water) .....	0016-454-CS	100
Quick disconnect w/flow control 3/8" (for air) .....	0016-456-CS	100
► Storage organizer w/bins & frame .....	1110-179-CS	420
► Swivel tray .....	1110-909-CS	335
Twin wheel carpet casters .....	0010-401-CS	N/C
Touchpad for 3900 chair (25 ft. Cable) .....	1700-081-1-CS	500
Chrome accents .....	9000-998-CS	135
Monitor assembly (monitor screen not included) .....	1300-025-CS	395
Standard color laminate work surface in lieu of solid surface .....	0023-776-LAM	-525
Special color solid or laminate work surface .....	0023-776-SSV / 0023-776-X	(call for price)
Additional umbilical (per foot) .....	1115-999-CS	30
Less self-contained clean water system, w/quick disconnect, single bottle, add water shutoff and regulator .....	1700-355-CS	N/C
Routing valve kit w/water shutoff & regulator/factory installation .....	1115-315-CS	235
Water shut-off valve & regulator kit .....	1115-298	160
► Ancillaries, uses existing handpiece position (see pages 62-63)		
Write in selection here .....		
Write in selection here .....		

**TOTAL**

► Often Overlooked

**DESIGNER FRIENDLY™ SURFACE TREATMENTS**  
*standard unless noted*

- |                             |                                   |                              |                              |                                     |                                |
|-----------------------------|-----------------------------------|------------------------------|------------------------------|-------------------------------------|--------------------------------|
| Powder coat paint finishes: | <input type="radio"/> Forest Grey | <input type="radio"/> Shadow | <input type="radio"/> White  | <input type="radio"/> Pearl         | <input type="radio"/> Charcoal |
|                             | <input type="radio"/> Pewter      | <input type="radio"/> Sage   | <input type="radio"/> Bronze | <input type="radio"/> Arctic Silver |                                |

# carts

## mobile carts

**MODEL 5386  
\$ 3380**



**MODEL 5886  
\$ 4950**



### STANDARD FEATURES

#### ASSISTANT'S CART W/INSTRUMENTATION 5386

Central vacuum accessory package, swing mount arm mounted central vacuum canister w/disposable solids collector  
asepsis 3-way air/water syringe w/tubing, 5' autoclavable standard HVE lever valve w/tubing, 5' autoclavable saliva ejector lever valve w/tubing, 5' self-contained clean water system w/quick disconnect, single bottle 1½ liter

Assistant's holder quad assembly

Laminate work surface, 19" x 24"

Adjustable height frame, 30" to 42"

Twin wheel hard floor casters

Electrical outlet/UL

External umbilical, 7'

Assistant's utility center with housing and cover  
air pressure regulator with filter  
master shut-off valve

#### DUO SWING CART 5886

Same as Unit above, with additional features:

Standard IC Control Head

automatic control for three handpieces  
asepsis 3-way air/water syringe w/tubing, 4'  
handpiece tubing, 4-hole Midwest, asepsis style, 4'  
pneumatic master control block w/diaphragm  
master on/off

individual handpiece pressure adjustment  
individual water coolant adjustment

non-retraction water coolant  
air coolant adjustment

handpiece flush

handpiece exhaust filter/oil collector  
wet/dry foot control

Standard utility center housing and cover

Air pressure regulator with filter

Master air shut-off valve

Shipping weight:

5386 100 lb / 46 kg

5886 105 lb / 48 kg

### FACTORY OPTIONS

### PART NO

### RETAIL PRICE

DELIVERY SYSTEM .....	MODEL CHOICE _____	PRICE: _____
Handpiece Illumination System (5886)		
one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595
add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265
one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595
add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190
Additional handpiece position, wet (5886) .....	1100-189-CS	405
Additional handpiece position, air drive only (5886) .....	1100-190-CS	195
Less handpiece position (5886) .....	1700-099-CS	-100
Additional standard autoclavable HVE		
w/tubing, 5' & autoclavable oral dry cup .....	1400-113-CS	85
Heated syringe conversion kit (l) per syringe .....	1107-450-CS	320
Additional assistant's quad holder bar .....	0024-149	150
Quick disconnect w/flow control 1/4" (for water) .....	0016-454-CS	100
Quick disconnect w/flow control 3/8" (for air) .....	0016-456-CS	100
Storage organizer w/bins & frame .....	1110-179-CS	420
Swivel tray .....	1110-909-CS	335
Touchpad for 3900 chair (25 ft. Cable) .....	1700-081-1-CS	500
Chrome accents .....	9000-998-CS	135
Twin wheel carpet casters .....	0010-401-CS	N/C
Monitor mount assembly (monitor screen not included) .....	1300-025-CS	395
Less standard utility center housing & cover .....	1115-020-CS	-50
Additional umbilical (per foot) (5386) .....	1115-998-CS	20
Additional umbilical (per foot) (5886) .....	1115-999-CS	30
Less self-contained clean water system, w/quick disconnect, single bottle, add water shutoff and regulator .....	1700-355-CS	N/C
Routing valve kit w/water shutoff & regulator/factory installation .....	1115-315-CS	235
Water shutoff valve & regulator kit .....	1115-298	160
Standard color solid work surface upgrade .....	0023-003-SS	525
Special color laminate or solid work surface upgrade . 0023-003-X / 0023-003-SSV	(call for price)	
Ancillaries, uses existing handpiece position (see pages 62-63)		
Write in selection here .....		
Write in selection here .....		
Often Overlooked		
		<b>TOTAL</b> _____

### DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Unit plastic finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> White	<input type="radio"/> Lemans	<input type="radio"/> Greywood	<input type="radio"/> Laser Leather
Powder coat paint finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> Shadow	<input type="radio"/> White	<input type="radio"/> Pearl	<input type="radio"/> Charcoal
	<input type="radio"/> Pewter	<input type="radio"/> Sage	<input type="radio"/> Bronze	<input type="radio"/> Arctic Silver	
Laminate work surface colors:	<input type="radio"/> Almond	<input type="radio"/> Graphite	<input type="radio"/> Surf	<input type="radio"/> Folkstone	<input type="radio"/> Oyster Gray
	<input type="radio"/> White Blanc	<input type="radio"/> Shadow			
Wilsonart:	<input type="radio"/> Shadow				



MODEL 5910  
\$ 2220



MODEL 5915  
\$ 2595

#### STANDARD FEATURES

##### DOCTOR CART 5910

Value Control Head  
automatic control for three handpieces  
asepsis 3-way air/water syringe w/tubing, 4'  
handpiece tubing, 4-hole Midwest, asepsis style, 4'  
individual pinch valve control blocks  
master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air coolant adjustment  
handpiece flush  
handpiece exhaust filter/oil collector  
self-contained clean water system w/quick disconnect,  
single bottle 1½ liter  
wet/dry foot control  
Adjustable height frame, 26" to 35"  
Twin wheel hard floor casters  
External umbilical, 12' with 1/4" male QD for air supply  
Air pressure regulator with filter inside the head

##### DUO CART 5915

*Same as Unit above, with additional features:*

Central vacuum canister w/disposable solids collector  
autoclavable standard HVE lever valve w/tubing, 5'  
autoclavable saliva ejector lever valve w/tubing, 5'

#### FACTORY OPTIONS

CART .....	MODEL CHOICE .....	PRICE: .....
Handpiece Illumination System		
one handpiece (5-hole), accommodates 1 add'l tubing .....	1700-021-CS	595
add'l handpiece illumination tubing (5-hole) .....	0017-615-CS	265
one handpiece (6-pin), accommodates 1 add'l tubing .....	1700-019-CS	595
add'l handpiece illumination tubing (6-pin) .....	0017-603-CS	190
Additional handpiece position, wet .....	1100-189-CS	405
Additional handpiece position, air drive only .....	1100-190-CS	195
Less handpiece position .....	1700-099-CS	- 100
Additional standard autoclavable HVE		
w/tubing, 5' & autoclavable oral dry cup (5915) .....	1400-113-CS	85
Heated syringe conversion kit (1) per syringe .....	1107-450-CS	320
Quick disconnect with flow control, 1/4" (specify air or water) .....	0016-455-CS	100
Foot control support shelf, clamps to 2" base .....	0023-665	180
Touchpad for 3900 chair (25 ft. Cable) .....	1700-081-1-CS	500
Chrome accents .....	9000-998-CS	135
Twin wheel carpet casters .....	0010-401-CS	N/C
Utility shelf assembly, 12" x 12" .....	1110-655	280
Standard utility center housing & cover only .....	1115-020	130
Routing valve kit w/water shutoff & regulator/factory installation .....	1115-315-CS	235
Water shut-off valve & regulator kit .....	1115-298	160
► Ancillaries, uses existing handpiece position (see pages 62-63)		
Write in selection here .....		
Write in selection here .....		
► Often Overlooked		
		<b>TOTAL</b> .....

#### ► DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Shipping weight:  
5910    35 lb / 16 kg  
5915    40 lb / 18 kg

Unit plastic finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> White	<input type="radio"/> Lemans	<input type="radio"/> Laser Leather	<input type="radio"/> Greywood
Powder coat paint finishes:	<input type="radio"/> Forest Grey	<input type="radio"/> Shadow	<input type="radio"/> White	<input type="radio"/> Pearl	<input type="radio"/> Charcoal
	<input type="radio"/> Pewter	<input type="radio"/> Sage	<input type="radio"/> Bronze	<input type="radio"/> Arctic Silver	





MODEL 5210  
\$ 3590



OPTIONAL MUFFLER  
PART NO. 3400-005



MODEL 5410  
\$ 4575

#### STANDARD FEATURES

##### **DOCTOR CART WITHOUT COMPRESSOR 5210**

Automatic control for three handpieces  
asepsis 3-way air/water syringe w/tubing, 4' handpiece tubing, 4-hole Midwest, asepsis style, 4' pneumatic master control block w/diaphragm master on/off  
individual handpiece pressure adjustment  
individual water coolant adjustment  
non-retraction water coolant  
air coolant adjustment  
handpiece flush  
handpiece exhaust filter/oil collector  
self-contained clean water system w/quick disconnect, single bottle 1½ liter  
wet/dry foot control  
Adjustable height frame, 26" to 35"  
Twin wheel hard floor casters  
External air line, 12' with 3/8 male QD for air supply (for remote compressor or nitrogen supply, existing compressor, etc.)

##### **DUO CART 5410**

*Same as Unit above, with additional features:*

Self-contained air-operated vacuum package with collection jar  
automatic autoclavable standard HVE lever valve w/ tubing, 5'  
automatic autoclavable saliva ejector lever valve w/ tubing, 5'

\* Nitrogen fitting not included

Shipping weight:  
5210 50 lb / 23 kg  
5410 55 lb / 25 kg

#### FACTORY OPTIONS

CART	MODEL CHOICE	PRICE:
Handpiece Illumination System		
one handpiece (5-hole), accommodates 1 add'l tubing	1700-021-CS	595
add'l handpiece illumination tubing (5-hole)	0017-615-CS	265
one handpiece (6-pin), accommodates 1 add'l tubing	1700-019-CS	595
add'l handpiece illumination tubing (6-pin)	0017-603-CS	190
Additional handpiece position, wet	1100-189-CS	405
Additional handpiece position, air drive only	1100-190-CS	195
Less handpiece position	1700-099-CS	-100
Heated syringe conversion kit	1107-450-CS	320
Quick disconnect w/flow control 1/4" (for water)	0016-454-CS	100
Quick disconnect w/flow control 3/8" (for air)	0016-456-CS	100
Foot control support shelf, clamps to 2" base	0023-665	180
Twin wheel carpet floor casters	0010-401-CS	N/C
Utility shelf assembly, 12" x 12"	1110-655	280
Standard utility center housing & cover	1115-020-CS	-50
Additional air supply umbilical (per foot)	0017-860	3
► Ancillaries, uses existing handpiece position (see pages 62-63)		
Write in selection here		
Write in selection here		
Muffler for air venturi (5410)	3400-005	180

**TOTAL**

► Often Overlooked

#### DESIGNER FRIENDLY™ SURFACE TREATMENTS standard unless noted

Plastic finishes:  Forest Grey  White

Powder coat paint finishes:  Forest Grey  White

# factory options + parts

## chair adapters, monitor posts + light posts

STANDARD CHAIR ADAPTERS - for 2" post	PART NO.	RETAIL
3900 FIXED accommodates external umbilical	1115-195	\$ 430
3900 FIXED accommodates internal umbilical	1115-530	565
3900 FIXED FOR SIDEBOX accommodates internal umbilical	1115-531	565
3900 PIVOT accommodates internal umbilical	1115-528	565

MONITOR POSTS - 2" diameter with oilite bushing (used with 1300TPM and 1800TPM)	PART NO.	RETAIL
30" W/DUAL SLOT, FIXED SIDEBOX MOUNT	1331-998	\$ 200
30" W/SINGLE SLOT, FIXED (INCLUDES POST CLAMP)	1331-997	260

LIGHT POSTS - 2" with light bushing	PART NO.	RETAIL
24" EXTENSION (INCLUDES POST CLAMP) FOREST PELTON & CRANE	1110-518	\$ 260
35" LIGHT POST (INCLUDES POST CLAMP) FOREST PELTON & CRANE	1110-820	260
35" LIGHT POST W/SLOT (FOR MONITOR CORD) (INCLUDES POST CLAMP) FOREST	1110-822	260
35" LIGHT POST W/SLOT (FOR MONITOR CORD) FOREST	1110-045	260
44" LIGHT POST FOREST PELTON & CRANE	1110-794	200
44" LIGHT POST W/SLOT (FOR MONITOR CORD) FOREST	1110-785	200
44" LIGHT POST W/SLOT (SIDEBOX OPTION) FOREST	1110-050	200
44" LIGHT POST W/SLOT (FOR MONITOR CORD & SIDEBOX ACCESS) FOREST	1110-900	200
	1110-890	200
	1110-170	200

# factory options + parts

## chair options



PEDO BOOSTER SEAT to accommodate smaller patients  
specify upholstery color ..... 3907-180 \$ 280

- A cushion that takes the lumbar out of the chair, essentially turning it into a bed, ideal for children.



IV BOARD ..... 3900-200 \$ 640  
 • Adjustable angle to accommodate any size patient  
 • Double articulating scissor arm for industry leading range of placement

### 3900 REMOTE MOUNT TOUCHPAD WITH 25' CABLE



#### FIXED HORIZONTAL MOUNT

Forest Grey .....	1700-079-G	\$ 625
Designer Friendly .....	1700-079-C	625



#### SWING HORIZONTAL MOUNT

Forest Grey .....	1700-089-HM-G	625
Designer Friendly .....	1700-089-HM-C	625



#### SWING VERTICAL MOUNT

Forest Grey .....	1700-089-VM-G	625
Designer Friendly .....	1700-089-VM-C	625

# factory options + parts

## ancillary options

PART NO RETAIL



### BIEN-AIR OPTIMA INT ELECTRIC HANDPIECE CONTROL & HANDPIECE CONTROL W/ENDO, uses existing handpiece position

Factory option ..... BIEN-AIR-MCX / BIEN-AIR-MX2 \$ 2195 / 4095

Please call 800-433-2436 for warranty and technical support of the Handpiece

#### BIEN-AIR OPTIMA INT ELECTRIC HANDPIECE CONTROL (BIEN-AIR-MCX)

- Attachments purchased separately from Bien-air
- Highest torque output on the market
- Swivel tubing reduces arm fatigue
- Quick-disconnect, autoclavable motor with type E connection compliant with ISO 3964
- High and low speed applications using one attachment
- MX motor speed from 100 to 40,000 rpm
- 24 volt power supply
- Ancillary shelf
- Bien-Air warranty 3 years from date of purchase

#### BIEN-AIR OPTIMA INT ELECTRIC HANDPIECE CONTROL W/ENDO (BIEN-AIR-MX2)

Same as above with additional features

- Auto-reverse and torque control for endo procedures
- 40 memory position with pre-programmed procedures



### SATELEC SP NEWTRON PIEZO ULTRASONIC SCALER, w/power supply, uses existing handpiece position

Factory option ..... SATELEC-S \$ 1845

### SATELEC SP NEWTRON LED LIT PIEZO ULTRASONIC SCALER, w/power supply, uses existing handpiece position

Factory option ..... SATELEC-S-LED 2320

Please call 800-289-6367 for warranty and technical support of the Scaler

- Color coding system ensures proper tip and power setting
- Three tips: #1, #0Z, H3 (Intro Scaling Kit)
- Tips available for periodontics, endodontics, scaling, restorative, and cavity preparation
- Cruise Control™ automatically adjusts frequency amplitude
- On/off controlled by rheostat
- First ever autoclavable LED lit handpiece – 100,000 LUX
- Newtron technology means power, control and reliability
- Satelec warranty 1 year / tips 30 days from date of purchase



### DENTSPLY CAVITRON MAGNETOSTRICTIVE ULTRASONIC SCALER

Factory option ..... CAVITRON \$ 3250

Please call 800-800-7202 for warranty and technical support of the Scaler

- Blue Zone™ for patient comfort
- SPS® (Sustained Performance Technology) Constant balance between scaling efficiency and patient comfort
- Handpiece with 330 degree swivel and lavage control (reduces cord drag)
- Autoclavable Steri-mate™ (infection control)
- PowerLINE 1000 30K insert
- Frequency 30 kHz
- Cavitron warranty 2 years from date of purchase



### ANCILLARY SHELF

Factory option ..... 0010-310 \$ 95

- Accommodates additional ancillary products
- Designed for easy access to installed devices
- Powder coated aluminum
- Installed on all control heads except Euro and MT compact controls

# factory options + parts

## ancillary options

	PART NO	RETAIL
<b>AIR TECHNIQUES POLARIS &amp; SPECTRA CAMERAS</b>		
Factory option .....	AIRTECH-P / AIRTECH-S	\$ 4500 / 4800
Please call 800-247-8324 for warranty and technical support of Camera.		
<ul style="list-style-type: none"> <li>• New Direct USB connection fits in standard chair side delivery unit</li> <li>• TWAIN interface for easy "Play &amp; Play" software integrations</li> <li>• Hot swappable USB umbilical between Polaris &amp; Spectra</li> <li>• New convenient On/Off rubber button</li> <li>• Revolutionary 120 degree rubber handgrip button ring allows for capture of all angles with ease and comfort</li> <li>• Top button capture on BUTTON RELEASE</li> <li>• Sapphire crystal lens (no scratching)</li> <li>• USB cable, 25' with signal booster</li> <li>• Air Techniques warranty 3 years from date of purchase</li> </ul>		
		
<b>POLARIS (AIRTECH-P)</b>		
<ul style="list-style-type: none"> <li>• Fixed focal distance: 6 to 40 mm</li> <li>• No manual focusing required</li> </ul>		
<b>SPECTRA (AIRTECH-S)</b>		
<ul style="list-style-type: none"> <li>• Uses fluorescence to detect caries within fissures &amp; on smooth surfaces that often go undetected by traditional means</li> <li>• Detects decay around the margins of existing composite and amalgam restorations</li> <li>• Can be used during the restorative phase to verify that all caries have been removed</li> <li>• Promotes minimally invasive treatment techniques</li> <li>• Displays a tangible color/numerical value to show the extent of decay present, providing means of tracking progression and regression of treatment plans over time</li> </ul>		
		
<b>USB CABLE, 16'</b>		
Factory installed in delivery system .....	0014-286	\$ 85
		
<b>USB BOOSTER CONNECTOR CABLE, 16'</b>		
Factory option .....	0014-355	100
<b>250 WATT (UNIVERSAL POWER SUPPLY)</b>		
Factory option .....	0014-376	\$ 435
<ul style="list-style-type: none"> <li>• Required when using two or more ancillary devices, or when using a Forest Operatory light and any Forest ancillary product (excluding the 3900 soft-touch massager power source)</li> <li>• Provides power for MULTIPLE electrical accessories such as: fiber-optics, heated syringe tubing, operatory lights, scalers, curing lights, and electric micro motors</li> <li>• Size: 7.5" - 5.5" - 2.5". Easily fits in all Forest floor boxes</li> <li>• Quick disconnect connectors for all of Forest ancillary products</li> <li>• Built in air-electric switch, enables power source to turn on/off with pilot air</li> <li>• Available in both 230 volt and 110 volt</li> <li>• Auto-resettable Klixon switch: inbuilt overcurrent and over temperature protection, eliminating the need for manually changing fuses or circuit breakers</li> </ul>		
		
<b>ASSISTANT'S INSTRUMENTATION (CONTROL HEAD MOUNT)</b>		
Factory option .....	1408-CS	\$ 320
Factory option .....	1409-CS	400
<ul style="list-style-type: none"> <li>• Canister with 1/2" drain</li> <li>• Mounting bracket</li> <li>• Drain tubing</li> <li>• Autoclavable standard HVE w/tubing, 5'</li> <li>• Autoclavable saliva ejector w/tubing, 5'</li> <li>• Standard vacuum holder</li> <li>• Universal holder</li> <li>• Installed on all control heads except Euro and MT compact controls</li> </ul>		

# factory options + parts



	PART NO	RETAIL
<b>HANDPIECE ILLUMINATION SYSTEM, 5 HOLE ADAPTER WITH KRYPTON LAMP</b>		
Factory option .....	1700-019-CS	\$ 595
with Sterling tubing .....	1700-052-STG	595
with Shadow tubing .....	1700-052-SHW	595
with Black tubing .....	1700-052-BLK	595
<b>HANDPIECE ILLUMINATION SYSTEM, 6 PIN</b>		
Factory option .....	1700-021-CS	595
with Sterling tubing .....	1700-051-STG	595
with Shadow tubing .....	1700-051-SHW	595
with Black tubing .....	1700-051-BLK	595
<ul style="list-style-type: none"> <li>• Dual channel outputs individually adjustable from 2VDC to 4VDC</li> <li>• Accommodates one additional handpiece tubing</li> <li>• Illumination delay (adjustable between 0-25 seconds) can be set individually for two handpieces</li> <li>• Automatic activation of handpiece illumination with the dental unit foot control</li> <li>• Simple installation into any new or existing dental unit</li> <li>• Indicator LEDs confirms correct operation of the control module</li> <li>• Powers all brands of optic handpieces and air-motors</li> <li>• ISO standard handpiece terminal with stainless steel nut</li> <li>• Asepsis style handpiece tubing, 7' long (4' to 5' installed)</li> <li>• Power transformer, output 12VAC, 1.5 Amps</li> <li>• Installation hardware and instructions</li> </ul>		




	PART NO	RETAIL
<b>HANDPIECE ILLUMINATION TUBING, 5 HOLE ADAPTER WITH KRYPTON LAMP</b>		
Factory option .....	0017-615-CS	\$ 265
Sterling .....	0017-615-STG	265
Shadow .....	0017-615-SHW	265
Black .....	0017-615-BLK	265
Used to replace power optic handpiece tubing		
• ISO standard power optic handpiece terminal with stainless steel nut		
• 5 hole adapter with krypton lamp		
• Asepsis style handpiece tubing, 7' (4' to 5' factory installed)		



	PART NO	RETAIL
<b>HANDPIECE ILLUMINATION TUBING, 6 PIN</b>		
Factory option .....	0017-603-CS	\$ 190
Factory option (in lieu of 5 hole) .....	0017-603-1-CS	N/C
Sterling .....	0017-603-STG	190
Shadow .....	0017-603-SHW	190
Black .....	0017-603-BLK	190
Used to replace power optic handpiece tubing		
• ISO standard power optic handpiece terminal with stainless steel adapter and nut (5-hole Midwest)		
• Asepsis style handpiece tubing, 7' (4' to 5' factory installed)		



KRYPTON LAMP .....	0014-067	\$ 80
For use in handpiece illumination system, 3.5VAC, 720mA		
Used with 0017-692, 0017-592 and 0017-675 tubing		
<b>ADAPTER WITH KRYPTON LAMP .....</b>		
For use in handpiece illumination system, 3.5VAC, 720mA		
Used with 1700-019 or 1700-052 system and 0017-615 tubing		

	PART NO	RETAIL
	HANDPIECE TUBING, ASEPSIS STYLE (VINYL) Sterling ..... Shadow ..... Black ..... 4-Hole, straight midwest terminal	0017-810 \$ 35 0017-510 35 0017-398 35
	DELUXE AIR/WATER SYRINGE Standard with quick-release tip ..... with sheathed tubing assembly, 5' Sterling ..... Shadow ..... Black .....  Straight with quick-release tip ..... with sheathed tubing assembly, 4' Sterling ..... Shadow ..... Black .....	1109-010 \$ 85 1107-755 135 1107-755-SH 135 1107-755-BLK 135  1109-025 105 1107-585-STG 155 1107-585-SHW 155 1107-585-BLK 155
	SYRINGE TUBING ASSEMBLY (POLYURETHANE), SHEATHED Sterling ..... Shadow ..... Black .....	0017-865 \$ 30 0017-565 30 0017-368 30
	HEATED SYRINGE CONVERSION KIT FACTORY OPTION ..... The heated syringe conversion kit adapts to existing 3-way syringes to provide instant warm water during normal operations. Includes: <ul style="list-style-type: none"><li>• Power source with variable temperature and on/off controls (1/2 Amp fuse, Volts 10 1/2 - 12 w/ load applied)</li><li>• Size, 2 3/4" x 5" x 2 1/2"</li><li>• Tubing</li><li>• Low voltage wire, 10'</li></ul>	1107-450-CS \$ 320
	HEATED SYRINGE TUBING ASSEMBLY, ASEPSIS STYLE Sterling ..... Shadow ..... Black .....	0017-122 \$ 130 0017-113 130 0017-522 130
	DELUXE STANDARD HVE VACUUM VALVE ..... Autoclavable with metal lever	1100-342 \$ 45
	DELUXE UNIVERSAL VACUUM VALVE ..... Autoclavable with metal lever	1100-343 \$ 45
	DELUXE SALIVA EJECTOR VALVE Autoclavable with metal lever with grey tip ..... with black tip .....	1100-341 \$ 40 1100-341-BLK 40

# factory options + parts

PART NO RETAIL



ADDITIONAL HVE with tubing 5' and autoclavable oral dry cup  
Standard vacuum valve

Factory option .....	1400-113-CS	\$ 85
Sterling .....	1400-113-G	85
Shadow .....	1400-113-SH	85
Black .....	1400-113-BLK	85
Universal vacuum valve		
Factory option .....	1400-119-CS	85
Sterling .....	1400-119-G	85
Shadow .....	1400-119-SH	85
Black .....	1400-119-BLK	85



ORAL DRY CUP ..... 1400-010 \$ 15  
Autoclavable polysulfone (beige). Fits all vacuum valves with 7/16" dia. tips.



SMOOTH VACUUM TUBING, 1/2" ID

Factory option .....	0017-933-CS	\$ 4.25 / foot
Sterling .....	0017-933	4.25 / foot
Shadow .....	0017-533	4.25 / foot
Black .....	0017-333	4.25 / foot



SMOOTH VACUUM TUBING, 1/2"HEX ID

Sterling .....	0017-932	3.25 / foot
Shadow .....	0017-532	3.25 / foot
Black .....	0017-332	3.25 / foot



SALIVA EJECTOR TUBING (VINYL)

Sterling .....	0017-930	\$.80 / foot
Shadow .....	0017-330	.80 / foot
Black .....	0017-530	.80 / foot




QUAD HOLDER ASSEMBLY, UNIVERSAL

White .....	0024-114	\$ 75
Black .....	0024-114-BLK	75



QUAD HOLDER ASSEMBLY, DOUBLE HVE

White .....	0024-127	75
Black .....	0024-127-BLK	75



DUAL ARTICULATING QUAD HOLDER ASSEMBLY

Factory option .....	0024-227-CS	\$ 235
4 standard, 4 universal, white .....	0024-227	310
4 standard, 4 universal, black .....	0024-227-BLK	310
Factory option .....	0024-224-CS	235
8 universal, white .....	0024-224	310
8 universal, black .....	0024-224-BLK	310

PART NO RETAIL



#### WET/DRY FOOT CONTROL

Forest Grey .....	1105-174-G	\$ 170
Designer Friendly (specify color) .....	1105-174-C	170



#### CONTOURED UTILITY CENTER HOUSING AND COVER (8 X 17 3/16 X 17 3/16)

Factory option .....	3900-101	\$ 420
Forest Grey .....	3900-101-G	420
Shadow .....	3900-101-SH	420
Designer Friendly .....	3900-101-C	520

#### COMPLETE CONTOURED UTILITY CENTER

• air and water pressure regulator with filter		
• master air and water shut-off valves		
Factory option .....	3900-102	\$ 765
Forest Grey .....	3900-102-G	765
Shadow .....	3900-102-SH	765
Designer Friendly .....	3900-102-C	865



#### DELUXE UTILITY CENTER HOUSING AND COVER (7 X 13 5/16 X 12 11/16)

• rear umbilical access		
Factory option .....	1115-007-CS	-\$ 200
Forest Grey .....	1115-007-G	225
Shadow .....	1115-007-SH	225
Black .....	1115-007-BLK	225
White .....	1115-007-WH	225
• right umbilical access		
Forest Grey .....	1115-009-G	225
Shadow .....	1115-009-SH	225
Black .....	1115-009-BLK	225
White .....	1115-009-WH	225

#### COMPLETE DELUXE UTILITY CENTER

• right umbilical access		
• air and water pressure regulator with filter		
• master air and water shut-off valves		
Forest Grey .....	1115-005-G	640
Shadow .....	1115-005-SH	640
Black .....	1115-005-BLK	640
White .....	1115-005-WH	640



#### STANDARD UTILITY CENTER HOUSING AND COVER (4 1/2 X 9 X 14)

Factory option .....	1115-020	\$ 130
Factory option .....	1115-020-CS	-\$ 50
Forest Grey .....	1115-020-G	130
Designer Friendly (specify color) .....	1115-020-C	130

#### COMPLETE STANDARD UTILITY CENTER

• air and water pressure regulator with filter		
• master air and water shut-off valves		
Factory option .....	1115-001	\$ 475
Factory option .....	1115-001-CS	-\$ 100
Forest Grey .....	1115-001-G	475
Designer Friendly (specify color) .....	1115-001-C	475

# factory options + parts

	PART NO	RETAIL
		
TRAY HOLDER ASSEMBLY, use with 1095 rear pivot assistant's instrumentation		
Factory option .....	1110-680	\$ 560
Forest Grey .....	1110-680-G	560
Designer Friendly (specify color) .....	1110-680-C	560
includes:		
• pivoting powder coated tray holder, 9 3/4" x 13 1/2"		
• clamp assembly adapts to 1 1/2" horizontal or vertical post		
		
SLIDE WITH SWIVEL 7 5/8" EXTENSION, use with under surface mount units 2400, 2700 & MT's		
Factory option .....	1700-140	\$ 280
Forest Grey .....	1700-140-G	280
Designer Friendly (specify color) .....	1700-140-C	280
SLIDE WITH SWIVEL 7 5/8" EXTENSION, use with under surface mount units 1280, 2180		
Factory option .....	1700-141	280
includes:		
• 4" x 4" horizontal mounting plate with smooth slide glides		
• 7 5/8" slide plate		
• swivel control bracket		
		
FOLDING ARM ASSEMBLY 8 X 12, use with under surface mount units 2400, 2745, MT's, 2180, 1280, 1700-141		
Factory option .....	1110-115-CS	\$ 335
Forest Grey .....	1110-115-G	335
Designer Friendly (specify color) .....	1110-115-C	335
Includes:		
• 8" folding arm with 1/2" pin and umbilical retainer		
• 12" folding arm with 1/2" pin and umbilical retainer		
• 4" x 4" wall mounting plate which accepts 1/2" pin		
• Control mount with 1/2" pin		
		
UNIVERSAL HEAVY DUTY MOUNTING PLATE		
10" X 24" beneath counter .....	0010-084-BLK	\$ 320
11" X 14 1/2" beneath counter .....	0010-086-BLK	195
10" X 24" in lieu of 11" X 14 1/2" .....	0010-084-CS	125
		
FLOATING WORK SURFACE ARM ASSEMBLY, 20"		
Factory option 1280/2180 .....	0023-806	\$ 730
with beneath counter and work surface support bracket		
Forest Grey .....	0023-815-G	730
Designer Friendly (specify color) .....	0023-815-C	730
		
OPTIONAL WALL MOUNT CONVERSION, allows 1055 to be mounted on vertical surface		
in lieu of under cabinet mount factory option .....	1110-349-CS	\$ 115
Includes:		
• 8" rigid arm with 1/2" pin and control mount adapter		
• 4" x 4" wall mounting plate which accepts 1/2" pin		
		
SUPPORT SHELF		
Factory option .....	0023-665	\$ 180
Forest Grey .....	0023-665-G	180
Designer Friendly (specify color) .....	0023-665-C	180
Provides adequate support for use of foot control while allowing cart mobility.		

	PART NO	RETAIL
<b>26" RIGID W/24" FLEXIBLE ARM</b>		
in lieu of top post mount factory option .....	1110-771-RM	N/C
22" rigid w/18" flexible arm in lieu of top post mount, factory option .....	1110-536-RM	N/C
26" rigid w/18" flexible arm in lieu of top post mount, factory option .....	1110-537-RM	N/C
Designed for use with 2" diameter post.		
<b>VACUUM DRAIN CUSPIDOR KIT</b> .....	1400-115	\$ 420
May be used to adapt all models of gravity drain cuspidor to central vacuum system. Includes:		
• Vacuum drain valve		
• Flush nipple		
• Filter assembly		
• Tubing		
• Installation hardware		
<b>AIR AND WATER ACCESSORY PANEL</b> Forest Grey .....	1300-030-G	\$ 165
Designer Friendly (specify color) .....	1300-030-C	165
• 3/8" female air quick disconnect		
• Toggled shut off for air quick disconnect		
• 1/4" female water quick disconnect		
• Flow adjustment for water quick disconnect		
• 4" x 3" flush mount panel		
<td></td> <td></td>		
<b>QUICK DISCONNECTS</b>		
1/4" QD with flow control, factory installed .....	0016-454-CS	\$ 100
1/4" QD with flow control, black / grey .....	0016-454-BLK / -G	100
1/4" QD with flow control, arm and control head mount, grey, factory installed .....	0016-455-CS	100
1/4" QD with flow control, arm and control head mount, black / grey .....	0016-455-BLK / -G	100
3/8" QD with flow control, grey, factory installed .....	0016-456-CS	100
3/8" QD with flow control, black / grey .....	0016-456-BLK / -G	100
<td></td> <td></td>		
<b>ROUTING VALVE KITS</b> — Check with local authorities about installation requirements		
Factory installed with water shut-off and regulator .....	1115-315-CS	\$ 235
Factory installed without water shut-off and regulator .....	1115-316-CS	60
Dealer installed includes; toggle valve, bracket, tubing,		
water shut-off and regulator, black / grey .....	1115-310-BLK / 1115-310	225
Dealer installed includes; toggle valve, bracket, and tubing		
black / grey .....	1115-312-BLK / 1115-312	225
<td></td> <td></td>		
<b>HIGH DENSITY BOTTLE</b>		
1½ liter .....	0010-464	\$ 30
1 liter .....	0010-364	20
For use with clean water system		
• Opaque bottle		
• High Density Polyethylene (HDPE)		
<td></td> <td></td>		
<b>AUTOMATIC HANDPIECE CONTROL BLOCK DIAPHRAGM</b>		
Tri-block .....	0022-020	\$ 3.25
Quad-block .....	0022-024	7.35

N/C lifetime warranty. Serial number and purchase order required. Freight and labor paid by dealer.

# factory options + parts

PART NO RETAIL



FOREST ALLEN WRENCH KIT .....	0008-019	\$ 18
Includes one each of the following sizes:		
.050, 1/16, 5/64, 3/32, 7/64, 1/8, 9/64, 5/32, 3/16, 7/32, 1/4, 5/16		



HOLDER for Automatic handpiece, 1/2" holder bar .....		
Forest Grey/Sterling .....	0010-836-G	\$ 70
Shadow .....	0010-836-SH	70
Black .....	0010-836-BLK	70
White .....	0010-836-WH	70



HOLDER, NOTCHED for automatic handpiece, fixed holder bar .....		
Forest Grey/Sterling .....	0010-835-G	\$ 70
Shadow .....	0010-835-SH	70
Black .....	0010-835-BLK	70
White .....	0010-835-WH	70



HOLDER, ELECTRIC automatic handpiece/universal vacuum/saliva ejector, normally open .....		
Forest Grey/Sterling .....	0010-856-G	\$ 75
Shadow .....	0010-856-SH	75
Black .....	0010-856-BLK	75
White .....	0010-856-WH	75



HOLDER, NOTCHED ELECTRIC automatic handpiece/universal vacuum/saliva ejector, normally open, fixed holder bar .....		
Forest Grey/Sterling .....	0010-857-G	\$ 75
Shadow .....	0010-857-SH	75
Black .....	0010-857-BLK	75
White .....	0010-857-WH	75



HOLDER for universal vacuum/syringe/saliva ejector .....		
Forest Grey/Sterling .....	0010-886-G	\$ 35
Shadow .....	0010-886-SH	35
Black .....	0010-886-BLK	35
White .....	0010-886-WH	35







HOLDER, NOTCHED for universal vacuum/syringe/saliva ejector, fixed holder bar .....		
Forest Grey/Sterling .....	0010-885-G	\$ 35
Shadow .....	0010-885-SH	35
Black .....	0010-885-BLK	35
White .....	0010-885-WH	35

	PART NO	RETAIL
	HOLDER for standard vacuum Forest Grey/Sterling ..... 0010-887-G Shadow ..... 0010-887-SH Black ..... 0010-887-BLK White ..... 0010-887-WH	\$ 35 35 35 35
	HOLDER, NOTCHED for standard vacuum, fixed holder bar Forest Grey/Sterling ..... 0010-888-G Shadow ..... 0010-888-SH Black ..... 0010-888-BLK White ..... 0010-888-WH	\$ 35 35 35 35
	HOLDER for automatic standard vacuum Forest Grey/Sterling ..... 0010-837-G Shadow ..... 0010-837-SH Black ..... 0010-837-BLK White ..... 0010-837-WH	\$ 70 70 70 70
	HOLDER, NOTCHED for automatic standard vacuum, for fixed holder bar Forest Grey/Sterling ..... 0010-838-G Shadow ..... 0010-838-SH Black ..... 0010-838-BLK White ..... 0010-838-WH	70 70 70 70
	HOLDER for automatic accessory/universal vacuum Forest Grey/Sterling ..... 0010-846-G Shadow ..... 0010-846-SH Black ..... 0010-846-BLK White ..... 0010-846-WH	\$ 70 70 70 70
	HOLDER, NOTCHED for automatic accessory/universal vacuum, for fixed holder bar Forest Grey/Sterling ..... 0010-847-G Shadow ..... 0010-847-SH Black ..... 0010-847-BLK White ..... 0010-847-WH	70 70 70 70
	HOLDER, ELECTRIC for automatic standard vacuum, normally open Forest Grey/Sterling ..... 0010-858-G Shadow ..... 0010-858-SH Black ..... 0010-858-BLK White ..... 0010-858-WH	\$ 75 75 75 75
	HOLDER, ELECTRIC for automatic standard vacuum, normally closed Forest Grey/Sterling ..... 0010-868-G Shadow ..... 0010-868-SH Black ..... 0010-868-BLK White ..... 0010-868-WH	75 75 75 75
	HOLDER, NOTCHED ELECTRIC for automatic standard vacuum, normally open for fixed holder bar Forest Grey/Sterling ..... 0010-859-G Shadow ..... 0010-859-SH Black ..... 0010-859-BLK White ..... 0010-859-WH	\$ 75 75 75 75





# general information

## dimensions

### PIVOT DELIVERY SYSTEMS







### 1495PI





# general information

## dimensions



### FIXED WITH SIDEBOX DELIVERY SYSTEMS



### FIXED WITHOUT SIDEBOX DELIVERY SYSTEMS




### 1636 / 1637 BASIC




# general information

## dimensions

### 1800RM/TPM FLEX TRAY ASSEMBLY




### 3900 CHAIR




### UTILITY CENTERS


#### Contoured




#### Deluxe





#### Standard




#### Recessed





### 6198 OFFICE STOOL




### 6176 DOCTOR'S STOOL





### 6187 ASSISTANT'S STOOL



### 6126 DOCTOR'S STOOL




### 6137C / 6147C ASSISTANT'S STOOL




# general information

## dimensions


### CHAIR MOUNT LIGHT




### CEILING MOUNT LIGHT




**1300RM**




**1300TPM**



### 1300-C CEILING MOUNT MONITOR



### TRACK MOUNT LIGHT/MONITOR




8 MOUNTING HEIGHT REQUIRES 12" MOUNTING POST  
9 MOUNTING HEIGHT REQUIRES 24" MOUNTING POST  
10 MOUNTING HEIGHT REQUIRES 36" MOUNTING POST  
11 MOUNTING HEIGHT REQUIRES 48" MOUNTING POST

8 MOUNTING HEIGHT REQUIRES 12" MOUNTING POST  
9 MOUNTING HEIGHT REQUIRES 24" MOUNTING POST  
10 MOUNTING HEIGHT REQUIRES 36" MOUNTING POST  
11 MOUNTING HEIGHT REQUIRES 48" MOUNTING POST


# general information

## dimensions


### STANDARD WALL & CABINET MOUNT LIGHT




### TALL WALL & CABINET MOUNT LIGHT




### CABINET MOUNT




### STANDARD WALL & CABINET MOUNT MONITOR



### OPTIONAL CABINET MOUNT





### TALL WALL & CABINET MOUNT MONITOR



# general information


## dimensions

### 7000BC BENEATH COUNTER




- \* work surface 9/16" thick
- \* shown with optional storage organizer, tray and additional assistant's holder bar
- \* Height adjustable work surface, 5"


### 7000BC Mounting Plates





### 0010-086



### 0010-084-CS Option




### 7000WM WALL MOUNT




- \* work surface 9/16" thick
- \* shown with optional storage organizer, tray and additional assistant's holder bar
- \* Height adjustable work surface, 5"


### 7000WM Mounting Plate



### 7000WM Leveling/Mounting Board




### 1300-025-CS Monitor Assembly Option




# general information

## dimensions


### 7000UC BENEATH CABINET MOUNT




7000UC Under Cabinet Leveling/Mounting Plate




1300-025-CS Monitor Assembly Option




### 7000FM FLOOR MOUNT



7000FM FLOOR MOUNT Leveling/Mounting Plate




### 3386BC/WM 3486BC/WM 3686BC/WM WALL OR CABINET MOUNT




- \* work surface 9/16" thick
- \* shown with optional storage organizer, tray and additional assistant's holder bar


### 3386/3486/3686WM Leveling/Mounting Board




### 3386/3486/3686WM 3686BC/WM Rear Beneath Duo




### 3486BC/WM REAR BENEATH COUNTER MOUNT ASSISTANT'S




### 3686BC/WM REAR BENEATH DUO



### 3386/3486/3686BC Heavy Duty Mounting Plate



### 3386/3486/3686BC Mounting Arm




Work surface

Counter


# general information

## dimensions


### 2400 / 2700 REAR MOUNT BEHIND COUNTER




### 2400/2700 Leveling/Mounting Plate




### 1110-115-CS Folding Arm System Option



### 1700-140 Slide Mount with Swivel Option




### 3235 / 3735 CABINET SIDE DELIVERY



### 3230 / 3730 WALL SIDE DELIVERY



Shown: Standard IC  
Not shown: Value = 21"




# general information

## dimensions


### 2180 REAR, BENEATH WORK SURFACE




### 2180 Leveling/Mounting Plate




### 2180 WITH OPTIONAL 0010-652 WORK SURFACE AND 0023-806 ARM




### 1280 REAR, BENEATH WORK SURFACE



### 1280 Leveling/Mounting Plate




### 0023-806 Option




Work Surface

### 0023-815



Counter


### Cabinet Mounting Plate



# general information

## dimensions


### 5386 ASSISTANT'S CART



\* work surface 9/16" thick

\* shown with optional storage organizer, tray and additional assistant's holder bar


### 5886 DUO CART



\* work surface 9/16" thick

\* shown with optional storage organizer, tray and additional assistant's holder bar


### 7000MC DUO CART




\* work surface 9/16" thick

\* shown with optional storage organizer, tray and additional assistant's holder bar


### 5915 / 5910 VALUE CART




### 5200 / 5210 / 5400 / 5410 IC CART



### 5200 / 5210 / 5400 / 5410 Work Surface




### 5915 / 5910 VALUE Work Surface




# general information

## dimensions


**2500 PANEL MOUNT**




**2505 PANEL MOUNT**




**2180 / 2500 / 2505 Wall cut-out**




**2510 PANEL MOUNT**




**2515 PANEL MOUNT**




**2510 / 2515 Wall cut-out**




**MT600 / MT610 COMPACT CONTROLS**




**1055 ASSISTANT'S INSTRUMENTATION**




**1055 Leveling/Mounting Plate**




**1110-115-CS Folding Arm System Option**




**MT600/610BS Leveling/Mounting Plate**



**1700-140 Slide Mount with Swivel Option**



**1110-115-CS Leveling/Mounting Plate**



# general information

## ordering guide

ORDERING: When ordering Forest products, the following information should be included:

QUANTITY	ACCESSORIES OR OPTIONS
MODEL/PART #	DATE REQUIRED
DESCRIPTION	PURCHASE ORDER NUMBER
COLOR	METHOD OF SHIPPING

*Deviation from our standard colors and Designer Friendly options are subject to an upcharge.*

### STANDARD LAMINATE WORK SURFACE COLORS

(special laminates: call customer service for pricing and quantity discounts.)

FORMICA:	Almond 920-58 Folkstone 927-58	Graphite 837-58 Oyster Gray 929-58	Surf 923-58 White Blanc 949-58
WILSONART:	Shadow D96-60		

### STANDARD SOLID WORK SURFACE COLORS (same as Modular and Custom Cabinetry)

(Special non-stock solid work surface colors: call customer service for pricing and quantity discounts.)

CORIAN:	Antarctica Doeskin	Fossil Savannah	Blue Pebble
GIBRALTAR:	Arctic Melange Northern Melange	Frosty White Mirage Light Beige Melange	Baja Melange Midnight Melange
MEGANITE™:	Moonstone	Goldrush Gemstone	

### STANDARD UPHOLSTERY COLORS (Pearlized colors tend to wear similarly to soft leather)

Ultraleather™®:	Baltic Charcoal Diplomat Blue Doe	Dove Grey Grenadine Papyrus Parrot	Raven Wing Sassy Stone Walnut	Pearlized Green Tea Pearlized Mica Pearlized Midnight Pearlized Pewter
Ultraleather™® "RAFFIA" (upcharge):				Durango Dust      Noche

### DESIGNER FRIENDLY PAINT COLORS

Arctic Silver	Forest Grey	Sage
Bronze	Pearl	Shadow
Charcoal	Pewter	White

### DESIGNER FRIENDLY PLASTIC COLORS

Solid Colors:	Forest Grey	White
Carbon Fiber Patterns:	Greywood	Lemans      Laser Leather

### TUBING COLORS (no upcharge - colors configured by Forest)

Black	Shadow
Sterling (standard on Forest Grey and White)	

### STOOL PLASTIC COLORS (colors configured by Forest)

Black	Beige
-------	-------

### CHROME ACCENTS (\$135 upcharge)

### CUSTOM COLORS DEVIATING FROM STANDARD SUBJECT TO A 10%-15% UPCHARGE AND EXTENDED DELIVERY TIME

# general information

## important details regarding your NEW FOREST 2013 catalog

### DOCTOR'S CONTROL HEADS

4195 = Value, 4485 = Standard IC, 4987 = Euro

*Contoured IC Doctor's Control Head* (example: 4585, 8350, 8138, 3530, 3535) have been discontinued.

*Hygienist Control Heads* no longer have a unique model number. Now simply add the assistant's instrumentation option part # 1409-CS (which will add an HVE, saliva ejector and canister) to any Standard IC or Value Chair Mounted Control Head (not available on Euro Control Head). Select 1700-099-CS to delete a handpiece position if desired. Please note this does NOT apply to the "SI/SE" sidebox chair mounted deliveries since they come standard with assistant's instrumentation on a telescoping arm.

### CHAIR MOUNT CONFIGURATIONS

*For 3900 chair or can adapt to existing 3900 chair:*

PI = Pivot Chair Mount Delivery w/internal umbilical, cast chair adapter and DCI bottle Quick Disconnect.

NEW Unit part number:

- 4195PI Value Control (previously was model 4195), 4987PI Euro Control (previously was model 4987), 4485PI Standard IC Control (previously was model 4485). For Hygienist Control on 4485PI or 4195PI add 1409-CS Assistant's instrumentation/control head mount and 1700-099-CS to delete a handpiece position.

SI = Fixed Chair Mount w/Sidebox Delivery w/internal umbilical, cast chair adapter and DCI bottle Quick Disconnect.

NEW Unit part number:

- 4195SI Value Control (previously was model 4118 w/external umbilical), 4485SI Standard IC Control (previously was model 8300 w/external umbilical), 4987SI Euro Control (previously was model 8807 w/external umbilical).

FI = Fixed Chair Mount w/o Sidebox Delivery w/internal umbilical, cast chair adapter and DCI bottle Quick Disconnect.

NEW Unit part number:

- 4195FI Value Control (previously was model 4124 w/external umbilical), 4485FI Standard IC Control (previously was model 4338 w/external umbilical), Euro control head not offered in this configuration (previously was model 4167). For Hygienist Control add 1409-CS Assistant's instrumentation/control head mount and 1700-099-CS to delete a handpiece position.

*Delivery Systems for other manufacturer's chairs:*

FE = Fixed Chair Mount w/o Sidebox Delivery w/external umbilical with flat steel adapter and DCI bottle Quick Disconnect.

NEW Unit part number:

- 4195FE Value Control (previously was model 4124), 4485FE Standard IC (previously was model 4338), Euro control head not offered in this configuration (previously was model 4167). For Hygienist Control add 1409-CS Assistant's instrumentation/control head mount and 1700-099-CS to delete a handpiece position.

SE = Fixed Chair Mount w/Sidebox Delivery w/external umbilical with flat steel adapter and DCI bottle Quick Disconnect.

NEW Unit part number:

- 4195SE Value Control (previously was model 4118), 4485SE Standard IC Control (previously was model 8300), 4987SE Euro Control (previously was model 8807).

### DELIVERY SYSTEMS: REAR/SIDE DELIVERIES AND CARTS

- All rear/side deliveries and carts come standard with DCI bottle Quick Disconnect (except 1055 self-contained clean water system is optional).
- No longer available: 7380, 7381, 7382, 7383, 7384 Assistant's Deliveries, 7374 Work Surface Only, 3680/3480 Beneath Counter Deliveries, MT600/610WM (Wall Mount) Compact Control, 5380/5880/5370 Carts, 1070 Assistant's Telescoping Arm, 1570, 1575, 1650, 1655, 1690 & 1695 Cuspidors.
- NEW part number: 7000BC, 7000FM, 7000UC, 7000WM, 7000MC Beneath Counter, Floor Mount, Under Cabinet, Wall Mount and Mobile Cart delivery; storage organizer, tray and additional assistant's holder optional (previously was model 7681, 7682, 7683, 7680, 7684).
- NEW part number: 3386WM/3486WM/3686WM Wall Mount, 3386BC/3486BC/3686BC Beneath Counter (previously was model 3485/3685), 5386/5686 Assistant's and Duo carts (previously was model 5380/5880) with semi-round top; storage organizer, tray and additional assistant's holder optional.
- Doctor's Standard IC Delivery 3235/3230: previously offered in Hygiene control 3245/3240.  
Doctor's Value Delivery 3735/3730: previously offered in Hygiene control 3745/3740.  
Doctor's Standard IC Delivery 2400: now available in a Hygiene control.  
NEW Doctor's Value Delivery 2700: previously offered in Hygiene control 2745.  
For Hygienist control, order Dr's Control (listed above) with 1408-CS Assistant's instrumentation/control head mount and 1700-099-CS to delete one of three handpiece positions.
- MT600/610BS and TS (Beneath shelf and Top shelf Mount) Compact Controls available only in Forest Grey and White powder coat and plastic.
- 1055 Beneath Counter Assistant's Instrumentation available only in Forest Grey and White powder coat.
- NEW part number: 5210/5410 Doctor's and Duo hospital carts ready for nitrogen hook-up. Available only in Forest Grey and White.

# general information

## terms + warranty

### TERMS

All prices unless otherwise listed are F.O.B. Hillsboro, Oregon U.S.A. Payments are net 30 days of invoice date. A service charge of 1-1/2% per month (annual rate 18%) will be charged on past due accounts, \$1.00 minimum. All prices are subject to change without notice.

### RETURN MERCHANDISE

Authorized Forest dealers wishing to return merchandise for credit consideration must:

Obtain prior factory authorization.

Provide basic reason for return.

Provide serial number, purchase order or invoice number of original order.

Provide a purchase order for warranty items when necessary. Credit will be issued by Forest upon determination of warranty coverage in accordance with the stated Warranty Policy.

Return product within a period of 90 days—if product is not returned within time period, authorization becomes invalid.

Forest will approve return and credit amount after evaluation of product. If product is used or modified, the reason for return will be used to determine if replacement product or credit will be allowed.

#### Handling Charges:

25% or \$15.00 minimum restock fee will be charged on returned equipment approved for credit.

10% or \$15.00 minimum restock fee will be charged on new and unopened parts/components.

Returned items must be shipped freight prepaid.

### WARRANTY

FOREST warrants its products to be free from defects in material and workmanship only. No other warranties are expressed or implied. This warranty shall extend for five years unless otherwise stated. Written notice of breach must be given to FOREST within this period. Buyer's remedy for breach of this warranty is limited to repair parts or replacement of items by FOREST. Factory repairs are covered for a period of 90 days. No claim for labor or consequential damages will be allowed. Warranty is valid only if Warranty Registration Card provided at time of purchase is returned to FOREST within 90 days of installation. FOREST is not responsible for ensuring equipment or installation conformity to local code or regulatory requirements. Check with local code or regulatory requirements before purchase and installation of all FOREST product.

#### Warranty is void if:

Items are carelessly used or improperly maintained or installed.

Damage results from using cleaning, disinfecting or sterilization chemicals and processes.

Damage is the result of freight damage and as such claims must be filed with the carrier.

#### Warranty term exceptions cover items that are normally expected service items and items not covered by Forest's warranty:

Handpiece illumination tubings, heated syringe tubings & assemblies and light bulbs are warranted for six months.

Power supplies, syringe buttons, HVE & SE valves are warranted for one year.

Water bottles, gas springs and handpiece tubings are warranted for two years.

Air & Water tubing connectors (T's), Distribution Block, Master Control Block, Flow Valves, 7000 Chassis are warranted for lifetime of dental unit (maximum of 20 years).

Bien Air handpiece motors, Satelec scaler, Cavitron and Air Techniques intraoral camera are covered under manufacturer's warranty. Installation of ancillary devices are warranted for two years.

### CANCELLATION POLICY

#### Standard orders:

- Fax cancellation request required 10 working days prior to promised ship date, otherwise a 25% restock fee will be charged.
- Expedited ship date requests are non-cancellable.

### CUSTOM ORDERS:

- Custom configurations are those that deviate from Forest's published product offering:

Orders are non-cancellable and non-returnable.

May be subject to an additional "Deviation" charge and/or a non-refundable engineering charge.

- Custom powder coat, custom plastic surface treatments and/or custom treatments such as multiple color combinations or custom colors may result in an upcharge of 10-15% per model and are non-cancellable and non-returnable.

# general information

## professional sales associates (PSA)

### MARK ASPLUND

Cell: 206-954-9421  
Email: masplund@profsales.com  
OR, WA, Northern ID, MT, Northern WY, AK

### MILT BAISLEY

Cell: 203-645-7431  
Fax: 203-248-0427  
Email: mbaisley@profsales.com  
Upstate NY-Albany, Syracuse, Newburgh

### MICHEL CABANIE

Cell: 514-793-3368  
EFax: 514-221-3368  
Email: mcabanie@profsales.com  
Quebec, Nova Scotia, New Brunswick,  
Ottawa Canada

### STEVE CARLIN

Cell: 952-491-1223  
Email: scarlin@profsales.com  
MN, SD, ND, IA, MB Canada

### CHARLIE CLARKE

859-268-8856  
Cell: 859-509-6298  
Fax: 859-263-0302  
Email: cclarke@profsales.com  
KY, Southern OH, IN, Western WV

### JIM DOMITRZ

Cell: 617-921-4580  
EFax: 603-218-7400  
Email: jdomitrz@profsales.com  
CT, MA, ME, NH, RI, VT

### RICH ELLIS

Cell: 301-502-7424  
EFax: 888-814-8630  
Email: rellis@profsales.com  
MD, VA, Southern WV, DC

### BILL HIGHLAND

Cell: 810-623-7297  
Fax: 517-540-9690  
Email: bhigland@profsales.com  
MI, OH, IN

### MATT HIGHLAND

Cell: 901-606-3633  
Fax: 901-756-0564  
Email: mdhighland@profsales.com  
New Orleans, MS, Memphis, Nashville,  
Little Rock

### DAN KNOX

Cell: 303-883-8544  
Fax: 623-551-9426  
Email: dknox@profsales.com  
AZ, NM, Southern WY, CO, El Paso TX

### MATTHEW LANGDON

Cell: 516-316-5913  
Fax: 516-706-2844  
Email: mlangdon@profsales.com  
Greater NYC, Westchester & Long Island,  
Northern & Central NJ

### MILTON LORRING

Cell: 216-396-2341  
Fax: 216-803-9910  
Email: mlorring@profsales.com  
Central and northern OH, Western PA; Buffalo  
& Rochester NY, Northern WV

### JIM McGONIGAL

Cell: 916-759-2939  
Fax: 916-405-4001  
Email: jmcgonigal@profsales.com  
President and Western Regional Manager

### ALEX MARSDEN

508-748-6789  
Cell: 508-728-6395  
Fax: 508-355-0539  
Email: amarsden@profsales.com  
Eastern Regional Manager

### CHARLIE MARTIN

Cell: 330-421-4959  
Fax: 330-769-2853  
Email: cmartin@profsales.com  
Central Regional Manager

### PAT MOYLAN

Cell: 816-729-2900  
Fax: 208-730-2759  
Email: pmoylan@profsales.com  
KS, NE, MO, Southern IL

### PAUL MYERS

Cell: 214-704-5100  
Fax: 253-369-1308  
Email: pmyers@profsales.com  
Northern TX, OK, Fayetteville AR,  
Shreveport LA

### KEITH O'BERRY

Cell: 336-210-1429  
Fax: 336-855-0184  
Email: koberry@profsales.com  
NC, SC, Savannah GA

### TOM O'BERRY

Cell: 404-307-4931  
Email: toberry@profsales.com  
GA (except Savannah & Valdosta), AL, East TN,  
Pensacola FL

### RICK ORME

Cell: 916-837-3839  
Fax: 916-984-8407  
Email: rorme@profsales.com  
Northern CA, NV (except L.V.) HI

### SAM PAOLA

Cell: 310-569-1178  
Fax: 206-770-7477  
Email: spaola@profsales.com  
Southern CA, Southern ID, UT, Las Vegas

### MICHAEL PARK

281-251-9155  
Cell: 281-702-3992  
Fax: 281-251-9125  
Email: mpark@profsales.com  
Southern TX

### BILL REINHARDT

847-382-2364  
Cell: 847-840-3330  
Email: breinhardt@profsales.com  
Northern & Central IL, WI

### AARON RENNER

Cell: 519-319-7457  
Fax: 519-637-2466  
Email: arenner@profsales.com  
Ontario Canada (except Ottawa)

### ROBERT RODRIGUEZ

Cell: 305-458-6538  
Fax: 305-723-5909  
Email: rrodriguez@profsales.com  
FL, Valdosta GA, Puerto Rico

### DON ROGERS, III

Cell: 717-201-7332  
Fax: 866-764-5980  
Email: drogers@profsales.com  
Eastern PA, Southern NJ, DE

### SCOTT RUMBLE

Cell: 403-801-5390  
Email: srumble@profsales.com  
BC, AB, SK Canada

### MEL SEAGLE

Cell: 949-678-4550  
Fax: 949-258-5800  
Email: mseagle@profsales.com  
Southern CA, Southern ID, UT, Las Vegas

### VICTOR YOON

Cell: 778-688-8429  
Email: vyoon@profsales.com  
BC, AB, SK Canada

# general information

## Index

### MODEL NUMBER + PAGE NUMBER INDEX

1055 .....	53	2505 .....	51	4987SI .....	2, 14
1075 .....	36, 38	2510 .....	50	5200 .....	58
1075FI .....	2, 18, 20	2515 .....	51	5210 .....	59
1095 .....	2, 4, 6, 8	2700 .....	41, 48	5386 .....	56
1280 .....	41, 49	3230 .....	46	5400 .....	58
1300 .....	3, 28, 32	3235 .....	41, 46	5410 .....	59
1310 .....	31	3386 .....	44, 45	5886 .....	56
1495PI .....	2, 10	3486 .....	40, 44, 45	5910 .....	57
1576 .....	2, 14, 16, 34	3686 .....	40, 44, 45	5915 .....	57
1577 .....	2, 12, 14, 16, 34	3730 .....	47	6126C .....	3, 23
1580 .....	2, 4, 6, 8	3735 .....	41, 47	6137C .....	3, 23
1590 .....	2, 4, 6, 8	3900 .....	24, 25	6147C .....	3, 23
1636 .....	38	4195FE .....	36	6176 .....	3, 23
1636FE .....	36	4195FI .....	2, 20	6187 .....	3, 23
1636FI .....	2, 18, 20	4195PI .....	2, 9	6198 .....	26
1637 .....	38	4195SE .....	34	7000 .....	40, 42, 43, 55
1637FE .....	36	4195SI .....	2, 16	9045 .....	3, 27, 28, 29
1637FI .....	2, 18, 20	4485FE .....	36	9065 .....	30
1678 .....	2, 14, 34	4485FI .....	2, 18	9073 .....	3, 27, 28, 29
1679 .....	2, 12, 14, 16, 34	4485PI .....	2, 5	9074 .....	30
1800 .....	39	4485SE .....	34	9080 .....	3, 27, 28, 29
2180 .....	41, 49	4485SI .....	2, 12	9090 .....	30
2400 .....	41, 48	4987PI .....	2, 7	MT600 .....	41, 52
2500 .....	50	4987SE .....	34	MT610 .....	41, 52

### FACTORY OPTION + PART NUMBER INDEX

0010-084 .....	68	1110-537-RM .....	69	1700-021-CS .....	64
0010-086 .....	68	1110-680 .....	68	1700-140 .....	41, 68
0010-652 .....	49	1110-771-RM .....	69	1700-141 .....	41, 68
0016-454 .....	69	1110-820 .....	3	3400-005 .....	59
0016-455 .....	69	1110-900 .....	3	3900-101 .....	67
0016-456 .....	69	1112-644 .....	3	3900-102 .....	67
0017-603 .....	64	1115-001 .....	67	3900-200 .....	61
0017-615 .....	64	1115-007 .....	67	3907-180 .....	61
0023-665 .....	68	1115-020 .....	67	3907-999-MC .....	25
0023-806 .....	68	1115-305 .....	74	3907-999-NB .....	25
0024-054 .....	10	1115-315 .....	69	3907-999-NS .....	25
0024-224 .....	66	1115-316 .....	69	AIRTECH .....	63
1107-450-CS .....	65	1400-113 .....	3, 66	BIEN-AIR .....	62
1110-115 .....	41, 68	1400-115 .....	69	CAVITRON .....	62
1110-349-CS .....	68	1408-CS .....	63	SATELEC .....	3, 62
1110-468 .....	3, 4, 8, 12, 18	1409-CS .....	3, 63		
1110-536-RM .....	69	1700-019-CS .....	3, 64		

## Surface Treatments

---

Forest is the only dental equipment with a WOW factor. It really is a kick watching a dentist react to Forest equipment when they first see it. A common reaction is “wow, that is beautiful.”

---

### Carbon Fiber

---

#### Control Heads

- Standard IC - chassis
- Euro - lid
- Value - lid

#### Sidebox - cover

#### Basic Cuspidor - cover

#### Assistant's Instrumentation - arm cover

---

---

### Powder Coat

---

Dentists and interior designers love being able to put something besides “computer monitor white” in their beautiful patient friendly decors.

#### Control Heads

- Standard IC - lid, brake handles & holder bar
- Euro - chassis, whip arms, brake handles & swivel tray
- Value - chassis, brake handles & holder bar

#### Sidebox - housing

#### Basic Cuspidor - housing & holder bar

#### Assistant's Instrumentation - arm

---

---

### Ultraleather™ (standard)

---

What if you could only order a white car with your sole color being your upholstery? You could live with it but if given the choice, why would you? Once you see a Forest chair, delivery system and light package with color coordinated paint and upholstery in a well-designed treatment room, you wonder why it's taken so long for someone to offer dental equipment with a finish other than white.

#### 3900 Chair

- plush upholstery
- asepsis upholstery
- memory comfort upholstery
- pedo wedge
- IV board

#### 6126C Doctor's Stool

#### 6137C/6147C Assistant's Stools

---

---

### Solid Surface

---

- 7000 Series rear mount, side delivery & mobile cart (standard)
- 3386/3486/3686 rear mount, side delivery (option)
- 5386/5886 mobile cart (option)


---

### Laminate


---

- 7000 Series rear mount, side delivery & mobile cart (option)
- 3386/3486/3686 rear mount, side delivery (standard)
- 5386/5886 mobile cart (standard)

**Formica™**




graphite



surf



folkstone



oyster grey



white blanc



almond

**Gibraltar™**



arctic melange



frosty white mirage



northern melange



baja melange




light beige mirage



midnight melange

**Corian™**



antarctica



fossil



savannah



blue pebble




doeskin


**Powder Coat**




arctic silver



bronze




charcoal




forest grey




pearl




pewter



sage



shadow




white

**Wilson Art™**



shadow

**Meganite™**



goldrush gemstone



moonstone

## Ultraleather™



baltic



charcoal



diplomat blue



doe



dove grey



stone




grenadine



papyrus



pearlized midnight



parrot



raven wing



pearlized pewter




sassy



walnut



pearlized green tea



raffia noche



raffia durango



pearlized mica


## Designer Friendly Surfaces



forest dental

### Carbon Fiber

laser leather



lemans



greywood



<b>4</b>	Pivot Chair Mount Packages
<b>12</b>	Fixed Chair Mount w/ Sidebox Packages
<b>18</b>	Fixed Chair Mount w/o Sidebox Packages
<b>23</b>	Chairs & Stools <b>23</b> Doctor & Assistant's Stools <b>24</b> 3900 Hydraulic Patient Chair <b>26</b> Office Stool
<b>27</b>	Lights & Monitor Mount Assemblies <b>27</b> 9080 / 9045 / 9073 Lights <b>31</b> Monitor Mount Assemblies
<b>33</b>	Delivery System w/o Chair <b>34</b> Fixed Chair Mount w/ Sidebox <b>36</b> Fixed Chair Mount w/o Sidebox <b>38</b> Post Mount Cuspidors & Assistant's Instrumentation <b>39</b> Support Systems
<b>40</b>	Rear & Side Delivery Systems <b>40</b> Cabinetry Mounting Configurations <b>42</b> 7000 Series <b>44</b> Rear & Side Mount Deliveries
<b>50</b>	Optional Delivery Systems <b>50</b> Panel Mounts <b>52</b> Compact Controls <b>53</b> Assistant's Instrumentation
<b>55</b>	Carts
<b>60</b>	Factory Options & Parts <b>60</b> Chair Adapters, Monitor Posts, Light Posts & Chair Options <b>62</b> Ancillary Options <b>64</b> Factory Options & Parts
<b>72</b>	General Information <b>72</b> Dimensions <b>84</b> Ordering Guide, Terms & Warranty <b>87</b> PSA, Index

Affordable Excellence + Designer Friendly™

Phone

- 800 . 423 . 3555
- 503 . 640 . 3012

Fax

- 503 . 693 . 9715

Digital

- [sales@forestdental.com](mailto:sales@forestdental.com)
- [forestdental.com](http://forestdental.com)

Address

- 6200 ne campus court hillsboro oregon 97124



\* sales videos available online at [forestdental.com](http://forestdental.com)

