ULTRA CLEAN OILLESS AIR COMPRESSOR



- COMPRESSOR DESIGN OFFERS THE QUIETEST COMPRESSORS IN THE INDUSTRY
- ALL AIR RECEIVER TANKS HAVE INTERNAL EPOXY COATING APPLIED
- THE AIR IS DRIED BEFORE THE TANK TO PREVENT BACTERIA BUILD UP IN THE TANK
- ALL MODELS HAVE DESICCANT AIR DRYING SYSTEMS
- ENERGY EFFICIENT MOTORS REDUCE THE NEED FOR COSTLY BOOST TRANSFORMERS
- COMPRESSORS ARE AVAILABLE IN SINGLE, DUAL OR TRIPLE HEAD CONFIGURATIONS
- TEFLON COATED PISTONS AND SLEEVES ENSURE LONG LIFE OF THE COMPRESSOR HEADS
- AIR INTAKE FILTERS TRAP DIRT, DUST AND AIR BORN POLLUTANTS
- ENERGY EFFICIENT THERMALLY PROTECTED MOTORS
- MOISTURE INDICATORS STANDARD ON ALL MODELS



Manufacturers of Dental Vacuum and Air Systems

2625 N. Argyle Ave. • Fresno, CA 93727 (559) 291-1650 • (800) 428-7139 • FAX (559) 348-9677

Tech West Ultra Clean Oilless Air Compressor provides 100% oil free air to the dental office. With our high performance drying system we guarantee the clean dry air you expect. The design of the Tech West Oilless Compressor provides smooth vibration free running at 72 DB saving the expense and hassle of a sound cover. They are offered in a variety of configurations; singles, duals or triples to meet any need.

Tech West will provide full plumbing and electrical design upon request.



OILLESS AIR COMPRESSOR

FEATURES

Redesigned Air Dryer makes the unit more compact and reliable.

BENEFITS

This allows the compressor to be installed under cabinets and in closets, saving space.

American made electrical motor with capacitor start and run windings.

Faster starting and lower running amps saves electricity.

Motor will start under full compressor load with a voltage range of 208/248. No Italian or European motor can.

Compressors can be installed where others would need a buck/boost transformer saving hundreds of dollars.

One handle dryer shutdown.

Quick shutdown and ability to continue working until a service call can be made.

Moisture indicators on the outlet air line.

Early indication of moisture in the air before any condensation occurs.

Integrated carrying handles on the tank.

Makes for easy moving and installation.

OILLESS AIR COMPRESSOR

NOTE: NUMBER SEPTE STAFF CHEFT SEPTE																				
ACO2S1	2	YES	NO	1	1	4.8	10	115	17	20	70	159	27	26	19	189	37	35	30	
ACO2S2	2	YES	NO	1	1	4.8	10	208/230	10	20	70	159	27	26	19	189	37	35	30	
ACO4D2	4	YES	NO	2	2	9.6	20	208/230	20	30	74	228	31	36	21	278	43	42	31	
ACO6T2	6	YES	NO	3	3	14.6	30	208/230	30	40	77	298	33	41	24	348	43	42	31	

Revised 4/14

ALL OILLESS COMPRESSORS CARRY A FULL THREE YEAR WARRANTY ON OILLESS HEADS AND A FULL TWO YEAR WARRANTY ON ALL OTHER COMPONENTS







